
**REF**  
**410 233**

# PROTECTOR III

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**PL**

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EN

Please read these safety instructions carefully before use, and keep them for future question.

#### 1. Safety indications

- ♥ ALWAYS consult with your physician or supplier to determine proper adjustment and use of the device.
- ♥ Don't store outside and don't let under the sun (Beware: don't let the product in a car in case of high outside temperature!!), store in a dry room in a flat or rolled position to prevent deforming, stretching and tearing of the mattress overlay.
- ♥ Don't use in extreme temperature (above 38°C or 100°F or below 0°C or 32°F). Indoor use only.
- ♥ Store away from sharp objects and pointed edges to prevent tears and punctures.
- ♥ Use in a flat unfolded position.
- ♥ Look for cracks in mattress overlay and compressor before use.
- ♥ Transport rolled in its box to avoid tears and punctures.
- ♥ No modification must be made on this device.
- ♥ ALWAYS test to see that the mattress overlay is properly and securely locked in place before using. End flaps must be tucked under the foam mattress.
- ♥ In case you didn't use your mattress overlay and compressor for a long time, let it checked by the technical staff, before using.
- ♥ ALWAYS observe the weight limit on the labelling of your device.
- ♥ Use only as an overlay mattress and in proper way: don't kneel or step on the product.
- ♥ Follow the advise of the professional who prescribed this medical device, respect duration and frequency of use.
- ♥ If your skin develops redness, discontinue the use of the mattress overlay immediately and see your doctor or therapist.
- ♥ In case of discomfort consult your professional.

- ♥ Avoid contact with any sharp objects or pointed edges that may tear or puncture the mattress overlay.
- ♥ The compressor shall be connected directly to the electrical power. Do not use electrical outlets and extension cords. Connect this product to a properly grounded outlet only.
- ♥ Make sure that the use of other medical or electronic devices is possible and ensure the safe operation of each device.
- ♥ In case of doubt, if the mattress overlay or the compressor appear to be damaged or if you suspect a malfunction, immediately unplug the compressor from wall outlet, mark it as „DEFECTIVE“ and contact your technical staff.
- ♥ Unplug the power cable of the compressor from wall outlet, before moving the mattress, the compressor or the bed (on which the mattress overlay is placed) around.
- ♥ Check the power cable regularly for damages (crushed, bare, sheared, porous..) and particularly:
  - 1/after every mechanical strains (for example: rolling over the cables).
  - 2/after having moved the Home Care bed and before plugging in the power cable.
  - 3/every month in a normal use.
- ♥ To reduce the risk of electrocution: always unplug this product immediately after use.
- ♥ Do not place in or drop into water or other liquids.
- ♥ Product should never be left unattended when plugged in.
- ♥ Never operate this product if it has a damaged cord or plug, if it has been dropped, damaged, or dropped into water.
- ♥ Keep the cord away from heated surfaces.
- ♥ Never drop or insert any object into any opening or hose.
- ♥ Do not operate where aerosol products are being used or where oxygen is being administered.
- ♥ Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## 2. Intended use, indications & contraindications

### Intended use:

It is easier to prevent pressure ulcers than to treat them. Bedsores can be avoided by changing frequently lying or seating positions to avoid stress on the skin and to minimize the risk of pressure on vulnerable bony areas. Other strategies include skin care, regular skin inspections and good nutrition.

Using an alternating overlay mattress can help relieve pressure. Ideal to prevent stage I bedsores.

This air alternating mattress system is designed for patients who are at risk of pressure of ulcers or require treatment for existing pressure ulcers. This device is intended to treat and prevent pressure ulcers by facilitating blood circulation, reactive hyperaemia and decreasing interface pressures.

The use of this system does not replace the regular repositioning, monitoring, and nursing of the patient.

### Indications:

- Prevention of pressure ulcers in high-risk patients
- Pain management
- Treatment of existing pressure ulcers
- Prevention of skin breakdown, bedridden patient

### Contraindications:

- Unstable fractures of vertebrae and conditions of the vertebrae.
- Gross oedema, burns or an intolerance to motion

The PROTECTOR III set contains a pump equipped with quick connector and cap, visual low pressure alarm and correct pressure

indicator. The pad is equipped with 17 cells, where cells from 6 to 11 feature a microventilation function. The pad is equipped with a CPR strap and enables a rapid air deflation in case of an emergency. The quick connector features a transport function : the cap allows to stop air loss when necessary to transport the patient in the bed to another room.

This set includes one free air cell for repairs, easy to replace by the user in case of damage.

## 3. Description and Fitting instructions

**Beware : always place a bed sheet between patient and mattress.**

1. Position pump on horizontal surface or hang it by means of the hooks.
2. Open mattress on the bed. Secure it with end flaps or straps. Air tubes must be at foot end of bed.
3. Connect quick connector with air feed tubes to air outlet of the pump, ensuring that the tubes are not “kinked” or twisted. Allow space for bed rails to drop freely.
4. Plug the pump into a suitable electrical power source and switch on. The switch indicator should now illuminate, and pressure must be adjusted to MAX position.
5. Allow 20 minutes for mattress to inflate fully.
6. Cover with a sheet and tuck loosely.
7. Please use pressure adjustable knob to give maximum patient comfort.
8. Correct working pressure is indicated by a green lamp „PRESSURE NORMAL“.
9. In case that yellow indicator „PRESSURE LOW“ appears, it might be, that there is damage somewhere in mattress pad or air tubes. Stop using pad, and locate the damage. If you don't find any damage - please contact service.
10. Pump must work continuously, 24 hours a day.

## 4. Washing instructions

Pad can be cleaned by disinfectant solution.

Do not immerse the pump in water.

Do not heat or steam in autoclave.

Regularly check the condition of air filter - in case it is dirty or if you notice lower efficiency of the pump - replace filter.

## 5. Disposal

Dispose of this product in accordance with your regional waste disposal regulations.

## 6. Technical specifications

This product has a 2 years warranty, from the date of purchase, against any manufacturing defect. Duration of use : 2 years  
This product is biocompatible and CE certified.

### Compressor:

Power supply: 230V/50Hz . Air output: 8 l/min.

Pressure range: 40 mmHg – 100 mmHg

Noise level: below 35dBA

Cycle time: 10 to 12 min. Power consumption: 7W. Fuse: 1A

Size: 28cm(L) x 15cm(W) x 9cm(H). Weight: 1.4 kg

### Mattress overlay:

17 inflatable cell rolls, including 6 with microventilation function

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serial No.:

date of purchase:

This product has one year guarantee, from the date of purchase, against any manufacturing defects, which appears during using in normal conditions, according to this instruction.

### Guarantee does not include:

- mechanical damage, thermal, chemical and any other provided by user or third persons
- damage coming from using the product improperly or using damaged product
- natural usage of product or it's parts, that becomes used according to it's properties or destination
- service an maintenance, such as cleaning or replacing dirty air filter

### Guarantee is not valid:

- in case that the product was modified by third person, as well as when not original parts have been used
- in case of discovered modifications in this guarantee document or purchasing document
- in case of improper storage of the product or not correct maintenance

Guarantee is valid together with purchasing document, so please keep those documents carefully.

Complain should be placed by your distributor, where all conditions, terms and way of repair will be established, and in case that repair is not possible - delay to exchange the product to a new one.

Delay wil not exceed 14 work days, unless it is necessary to import replacement parts - then the delay can by prolonged with time necessary to get repacement parts.

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ul. Powstania Wlkp. 45, 88-320 Strzelno

tel. 52 318 38 57, fax 52 318 37 01, [www.aston.pl](http://www.aston.pl)

HERDEGEN SNC

2, Rue Freyssinet

77500 CHELLES - FRANCE