

7







# Bristol Maid Hospital Metalcraft<sup>™</sup> Limited Product Catalogue 2022









## MEDICAL FURNITURE & EQUIPMENT 2022

Founded in 1953, Bristol Maid has been manufacturing and supplying medical furniture & equipment for over 65 years. From design and manufacture to supply and installation, Bristol Maid has a reputation for providing high product quality, excellent customer service and value for money. We take pride in working in close partnership with our customers to deliver effective solutions based on innovative design and technology that will meet the demands of today while anticipating the solutions for tomorrow.

Our highly-trained staff are experienced in every aspect of customer support, from simple product selection to bespoke technical solutions, from on-site audits to complex project management arrangements, according to our customers' needs. Bristol Maid's commitment to quality is based on years of experience working with healthcare professionals in an exacting and regulated environment, and is supported by an ISO 9001:2015 quality management system that is registered with a UKAS accredited certification body.

If you require further assistance, please do not hesitate to contact our Customer Service team who are here to help; or alternatively can arrange for your Regional Sales Representative to visit and help find a solution that best suits your requirements.









#### **Family business**

Hospital Metalcraft Ltd is a family company. We own our factory, production lines and delivery lorries. We continually invest into our staff and the development of new product ranges.



### Healthcare focused

With an extensive portfolio that covers almost every discipline, we have solutions to cater for the demands of the healthcare environment.



### Manufactured in the UK

Hospital Metalcraft is proud to be a British manufacturing success story. Strict quality criteria are applied to both the Bristol Maid product range designed and produced within our own manufacturing facility, and complementary equipment from carefully selected partners.



#### Delivered using our own lorries

Wherever possible we use our own fleet of lorries to ensure that your equipment arrives safe and secure.



## The best quality

With over 65 years' experience, we are focussed on the design and manufacture of high-quality products which will create value for our customers. With a worldwide customer base, the Bristol Maid brand is recognised by many thousands of satisfied customers.



#### **Approved supplier**

We are an accredited NHS & public sector accredited supplier. This allows you to buy with confidence knowing your requirements for compliant and competitive tendering has already been met.

## Support you can rely on

Our external sales representatives, sales support engineers and our customer service team are here to help you at every stage, from advice with product selection, managing the delivery of your order through to full service and repair support.

## **HOW TO ORDER**

#### How to Order

- Orders are required in writing. Should you wish to order over the phone, payment is required in full (debit or credit card). Alternatively, we can issue a pro-forma invoice, payment is required in full prior to manufacture
- Orders should be submitted as an 'official purchase order' on company headed paper with all information clearly displayed
- Your order should contain the following information
  - Item No.
  - Description
  - Quantity
  - Delivery address
  - Invoice address
  - Full contact details, including name, email address, phone & fax numbers
- · Where a quotation has been issued, please reference this to ensure any agreed discounts are applied
- Item codes & descriptions in this catalogue supersede all previous editions
- · A full price list can be requested from our website
- Payments can be made by cheque, bank transfer, credit or debit card
- If you are a new customer, the first three orders will be processed on a pro-forma basis, paid in full prior to manufacture & delivery, after which a credit account may be applied for
- · For our terms and conditions refer to www.bristolmaid.com
- For further information and assistance, please contact our customer services team on +44 (0)1258 484455

#### **Payment Options**

#### BACS

- · Payments to: HSBC Bank PLC
- Sort Code: 40-12-18
- Account Number: 11044052

#### Credit or Debit Card

• Major credit cards accepted (with exception of American Express)





## **HOW TO ORDER**

#### **NHS & Public Sector Accredited Supplier**

Public sector organisations can take advantage of our nationally agreed discounts through the NHS Supply Chain and Shared Business Services. Buy with confidence knowing your requirement for compliant and competitive tendering has already been met.

#### **NHS Supply Chain**

Medical Healthcare Furniture - Ref: 2020/S 252-638069

#### **Shared Business Services**

Design, Furniture & Appliances - Ref: SBS/19/MA/TYX/9363/23

#### **National Services Scotland**

Medical Ward Equipment - Ref: NP846/21

#### **HSC Business Services**

Medical Ward Equipment - Ref: NP846/21

#### Your Partner for Healthcare Equipment & Furniture

Through many years' experience we can support you with the facilitation and project management of tenders and contracts.

- · Specifying Requirements We will analyse your needs to identify products that meet your requirements
- Product Demonstrations We will allow you the opportunity to gain 'hands on' experience of our products to ensure that it is suitable for your needs
- Project Planning We will work with you and your team to identify any constraints and agree delivery phasing to ensure that we are able to fully meet your expectations
- Scheduled Deliveries We will deliver your order using our own transport to an agreed delivery schedule and, if agreed place the product and remove any packaging
- Installation & Commissioning The majority of the products that we supply requires no installation or commissioning, however we
  can provide this where it is required
- Training We can supply training to ensure that the products are used and maintained in the correct manner
- Full Service & Repair Support Bristol Maid branded products are designed and manufactured to meet the demands of everyday use, which under normal conditions should provide many years of trouble-free use. However, there are times where attention is required due to wear and tear or misuse
- Maintenance Contracts Peace of mind service contracts to ensure that your product is periodically checked, and where necessary adjustments and repairs are made
- · Site Audit & Equipment Recommendations We will support you with site and equipment audits

#### **Custom Design & Manufacture**

Whether it is a trolley that is slightly lower, a work surface that is slightly larger, a cabinet that must fit into a specific space or a unique colour scheme on a bedside cabinet, we recognise

that our standard product range may not always to be able to provide a product that meets your specific requirements.

If you require a unique product, please contact us so we can work together to develop a solution that will exactly match your needs. Our design capability allows us to provide everything from minor product specification changes through to a full design & development service.

#### **Bariatric Equipment**

| Folding Wheelchairs            |  |
|--------------------------------|--|
| Portering Chairs               |  |
| Mobile Shower / Commode Chairs |  |
| Shower / Commode Chairs        |  |
| Chairs                         |  |

#### **Bedside Cabinets & Overbed Tables**

| Polymer Bedside Cabinets                |
|---|
| Bedside Cabinets                        |
| Bedside Cabinet / Wardrobe Combinations |
| Polymer Overbed & Overchair Tables      |
| Laminate Overbed & Overchair Tables     |

#### **Caretray Trolleys & Racks**

| Features                                 |
|--|
| Disassembly & Cleaning                   |
| Stainless Steel, Low Level Trolleys      |
| Stainless Steel, Standard Level Trolleys |
| Stainless Steel, High Level Trolleys     |
| Stainless Steel, Mobile Racks            |
| Stainless Steel, Static Racks            |
| Mild Steel, Low Level Trolleys           |
| Mild Steel, Standard Level Trolleys      |
| Mild Steel, High Level Trolleys          |
| Mild Steel, Mobile Racks                 |
| Mild Steel, Static Racks                 |
| Accessories                              |

#### **Chart Trolleys, Workstations Room Dividers & Screens**

| Fixed Height Trolleys        |
|------------------------------|
| Fixed Height Workstations    |
| Variable Height Workstations |
| Chart Boards                 |
| Document Holders             |
| Room Dividers                |
| Curtain Screens              |
| Solid Screens                |
|                              |

#### **Concealment Trolleys**

| Variable Height Trolleys |  |
|--------------------------|--|
| Body Trays & Covers      |  |

#### **Cylinder Storage & Transportation**

| Cylinder Trolleys - F/G Size |  |
|------------------------------|--|
| Cylinder Trolleys - J Size   |  |

#### **Cylinder Storage & Transportation**

| Bulk Transportation Trolleys     |
|----------------------------------|
| Cylinder Trolleys - E, CD/D Size |
| Wall Brackets & Cylinder Holders |

#### **Dispensers, Holders & Racks**

| Apron Dispensers, Pocket Mask Holders            |
|--|
| Glove / Mask Box Holders                         |
| Alcohol Gel Holders                              |
| Sharps Box Holder, Urine Bottle Holder, Catheter |
| Bag Holders                                      |
| Catheter Racks                                   |
| Drainage Racks                                   |
|  |

#### **Dressing / Instrument Trolleys**

| Removable Shelves - Stainless Steel         |
|---|
| Removable Shelves - Mild Steel              |
| Fixed Shelves, Sides Up - Stainless Steel   |
| Fixed Shelves, Sides Down - Stainless Steel |
| Drawer Units                                |
| Options & Accessories                       |
|   |

#### **Drug & Medicine Dispensing**

| Dispensing Trolleys                       | 124 |
|---|-----|
| Laptop Holders                            | 129 |
| CLIQ <sup>®</sup> & RFID Solutions        | 131 |
| Patient Administration Trolleys           | 132 |
| Drawer Trolleys                           | 134 |
| Tray Trolleys                             | 135 |
| Blister Pack (MDS) Trolleys               | 136 |
| Nomad <sup>®</sup> Cassette Trolleys      | 137 |
| Biodose <sup>®</sup> & MultiMeds Trolleys | 138 |
| Medinoxx Trolleys                         | 139 |
| Security Clamp                            | 142 |
| Distribution Trolleys                     | 143 |
|   |     |

#### **Drug & Medicine Storage**

| Drug & Medicine Cabinets - Shelves & Door Shelves      | 146 |
|--|-----|
| Drug & Medicine Cabinets - Deep Shelves                | 148 |
| Drug & Medicine Cabinets - Blister Pack (MDS)          | 150 |
| Controlled Drug Cabinets                               | 152 |
| Auditable Access Control for Medicines Management      | 154 |
| Drug & Medicine Cupboards                              | 156 |
| Options  | 158 |
| Patients / Residents Own Medication Cabinets - Self    | 160 |
| Administration   |     |
| Patients / Residents Own Medication Cabinets - Blister | 163 |
| Pack (MDS)   |     |
|  |     |

#### **Drug & Medicine Storage**

| Tablet Counter & Storage Boxes   | 164 |
|----------------------------------|-----|
| Pharmacy & Vaccine Refrigerators | 165 |

#### **Examination / Treatment Couches & Chairs**

| Examination & Treatment Couches - Variable Height  | 168 |
|--|-----|
| Examination & Treatment Couches - Fixed Height     | 172 |
| Couch Steps  | 173 |
| Leg Rests  | 174 |
| Examination & Treatment Couches - Variable Height, | 176 |
| Bariatric  |     |
| Examination & Treatment Chairs - Variable Height   | 180 |
| Examination & Treatment Chairs - Variable Height,  | 181 |
| Bariatric  |     |
| Phlebotomy Chairs - Variable Height                | 182 |
| Phlebotomy Chairs - Variable Height Bariatric      | 183 |
| Phlebotomy Chairs - Fixed Height                   | 184 |
| Options  | 185 |
|  |     |

#### **General Purpose Trolleys**

| Stainless Steel - Heavy Duty     | 190 |
|----------------------------------|-----|
| Stainless Steel - Light Duty     | 191 |
| Polymer - Shallow & Deep Shelves | 192 |
| Mild Steel - Laminate Shelves    | 193 |

#### **HTM71 Storage Solutions**

| ScanCell <sup>®</sup> Carts            | 196 |
|--|-----|
| Distribution Trolleys                  | 201 |
| Open Transfer / Transport Trolleys     | 203 |
| HTM71 Racking                          | 204 |
| Trays, Baskets, Dividers & Accessories | 206 |
|  |     |

#### Infusion / Pump Stands & Accessories

| Mobile Pump Stands            | 212 |
|-------------------------------|-----|
| Mobile Infusion Stands        | 216 |
| Infusion Poles                | 220 |
| Clamping Infusion Stands      | 221 |
| Wall Mounted Infusion Holders | 221 |

#### **Linen Trolleys & Bags**

| Soiled Linen Trolleys | 224 |
|-----------------------|-----|
| Linen Trolleys        | 227 |
| Hamper Trolleys       | 228 |
| Soiled Linen Bags     | 229 |
| Clean Linen Trolleys  | 230 |

#### **Maternity Equipment**

| Baby Cribs - Variable Height, Battery Operation   |  |
|---|--|
| Baby Cribs - Fixed Height                         |  |
| Baby Cribs - Variable Height, Gas Strut Operation |  |

#### **Medical Records**

| Medical Records Trolleys - Locking     |
|--|
| Medical Records Workstations - Locking |
| Laptop Holders                         |
| Carefile Filing Pockets                |
| Medical Records Trolleys - Filing      |
| Medical Records Transfer Trolleys      |

#### Phlebotomy

| Caretray Trolley - Phlebotomy                  |
|--|
| Phlebotomy Chairs - Variable Height            |
| Phlebotomy Chairs - Variable Height, Bariatric |
| Phlebotomy Chairs - Fixed Height               |

#### **Portering Chairs**

| Portering Chair - Rear Steer            |
|---|
| Portering Chair - Front Steer           |
| Portering Chair - Rear Steer, Bariatric |
| Options & Accessories                   |

#### **Procedure Trolleys & Carts**

| Caretray Trolley - Anaesthesia             |
|--|
| Caretray Trolley - Cast                    |
| Caretray Trolley - Critical Care           |
| Caretray Trolley - Infection               |
| Caretray Trolley - IV Preparation          |
| Caretray Trolley - Medical / Surgical      |
| Caretray Trolley - Plaster                 |
| Caretray Trolley - Vaccination             |
| Caretray Trolley - Basics                  |
| Mobile Workstation - Mild Steel, Low Level |
| Procedure Carts - Polymer, Metro Flexline™ |
| Procedure Carts - Polymer, Metro Starsys™  |
| Procedure Carts - Polymer, Metro Basics™   |

#### **Resuscitation Trolleys & Carts**

| Resuscitation Trolleys - Mild Steel                 |
|---|
| Caretray Trolley - Resuscitation                    |
| Resuscitation Carts - Polymer, Metro Lifeline™      |
| Resuscitation Carts - Polymer, Metro Flexline™      |
| Airways / Difficult Intubation Trolley - Mild Steel |
| Caretray Trolley - Airways / Difficult Intubation   |

#### **Sackholders & Waste Bins**

| Removable Body Bins                  |  |
|--------------------------------------|--|
| Free Standing Sackholders            |  |
| Fixed Body Sackholders - Plastic     |  |
| Removable Body Sackholders - Plastic |  |

#### **Seating Solutions**

| Office, Task & Operators Chairs            |
|--|
| 24 Hour Task & Operator Chairs             |
| Bariatric Chairs                           |
| Saddle & Perch Stools                      |
| Low Back Patient Chairs                    |
| Medium & High Back Patient Armchairs       |
| Reclining Armchairs                        |
| Reclining Armchairs, Bariatric             |
| Reclining Porters Armchairs                |
| Rise & Recline Porters Armchairs           |
| Seat, Back Cushion & Headrest Options      |
| Overnight Chair                            |
| Soft Seating                               |
| Visitor & Waiting Room Seating             |
| Waiting Room Seating                       |
| Tech Chairs & Stools                       |
| PU Tech Chairs & Stools                    |
| Static Safe Tech Chairs & Stools           |
| Static Safe & Sterile Tech Chairs & Stools |
| Medical Grade Clean Room Tech Chairs       |
|  |

#### Sepsis Trolleys

#### **Shelving & Racking Solutions**

| 390 |
|-----|
| 392 |
| 394 |
| 396 |
| 398 |
| 399 |
|     |

#### **Sterile Services Equipment**

| Preparation Tables - Variable Height          | 402 |
|---|-----|
| Preparation Tables - Fixed Height             | 403 |
| Sterile Supplies Trolleys - UN3291 Certified  | 404 |
| Sterile Supplies Trolleys - UN3291 Certified, | 406 |
| Trolley Washer Compatible                     |     |

#### Storage Cupboards, Cabinets & Lockers

| Storage Cupboards  | 410 |
|--|-----|
| Storage Cabinets - Single & Double Door                  | 412 |
| Storage Cabinets - First Aid & Clinical Reagents Storage | 413 |
| Cabinets - Traka 21                                      | 414 |
| Storage Cabinets - Keys                                  | 415 |
| Storage Lockers - Personal Effects                       | 416 |

#### **Theatre Equipment**

| 420 |
|-----|
| 423 |
| 424 |
| 424 |
| 120 |
| 426 |
| 427 |
| 428 |
|     |

#### Wireless Mobile Carts & Workstations

| Mobile Workstations - Non-Powered | 432 |
|-----------------------------------|-----|
| Options & Accessories             | 434 |
| Caretray Trolley - Laptop         | 436 |

#### Workstations & Speciality Storage

| Workstations           | 440 |
|------------------------|-----|
| Storage Units          | 441 |
| Catheter Storage Units | 442 |



# **Bariatric Equipment**

- Folding Wheelchairs 8
  - Portering Chairs 10
- Mobile Shower / Commode Chairs 11
  - Shower / Commode Chairs 12
    - Chairs 13



## **Bariatric Folding Wheelchair - Minimaxx & Disc**

#### **Features**

- Maximum capacity 325kg
- · Powder coated chromium molybdenum steel frame
- Rigid folding mechanism
- Low seat height to assist access
- Viscoelastic foam cushion with two way stretch vapour permeable incontinence cover and non-slip base
- · Padded push handles to assist the attendant/porter
- Puncture proof wheels & castors
- Parking brakes as standard/optional disc brake for improved control
- Side/width adjustable polyurethane moulded footrests, 100kg capacity

#### Dimensions (w x d x h)

- Seat Width 610/710mm
- Seat Height 430mm
- Seat Depth 460mm
- Backrest Height 510mm

#### **Options / Accessories**

- 5X-0100-099-005 Butterfly Armrests, Pair
- 5X-0100-099-006 Angled Armrests, Pair
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0100-099-001 Universal Table
- 5X-0100-099-020 IV & O2 Holder, 100mm Dia

#### **Spares**

- 5X-0230-061-045 Viscoelastic Cushion with Incontinence Cover, 610mm
- 5X-0230-071-045 Viscoelastic Cushion with Incontinence Cover, 710mm





IV & O<sup>2</sup> Holder

Universal Table

Elevating Calf Support

sales@bristolmaid.com

| ltem              | Description          | Overall Dimensions (w x d x h) | Weight |
|-------------------|----------------------|--------------------------------|--------|
| 5X-0100-061-000/K | Minimaxx, 610mm      | 850 x 1100 x 990mm             | 27.0kg |
| 5X-0100-071-000/K | Minimaxx, 710mm      | 950 x 1100 x 990mm             | 27.5kg |
| 5X-0109-061-000/K | Minimaxx Disc, 610mm | 850 x 1100 x 990mm             | 28.0kg |
| 5X-0109-071-000/K | Minimaxx Disc, 710mm | 950 x 1100 x 990mm             | 28.5kg |

Bariatric Equipment

## **Bariatric Folding Wheelchair - Minimaxx Push Motor**

#### Features

- Maximum capacity 325kg
- · Powder coated chromium molybdenum steel frame
- · Rigid folding mechanism
- · Low seat height to assist access
- Viscoelastic foam cushion with two way stretch vapour permeable incontinence cover and non-slip base
- · Padded push handles to assist the attendant/porter
- · Puncture proof wheels & castors
- Parking brakes as standard/optional disc brake for improved control
- Side/width adjustable polyurethane moulded footrests, 100kg capacity
  Push motor to assist the attendant/porter, reduces the risk of sore
- shoulders & stress on wrists
- Compatible on flat surfaces up to 5 degrees
- Operated using a control panel fitted to the handle, offering drive & reverse modes with speed control from 0.6-3.1mph
- · One charge provides a range of 2.75miles
- Supplied with one removable battery and charger (battery has to be removed from the unit when charging)

#### Dimensions (w x d x h)

- Seat Width 610/710mm
- · Seat Height 430mm
- · Seat Depth 460mm
- Backrest Height 510mm

#### **Options / Accessories**

- 5X-0100-099-005 Butterfly Armrests, Pair
- 5X-0100-099-006 Angled Armrests, Pair
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0100-099-001 Universal Table
- 5X-0100-099-020 IV & O2 Holder, 100mm Dia

#### **Spares**

- 5X-0230-061-045 Viscoelastic Cushion with Incontinence Cover, 610mm
- 5X-0230-071-045 Viscoelastic Cushion with Incontinence Cover, 710mm

| Item              | Description                 | Overall Dimensions (w x d x h) | Weight |
|-------------------|-----------------------------|--------------------------------|--------|
| 5X-0110-061-000/K | Minimaxx, Push Motor, 610mm | 850 x 1100 x 990mm             | 50.0kg |
| 5X-0110-071-000/K | Minimaxx, Push Motor, 710mm | 950 x 1100 x 990mm             | 50.5kg |







Push Motor

Handle Mounted Control Panel

## Portering Chair - Bariatric, Rear Steer, Sliding Foot Rest





**Features** 

- · Maximum capacity 380kg
- Seamless & contoured upholstery using flame retardant vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- · X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles mounted to rear uprights
- Retractable diagonal seatbelt
- 200mm fixed wheels/200mm total stop braking castors (cushioned)
- We recommend that two members of staff are in attendance when used with a patient in excess of 190kg

#### Dimensions (w x d x h)

- Overall Dimensions 930 x 1000 x 1090mm
- Seat Width 740mm
- Seat Height 530mm
- Seat Depth 670mm
- Backrest Height 480mm

#### **Standard Colours**

- · Paint Grey
- Upholstery Bristol Blue

#### **Options / Accessories**

- G/500/CYL Cylinder Holder, CD/D
- G/CLS/500 Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/CARD/STATION Head of Station
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- 2LOKE011 Release Card

Alternative colours are available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Rear Steer Portering Chairs - Bariatric, Sliding Foot Rest in Aqua Marine - G/500/RS/VAQ

Vinyl - Ambla New Richmond in Colour



| Item         | Description                                 | Weight |
|--------------|---|--------|
| G/500/RS/VBB | Rear Steer, Sliding Foot Rest, Bristol Blue | 43.8kg |

+44 (0)1258 484455

## **Bariatric Mobile Shower / Commode Chairs**

#### **Features**

- Suitable for showering or toileting (supplied with a bucket & lid)
- Powder coated chromium molybdenum steel frame
- Available in two widths 610/710mm
- Rotating & removable straight or butterfly armrests to assist access
- Open back to allow easy sling access & assist positioning
- Polyurethane moulded seat & back
- Seat supplied with a removable lid
- Seat panel removable for cleaning
- Long open hole in seat (220mm at the widest point)
- Deep bucket with lid is removable from the front or back of the chair
- Side/width adjustable polyurethane moulded footrests, 100kg capacity
- 125mm swivel castors, 2 off braking
- Maximum capacity 325kg

#### Dimensions (w x d x h)

- Seat Width 610/710mm
- · Seat Height 550mm
- · Seat Depth 560mm
- Backrest Height 480mm

#### **Standard Colours**

- Paint White
- Polyurethane Blue

#### **Options / Accessories**

- 5X-0150-999-022 Seat Panel 610mm
- 5X-0150-999-023 Seat Panel 660/710mm
- 5X-0140-999-001 Commode Bucket & Lid
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0156-099-010 Butterfly Armrests, Pair
- 5X-0150-099-010/K Straight Armrests, Pair

| Item               | Description   | Overall Dimensions (w x d x h) | Weight |
|--------------------|---|--------------------------------|--------|
| Standard Armrest   |   |                                |        |
| 5X-0150-061-000    | Bariatric Mobile Shower / Commode Chair - Standard Armrests, 610mm  | 670 x 1120 x 990mm             | 27.0kg |
| 5X-0150-071-000    | Bariatric Mobile Shower / Commode Chair - Standard Armrests, 710mm  | 770 x 1120 x 990mm             | 27.5kg |
| Butterfly Armrests |   |                                |        |
| 5X-0151-061-000    | Bariatric Mobile Shower / Commode Chair - Butterfly Armrests, 610mm | 840 x 1120 x 990mm             | 27.5kg |
| 5X-0151-071-000    | Bariatric Mobile Shower / Commode Chair - Butterfly Armrests, 710mm | 940 x 1120 x 990mm             | 28.0kg |



+44 (0)1258 484455

CE



#### 5X-0142-061-000

## **Bariatric Shower / Commode Chair**

#### **Features**

- Suitable for showering or toileting (supplied with a bucket & lid) Powder coated chromium molybdenum steel frame
- Available in two widths 610/710mm
- Seat height adjustable 420-550mm
- · Long armrests to assist access
- Open back to allow easy sling access & assist positioning
- Polyurethane moulded seat & back
- Seat supplied with a removable lid
- Seat panel removable for cleaning
- Long open hole in seat (220mm at the widest point)
- Optional deep bucket is removable from the front or back of the chair
- Maximum capacity 325kg

#### Dimensions (w x d x h) Standard Colours

- Seat Width 610/710mmSeat Height 420-550mm
- Paint WhitePolyurethane Blue
- Seat Depth 400-520mm

#### **Options / Accessories**

- 5X-0150-999-022 Seat Panel 610mm
- 5X-0150-999-023 Seat Panel 660/710mm
- 5X-0140-999-001 Commode Bucket & Lid



| Item            | Description                              | Overall Dimensions (w x d x h) | Weight |
|-----------------|--|--------------------------------|--------|
| 5X-0142-061-000 | Bariatric Shower / Commode Chair - 610mm | 730 x 1010 x 790 - 940mm       | 15.0kg |
| 5X-0142-071-000 | Bariatric Shower / Commode Chair - 710mm | 830 x 1010 x 790 - 940mm       | 16.0kg |

## **Bariatric High Back Chairs**

#### **Features**

- Deep CMHR foam
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- · Fabric Fire retardant
- Black heavy duty steel five star base
- · Black heavy duty gas lift
- · Black height adjustable soft touch armrests
- Easy lift ratchet height adjustable back
- Independent seat & back mechanism with side tension control
- Maximum capacity 315kg

Replace suffix on the product code with any of the following to select your colour i.e. Bariatric High Back in Lime - 5660XL/VLI

#### Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage









## **Bariatric Chairs**

#### Features

- Deep CMHR foam
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- · Beech, varnished four leg frame
- Maximum capacity 315kg

Alternative colours are available at extra cost. Replace suffix on the product code with any of the following to select your colour i.e. Bariatric, Beech Frame in Lime - 5671/VLI

Vinyl - Ambla New Richmond in Colour



Fabric - Camira Advantage



| ltem | Description            | Overall Dimensions (w x d x h) | Seat Height | Weight |
|------|------------------------|--------------------------------|-------------|--------|
| 5671 | Bariatric, Beech Frame | 844 x 580 x 890mm              | 490mm       | 16.0kg |

## Other Products you may be interested in...

See Pages 176 & 178

Examination & Treatment Couch Two Section, Bariatric Examination & Treatment Couch Three Section, Bariatric







# **Bedside Cabinets** & Overbed Tables

- Polymer Bedside Cabinets 18
  - Bedside Cabinets 20
- Bedside Cabinet / Wardrobe Combinations 26
  - Polymer Overbed & Overchair Tables 30
  - Laminate Overbed & Overchair Tables 32

# eClean™



## **Polymer Bedside Cabinet - Easy Clean**

#### **Features**

- BS 3962 Part 6:1980 Assessment of resistance to mechanical damage
- BS EN 12720:2009 Furniture assessment of surface resistance to cold liquids
- BS EN 12721:2009 Furniture assessment of surface resistance to wet heat
- BS EN 12722:2009 Furniture assessment of surface resistance to dry heat
- BS 5852:2006 Methods of the test for assessment of the ignitability
- BS 4875-7:2006 Domestic & contract storage furniture
- Developed as part of Design Bugs Out & evaluated by the Department of Health
- Manufactured from recyclable polyethylene
- Open backless design ensures that all areas are accessible, whilst surfaces allow for easy & effective cleaning
- Distinct upper, middle & lower zones
  - Upper Nurse zone with integral rim to contain any spillages. Ideally positioned for staff to use while standing
  - Middle Patient zone provides storage area for patient use
  - Lower Patient zone provides shelf for washbowl, footwear etc.
- Drawers can be inserted from either side, allowing the cabinet to be positioned either side of the bed
- Upper & middle drawers can be fitted with an L15 125KHz RFID lock, providing secure storage for patient medicines / property\*
- 100mm polymer twin-wheel braking castors

\*cabinets with locks are supplied with one patient card - a programming card will be required per ward / directorate as per hospital policy

| Item   | Description                        | Weight |
|--------|------------------------------------|--------|
| EBC015 | No Locks                           | 29.0kg |
| EBC020 | Upper Drawer with Lock*            | 29.1kg |
| EBC025 | Middle Drawer with Lock*           | 29.1kg |
| EBC030 | Upper & Middle Drawers with Locks* | 29.2kg |

#### Dimensions (w x d x h)

- Overall Dimensions 480 x 665 x 1237mm
- Drawers (internal)
  - Upper 149 x 222 x 217mm
    Middle 338 x 491 x 197mm

  - Lower 338 x 491 x 220mm

### **Standard Colours**

- Body Grey White
- Drawers Spearmint

#### **Options / Accessories**

• RFID Lock Cards - page 25





BE

CC

ΒA

## **Bedside Cabinet - Single Upper Drawer**

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper & personal effects drawer
- Lower storage offered through...
  - Cupboard with double doors (adjustable shelf model dependant)
  - Drawer (recessed shelf model dependant)
- Available as either two-tone / single colour
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions: 485 x 485 x 905mm
- Internal Dimensions...
  - Upper Drawer 413 x 381 x 182mm
  - Personal Drawer 413 x 381 x 102mm
  - Lower Drawer 413 x 381 x 505 / 297mm
  - Cupboard 445 x 443 x 537 / 394mm

#### **Standard Colours**

- Two-Tone
  - · Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour Beech

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White | Beech     | Description  | Weight |
|------------|-----------|--|--------|
| BC1C/GW    | BC1C/BE   | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 34.0kg |
| BC1CPD/GW  | BC1CPD/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors  | 35.8kg |
| BC1D/GW    | BC1D/BE   | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 40.4kg |
| BC1DPD/GW  | BC1DPD/BE | Single Upper Drawer, Personal Drawer & Lower Drawer                      | 42.1kg |

NO

NC

sales@bristolmaid.com +

## **Bedside Cabinet - Double Upper Drawer**

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawers suitable for the storage of patient's medication or personal effects (1 x fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for either drawer
- Lower storage offered through...
  - Cupboard with double doors (adjustable shelf)
  - Drawer (recessed shelf)
- Available as either two-tone / single colour
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions 485 x 485 x 905mm
- Internal Dimensions...
  - Upper Drawer 180 x 381 x 182mm
  - Lower Drawer 413 x 381 x 505mm
  - Cupboard 445 x 443 x 537mm

#### **Standard Colours**

- Two-Tone
  - Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- · Single Colour Beech

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23



| <b>Grey White</b> | Beech   | Description  | Weight |
|-------------------|---------|--|--------|
| BC2C/GW           | BC2C/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 34.2kg |
| BC2D/GW           | BC2D/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf                    | 40.5kg |



Alternative Colours Available



+44 (0)1258 484455







BC3C/GW

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper section suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- · Rear side doors allow access to contents
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper section & personal effects drawer
- Lower storage offered through...
  - Cupboard with double doors (adjustable shelf model dependant)
  - Drawer (recessed shelf model dependant)
- · Available as either two-tone / single colour
- 50mm twin-wheel polypropylene castors

#### **Dimensions** (w x d x h)

- Overall Dimensions 485 x 485 x 905mm
- Internal Dimensions...
  - Upper Section 445 x 135 x 249mm
  - Rear Side Section 290 x 445 x 249mm
  - Personal Drawer 413 x 381 x 102mm
  - Lower Drawer 413 x 381 x 505 / 297mm
  - Cupboard 445 x 443 x 537 / 394mm

#### **Standard Colours**

Two-Tone

www.bristolmaid.com

- · Body Grey White
- Top / Drawer / Door Fronts Spearmint
- · Single Colour Beech

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White | Beech     | Description  | Weight |
|------------|-----------|--|--------|
| BC3C/GW    | BC3C/BE   | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 31.8kg |
| BC3CPD/GW  | BC3CPD/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors  | 33.5kg |
| BC3D/GW    | BC3D/BE   | Upper Section, Lower Drawer with Recessed Shelf                    | 38.1kg |
| BC3DPD/GW  | BC3DPD/BE | Upper Section, Personal Drawer & Lower Drawer                      | 39.8kg |



**Overbed Tables** 

**Bedside Cabinets &** 

## **Towel Holder**

#### Features

- · Self-adhesive hook
- Supplied for customer to fit



| Item | Description  | Weight |
|------|--------------|--------|
| TH   | Towel Holder | 0.1kg  |

## **Removable Plastic Top**

#### **Features**

- Raised edges on three sides reduces misplacement of objects
- Clips to the top of the bedside cabinet, & is removable to assist cleaning
- Not compatible with bedside cabinets fitted with wardrobes



RPT shown with BC1C/GW

| Item | Description           | Weight |
|------|-----------------------|--------|
| RPT  | Removable Plastic Top | 0.8kg  |

## **Electronic Push Button & RFID Locks Solutions**

Through either electronic push button or contactless RFID, we can offer a secure locking solution to suit your specific requirements.

#### **RFID Locks**

- Available as 125KHz or 13.56Mhz frequencies
- Easy programming using a dedicated card
- Three / four user levels offer differing levels of authority Administration / Master / Sub-Master / User
- Available with clamshell cards or fobs
- Compatible with Abloy CLIQ RFID chip inserted into the key







EBC Lock (L15)

PS Solo RFID Lock (PL)

ABLOY RFID Lock (AL)

|                        | EBC Lock (L15)          | PS Solo Lock (PL)       | ABLOY RFID Lock (AL)   |
|------------------------|-------------------------|-------------------------|--|
| Frequency              | 125KHz                  | 125KHz                  | 13,56Mhz   |
| Programming            | Programming Card        | Programming Card        | Programming Card   |
| Modes                  | Public & Private        | Public & Private        | Public & Private   |
| Audit Trail            | No                      | No                      | No   |
| User Levels            | 3 - Master, Staff, User | 3 - Master, Staff, User | 4 - Administration, Master, Staff, User                        |
| No. of Cards / Fobs    | 50 x Cards              | 50 x Cards              | 50 x Cards / Fobs  |
| Batteries              | 1 x CR123               | 1 x CR123               | 3 x AAA  |
| Openings (Battery Set) | 20,000                  | 20,000                  | 4,000  |
| Low Battery Warning    | Audible Alarm           | Audible Alarm           | Audible Alarm  |
| Low Battery Override   | Fails Unlocked          | Fails Unlocked          | Micro USB - External Power Source (power bank / phone charger) |
| Guarantee              | 1 Year                  | 2 Years                 | 1 Year   |
| Abloy CLIQ Compatible* | Yes                     | Yes                     | Yes  |

\*RFID chip inserted into the CLIQ key

#### **RFID Cards & Fobs**

- Compatible with 125KHz & 13,56Mhz frequencies
- Available as clamshell cards or fobs
- Can be fitted to a key ring
- A programming card is required for each hospital or ward as per hospital procedure

#### Dimensions (w x d x h)

- Clamshell Card 86 x 54 x 2mm
- Key Fob 27 x 35 x 5mm



#### Frequency - 13.56 mhz

| Item        | Description                 | Lock Compatibility |
|-------------|-----------------------------|--------------------|
| RFID1356/PC | 13,56KHz Programming Card   | AL                 |
| RFID1356/SC | 13,56KHz Staff Card (Green) | AL                 |
| RFID1356/UC | 13,56KHz User Card (Blue)   | AL                 |
| RFID1356/SF | 13,56KHz Staff Fob (Green)  | AL                 |
| RFID1356/UF | 13,56KHz User Fob (Blue)    | AL                 |

#### Frequency - 125KHz

| Item       | Description                                     | Lock Compatibility |
|------------|---|--------------------|
| RFID125/PC | 125KHz Programming Card - Staff<br>Pair (Black) | PL, L15            |
| RFID125/SC | 125KHz Staff Card (Green)                       | PL, L15            |
| RFID125/MC | 125KHz Master Card (Red)                        | PL, L15            |
| RFID125/UC | 125KHz User Card (Blue)                         | PL, L15            |

#### **Electronic Push Button Lock**

- Uses a 10-button keypad & # button
- Easy programming & code change
- Four user levels offer differing levels of authority...
  - Master / Technician 8 digits
    Staff / User 4 digits
- Available with a key override, allowing the lock to be opened immediately in the event of a low battery





CodeLock Key Overide (CLOR)

Electronic Push Button Lock (CL)

|                        | CodeLock Key Overide (CLOR)         | Electronic Push Button Lock (CL)    |
|------------------------|-------------------------------------|-------------------------------------|
| Programming            | 10-Button Keypad & # Button         | 10-Button Keypad & # Button         |
| Modes                  | Public & Private                    | Public & Private                    |
| Audit Trail            | No                                  | No                                  |
| User Levels            | 4 - Master, Technician, Staff, User | 4 - Master, Technician, Staff, User |
| User / Staff Codes     | Twenty                              | Single                              |
| Batteries              | 2 x AAA                             | 2 x AAA                             |
| Openings (Battery Set) | 100,000                             | 15,000                              |
| Low Battery Warning    | Red Flashing LED                    | Red Flashing LED                    |
| User / Staff Codes     | 9V PP3 Battery + Code               | 9V PP3 Battery + Code               |
| Guarantee              | 2 Years                             | 2 Years                             |
| Abloy CLIQ Compatible  | N/A                                 | N/A                                 |

## **Compact Bedside Cabinet / Wardrobe Combination - Single Upper Drawer**





Alternative Colours Available





BCC1CWL/GW

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper drawer
- Lower storage offered through...
  - · Cupboard with double doors (adjustable shelf)
  - Drawer (recessed shelf)
- · Available as either two-tone / single colour
- · Wardrobe supplied with a hanging rail, available right / left hand opening
- 50mm twin-wheel polypropylene castors

#### **Dimensions** (w x d x h)

- Overall Dimensions 633 x 489 x 1420mm
- Internal Dimensions...
  - Upper Drawer 328 x 381 x 182mm
  - Lower Drawer 328 x 381 x 505mm
  - Cupboard 360 x 443 x 537mm
  - Wardrobe 215 x 443 x 1420mm

#### **Standard Colours**

- Two Tone
  - Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour Beech

#### **Options / Accessories**

www.bristolmaid.com

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White | Beech      | Description  | Weight |
|------------|------------|--|--------|
| Right Hand |            |  |        |
| BCC1CWR/GW | BCC1CWR/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 50.6kg |
| BCC1DWR/GW | BCC1DWR/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 56.9kg |
| Left Hand  |            |  |        |
| BCC1CWL/GW | BCC1CWL/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 50.6kg |
| BCC1DWL/GW | BCC1DWL/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 56.9kg |

sales@bristolmaid.com

## **Bedside Cabinet / Wardrobe Combination - Single Upper Drawer**

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper drawer & personal effects drawer
- Lower storage offered through...
  - Cupboard with double doors (adjustable shelf)
- Drawer (recessed shelf model dependant)
- Available as either two-tone / single colour
- · Wardrobe supplied with a hanging rail, available right / left hand opening
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions: 718 x 489 x 1420mm
- Internal Dimensions...
  - Upper Drawer: 413 x 381 x 182mm
  - Personal Drawer: 413 x 381 x 102mm
  - Lower Drawer: 413 x 381 x 505 / 297mm
  - Cupboard: 445 x 443 x 537 / 394mm
  - Wardrobe: 215 x 443 x 1420mm

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White  | Beech       | Description  | Weight |
|-------------|-------------|--|--------|
| Right Hand  |             |  |        |
| BC1CWR/GW   | BC1CWR/BE   | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 54.5kg |
| BC1CPDWR/GW | BC1CPDWR/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors  | 56.2kg |
| BC1DWR/GW   | BC1DWR/BE   | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 60.8kg |
| BC1DPDWR/GW | BC1DPDWR/BE | Single Upper Drawer, Personal Drawer & Lower Drawer                      |        |
| Left Hand   |             |  |        |
| BC1CWL/GW   | BC1CWL/BE   | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 54.5kg |
| BC1CPDWL/GW | BC1CPDWL/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors  | 56.2kg |
| BC1DWL/GW   | BC1DWL/BE   | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 60.8kg |
| BC1DPDWL/GW | BC1DPDWL/BE | Single Upper Drawer, Personal Drawer & Lower Drawer                      | 62.5kg |

#### **Standard Colours**

- Two Tone
  - Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour- Beech



Alternative Colours Available CP GE RB LA VI BI GW GP BF CA NO NC BA CC BE





BC2CWL/GW

Alternative Colours Available



## **Bedside Cabinet / Wardrobe Combination - Double Upper Drawer**

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawers suitable for the storage of patient's medication or personal effects (1 x fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for either drawer
- Lower storage offered through...
  - · Cupboard with double doors (adjustable shelf)
  - Drawer (recessed shelf)
- Available as either two-tone / single colour
- Wardrobe supplied with a hanging rail, available right / left hand opening
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions: 718 x 489 x 1420mm
- Internal Dimensions...
  - Upper Drawer: 180 x 381 x 182mm
  - Lower Drawer: 413 x 381 x 505mm
  - Cupboard: 445 x 443 x 537mm
  - Wardrobe: 215 x 443 x 1420mm

#### **Standard Colours**

- Two Tone
  - Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour Beech

#### **Options / Accesssories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White | Beech     | Description  | Weight |
|------------|-----------|--|--------|
| Right Hand |           |  |        |
| BC2CWR/GW  | BC2CWR/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 55.0kg |
| BC2DWR/GW  | BC2DWR/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf                    | 61.3kg |
| Left Hand  |           |  |        |
| BC2CWL/GW  | BC2CWL/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 55.0kg |
| BC2DWL/GW  | BC2DWL/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf                    | 61.3kg |

sales@bristolmaid.com +44

## Bedside Cabinet / Wardrobe Combination - Upper Section & Rear Side Doors

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper section suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- · Rear side doors allow access to contents
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper section & personal effects drawer
- Lower storage offered through...
  - · Cupboard with double doors (adjustable shelf model dependant)
  - Drawer (recessed shelf model dependant)
- · Available as either two-tone / single colour
- Wardrobe supplied with a hanging rail, available right / left hand opening
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions: 718 x 489 x 1420mm
- Internal Dimensions...
  - Upper Section: 445 x 135 x 249mm
  - Rear Side Section: 290 x 445 x 249mm
  - Personal Drawer: 413 x 381 x 102mm
  - Lower Drawer: 413 x 381 x 505 / 297mm
  - Cupboard: 445 x 443 x 537 / 394mm
  - Wardrobe: 215 x 443 x 1420mm

#### **Standard Colours**

- Two Tone
  - · Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour Beech

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White  | Beech       | Description  | Weight |
|-------------|-------------|--|--------|
| Right Hand  |             |  |        |
| BC3CWR/GW   | BC3CWR/BE   | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 52.1kg |
| BC3CPDWR/GW | BC3CPDWR/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors  | 53.8kg |
| BC3DWR/GW   | BC3DWR/BE   | Upper Section, Lower Drawer with Recessed Shelf                    | 58.4kg |
| BC3DPDWR/GW | BC3DPDWR/BE | Upper Section, Personal Drawer & Lower Drawer                      | 60.1kg |
| Left Hand   |             |  |        |
| BC3CWL/GW   | BC3CWL/BE   | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 52.1kg |
| BC3CPDWL/GW | BC3CPDWL/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors  | 53.8kg |
| BC3DWL/GW   | BC3DWL/BE   | Upper Section, Lower Drawer with Recessed Shelf                    | 58.4kg |
| BC3DPDWL/GW | BC3DPDWL/BE | Upper Section, Personal Drawer & Lower Drawer                      | 60.1kg |



BC3CPDWR/GW

#### Alternative Colours Available





## **Overbed & Overchair Table - Polymer, Easy Clean**

#### **Features**

- BS 2483: 1977 Specification for overbed tables (selected tests)
- Finish Performance Tests (FIRA Standard 6250 2005)
  - BS 3962-6: 1980 Methods of test for finishes for wooden furniture
  - BS EN 12720: 2009 + A1: 2013 Assessment of surface resistance to cold liquids
  - BS EN 12721: 2009 + A1: 2013 Assessment of surface resistance to wet heat
  - BS EN 12722: 2009 + A1: 2013 Assessment of surface resistance to dry heat
- Height adjustable using a gas damper (760 1070mm)
- Safety feature prevents trapping of the patient between the bed / chair & table
- · Extruded aluminium column provides a smooth linear height adjustment
- Cantilever tubular base available as...
  - C shaped suitable for use with beds & chairs
  - · Straight leg suitable for use with beds
- Recyclable polypropylene one-piece top, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
- · Raised edges reduces misplacement of objects & contains spillages
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors
- Maximum load 25kg

#### Dimensions (w x d x h)

• Table Top - 880 x 440mm

#### **Standard Colours**

Top - Grey White

- Upper Extrusion Spearmint
- Lower Extrusion / Base Grey White

#### **Options / Accessories**

• B - Braking Castors (x 2)

| Item     | Description  | Overall Dimensions (w x d x h) | Weight |
|----------|--|--------------------------------|--------|
| OBP/SP/S | Overbed Table - Height Adjustable, Polymer, Easy Clean   | 895 x 440 x 760 - 1070mm       | 14.0kg |
| OCP/SP/S | Overchair Table - Height Adjustable, Polymer, Easy Clean | 918 x 440 x 760 - 1070mm       | 13.5kg |

sales@bristolmaid.com

+44 (0)1258 484455



## **Overbed & Overchair Table - Laminate**

#### **Features**

• BS 2483: 1977 Specification for overbed tables (selected tests)

 BS EN 438-2: 2016 High-pressure decorative laminates (HPL). Sheets based on thermosetting resins (usually called laminates). Determination of properties

- Height adjustable using a gas damper (735 1045mm)
- Safety feature prevents trapping of the patient between the bed / chair & table
- Extruded aluminium column provides a smooth linear height adjustment
- Cantilever tubular base available as...
  - · C shaped suitable for use with beds & chairs
  - · Straight leg suitable for use with beds
- · Laminate top with varnished edges available in three configurations...
  - Flat top
  - One varnished beech lip
  - Three varnished beech lips
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors
- · Maximum load 25kg

#### Dimensions (w x d x h)

- Table Top 880 x 395mm
- Overbed Table 895 x 440 x 735 -1045mm
- Overchair Table 918 x 440 x 735 -1045mm

#### **Standard Colours**

- Top Grey White or Beech
- Upper Extrusion Spearmint or Grey White
- · Lower Extrusion / Base Grey White

#### **Options / Accessories**

• B - Braking Castors (x 2)

#### **Overbed Tables**

| Grey White   | Beech        | Description | Weight |
|--------------|--------------|-------------|--------|
| OBL/GW0/SP/S | OBL/BE0/GW/S | Flat Top    | 14.5kg |
| OBL/GW1/SP/S | OBL/BE1/GW/S | One Lip     | 15.4kg |
| OBL/GW3/SP/S | OBL/BE3/GW/S | Three Lips  | 15.8kg |



sales@bristolmaid.com

OBL/GW1/SP/S

+44 (0)1258 484455

HEIGHT



www.bristolmaid.com



# **Caretray Trolleys & Racks**

- Mild Steel, Low Level Trolleys 50
- Mild Steel, Standard Level Trolleys 52
  - Mild Steel, High Level Trolleys 54
    - Mild Steel, Mobile Racks 56
    - Mild Steel, Static Racks 58
      - Accessories 60

- Features 36
- Disassembly & Cleaning 38
- Stainless Steel, Low Level Trolleys 40
- Stainless Steel, Standard Level Trolleys 42
  - Stainless Steel, High Level Trolleys 44
    - Stainless Steel, Mobile Racks 46
    - Stainless Steel, Static Racks 48

### **Features**

### **Carerails**

 Fitted to both sides of the trolley, compatible with a wide range of accessories

### Frame

- · Fully welded eliminating gaps
- 304 grade stainless steel or mild steel that is powder coated in a durable epoxy polyester paint

### **Side Panels**

- One-piece moulding
- Smooth surfaces & rounded edges assist cleaning (compatible with steam cleaning)
- Runners with positive stops prevent trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)

### **Easy Clean**

- One-piece mouldings
- Smooth surfaces & rounded edges assist cleaning (compatible with steam cleaning)
- Quick & easy tool-less disassembly assists thorough or deep cleaning programmes



### Тор

- One-piece moulding
- Smooth surfaces & rounded edges assist cleaning (compatible with steam cleaning)
- Raised edges on three sides provides defined grab points to act as push handles, whilst an integral rim at the front contains spillages
- Maximum load 30kg

### **Trays & Dividers**

- Translucent & impact-resistant (100 / 150mm depths)
- Smooth surfaces & rounded edges assist cleaning (compatible with steam cleaning & able to withstand high temperature machine washing)
- Trays feature an integral handle with plastic clip on label cover, (suitable for use with our perforated coloured label strips)
- Trays can be subdivided in both length & width using simple slot in dividers
- Maximum load 5kg

### **Castors & Ring Buffers**

- 100mm pressed steel or polymer castors (model dependant)
- Two braking castors fitted as standard
- Ring buffers offer protection to both the trolley & fabric of the building
- Smooth surfaces assist cleaning

Racks

øð

### Accessories

The Caretray Trolley has been developed specifically to meet the needs of Healthcare Professionals across a wide variety of disciplines including Minor Treatment, A&E. Theatre. Phlebotomy, ITU & Outpatients. Through a wide range of accessories, you can configure your trolley to suit your specific needs.

### Locking & Security - Page 60

- · Solid grade laminate / clear acrylic doors
- · Range of locking options
- Security clamp & seals to safeguard contents

### Laptop Holders - Page 74

- · Suitable for use with Electronic Prescribing & Medicines Administration (EPMA) & Electronic Patient Records (EPR)
- Offers height, rotation & angle adjustment
- Laptop can be secured using a security clamp



### **Overbridge & Tilt Bins - Page 66**

- Provides high level visible storage & assists organisation
- Tilt-bins available with 4, 5 or 6 clear compartments
- Compartments tilt to 45 degrees, to assist access of contents



### **Defibrillator Platform - Page 67**

- · Compatible with all popular defibrillators
- · Adjustable, allowing compatibility with all popular models
- Platform rotates 360 degrees, whilst the arm swings 90 degrees from the worksurface

### Castors & Ring Buffers - Page 76

- Range of polymer castors
- Available in various braking & non-braking configurations
- Ring buffers offer protection to both the trolley & fabric of the building



### **Carerail Mounted - Page 68**

- · An extensive range of accessories can be attached/suspended from the Carerails
- Allows the trolley to be configured to meet specific user requirements
- Accessories range from sharps & glove box holders, catheter racks & waste bins, through to work flaps & suction pump shelves



### Trays, Shelves, Dividers & Labels - Page 77

- · Additional dividers assist segregation
- · Coloured labels ensure that trays can be easily identified
- · Flat shelves provide an alternative to trays



sales@bristolmaid.com

## **Disassembly & Cleaning**

The design of the trolley allows for quick & easy tool-less disassembly providing access for all areas when conducting thorough or deep cleaning programmes. The simple guide below illustrates the disassembly & re-assembly process...






Remove the trays from the unit



### **Remove Top**

The top is held in place with two thumb screws, located on the underside of the top. Once removed, the top can be removed with a firm tug



### **Remove Back Panel** (if fitted)

Once the top has been removed, the back panel can be lifted clear of the cross members on the rear of the frame



### **Remove Side Panels**

The side panels are removed by lifting off the positioning inserts. Double height trolleys require thumb screws at the top of the lower panel to be removed.



### Clean

The frame, side panels & top can now be cleaned



### Reassembly

The trolley is re-assembled by repeating the steps in reverse

### Caretray Trolleys - Stainless Steel, Low Level, Single & Double Column





CTS1/08/3S1D - shown with optional CT/WF, CT/ALC & GMBH/15



Easy disassembly & cleaning - page 38

#### **Features**

- · Low level, single & double column easy clean procedure trolleys
- Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- · 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 807mm
  - Double Column 1005 x 550 x 807mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 782mm
- Double Column 915 x 455 (internal) at height 782mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

 Item
 Description
 Weight

 CTS1/08/1S2D
 1 x 100mm, 2 x 150mm Trays
 21.2kg

 CTS1/08/3D
 3 x 150mm Trays
 21.5kg

 CTS1/08/3D1
 3 x 100mm, 1 x 150mm Trays
 22.3kg

 CTS1/08/4S
 4 x 100mm Trays
 21.9kg

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67

- Carerail Mounted page 68
  Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77



CTS2/08/4S3D

| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTS2/08/3S4D | 3 x 100mm, 4 x 150mm Trays | 37.3kg |
| CTS2/08/4S3D | 4 x 100mm, 3 x 150mm Trays | 36.9kg |
| CTS2/08/6D   | 6 x 150mm Trays            | 36.5kg |
| CTS2/08/6S2D | 6 x 100mm, 2 x 150mm Trays | 38.1kg |
| CTS2/08/7S1D | 7 x 100mm, 1 x 150mm Trays | 37.7kg |
| CTS2/08/8S   | 8 x 100mm Trays            | 37.3kg |

CTS2/08/4S3D - shown with optional CT/IV, GMBH/15 & CT/ALC

### Caretray Trolleys - Stainless Steel, Standard Level, Single & Double Column



| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTS1/10/1S3D | 1 x 100mm, 3 x 150mm Trays | 27.9kg |
| CTS1/10/3S2D | 3 x 100mm, 2 x 150mm Trays | 29.0kg |
| CTS1/10/4D   | 4 x 150mm Trays            | 28.3kg |
| CTS1/10/4S1D | 4 x 100mm, 1 x 150mm Trays | 28.7kg |
| CTS1/10/6S   | 6 x 100mm                  | 29.8kg |



CTS1/10/3S2D - shown with optional CT/IV, CT/CR2, GMBH/15 & CT/ALC



Easy disassembly & cleaning - page 38

#### **Features**

- Standard level, single & double column easy clean procedure trolleys
- Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 993mm
  - Double Column 1005 x 550 x 993mm
- Plastic Top
  - Single Column 455 x 455 (internal) at height 968mm
  - Double Column 915 x 455 (internal) at height 968mm

#### Trays

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

aid.com +44 (0)1258 484455

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
- Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77



CTS2/10/6S4D

| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTS2/10/3S6D | 3 x 100mm, 6 x 150mm Trays | 46.6kg |
| CTS2/10/4S5D | 4 x 100mm, 5 x 150mm Trays | 46.2kg |
| CTS2/10/5S4D | 5 x 100mm, 4 x 150mm Trays | 45.9kg |
| CTS2/10/6S4D | 6 x 100mm, 4 x 150mm Trays | 47.3kg |
| CTS2/10/8D   | 8 x 150mm Trays            | 45.8kg |
| CTS2/10/8S2D | 8 x 100mm, 2 x 150mm Trays | 46.6kg |



CTS2/10/6S4D - shown with optional CT/IV, CT/CR2, GMBH/15 & CT/ALC



#### CTS1/14/4S3D

| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTS1/14/2S5D | 2 x 100mm, 5 x 150mm Trays | 36.2kg |
| CTS1/14/3S4D | 3 x 100mm, 4 x 150mm Trays | 35.8kg |
| CTS1/14/4S3D | 4 x 100mm, 3 x 150mm Trays | 35.4kg |
| CTS1/14/5S3D | 5 x 100mm, 3 x 150mm Trays | 37.8kg |
| CTS1/14/6D   | 6 x 150mm Trays            | 35.0kg |
| CTS1/14/8S1D | 8 x 100mm, 1 x 150mm Trays | 37.8kg |



CTS1/14/4S3D - shown with optional CT/DH/A4



Easy disassembly & cleaning - page 38

#### Features

Caretray Trolleys - Stainless Steel, High Level, Single & Double Column

- High level, single & double column easy clean procedure trolleys
- Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 1410mm
  - Double Column 1005 x 550 x 1410mm
- Plastic Top
  - Single Column 455 x 455 (internal) at height 1385mm
- Double Column 915 x 455 (internal) at height 1385mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

- Locking & Security page 60
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



CTS2/14/10S6D - shown with optional CT/DH/A4

| Item          | Description                 | Weight |
|---------------|-----------------------------|--------|
| CTS2/14/4S10D | 2 x 100mm, 10 x 150mm Trays | 60.6kg |
| CTS2/14/6S8D  | 6 x 100mm, 8 x 150mm Trays  | 62.8kg |
| CTS2/14/8S6D  | 8 x 100mm, 6 x 150mm Trays  | 63.0kg |
| CTS2/14/10S6D | 10 x 100mm, 6 x 150mm Trays | 65.0kg |
| CTS2/14/12D   | 12 x 150mm Trays            | 61.3kg |
| CTS2/14/16S2D | 16 x 100mm, 2 x 150mm Trays | 66.6kg |



CTS2/14/10S6D





#### CRSM1/18/6S4D

| Item          | Description                | Weight |
|---------------|----------------------------|--------|
| CRSM1/18/2S6D | 2 x 100mm, 6 x 150mm Trays | 41.7kg |
| CRSM1/18/3S6D | 3 x 100mm, 6 x 150mm Trays | 44.1kg |
| CRSM1/18/4S5D | 4 x 100mm, 5 x 150mm Trays | 42.8kg |
| CRSM1/18/5S5D | 5 x 100mm, 5 x 150mm Trays | 44.3kg |
| CRSM1/18/6S4D | 6 x 100mm, 4 x 150mm Trays | 43.9kg |
| CRSM1/18/8D   | 8 x 150mm Trays            | 43.3kg |

Easy disassembly & cleaning - page 38

### Caretray Racks - Stainless Steel, Mobile, Single & Double Column

#### **Features**

- Mobile, single & double column easy clean racks
- · Quick & easy disassembly, assists thorough & deep cleaning
- · Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recvclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each trav
- · Polyurethane ring buffers offer protection to both the rack & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 1782mm
  - Double Column 1005 x 550 x 1782mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 1757mm
- Double Column 915 x 455 (internal) at height 1757mm
- Travs
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

sales@bristolmaid.com

- Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



CRSM2/18/8S10D -shown with optional C004

| Item            | Description                  | Weight |
|-----------------|------------------------------|--------|
| CRSM2/18/4S12D  | 4 x 100mm, 12 x 150mm Trays  | 72.1kg |
| CRSM2/18/6S12D  | 6 x 100mm, 12 x 150mm Trays  | 77.1kg |
| CRSM2/18/8S10D  | 8 x 100mm, 10 x 150mm Trays  | 75.4kg |
| CRSM2/18/10S10D | 10 x 100mm, 10 x 150mm Trays | 78.4kg |
| CRSM2/18/12S8D  | 12 x 100mm, 8 x 150mm Trays  | 77.7kg |
| CRSM2/18/16D    | 16 x 150mm Trays             | 74.6kg |







#### Features

Caretray Racks - Stainless Steel, Static, Single & Double Column

- Static, single & double column easy clean racks
- Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- · Adjustable feet allow levelling on uneven floors

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 536 x 536 x 1685mm
  - Double Column 995 x 536 x 1685mm
- Plastic Top
  - Single Column 455 x 455 (internal) at height 1645mm
- Double Column 915 x 455 (internal) at height 1645mm

#### Trays

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray -355 x 405 x 150mm (internal)

CRSS1/16/6S4D

| Item          | Description                | Weight |
|---------------|----------------------------|--------|
| CRSS1/16/2S6D | 2 x 100mm, 6 x 150mm Trays | 38.4kg |
| CRSS1/16/3S6D | 3 x 100mm, 6 x 150mm Trays | 40.8kg |
| CRSS1/16/4S5D | 4 x 100mm, 5 x 150mm Trays | 39.5kg |
| CRSS1/16/5S5D | 5 x 100mm, 5 x 150mm Trays | 41.0kg |
| CRSS1/16/6S4D | 6 x 100mm, 4 x 150mm Trays | 40.6kg |
| CRSS1/16/8D   | 8 x 150mm Trays            | 40.0kg |

Easy disassembly & cleaning - page 38

www.bristolmaid.com

• Trays, Shelves, Dividers & Labels - page 77

CRSS2/16/8S10D

| Item            | Description                  | Weight |
|-----------------|------------------------------|--------|
| CRSS2/16/4S12D  | 4 x 100mm, 12 x 150mm Trays  | 71.1kg |
| CRSS2/16/6S12D  | 6 x 100mm, 12 x 150mm Trays  | 74.0kg |
| CRSS2/16/8S10D  | 8 x 100mm, 10 x 150mm Trays  | 72.3kg |
| CRSS2/16/10S10D | 10 x 100mm, 10 x 150mm Trays | 75.3kg |
| CRSS2/16/12S8D  | 12 x 100mm, 8 x 150mm Trays  | 74.6kg |
| CRSS2/16/16D    | 16 x 150mm Trays             | 71.5kg |



### Caretray Trolleys - Mild Steel, Low Level, Single & Double Column





CTM1/08/3S1D - shown with optional CT/WF, SBH/15 & CT/ALC



Easy disassembly & cleaning - page 38

#### **Features**

- Low level, single & double column easy clean procedure trollevs
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows trav insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant travs & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- 100 / 150mm deep travs, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each trav
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 807mm
  - Double Column 1005 x 550 x 807mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 782mm
- Double Column 915 x 455 (internal) at height 782mm
- Travs
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### Standard Colours

Paint - Spearmint

Item Description Weight CTM1/08/1S2D 1 x 100mm, 2 x 150mm Trays 20.9kg CTM1/08/3D 3 x 150mm Trays 21.2kg CTM1/08/3S1D 3 x 100mm, 1 x 150mm Trays 22.0kg CTM1/08/4S 4 x 100mm Trays 21.6kg

sales@bristolmaid.com

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67

- Carerail Mounted page 68
  Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77

CTM2/08/4S3D

| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTM2/08/3S4D | 3 x 100mm, 4 x 150mm Trays | 36.8kg |
| CTM2/08/4S3D | 4 x 100mm, 3 x 150mm Trays | 36.4kg |
| CTM2/08/6D   | 6 x 150mm Trays            | 36.0kg |
| CTM2/08/6S2D | 6 x 100mm, 2 x 150mm Trays | 37.6kg |
| CTM2/08/7S1D | 7 x 100mm, 1 x 150mm Trays | 37.2kg |
| CTM2/08/8S   | 8 x 100mm Trays            | 36.8kg |

CTM2/08/4S3D - shown with optional CT/WF, SBH/15, CT/ALC

### Caretray Trolleys - Mild Steel, Standard Level, Single & Double Column





#### Easy disassembly & cleaning - page 38

#### **Features**

- Standard level, single & double column easy clean procedure trolleys
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 993mm
  - Double Column 1005 x 550 x 993mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 968mm
- Double Column 915 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

Paint - Spearmint

Description Weight Item CTM1/10/1S3D 1 x 100mm, 3 x 150mm Trays 27.5kg CTM1/10/3S2D 3 x 100mm, 2 x 150mm Trays 28.6kg CTM1/10/4D 4 x 150mm Travs 27.9kg CTM1/10/4S1D 4 x 100mm, 1 x 150mm Trays 28.3kg CTM1/10/6S 6 x 100mm 29.4kg

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
- Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77

CTM2/10/6S4D

| Item   | Description   | Weight                           |
|--|---|----------------------------------|
| CTM2/10/3S6D   | 3 x 100mm, 6 x 150mm Trays  | 46.0kg                           |
| CTM2/10/4S5D   | 4 x 100mm, 5 x 150mm Trays  | 45.6kg                           |
| CTM2/10/5S4D   | 5 x 100mm, 4 x 150mm Trays  | 45.3kg                           |
| CTM2/10/6S4D   | 6 x 100mm, 4 x 150mm Trays  | 46.7kg                           |
| CTM2/10/8D   | 8 x 150mm Trays   | 45.2kg                           |
| CTM2/10/8S2D   | 8 x 100mm, 2 x 150mm Trays  | 46.0kg                           |
| CTM2/10/4S5D<br>CTM2/10/5S4D<br>CTM2/10/6S4D<br>CTM2/10/8D | 4 x 100mm, 5 x 150mm Trays<br>5 x 100mm, 4 x 150mm Trays<br>6 x 100mm, 4 x 150mm Trays<br>8 x 150mm Trays | 45.6k<br>45.3k<br>46.7k<br>45.2k |

and the second

CTM2/10/6S4D - shown with optional CT/IV, CT/CR2, GMBH/15 & CT/ALC



sales@bristolmaid.com

### High level, single & double column easy clean procedure trolleys Quick & easy disassembly, assists thorough & deep cleaning

- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly) Recvclable polypropylene one-piece top & side panels, smooth
- surfaces & radiused edges assists cleaning (compatible with steam cleaning)
- Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
- Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant travs & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- 100 / 150mm deep travs, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each trav
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 1410mm
  - Double Column 1005 x 550 x 1410mm
- Plastic Top

**Features** 

- Single Column 455 x 455 (internal) at height 1385mm
- Double Column 915 x 455 (internal) at height 1385mm
- Travs
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### Standard Colours

Paint - Spearmint

Item Description CTM1/14/2S5D 2 x 100mm, 5 x 150mm Trays 36.2kg CTM1/14/3S4D 3 x 100mm, 4 x 150mm Trays 35.8kg 4 x 100mm, 3 x 150mm Trays CTM1/14/4S3D 35.4kg CTM1/14/5S3D 5 x 100mm. 3 x 150mm Travs 37.3kg CTM1/14/6D 6 x 150mm Trays 35.0kg CTM1/14/8S1D 8 x 100mm, 1 x 150mm Trays 37.2kg Easy disassembly & cleaning - page 38

- Locking & Security page 60

- Carerail Mounted Accessories page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



CTM2/14/10S6D - shown with optional CTD2/14SS (Double door fitted with Security Seal Lock)

| Item          | Description                 | Weight |
|---------------|-----------------------------|--------|
| CTM2/14/4S10D | 2 x 100mm, 10 x 150mm Trays | 60.6kg |
| CTM2/14/6S8D  | 6 x 100mm, 8 x 150mm Trays  | 62.8kg |
| CTM2/14/8S6D  | 8 x 100mm, 6 x 150mm Trays  | 62.3kg |
| CTM2/14/10S6D | 10 x 100mm, 6 x 150mm Trays | 65.0kg |
| CTM2/14/12D   | 12 x 150mm Trays            | 61.3kg |
| CTM2/14/16S2D | 16 x 100mm, 2 x 150mm Trays | 66.6kg |







CRMM1/18/6S4D

Easy disassembly & cleaning - page 38

|  | Item          | Description                | Weight |
|--|---------------|----------------------------|--------|
|  | CRMM1/18/2S6D | 2 x 100mm, 6 x 150mm Trays | 41.2kg |
|  | CRMM1/18/3S6D | 3 x 100mm, 6 x 150mm Trays | 43.6kg |
|  | CRMM1/18/4S5D | 4 x 100mm, 5 x 150mm Trays | 42.3kg |
|  | CRMM1/18/5S5D | 5 x 100mm, 5 x 150mm Trays | 43.8kg |
|  | CRMM1/18/6S4D | 6 x 100mm, 4 x 150mm Trays | 43.4kg |
|  | CRMM1/18/8D   | 8 x 150mm Trays            | 42.8kg |

#### **Features**

- Mobile, single & double column easy clean racks
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- · 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the rack & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 1782mm
  - Double Column 1005 x 550 x 1782mm
- · Plastic Top
  - Single Column 455 x 455 (internal) at height 1757mm
- Double Column 915 x 455 (internal) at height 1757mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

· Paint - Spearmint

sales@bristolmaid.com +

- Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



CRMM2/18/8S10D

| Item            | Description                  | Weight |
|-----------------|------------------------------|--------|
| CRMM2/18/4S12D  | 4 x 100mm, 12 x 150mm Trays  | 73.3kg |
| CRMM2/18/6S12D  | 6 x 100mm, 12 x 150mm Trays  | 76.3kg |
| CRMM2/18/8S10D  | 8 x 100mm, 10 x 150mm Trays  | 74.6kg |
| CRMM2/18/10S10D | 10 x 100mm, 10 x 150mm Trays | 77.6kg |
| CRMM2/18/12S8D  | 12 x 100mm, 8 x 150mm Trays  | 76.9kg |
| CRMM2/18/16D    | 16 x 150mm Trays             | 73.8kg |





| Item          | Description                | Weight |
|---------------|----------------------------|--------|
| CRMS1/16/2S6D | 2 x 100mm, 6 x 150mm Trays | 38.4kg |
| CRMS1/16/3S6D | 3 x 100mm, 6 x 150mm Trays | 39.5kg |
| CRMS1/16/4S5D | 4 x 100mm, 5 x 150mm Trays | 39.5kg |
| CRMS1/16/5S5D | 5 x 100mm, 5 x 150mm Trays | 41.0kg |
| CRMS1/16/6S4D | 6 x 100mm, 4 x 150mm Trays | 40.6kg |
| CRMS1/16/8D   | 8 x 150mm Trays            | 39.6kg |
|               |                            |        |



Easy disassembly & cleaning - page 38

### Caretray Racks - Mild Steel, Static, Single & Double Column

#### **Features**

- Mobile, single & double column easy clean racks
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- · Adjustable feet allow levelling on uneven floors

#### **Dimensions**

- Overall Dimensions
  - Single Column 536 x 536 x 1670-1685mm
  - Double Column 995 x 536 x 1670-1685mm
- Plastic Top
  - Single Column 455 x 455 (internal) at height 1645mm
  - Double Column 915 x 455 (internal) at height 1645mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

· Paint - Spearmint

sales@bristolmaid.com

www.bristolmaid.com

• Trays, Shelves, Dividers & Labels - page 77



| Item            | Description                  | Weight |
|-----------------|------------------------------|--------|
| CRMS2/16/4S12D  | 4 x 100mm, 12 x 150mm Trays  | 70.2kg |
| CRMS2/16/6S12D  | 6 x 100mm, 12 x 150mm Trays  | 73.2kg |
| CRMS2/16/8S10D  | 8 x 100mm, 10 x 150mm Trays  | 72.3kg |
| CRMS2/16/10S10D | 10 x 100mm, 10 x 150mm Trays | 75.3kg |
| CRMS2/16/12S8D  | 12 x 100mm, 8 x 150mm Trays  | 74.6kg |
| CRMS2/16/16D    | 16 x 150mm Trays             | 71.5kg |
|                 |                              |        |





### Locking

#### **Features**

- · Doors available as either
  - Solid Grade Laminate Single / double column units
  - Clear Acrylic Single column units only
- Back panels are solid grade laminate, secures the contents of the trolley (removable to assist cleaning)
- Available for both single & double column units (mild & stainless steel)
  - Single Column Units Single door
  - Double Column Units Double door, features a return to ensure that only one lock is required to secure both doors
- 270-degree opening hinges allow the doors to fold flat against the sides of the trolley
- Four locking options
  - CAM Lock with Wing Handle Solid grade laminate / clear acrylic
  - Electronic Push Button Lock Solid grade laminate / clear acrylic
  - High Security Bolt Lock Solid grade laminate only
  - Security Seal Lock Solid grade laminate / clear acrylic
- Factory fitted



Door Open - CTM1/10/3S2D/CTD1/10EPB

SS

EPB

### **CAM Lock with Wing Handle**

#### Features

- Supplied with two keys
  Wing handle can operate independently of the key
  Locks keyed to differ, alike on request



CTM1/08/3S1D/CTD1/08CA

| Description                   | Solid Grade Laminate | Clear Acrylic | Weight |
|-------------------------------|----------------------|---------------|--------|
| Low Level, Single Column      | CTD1/08CA            | CTDC1/08CA    | 2.7kg  |
| Low Level, Double Column      | CTD2/08CA            | -             | 5.3kg  |
| Standard Level, Single Column | CTD1/10CA            | CTDC1/10CA    | 3.5kg  |
| Standard Level, Double Column | CTD2/10CA            | -             | 6.8kg  |
| High Level, Single Column     | CTD1/14CA            | CTDC1/14CA    | 5.3kg  |
| High Level, Double Column     | CTD2/14CA            | -             | 10.3kg |



CTM1/10/3S2D/CTDC1/10CA



### **Electronic Push Button Lock**

#### **Features**

- Uses a 10-button keypad
- Easy programming & code change
- Four operating codes with differing levels of authority Master / Sub-Master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code



CTM1/08/3S1D/CTD1/08EPB

| Description                   | Solid Grade Laminate | Clear Acrylic | Weight |
|-------------------------------|----------------------|---------------|--------|
| Low Level, Single Column      | CTD1/08EPB           | CTDC1/08EPB   | 2.9kg  |
| Low Level, Double Column      | CTD2/08EPB           | -             | 5.5kg  |
| Standard Level, Single Column | CTD1/10EPB           | CTDC1/10EPB   | 3.7kg  |
| Standard Level, Double Column | CTD2/10EPB           | -             | 7.0kg  |
| High Level, Single Column     | CTD1/14EPB           | CTDC1/14EPB   | 5.5kg  |
| High Level, Double Column     | CTD2/14EPB           | -             | 10.5kg |



www.bristolmaid.com

KitLock

### **High Security Bolt-Lock**

#### Features

Supplied with two keys



CTD1/10HSB shown on CTM1/10/3S2D





| Description                   | Solid Grade Laminate | Weight |
|-------------------------------|----------------------|--------|
| Low Level, Single Column      | CTD1/08HSB           | 2.9kg  |
| Low Level, Double Column      | CTD2/08HSB           | 5.6kg  |
| Standard Level, Single Column | CTD1/10HSB           | 3.6kg  |
| Standard Level, Double Column | CTD2/10HSB           | 7.1kg  |
| High Level, Single Column     | CTD1/14HSB           | 5.3kg  |
| High Level, Double Column     | CTD2/14HSB           | 10.6kg |

www.bristolmaid.com

### **Security Seal Lock**

#### **Features**

• Accepts standard plastic security seals (order separately)



| Description                   | Solid Grade Laminate | Clear Acrylic | Weight |
|-------------------------------|----------------------|---------------|--------|
| Low Level, Single Column      | CTD1/08SS            | CTDC1/08SS    | 2.6kg  |
| Low Level, Double Column      | CTD2/08SS            |               | 5.2kg  |
| Standard Level, Single Column | CTD1/10SS            | CTDC1/10SS    | 3.7kg  |
| Standard Level, Double Column | CTD2/10SS            |               | 6.7kg  |
| High Level, Single Column     | CTD1/14SS            | CTDC1/14SS    | 5.2kg  |
| High Level, Double Column     | CTD2/14SS            |               | 10.2kg |

CTM1/10/3S2D/CTDC1/10SS

### **Security Clamp**


#### **Features**

- Allows the trolley to be secured to the 'fabric of the building'
- Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door of the trolley when closed & locked
- Recommended wall fitting height approx. 1000mm

### **Security Seals**

#### **Features**

- Snap on plastic security seals
- Pack of 100





#### Important Note - Security Clamps

It is a requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area & broken into, but also encourages good nursing practice, giving the trolley a "permanent home".

'Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB-March 2005"

| Item    | Description    | Cable Length | Weight |
|---------|----------------|--------------|--------|
| MT/1/6F | Security Clamp | 950mm        | 0.75kg |

| ltem     | Description    |
|----------|----------------|
| 2PLCO026 | Security Seals |

www.bristolmaid.com

sales@bristolmaid.com

+44 (0)1258 484455

### Overbridge

#### **Features**

- Attaches to the rear uprights of the trolley (increases overall height by 625mm)
- Available as either painted mild steel or self-colour stainless steel
- Provides a fixing location for tilt bins (order separately)
- Compatible with single column trolleys only
- Factory fitted only

#### **Standard Colours**

- Paint Grey White
- Stainless Steel Self Colour



| Item     | Description                  | Weight |
|----------|------------------------------|--------|
| CT/OB/MS | Overbridge - Mild Steel      | 2.0kg  |
| CT/OB/SS | Overbridge - Stainless Steel | 2.0kg  |

# Tilt Bins

- Available with 4, 5 or 6 compartments
- Clear bins provide visible storage & assists organisation
- Bins tilt to 45° allowing easy access of contents, & can be removed to assist refiling / cleaning
- Label holder for easy identification
- Retaining bar secures bins in place, preventing from tipping forward during transit
- Supplied with a mild / stainless steel mounting

#### Dimensions (w x d x h)

- TBIN4 143 x 127 x 190mm
- TBIN5 108 x 95 x 146mm
- TBIN6 89 x 67 x 95mm

#### **Standard Colours**

- Plastic Cream
- Paint Grey White
- Stainless Steel Self Colour



Label Holder

Tilts to 45°

| Mild Steel | Stainless Steel | Description      | Dimensions (w x d x h) | Weight |
|------------|-----------------|------------------|------------------------|--------|
| TBIN4/MS   | TBIN4/SS        | Four Compartment | 600 x 171 x 207mm      | 2.3kg  |
| TBIN5/MS   | TBIN5/SS        | Five Compartment | 600 x 135 x 164mm      | 1.6kg  |
| TBIN6/MS   | TBIN6/SS        | Six Compartment  | 600 x 93 x 113mm       | 0.9kg  |

Racks

+44 (0)1258 484455 www.br

## **Adjustable Defibrillator Platform**

#### **Features**

- Compatible with all popular defibrillators
- Defibrillator supported on a 365mm wide platform & is held securely in place using a clamp (adjusts using a handwheel 205 380mm)
- Platform mounted on an arm that swings 90 degrees from the worksurface & can be locked in position using a tensioning knob (rotates 360 degrees)
- Base of platform sits 130mm above the top of the trolley
- Fitted to the rear left hand upright of the trolley
- Maximum capacity 6kg
- Factory fitted

### **Standard Colours**

• Paint - Grey White

### Dimensions (w x d x h)

• 477 x 365 x 642mm (increases overall height by 170mm)

| Item     | Description                       | Weight |
|----------|-----------------------------------|--------|
| CT/DEFIB | Adjustable Defibrillator Platform | 5.5kg  |

MAXIMUM SAFE

### **Carerail Mounted**

### **Alcohol Gel Holder**

#### **Features**

- Accepts bottles with 90 x 90mm square base
- Fitted with a 15mm hook

#### **Standard Colours**

• Paint - Grey White

#### Dimensions (w x d x h)

• 95 x 110 x 115mm



| Item   | Description        | Weight |
|--------|--------------------|--------|
| CT/ALC | Alcohol Gel Holder | 0.3kg  |

### **Airway Tube Holder**

#### **Features**

• Internal dimensions of holder 95 x 50 x 510mm

ol Maid 🙎

- Outer sleeve removable for cleaning
- Fitted with a 15mm hook

#### **Standard Colours**

Paint - Grey White

#### **Dimensions** (w x d x h)

• 140 x 70 x 555mm

### **Debris Bin & Bracket**

#### **Features**

- Plastic debris bin & bracket
- Fitted with a 15mm hook

#### **Standard Colours**

• Paint - Grey White

#### Dimensions (w x d x h)

- Overall Dimensions 210 x 164 x 109mm
- Internal Dimensions 195 x 130 x 95mm





| Item   | Description          | Weight |
|--------|----------------------|--------|
| CT/BIN | Debris Bin & Bracket | 0.4kg  |

sales@bristolmaid.com +44

Item CT/AIR

+44 (0)1258 484455

Airway Tube Holder

Weight

1.4kg

### **Bowl Support & Bracket**

#### **Features**

- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- Bowl retaining ring fitted with removable anti-static spiral buffers to support bowls

#### **Standard Colours**

- Paint Grey White
- Stainless Steel

#### Dimensions (w x d x h)

• 510 x 360 x 325mm

#### Accessories

- 2MECO079 360mm Diameter, Stainless Steel Bowl
- 2PLCO054 360mm Diameter, Polythene Bowl



| Item   | Description            | Weight |
|--------|------------------------|--------|
| CT/BOW | Bowl Support & Bracket | 2.5kg  |

### **Bucket Support & Bracket**

#### **Features**

- Suitable for use with 15-litre graduated stainless steel buckets (order separately)
- Bucket retaining ring fitted with removable anti-static spiral buffers to support the bucket

#### **Standard Colours**

- Paint Grey White
- Stainless Steel

#### Dimensions (w x d x h)

• 510 x 360 x 335mm

#### Accessories

 2MECO081 - 15-litre Graduated Stainless Steel Bucket



| Item   | Description              | Weight |
|--------|--------------------------|--------|
| CT/BUC | Bucket Support & Bracket | 3.9kg  |

### **Document Holder**

### Features

- Suitable for A4 ring binders or foolscap notes (pocket size 274 x 63mm)
- Fitted with a 15mm hook

#### **Standard Colours**

· Paint - Grey White

#### **Dimensions** (w x d x h)

• 275 x 82 x 402mm



| Item     | Description     | Weight |
|----------|-----------------|--------|
| CT/DH/A4 | Document Holder | 1.7kg  |

### **Catheter Racks**

#### **Features**

- Available with 2 or 5 sections
- Clear perspex panel provides visibility of contents, & can be removed to assist cleaning
- Fitted with a 15mm hook

#### Standard Colours

· Paint - Grey White

#### **Dimensions** (w x d x h)

• Sections - 64 x 75 x 400mm (internal)



| Item   | Description  | Dimensions (w x d x h ) | Weight |
|--------|--------------|-------------------------|--------|
| CT/CR2 | Two Section  | 130 x 108 x 405mm       | 1.4kg  |
| CT/CR5 | Five Section | 325 x 108 x 405mm       | 3.0kg  |

### Cylinder Holder, C/CD/D Size

#### **Features**

- Compatible with C/CD/D size cylinders (max diameter 100mm)
- Cylinder held in position using a chain, whilst a rubber pad protects wear points & reduces noise

#### **Standard Colours**

· Paint - Grey White

#### Dimensions (w x d x h)

• 145 x 135 x 400mm



#### Oxygen cylinder not supplied

| Item   | Description                  | Weight |
|--------|------------------------------|--------|
| CT/CYD | Cylinder Holder, C/CD/D Size | 1.4kg  |

### **Infusion Pole & Bracket**

#### **Features**

- Two hook infusion pole & bracket (pole length 1,160mm)
- Attaches to upper carerail
- Maximum capacity on hooks 2kg

#### **Standard Colours**

- Paint Grey White
- Stainless Steel

#### **Dimensions** (w x d x h)

• 238 x 104 x 1160mm



| Item  | Description             | Weight |
|-------|-------------------------|--------|
| CT/IV | Infusion Pole & Bracket | 1.6kg  |

Racks

ංඊ

**Caretray Trolleys** 

## **Phlebotomy Armrest**

#### **Features**

- Fully adjustable easy clean moulded polyurethane armrest
- Attaches to lower carerail

#### Standard Colours

- · Paint Grey White
- Stainless Steel

#### Dimensions (w x d x h)

Armrest Pad - 350 x 120 x 40mm



| Item    | Description        | Weight |
|---------|--------------------|--------|
| CT/PARM | Phlebotomy Armrest | 2.1kg  |

## **Sharps Box Bracket, Universal**

#### **Features**

- · Compatible with the following brackets (customer to supply & fit the plastic brackets - fixings provided)
- Daniels Healthcare
  - POUDS® Wall-Bracket (up to 5L) SHARPSGUARD® mini. 1. 2.5. 3.75. 5 & com plus
  - POUDS® Wall-Bracket (7L & above) SHARPSGUARD® 7.11.5 & 22
- Frontier (Sharpsafe) 2. 3. 4 & 7
- Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS
- Fitted with a 15mm hook

#### Standard Colours

· Paint - Grey White

#### **Dimensions** (w x d x h)

• 75 x 32 x 152mm





Plastic brackets not supplied

Sharps Box Bracket, Universal

Description



Weiaht

0.2kg





| Item   | Description        | Weight |
|--------|--------------------|--------|
| CT/SPS | Suction Pump Shelf | 5.0kg  |

## **Suction Pump Shelf**

#### **Features**

- · Compatible with all popular suction pumps
- Shelf provides an area 415 x 300mm (w x d)
- · 30mm upstand at front of shelf
- Corner bumpers offer protection to both the trollev & fabric of the building
- Floor to base of shelf at height 180mm
- Maximum capacity 15kg

#### Standard Colours

· Paint - Grey White

#### **Dimensions** (w x d x h)

• 460 x 245 x 390mm

Item CT/SBB/15

## Waste Bin & Bracket

#### **Features**

- 8 litre plastic swing bin & metal support frame
- Fitted with a 15mm hook

#### Standard Colours

- · Paint Grey White
- Bin Cream

#### Dimensions (w x d x h)

• 8 Litre Plastic Swing Bin - 220 x 190 x 370mm



| Item  | Description         | Weight |
|-------|---------------------|--------|
| CT/WB | Waste Bin & Bracket | 0.6kg  |

## **Work Flap**

#### **Features**

- Fold down flap (provides work area 300 x 455mm)
- Machined from solid grade laminate allowing for easy cleaning
- Maximum capacity 2kg

#### **Standard Colours**

- Paint Grey White
- Solid Grade Laminate Grey White

#### **Dimensions** (w x d x h)

• 325 x 465 x 70mm



#### **Features**

- Suitable for use with roll type aprons, maximum width 415mm
- Fitted with a 15mm hook

#### **Standard Colours**

• Paint - Grey White

#### Dimensions (w x d x h)

• 434 x 187 x 222mm





| Item     | Description           | Weight |
|----------|-----------------------|--------|
| DAR/1/15 | Apron Dispenser, Roll | 2.2kg  |

Item CT/WF

Work Flap

Weight

3.2kg

## Glove, Mask Box Holder, One Box Glove, Mask Box Holder, Three Box

· Plastic coated wire frame

Standard Colours

Dimensions (w x d x h)

• DG/3/15 - 290 x 141 x 405mm

· Paint - White

• Fitted with a 15mm square hook

· Boxes presented in landscape orientation

• Accommodates glove & mask boxes up to 150mm in width

4. 5. - 2 4

**Features** 

#### **Features**

- Accommodates glove & mask boxes up to 150mm in width
- · Box presented in portrait orientation
- Fitted with a 15mm square hook

#### **Standard Colours**

· Paint - Grey White

#### Dimensions (w x d x h)

• GMBH/15 - 143 x 122 x 177mm



| Item    | Description               | Weight |
|---------|---------------------------|--------|
| GMBH/15 | One Box, 15mm Square Hook | 0.6kg  |

| Item    | Description                 | Weight |
|---------|-----------------------------|--------|
| DG/3/15 | Three Box, 15mm Square Hook | 0.8kg  |

## **Sharps Box Holder**

#### **Features**

- · Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a toggle
- Fitted with a 15mm hook

#### **Standard Colours**

· Paint - Grey White

#### Dimensions (w x d x h)

• 127 x 118 x 230mm



| Item   | Description       | Weight |
|--------|-------------------|--------|
| SBH/15 | Sharps Box Holder | 0.5kg  |



## Laptop Holder, Removeable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with all standard height trolleys
- Fitted to the rear right upright of the trolley
- Laptop supported in a cradle, fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

\*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**

• Paint - Grey White



| Item     | Description          | Weight |
|----------|----------------------|--------|
| CT/LHR/S | Removable, Swing Arm | 5.4kg  |

## Laptop Holder, Non-Removeable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with all standard height trolleys
- Fitted to the rear right upright of the trolley
- Laptop supported in a cradle fitted with a steel mounting bracket, fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

\*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**

· Paint - Grey White



www.bristolmaid.com

| Item      | Description              | Weight |
|-----------|--------------------------|--------|
| CT/LHNR/S | Non-Removable, Swing Arm | 5.7kg  |

Caretray Trolleys



C004



## **Polymer Castors**

#### **Features**

- Offers reduced rolling resistance making the trolley easier to manoeuvre & lowers noise levels
- Smooth surfaces assist cleaning
- Supplied as a set of 4 castors
  - 2 x braking / 2 x non-braking
  - 4 x braking
  - 4 x non-braking
- Factory fitted only

| Item | Description   | Weight |
|------|---|--------|
| C004 | 100mm, Grey White Polymer, 2 x Braking / 2 x Non-Braking (set of 4) | 2.0kg  |
| C005 | 100mm, Grey White Polymer, 4 x Braking (set of 4)                   | 2.0kg  |
| C006 | 100mm, Grey White Polymer, 4 x Non-Braking (set of 4)               | 2.0kg  |

## **Ring Buffers**

#### **Features**

- Set of four ring buffers
- Offers protection to both the trolley & fabric of the building
- Increases the length & width of trolley by approx. 80mm

www.bristolmaid.com

Factory fitted



| Item  | Description                | Weight |
|-------|----------------------------|--------|
| ZP010 | Ring Buffers (set of four) | 0.3kg  |

## Trays

#### **Features**

- 100/150mm deep translucent high impact plastic trays with integral handle (max load per tray 5kg)
- Each tray is supplied with plastic clip on label cover & four removable dividers (2 x long & 2 x short)
- Trays are supplied without coloured label strips (see CT/LAB)

#### Internal Dimensions (w x d x h)

- CT/TRAY/S 355 x 405 x 100mm
- CT/TRAY/D 355 x 405 x 150mm





| Item      | Description        | Weight |
|-----------|--------------------|--------|
| CT/TRAY/S | 100mm Shallow Tray | 1.5kg  |
| CT/TRAY/D | 150mm Deep Tray    | 1.8kg  |

## **Dividers**

#### **Features**

- 100/150mm deep translucent dividers
- Supplied as a set of 2 x long & 2 x short



www.bristolmaid.com



| Item     | Description                              | Weight |
|----------|--|--------|
| CT/DIV/S | Set of 4 dividers for 100mm Shallow Tray | 0.6kg  |
| CT/DIV/D | Set of 4 dividers for 150mm Deep Tray    | 0.4kg  |

## Labels, Coloured

#### **Features**


- Coloured labels Red / Yellow / Green / Blue / White
- Labels position beneath the clip-on label cover on the plastic trays
- Supplied as a set of 6 sheets



## Labels, Sepsis Six

#### Features

- · Labels to support the delivery of 'The Sepsis Six'
- Printed Oxygen / Blood Cultures / Intravenous Antibiotics / Intravenous Fluid / Haemoglobin & Lactates / Urine Measurement
- Labels position beneath the clip-on label cover on the plastic trays
- Supplied as a set of 6 labels



| ltem   | Description       |
|--------|-------------------|
| CT/LAB | Labels - Coloured |

| Item       | Description         |
|------------|---------------------|
| CT/LAB/SEP | Labels - Sepsis Six |

sales@bristolmaid.com +44 (0)1258 484455

## Flat Shelf, Removable

#### Features

- Solid grade laminate (max load per shelf 5kg)
- Supported on the runner guides (non-pull out)

#### **Standard Colours**

• Solid Grade Laminate - Grey White

#### Dimensions (w x d x h)

• 413 x 457 x 4mm



CT/FS Flat Shelf - Removable 1.1kg

## Full Width Shelf, Removable

#### **Features**

- Solid grade laminate (max load per shelf 5kg)
- Supported on the trolley frame

#### **Standard Colours**

• Solid Grade Laminate - Grey White

#### Dimensions (w x d x h)

• Overall Dimensions - 480 x 480 x 4mm





#### www.bristolmaid.com



# Chart Trolleys Workstations Room Dividers & Screens

- Fixed Height Trolleys 82
- Fixed Height Workstations 84
- Variable Height Workstations 86
  - Chart Boards 88
  - Document Holders 89
    - Room Dividers 90
    - Curtain Screens 92
      - Solid Screens 93





## **Caretray Trolley - Chart**

#### Features

- Standard level, single & double column easy clean fixed height chart trolleys
- Angle adjustable chartboard fitted with a spring clip (adjustable 0 50 degrees)
- Usable chart area...
  - Single Column 500 x 530mm (w x h)
  - Double Column 800 x 530mm (w x h)
- · Lower lip on chartboard to support pens etc.
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side (2 x L/H, 1 x R/H) support a wide range of accessories
- Drop down work flap fitted to R/H side
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- 100mm pressed steel castors, 2 off braking

### Dimensions (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 620 x 1040 1450mm
  - Double Column 1010 x 550 620 x 1040 1450mm
- Usable Chart Area
  - Single Column 500 x 530mm
  - Double Column 800 x 530mm
- Work Flap 420 x 290mm (w x d)
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

- Paint White
- Chartboard / Work Flap White

| Item             | Description                | Weight |
|------------------|----------------------------|--------|
| CTCHARM1/10/3S2D | 3 x 100mm, 2 x 150mm Trays | 33.2kg |
| CTCHARM1/10/4S1D | 4 x 100mm, 1 x 150mm Trays | 32.9kg |

#### **Options / Accessories**

- Carerail Mounted page 68
- Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77



Angle adjustable chartboard (adjustable 0 - 50 degrees)

| Item             | Description                | Weight |
|------------------|----------------------------|--------|
| CTCHARM2/10/6S4D | 6 x 100mm, 4 x 150mm Trays | 53.2kg |
| CTCHARM2/10/8S2D | 8 x 100mm, 2 x 150mm Trays | 51.2kg |





## **Chart Workstations - Fixed Height**

#### **Features**

- Fixed height chart workstation
- Angle adjustable chartboard fitted with a spring clip (0 50 degrees)
- Usable chart area 800 x 530mm (w x h)
- Lower lip on chartboard to support pens etc.
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- · Models available with
  - Chartboard only
  - Storage pocket
- Carerails fitted to each side of the trolley (2 x L/H, 1 x R/H) support a wide range of accessories
- Drop down work flap fitted to R/H side
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions 1010 x 550 620 x 1040 1450mm
- Usable Chart Area 800 x 530mm
- Storage Pocket 825 x 450 x 150mm (internal)
- Work Flap 420 x 290mm (w x d)

#### **Standard Colours**

- Paint White
- · Chartboard / Work Flap White

| Item      | Description                                      | Weight |
|-----------|--|--------|
| CWF/80    | Chart Workstation - Fixed Height                 | 23.0kg |
| CWF/80/SP | Chart Workstation - Fixed Height, Storage Pocket | 33.0kg |

#### **Options / Accessories**

- Carerail Mounted page 68
- Castors & Ring Buffers page 76

Angle adjustable chartboard (0 - 50 degrees)





## **Chart Workstations - Variable Height**

#### **Features**

- Height adjustable using a gas damper (795 1075mm)
- Angle adjustable chartboard fitted with spring clips (adjustable 0 - 90 degrees)
- Chartboard available in two sizes (w x h)
  - Small 790 x 700mm
  - Large 790 x 800mm
- Lower lip on chartboard to support pens etc
- · Slide out shelf at front
- Document storage pocket on side to accept patients notes, ring binders, etc
- 80mm swivel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions 800 x 700 / 800 x 1095mm
- Document Storage Pocket 75 x 255 x 360mm (internal)

#### **Standard Colours**

- Paint White
- · Chartboard White

#### **Options / Accessories**

• ZP010 - 100mm Ring Buffers

| Item     | Description  | Weight |
|----------|--|--------|
| CB/100/A | Chart Workstation - Variable Height, Small, 700mm High Board | 40.0kg |
| CB/100/B | Chart Workstation - Variable Height, Large, 800mm High Board | 41.0kg |

## Other products that you may be interested in...

See Pages 374 & 375

### **Static Safe Tech Chairs**



+44 (0)1258 484455

www.bristolmaid.com

sales@bristolmaid.com

87



## Chart Boards - A4 & A3

#### Features

- A4 one spring clip (holds one chart 220 x 280mm)
  A3 two spring clips (holds one chart 450 x 280mm or two charts 220 x 280mm)
- Available with a 30 or 50mm mild steel square hook

### **Standard Colours**

White

| Item    | Description   | Overall Dimensions (w x d x h) | Weight |
|---------|---------------|--------------------------------|--------|
| A4 Char | t Board       |                                |        |
| CB005   | A4, 30mm Hook | 230 x 80 x 370mm               | 0.5kg  |
| CB010   | A4, 50mm Hook | 230 x 100 x 370mm              | 0.5kg  |
| A3 Char | t Board       |                                |        |
| CB015   | A3, 30mm Hook | 460 x 80 x 370mm               | 0.8kg  |
| CB020   | A3, 50mm Hook | 460 x 100 x 370mm              | 0.8kg  |

sales@bristolmaid.com



| Item      | Description      | Overall Dimensions (w x d x h) | Weight |
|-----------|------------------|--------------------------------|--------|
| A4 Docume | ent Holder       |                                |        |
| DH/A4/WF  | A4, Wall Fitting | 280 x 70 x 400mm               | 1.4kg  |
| DH/A4/50  | A4, 50mm Hook    | 280 x 125 x 405mm              | 1.4kg  |
| A3 Docume | ent Holder       |                                |        |
| DH/A3/WF  | A3, Wall Fitting | 405 x 105 x 400mm              | 2.9kg  |
| DH/A3/50  | A3, 50mm Hook    | 405 x 160 x 405mm              | 3.0kg  |

## **Document Holders - A4 & A3**

#### Features

- A4 size suitable for A4 ring binders or foolscap notes (pocket size 260 x 65mm)
- A3 size suitable for A3 notes & x-rays (pocket size 398 x 100mm)
- Open front design allows for easy patient identification
- Available for wall fitting or with a 50mm square hook

#### **Standard Colours**

• Paint - Grey White





## **Room Dividers**

#### Features

- · Fixed height mobile room divider/workstation
- Ideal for use in critical care areas A&E. ICU. ICT. HDU
- Modular design allows the unit to be configured to suit your specific requirements
- Frame, worktop & upstand/panel available as...
  - Full width upstand
  - Full width upstand & right hand mounted panel
  - Full width upstand & left hand mounted panel
- Full width panel
- · Four configurations of lower units...
  - Three drawer unit, upper drawer fitted with an electronic push button lock (2 x 170, 1 x 285mm drawers)
  - Five drawer unit, upper drawer fitted with an electronic push button lock (5 x 125mm drawers)
  - · Cupboard with double doors & removable shelf
  - Cover panel (creates a space within the unit)
- Easy clean solid grade laminate work surface & up-stand/panel
- · Upper panel fitted with
  - Triple glove box holder
  - Alcohol gel holder (compatible with 500ml bottles, 70 x 70 x 225mm)
  - Apron roll dispenser
- 100mm ring buffers to protect the fabric of the building
- 100mm swivel castors, 2 off braking

#### Dimensions (w x d x h)

- RD1 1768 x 698 x 1000mm
- RD2R/RD2L 1768 x 698 x 1588mm
- RD3 1768 x 698 x 1350mm


#### **Standard Colours**

- · Paint White
- · Worksurface / Up-Stand / Panel White

#### **Options / Accessories**

- CC Colour coded drawers available on request
- DTSF/RD/450/2 Dressing Trolley
  - · Seamlessly welded 304 grade stainless steel frame
  - · Fixed shelves with down turned sides
  - Shelf Heights 375 & 840mm
  - Max load per shelf 40kg

## Configure Your Workstation in Four Easy Steps...







MS/C4

+44 (0)1258 484455

sales@bristolmaid.com

## Screen - Curtain, Four Section, Folding

#### Features

- PVC fire retardant curtain (machine washable)
- Available as four section
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors

#### Dimensions (w x d x h)

- Overall Screen Height 1708mm
- Lowest Screen Height 288mm
- Folded 650 x 420 x 1790mm

#### **Standard Colours**

- Paint White
- Curtain Grey

#### **Spares**

• WS/SPARES/001 - Set of 4 replacement curtains

www.bristolmaid.com

Overall Dimensions (w x d x h)

Weight

14.6kg

## Mobile Screen - Solid, Three & Four Section, Folding

#### Features

- Fire retardant folding compact laminate solid screens
- Available as either three or four section
- Smooth surface assists cleaning
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors

### Dimensions (w x d x h)

- Overall Screen Height 1680mm
- Lowest Screen Height 180mm
- Folded 600 x 435 x 1680mm

#### **Standard Colours**

- · Paint Spearmint
- · Screen White



| Item  | Description                                   | Overall Dimensions (w x d x h) | Weight |
|-------|---|--------------------------------|--------|
| MS/S3 | Mobile Screen - Solid, Three Section, Folding | 1802 x 435 x 1680mm            | 20.8kg |
| MS/S4 | Mobile Screen - Solid, Four Section, Folding  | 2403 x 435 x 1680mm            | 27.8kg |





# **Concealment Trolleys**

- Variable Height Trolleys 96
  - Body Trays & Covers 97



## **Variable Height**

#### **Features**

- · Height adjustable
  - Fixed Tray 580 1280mm
  - Loose Tray 610 1310mm
- Scissor action provides vertical lift, ideal for use in confined areas where space is limited
- Battery powered actuator
  - Rechargeable battery is charged through the handset using a UK 3 pin plug
  - Audible warning when battery power is at 50%
  - Controlled using a handset
- Liquid retaining tray
  - Fixed Tray 1900 x 540 x 25mm
  - Loose Tray 2100 x 560 x 70mm
- Drop down side-rails, distance between rails 660mm
- Bumper bar / push handle at either end
- · Heavy duty reinforced PVC cover (locates over side rails)
- 200mm heavy duty castors, 2 x braking, 2 x directional locking
- Maximum capacity 190kg

#### **Standard Colours**

• Paint - White



**Concealment Trolleys** 

2150 x 770 x 580 - 1280mm

2270 x 770 x 610 - 1310mm

www.bristolmaid.com

Variable Height, Fixed Body Tray

Variable Height, Loose Body Tray

VH/950

VH/950/LBT

## **Body Tray**

#### Features

- · Heavy duty body tray, mild steel frame & stainless steel tray
- Integral runners
- Fits most standard mortuary racks (check prior to ordering)
- Maximum capacity 190kg

#### Dimensions (w x d x h)

• Overall Dimensions - 2100 x 550 x 70mm

#### Standard Colours

· Paint - White

Item



#### Features

· Heavy duty reinforced PVC cover (locates over side rails)





| Item     | Description                       | Weight |
|----------|-----------------------------------|--------|
| 2FACO024 | Cover - PVC Heavy Duty Reinforced | 5.0kg  |

www.bristolmaid.com

Shown with VH/950



# Cylinder Storage & Transportation

- Cylinder Trolleys F/G Size 100
  - Cylinder Trolleys J Size 101
- Bulk Transportation Trolleys 102
- Cylinder Trolleys E, CD/D, Size 102
- Wall Brackets & Cylinder Holders 103

## Cylinder Trolley F/G Size

#### **Features**

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinder
- 160mm anti-static wheels

#### Dimensions (w x d x h)

• 510 x 405 x 1145mm

#### **Standard Colours**

Paint - Grey



| Item    | Description | Weight |
|---------|-------------|--------|
| CYT/FG1 | F/G Size    | 10.0kg |

## Cylinder Trolley with Stabilisers F/G Size

#### **Features**

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinder
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 160mm anti-static wheels, 100mm swivel castors

#### Dimensions (w x d x h)

• 520 x 700 x 1145mm

#### **Standard Colours**

Paint - Grey



| Item       | Description | Weight |
|------------|-------------|--------|
| CYT/FG1/ST | F/G Size    | 14.0kg |

www.bristolmaid.com

## **Cylinder Trolley J** Size

#### **Features**

- Compatible with J (240 cu ft), AK, L & HL cylinders, max diameter 245mm, height 1550mm
- One piece construction
- Drop over chains secure the cylinder
- 160mm anti-static wheels

#### Dimensions (w x d x h)

• 510 x 475 x 1290mm

#### Standard Colours

Paint - Grey



| Item   | Description |
|--------|-------------|
| CYT/J1 | J Size      |

## **J** Size

- **Features**
- Compatible with J (240 cu ft), AK, L & HL cylinders, max diameter, max diameter 245mm, height 1550mm

**Cylinder Trolley with Stabilisers** 

- One piece construction
- Drop over chains secure the cylinder
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 160mm anti-static wheels, 100mm swivel castors

#### Dimensions (w x d x h)

• 510 x 740 x 1290mm

#### **Standard Colours**

· Paint - Grey



| Item      | Description | Weight |
|-----------|-------------|--------|
| CYT/J1/ST | J Size      | 14.8kg |

+44 (0)1258 484455

Weight 10.5kg

## **Bulk Transportation Trolley** D & E Cylinders, 12 x Capacity

#### **Features**

- Compatible with D, AD, CD, DD, PD, RD, ZD & E (24 cu ft) cylinders (12 x capacity) max diameter 105mm, height 870mm
- One piece construction
- Cylinders locate into individual sections, supported on a non slip rubber base
- 125mm swivel castors, 2 off braking

#### Dimensions (w x d x h)

• 460 x 590 x 1005mm

#### Standard Colours

· Paint - Grey



| ltem     | Description                    | Weight |
|----------|--------------------------------|--------|
| CYT/DE12 | D & E Cylinders, 12 x Capacity | 23kg   |

## **Cylinder Trolley** F Size

#### **Features**

- · Compatible with E (24 cu ft) cylinders, max diameter 115mm, height 900mm
- One piece construction
- · Recessed base plate and drop over chain to secure the cylinder
- 40mm hook allows trolley to hang at bed end, a useful feature when transporting patients
- 100mm anti-static wheels

#### Dimensions (w x d x h)

• 320 x 235 x 980mm

#### **Standard Colours**

· Paint - Grey



| Item   | Description | Weight |
|--------|-------------|--------|
| CYT/E1 | E Size      | 4.3kg  |

## **Cylinder Trolley CD/D Size**

#### **Features**

- Compatible with CD & D cylinders, max diameter 115mm, height 540mm
- One piece construction
- Recessed base plate and drop over chain to secure the cvlinder
- 40mm hook allows trolley to hang at bed end, a useful feature when transporting patients
- 100mm anti-static wheels

#### Dimensions (w x d x h)

• 320 x 235 x 980mm

#### Standard Colours

· Paint - Grey



| Item    | Description | Weight |
|---------|-------------|--------|
| CYT/CD1 | CD/D Size   | 4.3kg  |

sales@bristolmaid.com


+44 (0)1258 484455

## Wall Brackets Cylinder

#### Features

- Max cylinder diameters
- D/E 112mm, F 148mm, G 185mm, J 235mm
- One piece construction
- Drop over chains to secure the cylinder
- · Holes for wall mounting (fixings not provided)
- Should be positioned on the wall at approx half cylinder height

#### **Standard Colours**



| Item    | Description       | Dimensions (w x d x h) | Weight |
|---------|-------------------|------------------------|--------|
| CYB/DE2 | 2 x D/E Cylinders | 359 x 107 x 88mm       | 1.1kg  |
| CYB/F2  | 2 x F Cylinders   | 435 x 126 x 88mm       | 1.4kg  |
| CYB/G2  | 2 x G Cylinders   | 511 x 145 x 88mm       | 1.7kg  |
| CYB/J2  | 2 x J Cylinders   | 613 x 170 x 88mm       | 2.1kg  |

## Cylinder Holder D & E Size, Vertical

#### Features

- Compatible with...
  - D Size D, AD, CD, DD, PD, RD & ZD cylinders (max 112mm diameter)
  - E Size E (max 112mm diameter)
- One piece construction
- Carry handle and extended hook
- Suitable for beds with 25mm wide side/end members
- Holder can be hung and/or bolted to 25mm section bed-frame

#### Dimensions (w x d x h)

- CYH/D1 200 x 160 x 290mm
- CYH/E1 200 x 160 x 595mm

CYH/D1

+44 (0)1258 484455

Description

D Size, Vertical

E Size, Vertical

#### **Standard Colours**

• Paint - Grey

Item

CYH/D1

CYH/E1

## Cylinder Holder Wall Plate

#### **Features**

- Painted stainless steel wall plate to support cylinder holders
- Compatible with CYH/D1, CYH/E1, CYHP/D1 & CYHP/E1

#### Dimensions (w x d x h)

• 240 x 47 x 254mm

#### **Standard Colours**

• Paint - Grey





CYH/F1

Weight

1.4kg

2.5kg



## **Dispensers, Holders & Racks**

- Apron Dispensers, Pocket Mask Holders 106
  - Glove / Mask Box Holders 107
    - Alcohol Gel Holders 108
- Sharps Box Holders, Urine Bottle Holder, Catheter Bag Holders 109
  - Catheter Racks 110
  - Drainage Racks 111

## **Apron Dispensers - Roll**

#### **Features**

- Apron roll holders, maximum width 415mm
- Holes for wall mounting (fixings not provided)

#### Dimensions (w x d x h)

- One Roll 420 x 168 x 268mm
- Three Roll 420 x 168 x 615mm

#### Standard Colours

· Paint - White



| Item  | Description | Weight |
|-------|-------------|--------|
| DAR/1 | One Roll    | 2.0kg  |
| DAR/3 | Three Rolls | 8.2kg  |

## **Apron Dispenser - Flat Pack**

#### Features

- Holder to accept flat packed aprons, maximum width 400mm
- Packs of aprons fold over the centrally positioned divider
- Holes for wall mounting (fixings not provided)

#### Dimensions (w x d x h)

• 410 x 105 x 290mm

#### Standard Colours

· Paint - White



| Item  | Description | Weight |
|-------|-------------|--------|
| DAF/1 | Flat Pack   | 2.0kg  |

## Pocket Mask Holder

#### **Features**

- Holder to accept Laerdal & Ambu Pocket Masks (hard cases)
- Allows the case to be guickly removed from the holder
- Supplied with self adhesive tape/holes for wall mounting (fixings not provided)

#### Dimensions (w x d x h)

• 110 x 40 x 110mm

#### Standard Colours

· Paint - White



Pocket mask not supplied



| Item   | Description        | Weight |
|--------|--------------------|--------|
| PMH005 | Pocket Mask Holder | 0.1kg  |

& Racks

**Dispensers, Holders** 

sales@bristolmaid.com

# Glove / Mask Box Holders Three Box

#### **Features**

- Plastic coated wire frame
- Accommodates glove & mask boxes up to 150mm in width
- Boxes presented in landscape orientation
- Available for wall mounting or fitted with a 15mm square hook
   (fixings not provided)

#### Dimensions (w x d x h)

- DG/3 290 x 115 x 405mm
- DG/3/15 290 x 141 x 405mm

#### **Standard Colours**

· Paint - White



| Item    | Description                 | Weight |
|---------|-----------------------------|--------|
| DG/3    | Three Box, Wall Mounting    | 0.7kg  |
| DG/3/15 | Three Box, 15mm Square Hook | 0.8kg  |

# Glove / Mask Box Holders One Box

#### Features

- Accommodates glove & mask boxes up to 150mm in width
- · Box presented in portrait orientation
- Available for wall mounting or fitted with a 15mm square hook (fixings not provided)

#### Dimensions (w x d x h)

- GMBH/WM 143 x 100 x 152mm
- GMBH/15 143 x 122 x 177mm

#### **Standard Colours**

· Paint - Grey White



#### **Features**

- Packs of sterile gloves separated by two
   adjustable dividers
- Full width card holder for identification
- Clear front panel assists with stock level inspection
- Hinged lid assists replenishment, whilst top loading ensures stock rotation
- Holes for wall mounting (fixings not provided)

#### Dimensions (w x d x h)

• 565 x 260 x 315mm

#### **Standard Colours**

• Paint - White



| Item    | Description              | Weight |
|---------|--------------------------|--------|
| SS/48/3 | Three Column, Adjustable | 7.0kg  |

Item

GMBH/WM

GMBH/15

Weight

0.5ka

Description

One Box. Wall Mounted

One Box, 15mm Square Hook

## **Alcohol Gel Holders - Wire Frame**

#### **Features**

- Compatible with popular 500ml bottles, 70mm x 70mm x 225mm (not supplied)
- Stainless steel construction
- · Plastic drip tray
- Available with holes for wall mounting (fixings not provided) or with a 30mm square hook

#### **Standard Colours**

Stainless Steel



# Item Description Overall Dimensions (w x d x h) Weight GH005 Lever Operated, Wall Mounted, Angled 77 x 300 x 278mm 0.4kg GH015 30mm Hook 77 x 117 x 180mm 0.3kg

sales@bristolmaid.com

# Syringe Dispenser - Five Column, Adjustable

#### **Features**

- Syringes separated by four adjustable dividers
- Full width card holder for identification
- Clear front panel assists with stock level inspection
- Hinged lid assists replenishment, whilst top loading ensures stock rotation
- Holes for wall mounting (fixings not provided)

#### Dimensions (w x d x h)

• 565 x 317 x 320mm

#### **Standard Colours**

• Paint - White



| Item  | Description             | Weight |
|-------|-------------------------|--------|
| SS/36 | Five Column, Adjustable | 9.2kg  |

# **Sharps Box Holder One Box, Wall Mounted**

#### Features

- · Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a togale
- Holes for wall mounting (fixings not provided

#### Dimensions (w x d x h)

• 127 x 100 x 190mm

#### Standard Colours

· Paint - White



# **Urine Bottle Holder** Hanging

#### **Features**

- · Compatible with bottles up to 115mm Diameter
- Plastic coated wire frame
- 30mm square hook

#### Dimensions (w x d x h)

• 150 x 182 x 282mm

#### Standard Colours

· Plastic - White



| Item | Description  | Weight |
|------|--------------|--------|
| UB/9 | Urine Bottle | 0.3kg  |

# **Catheter Bag Holders Free Standing/Hanging**

#### **Features**

- · Plastic coated wire frame
- · Accepts most widely used bags
- Available freestanding or hanging (30mm square hook)
- · Upper clip to hold the inlet tube
- Lower clip angles the drainage tube upwards to prevent drips should the tap leak

#### Standard Colours

· Plastic - White

www.bristolmaid.com



| Item  | Description  | Overall Dimensions (w x d x h) | Weight |
|-------|--------------|--------------------------------|--------|
| UB005 | Freestanding | 250 x 250 x 550mm              | 0.6kg  |
| UB010 | Hanging      | 250 x 77 x 410mm               | 0.4kg  |

# **Catheter Racks - Two & Five Section**

#### Features

• Internal dimensions of each section 65 x 75 x 405mm

RC/2

- Removable clear perspex panel to aid cleaning
- Available with holes for wall mounting (fixings not provided) or with a 15mm square hook

#### **Standard Colours**

• Paint - White



| Item Description |              | Overall Dimensions (w x d x h) | Weight |
|------------------|--------------|--------------------------------|--------|
| Two Section      |              |                                |        |
| RC/2             | Wall Fitting | 135 x 85 x 410mm               | 1.3kg  |
| RC/2/15          | 15mm Hook    | 135 x 115 x 410mm              | 1.4kg  |
| Five Section     |              |                                |        |
| RC/5             | Wall Fitting | 320 x 85 x 410mm               | 3.3kg  |
| RC/5/15          | 15mm Hook    | 320 x 115 x 410mm              | 3.4kg  |

sales@bristolmaid.com

# **Catheter Racks - One Section**

#### Features

- Outer sleeve removable for cleaning
- 20mm square hook

#### Dimensions (w x d x h)

• 185 x 55 x 555mm

#### Standard Colours

· Paint - White



| Item     | Description | Weight |
|----------|-------------|--------|
| SS/40/MS | One Section | 2.0kg  |


# **Drainage Racks - Urinal**

#### Features

- Holds five urine bottles
- Stainless steel construction
- Supplied with a stainless steel drip tray
- Holes for wall mounting (fixings not provided)

#### Dimensions (w x d x h)

• 535 x 205 x 415mm



# Drainage Racks - Bedpan

#### Features

- Holds five bedpans
- Stainless steel construction
- Supplied with a stainless steel drip tray
- Holes for wall mounting (fixings not provided)

#### Dimensions (w x d x h)

• 530 x 295 x 375mm



| Item     | Description             | Weight |
|----------|-------------------------|--------|
| RD/U5/SS | Urinal, Stainless Steel | 1.7kg  |

| Description             | Weight |
|-------------------------|--------|
| Bedpan, Stainless Steel | 3.3kg  |

Item RD/B5/SS



# **Dressing / Instrument Trolleys**

- Removable Shelves Stainless Steel 114
  - Removable Shelves Mild Steel 115
- Fixed Shelves, Sides Up Stainless Steel 116
- Fixed Shelves, Sides Down Stainless Steel 117
  - Drawer Units 119
  - Options & Accessories 120



# Removable Shelves, Stainless Steel, 450mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Flat or reversible 304 grade stainless steel shelves, cushioned on anti-static spiral buffers (removable for cleaning)
- Reversible shelves allow the end user to choose between a flat shelf or one with three upturned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

#### **Options / Accessories**

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C013 75mm, Pressed Steel Castors (4 x braking) refer to page 121
- Drawer Units page 119 Factory Fitted Only
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only

| Flat Shelves | Reversible Shelves | Size    | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|--------------|--------------------|---------|--------------------|--------------------------------|--------|
| DTSR/450     | DTSR/450/SR        | Small   | 450 x 450mm        | 470 x 470 x 890mm              | 11.5kg |
| DTSR/600     | DTSR/600/SR        | Medium  | 600 x 450mm        | 620 x 470 x 890mm              | 16.0kg |
| DTSR/750     | DTSR/750/SR        | Large   | 750 x 450mm        | 770 x 470 x 890mm              | 18.0kg |
| DTSR/900     | DTSR/900/SR        | X-Large | 900 x 450mm        | 920 x 470 x 890mm              | 20.0kg |

sales@bristolmaid.com

```
www.bristolmaid.com
```



# Removable Shelves, Mild Steel, 450mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- · Seamlessly welded mild steel frame
- Flat or reversible 304 grade stainless steel shelves, cushioned on anti-static spiral buffers (removable for cleaning)
- Reversible shelves allow the end user to choose between a flat shelf or a shelf with three upturned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

#### **Standard Colours**

Paint - White

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C013 75mm, Pressed Steel Castors (4 x braking) refer to page 121
- Drawer Units page 119 Factory Fitted Only
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only

| Flat Shelves | Reversible Shelves | Size   | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|--------------|--------------------|--------|--------------------|--------------------------------|--------|
| DTMR/450     | DTMR/450/SR        | Small  | 450 x 450mm        | 470 x 470 x 890mm              | 11.5kg |
| DTMR/600     | DTMR/600/SR        | Medium | 600 x 450mm        | 620 x 470 x 890mm              | 16.0kg |



### Fixed Shelves, Sides Up, Stainless Steel, 450mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with three sides up & one down
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

#### **Options / Accessories**

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C013 75mm, Pressed Steel Castors (4 x braking) refer to page 121
- Drawer Units page 119 Factory Fitted Only
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only

| Item     | Size    | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|----------|---------|--------------------|--------------------------------|--------|
| DTSF/450 | Small   | 450 x 450mm        | 470 x 470 x 890mm              | 11.5kg |
| DTSF/600 | Medium  | 600 x 450mm        | 620 x 470 x 890mm              | 16.0kg |
| DTSF/750 | Large   | 750 x 450mm        | 770 x 470 x 890mm              | 18.0kg |
| DTSF/900 | X-Large | 900 x 450mm        | 920 x 470 x 890mm              | 20.0kg |

sales@bristolmaid.com

DTSF/450



# Fixed Shelves, Sides Down, Stainless Steel, 450mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trollevs for instruments. dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- · Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

#### **Options / Accessories**

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C013 75mm, Pressed Steel Castors (4 x braking) refer to page 121
- Drawer Units page 119 Factory Fitted Only
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only

| sales@bristolmaid.com |
|-----------------------|
|                       |

18.0kg

21.0kg

770 x 470 x 870mm

920 x 470 x 870mm

750 x 450mm

DTSF/750/2/SD

Large DTSF/900/2/SD X-Large 900 x 450mm



# Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel anti-static non-marking castors
- \* Drawers not suitable for this product

- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C015 100mm, Pressed Steel Castors (4 x braking) Factory Fitted
   Only refer to page 121
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only refer to page 121

| Item          | Size   | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|---------------|--------|--------------------|--------------------------------|--------|
| DTSF/66/2/SD  | Small  | 600 x 600mm        | 620 x 620 x 870mm              | 17.5kg |
| DTSF/96/2/SD  | Medium | 900 x 600mm        | 920 x 620 x 870mm              | 24.0kg |
| DTSF/126/2/SD | Large  | 1200 x 600mm       | 1220 x 620 x 870mm             | 29.6kg |





2HH-7011100

#### DTSR/450/DUS/R

| Item                 | Description                                   | Trolley Compatibility | Weight |
|----------------------|---|-----------------------|--------|
| Removal              | ble Shelf Trolleys                            |                       |        |
| DUS/R                | Drawer Unit - Small, Removable Shelf Trolleys | 450mm & 600mm         | 5.1kg  |
| DUL/R                | Drawer Unit - Large, Removable Shelf Trolleys | 750mm & 900mm         | 6.9kg  |
| Fixed Shelf Trolleys |   |                       |        |
| DUS/F                | Drawer Unit - Small, Fixed Shelf Trolleys     | 450mm & 600mm         | 4.7kg  |
| DUL/F                | Drawer Unit - Large, Fixed Shelf Trolleys     | 450mm & 600mm         | 6.5kg  |

# **Drawer Units**

#### Features

- Compatible with removable & fixed shelf ranges, with exception of 600mm deep trolleys
- · Available in two sizes
  - Small 300 x 400mm
  - Large 600 x 400mm
- Each drawer unit is supplied with a 100mm deep HTM71 modular basket
- Stainless steel basket support & drawer carriage
- Precision telescopic runners
- Solid grade laminate drawer front
- Maximum load 5kg
- Trolleys accept only one tier of drawers
- Drawers are not lockable
- Factory fitted & must be ordered at the same time as the trolley

#### Internal Dimensions (w x d x h)

- Small 235 x 335mm (internal)
- Large 535 x 335mm (internal)

#### **Options / Accessories**

- 2HH-7014100 Basket, 300 x 400 x 100mm
- 2HH-7012100 Basket, 600 x 400 x 100mm
- 2HH-7811206N Divider, 300 x 100mm
- 2HH-7811205N Divider, 400 x 100mm
- 2HH-7811204N Divider, 600 x 100mm
- Additional Trays, Baskets, Dividers & Accessories page 206

#### sales@bristolmaid.com

# **Reversible Shelves**

#### Features

- Allows the end user to choose between a flat shelf or one with three upturned sides
- Removable for cleaning
- Compatible with DTMR & DTSR ranges only



Description

**Reversible Shelves** 

# **Clamping Infusion Stands - Trolley**

#### **Features**

- Clamp secures to 25mm vertical round tube (designed to fit our range of dressing trolleys)
- When fitting to other trolleys, we strongly recommend that a risk assessment is carried out to assess stability
- When fitted protrudes 45mm
- Supplied with a two pigtail hook infusion pole
- Maximum capacity on hooks 2kg

#### Dimensions (w x d x h)

• 55 x 220 x 1200mm





| Item | Description                        | Weight |
|------|------------------------------------|--------|
| DP/C | Clamping Infusion Stands - Trolley | 1.8kg  |

Item

SR

Weight

5-10kg

# Castors

#### Features

- Polymer & pressed steel castors
- Anti-static & non-marking
- Supplied as a set of 4 castors
- 80mm polymer castors increase the height of the trolley by approx. 5mm
- 100mm polymer castors increase the height of the trolley by approx. 30mm
- 100mm steel castors increase the height of the trolley by approx. 30mm
- · Factory fitted



#### Item Description

- C001 Castor 80mm, Grey Polymer, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C002 Castor 100mm, Grey Polymer, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C003 Castor 100mm, Pressed Steel, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C013 Castor 75mm, Pressed Steel (4 x braking) Factory Fitted Only
- C015 Castor 100mm, Pressed Steel (4 x braking) Factory Fitted Only

| <b>Ring B</b> | uffers |
|---------------|--------|
|---------------|--------|

#### Features

- · Set of four ring buffers
- Offers protection to both the trolley & fabric of the building
- Increases the length & width of trolley by approx. 80mm
- Factory fitted





| Item  | Description                                      | Weight |
|-------|--|--------|
| ZP005 | Ring Buffers (set of four) - Factory Fitted Only | 0.3kg  |



# **Drug & Medicine Dispensing**

- Blister Pack (MDS) Trolleys 136
- Nomad<sup>®</sup> Cassette Trolleys 137
- Biodose<sup>®</sup> & MultiMeds Trolleys 138
  - Medinoxx Trolleys 139
    - Security Clamp 142
  - Distribution Trolleys 143

- Dispensing Trolleys 124
  - Laptop Holders 129
- CLIQ<sup>®</sup> & RFID Solutions 131
- Patient Administration Trolleys 132
  - Drawer Trolleys 134
    - Tray Trolleys 135





CTDISPM1/10/3S2D/HSB

| Description   | Weight   |
|---|--|
|   |  |
| 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 35.5kg   |
| 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 35.2kg   |
| 6 x 100mm Trays, High Security Bolt Lock            | 36.3kg   |
|   | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock<br>4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock |

# **Caretray Trolley - Dispensing**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines
- Suitable for use with Electronic Prescribing & Medicines Administration (EPMA) & Electronic Patient Records (EPR)
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Waste Bin & Bracket, 8 Litre
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking
- When not in use, trolleys storing medicines should be secured to the wall using a Security Clamp (see locking & security accessories see page 142)

#### Dimensions (w x d x h)

#### **Standard Colours**

• Paint - Spearmint

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Carerail Mounted page 68

- Laptop Holders page 74
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77





| CTDISPM1/ | /10/3S2D/EI | B |
|-----------|-------------|---|
|           |             |   |

| Item                     | Description   | Weight |
|--------------------------|---|--------|
| Electronic Push Button L | ock   |        |
| CTDISPM1/10/3S2D/EPB     | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 35.5kg |
| CTDISPM1/10/4S1D/EPB     | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 35.2kg |
| CTDISPM1/10/6S/EPB       | 6 x 100mm Trays, Electronic Push Button Lock            | 36.3kg |

sales@bristolmaid.com



| Item                 | Description                         | Weight |
|----------------------|-------------------------------------|--------|
| <b>High Security</b> | Bolt Lock                           |        |
| MT100                | Small, High Security Bolt Lock      | 26.0kg |
| MT110                | Medium, High Security Bolt Lock     | 35.0kg |
| MT120                | Large, High Security Bolt Lock      | 42.5kg |
| Electronic Pus       | sh Button Lock                      |        |
| MT100CL              | Small, Electronic Push Button Lock  | 26.0kg |
| MT110CL              | Medium, Electronic Push Button Lock | 35.0kg |
| MT120CL              | Large, Electronic Push Button Lock  | 42.5kg |

# **Dispensing Trolleys - Wooden, Standard Capacity**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines
  Lockable hinged lids...
  - Upper Assisted using a gas strut, providing document support
  - Lower Provides a work area, at height 980mm (max load 4.0kg)
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Hinged work flap provides a work area at height 930mm (max load 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Two slide out storage / waste bins
- 100mm pressed steel castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

- Overall...
  - Small 750 x 490 x 1130mm
  - Medium 890 x 490 x 1130mm
  - Large 1030 x 525 x 1130mm
- Internal...
  - Small 530 x 427 x 263mm sloping to 112mm at front
  - Medium 683 x 427 x 263mm sloping to 112mm at front
  - Large 805 x 461 x 263mm sloping to 112mm at front

#### **Standard Colours**

- Frame Grey
- Shelves / Dividers White
- Cabinet Spearmint
- Shelf / Flap Grey

- MT/1/6F Security Wall Clamp
- ZP005 Ring Buffers
- MT/LHR/S Laptop Holder Removable, Swing Arm
- MT/LHNR/S Laptop Holder Non-Removable, Swing Arm



| Item       | Description                         | Weight |
|------------|-------------------------------------|--------|
| High Secur | ity Bolt Lock                       |        |
| MT220      | Medium, High Security Bolt Lock     | 42.0kg |
| MT230      | Large, High Security Bolt Lock      | 50.0kg |
| Electronic | Push Button Lock                    |        |
| MT220CL    | Medium, Electronic Push Button Lock | 42.0kg |
| MT230CL    | Large, Electronic Push Button Lock  | 50.0kg |

# **Dispensing Trolleys - Wooden, High Capacity**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines
  Lockable hinged lids...
  - Upper Assisted using a gas strut, providing document support
  - Lower Provides a work area, at height 921mm (max load 4.0kg)
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Hinged work flap provides a work area at height 850mm (max load 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Two slide out storage / waste bins
- 125mm pressed steel castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

- Overall...
  - Medium 760 x 560 x 1290mm
  - Large 905 x 560 x 1290mm
- Internal...
  - Medium 530 x 461 x 513mm sloping to 137mm at the front
  - Large 683 x 461 x 513mm sloping to 137mm at the front

#### **Standard Colours**

- Frame Grey
- · Shelves / Dividers White
- Cabinet Spearmint
- Shelf / Flap Grey

- MT/1/6F Security Wall Clamp
- MT/LHR/S Laptop Holder Removable, Swing Arm
- MT/LHNR/S Laptop Holder Non-Removable, Swing Arm



# Features

- Trolleys suitable for the storage & dispensing of drugs & medicines
- Lockable hinged lid provides document support
- Fitted with a high security bolt lock (supplied with two keys)

**Dispensing Trolleys - Metal, Standard Capacity** 

- Hinged writing flap provides a work area at height 907mm (max capacity 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Two slide out storage / waste bins
- 100mm pressed steel castors
- When not in use, trolley should be secured to the building see optional  $\ensuremath{\mathsf{MT}}\xspace/1/6\ensuremath{\mathsf{F}}\xspace$

#### Dimensions (w x d x h)

- Overall...
  - Small 750 x 459 x 1112mm
  - Medium 891 x 459 x 1112mm
  - Large 1030 x 492 x 1112mm
- Internal...
  - Small 489 x 457 x 292mm sloping to 127mm at the front
  - Medium 706 x 457 x 292mm sloping to 127mm at the front
  - Large 832 x 490 x 292mm sloping to 127mm at the front

#### **Standard Colours**

- Frame Grey
- Cabinet White
- Shelf / Flap Spearmint

- MT/1/6F Security Wall Clamp
- ZP005 Ring Buffers
- MT/LHR/S Laptop Holder Removable, Swing Arm
- MT/LHNR/S Laptop Holder Non-Removable, Swing Arm

| Item  | Description                     | Weight |
|-------|---------------------------------|--------|
| MT130 | Small, High Security Bolt Lock  | 22.0kg |
| MT140 | Medium, High Security Bolt Lock | 25.0kg |
| MT150 | Large, High Security Bolt Lock  | 31.0kg |

# Laptop Holder - Removable, Swing Arm

#### Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MT range of dispensing trolleys
- Fitted to the rear right hand upright of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle, fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports
- \*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**

Paint - Grey



| Item     | Description                          | Weight |
|----------|--------------------------------------|--------|
| MT/LHR/S | Laptop Holder - Removable, Swing Arm | 3.7kg  |

# Laptop Holder - Non-Removable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MT range of dispensing trolleys
- Fitted to the rear right hand upright of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle fitted with a steel mounting bracket, fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports
- \*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**


Paint - Grey

| Item      | Description                              | Weight |
|-----------|--|--------|
| MT/LHNR/S | Laptop Holder - Non-Removable, Swing Arm | 5.0kg  |

fixed using either nylock or wing nuts

Laptop supported in a cradle fitted with a steel mounting bracket that can be

**Drug & Medicine Dispensing** 





#### **Single Door**

# High Security Bolt LockElectronic Push Button LockDescriptionWeightPTS/HSB/PA24PTS/EPB/PA24Patient Administration, 24 Compartments51.4kg

# **Pharmacy Trolleys - Patient Administration**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using individual compartments dedicated for each patient
- Two locking options...
  - · High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Supplied with 4 / 8 trays subdivided into equal compartments...
  - 24 Compartment 4 trays subdivided into 6 equal compartments, each compartment - 190mm x 132mm x 220mm (5.5 Litres). Label holders provided for each compartment
  - 32 Compartment 8 trays subdivided into 4 equal compartments, each compartment - 192mm x 143mm x 220mm (6.1 Litres). Label holders provided for each compartment
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional  $\rm MT/1/6F$

#### Dimensions (w x d x h)

#### **Standard Colours**

• Work Surface - Spearmint

· Paint - Grev

- Single Door 627 x 579 x 1165mm
- Double Door 820 x 579 x 1165mm

#### **Options / Accessories**

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm



PTS/EPB/PA24



#### **Double Door**

| High Security Bolt Lock | Electronic Push Button Lock | Description                             | Weight |
|-------------------------|-----------------------------|---|--------|
| PTD/HSB/PA32            | PTD/EPB/PA32                | Patient Administration, 32 Compartments | 68.8kg |



Trays tilt to assist access to contents



PTD/HSB/PA32



| Item        | Description                              | Weight |
|-------------|--|--------|
| PTS/HSB/DR6 | High Security Bolt Lock, Six Drawers     | 53.0kg |
| PTS/EPB/DR6 | Electronic Push Button Lock, Six Drawers | 53.0kg |

# **Pharmacy Trolleys - Drawers**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using drawers
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- 3 adjustable door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

- Overall Dimensions 511 x 542 x 1072mm
- Drawer 335 x 295 x 113mm (internal)

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint

#### **Options / Accessories**

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- SEP/DIS Drawer Inserts, single door trolley refer to page 384



PTS/EPB/DR6



| ecurity Bolt Lock Electro | onic Push Button Lock | Description                 | Weight |
|---------------------------|-----------------------|-----------------------------|--------|
| Door                      |                       |                             |        |
| SB/LA16 PTS/E             | PB/LA16               | Sixteen LA Trays            | 39.0kg |
| e Door                    |                       |                             |        |
| SB/A18C12 PTD/E           | PB/A18C12             | Eighteen A & Twelve C Trays | 57.0kg |
| SB/LP24 PTD/E             | PB/LP24               | Twenty Four LP Trays        | 57.0kg |
| SB/A18C12 PTD/E           |                       | б ў                         |        |

# **Pharmacy Trolleys - Trays**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines
   using individual patient trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Trays supported on adjustable shelves
- Tray Sizes (w x d x h)...
  - LA 86 x 268 x 57mm (aluminium)
  - A 85 x 181 x 54mm (plastic)
  - C 129 x 181 x 54mm (plastic)
  - LP 115 x 266 x 95mm (plastic)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- Adjustable door shelves...
  - Single Door 3 door shelves
  - Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional  $\ensuremath{\mathsf{MT}}\xspace{-1.5mu}$

#### Dimensions (w x d x h)

- Single Door 511 x 542 x 1072mm
- Double Door 711 x 568 x 1072mm

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm





PTS/HSB/BP6 (frames not included)

# Pharmacy Trolleys - Blister Pack (MDS)

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using the blister packed monitored dosage system (MDS)
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Frames supported on three adjustable slide out shelves (frames not supplied)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- Adjustable door shelves...
  - Single Door 3 door shelves
  - Double Door 6 door shelves
- Fold over MDS dispensing surface
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

#### **Options / Accessories**

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- Single Door 511 x 542 x 1072mm
  Double Door 711 x 568 x 1072mm

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint



Tab prevents the shelf from being removed

| High Security Bolt Lock | Electronic Push Button Lock | Description               | Weight |
|-------------------------|-----------------------------|---------------------------|--------|
| Single Door             |                             |                           |        |
| PTS/HSB/BP6             | PTS/EPB/BP6                 | Blister Pack, Six Frames  | 38.0kg |
| Double Door             |                             |                           |        |
| PTD/HSB/BP9             | PTD/EPB/BP9                 | Blister Pack, Nine Frames | 42.0kg |

# Pharmacy Trolleys - Nomad<sup>®</sup> Cassettes

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using the Nomad® cassette monitored dosage system
- Two locking options...
  - · High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Cassettes supported on three adjustable shelves (cassettes not supplied)
- · Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- · Six adjustable door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

• Overall Dimensions 711 x 568 x 1072mm



- · Paint Grey
- Work Surface Spearmint

#### **Options / Accessories**

High Security Bolt Lock

PTD/HSB/NC56

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm

**Electronic Push Button Lock** 

PTD/EPB/NC56



sales@bristolmaid.com +44 (0)1258 484455

PTD/EPB/NC56

Weight

48.0kg

(-REER

16

Description

Nomad<sup>®</sup>, Fifty Six Cassettes

# Pharmacy Trolleys - Biodose® & MultiMeds



| High Security Bolt Lock | Electronic Push Button Lock | Description                                | Weight |
|-------------------------|-----------------------------|--|--------|
| Single Door             |                             |  |        |
| PTS/HSB/BD46            | PTS/EPB/BD46                | Biodose <sup>®</sup> & MultiMeds, 46 Trays | 37.3kg |
| Double Door             |                             |  |        |
| PTD/HSB/BD69            | PTD/EPB/BD69                | Biodose <sup>®</sup> & MultiMeds, 69 Trays | 52.0kg |

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using the Biodose<sup>®</sup> & MultiMeds monitored dosage system
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Trays supported on fixed runners (trays not supplied)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- · Adjustable door shelves...
  - Single Door 3 door shelves
  - Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional  $\ensuremath{\mathsf{MT}}\xspace{1/6}$

#### Dimensions (w x d x h)

- Single Door 581 x 568 x 1072mm
- Double Door 812 x 568 x 1072mm

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint

#### **Options / Accessories**

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm

sales@bristolmaid.com +44 (0)1258 484455



| High Security Bolt Lock | Electronic Push Button Lock | Description        | Weight |
|-------------------------|-----------------------------|--------------------|--------|
| Single Door             |                             |                    |        |
| PTS/HSB/MX40            | PTS/EPB/MX40                | Medinoxx, 40 Trays | 41.4kg |
| Double Door             |                             |                    |        |
| PTD/HSB/MX60            | PTD/EPB/MX60                | Medinoxx, 60 Trays | 55.2kg |

# **Pharmacy Trolleys - Medinoxx**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using the Medinoxx monitored dosage system (compatible with 7 x 4 medi-trays)
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Trays supported on fixed runners (trays not supplied)
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the building
- · Adjustable door shelves
  - Single Door 3 door shelves
  - Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

- Single Door 581 x 568 x 1072mm
- Double Door 812 x 568 x 1072mm

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm



PTS/EPB/MX40

# Laptop Holder - Removable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the PT range of single & double door pharmacy trolleys
- Fitted to an upright bolted to the rear left-hand corner of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle, fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

#### Standard Colours

• Paint - Grey



| Item     | Description          | Weight |
|----------|----------------------|--------|
| PT/LHR/S | Removable, Swing Arm | 4.7kg  |

# Laptop Holder - Non-Removable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the PT range of single & double door pharmacy trolleys
- Fitted to an upright bolted to the rear left-hand corner of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle fitted with a steel mounting bracket, fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg

**Standard Colours** 

· Paint - Grey

• Design allows easy access to laptop ports



Laptop supported in a cradle fitted with a steel mounting bracket that can be fixed using either nylock or wing nuts

| ço         |
|------------|
| S          |
| Medicine   |
| <u>0</u>   |
| ne         |
|            |
| Dispensing |
| ĕ          |
| <u></u>    |
| J.         |
|            |
|            |


Drug

| Item      | Description              | Weight |
|-----------|--------------------------|--------|
| PT/LHNR/S | Non-Removable, Swing Arm | 5.0kg  |

# **Fitting Your Security Clamp**



Securely fasten one end of the clamp to the wall



The 'hook shaped' end is trapped by the door when closed & locked



sales@bristolmaid.com

# **Security Clamp**

#### **Features**

- · Allows the trolley to be secured to the 'fabric of the building'
- · Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door/lid of the trolley when closed & locked
- Recommended wall fitting height approx 1000mm



**Important Note - Security Clamps** 

It is a requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area & broken into, but also encourages good nursing practice, giving the trolley a "permanent home". "Safe & Secure Handling of Medicines - Royal Pharmaceutical Society,

GB-March 2005"

| Item    | Description    | Cable Length | Weight |
|---------|----------------|--------------|--------|
| MT/1/6F | Security Clamp | 950mm        | 0.75kg |



| Item  | Description                       | Weight |
|-------|-----------------------------------|--------|
| SA005 | Distribution Trolley - 15 Drawers | 49.0kg |

Distribution Trolley - 18 Drawers

SA025

| rawers | 56.0kg |  |
|--------|--------|--|
|        |        |  |
|        |        |  |

#### www.bristolmaid.com

# **Distribution Trolleys - Lockable Drawers**

#### **Features**

- Trolleys suitable for the distribution of drugs & medicines from pharmacy to the ward & also with self-administration / patient own drug schemes
- Provides secure drug delivery using drawers which when locked are securely fastened to the trolley
- Drawers are fitted with a CAM lock (supplied with two keys)
- Max load per drawer 3kg
- Supplied as standard with locks to differ, allowing each ward to have its own key & be responsible for removing the drawer from the trolley (remaining key should be kept at the pharmacy)
- Master keys available (additional cost)
- Upper work surface
- · Push handles at either end
- · Document pocket fitted to end of trolley
- Low level corner bumpers offer protection to both the trolley & fabric of the building
- 125mm pressed steel castors. 2 off braking

#### Dimensions (w x d x h)

- SA005 1026 x 523 x 1002mm
- SA025 1180 x 523 x 1002mm
- Drawer 145 x 423 x 230mm (Internal)

#### **Standard Colours**

- Paint Grev
- · Work Surface Spearmint

#### **Options / Accessories**

· Master Keys - contact us for details



# **Drug & Medicine Storage**

- Drug & Medicine Cabinets Shelves & Door Shelves 146
  - Drug & Medicine Cabinets Deep Shelves 148
  - Drug & Medicine Cabinets Blister Pack (MDS) 150
    - Controlled Drug Cabinets 152
- Auditable Access Control for Medicines Management 154
  - Drug & Medicine Cupboards 156
    - Options 158
- Patients / Residents Own Medication Cabinets Self-Administration 160
- Patients / Residents Own Medication Cabinets Blister Pack (MDS) 163
  - Tablet Counter & Storage Boxes 164
  - Pharmacy & Vaccine Refrigerators 165



# Drug & Medicine Cabinets - Shelves & Door Shelves

#### **Features**

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Adjustable & removable shelves (190mm deep) & door shelves (75mm deep)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems
- Compatible with our range of CD cabinets for the storage of Controlled Drugs (see page 152)

#### **Options / Accessories**

- DGL Digital Push Button Lock (Not available on independent locking cabinets)
- SL Suited Locks (contact us for more information)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/800/ST Sloping Top (800mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

#### **Standard Colours**

• Paint - White

| Right I | Hand Opening | Left Ha  | nd Opening  |                                |                        | * Transformer casing adds 105mm to | overall height |
|---------|--------------|----------|-------------|--------------------------------|------------------------|------------------------------------|----------------|
| Light*  | No Light     | Light*   | No Light    | Overall Dimensions (w x d x h) | Shelves / Door Shelves | Compatible CD Cabinets             | Weight         |
| PC536   | PC536/LX     | PC536/LH | PC536/LH/LX | 500 x 300 x 600mm              | 3/3                    | CD005, CD010                       | 29.9kg         |
| PC539   | PC539/LX     | PC539/LH | PC539/LH/LX | 500 x 300 x 900mm              | 4 / 4                  | CD005, CD010, CD015                | 34.0kg         |
| PC636   | PC636/LX     | PC636/LH | PC636/LH/LX | 600 x 300 x 600mm              | 3/3                    | CD005, CD010                       | 33.9kg         |
| PC639   | PC639/LX     | PC639/LH | PC639/LH/LX | 600 x 300 x 900mm              | 4 / 4                  | CD005, CD010, CD015                | 38.0kg         |

#### Single Door Cabinets

+44 (0)1258 484455

sales@bristolmaid.com

www.bristolmaid.com



PC836/LX

#### \*Transformer casing adds 105mm to overall height \*\*CD010 not compatible with PC836/IND & PC836/IND/LX **Double Door Cabinets** Single Lock (Secures Both Doors) Each Door Fitted With A Lock \*\*\*CD015 not compatible with PC839/IND & PC839/IND/LX Overall Dimensions (w x d x h) Shelves / Door Shelves Compatible CD Cabinets Weight PC836/LX PC836 PC836/IND PC836/IND/LX 800 x 300 x 600mm 6/5 CD005, CD010\*\* 40.9kg PC839 PC839/LX PC839/IND PC839/IND/LX 800 x 300 x 900mm 8/8 CD005, CD010, CD015\*\*\* 45.9kg PC1036 PC1036/LX PC1036/IND PC1036/IND/LX 1000 x 300 x 600mm 6/6 CD005, CD010 48.9kg PC1039 PC1039/LX PC1039/IND PC1039/IND/LX 1000 x 300 x 900mm 8/8 CD005, CD010, CD015 53.9kg

+44 (0)1258 484455







sales@bristolmaid.com

#### **Single Door Cabinets**

# **Drug & Medicine Cabinets - Deep Shelves**

#### **Features**

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Adjustable & removable deep shelves (250mm deep)
- Key profile unique to Hospital Metalcraft Ltd. offering 250,000 differs
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- · Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems
- · Compatible with our range of CD cabinets for the storage of Controlled Drugs (see page 152)

#### **Options / Accessories**

- DGL Digital Push Button Lock (Not available on independent locking cabinets)
- SL Suited Locks (contact us for more information)
- PC/300/ST Sloping Top (300mm Wide Cabinets)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/800/ST Sloping Top (800mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

#### Standard Colours

· Paint - White

| Right Ha | and Opening | Left Ha     | nd Opening     |                                |         | * Trans                | former casing adds 105mm to | o overall height |
|----------|-------------|-------------|----------------|--------------------------------|---------|------------------------|-----------------------------|------------------|
| Light*   | No Light    | Light*      | No Light       | Overall Dimensions (w x d x h) | Shelves | Compatible CD Cabinets | Nomad Capacity              | Weight           |
| PC336/FS | PC336/FS/LX | PC336/FS/LH | PC336/FS/LH/LX | 300 x 300 x 600mm              | 3       | CD005                  | 13 x Cassettes              | 16.9kg           |
| PC536/FS | PC536/FS/LX | PC536/FS/LH | PC536/FS/LH/LX | 500 x 300 x 600mm              | 3       | CD005, CD010           | 24 x Cassettes              | 29.9kg           |
| PC636/FS | PC636/FS/LX | PC636/FS/LH | PC636/FS/LH/LX | 600 x 300 x 600mm              | 3       | CD005, CD010           | 32 x Cassettes              | 33.9kg           |
| PC539/FS | PC539/FS/LX | PC539/FS/LH | PC539/FS/LH/LX | 500 x 300 x 900mm              | 4       | CD005, CD010, CD015    | 35 x Cassettes              | 34.9kg           |
| PC639/FS | PC639/FS/LX | PC639/FS/LH | PC639/FS/LH/LX | 600 x 300 x 900mm              | 4       | CD005, CD010, CD015    | 44 x Cassettes              | 38.9kg           |



#### **Double Door Cabinets**

| Doub   | le Door Cabin       | ets             |                       | _                              |         | *Transformer<br>**CD010 not compatible with | casing adds 105mm to or<br>PC836/FS/IND & PC836 |            |
|--------|---------------------|-----------------|-----------------------|--------------------------------|---------|---|---|------------|
| Single | e Lock (Secures Doo | ors) Each Doc   | or Fitted With A Lock |                                |         | *** CD015 not compatible with               | n PC839/FS/IND & PC839                          | /FS/IND/LX |
| Light* | No Light            | Light*          | No Light              | Overall Dimensions (w x d x h) | Shelves | Compatible CD Cabinets                      | Nomad Capacity                                  | Weight     |
| PC836  | S/FS PC836/FS/L     | X PC836/FS/IND  | PC836/FS/IND/LX**     | 800 x 300 x 600mm              | 6       | CD005, CD010**                              | -   | 40.9kg     |
| PC839  | /FS PC839/FS/L      | X PC839/FS/IND  | PC839/FS/IND/LX***    | 800 x 300 x 900mm              | 8       | CD005, CD010, CD015***                      | -   | 45.9kg     |
| PC103  | 86/FS PC1036/FS/    | LX PC1036/FS/IN | D PC1036/FS/IND/LX    | 1000 x 300 x 600mm             | 6       | CD005, CD010                                | 58 Cassettes                                    | 48.9kg     |
| PC103  | 39/FS PC1039/FS/    | LX PC1039/FS/IN | D PC1039/FS/IND/LX    | 1000 x 300 x 900mm             | 8       | CD005, CD010, CD015                         | 80 Cassettes                                    | 53.9kg     |

sales@bristolmaid.com



PC536/MDS (Frames not included)

## **Drug & Medicine Cabinets - Blister Pack (MDS)**

#### **Features**

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Hooks at rear of cabinet support frames for the popular blister pack Monitored Dosage System (MDS)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/ EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems

#### **Standard Colours**

· Paint - White

| Right Hand Opening Left Hand Opening |              | * Transform  | er casing adds 105mm to | o overall height               |              |        |
|--------------------------------------|--------------|--------------|-------------------------|--------------------------------|--------------|--------|
| Right Ha                             | and Opening  | Left Har     | nd Opening              |                                |              |        |
| Light*                               | No Light     | Light*       | No Light                | Overall Dimensions (w x d x h) | MDS Capacity | Weight |
| PC536/MDS                            | PC536/MDS/LX | PC536/LH/MDS | PC536/LH/MDS/LX         | 500 x 300 x 600mm              | 2 x Frames   | 26.9kg |
| PC539/MDS                            | PC539/MDS/LX | PC539/LH/MDS | PC539/LH/MDS/LX         | 500 x 300 x 900mm              | 4 x Frames   | 31.4kg |
| PC639/MDS                            | PC639/MDS/LX | PC639/LH/MDS | PC639/LH/MDS/LX         | 600 x 300 x 900mm              | 6 x Frames   | 35.2kg |

#### **Single Door Cabinets**

#### **Options / Accessories**

Double Door Cabinets
Single Lock (Secures Both Doors)

• DGL - Digital Push Button Lock (Not available on independent locking cabinets)

PC1039/MDS PC1039/MDS/LX PC1039/IND/MDS PC1039/IND/MDS/LX 1000 x 300 x 900mm

- SL Suited Locks (contact us for more information)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)



52.1kg

8 x Frames



# Features Exceeds the statutory requires of Drugs (Safe Custody) Reg CD140 & CD145) BS 3621 anti-pick & drill, key identification card) Single-point locking Key profile unique to Hospita Compatible with CLIQ®, integ cylinders & keys, combining Can be supplied with a microsystem by others (CD905), ow with volt-free contact to allow -240y ~ 50Hz -60Hz 0.3A

# **Controlled Drug Cabinets - Single Point Locking**

- Exceeds the statutory requirement for the storage of controlled drugs "The Misuse of Drugs (Safe Custody) Regulations 1973" (with the exception of CD040, CD045, CD140 & CD145)
- BS 3621 anti-pick & drill, key retaining six pin cylinder lock (supplied with 2 keys & identification card)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs per profile
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Can be supplied with a micro-switch to allow connection to remote lighting/nurse call system by others (CD905), or supplied with a 12v DC LED lighting/warning system with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer casing adds 105mm to the overall height
- 1.5mm (16 gauge) mild steel, seamless welded construction
- Mounting holes provided in back (fixings not included)
- Adjustable & removable shelves
- Security of models CD005, CD010 & CD015 can be enhanced by concealing inside one of our 'Drug & Medicine Cabinets'
- Modular design allows for easy interchange with HTM 63 storage systems
- Supplied right hand hinging as standard, left hand available on request

#### **Standard Colours**

• Paint - White

#### **Options / Accessories**

- CD905 Micro-switch to connect to remote light box / alarm
- CD910 / CD911 / CD913 12v DC LED Lighting/Warning System (for compatibility see tables below)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

R/H Door L/H Door Overall Dimensions (w x d x h) Shelves Lighting / Warning System (Optional) Weight CD005 210 x 270 x 300mm CD910 -1 8.0kg CD010 CD110 CD910 335 x 270 x 300mm 1 9.0kg CD015 CD115 335 x 270 x 550mm 2 CD910 17.0kg CD020 CD120 500 x 300 x 550mm 2 CD910 20.0kg CD025 CD125 500 x 300 x 850mm 3 CD910 37.0kg

CD125







CD135

| R/H Door | L/H Door | Overall Dimensions (w x d x h) | Shelves | Lighting / Warning System (Optional) | Weight |
|----------|----------|--------------------------------|---------|--------------------------------------|--------|
| CD030    | CD130    | 500 x 450 x 850mm              | 3       | CD911                                | 42.0kg |
| CD035    | CD135    | 500 x 450 x 1250mm             | 4       | CD911                                | 59.0kg |
| CD045    | CD145    | 500 x 465 x 1900mm             | 6       | CD913                                | 91.0kg |



sales@bristolmaid.com

# Auditable Access Control for Medicines Management



#### What Is CLIQ<sup>®</sup>?

CLIQ<sup>®</sup> is the very latest in user-friendly, web-based security management. Part of the Protec<sup>2</sup> CLIQ<sup>®</sup> from Abloy range, CLIQ<sup>®</sup> is designed for remote sites, enabling the commissioning, programming and full administration of security management centrally via web-based software.

#### How Does It Work?

CLIQ<sup>®</sup> integrates electronics into mechanical cylinders and keys, providing a security solution that can be managed remotely. In short, this works like an electronic handshake: when the key is inserted it engages with the cylinder, sending a random number back to the key. The key replies with an encrypted response which is decrypted by the cylinder and if valid the cylinder then allows access. It is important to remember that the mechanical combination must match as well as the digital combination. CLIQ<sup>®</sup> allows a single key to be used with multiple cabinets or trolleys.

#### **System Features**

- Secure Provides dual security through the patented Protec<sup>2</sup> CLIQ<sup>®</sup> from Abloy mechanical disc cylinders and keys along with DES encryption, a recognised electronic standard
- · Easily Programmable Keys can be added quickly and programmed to specific time and date schedules
- · Fully Auditable- Audit trails can be accessed from both the cylinder and key
- · Reduces Keys One key can be programmed to open multiple cabinets or trolleys
- The CLIQ® environment works in an Amazon cloud hosted server
- · maintained and supported by Abloy allowing remote user access via a secure
- · connection from anywhere in the world
- Cost Effective Low installation & maintenance costs along with a long lifecycle (requires no wires or system updates, power for the cylinder is provided by the key)
- Easily Retrofitable CLIQ<sup>®</sup> can be used with a variety of locking solutions including euro lock profile cases, padlocks, camlocks & furniture locks

#### **System Benefits**

- · Save valuable nursing time which can be used to deliver patient care
- · Reduces delay in administering medication
- · Improves drug & medicine security
- · No more unlocked drug & medicine cabinets
- · Eliminates the need to replace locks as a result of lost keys
- · One key can be programmed to open multiple cabinets or trolleys
- The key is made personal to that member of staff, resulting in increased ownership and a reduction in the number of keys lost
- Easy to apply remote management, saves money and has the ability to cover several hospitals
- Time and date based access rights, with a full audit trail allow for full visibility of when visits were made
- New keys can be added to the system quickly and easily, removal of lost ones is hassle-free
- · Reduction in time spent looking for the person holding the key

### **Process Flow**

A JA



#### **During The Shift**

#### Operation

The key will only open cabinets where access permissions have been granted. Each time the key is used in a cabinet fitted with an PROTEC<sup>2</sup> cylinder an audit trail is created on both the key and cylinder



#### **Start of Shift**

#### Key Validation & Authorisation

Recurring validation forces users to place their key into the Wall PD. User authorisations are uploaded to the key



#### End of Shift

#### Automated Data Transfer

Key is placed into the Wall PD. Previous key transactions are automatically downloaded from the key and stored in the audit trail archive

# **System Components**



#### **CLIQ®** Kev

- · Provides dual security through the patented PROTEC<sup>2</sup> mechanical disc cylinders and keys along with DES encryption, a recognised electronic standard
- Thin, light but incredibly strong key

#### **PROTEC<sup>2</sup>** Lock Barrel

A bump proof mechanism

as well as padlocks

· Exceptionally smooth key entry and operation




· Can be master keyed and keyed alike with cylinders, cam and furniture locks



#### Wall Programming Device

- · Updates keys remotely (linking the server and users keys)
- Offers high level data encryption
- Allows the use of recurring validation, forcing users to validate their keys at the start of their shift whilst automatically downloading audit trails which are relayed to the server

#### Traka21



- · Keys are accessed using PIN number technology
- Holds up to 21 keys or keysets



#### **Questions & Answers** Answers to common questions about CLIQ®

#### O How secure is it?

- A CLIQ<sup>®</sup> uses a DES encryption (Data Encryption Standard) so it is just as safe as online banking
- Q Can I use this with an existing Ablov master key system?
- A Yes: CLIQ<sup>®</sup> can be integrated with an existing PROTEC<sup>2</sup> master key system
- Q What is the benefit?
- A CLIQ<sup>®</sup> will give you some of the benefits of access control, such as adding, deleting and amending access rights, take audit trails from keys and cylinders as well as setting the key to expire on a specific date at a specific time. Also, lost or stolen key can be easily be deleted electronically

#### Q Is the system expensive to install?

A The system has a low installation & maintenance cost along with a long life cycle. CLIQ<sup>®</sup> cylinders are designed to fit into Euro Profile mechanical lock cases. No wiring is required to the cabinet. The wall programming device requires an Ethernet connection and power through a 12-24V adapter or alternatively power over Ethernet

#### Q Who maintains the system?

- A System maintenance is normally assigned to a member of the hospital staff who will be responsible for programming/issuing of keys along with interrogation of audit trails
- Q How many audit trails can be collected from a single key?
- A The programmable user key stores the last 1.800 transactions and the last 20 foreign transactions
- Q How large is the cylinder memory?
- A 200 key groups; each key group holding 65,000 individual unique key combinations. This provides up to 13 million unique possible key combinations in one cylinder and can retain a black list of 3.000 keys
- Q Can I take an audit trail from a cylinder?
- A Yes, the cylinder stores an audit trail of the last 1,000 key insertions

#### Q What happens if a key is lost?

A We recommend that time schedules are set to each key, so if the key is not returned it will expire. Keys can be removed from the cylinder with the programming key or directly from the system. If subsequently found keys can be added back and re-issued. In the event that an audit trail is required for a lost key. lock cylinders can be interrogated using a programming key (each cylinder stores an audit trail of the last 1,000 key insertions)

#### Q How long do the batteries last?

- A The battery is located in the key, not the cylinder. Under normal usage a battery will operate up to 20.000 cvcles or 2 vears
- Q How is the information sent?
- A Via the internet to a wall mounted PD
- Q Can the information be taken off the Wall mounted PD by anyone?
- A No, there is no information stored on the PD, it is stored on a server until the correct key is inserted and only then is the information securely transferred
- Q Is it possible to access the system remotely?
- A The CLIQ<sup>®</sup> environment works in an Amazon cloud hosted server environment ISO27001 Accredited. maintained and supported by Abloy allowing remote user access via a secure connection from anywhere in the world.
- Q Where is the server?
  - A Amazon Republic Of Ireland
- Q Can I host this myself?
- A Yes, but there would be additional costs
- Q What if there is no internet connection on a particular site?
- A mobile PD is available which will work on the GSM network
- Q How can I better control the issuing of keys?
- A Through the use of Traka21 (key safe) keys can be automatically assigned to specific users. Keys are issued using either the staff members own security card or via pin number technology.
- Q Is it possible to integrate the CLIQ<sup>®</sup> key with other systems?
- A Unique RFID tag identifiers can be inserted into the key allowing it to be used with other MIFARE card based RFID systems
  - Self-Administration Schemes
  - Medicines Management
  - Patient Records



# **Drug & Medicine Cupboards - Single Door**

#### **Features**

- Independently tested to BS2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Single door cupboards available as either 870 / 1870mm high
- Available with either...
  - Plastic Trays
  - Shelves
  - Shelves & Dividers
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- Worksurface supplied with 870mm cupboards (removable to assist with positioning under existing worktops)
- 0.9mm (20 Gauge) mild steel welded construction
- Modular design allows for easy interchange with HTM 63 storage systems

#### Dimensions (w x d x h)

- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal

#### **Standard Colours**

· Paint - White

#### **Options / Accessories**

- SL Suited Locks (contact us for more information)
- DGL Digital Push Button Lock
- MC/JK/870 Joining Kit & Full Width Worktop (required when connecting 870mm cupboards)
- MC/JK/1870 Joining Kit (required when connecting 1870mm cupboards)
- MC/ST Sloping Top (1870mm cupboards)
- MC/DIV Set of 16 Dividers

| R/H Hinged    | L/H Hinged    | Overall Dimensions (w x d x h) | Plastic Trays / Shelves | Weight |
|---------------|---------------|--------------------------------|-------------------------|--------|
| MC/R/870/3S2D | MC/L/870/3S2D | 500 x 536 x 870mm              | Plastic Trays           | 45.6kg |
| MC/R/870/FS   | MC/L/870/FS   | 500 x 536 x 870mm              | Shelves                 | 37.0kg |



MC/R/1870/7S5D MC/L/1870/7S5D

MC/L/1870/FS

MC/L/1870/FSD

MC/R/1870/FS

MC/R/1870/FSD





#### **Plastic Trays**

- Interchangeable with our popular Caretray range
  - 870mm Cupboards 3 x 100mm, 2 x 150mm trays
  - 1870mm Cupboards 7 x 100mm, 5 x 150mm trays
  - Supplied with 7 x 100mm & 5 x 150mm trays
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray



- Adjustable at 50mm increments
  - 870mm Cupboards 3 x shelves
  - 1870mm Cupboards 8 x shelves
- Suitable for the storage of MDS frames (if required)
  - 870mm Cupboards 6 x frames
  - 1870mm Cupboards 16 x frames



#### **Shelves & Dividers**

- Adjustable at 50mm increments
- Removable dividers 100mm high adjustable at 25mm increments
  - 870mm Cupboards N/A
  - 1870mm Cupboards 8 x shelves & 16 x dividers

sales@bristolmaid.com

**Plastic Trays** 

Shelves & Dividers

Shelves

Overall Dimensions (w x d x h)

500 x 536 x 1870mm

500 x 536 x 1870mm

500 x 536 x 1870mm

88.6kg

69.0kg

79.0kg

Plastic Trays / Shelves / Shelves & Dividers

# **Sloping Top**

#### Features

- Compatible with 300/500/600/800 & 1000mm wide PC cabinets & 500mm wide MC cupboards
- · Prevents items being stored on the top of the cabinet & assists cleaning
- Increases the overall height of the cabinet by 193mm
- Void within the top provides a space to locate the transformer for the lighting/warning system
- · Access to the transformer provided via a removable plate

#### **Standard Colours**

· Paint - White


# Warning Light System

#### **Features**

- Suitable for use with Bristol Maid CD Controlled Drug cabinets
- · Provides a visual indication when the cabinet is open
- LED illuminates...
  - Green when the door is closed
  - Red when the door is open
  - Flashes red when open for 2 1/2 minutes
- Operates using a door activated micro-switch connected to a 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height



| Item     | Description  |
|----------|--|
| CD Range | Compatible   |
| CD910    | 12v DC LED Lighting/Warning System - Refer To Page 152 For Compatibility |
| CD911    | 12v DC LED Lighting/Warning System - Refer To Page 153 For Compatibility |
| CD913    | 12v DC LED Lighting/Warning System - Refer To Page 153 For Compatibility |



| Item       | Description             | Weight |
|------------|-------------------------|--------|
| ntem       | Description             | weight |
| PC/300/ST  | 300mm Wide PC Cabinets  | 1.3kg  |
| PC/500/ST  | 500mm Wide PC Cabinets  | 1.5kg  |
| PC/600/ST  | 600mm Wide PC Cabinets  | 1.6kg  |
| PC/800/ST  | 800mm Wide PC Cabinets  | 1.7kg  |
| PC/1000/ST | 1000mm Wide PC Cabinets | 1.8kg  |
| MC/ST      | 500mm Wide MC Cupboards | 2.9kg  |

# **Joining Kits - MC Range of Cabinets**

#### **Features**

- · Allows two units to be joined together
- 870mm joining kit supplied with a full width worktop (2 units)



| Item       | Description                      | Compatibility     |
|------------|----------------------------------|-------------------|
| MC/JK/870  | Joining Kit & Full Width Worktop | 870mm High Units  |
| MC/JK/1870 | Joining Kit                      | 1870mm High Units |

# **Remote Light Box**

#### **Features**

- Connects to all Bristol Maid PC, MS & CD cabinets fitted with a light)
- Provides a remote visual indication when the cabinet is open
- · Models available for surface or flush mounting
- Units fitted without a light require a separate micro-switch & a 12v DC supply



| Item    | Description                        | <b>Overall Dimensions</b> |
|---------|------------------------------------|---------------------------|
| PC/RLBS | Remote Light Box - Surface Mounted | 100 x 50 x 150mm          |
| PC/RLBF | Remote Light Plate - Flush Mounted | 100 x 50 x 150mm          |

# **Digital Push Button Lock**

#### **Features**

- Requires a five digit code to operate (set by customer)
- Cabinets fitted with this lock will not conform to BS2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises



Digital Push Button Lock

# Dividers - MC Range of Cabinets

#### **Features**

DGL

- Spare / replacement dividers for MC cupboards supplied with shelves & dividers
- 100mm high dividers locate into slots on the shelves, adjustable at 25mm increments
- · Compatible with 870mm & 1870mm high units



#### Shelf not supplied

| Item   | Descrip  |
|--------|----------|
| MC/DIV | Set of 1 |

| Description        |  |
|--------------------|--|
| Set of 16 Dividers |  |

# **Micro Switch**

#### **Features**

• Micro-switch to allow connection to



CD905

#### Item Description

CD905Micro-switch (CD Range Compatible)PC905Micro-switch (PC Range Compatible)

# **Buzzer Alarm**

#### Features

- · Replacement circuit board with buzzer
- Compatible for PC, MS & CD cabinets fitted with a light unit
- Buzzer pulses on/off after the door has been open for 2 ½ minutes
- Factory fit only to cabinets fitted with lights

| Item | Description  |
|------|--------------|
| PCB  | Buzzer Alarm |









SAC/313/S/VP

#### Features

- Suitable for use with...
  - Drugs retained in original packs
  - Individual MDS packs with dividers & name cards
- Nomad cassettes
- Roll type blister packs
- Available as side or bottom hinged fitted with either a CAM or electronic push button lock
  - CAM Supplied with two keys & uses 40 different master suites (for security we recommend that each ward/wing operates a different series, master keys available at additional cost)
  - Electronic Push Button Battery operated
- Available with a solid door or with an inset perspex visibility panel
- 0.9mm (20 gauge) mild steel, spot welded construction
- Mounting holes in back/bottom (fixings not included)
- Range of mounting plates (optional)
- Shelves not included (optional)
- Supplied with a self-adhesive name card holder (length 60mm)

#### **Standard Colours**

• Paint - White

#### **Electronic Push Button Lock**

- Uses a 10-button keypad
- Easy programming & code change
- Four operating codes with differing levels of authority - Master / Sub-Master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code



|          |            | yea                      |                  |                          |                                |                     |        |
|----------|------------|--------------------------|------------------|--------------------------|--------------------------------|---------------------|--------|
| CAM Lock |            | Electric                 | Push Button Lock |                          |                                |                     |        |
|          | Solid Door | Perspex Visibility Panel | Solid Door       | Perspex Visibility Panel | Overall Dimensions (w x d x h) | Shelf Compatibility | Weight |
|          | SAC/213/S  | SAC/213/S/VP             | SAC/213/S/CL     | SAC/213/S/CL/VP          | 210 x 155 x 315mm              | SAC/FS/2            | 2.6kg  |
|          | SAC/313/S  | SAC/313/S/VP             | SAC/313/S/CL     | SAC/313/S/CL/VP          | 310 x 155 x 315mm              | SAC/FS/3            | 3.4kg  |
|          | SAC/413/S  | SAC/413/S/VP             | SAC/413/S/CL     | SAC/413/S/CL/VP          | 410 x 155 x 315mm              | SAC/FS/4            | 4.1kg  |

#### **Side Hinged**

| sales@bristolmaid.com | +44 (0)1258 484455 |
|-----------------------|--------------------|

4455 www.bristolmaid.com

#### Options

- SAC/FS/2 Shelf, 200mm Wide
- SAC/FS/3 Shelf, 300mm Wide
- SAC/FS/4 Shelf, 400mm Wide
- SAC/WP Mounting Plate, Wall
- SAC/LT Mounting Plate, Locker Top



#### **Bottom Hinged**

CAM Lock **Electric Push Button Lock** SAC/FS/2 2.6kg SAC/213/B SAC/213/B/VP SAC/213/B/CL SAC/213/B/CL/VP 210 x 155 x 315mm SAC/313/B SAC/313/B/VP SAC/313/B/CL SAC/313/B/CL/VP 310 x 155 x 315mm SAC/FS/3 3.4kg SAC/413/B SAC/413/B/VP SAC/413/B/CL SAC/413/B/CL/VP 410 x 155 x 315mm SAC/FS/4 4.1kg SAC/213/B/VP

sales@bristolmaid.com

+44 (0)1258 484455

www.bristolmaid.com

# **Shelves**

#### Features

- Flat shelves (100mm deep)
- Removable to assist cleaning, allow the storage of larger items

#### **Standard Colours**

Paint - White

# **Mounting Plates**

#### **Features**


- Compatible with all SAC cabinets using existing holes (fixings not provided)
- SAC/WP Cabinet secured to wall plate using a small thumbscrew. Allows the cabinet to be repositioned
- SAC/LT Mounts to upper surface of traditional bedside cabinets

#### Standard Colours



| Item     | Description                                 | Weight |
|----------|---|--------|
| SAC/FS/2 | Shelf - 200mm Wide Cabinets (SAC/213 Range) | 0.2kg  |
| SAC/FS/3 | Shelf - 300mm Wide Cabinets (SAC/313 Range) | 0.2kg  |
| SAC/FS/4 | Shelf - 400mm Wide Cabinets (SAC/413 Range) | 0.3kg  |

| Item   | Description | Overall Dimensions (w x d x h) | Weight |
|--------|-------------|--------------------------------|--------|
| SAC/WP | Wall        | 190 x 12 x 162mm               | 0.5kg  |
| SAC/LT | Locker Top  | 190 x 150 x 150mm              | 0.9kg  |



SAC/FS/4

sales@bristolmaid.com

```
www.bristolmaid.com
```

# Patients / Residents Own Medication Cabinets - Monitored Dosage, MDS Frames



#### **Features**

- Suitable for the storage of MDS frames with dividers & name cards (max 100mm depth)
- · Models available to hold one or two frames (not included)
  - Small One frame
  - Large Two frames
- Available with either a CAM or electronic push button lock
  - CAM Supplied with two keys, master keys available at additional cost
  - · Electronic Push Button Battery operated
- Mounting holes provided in back (fixings not included)
- Supplied with an adjustable shelf (SA045 & SA045/CL only)
- Constructed from melamine faced chipboard offering an extremely strong & durable, yet aesthetically pleasing appearance
- MDS frames not supplied

#### **Standard Colours**

• Body / Door - Beech Woodgrain

#### **Electronic Push Button Lock**

- Uses a 10-button keypad
- Easy programming & code change
- Four operating codes with differing levels of authority - Master / Sub-Master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

02

34

56

28

00

SA040 - frames not included

| CAM Lock | Digital Push Button Lock | Overall Dimensions (w x d x h) | MDS Capacity | Weight |
|----------|--------------------------|--------------------------------|--------------|--------|
| SA040    | SA040/CL                 | 325 x 245 x 420mm              | 1 x Frame    | 7.5kg  |
| SA045    | SA045/CL                 | 325 x 245 x 775mm              | 2 x Frames   | 12.1kg |

+44 (0)1258 484455

# **Tablet Counter**

#### **Features**

- Triangular stainless steel tablet counter
- Etched legend to interpret rows (max 22) into number of tablets

# **Storage Boxes - Lightweight**


#### **Features**

- Lightweight aluminium riveted construction
- Fitted with handle
- Single-point locking using a CAM lock (supplied with 2 keys)
- Cannot be wall mounted

#### **Standard Colours**

• Paint - White







| Item  | Description    | Dimensions (w x d x h) |
|-------|----------------|------------------------|
| PE/30 | Tablet Counter | 205 x 195 x 10mm       |

+44 (0)1258 484455

sales@bristolmaid.com

www.bristolmaid.com





Powered by

5RLDF0219

| CAM Lock          | Digital Push Button | Overall Dimensions (w x d x h) | Capacity   | Description          | Shelves | Weight |
|-------------------|---------------------|--------------------------------|------------|----------------------|---------|--------|
| Solid Door        |                     |                                |            |                      |         |        |
| 5RLDF0119         | 5RLDF0119DIGILOCK   | 450 x 510 x 538mm              | 36 Litres  | Bench/Wall Mounting* | 2       | 30.0kg |
| 5RLDF0219         | 5RLDF0219DIGILOCK   | 450 x 510 x 740mm              | 66 Litres  | Bench/Wall Mounting* | 3       | 34.0kg |
| 5RLDF0519         | 5RLDF0519DIGILOCK   | 600 x 600 x 835mm              | 150 Litres | Under Bench          | 3       | 51.0kg |
| Glass Door        |                     |                                |            |                      |         |        |
| 5RLDG0519         | 5RLDG0519DIGILOCK   | 600 x 600 x 835mm              | 150 Litres | Under Bench          | 3       | 54.0kg |
| *Integral Wall Mo | ounting Brackets    |                                |            |                      |         |        |
|                   |                     |                                |            | 44 (0) 4050 404455   |         |        |

#### Features

- Compliant with the Department of Health Green Book
- Recommendations in the General Pharmaceutical Council Standards
- Royal College of Veterinary Surgeons Guidelines
- Designed specifically for the safe storage of vaccines and pharmaceuticals
- Available as bench top, under bench, wall fitting or freestanding
- Temperature range +2°C to +8°C
- Second wandering temperature probe (can be placed anywhere within the refrigerator, ideal for pharmacists needing to comply with MHRA regulations)
- IntelliCold<sup>®</sup> touch screen controller
  - Automatic 24/7 temperature & fridge operation data recording (micro SD card required)
  - Min/max temperature control, display & recording
  - Displays temperature from both probes
  - High/low temperature alarm
  - Door open alarm
  - Mains failure alarm (requires batteries)
  - Alarm battery back-up
- Internal LED light
- · Access port & remote alarm contacts
- · Fan assisted cooling
- · Off cycle auto defrost
- Low GWP environmentally friendly refrigerant
- · Self-closing doors with door open alarm
- Available with either a CAM or digital push button lock
  - CAM Supplied with two keys
  - Digital Push Button Mains operated (uses a 4-digit pin)
- · Adjustable shelves
- Micro SD card (optional)
- UKAS traceable calibration service (optional)
- · Repeat alarm (optional)

#### **Standard Colours**

- Exterior White / Blue
- · Interior White



# Examination / Treatment Couches & Chairs

- Examination & Treatment Couches Variable Height 168
  - Examination & Treatment Couches Fixed Height 172
    - Couch Steps 173
      - Leg Rests 174
- Examination & Treatment Couches Variable Height, Bariatric 176
  - Examination & Treatment Chairs Variable Height 180
  - Examination & Treatment Chairs Variable Height, Bariatric 181
    - Bariatric Examination & Treatment Couches 192
      - Phlebotomy Chairs Variable Height 182
    - Phlebotomy Chairs Variable Height, Bariatric 183
      - Phlebotomy Chairs Fixed Height 184
        - Options 185



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Hydraulic in Aquamarine - EC1205/VAQ



|            |   |        | 1170mm  | 615mi  |
|------------|---|--------|---------|--------|
| Item       | Description                                   | Weight | TT7 Omm | 015111 |
| EC1205/VBB | Hydraulic                                     | 68.8kg |         |        |
| EC1210/VBB | Electric, Hand Switch                         | 68.8kg |         |        |
| EC1215/VBB | Electric, Foot Switch                         | 68.8kg |         |        |
| EC1220/VBB | Electric with Electric Backrest - Hand Switch | 71.1kg |         |        |
| EC1225/VBB | Electric with Electric Backrest - Foot Switch | 71.1kg |         |        |
|            |   |        |         |        |

# Examination & Treatment Couch Two Section

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

#### Dimensions (w x d x h)

• 630 x 1790 x 490 - 950mm

#### **Standard Colours**

- · Paint White
- Vinyl Bristol Blue

#### **Options**

- E02 Paper Roll Holder
- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E11 Lithotomy Stirrups (Pair)
- E12 Leg Supports (Pair)

630mm

# **Examination & Treatment Couch Two Section, Mobile**



#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- · Height adjustment from 490 950mm
- · Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- · Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinvl
- · Swing down side support rails fitted to either side
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

#### Dimensions (w x d x h)

• 630 x 1790 x 490 - 950mm

#### Standard Colours

- · Paint White
- · Vinyl Bristol Blue

#### Options

- · E02 Paper Roll Holder
- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- · E05 Breather Hole & Plug
- E06 Battery Back Up
- E11 Lithotomy Stirrups (Pair)
- E12 Leg Supports (Pair)

| Foot Switch | Alternative | E Colours avai | ilable at ex<br>gi.e. Electr | tra cost. Changic, Hand Switc                | ge suffix 'V   | EC1235/V                                      | oduct code |             | HI        | EIGHT                      |
|-------------|-------------|----------------|------------------------------|--|----------------|---|------------|-------------|-----------|----------------------------|
| Hand Switch | VBB         | VAQ            | VFR                          | VGH  | VSD            | VOB   | VLI        | VCH         | VCB       | VPE                        |
| 1170mr      |             | 615mm          |                              | Item<br>EC1230/V                             | D              | escription                                    |            |             |           | Weight<br>75.0kg           |
|             |             |                | 630mm                        | EC1235/V<br>EC1235/V<br>EC1240/V<br>EC1245/V | BB EI<br>BB EI | ectric, Hanc<br>ectric, Foot<br>ectric with E | Switch     | krest - Har | nd Switch | 75.0kg<br>75.0kg<br>75.3kg |

EC1250/VBB

Electric with Electric Backrest - Foot Switch

77.3kg



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric, Foot Switch in Onyx Black - EC1315/VOB



| Item       | Description                                   | Weight | 690mm | 470mm |
|------------|---|--------|-------|-------|
| EC1305/VBB | Hydraulic                                     | 70.2kg |       |       |
| EC1310/VBB | Electric, Hand Switch                         | 70.2kg |       |       |
| EC1315/VBB | Electric, Foot Switch                         | 70.2kg |       |       |
| EC1320/VBB | Electric with Electric Backrest - Hand Switch | 72.5kg |       |       |
| EC1325/VBB | Electric with Electric Backrest - Foot Switch | 72.5kg |       |       |

# Examination & Treatment Couch Three Section

#### Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch

Planned

Preventative

Maintenance Contact us for

details

- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

#### Dimensions (w x d x h)

• 630 x 1790 x 490 - 950mm

#### **Standard Colours**

- · Paint White
- Vinyl Bristol Blue

#### Options

- E02 Paper Roll Holder
- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole

VCB

VPE

630mm

Hand Switch

615mm

# Examination & Treatment Couch Three Section, Mobile

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- · Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- 125mm swivel castors featuring a central locking & steering facility
- · Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless & contoured upholstery using flame retardant & antimicrobial vinyl
- · Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

#### Dimensions (w x d x h)

• 630 x 1790 x 490 - 950mm

#### **Standard Colours**

- · Paint White
- Vinyl Bristol Blue

#### **Options**

- E02 Paper Roll Holder
- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E05 Breather Hole & Plug
- · E06 Battery Back Up



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Coastal Blue - EC1345/VCB



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC1330/VBB | Hydraulic                                     | 76.4kg |
| EC1335/VBB | Electric, Hand Switch                         | 76.4kg |
| EC1340/VBB | Electric, Foot Switch                         | 76.4kg |
| EC1345/VBB | Electric with Electric Backrest - Hand Switch | 78.7kg |
| EC1350/VBB | Electric with Electric Backrest - Foot Switch | 78.7kg |
|            |   |        |

Hand Switch

Foot Switch

630mm



#### **Features**

- Fixed height approx. 786mm
- Adjustable angle backrest 0 to +90 degrees using a scissor jack
- · Adjustable feet allow levelling on uneven floors
- Seamless upholstery using flame retardant anti-microbial vinyl
- Paper roll holder fitted to foot end (suitable for rolls max 510mm wide x 200mm diameter)
- Supplied boxed for self-assembly
- Maximum capacity 125kg

#### Dimensions (w x d x h)

• 610 x 1860 x 795mm

#### **Standard Colours**

- Paint White
- · Legs Chrome
- Vinyl Bristol Blue





sales@bristolmaid.com

Chairs లర **Examination / Treatment Couches** 

CE

PETAPDAN

www.bristolmaid.com

# Couch Steps One & Two Tread

#### Features

- Non-slip rubber feet & tread
- Tread dimension 460 x 250mm
- Tread Heights
  - One Tread 180mm
  - Two Tread 185 & 365mm
- One Step available with a handle for additional safety
- Handle height max. 1005mm
- Maximum capacity 125kg

#### **Standard Colours**

• Paint - White





| ltem  | Description                    | Overall Dimensions (w x d x h) | Weight |
|-------|--------------------------------|--------------------------------|--------|
| CS005 | Couch Step - One Tread         | 519 x 284 x 228mm              | 4.4kg  |
| CS010 | Couch Step - One Tread, Handle | 519 x 284 x 1005mm             | 6.0kg  |
| CS015 | Couch Step - Two Tread         | 519 x 599 x 405mm              | 9.6kg  |

www.bristolmaid.com





# Leg Rest - Adjustable

#### **Features**

- Independently adjustable legs
- Height & angle adjustment from 355 510mm
- Seamless upholstery using flame retardant anti-microbial vinyl
- Maximum capacity 25kg

#### Dimensions (w x d x h)

- Overall Dimensions 460 x 305 x 350 510mm
- Upholstered Pad 460 x 305mm

#### **Standard Colours**

- Paint White
- Vinyl Bristol Blue



| Item      | Description           | Weight |
|-----------|-----------------------|--------|
| CS/28/VBB | Leg Rest - Adjustable | 3.5kg  |

VOB

VLI

VCH

VCB

VSD

VPE

VBB

VAQ

VFR

VGH

# Other products that you may be interested in...

See Pages 52 & 53

Caretray Trolleys - Mild Steel, Standard Level, Single & Double Column







Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Peach - EC2215/VPE



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC2205/VBB | Electric, Hand Switch                         | 79.0kg |
| EC2210/VBB | Electric, Foot Switch                         | 79.0kg |
| EC2215/VBB | Electric with Electric Backrest - Hand Switch | 81.3kg |
| EC2220/VBB | Electric with Electric Backrest - Foot Switch | 81.3kg |

# **Examination & Treatment Couch Two Section, Bariatric**

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- · Height adjustment from 520 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

#### Dimensions (w x d x h)

• 850 x 1790 x 520 - 950mm

#### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

#### **Options**

- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E18 Battery Back Up
- · E19 Paper Roll Holder
- E21 Breather Hole & Plug

#### Planned Preventative Maintenance Contact us for details

1170mm

850mm

Hand Switch

615mm

VPF

Chairs

ංඊ

**Examination / Treatment Couches** 

# Examination & Treatment Couch Two Section, Mobile, Bariatric

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- · Height adjustment from 520 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- 230V  $\sim,$  50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinyl
- · Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

#### Dimensions (w x d x h)

• 850 x 1790 x 520 - 950mm

#### **Standard Colours**

- Paint White
- · Vinyl Bristol Blue

#### Options

- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

|   |  | CE                   |
|---|--|----------------------|
|   |  | /                    |
|   |  | HEIGHT<br>ADJUSTABLE |
|   | EC2225/VBB   |                      |
| ¢ | Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Foot Switch in Bristol Blue - EC2240/VBB |                      |
|   |  |                      |



| 1170mm |
|--------|
|        |

| 615mm |       |            |   |        |
|-------|-------|------------|---|--------|
|       |       | ltem       | Description                                   | Weight |
|       | 8     | EC2225/VBB | Electric, Hand Switch                         | 85.2kg |
|       | 850mm | EC2230/VBB | Electric, Foot Switch                         | 85.2kg |
|       | З     | EC2235/VBB | Electric with Electric Backrest - Hand Switch | 87.5kg |
|       |       | EC2240/VBB | Electric with Electric Backrest - Foot Switch | 87.5kg |

Foot Switch

Hand Switch



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC2305/VBB | Electric, Hand Switch                         | 81.8kg |
| EC2310/VBB | Electric, Foot Switch                         | 81.8kg |
| EC2315/VBB | Electric with Electric Backrest - Hand Switch | 84.1kg |
| EC2320/VBB | Electric with Electric Backrest - Foot Switch | 84.1kg |
|            |   |        |

# **Examination & Treatment Couch Three Section, Bariatric**

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- · Height adjustment from 520 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinvl
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

#### Dimensions (w x d x h)

850 x 1790 x 520 - 950mm

#### Standard Colours

- Paint White
- Vinvl Bristol Blue

#### Options

- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- · E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

Planned Preventative Maintenance Contact us for details

600mm

470mm

850mm

615mm

Chairs ංඊ **Examination / Treatment Couches** 

# **Examination & Treatment Couch Three Section, Mobile, Bariatric**

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- · Height adjustment from 520 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections (optional electric backrest)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- 230V  $\sim,$  50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinyl
- · Swing down side support rails fitted to either side
- · Frames prepared to accept optional accessories
- Maximum capacity 325kg

#### Dimensions (w x d x h)

• 850 x 1790 x 520 - 950mm

#### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

#### Options

- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug



- Foot Switch
- and the second s

| Hand | Switch |
|------|--------|
|------|--------|

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric, Foot Switch in Flame Red - EC2330/VFR

www.bristolmaid.com



| 470mm | 615mm |       |            |   |        |
|-------|-------|-------|------------|---|--------|
|       |       |       | Item       | Description                                   | Weight |
|       |       | 8     | EC2325/VBB | Electric, Hand Switch                         | 88.0kg |
|       |       | 850mm | EC2330/VBB | Electric, Foot Switch                         | 88.0kg |
|       |       | З     | EC2335/VBB | Electric with Electric Backrest - Hand Switch | 90.3kg |
|       |       |       | EC2340/VBB | Electric with Electric Backrest - Foot Switch | 90.3kg |

EC2340/VBB



HEIGHT

690mm



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric - Foot Switch in Peach - EC4315/VPE

| VBB     | VAQ        | VFR           | VGH         | VSD          | VOB        | VLI    | VCH  | VCB    | VPE |       |
|---------|------------|---------------|-------------|--------------|------------|--------|------|--------|-----|-------|
| Item    |            | Description   | 1           |              |            | Weight |      |        |     |       |
| EC4305/ | VBB        | Hydraulic     |             |              |            | 66.5kg | 360m | m 470m | m   | 820mm |
| EC4310/ | VBB        | Electric - Ha | ind Switch  |              |            | 66.5kg |      |        |     |       |
| EC4315/ | <b>VBB</b> | Electric - Fo | ot Switch   |              |            | 66.5kg |      |        |     |       |
| EC4320/ | VBB        | Electric with | Electric Ba | ackrest - Ha | and Switch | 68.8kg |      |        |     |       |
| EC4325/ | VBB        | Electric with | Electric Ba | ackrest - Fo | oot Switch | 68.8kg |      |        |     |       |
|         |            |               |             |              |            |        |      |        |     |       |

# Examination & Treatment Chair Drop End, Three Section

#### Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

#### Dimensions (w x d x h)

• 630 x 1710 x 490 - 950mm

#### **Standard Colours**

- Paint White
- Vinyl Bristol Blue

#### Options

- E02 Paper Roll Holder
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E09 Bowl

www.bristolmaid.com

630mm

- E10 Bowl Support
- E13 Phlebotomy Armrests (Pair)
- E14 Side Support Loops (Pair)



Foot Switch



Hand Switch

# Examination & Treatment Chair Drop End, Three Section, Bariatric

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- · Height adjustment from 520 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- · Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant antimicrobial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

#### Dimensions (w x d x h)

• 850 x 1710 x 520 - 950mm

#### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

#### **Options**

- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E09 Bowl
- E10 Bowl Support
- E13 Phlebotomy Armrests (Pair)
- · E14 Side Support Loops (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- · E21 Breather Hole & Plug



Foot Switch

Hand Switch







| 820mm |       |            |   |        |
|-------|-------|------------|---|--------|
|       |       | Item       | Description                                   | Weight |
|       | 85    | EC4330/VBB | Electric - Hand Switch                        | 76.7kg |
|       | 850mm | EC4335/VBB | Electric - Foot Switch                        | 76.7kg |
|       | 2     | EC4340/VBB | Electric with Electric Backrest - Hand Switch | 79.0kg |
|       |       | EC4345/VBB | Electric with Electric Backrest - Foot Switch | 79.0kg |
|       |       |            |   |        |

sales@bristolmaid.com +44

470mm

HEIGHT



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric - Foot Switch in Graphite - EC5315/VGH



| Item       | Description  | Weight   | 360mm  |
|------------|--|--|--|
| EC5305/VBB | Hydraulic  | 67.0kg   |  |
| EC5310/VBB | Electric - Hand Switch                               | 67.0kg   |  |
| EC5315/VBB | Electric - Foot Switch                               | 67.0kg   |  |
| EC5320/VBB | Electric with Electric Backrest - Hand Switch        | 69.3kg   |  |
| EC5325/VBB | Electric with Electric Backrest - Foot Switch        | 69.3kg   |  |
|            | EC5305/VBB<br>EC5310/VBB<br>EC5315/VBB<br>EC5320/VBB | EC5305/VBBHydraulicEC5310/VBBElectric - Hand SwitchEC5315/VBBElectric - Foot SwitchEC5320/VBBElectric with Electric Backrest - Hand Switch | EC5305/VBBHydraulic67.0kgEC5310/VBBElectric - Hand Switch67.0kgEC5315/VBBElectric - Foot Switch67.0kgEC5320/VBBElectric with Electric Backrest - Hand Switch69.3kg |

# Phlebotomy Chair Variable Height

#### Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (selfassembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

### Dimensions (w x d x h)

• 770 x 1710 x 490 - 950mm

#### **Standard Colours**

- Paint White
- Vinyl Bristol Blue

#### Options

E06 - Battery Back Up

470mm

820mm

630mm

# Phlebotomy Chair Variable Height, Bariatric

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- Height adjustment from 490 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- · Clearance beneath frame allows access for mobile hoists
- 230V  $\sim,$  50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (selfassembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- · Frames prepared to accept optional accessories
- Maximum capacity 325kg

#### Dimensions (w x d x h)

• 990 x 1710 x 490 - 950mm

#### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

#### Options


• E18 - Battery Back Up



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Onyx Black - EC5340/VOB



360mm 470mm 820mm



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC5330/VBB | Electric - Hand Switch                        | 77.2kg |
| EC5335/VBB | Electric - Foot Switch                        | 77.2kg |
| EC5340/VBB | Electric with Electric Backrest - Hand Switch | 79.5kg |
| EC5345/VBB | Electric with Electric Backrest - Foot Switch | 79.5kg |

Planned Preventative

Maintenance

Contact us for

details



| Item       | Description  | Seat Height | Weight |
|------------|--------------|-------------|--------|
| EC5000/VBB | Fixed Height | 485mm       | 15.0kg |

# **Phlebotomy Chair Fixed Height**

#### **Features**

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Certified for use up to 32 stone (203kg)
- Stable four leg design
- · Black finish tubular frame
- · Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrest
- · Armrest can be positioned on either side of the chair (additional armrest available if required)
- Maximum capacity 203kg

#### Dimensions (w x d x h)

• 620 x 580 x 800mm

#### **Standard Colours**

- · Paint Black
- · Vinyl Bristol Blue

#### **Options**

• E15 - Phlebotomy Armrest (Single)

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Fixed Height Phlebotomy Chair in Aquamarine - EC5000/VAQ



www.bristolmaid.com

Chairs

Examination / Treatment Couches &



# Paper Roll Holder

· Fitted to the head section

E02

 Suitable for rolls max 510mm wide x 200mm diameter



# Side Support Rails (Pair)

#### **E03**

- Swing down-side support rails
- Can be added to either one or both sides of the couch to provide additional safety for patients



# Cushioned Side Support Rails (Pair)

#### E04 - E04/U

- Vinyl upholstered swing down-side support rails
- Can be added to either one or both sides of the couch to provide additional safety for patients
- Upgrade option for factory fitted side support rails (E04/U)



# Breather Hole & Plug E05

- Breather hole located in backrest section to improve patient comfort when in the prone position
- Supplied with a plug / insert



# **Battery Back Up**

#### **E06**

- Allows the couch to be lowered in a controlled manner in the event of a power failure
- Provides enough charge for 1-3 cycles



# Height Adjustable Armrests (Pair)

#### **E07**

 Height adjustable easy clean moulded polyurethane armrests





### **Infusion Pole**

#### E08

- Height adjustable two hook infusion poleCannot be fitted in conjunction with side
- support rails E03, E04, E04/U



# **Bowl**

 360mm diameter 304 Grade stainless steel bowl



# **Bowl Support**

#### E10

- Support to hold 360mm diameter 304 Grade stainless steel bowl
- Order bowl separately



# Lithotomy Stirrups (Pair) E11

• Traditional lithotomy stirrups



# Leg Supports (Pair) E12

• Fully adjustable easy clean moulded polyurethane leg supports



# Phlebotomy **Armrests (Pair)** E13

· Fully adjustable easy clean moulded polyurethane armrests



# Side Support Loops (Pair)

 Side support loops to assist mount / dismount & provide patient security



# Phlebotomy Armrest (Single)

#### E15

- Fully adjustable easy clean moulded polyurethane armrest
- Compatible with EC5000



# **Battery Back Up**

#### E18 (Bariatric Couch)

- Allows the couch to be lowered in a controlled manner in the event of a power failure
- Provides enough charge for 1-3 cycles



# Paper Roll Holder

#### E19 (Bariatric Couch)

- Fitted to the head section
- Suitable for rolls max 510mm wide x 200mm diameter



# Breather Hole & Plug

#### E21 (Bariatric Couch)

- Breather hole located in backrest section to improve patient comfort when in the prone position
- Supplied with a plug / insert



# **General Purpose Trolleys**

- Stainless Steel Trolleys Heavy Duty 190
  - Stainless Steel Light Duty 191
  - Polymer Shallow & Deep Shelves 192
    - Mild Steel Laminate Shelves 193





# General Purpose Trolley Stainless Steel, Heavy Duty

#### Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with three sides up & one down (shelf size 920 x 535mm)
- Available with two or three shelves, shelf heights
  - Two shelf 280 & 860mm
  - Three shelf 220, 540 & 860mm
- Maximum capacity per shelf 75kg
- 125mm pressed steel castors

#### Dimensions (w x d x h)

- No Ring Buffers 1000 x 570 x 1020mm
- Ring Buffers 1030 x 600 x 1020mm

#### **Standard Colours**

Stainless Steel

#### Options

• C008 - 125mm Castors, Heavy Duty Pressed Steel, 2 x Braking, 2 x Non-Braking (set of 4) - Factory fitted

| Item          | Description   | Weight |
|---------------|---|--------|
| WT/37/2       | Stainless Steel - Heavy Duty, Two Shelves                 | 28.0kg |
| WT/37/3       | Stainless Steel - Heavy Duty, Three Shelves               | 37.5kg |
| WT/37/2/ZP010 | Stainless Steel - Heavy Duty, Two Shelves, Ring Buffers   | 28.2kg |
| WT/37/3/ZP010 | Stainless Steel - Heavy Duty, Three Shelves, Ring Buffers | 37.7kg |

# General Purpose Trolley Stainless Steel, Light Duty

#### Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed liquid retaining sunken shelves (shelf size 825 x 430mm)
- Available with two or three shelves, shelf heights
  - Two shelf 250 & 870mm
  - Three shelf 250, 560 & 870mm
- Maximum capacity per shelf 50kg
- 100mm pressed steel castors

#### Dimensions (w x d x h)

- No Ring Buffers 1100 x 540 x 1040mm
- Ring Buffers 1100 x 580 x 1040mm

#### **Standard Colours**

Stainless Steel



www.bristolmaid.com

| Item             | Description   | Weight |
|------------------|---|--------|
| WT/33/SS/2       | Stainless Steel - Light Duty, Two Shelves                 | 16.0kg |
| WT/33/SS/3       | Stainless Steel - Light Duty, Three Shelves               | 22.0kg |
| WT/33/SS/2/ZP005 | Stainless Steel - Light Duty, Two Shelves, Ring Buffers   | 16.2kg |
| WT/33/SS/3/ZP005 | Stainless Steel - Light Duty, Three Shelves, Ring Buffers | 22.2kg |

# **General Purpose Trollev Polymer, Shallow Shelves**

#### **Features**

- Easy to clean corrosion resistant high-density polyethylene shelves (not suitable for use with automated washers)
- 11mm raised edge retains spills & prevents items from sliding off during transport
- Microban® product protection inhibits growth of bacteria & aids infection control
- Ergonomic handle & recessed shelves provide additional legroom allowing you to push the cart without hitting your shins & feet
- Chrome plated posts offers shelf adjustment in 25mm increments
- Available with two / three shelves in three sizes
- Small 620 x 415mm (internal) at heights 185, 515 & 845mm (lower shelves) adiustable)
  - · Easy grip handle & utility trav
  - Maximum capacity per shelf 68kg / per cart 136kg
  - 100mm pressed steel castors
- Medium 693 x 541mm (internal) at heights 185, 515 & 845mm (lower shelves) adjustable)
  - Handle has built in recessed support channel to manage odd shaped items
  - Maximum capacity per shelf 68kg / per cart 136kg
  - 100mm pressed steel castors
- Large 833 x 652mm (internal) at heights & 218, 548 & 878mm (lower shelves) adiustable)
  - Handle has built in slot to hold scanners & labellers
  - Maximum capacity per shelf 90kg / per cart 180kg
  - 125mm pressed steel castors
- Self-assembly (supplied boxed)

# **General Purpose Trolley Polymer, Deep Shelves**

#### Features

- Easy to clean corrosion resistant high-density polyethylene shelves (not suitable for use with automated washers)
- 70mm deep shelves retains spills & prevents items from sliding off during transport
- Microban® product protection inhibits growth of bacteria & aids infection control
- · Chrome plated posts offers shelf adjustment in 25mm increments
- Available with two / three shelves in two sizes
  - Small 730 x 469mm (internal) at heights 140, 470 & 800mm (lower shelves) adiustable)
  - Large 877 x 616mm (internal) at heights 140, 470 & 800mm (lower shelves adjustable)
- Max capacity per shelf 68kg / per cart 136kg
- 100mm pressed steel castors
- Self-assembly (supplied boxed)



#### Standard Colours

Standard Colours

Small - 850 x 547 x 1067mm

Large - 997 x 691 x 1067mm

· Uprights - Chrome Shelves - Blue

8BC2030-2DMB

8BC2030-3DMB

8BC2636-2DMB

8BC2636-3DMB

Item

- · Uprights Chrome
- · Shelves Blue

| ltem              | Shelves     | Weight |
|-------------------|-------------|--------|
| Small - 805 x 470 | x 902mm     |        |
| 8MY1627-24BU      | 2           | 12.9kg |
| 8MY1627-34BU      | 3           | 15.4kg |
| Medium - 877 x 6  | 600 x 902mi | m      |
| 8MY2030-24BU      | 2           | 14.9kg |
| 8MY2030-34BU      | 3           | 18.8kg |
| Large - 1023 x 70 | 04 x 937mm  | 1      |
| 8MY2636-25BU      | 2           | 18.1kg |
| 8MY2636-35BU      | 3           | 23.1kg |





sales@bristolmaid.com +44 (0)1258 484455 Shelves

2

3

2

3

www.bristolmaid.com

Weight

13.0ka

17.0kg

15.9kg

21.1kg

# General Purpose Trolley Mild Steel, Laminate Shelves

#### Features

- Laminate shelves
- Available with two or three shelves in two sizes
- Shelf heights
  - Two shelf 210 & 845mm
  - Three shelf 210, 525 & 845mm
- Shelf size
  - Small 705 x 445mm
  - Large 855 x 445mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel castors

#### **Standard Colours**

- Paint Grey
- Shelf Spearmint

#### Options

• ZP005 - Ring Buffers





| Item         | Description                                | Weight |
|--------------|--|--------|
| Small - 972  | x 545 x 1035mm                             |        |
| GP005        | Mild Steel - Two Laminate Shelves, Small   | 18.0kg |
| GP015        | Mild Steel - Three Laminate Shelves, Small | 23.0kg |
| Large - 1122 | 2 x 545 x 1035mm                           |        |
| GP025        | Mild Steel - Two Laminate Shelves, Large   | 20.0kg |
| GP035        | Mild Steel - Three Laminate Shelves, Large | 26.5kg |



# **HTM71 Storage Solutions**

- ScanCell<sup>®</sup> Carts 196
- Distribution Trolleys 201
- Open Transfer / Transport Trolleys 203
  - HTM71 Racking 204
- Trays, Baskets, Dividers & Accessories 206





### ScanCell<sup>®</sup> Carts - Open Fronts



#### **Features**

- · Mobile open front storage carts for use with 600 x 400mm trays & baskets
- Available as single, double & treble columns in a range of heights
- Lightweight easyclean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides. & back
- Fitted with angled / horizontal side panels to suit 400 or 600mm insertions
  - Type A 600mm Insertions
  - Type B 400mm Insertions
- 125mm polymer castors, 2 off braking

#### Standard Colours

Trays, Baskets, Dividers & Accessories Page 206 - 209

· Body - Grey

#### Options

- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
- 4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)
- 4PR1A Push Rail, Type A (Cart Heights 1503mm, 1803mm)
- 4PR1B Push Rail, Type B (Cart Heights 1503mm, 1803mm)



(trays & baskets not included)

(trays & baskets not included)

# ScanCell<sup>®</sup> Carts - Roller Shutters



#### **Features**

- Mobile storage carts fitted with roller shutters for use with 600 x 400mm trays & baskets
- · Available as single, double & treble columns in a range of heights
- Lightweight easyclean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back, polypropylene roller shutter
- · Fitted with angled / horizontal side panels to suit 400 or 600mm insertions
  - Type A 600mm Insertions
  - Type B 400mm Insertions
- 125mm polymer castors, 2 off braking

Trays, Baskets, Dividers & Accessories Page 206 - 209

# Body - Grey Options

Standard Colours

- 411P1A Litility Pail Type A (Cart Heights 005mm 1055
- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
   4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)









### **ScanCell® Carts - Doors**



Trays, Baskets, Dividers

& Accessories

Page 206 - 209

#### **Features**

- Mobile storage carts fitted with doors for use with 600 x 400mm trays & baskets
- Available as single, double & treble columns in a range of heights
- Lightweight easy clean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back, post form / plastic finish doors
- Fitted with angled / horizontal side panels to suit 400 or 600mm insertions
- Type A 600mm Insertion Single door for each column
  - Single Column R/H Hinged
  - Double Column R/H & L/H Hinged
  - Treble Column 1 x R/H & 2 x L/H Hinged
- Type B 400mm Insertion Double doors for each column
- 125mm polymer castors, 2 off braking

#### **Standard Colours**

Body - Grey

#### Options

- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
- 4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)
- 4PR1A Push Rail, Type A (Cart Heights 1503mm, 1803mm)
- 4PR1B Push Rail, Type B (Cart Heights 1503mm, 1803mm)



HTM71 Storage Solutions

sales@bristolmaid.com +44 (0)1258 484455

www.bristolmaid.com

# ScanCell® Carts - Vario Drawers



#### **Features**

- Mobile storage carts fitted with vario drawers & supplied with 600 x 400mm trays
- Available as single, double & treble columns in a range of heights
- Lightweight easy clean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back
- Type B 400mm Insertions
- 125mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

#### **Standard Colours**

- Single Column 662 x 494 x Height
- Body / Drawers Grey
  Drawer Handle Blue
- Double Column 1324 x 494 x Height
- Treble Column 1986 x 494 x Height

#### Options

- 4UR1B Utility Rail (Cart Height 905mm, 1055mm, 1205mm)
- 4PR1B Push Rail (Cart Height 1503mm)

| Item                    | 50mm Tray | 100mm Tray | 200mm Tray | 300mm Tray |
|-------------------------|-----------|------------|------------|------------|
| Height - 905mm          |           |            |            |            |
| 4B1VD09/A               | 1         | 2          | 1          |            |
| 4B1VD09/B               |           | 4          |            |            |
| 4B1VD09/C               | 2         | 3          |            |            |
| 4B2VD09/AA*             | 2         | 4          | 2          |            |
| 4B3VD09/AAA**           | 3         | 6          | 3          |            |
| Height - 1055mm         |           |            |            |            |
| 4B1VD10/D               | 1         | 3          | 1          |            |
| 4B1VD10/E               |           | 5          |            |            |
| 4B1VD10/F               | 2         | 4          |            |            |
| 4B2VD10/DD*             | 2         | 6          | 2          |            |
| 4B3VD10/DDD**           | 3         | 9          | 3          |            |
| Height - 1205mm         |           |            |            |            |
| 4B1VD12/G               | 1         | 4          | 1          |            |
| 4B1VD12/H               | 2         | 2          | 2          |            |
| Height - 1503mm         |           |            |            |            |
| 4B1VD15/I               | 1         | 6          | 1          |            |
| 4B1VD15/J               | 1         | 4          | 1          | 1          |
| *double column **treble | e column  |            |            |            |



**TYPE B - 400mm Insertion** 







#### sales@bristolmaid.com



# **ScanCell® Carts - Telescopic Drawers**

#### **Features**

- Mobile storage carts with telescopic drawers supplied with 600 x 400mm trays & baskets
- Available as single, double & treble columns in a range of heights
- Lightweight easy clean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back
- Type B 400mm Insertions
- Fitted with horizontal side panels to suit 400 insertions
- 125mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

- Single Column 662 x 494 x Height
- Double Column 1324 x 494 x Height

#### **Standard Colours**

• Body / Drawers - Grey

#### Options

- Trays, Baskets, Dividers & Accessories Page 206 - 209
- lone
- 4UR1B Utility Rail (Cart Height 905mm, 1055mm, 1205mm)

| Item            | 50mm Tray | 100mm Tray | 200mm Basket |
|-----------------|-----------|------------|--------------|
| Height - 905mm  |           |            |              |
| 4B1TD09/A       | 1         | 2          | 1            |
| 4B1TD09/B       |           | 4          |              |
| 4B1TD09/C       | 2         | 3          |              |
| 4B2TD09/AA *    | 2         | 4          | 2            |
| Height - 1055mm |           |            |              |
| 4B1TD10/D       | 1         | 3          | 1            |
| 4B1TD10/E       |           | 5          |              |
| 4B1TD10/F       | 2         | 4          |              |
| 4B1TD10/G       | 10        |            |              |
| Height - 1205mm |           |            |              |
| 4B1TD12/H       | 1         | 4          | 1            |
| 4B1TD12/I       | 2         | 5          |              |
| 4B1TD12/J       | 6         | 3          |              |

\*double column

# **ScanCell® Carts - Compact**



#### **Features**

- · Mobile storage carts with telescopic drawers supplied with 300 x 400mm trays & baskets
- Lightweight easyclean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back
- Type D 400mm Insertions
- 125mm polymer castors, 2 off braking

#### Dimensions (w x d x h) Options

- 462 x 494 x 905mm
- 4UR1B Utility Rail (Type B)

#### Standard Colours

· Body / Drawers - Grev



4D1TD09/A shown with optional utility rail

| Item      | 50mm Tray | 100mm Tray | 200mm Basket |
|-----------|-----------|------------|--------------|
| 4D1TD09/A | 1         | 2          | 1            |
| 4D1TD09/B |           | 4          |              |

# **Distribution Trolleys - Compact**

#### **Features**

- Mobile distribution trollevs for use with 300 x 400mm travs & baskets
- Available in mild & 304 grade stainless steel
  - · Mild Steel Welded & powder coated
  - · Stainless Steel Seamlessly welded for ease of cleaning
- Runners allow travs to be accessed from either side (4 stops required per trav / basket)
  - Type D 400mm insertion
- Trays & baskets can be supported on the upper & lower frame (optional solid grade laminate worktop)
- 75mm pressed steel castors

# Dimensions (w x d x h)

• 360 x 460 x 850mm

#### Standard Colours

· Paint - Grey

#### **Options**

- ZP010 Ring Buffers
- 4D1SLWT- Solid Grade Laminate Worktop







4D1DSTCM (trays & baskets not included)

| Item     | Description                             |
|----------|---|
| 4D1DSTCM | Compact, Single Column, Mild Steel      |
| 4D1DSTCS | Compact, Single Column, Stainless Steel |



# **Distribution Trolleys**

#### **Features**

- Mobile distribution trolleys for use with 600 x 400mm trays & baskets
- · Available in mild & 304 grade stainless steel
  - Mild Steel Welded & powder coated
  - · Stainless Steel Seamlessly welded for ease of cleaning
- Runners allow trays to be accessed from either side (4 stops required per tray / basket)
  - Type A 600mm insertion
  - Type B 400mm insertion
- Trays & baskets can be supported on the upper & lower frame (optional solid grade laminate worktop)
- 75mm pressed steel castors

#### **Standard Colours**

• Paint - Grey

#### Options

- ZP010 Ring Buffers
- 4A1SLWT Solid Grade Laminate Worktop, Single Column (Type A)
- 4A2SLWT Solid Grade Laminate Worktop, Double Column (Type A)
- 4B1SLWT Solid Grade Laminate Worktop, Single Column (Type B)

#### TYPE A - 600mm Insertion

#### **TYPE B - 400mm Insertion**

Trays, Baskets, Dividers

& Accessories

Page 206 - 209



| ltem            | Description                            | Overall Dimensions (w x d x h) |
|-----------------|--|--------------------------------|
| Mild Steel      |  |                                |
| 4A1DSTM         | Single Column, Type A, Mild Steel      | 460 x 655 x 850mm              |
| 4A2DSTM         | Double Column, Type A, Mild Steel      | 895 x 655 x 850mm              |
| 4B1DSTM         | Single Column, Type B, Mild Steel      | 660 x 455 x 850mm              |
| Stainless Steel |  |                                |
| 4A1DSTS         | Single Column, Type A, Stainless Steel | 460 x 655 x 850mm              |
| 4A2DSTS         | Double Column, Type A, Stainless Steel | 895 x 655 x 850mm              |
| 4B1DSTS         | Single Column, Type B, Stainless Steel | 660 x 455 x 850mm              |
|                 |  |                                |

sales@bristolmaid.com

665 x 725 x 1555mm

1100 x 725 x 1555mm

1535 x 725 x 1635mm

665 x 725 x 1555mm

1100 x 725 x 1555mm

1535 x 725 x 1635mm

Solutions

# **Open Transfer / Transport Trolleys**

#### Features

- Mobile transfer / transport trolleys for use with 600 x 400mm trays & baskets
- Available in mild & 304 grade stainless steel
  - Mild Steel Welded & powder coated
  - Stainless Steel Seamlessly welded for ease of cleaning
- Runners allow trays to be accessed from either side (4 stops required per tray / basket)
  - Type A 600mm insertion
- Push handles at either end
- · Ring buffers to protect the trolley / fabric of the building
- Castors
  - Single / Double Column 125mm pressed steel castors, 2 off braking
  - Treble Column 200mm heavy duty pressed steel castors, 2
     off fixed, 2 off swivel braking

#### **Standard Colours**

· Paint - Grey

#### Options

Item

Mild Steel 4A10TTM

4A2OTTM

4A3OTTM

4A2OTTS

4A3OTTS

Stainless Steel 4A10TTS

- 4A1OTTNC Nylon Cover (Single Column)
- 4A2OTTNC Nylon Cover (Double Column)

Description

Single Column, Type A, Mild Steel

Double Column, Type A, Mild Steel

Treble Column, Type A, Mild Steel

Single Column, Type A, Stainless Steel

Double Column, Type A, Stainless Steel

Treble Column, Type A, Stainless Steel

• 4A3OTTNC - Nylon Cover (Treble Column)













2HH-7532293 Wall Mounting Bracket



2HH-7532277 10mm Plastic Cover Cap

sales@bristolmaid.com

4HH-HTM71/S1/400 Runners, trays, baskets, dividers & angle brackets not included (order separately)

#### Components

- 2HH-7532301-HH Base Unit, Type A, 600mm Insertion (x 1)
- 2HH-7532311 Extension Unit, Type A, 600mm Insertion (x 1)
- 2HH-7532356 Runners, Type A, 600mm Insertion (pair)
- 2HH-7532329 Angle Adaptor Plate, Type A, 600mm Insertion (pair)
- 2HH-7532304-HH Base Unit, Type B, 400mm Insertion (x 1)
- 2HH-7532314 Extension Unit, Type B, 400mm Insertion (x 1)
- 2HH-7532355 Runners, Type B, 400mm Insertion (pair)
- 2HH-7532328 Angle Adaptor Plate, Type B, 400mm Insertion (pair)
- 2HH-7532320 Adjustable Foot Set (x 4)
- 2HH-7532321 Adjustable Foot Set (x 2)
- 2HH-7532322 50mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532323 50mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)
- 2HH-7532324 80mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532325 80mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)

# HTM71 Racking

#### • Fully configurable freestanding modular racking for use with 600 x 400mm modular trays & baskets

- Available to suit 600 or 400mm insertions
  - Type A 600mm insertion
  - Type B 400mm insertion
- · Trays & baskets supported on adjustable moulded plastic/glass fibre runners
- Runners can be positioned horizontal or angled (approx. 25 degrees requires angle adaptor brackets)
- · Integral stops on the runners allow trays/baskets to be accessible in two positions
  - Straight (approx. 5 degrees)
  - Angled (approx. 25 degrees)
- · Runners offer pass through facility allowing access on either side
- · Anodised extruded aluminium uprights & cross members
- Bolted construction using high quality moulded fixings provides rigidity & allows the system to be altered easily
- · Available as either static or mobile units
  - · Static Adjustable feet allow levelling on uneven floors
  - Mobile 50/80mm pressed steel castors
- Optional wall mounting brackets allow static units to be fixed securely to the wall, (when units are loaded over 100kg
- Supplied for self-assembly

#### Dimensions (w x d x h)

#### Static

Type A - 600mm Insertion

- Single Bay 494 x 642 x 1940mm
- Double Bay 937 x 642 x 1940mm

Type B - 400mm Insertion

- Single Bay 694 x 442 x 1940mm
- Double Bay 1337 x 442 x 1940mm

#### **Standard Colours**

- Uprights & Cross Members Anodised Aluminium
- Plastic Mouldings Green/Grey

#### **Static - Adjustable Feet**

Item Description Weight Type A - 600mm Insertion 4HH-HTM71/S1/600 Static, Single Bay, 600mm Insertion 10.3kg 4HH-HTM71/S2/600 Static, Double Bay, 600mm Insertion 16.1kg Type B - 400mm Insertion 4HH-HTM71/S1/400 Static, Single Bay, 400mm Insertion 10.2kg 4HH-HTM71/S2/400 Static, Double Bay, 400mm Insertion 16.3kg

- Mobile
  - Type A 600mm Insertion
    - Single Bay 494 x 642 x 1996mm
    - Double Bay 937 x 642 x 1996mm

#### Type B - 400mm Insertion

- Single Bay 694 x 442 x 1996mm
- Double Bay 1337 x 442 x 1996mm

#### Options

- 2HH-7532293 Wall Mounting Bracket
- 2HH-7532277 10mm Plastic Cover Cap

| +44 ( | (0)1258 484455 |
|-------|----------------|
|-------|----------------|



#### Mobile - 50mm Castors

| Item                     | Description                         | Weight |
|--------------------------|-------------------------------------|--------|
| Type A - 600mm Insertion |                                     |        |
| 4HH-HTM71/M51/600        | Mobile, Single Bay, 600mm Insertion | 11.1kg |
| 4HH-HTM71/M52/600        | Mobile, Double Bay, 600mm Insertion | 17.3kg |
| Type B - 400mm Insertion |                                     |        |
| 4HH-HTM71/M51/400        | Mobile, Single Bay, 400mm Insertion | 11.0kg |
| 4HH-HTM71/M52/400        | Mobile, Double Bay, 400mm Insertion | 17.5kg |



Runners, travs, baskets, dividers not included (order separately)

V LITTLE CONTRACTOR 

utation and the second

THE OWNER WATCHING THE PARTY OF THE PARTY OF

NUMBER OF BRIDE

manunununun

Runners, trays, baskets, dividers

not included (order separately)

TANK AND REAL PROPERTY AND REAL PROPERTY.

annin mark 

4HH-HTM71/M52/400

AND THURSDAY

......

111111

OTT DATA

TTLLIT

ITING

HUHB TZ

TETHON

....

1111111111111111111

## **Trays & Baskets** Type D - Size 300 x 400mm

Item







2HH-7014200

Basket - ABS, 200mm Dimensions (w x d x h): 300 x 400 x 200mm Material: ABS Grey

Description



2HH-7012100



2HH-7012200

Basket - ABS, 100mm Dimensions (w x d x h): 600 x 400 x 100mm Material: ABS Grey

Description

Tray - ABS, 50mm

Material: ABS Grev

Tray - ABS, 100mm

Material: ABS Grey

Basket - ABS, 200mm Dimensions (w x d x h): 600 x 400 x 200mm Material: ABS Grey



Tray/Lid - ABS, 15mm Dimensions (w x d x h): 600 x 400 x 15mm Material: Plastic Grey

2HH-K561039D


# **Dividers**

## Type D - Size 300 x 400mm

| Item         | Description  | Tray/Basket Compatibility  | Item        | Description  | Tray/Basket Compatibility  |
|--------------|--|----------------------------|-------------|--|----------------------------|
| 2HH-7810103  | Divider - PC<br>Dimensions (w x h): 300 x 50mm<br>Material: Polycarbonate, Transparent | 2HH-7014050                |             | Divider - ABS<br>Dimensions (w x h): 400 x 200mm<br>Material: ABS Grey   | 2HH-7014200<br>2HH-7012200 |
| 2HH-7810102  | Divider - PC<br>Dimensions (w x h): 400 x 50mm<br>Material: Polycarbonate, Transparent | 2HH-7014050<br>2HH-7011050 | 2HH-7811305 | Sliding Divider - SAN<br>Dimensions (w x h): 170 x 45mm<br>Material: Plastic, Transparent                            | 2HH-7014050<br>2HH-7011050 |
| 2HH-7811206N | Divider - ABS<br>Dimensions (w x h): 300 x 100mm<br>Material: ABS Grey                 | 2HH-7014100                | 2HH-241302  | Sliding Divider - SAN<br>(including cardholder)<br>Dimensions (w x h): 170 x 45mm<br>Material: Plastic, Transparent  | 2HH-7014050<br>2HH-7011050 |
| 2HH-7811205N | Divider - ABS<br>Dimensions (w x h): 400 x 100mm<br>Material: ABS Grey                 | 2HH-7014100<br>2HH-7012100 | 2HH-241301  | Sliding Divider - SAN<br>Dimensions (w x h): 170 x 100mm<br>Material: Plastic, Transparent                           | 2HH-7014100                |
|              | Divider - ABS<br>Dimensions (w x h): 300 x 200mm<br>Material: ABS Grey                 | 2HH-7014200                | 2HH-240820  | Sliding Divider - SAN<br>(including cardholder)<br>Dimensions (w x h): 170 x 100mm<br>Material: Plastic, Transparent | 2HH-7014100                |
| 2HH-7811306  |  |                            | 2HH-240821  |  |                            |

# **Dividers**


#### Type A & B - Size - 600 x 400mm



# Accessories

| Item         | Description   | Tray/Basket Compatibility |
|--------------|---|---------------------------|
| 2HH-240501   | T/Clip  | All Trays & Baskets       |
| 2HH-7910204N | Safety Stop - Long<br>(Scan Modul Compatible)<br>Material: Plastic  | All Trays & Baskets       |
| 2HH-7910201N | Safety Stop - Short<br>(H+H System Compatible)<br>Material: Plastic | All Trays & Baskets       |
| 2HH-100072   | Label Holder<br>Credit Card Size: (85 x 54mm)<br>Material: Plastic  | All Trays & Baskets       |
|              | Label<br>Credit Card Size: (85 x 54mm)<br>Material: Plastic, White  | All Trays & Baskets       |

2HH-K100187



# Infusion / Pump Stands & Accessories

- Mobile Pump Stands 212
- Mobile Infusion Stands 216
- Clamping Infusion Stands 221
- Wall Mounted Infusion Holders 221



## Mobile Pump Stands - Polyurethane, Easy Clean

#### **Features**

- Range of stands designed specifically to hold blood warmers, pumps & syringe drivers
- · Five leg base with low centre of gravity to assist stability
- 32mm diameter 304 grade stainless steel upright & 25mm diameter pump support bows
- Polyurethane overmoulded base with smooth edges & contours provides...
  - Water-resistant finish, allowing for easy & effective cleaning (including steam cleaning)
  - · Shield for the castors
  - Impact protection to patients & staff, reducing the risk of skin abrasions & bruising
  - Protection to the fabric of the building
- Height adjustable upright using a handwheel (1,300 2,000mm)
- · Four U shaped hooks
- Maximum recommended load...
  - Hooks 12kg per stand plus
  - Fixed Upright/s 22kg per stand
- 75mm polypropylene / thermoplastic rubber twin-wheel castors (1 x anti-static)

#### Dimensions (w x d x h)

- Overall 620 x 620 x 1,300 2,000mm
- Vertical Clamping Area
  - PS/P/SC/4H 975mm
  - PS/P/HC/4H 771mm (upright) 2 x 646mm (support bows)

#### **Standard Colours**

- Base Grey
- Upright & Support Bows Stainless Steel
- Vertical Clamping Area
- PS/P/SC/4H 975mm
- PS/P/HC/4H 771mm (upright) 2 x 646mm (support bows)

ItemDescriptionWeightPS/P/SC/4HMobile Pump Stand - Polyurethane, Easy Clean, Standard Capacity, Four Hook11.0kgPS/P/HC/4HMobile Pump Stand - Polyurethane, Easy Clean, High Capacity, Four Hook13.0kg

sales@bristolmaid.com +44 (0)1258 484455





# Mobile Pump Stands - Mild Steel, Four Leg

#### **Features**

- Range of stands designed specifically to hold blood warmers, pumps & syringe drivers
- · Fully welded four leg base
- 32mm diameter 304 grade stainless steel upright & 25mm diameter pump support bows
- Height adjustable upright using a handwheel (1,270 2,000mm)
- Four U-shaped hooks
- Maximum recommended load...
  - Hooks 12kg per stand plus
  - Fixed Upright/s 22kg per stand
- 75mm polypropylene / thermoplastic rubber twin-wheel castors (1 x anti-static)

#### Dimensions (w x d x h)

- Overall 560 x 621 x 1,270 2,000mm
- Vertical Clamping Area
  - PS/M/SC/4H 975mm
  - PS/M/HC/4H 771mm (upright) 2 x 646mm (support bows)

#### **Standard Colours**

- Base Grey
- Upright & Support Bows Stainless Steel

| Item       | Description  | Weight |
|------------|--|--------|
| PS/M/SC/4H | Mobile Pump Stand - Mild Steel, Four Leg, Standard Capacity, Four Hook | 14.7kg |
| PS/M/HC/4H | Mobile Pump Stand - Mild Steel, Four Leg, High Capacity, Four Hook     | 16.7kg |



www.bristolmaid.com



IS/P/2H/H/SP IS/P/2H/H/BL IS/P/2H/H/GR IS/P/2H/H/RD IS/P/2H/H/YE Two Hook, Handle

IS/P/4H/H/SP IS/P/4H/H/BL IS/P/4H/H/GR IS/P/4H/H/RD IS/P/4H/H/YE Four Hook, Handle

sales@bristolmaid.com

#### Mobile Infusion Stands - Polyurethane, Easy Clean

#### **Features**

- Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile Infusion Stands
- · Five leg base with low centre of gravity to assist stability
- · 25mm diameter 304 grade stainless steel upright
- Polyurethane overmoulded base with smooth edges & contours provides...
  - Water-resistant finish, allowing for easy & effective cleaning (including steam cleaning)
  - Shield for the castors
  - Impact protection to patients & staff, reducing the risk of skin abrasions & bruising
  - Protection to the fabric of the building
- Coloured indicator caps provide department identification & segregation, reducing the risk of cross contamination
- Height adjustable upright using a handwheel (1,260 2,100mm)
- Available with two or four pigtail hooks
- Maximum recommended load...
  - Hooks 4kg per stand plus
  - Fixed Upright 4kg per stand
- Optional handle keeps feet clear of the base, assists manoeuvrability
   & improves patient safety
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors (1 x anti-static)

#### Dimensions (w x d x h)

- Overall 620 x 620 x 1,260 2,100mm
- Vertical Clamping Area
  - IS/P/2H & IS/P/4H 970mm
  - IS/P/2H/H & IS/P/4H/H 840mm

#### **Standard Colours**

- · Base Grey
- · Upright & Handle Stainless Steel
- · Coloured Indicator Cap Spearmint

#### Options

• Coloured Indicator Cap - Blue / Green / Red / Yellow



8.2kg

8.3kg





#### **Mobile Infusion Stands - Stainless Steel**

#### **Features**

- Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile
   Infusion Stands
- · Seamless welded 304 grade stainless steel construction
- Height adjustable upright using a handwheel (1,270 2,000mm)
- Available with two or four pigtail hooks
- 9kg plastic coated weighted base available for additional stability
- Maximum recommended load...
- · Hooks 4kg per stand plus
- Fixed Upright
  - Non-Weighted Base 2kg per stand
  - Weighted Base 4kg per stand
- 50mm polypropylene / thermoplastic rubber twin wheel castors (1 x anti-static)

#### Dimensions (w x d x h)

• 650 x 650 x 1,270 - 2,000mm

#### **Standard Colours**


- Base / Upright Stainless Steel
- Weighted Base Grey



sales@bristolmaid.com



www.bristolmaid.com



| Item      | Description                                  | Weight |
|-----------|--|--------|
| DPKS/1    | Mobile Infusion Stand - Knockdown, Two Hook  | 3.8kg  |
| DPKS/1/4H | Mobile Infusion Stand - Knockdown, Four Hook | 3.9kg  |

sales@bristolmaid.com

| Item  | Description             | Weight |
|-------|-------------------------|--------|
| DP/2H | Two Hook Infusion Pole  | 0.6kg  |
| DP/4H | Four Hook Infusion Pole | 0.8kg  |

+44 (0)1258 484455

www.bristolmaid.com

#### **Clamping Infusion Stands - Trolley**

#### **Features**

- Clamp secures to 25mm vertical round tube (designed to fit our range of dressing trolleys)
- When fitting to other trolleys, we strongly recommend that a risk assessment is carried out to assess stability
- When fitted protrudes 45mm
- Supplied with a two pigtail hook infusion pole
- Maximum load on hooks 2kg

#### Dimensions (w x d x h)

• 55 x 220 x 1200mm





| Item | Description                        | Weight |
|------|------------------------------------|--------|
| DP/C | Clamping Infusion Stands - Trolley | 1.8kg  |

#### Wall Mounted Infusion Holders

#### Features

CE

- Wall mounted infusion bag holder
- Available with one/two U shaped hooks
- 25mm diameter "D" shaped tube to accept pumps (max load 10kg)
- Maximum load per hook 2kg
- Fixings not included

#### Dimensions (w x d x h)

- IV001 344 x 75 x 500mm
- IV002 344 x 75 x 570mm



| Item  | Description                            | Weight |
|-------|--|--------|
| IV001 | One Hook, Wall Mounted Infusion Holder | 1.2kg  |
| IV002 | Two Hook, Wall Mounted Infusion Holder | 1.6kg  |

Infusion &

€



## **Linen Trolleys & Bags**

- Soiled Linen Trolleys 224
  - Linen Trolleys 227
  - Hamper Trolleys 228
  - Soiled Linen Bags 229
- Clean Linen Trolleys 230

#### Soiled Linen Trolleys - Stainless Steel, Cantilever Frame



#### **Features**

- · Seamlessly welded from from 304 grade stainless steel for strength & ease of cleaning
- Suitable for use with fabric linen bags using a tubular C shaped ring\*
- · Bags supported on an easy clean solid base
- · Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors

\*Recommended minimum bag circumference 1360mm, height 750mm

#### Dimensions (w x d x h)

- Single 440 x 460 x 820mm
- Double 835 x 460 x 820mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

#### **Options / Accessories**

Linen Bags:

- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grev
- R/BAG/ORANGE Orange

Linen Bags page 229



Pedal operated lid

WT/6/2/PL/SS

Weight

6.5kg

10.0kg

8.5kg

```
www.bristolmaid.com
```

WT/6/PL/SS

Description

Stainless Steel Cantilever, Single

Stainless Steel Cantilever, Double

Stainless Steel Cantilever, Single, Pedal Operated

Item

WT/6/SS

WT/6/2/SS

WT/6/PL/SS

#### Soiled Linen Trolleys - Mild Steel, Cantilever Frame

#### Features

- · Welded mild steel frame
- · Suitable for use with fabric linen bags using a tubular C shaped ring\*
- Bags supported on an easy clean solid base
- Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors

\*Recommended minimum bag circumference 1360mm, height 750mm

#### Dimensions (w x d x h)

- Single 440 x 460 x 820mm
- Double 835 x 460 x 820mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

#### **Standard Colours**

· Paint - Spearmint

#### **Options / Accessories**

Linen Bags:

- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grey
- R/BAG/ORANGE Orange

Linen Bags page 229



Mild Steel Cantilever, Single, Pedal Operated

Mild Steel Cantilever, Double, Pedal Operated

WT/6/PL/MS

WT/6/2/PL/MS

www.bristolmaid.com

8.5ka

14.5kg



#### Soiled Linen Trolleys - Mild Steel, Sprung Frame

#### **Features**

- Welded mild steel frame
- Suitable for use with either fabric or disposable polythene bags using a sprung rod frame\*
- · Bags supported on an easy clean solid base
- Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors

\*Factory set for standard 1360mm circumference bags. Can be manually adjusted to suit other bags. Recommended minimum bag height 750mm

#### Dimensions (w x d x h)

- Single 440 x 460 x 825mm
- Double 835 x 460 x 825mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

#### **Standard Colours**

Paint - Spearmint

#### **Options / Accessories**

- Linen Bags:
  - R/BAG/WHITE White
  - R/BAG/BLUE Blue
  - R/BAG/RED Red
  - R/BAG/GREEN Green
  - R/BAG/YELLOW Yellow
  - R/BAG/GREY Grey
  - R/BAG/ORANGE Orange

Linen Bags page 229

www.bristolmaid.com

10.0kg

14.5kg

WT/9/2/MS

Mild Steel, Sprung Frame, Double

WT/9/2/PL/MS Mild Steel, Sprung Frame, Double, Pedal Operated

#### Soiled Linen Trolley - Mild Steel, Folding

#### **Features**

- · Folds to minimize space
- Suitable for use with fabric linen bags
- 50mm polypropylene / thermoplastic rubber twin-wheel castors
- \*Recommended minimum bag circumference 1360mm, height 750mm

#### Dimensions (w x d x h)

- Overall Dimensions 430 x 430 x 1030mm
- Folded 430 x 200 x 1040mm

#### **Standard Colours**

· Paint - Grey

#### **Options / Accessories**

- Linen Bags:
  - R/BAG/WHITE White
  - R/BAG/BLUE Blue
  - R/BAG/RED Red
  - R/BAG/GREEN Green
  - R/BAG/YELLOW Yellow
  - R/BAG/GREY Grey
  - R/BAG/ORANGE Orange

#### Linen Bags page 229

#### Linen Trolley - Bed Changing

#### **Features**

- Welded mild steel frame
- Three removable, sunken, liquid retaining trays
- Tray heights 230, 590 & 960mm
- Dividing panel separates clean & soiled linen
- Suitable for use with fabric linen bags using a tubular C shaped ring\*
- Bag supported on a base frame
- 100mm pressed steel castors

\*Recommended minimum bag circumference 1360mm, height 750mm

#### Dimensions (w x d x h)

- Overall Dimensions 830 x 490 x 1040mm
- Tray 390 x 400 x 25mm

#### **Standard Colours**

· Paint - Grey

#### **Options / Accessories**

- ZP005 100mm Ring Buffers
- Linen Bags:
  - R/BAG/WHITE White
  - R/BAG/BLUE Blue
  - R/BAG/RED Red
  - R/BAG/GREEN Green
  - R/BAG/YELLOW Yellow
  - R/BAG/GREY Grey
  - R/BAG/ORANGE Orange



| Item    | Description         | Weight |
|---------|---------------------|--------|
| WT/7/MS | Folding, Mild Steel | 3.25kg |

| Item     | Description  | Weight |
|----------|--------------|--------|
| WT/13/MS | Bed Changing | 15.5kg |

WT/7/MS



#### Hamper Trolleys - Mild Steel

#### Features

- Welded mild steel construction
- Suitable for use with fabric hamper bags supported on 4 hooks
- Bag supported on an easy clean solid base
- 100mm pressed steel castors

#### Dimensions (w x d x h)

• Overall Dimensions - 815 x 514 x 1000mm

#### Standard Colours

Paint - Spearmint

#### **Options / Accessories**

- ZP005 100mm Ring Buffers 2FACO013 Hamper Bag (190 litre)

| Item     | Description                 | Weight |
|----------|-----------------------------|--------|
| WT/11/MS | Hamper Trolley - Mild Steel | 12.5kg |

www.bristolmaid.com

#### **Soiled Linen Bags**

#### Features

- Washable 7oz polyester soiled linen bag
  Circumference 1360mm, height 750mm, 120 litres
- Fitted with eyelets and draw cord with plastic toggle
  Available in a range of colours to allow compliance with local, national & international guidelines





| Item         | Description      | Weight |
|--------------|------------------|--------|
| R/BAG/WHITE  | White Linen Bag  | 0.3kg  |
| R/BAG/BLUE   | Blue Linen Bag   | 0.3kg  |
| R/BAG/RED    | Red Linen Bag    | 0.3kg  |
| R/BAG/GREEN  | Green Linen Bag  | 0.3kg  |
| R/BAG/YELLOW | Yellow Linen Bag | 0.3kg  |
| R/BAG/GREY   | Grey Linen Bag   | 0.3kg  |
| R/BAG/ORANGE | Orange Linen Bag | 0.3kg  |



R/BAG/BLUE



R/BAG/YELLOW



R/BAG/GREY



R/BAG/RED



R/BAG/GREEN



R/BAG/ORANGE



#### Clean Linen Trolleys - Small (390 Litre)

#### Features

- Suitable for ward clean linen storage
- Mild steel frame
- Lightweight aluminium cladding & shelves
- Powder coated for durability
- Two adjustable shelves plus base (max 45kg per shelf)
- Shelf size 1000 x 450mm (w x d)
- · Push handles at either end
- Available with a flexible cover or fold back doors, preventing dust ingress
- 125mm pressed steel castors
- Maximum capacity 135kg per trolley

#### Dimensions (w x d x h)

- Overall Dimensions 1175 x 530 x 1065mm
- Internal Dimensions 1000 x 500 x 810mm

#### **Standard Colours**

• Paint - Grey

#### Options

- C008 125mm Castors, Heavy Duty Pressed Steel, 2 x Braking, 2 x Non-Braking (set of 4) Factory fitted
- ZP010 100mm Ring Buffers

| Item        | Description                                  | Weight |
|-------------|--|--------|
| WT/42       | Clean Linen Trolley - Small                  | 39.7kg |
| WT/42/COVER | Clean Linen Trolley - Small, Flexible Cover  | 42.4kg |
| WT/42/FBD   | Clean Linen Trolley - Small, Fold Back Doors | 46.0kg |



#### **Flexible Cover**

- · Fire retardant cover
- Secured to the top of the trolley & drops down to enclose the front
- Weighted bottom hem with velcro style side fastenings



#### **Fold Back Doors**

- Fold back doors encloses the front
- When opened, doors fold against the push handles at either end of the trolley
- Drop down bolt secures base of doors
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)

# -inen Trolleys & Bags

#### AGV Compatible Trolleys & Carts Bristol Maid<sup>™</sup> is your ideal partner for transport trolleys & carts for Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy

distribution



WT/50

sales@bristolmaid.com

| Item             | Description  | Weight  |
|------------------|--|---------|
| Non-Towing       |  |         |
| WT/50            | Clean Linen Trolley - Medium                           | 87.0kg  |
| WT/50/COVER      | Clean Linen Trolley - Medium, Flexible Cover           | 90.0kg  |
| WT/50/FBD        | Clean Linen Trolley - Medium, Fold Back Doors          | 97.3kg  |
| Towing           |  |         |
| WT/50/TBAR       | Clean Linen Trolley - Medium, Tow Bar                  | 89.5kg  |
| WT/50/COVER/TBAR | Clean Linen Trolley - Medium, Flexible Cover, Tow Bar  | 92.5kg  |
| WT/50/FBD/TBAR   | Clean Linen Trolley - Medium, Fold Back Doors, Tow Bar | 100.8kg |

#### Clean Linen Trolleys - Medium (690 Litre)

#### **Features**

- Widely used for bulk linen delivery, part of a linen exchange system, or as an alternative to dedicated ward linen rooms
- Mild steel frame
- Lightweight aluminium cladding & shelves
- Powder coated for durability
- Three adjustable shelves plus base (max 45kg per shelf)
- Shelf size 1180 x 460mm (w x d)
- · Inset push handles at either end
- Low level bumper bar offers protection to both the trolley & fabric of the building
- Available with a flexible cover or fold back doors, preventing dust ingress
- Optional tow bar fitting
- 200mm heavy duty castors, 2 x fixed, 2 x swivel
- Maximum capacity 180kg per trolley

#### Dimensions (w x d x h)

- Overall Dimensions 1430 x 600 x 1495 mm
- Internal Dimensions 1180 x 500 x 1180 mm

#### **Standard Colours**

• Paint - Grey



#### **Tow Bar Fitting**

- Tow bar fitted to the front (swivel castor)
- Antiluce fastener fitted to rear (fixed castor)
- Tow bar position approx. 360 380mm above floor
- Recommended max speed 7kph for towed trolleys



#### **Flexible Cover**

- · Fire retardant cover
- Secured to the top of the trolley & drops down to enclose the front
- Weighted bottom hem with velcro style side fastenings



WT/50/FBD

#### **Fold Back Doors**

- · Fold back doors, encloses the front
- When opened, doors fold back flat at either end of the trolley\*
- Drop down bolt secures base of doors
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)

\*Doors will not fold back flat when a tow bar is fitted

#### AGV Compatible Trolleys & Carts

For Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy distribution



| Item             | Description   | Weight  |
|------------------|---|---------|
| Non-Towing       |   |         |
| WT/48            | Clean Linen Trolley - Large                           | 99.6kg  |
| WT/48/COVER      | Clean Linen Trolley - Large, Flexible Cover           | 102.9kg |
| WT/48/FBD        | Clean Linen Trolley - Large, Fold Back Doors          | 111.5kg |
| Towing           |   |         |
| WT/48/TBAR       | Clean Linen Trolley - Large, Tow Bar                  | 103.1kg |
| WT/48/COVER/TBAR | Clean Linen Trolley - Large, Flexible Cover, Tow Bar  | 106.4kg |
| WT/48/FBD/TBAR   | Clean Linen Trolley - Large, Fold Back Doors, Tow Bar | 115.0kg |

#### Clean Linen Trolleys - Large (1,100 Litre)

#### **Features**

- Widely used for bulk linen delivery, part of a linen exchange system, or as an alternative to dedicated ward linen rooms
- Mild steel frame
- Lightweight aluminium cladding & shelves
- Powder coated for durability
- Three adjustable shelves plus base (max 45kg per shelf)
- Shelf size 1430 x 555mm (w x d)
- · Inset push handles at either end
- Low level bumper bar offers protection to both the trolley & fabric of the building
- Available with a flexible cover or fold back doors, preventing dust ingress
- Optional tow bar fitting
- 200mm heavy duty castors, 2 x fixed, 2 x swivel
- Maximum capacity 180kg per trolley

#### Dimensions (w x d x h)

- Overall Dimensions 1710 x 700 x1620mm
- Internal Dimensions 1420 x 610 x 1295mm

#### **Standard Colours**

· Paint - Grey



#### **Tow Bar Fitting**

- Tow bar fitted to the front (swivel castor)
- Antiluce fastener fitted to rear (fixed castor)
- Tow bar position approx. 360 380mm above floor
- Recommended max speed 7kph for towed trolleys

www.bristolmaid.com



# 0 e 6 o

#### **Flexible Cover**

- · Fire retardant cover
- Secured to the top of the trolley & drops down to enclose the front
- Weighted bottom hem with velcro style side fastenings



#### **Fold Back Doors**

- Fold back doors, encloses the front
- When opened, doors fold back flat at either end of the trolley\*
- Drop down bolt secures base of doors
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)

\*Doors will not fold back flat when a tow bar is fitted



# **Maternity Equipment**

- Baby Cribs Variable Height, Battery Operated 238
  - Baby Cribs Fixed Height 240
- Baby Crib Variable Height, Gas Strut Operation 242



#### **Baby Crib - Variable Height, Battery Operation**

#### **Features**

- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Rechargeable battery powered height adjustment (700 910mm)
- Crib can be tilted at either end (horizontal to approx. 6 degrees)
- Push handles at either end
- Lower cupboard with fixed shelf
- Two lift up work flaps (max capacity 2.0kg)
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- Supplied with one rechargeable battery & handset\*
- 100mm pressed steel braking castors

\*Battery must be removed from the unit for charging (see options / accessories for additional batteries)

#### Dimensions (w x d x h)

- Overall Dimensions 940 x 510 x 700 910mm
- Mattress 600 x 330 x 60mm
- Cupboard 210 x 265 x 315mm (internal, either side of the column)
- Work Flap 415 x 305mm

#### **Standard Colours**

- · Paint Spearmint
- Cabinet / Lower Shelf / Work Flaps Spearmint / Grey

#### **Accessories / Spares**

- 2LIDE033 Additional Rechargeable Battery
- 2LIDE037 Charger for Rechargable Battery
- 2PLCO014 Clear Plastic Crib, Solid Side (Spare / Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare / Replacement)

| Item     | Description                       | Weight |
|----------|-----------------------------------|--------|
| BA/VH    | Battery Operated, With Charger    | 49.5kg |
| BA/VH/NC | Battery Operated, Without Charger | 49.5kg |

```
184455 v
```

```
www.bristolmaid.com
```





Suppled with a handset



Rechargeable battery located at rear of crib



#### **Baby Cribs - Fixed Height**

#### Features

- Available with...
  - Lower shelf
  - · Lower shelf & small cupboard (with adjustable shelf)
  - Large cupboard (with adjustable shelf)
- Crib can be tilted at either end (horizontal to approx. 6 degrees)
- Two lift up work flaps (max capacity 2.0kg)
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions 952 x 500 x 929mm
- Mattress 600 x 330 x 60mm
- Large Shelf 810 x 400mm at height 270mm (BA/3/B)
- Small Shelf 460 x 355mm at height 270mm (BA/6)
- Small Cupboard 355 x 420 x 345mm (internal) at height 270mm (BA/6)
- Large Cupboard 710 x 420 x 345mm (internal) at height 270mm (BA/6/A)
- Work Flap 415 x 305mm

#### **Standard Colours**

- · Paint Grey White
- Cabinet / Lower Shelf / Work Flaps Spearmint / Grey White

#### **Accessories / Spares**

- 2PLCO014 Clear Plastic Crib, Solid Side (Spare / Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare / Replacement)

| Item   | Description                  | Weight |
|--------|------------------------------|--------|
| BA/3/B | Fixed Height, Lower Shelf    | 23.5kg |
| BA/6   | Fixed Height, Small Cupboard | 34.0kg |
| BA/6/A | Fixed Height, Large Cupboard | 42.0kg |





#### **Baby Crib - Variable Height, Gas Strut Operation**

#### **Features**

- Gas damper powered height adjustment (710 1025mm)
- Crib can be tilted at either end (horizontal to approx. 6 degrees)
- · Push handles at either end
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- 75mm polymer twin-wheel braking castors

#### **Dimensions** (w x d x h)

- Overall Dimensions 930 x 476 x 775 1025mm
- Mattress 600 x 330 x 60mm

#### **Standard Colours**

- · Paint Spearmint / Grey
- · Work Flap Spearmint

#### **Accessories / Spares**

- BA902 Slide Out Debris Bin
- BA903 Wire Basket Internal Dimensions 420 x 170 x 220mm
- BA904 Lift Up Flap Dimensions 400 x 300mm
- 2PLCO014 Clear Plastic Crib (Spare / Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare / Replacement)





#### Other products that you may be interested in...



#### **Overnight Chair - Portland, Three Position**





### **Medical Records**

- Medical Records Trolleys Locking 246
- Medical Records Workstations Locking 248
  - Laptop Holders 250
  - Carefile Filing Pockets 252
  - Medical Records Trolleys Filing 254
  - Medical Records Transfer Trolleys 256





#### **Medical Records Trolleys - Locking**

#### **Features**

- Suitable for use with our popular narrow filing pockets (order separately)
- Filing capacity
  - Small 470mm
  - Medium 695mm
  - Large 940mm
- Locking lid/s act as a work surface when open & is fitted with either a...
  - CAM Lock Supplied with 2 keys
  - Electronic Push Button Lock Programmable four-digit code
- Maximum load on locking lid/s 2.0kg
- Document pocket/s fitted to end of trolley
- Optional ring buffers offer protection to both the trolley & fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

#### Dimensions (w x d x h)

- Small 671 x 522 x 1056mm (width open 1190mm)
- Medium 920 x 522 x 1056mm (width open 1420mm)
- Large 1240 x 522 x 1056mm (width open 2085mm)

#### **Standard Colours**

- Paint Spearmint / Grey
- Lid / Flap Spearmint

#### **Options / Accessories**

- CF/F/PK10 Foolscap Filing Pockets (x 10)
- R/C/T Paper Tabs (x 500)
- R/C/TH Tab Holders (x 25)
- ZP010 Ring Buffers Factory Fitted Only
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C005 100mm Polymer Castor, 4 x Braking (set of 4) Factory Fitted Only

| High Security Bolt Lock | Electronic Push Button Lock | Description | Weight |
|-------------------------|-----------------------------|-------------|--------|
| MR210                   | MR210CL                     | Small       | 27.2kg |
| MR310                   | MR310CL                     | Medium      | 34.4kg |
| MR410                   | MR410CL                     | Large       | 46.4kg |



02

56 78

00

\*\*\*\*\*\*\*\*\*\*\*\*

#### www.bristolmaid.com



#### **Electronic Push Button Lock**

Uses a 10-button keypad

11-12

- Easy programming & code change
- Four operating codes with differing levels of authority -Master / Sub-Master / User / Technician

MRW/S1DXP/CAM

- User codes are 4 digits
- · Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

#### **Medical Records Workstations - Locking**

#### Features

- Suitable for use with our filing pockets (order separately)
- Provides 375mm filing capacity per drawer
- Available as single & double column units with...
  - One / two drawers & x-ray pigeon hole (5 / 10 sections)
  - One / two drawers & open cupboard
  - Two / four drawers
- Two locking options...
  - · CAM Lock Supplied with two keys
- · Electronic Push Button Lock Programmable four-digit code
- Where relevant one lock secures all drawers
- Double drawer units fitted with an anti-tilt mechanism preventing two drawers from being opened at the same time
- Max load per drawer 16.0kg
- Push handle & document pocket fitted to end of trolley
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Single Column 604 x 538 x 1028mm
- Double Column 1056 x 538 x 1034mm

#### **Standard Colours**

- Paint Grey White
- Cabinet Spearmint / Grey White

MRW/S1DOC/CAM

#### Single Column

MRW/S2D/CAM

fitted with MRW/LHR/S

(laptop not included)

| High Security Bolt Lock | Electronic Push Button Lock | Description                   | Weight |
|-------------------------|-----------------------------|-------------------------------|--------|
| MRW/S1DXP/CAM           | MRW/S1DXP/CL                | One Drawer, X-Ray Pigeon Hole | 34.2kg |
| MRW/S1DOC/CAM           | MRW/S1DOC/CL                | One Drawer, Open Cupboard     | 34.0kg |
| MRW/S2D/CAM             | MRW/S2D/CL                  | Two Drawers                   | 40.6kg |



#### **Options / Accessories**

MRW/D2DXP/CAM

- MRW/LHR/S Laptop Holder, Removable, Swing Arm
   MRW/LHNR/S Laptop Holder, Non-Removable, Swing Arm
- CF/F/PK10 Foolscap Filing Pockets (x 10)
- R/C/T Paper Tabs (x 500)
- R/C/TH Tab Holders (x 25)
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
   C005 100mm Polymer Castor, 4 x Braking (set of 4) Factory Fitted Only

#### **Double Column**

| High Security Bolt Lock | Electronic Push Button Lock | Description                    | Weight |
|-------------------------|-----------------------------|--------------------------------|--------|
| MRW/D2DXP/CAM           | MRW/D2DXP/CL                | Two Drawers, X-Ray Pigeon Hole | 62.0kg |
| MRW/D2DOC/CAM           | MRW/D2DOC/CL                | Two Drawers, Open Cupboard     | 58.0kg |
| MRW/D4D/CAM             | MRW/D4D/CL                  | Four Drawers                   | 72.0kg |

sales@bristolmaid.com +44 (0)1258 484455

MRW/D2DOC/CL

MRW/D4D/CAM

-2

#### **Laptop Holders**



#### Removeable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MRW range of medical record workstations
- Fitted to an upright bolted to the rear left-hand corner of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle, fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

\*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**

• Paint - Grey White



| Item      | Description          | Weight |
|-----------|----------------------|--------|
| MRW/LHR/S | Removable, Swing Arm | 3.7kg  |



## Non-Removeable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MRW range of medical record workstations
- Fitted to an upright bolted to the rear left-hand corner of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle fitted with a steel mounting bracket, fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

\*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**

• Paint - Grey White



| Item       | Description              | Weight |
|------------|--------------------------|--------|
| MRW/LHNR/S | Non-Removable, Swing Arm | 5.0kg  |




## **Carefile Filing Pockets**

#### **Features**

- Filing pockets suitable with Bristol Maid trolleys / workstations & most standard filing suspension systems
- Provides a durable alternative to paper or card types
- Suitable for runners with 380mm spacings / centres
- Provides an internal pocket (w x h) 360 x 260mm
- Fire retardant PVC coated nylon pockets
- Unique 'stay-put' tab holders, with pink / blue reversible card tabs
- Supplied in packs of ten

#### **Standard Colours**

• Nylon Pockets - Grey



| Item      | Description                    |
|-----------|--------------------------------|
| CF/F/PK10 | Foolscap Filing Pockets (x 10) |
| R/C/T     | Paper Tabs (x 500)             |
| R/C/TH    | Tab Holders (x 25)             |

## Other products that you may be interested in...

See Page 339

#### Medium & High Back Patient Armchairs - Alderbury

5ND/MAR/ALDE/VVDEL





## **Medical Records Filing Trolleys - Tiered**

#### **Features**

- Suitable for use with ring binders, case note & lever arch files (max file height 350mm)
- Available in two sizes with either an open front or locking door/s
- Small 900mm filing capacity
- Large 1800mm filing capacity Adjustable dividers (small 4 / large 8)
- Doors fold against the sides of the trolley (not flat)
- Large trolleys use one lock to secure both doors
- Two locking options...
  - · CAM Lock Supplied with two keys
  - Electronic Push Button Lock Programmable four-digit code
- Optional ring buffers offer protection to both the trolley & fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

#### Dimensions (w x d x h)

- Small 546 x 440 x 1016mm
- Large 998 x 440 x 1016mm

#### **Standard Colours**

- · Paint Spearmint / Grey
- Top Spearmint

| Item         | Description  | Weight |
|--------------|--|--------|
| Small        |  |        |
| MRF/S/OC     | Small - Open Cupboard                              | 21.5kg |
| MRF/S/LD/CAM | Small - Locking Door, CAM Lock                     | 26.5kg |
| MRF/S/LD/CL  | Small - Locking Door, Electronic Push Button Lock  | 26.5kg |
| Large        |  |        |
| MRF/L/OC     | Large - Open Cupboard                              | 30.0kg |
| MRF/L/LD/CAM | Large - Locking Doors, CAM Lock                    | 41.0kg |
| MRF/L/LD/CL  | Large - Locking Doors, Electronic Push Button Lock | 41.0kg |
|              |  |        |

- MRF/LHR/S Laptop Holder Removable, Swing Arm
- MRF/LHNR/S Laptop Holder Non-Removable, Swing Arm
- ZP010 Ring Buffers Factory Fitted Only
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C005 100mm Polymer Castor, 4 x Braking (set of 4) Factory Fitted Only
- C006 100mm Polymer Castor, 4 x Non-Braking (set of 4) Factory Fitted Only

\$

• Low battery warning with failure override using 9v battery

across LED's (+/-) & master code



#### Capacity

|                | Case Note Files  | Ring Binders     | Lever Arch Files |  |
|----------------|------------------|------------------|------------------|--|
|                | 235 x 50 x 320mm | 250 x 40 x 315mm | 290 x 75 x 315mm |  |
| Small - 900mm  | 18               | 22               | 12               |  |
| Large - 1800mm | 36               | 44               | 24               |  |

## Medical Records Transfer Trolleys Single Box



#### **Features**

- Suitable for use the transportation of bulk case notes & x-rays
- Records stored in an upper box
- Lower tray
- Two locking options...
  - CAM lock supplied with 2 keys
  - Electronic Push Button lock Programmable four-digit code
- Trolley max capacity 50kg
- Optional ring buffers offer protection to both the trolley & fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

#### Dimensions (w x d x h)

- Overall Dimensions 817 x 523 x 975mm
- Box (internal) 600 x 450 x 250mm
- Top of Box (to floor) 805mm

#### **Standard Colours**

• Paint - Spearmint / Grey

#### **Options**

• ZP005 - Ring Buffers

| Item    | m Description                            |        |  |  |
|---------|--|--------|--|--|
| Open    |  |        |  |  |
| MR130   | Open Top                                 | 16.0kg |  |  |
| MR140   | Locking Lid, CAM Lock                    | 19.0kg |  |  |
| MR140CL | Locking Lid, Electronic Push Button Lock | 19.0kg |  |  |

## Medical Records Transfer Trolley Double Box, Extending

#### **Features**

- Trolley for the transportation of bulk case notes & x-rays
- · Records stored in two removable boxes
- Clearance between upper & lower boxes, 235mm
- Max capacity 25kg (per box)
- · 80mm ring buffers fitted to front of unit
- 125mm pressed steel castors & 200mm wheels (braking)

#### Dimensions (w x d x h)

- Overall Dimensions 875 x 598 x 1040mm
- Box (internal) 568 x 368 x 316mm
- Top of Box (to floor)
  - Lower 440mm
  - Upper 975mm

#### **Standard Colours**

· Paint - Spearmint

Alternative

Colours

Page 450

Description

Double Box, Extending

Item

MR110

#### **Options / Accessories**

• 2PLCO271 - Removable Box

## Medical Records Transfer Trolley Double Box, Duo

#### **Features**

- · Unique 'two' trolley system for the bulk transportation of case notes & x-rays
- Records stored in fixed boxes, each supplied with an adjustable divider
- · Clearance between upper & lower boxes, 87mm
- Trolleys are loaded independently, before being linked together & transported
- Trolleys link together using a simple clevis pin
- Max capacity 25kg (per box)
- 100mm pressed steel castors (2 off braking) & 200mm wheels (braking)

#### Dimensions (w x d x h)

- Overall Dimensions 866 x 522 x 970mm
- Box (internal) 600 x 420 x 255mm
- Top of Box (to floor)
- Lower 527mm
- Upper 875mm

#### **Standard Colours**

Paint - Spearmint / Grey



www.bristolmaid.com

| em   | Description     | Weight |
|------|-----------------|--------|
| R120 | Double Box, Duo | 34.6kg |

Weight

31.6kg

Ite

MF



# Phlebotomy

- Caretray Trolley Phlebotomy 260
- Phlebotomy Chairs Variable Height 262
- Phlebotomy Chairs Variable Height, Bariatric 263
  - Phlebotomy Chairs Fixed Height 264



## **Caretray Trolley - Phlebotomy**

#### **Features**

- Trollevs designed for use with Phlebotomy
- Quick & easy disassembly, assists thorough & deep cleaning
- Available in mild or stainless steel (requires no on-site assembly)
  - Fully welded mild steel frame, powder coated in a durable epoxy paint
  - · Seamlessly welded, one piece 304 grade stainless steel frame
- Ergonomic push handle fitted to one side of the trolley ensures correct posture & keeps feet away from the frame, ideal for mobile phlebotomy or in cases where the trolley is wheeled for long periods of time
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Available with a solid grade laminate door & back panel providing security for the contents of the travs
- Two locking options...
  - · CAM Lock with Wing Handle Supplied with two keys
  - · Electronic Push Button Lock Programmable four-digit code
- Trollev pre-configured with the following accessories...
  - Alcohol Gel Holder
  - Debris Bin & Bracket
  - Glove / Mask Box Holder, One Box
  - · Sharps Box Holder
- Translucent, impact-resistant travs & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each trav
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### CTPHLEM1/10/3S2D

#### Mild Steel

| Description                | No Door          | Door with CAM Lock  | Door with Electronic<br>Push Button Lock | Weight |
|----------------------------|------------------|---------------------|--|--------|
| 3 x 100mm, 2 x 150mm Trays | CTPHLEM1/10/3S2D | CTPHLEM1/10/3S2D/CA | CTPHLEM1/10/3S2D/EPB                     | 31.1kg |
| 4 x 100mm, 1 x 150mm Trays | CTPHLEM1/10/4S1D | CTPHLEM1/10/4S1D/CA | CTPHLEM1/10/4S1D/EPB                     | 30.8kg |
| 6 x 100mm Trays            | CTPHLEM1/10/6S   | CTPHLEM1/10/6S/CA   | CTPHLEM1/10/6S/EPB                       | 31.9kg |

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com

#### Dimensions (w x d x h)

- Overall Dimensions 610 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

- Mild Steel Paint Spearmint
- Stainless Steel Natural Finish

#### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
- Defibrillator Platform page 67
- Carerail Mounted Accessories page 68
- Laptop Holders page 74
- Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77



#### **Stainless Steel**

| Description                | No Door          | Door with CAM Lock  | Door with Electronic<br>Push Button Lock | Weight |
|----------------------------|------------------|---------------------|--|--------|
| 3 x 100mm, 2 x 150mm Trays | CTPHLES1/10/3S2D | CTPHLES1/10/3S2D/CA | CTPHLES1/10/3S2D/EPB                     | 31.1kg |
| 4 x 100mm, 1 x 150mm Trays | CTPHLES1/10/4S1D | CTPHLES1/10/4S1D/CA | CTPHLES1/10/4S1D/EPB                     | 30.8kg |
| 6 x 100mm Trays            | CTPHLES1/10/6S   | CTPHLES1/10/6S/CA   | CTPHLES1/10/6S/EPB                       | 31.9kg |

CTPHLES1/10/3S2D/CA



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric - Foot Switch in Graphite - EC5315/VGH



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC5305/VBB | Hydraulic                                     | 67.0kg |
| EC5310/VBB | Electric - Hand Switch                        | 67.0kg |
| EC5315/VBB | Electric - Foot Switch                        | 67.0kg |
| EC5320/VBB | Electric with Electric Backrest - Hand Switch | 69.3kg |
| EC5325/VBB | Electric with Electric Backrest - Foot Switch | 69.3kg |
|            |   |        |

## Phlebotomy Chair Three Section

#### Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (selfassembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

#### Dimensions (w x d x h)

• 770 x 1710 x 490 - 950mm

#### **Options**

· E06 - Battery Back Up

#### **Standard Colours**

- Paint White
- Vinyl Bristol Blue

sales@bristolmaid.com +44 (0)1258 484455

360mm

470mm

630mm

820mm

## Phlebotomy Chair Three Section, Bariatric

#### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- · Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- · Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (selfassembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

#### Dimensions (w x d x h)

• 990 x 1710 x 520 - 950mm

#### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

#### Options

· E18 - Battery Back Up



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Onyx Black - EC5340/VOB

www.bristolmaid.com





| Item       | Description                                   | Weight |
|------------|---|--------|
| EC5330/VBB | Electric - Hand Switch                        | 77.2kg |
| EC5335/VBB | Electric - Foot Switch                        | 77.2kg |
| EC5340/VBB | Electric with Electric Backrest - Hand Switch | 79.5kg |
| EC5345/VBB | Electric with Electric Backrest - Foot Switch | 79.5kg |



## Phlebotomy Chair Fixed Height

#### Features

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- BS 4875: Part1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Certified for use up to 32 stone (203kg)
- Stable four leg design
- Black finish tubular frame
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrest
- Armrest can be positioned on either side of the chair (additional armrest available if required)
- Maximum capacity 203kg

#### Dimensions (w x d x h)

• 620 x 580 x 800mm

#### **Standard Colours**

- Paint Black
- Vinyl Bristol Blue

#### Options

• E15 - Phlebotomy Armrest (Single)

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Fixed Height Phlebotomy Chair in Aquamarine - EC5000/VAQ





# **Portering Chairs**

- Portering Chair Rear Steer 268
- Portering Chair Front Steer 269
- Portering Chair Bariatric, Rear Steer 271
  - Options & Accessories 272



## Portering Chair Rear Steer, Nesting, Sliding Foot Rest

#### Features

- Nests to occupy reduced floor space & encourage tidy storage
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- Design ensures chair will not tip should the patient stand on the step
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- · Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 200mm fixed wheels/150mm total stop braking castors (cushioned)
- Maximum capacity 190.0kg

#### Dimensions (w x d x h)

- Overall Dimensions 660 x 840 x 1040mm
- Seat 540 x 420mm at Height 530mm
- Back 490 x 455mm

#### **Standard Colours**

- Paint Grey
- Upholstery Bristol Blue

#### **Options / Accessories**

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cylinder Holder, E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/SB Lap Seatbelt, Factory Fitted
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card

Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Rear Steer, Nesting, Sliding Foot Rest in Bristol Blue - G/200/RS/VBB



| ltem         | Description                            | Weight |
|--------------|--|--------|
| G/200/RS/VBB | Rear Steer, Nesting, Sliding Foot Rest | 29.0kg |

sales@bristolmaid.com

+44 (0)1258 484455

X-Ray/Notes Pocket

## Portering Chair Front Steer, Nesting, Sliding Foot Rest

#### **Features**

- · Nests to occupy reduced floor space & encourage tidy storage
- · Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- · Design ensures chair will not tip should patient stand on the step
- X-ray/notes pocket on rear of seat
- · Two hooks for patients belongings
- · Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 125mm swivel castors/200mm braking fixed wheels (cushioned)
- Maximum capacity 190.0kg
- Ergonomically designed over swept push handle gives greater control making chair easier to push/use

#### Dimensions (w x d x h)

- Overall Dimensions 660 x 1005 x 1090mm
- Seat 540 x 420mm at Height 600mm
- Back 490 x 455mm

#### **Standard Colours**

- · Paint Grey
- Upholstery Bristol Blue

#### **Options / Accessories**

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cylinder Holder, E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/SB Lap Seatbelt, Factory Fitted
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card



X-Ray/Notes Pocket



Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Front Steer, Nesting, Sliding Foot Rest in Lime - G/305/FS/VLI

|  | VGH |  |  |  |
|--|-----|--|--|--|

| Item         | Description                             | Weight |
|--------------|---|--------|
| G/305/FS/VBB | Front Steer, Nesting, Sliding Foot Rest | 29.0kg |

sales@bristolmaid.com

www.bristolmaid.com

+44 (0)1258 484455



## **Portering Chair** Front Steer, Nesting, Hinged Foot Rests

#### **Features**

- Nests to occupy reduced floor space & encourage tidy storage
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & helpers
- · Folding/hinging foot rests aids access as well as supporting feet during transit
- · X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 125mm swivel castors/200mm braking fixed wheels (cushioned)
- Maximum capacity 190.0kg

#### Dimensions (w x d x h)

- Overall Dimensions 660 x 710 x 1055mm
- Seat 540 x 420mm at Height 530mm
- Back 490 x 455mm

#### **Standard Colours**

- · Paint Grev
- Upholstery Bristol Blue

#### **Options / Accessories**

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cvlinder Holder, E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- · G/SB Lap Seatbelt, Factory Fitted
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card



| Item         | Description                             | Weight |
|--------------|---|--------|
| G/205/FS/VBB | Front Steer, Nesting, Hinged Foot Rests | 27.5kg |

+44 (0)1258 484455

X-Rav/Notes Pocket

## Portering Chair Bariatric, Rear Steer, Sliding Foot Rest

#### **Features**

- · Heavy duty bariatric portering chair
- · Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- · Design ensures chair will not tip should the patient stand on the step
- · X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- · Fittings to accept IV poles mounted to rear uprights
- · Retractable diagonal seatbelt
- 200mm fixed wheels/200mm total stop braking castors (cushioned)
- We recommend that two members of staff are in attendance when used with a patient in excess of 190kg
- Maximum capacity 380.0kg

#### Dimensions (w x d x h)

- Overall Dimensions 930 x 1000 x 1090mm
- Seat 740 x 670mm at Height 530mm
- Back 700 x 480mm

#### **Standard Colours**

- · Paint Grey
- Upholstery Bristol Blue

#### **Options / Accessories**

- G/500/RS/CYL Cylinder Holder, CD/D
- G/CLS/500 Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card



X-Ray/Notes Pocket





Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Front Steer, Nesting, Sliding Foot Rest in Lime - G/305/FS/VLI

|  | VGH |  |  |  |
|--|-----|--|--|--|

| Item         | Description                              | Weight |
|--------------|--|--------|
| G/500/RS/VBB | Bariatric, Rear Steer, Sliding Foot Rest | 43.8kg |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com

**Portering Chairs** 

## **Cylinder Holders**

#### **Features**

- Compatible with 100mm diameter cylinders (C/CD/D/E)
- . Mounts to IV support tube at rear of chair using a handwheel
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models

#### **Standard Colours**

Paint - Grey



| Item         | Description             | Weight |
|--------------|-------------------------|--------|
| CYHP/D1      | C/CD/D Size             | 1.0kg  |
| CYHP/E1      | E Size                  | 1.8kg  |
| G/500/RS/CYL | Bariatric - C/CD/D Size | 1.0kg  |

## Calf / Leg Support - Pair

#### **Features**

- Moulded calf/leg support on captive retractable support tube
- Retracts beneath the seat when not in use
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models
- Supplied as a pair

#### **Standard Colours**

Stainless Steel/Black



| Item      | Description                      | Weight |
|-----------|----------------------------------|--------|
| G/CLS     | Calf/Leg Support, Pair           | 2.3kg  |
| G/CLS/500 | Bariatric Calf/Leg Support, Pair | 2.3kg  |

## **Infusion Pole**

#### **Features**

- Stainless steel infusion pole with a U shaped hook (length of pole 980mm)
- Mounts to IV support tube at rear of chair using a handwheel
- · Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models



## Lap Seatbelt

#### **Features**

- Lap seat belt with snap buckle
- Compatible with G/200/RS, G/205/FS & G/305/FS models
- · Factory fitted

#### **Standard Colours**

Black





| Item | Description  | Weight |
|------|--------------|--------|
| G/SB | Lap Seatbelt | 0.4kg  |

## **Card/Coin Lock**

#### Features

- Compatible with old & new one pound coin (euro compatible)
- Coin/card remains in the lock until the chair is connected to another chair or the head of station
- Plastic release card available for use by porters
- Factory fitted
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models
- Replacement card/coin lock cassette available (requires cassette removal tool)



| Item        | Description                                     | Weight |
|-------------|---|--------|
| G/CARD/LOCK | Card/Coin Lock, Factory Fitted                  | 1.0kg  |
| 2LOKE011    | Plastic Release Card                            | 0.01kg |
| 2LOCK217    | Card/Coin Lock Cassette, New £1 Coin Compatible | 0.2kg  |
| 2LOCK218    | Cassette Removal Tool                           | 0.01kg |

## ItemDescriptionWeightG/CARD/STATIONHead of Station1.6kg

#### **Features**

- Head of station for chairs fitted with card/coin locks
- Requires fitting to the fabric of the building (wall fixings not supplied)
- Chain length 1000mm
- One required for each row of chairs, (G/200/RS, G/205/FS & G/305/FS)
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models





## **Procedure Trolleys & Carts**

- Caretray Trolley Anaesthesia 276
  - Caretray Trolley Cast 278
- Caretray Trolley Critical Care 280
- Caretray Trolley Infection 282
- Caretray Trolley IV Preparation 284
- Caretray Trolley Medical / Surgical 286
  - Caretray Trolley Plaster 288
- Caretray Trolley Vaccination290Caretray Trolley Basics292Mobile Workstation Mild Steel, Low Level294Procedure Carts Polymer, Metro Flexline™296Procedure Carts Polymer, Metro Starsys™298Procedure Carts Polymer, Metro Basics™300

#### CTANAEM1/10/3S2D/HSB

| Item                    | Description   | Weight |
|-------------------------|---|--------|
| High Security Bolt Lock |   |        |
| CTANAEM1/10/3S2D/HSB    | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 41.4kg |
| CTANAEM1/10/4S1D/HSB    | 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 41.1kg |
| CTANAEM1/10/6S/HSB      | 6 x 100mm Trays, High Security Bolt Lock            | 42.2kg |

## **Caretrav Trollev - Anaesthesia**

#### **Features**

- Trollevs suitable for use with Anaesthesia
- Quick & easy disassembly, assists thorough & deep cleaning
- · Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trollev pre-configured with the following accessories...
  - Overbridge, Mild Steel
  - · Four Compartment Tilt Bin, Mild Steel
  - · Five Compartment Tilt Bin, Mild Steel
  - Waste Bin & Bracket, 8 Litre
- · Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- · 100mm polymer castors, 2 off braking

#### Standard Colours · Paint - Grey White

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Travs
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Options / Accessories**

- Locking & Security page 60
  Overbridge & Tilt Bins page 66
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTANAEM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 41.4kg |
| CTANAEM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 41.1kg |
| CTANAEM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 42.2kg |







| Item                      | Description   | Weight |
|---------------------------|---|--------|
| CAM Lock with Wing Handle |   |        |
| CTCASTM1/10/3S2D/CA       | 3 x 100mm, 2 x 150mm Trays, CAM Lock with Wing Handle | 43.5kg |
| CTCASTM1/10/4S1D/CA       | 4 x 100mm, 1 x 150mm Trays, CAM Lock with Wing Handle | 43.2kg |
| CTCASTM1/10/6S/CA         | 6 x 100mm Trays, CAM Lock with Wing Handle            | 44.3kg |

## **Caretray Trolley - Cast**

#### **Features**

- Trolleys suitable for use within A&E / Plaster Room
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · CAM Lock with Wing Handle Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Overbridge, Mild Steel
  - Four Compartment Tilt Bin, Mild Steel
  - Bowl Support & Bracket
  - Stainless Steel Bowl, 7 Litre
  - Waste Bin & Bracket, 8 Litre
  - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### **Standard Colours**

· Paint - Grey White

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

sales@bristolmaid.com +44 (0)1258 484455

#### **Options / Accessories**

- Locking & Security page 60
  Overbridge & Tilt Bins page 66
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTCASTM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 43.7kg |
| CTCASTM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 43.4kg |
| CTCASTM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 44.5kg |



CTCASTM1/10/3S2D/EPB



| Item                      | Description   | Weight |
|---------------------------|---|--------|
| CAM Lock with Wing Handle |   |        |
| CTCRITM1/10/3S2D/CA       | 3 x 100mm, 2 x 150mm Trays, CAM Lock with Wing Handle | 44.6kg |
| CTCRITM1/10/4S1D/CA       | 4 x 100mm, 1 x 150mm Trays, CAM Lock with Wing Handle | 44.3kg |
| CTCRITM1/10/6S/CA         | 6 x 100mm Trays, CAM Lock with Wing Handle            | 45.4kg |

## **Caretray Trolley - Critical Care**

#### **Features**

- Trolleys suitable for use with Critical Care
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · CAM Lock with Winged Handle Supplied with two keys
  - Electric Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Overbridge, Mild Steel
  - Five Compartment Tilt Bin, Mild Steel
  - Six Compartment Tilt Bin, Mild Steel
  - Alcohol Gel Holder
  - Infusion Pole & Bracket
  - Waste Bin & Bracket, 8 Litre
  - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### **Standard Colours**

· Paint - Grey White

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Options / Accessories**

- Locking & Security page 60
  Overbridge & Tilt Bins page 66
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTCRITM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 44.8kg |
| CTCRITM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 44.5kg |
| CTCRITM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 45.6kg |



CTCRITM1/10/3S2D/EPB

σ



CTINFEM1/10/3S2D/HSB

| Description   | Weight   |
|---|--|
|   |  |
| 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 35.5kg   |
| 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 35.2kg   |
| 6 x 100mm Trays, High Security Bolt Lock            | 36.3kg   |
|   | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock<br>4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock |

sales@bristolmaid.com

## **Caretray Trolley - Infection**

#### **Features**

- Trolleys designed for use when treating patients that have been infected with multi-resistant pathogens
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Waste Bin & Bracket, 8 Litre
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### **Standard Colours**

- Dimensions (w x d x h)
- Paint Grey White
   Overall Dimensions 550 x 550 x 993mm
  - Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
  - Plastic Top 455 x 455 (internal) at height 968mm
  - Trays
    - Shallow Tray 355 x 405 x 100mm (internal)
    - Deep Tray 355 x 405 x 150mm (internal)

#### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
- Laptop Holders page 74
- Castors & Ring Buffers page 76
   Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTINFEM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 35.5kg |
| CTINFEM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 35.2kg |
| CTINFEM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 36.3kg |



CTINFEM1/10/3S2D/EPB



| Item                    | Description   | Weight |
|-------------------------|---|--------|
| High Security Bolt Lock |   |        |
| CTIVPRM1/10/3S2D/HSB    | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 39.8kg |
| CTIVPRM1/10/4S1D/HSB    | 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 39.5kg |
| CTIVPRM1/10/6S/HSB      | 6 x 100mm Trays, High Security Bolt Lock            | 40.6kg |

sales@bristolmaid.com

### **Caretray Trolley - IV Preparation**

#### **Features**

- Trolleys designed for use when preparing & administering intravenous infusions
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Infusion Pole & Bracket
  - Waste Bin & Bracket, 8 Litre
  - Glove / Mask Box Holder, One Box
  - · Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking
- When not in use, trolleys storing medicines should be secured to the wall using a Security Clamp (see locking & security accessories)

## • Paint - Grey White

#### **Dimensions** (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

+44 (0)1258 484455

#### **Options / Accessories**

- Locking & Security page 60
   Overbridge & Tilt Bins page 66
   Defibrillator Platform page 67
   Carerail Mounted page 68

- Laptop Holders page 74
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTIVPRM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 39.8kg |
| CTIVPRM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 39.5kg |
| CTIVPRM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 40.6kg |



+44 (0)1258 484455



| Item                      | Description   | Weight |
|---------------------------|---|--------|
| CAM Lock with Wing Handle |   |        |
| CTMEDSM1/10/3S2D/CA       | 3 x 100mm, 2 x 150mm Trays, CAM Lock with Wing Handle | 43.2kg |
| CTMEDSM1/10/4S1D/CA       | 4 x 100mm, 1 x 150mm Trays, CAM Lock with Wing Handle | 42.9kg |
| CTMEDSM1/10/6S/CA         | 6 x 100mm Trays, CAM Lock with Wing Handle            | 44.0kg |

CTMEDSM1/10/3S2D/CA

### **Caretray Trolley - Medical / Surgical**

#### **Features**

- Trolleys designed for use within the Medical / Surgical Environment
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · CAM Lock with Wing Handle Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Overbridge, Mild Steel
  - Four Compartment Tilt Bin, Mild Steel
  - Infusion Pole & Bracket
  - Waste Bin & Bracket, 8 Litre
  - Glove / Mask Box Holder, One Box
  - · Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### **Standard Colours**

· Paint - Grey White

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

sales@bristolmaid.com +44 (0)1258 484455

#### **Options / Accessories**

- Locking & Security page 60
  Overbridge & Tilt Bins page 66
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock | (   |        |
| CTMEDSM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 43.4kg |
| CTMEDSM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 43.1kg |
| CTMEDSM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 44.0kg |



CTMEDSM1/10/3S2D/EPB

O



# **Caretray Trolley - Plaster**

#### **Features**

- Trolleys suitable for use within A&E / Plaster Room
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Trolley pre-configured with the following accessories...
  - Bucket Support & Bracket
  - Stainless Steel Graduated Bucket, 15 Litre
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions 1005 x 550 x 993mm
- Plastic Top 915 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

## **Standard Colours**

• Paint - Grey White

| Item             | Description                | Weight |
|------------------|----------------------------|--------|
| CTPLASM2/10/6S4D | 6 x 100mm, 4 x 150mm Trays | 50.6kg |
| CTPLASM2/10/8S2D | 8 x 100mm, 2 x 150mm Trays | 49.9kg |

sales@bristolmaid.com +44 (0)1258 484455

#### **Options / Accessories**

- Carerail Mounted page 68
  Laptop Holders page 74
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77







## **Caretray Trolley - Vaccination**

#### **Features**

- Trolleys suitable for use with vaccination
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Trolley pre-configured with the following accessories...
  - Alcohol Gel Holder
  - Debris Bin & Bracket
  - Glove, Mask Box Holder, One Box
  - Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

## **Standard Colours**

### Dimensions (w x d x h)

- Paint Grey White
- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

| Item             | Description                | Weight |
|------------------|----------------------------|--------|
| CTVACIM1/10/3S2D | 3 x 100mm, 2 x 150mm Trays | 32.0kg |
| CTVACIM1/10/4S1D | 4 x 100mm, 1 x 150mm Trays | 31.7kg |
| CTVACIM1/10/6S   | 6 x 100mm Trays            | 32.8kg |

#### **Options / Accessories**

- Locking & Security page 60
   Overbridge & Tilt Bins page 66
   Defibrillator Platform page 67
   Carerail Mounted page 68

- Laptop Holders page 74
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77







|          | C |  |
|----------|---|--|
| <b>S</b> |   |  |

| Item             | Description            | Weight |
|------------------|------------------------|--------|
| CTBASIM1/10      | Top Panel              | 14.2kg |
| CTBASIM1/10/1FWS | 1 x Full Width Shelf   | 15.5kg |
| CTBASIM1/10/2FWS | 2 x Full Width Shelves | 16.8kg |

# **Caretray Trolley - Basics**

#### **Features**

- Standard level single column easy clean procedure trolleys
- Design allows for the trolley to be upgraded with side panels & travs
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top...
  - Smooth surfaces & radiused edges to assist cleaning (compatible with steam cleaning)
- Raised edges on three sides, providing defined grab points to act as push handles whilst an integral rim at the front contains spillages (max load 30kg)
- · Solid grade laminate full width shelves, removable & supported on the trolley frame (max load per shelf 5kg)
- · Carerails fitted to either side support a wide range of accessories
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Full Width Shelves 420 x 420 (internal) at heights 193 & 554mm

## **Standard Colours**

· Paint - Spearmint

#### **Options / Accessories**

- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
  Laptop Holders page 74

- Castors & Ring Buffers page 76
   Trays, Shelves, Dividers & Labels page 77





CTBASIM1/10/2FWS



## Mobile Workstation - Mild Steel, Low Level

#### Features

- Recessed liquid retaining top
- Single Point Locking Using a CAM lock Supplied with 2 keys
- Anti-tilt mechanism, allows only one drawer to be opened at any one time
- 85, 170 & 255mm deep drawers, fitted on full extension glides (max load 30kg per drawer)
- Two tone colour keyway, drawer fronts available in a range of colours
- Low level corner bumpers protect the fabric of the building
- 80mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall 615 x 500 x 810mm
- Drawers 480 x 375 x 85 / 170 / 255mm (internal)

#### **Standard Colours**

- Body Light Grey
- Drawers Dark Grey

| Item  | Description                            | Weight |
|-------|--|--------|
| MW005 | 2 x 85mm, 1 x 170mm, 1 x 255mm Drawers | 49.0kg |

```
www.bristolmaid.com
```





Single Point Locking Using a CAM lock - supplied with 2 keys)

Low level corner bumpers to protect the fabric of the building



## Procedure Carts - Polymer, Metro Flexline™

#### **Features**

- · Lightweight polymer construction with smooth rounded corners to assist cleaning
- Incorporates Microban<sup>®</sup> offering antimicrobial protection
- Upper work surface provides a writing surface / space for preparation work
- Model specific locking
  - CAM Lock Supplied with 2 keys
  - Keyless Entry Touch Pad Lock
  - Security Seal Lock Security seal compatible
- Full extension lightweight polymer drawers, depths 76, 152 & 229mm (max load 7 - 20kg per drawer)
- Two tone colour keyway, drawer handles available in a range of colours
- 127mm polymer castors, 2 off directional & 2 off total lock

#### **Dimensions** (w x d x h)

- Drawers
  - 76mm 514 x 432 x 57mm (internal)
  - 152mm 514 x 432 x 133mm (internal)
  - 229mm 514 x 432 x 210mm (internal)

#### Standard Colours

- Body Light Taupe
- Drawers Model Dependant

#### **Options / Accessories**

• 8FL190 - Drawer Label Holders (x 10)

| Item      | Description  | Overall Dimensions (w x d x h) | Weight |
|-----------|--|--------------------------------|--------|
| 8FLCRP4   | Emergency - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Code Blue Handles                  | 1056 x 607 x 1273mm            | 91.0kg |
| 8FLPED1   | Paediatric Emergency - 8 x 76mm, 1 x 152mm Drawers, Coloured Handles                   | 1056 x 607 x 1273mm            | 91.0kg |
| 8FLTMENT2 | Treatment - 2 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Slate Blue Handles                 | 819 x 568 x 1064mm             | 91.0kg |
| 8FLCOM-LT | Laptop - 1 x Keyboard Tray, 1 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Slate Blue Handles | 819 x 568 x 1064mm             | 83.0kg |
| 8FLXCLINE | Central Line - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles              | 819 x 568 x 1064mm             | 91.0kg |
| 8FLXCCU1  | Critical Care - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Red Handles                    | 840 x 569 x 1807mm             | 90.0kg |



8FLCRP4



|   | Emergency<br>8FLCRP4 | Paediatric Emergency<br>8FLPED1 | Treatment<br>8FLTMENT2 | Laptop<br>8FLCOM-LT | Central Line<br>8FLCLINE | Critical Care<br>8FLCCU1 |
|---|----------------------|---------------------------------|------------------------|---------------------|--------------------------|--------------------------|
| Locking & Security                      |                      |                                 |                        |                     |                          |                          |
| CAM Lock                                |                      |                                 | $\checkmark$           |                     |                          |                          |
| Security Seal Lock                      | $\checkmark$         | √                               |                        |                     |                          |                          |
| Keyless Entry Touchpad Lock             |                      |                                 |                        | $\checkmark$        | $\checkmark$             | $\checkmark$             |
| Drawer Seal Lock Bar                    |                      | √                               |                        |                     |                          |                          |
| Plastic Security Seals (x 100)          | $\checkmark$         | $\checkmark$                    |                        |                     |                          |                          |
| Overbridge & Tilt Bins                  |                      |                                 |                        |                     |                          |                          |
| Overbridge with Top Shelf & Hanger Rail |                      |                                 |                        |                     |                          | $\checkmark$             |
| Six Compartment Tilt Bin for Overbridge |                      |                                 |                        |                     |                          | $\checkmark$             |
| Defibrillator Platform                  |                      |                                 |                        |                     |                          |                          |
| Adjustable Defibrillator Platform       | $\checkmark$         | $\checkmark$                    |                        |                     |                          |                          |
| Drawer Trays & Dividers                 |                      |                                 |                        |                     |                          |                          |
| 76mm Drawer Divider Kit                 |                      |                                 | 2                      |                     | 2                        |                          |
| 152 / 229mm Drawer Divider Kit          |                      |                                 | 2                      |                     | 1                        |                          |
| 76mm Drawer Divider Tray & Dividers     | 3                    |                                 |                        |                     |                          |                          |
| 152mm Drawer Divider Tray & Dividers    | 2                    |                                 |                        |                     |                          |                          |
| Side Bins                               |                      |                                 |                        |                     |                          |                          |
| Side Bin, Locking                       | 4                    |                                 |                        |                     | 2                        | 4                        |
| Side Bin, Non-Locking                   |                      |                                 |                        | 3                   | 2                        | 4                        |
| Laptop Holders                          |                      |                                 |                        |                     |                          |                          |
| Laptop Arm                              |                      |                                 |                        | $\checkmark$        |                          |                          |
| Accessories                             |                      |                                 |                        |                     |                          |                          |
| Waste Basket (28 Litre)                 |                      |                                 |                        |                     | $\checkmark$             |                          |
| Oxygen Cylinder Holder                  | $\checkmark$         | √                               |                        |                     |                          |                          |
| Lockable Sharps Container               |                      |                                 |                        |                     | $\checkmark$             |                          |
| Cord Manager                            | $\checkmark$         |                                 |                        |                     |                          |                          |
| Hospital Grade 6, Outlet Strip & Holder | $\checkmark$         |                                 |                        |                     |                          |                          |
| Suction Pump Shelf                      | $\checkmark$         | √                               |                        |                     |                          |                          |
| Rear Fitting Backboard                  | $\checkmark$         | √                               |                        |                     |                          |                          |
| Pull Out-Side Shelf                     | $\checkmark$         | √                               |                        |                     |                          |                          |
| Two Hook Infusion Pole                  | $\checkmark$         | $\checkmark$                    |                        |                     |                          |                          |



## Procedure Carts - Polymer, Metro Starsys<sup>™</sup>

#### **Features**

- · Lightweight polymer construction with smooth rounded corners to assist cleaning
- Incorporates Microban® offering antimicrobial protection
- Upper work surface provides a writing surface / space for preparation work
- CAM lock Supplied with 2 keys\*
- Full extension lightweight polymer drawers, depths 76, 152 & 229mm (max load 7 - 20kg per drawer)
- Two tone colour keyway, drawer handles available in a range of colours
- 127mm polymer castors, 2 off braking
- \*model specific

#### Dimensions (w x d x h)

- Drawers
  - 76mm 412 x 400 x 60mm (internal)
  - 152mm 412 x 400 x 136mm (internal)
  - 229mm 412 x 400 x 225mm (internal)

#### **Standard Colours**

- Body Light Taupe
- Drawers Handles Model Dependant

#### **Options / Accessories**

• 8FL190 - Drawer Label Holders (x 10)



| Item       | Description   | Overall Dimensions (w x d x h) | Weight  |
|------------|---|--------------------------------|---------|
| 8SXRSANES  | Anaesthesia - 4 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock        | 882 x 633 x 1136mm             | 104.0kg |
| 8SXRSDRS   | Dressing - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock | 577 x 633 x 1060mm             | 83.0kg  |
| 8SXRSTRAMA | Trauma - 3 x 76mm, 4 x 152mm Drawers, Slate Blue Handles, CAM Lock              | 882 x 633 x 1136mm             | 116.0kg |
| 8SXRSCCU   | Critical Care - 1 x 76mm, 1 x 152mm, 3 x 229mm Drawers, Slate Blue Handles      | 882 x 633 x 1212mm             | 107.0kg |
| 8SXRSARTH  | Arthroscopic - 6 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock       | 882 x 633 x 1136mm             | 110.0kg |
| 8SXRSIV    | IV Therapy - 5 x 76mm, 2 x 229mm Drawers, Green Handles, CAM Lock               | 882 x 633 x 1136mm             | 106.0kg |
| 8SXRSCAST  | Cast - 2 x 76mm, 2 x 152mm, 2 x 229mm Drawers, Slate Blue Handles               | 882 x 633 x 1212mm             | 109.0kg |



|   | Anaesthesia<br>8SXRSANES | Dressing<br>8SXRSDRS | Trauma<br>8SXRSTRAMA | Critical Care<br>8SXRSCCU | Arthroscopic<br>8SXRSARTH | IV Therapy<br>8SXRSIV | Cast<br>8SXRSCAST |
|---|--------------------------|----------------------|----------------------|---------------------------|---------------------------|-----------------------|-------------------|
| Locking & Security  |                          |                      |                      |                           |                           |                       |                   |
| CAM Lock  |                          |                      |                      | $\checkmark$              |                           |                       | $\checkmark$      |
| Overbridge & Accessory Rails                              |                          |                      |                      |                           |                           |                       |                   |
| Overbridge with One Shelf & Accessory Rail                |                          |                      | ~                    |                           | $\checkmark$              |                       |                   |
| Overbridge with Two Accessory Rails                       |                          |                      | ~                    | $\checkmark$              | $\checkmark$              | $\checkmark$          | $\checkmark$      |
| Drawer Trays & Dividers                                   |                          |                      |                      |                           |                           |                       | $\checkmark$      |
| 76mm Drawer Divider Kit                                   |                          |                      |                      | $\checkmark$              |                           |                       | $\checkmark$      |
| Pods & Tilt-Bins  |                          |                      |                      |                           |                           |                       |                   |
| Left Swing Out, 3 x Tilt-Bins                             |                          |                      |                      | $\checkmark$              |                           |                       | $\checkmark$      |
| Left Swing Out, 3 x Locking Tilt-Bins (CAM Lock)          |                          |                      |                      |                           |                           | $\checkmark$          |                   |
| Left Swing Out Locking, 3 x Locking Tilt-Bins (CAM Lock)  | ~                        |                      |                      |                           |                           |                       |                   |
| Left Swing Out, Top Locking Tilt-Bin (CAM Lock)           |                          |                      | √                    |                           | √                         | $\checkmark$          |                   |
| Right Fixed Pod, 1 x Tilt-Bin                             |                          |                      |                      | $\checkmark$              |                           |                       | √                 |
| Right Swing Out   |                          |                      |                      |                           | √                         |                       |                   |
| Right Swing Out Locking, 3 x Locking Tilt-Bins (CAM Lock) | $\checkmark$             |                      |                      |                           |                           |                       |                   |
| Right Swing Out, Top Locking Tilt-Bin (CAM Lock)          |                          |                      | √                    |                           |                           |                       |                   |
| Middle Non-Locking Tilt-Bin                               |                          |                      | √                    |                           |                           |                       |                   |
| Pod Tray & Dividers                                       |                          |                      | ✓                    | $\checkmark$              | $\checkmark$              | $\checkmark$          | $\checkmark$      |
| Tilt-Bin  |                          |                      |                      |                           | √                         |                       |                   |
| Accessories   |                          |                      |                      |                           |                           |                       |                   |
| Waste Basket & Holder                                     |                          |                      | $\checkmark$         | $\checkmark$              | ✓                         |                       | $\checkmark$      |
| Universal Clamp   | ~                        |                      |                      |                           |                           |                       |                   |
| Laryngoscope Blade / Specimen Holder                      | ~                        |                      |                      |                           |                           |                       |                   |
| Half Size Utility Bin (Metal)                             | √                        |                      |                      |                           |                           |                       |                   |
| Full Size Utility Bin (Metal)                             | $\checkmark$             |                      |                      |                           |                           |                       |                   |
| Label / Tape Dispenser (x 3)                              | $\checkmark$             |                      |                      |                           |                           |                       |                   |
| Plastic Utility Bins with Clear Covers (x 6)              | $\checkmark$             |                      |                      | $\checkmark$              | $\checkmark$              |                       |                   |
| Sharps Container Bracket                                  |                          |                      | $\checkmark$         |                           | $\checkmark$              |                       |                   |
| Two Hook Infusion Pole                                    |                          |                      | $\checkmark$         |                           |                           | $\checkmark$          |                   |
| Wire Basket   |                          |                      | $\checkmark$         | $\checkmark$              | $\checkmark$              |                       | $\checkmark$      |



 $\overline{\mathbf{O}}$ 

## Procedure Carts - Metro Basix™

#### **Features**

σ

- · Recessed liquid retaining plastic work surface
- Model specific locking
  - CAM Lock Supplied with 2 keys
  - Push Button Lock
  - Security Seal Lock Security seal compatible
- Full extension drawers, depths 76, 152 & 230mm (max load 7 24kg per drawer)
- Two tone colour keyway, drawer handles available in a range of colours
- Low level bumpers fitted to Trauma & Anaesthesia model
- 127mm polymer castors, 2 off directional & 2 off total lock

#### Dimensions (w x d x h)

- Compact
  - Overall Dimensions 489 x 470 x 876mm
  - Drawers 381 x 432 x 76 / 152 / 230mm (internal)
- Trauma / Anaesthesia
  - Overall Dimensions 762 x 613 x 1440mm
  - Drawers 432 x 565 x 76 / 152 / 230mm (internal)

#### **Standard Colours**

- Body Light Taupe
- Drawers Dark Blue (see below for alternative colours)

#### **Options / Accessories**

• 8MBA172 - Label Holders (x 12)

Alternative drawer colours are available. Change suffix 'DB' on the product code with any of the following i.e. - Compact Basix Cart with red drawers 8MBC1210RE



| $\checkmark$ | Item             | Description  |        |
|--------------|------------------|--|--------|
| $\checkmark$ | 8MBC1210DB       | Compact - 1 x 76mm, 2 x 152mm, 1 x 230mm Drawers     | 63.0kg |
| $\checkmark$ | 8MBX3210TRAUMADB | Trauma - 3 x 76mm, 2 x 152mm, 1 x 230mm Drawers      | 78.0kg |
|              | 8MBX1310ANES4DB  | Anaesthesia - 1 x 76mm, 3 x 152mm, 1 x 230mm Drawers | 81.0kg |





# **Resuscitation Trolleys & Carts**

- Resuscitation Trolleys Mild Steel 316
  - Caretray Trolley Resuscitation 322
- Resuscitation Carts Polymer, Metro Lifeline<sup>™</sup> 324
- Resuscitation Carts Polymer, Metro Flexline<sup>™</sup> 326
- Airways / Difficult Intubation Trolley Mild Steel 327
- Caretray Trolley Airways / Difficult Intubation 328



Drop down cupboard door prevents drawers from being opened. Can be secured using security seals.



## **Resuscitation Trolleys - Mild Steel**

#### **Features**

- Popular resuscitation trolley
- Drawer contents are secured when the lower cupboard is closed (security seal compatible)
- Upper drawer supplied with a plastic drawer liner, with dividers
- Recessed liquid retaining stainless steel work surface (usable 630 x 445)
- Lift up work flap (usable 305 x 455mm)
- Two hook infusion pole mounted to rear of trolley
- Low level bumper bar
- · Push handles one end
- 125mm pressed steel castors

#### Dimensions (w x d x h)

- Drawers 600 x 430 x 70 / 96 / 155mm (internal)
- Lower Cupboard 625 x 475 x 295mm (internal)

#### **Standard Colours**

- · Paint Red
- Lift Up Flap White

| Item  | Description                                       | Overall Dimensions (w x d x h) | Weight |
|-------|---|--------------------------------|--------|
| WT200 | 2 x 70mm & 1 x 155mm Deep Drawers, Lower Cupboard | 825 x 575 x 999mm              | 74.0kg |
| WT210 | 3 x 96mm Deep Drawers, Lower Cupboard             | 825 x 575 x 999mm              | 75.1kg |
| WT250 | 4 x 70mm & 1 x 155mm Deep Drawers, Lower Cupboard | 825 x 575 x 1165mm             | 86.0kg |

www.bristolmaid.com



WT200



## **Adjustable Defibrillator Platform**

#### **Features**

- Compatible with all popular defibrillators
- Defibrillator supported on a 365mm wide platform & is held securely in place using a clamp (adjusts using a handwheel 205 380mm)
- Platform mounted on an arm that swings 90 degrees from the worksurface & can be locked in position using a tensioning knob (rotates 360 degrees)
- Base of platform sits 150mm above the top of the trolley
- Fitted to the rear left-hand side of the trolley
- Maximum capacity 6kg
- Factory fitted

#### Dimensions (w x d x h)

 477 x 365 x 642mm (increases overall height by 185mm)

#### **Standard Colours**

· Paint - Grey White



#### WT/DEFIB

 Item
 Description

 WT/DEFIB
 Adjustable Defibrillator Platform

# **Suction Pump Shelf**

#### **Features**

- Compatible with all popular suction pumps
- Shelf provides an area 415 x 200mm (w x d)
- 30mm upstand at front of shelf
- Floor to base of shelf at height 260mm
- Corner bumpers offer protection to both the trolley & fabric of the building
- Maximum capacity 15kg
- Factory fitted

#### Dimensions (w x d x h)

• Overall Dimensions - 420 x 205 x 380mm

#### **Standard Colours**

Paint - Red



WT/SPS shown on WT/200

| Item   | Description        | Weight |
|--------|--------------------|--------|
| WT/SPS | Suction Pump Shelf | 5.0kg  |

Weight

5.5kg

## **Cardiac Board**

#### Features

- Laminate cardiac board approx. 500 x 400 x 17mm
- Positioned at the rear of the trolley for easy access
- Factory fitted



#### Features

- WT/CYL/D/E compatible with AD, CD, D, DD, PD, RD, ZD & E cylinders, max diameter 100mm
- WT/CYL/F compatible with F, LF, VF, HX & ZX cylinders, max diameter 145mm
- Maximum two-cylinder holders per trolley
- Factory fitted



- Plastic drawer liner with removable dividers
- One tidy supplied as standard with each trolley





| • |  |
|---|--|
|   |  |
|   |  |

| Item      | Description   | Weight |
|-----------|---------------|--------|
| WT/CBOARD | Cardiac Board | 1.5kg  |

| Item       | Description                  | Weight |
|------------|------------------------------|--------|
| WT/CYL/D/E | D/E Size, Max Diameter 100mm | 3.0kg  |
| WT/CYL/F   | F Size, Max Diameter 145mm   | 3.4kg  |

| Item  | Description | Weight |
|-------|-------------|--------|
| WT/DT | Drawer Tidy | 1.6kg  |

| em | Description | Weight |
|----|-------------|--------|



# Four Gang Socket

#### **Features**

- Switched & fused four gang socket
- Supplied with 3 metres of beige cable & fitted with a UK 3 pin plug
- Fitted to the rear of the trolley, complete with a cleat to secure the cable when not in use
- Factory fitted

## **Guard Rail**

#### Features

- Stainless steel guard rail
- Increases the height of the trolley by 60mm
- Factory fitted

# **Security Seals**

#### **Features**

- Snap on plastic securityPack of 100







| Item       | Description      | Weight |
|------------|------------------|--------|
| E/SOCKET/4 | Four Gang Socket | 1.4kg  |

| Item     | Description | Weight |
|----------|-------------|--------|
| G/RAIL-C | Guard Rail  | 1.5kg  |

| Item     | Description    |
|----------|----------------|
| 2PLCO026 | Security Seals |

## **Heavy Duty Castors**

#### Features

- 125mm heavy duty pressed steel castorsOffers additional cushioning, lowering noise levels
- Factory fitted



6

www.bristolmaid.com



# **Caretray Trolley - Resuscitation**

#### **Features**

- Trolleys suitable for use with Resuscitation
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Door fitted with a security seal lock (security seal compatible)
- Trolley pre-configured with the following accessories...
  - Cylinder Holder, C/CD/D
  - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

#### **Options / Accessories**

• Paint - Red

Defibrillator Platform - page 67

| Item                 | Description                | Weight |
|----------------------|----------------------------|--------|
| CTRESUSM1/10/3S2D/SS | 3 x 100mm, 2 x 150mm Trays | 38.3kg |
| CTRESUSM1/10/4S1D/SS | 4 x 100mm, 1 x 150mm Trays | 38.0kg |
| CTRESUSM1/10/6S/SS   | 6 x 100mm Trays            | 39.1kg |

# **Overbridge**

#### Features

- · Fixes onto the rear uprights of the trolley
- Provides a fixing location for tilt bins (order separately)
- Factory fitted

#### Standard Colours

· Paint - Grey White



# **Tilt Bins**

## **Features**

- · Available with 4.5 or 6 compartments
- · Clear bins provide visible storage & assists organisation
- Bins tilt to 45° allowing easy access of contents, & can be removed to assist refilling / cleaning
- · Label holder for easy identification
- · Retaining bar secures bins in place, preventing from tipping forward during transit
- Supplied with a mild steel mounting

### Dimensions (w x d x h)

TBIN6/MS - 89 x 67 x 95mm

### Standard Colours

**Compartment Dimensions** 

45° Angle Tilt

Label Holder

+44 (0)1258 484455

- · Plastic Cream
- TBIN4/MS 143 x 127 x 190mm • TBIN5/MS - 108 x 95 x 146mm
- · Paint Grey White



Middle: TBIN5/MS Bottom: TBIN4/MS

| Item     | Description      | Overall Dimensions (w x d x h) | Weight |
|----------|------------------|--------------------------------|--------|
| TBIN4/MS | Four Compartment | 600 x 171 x 207mm              | 2.3kg  |
| TBIN5/MS | Five Compartment | 600 x 135 x 164mm              | 1.6kg  |
| TBIN6/MS | Six Compartment  | 600 x 93 x 113mm               | 0.9kg  |

Weight

2.0kg

www.bristolmaid.com



## Resuscitation Carts - Polymer, Metro Lifeline™

#### **Features**

- Range of pre-configured emergency carts
- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Incorporates Microban<sup>®</sup> offering antimicrobial protection
- Single seal mechanism secures the top compartment, drawers & tiltout side bins (security seal compatible)
- Top storage compartment features a clear removable cover, & provides instant access to first-line equipment
- Strapless defibrillator platform on swing arm adjusts to accommodate various models
- Full extension lightweight polymer drawers, depths 76, 152, 229 & 305mm, (max load 7 20kg per drawer)
- Ergonomic handle with fifth wheel steering provides maximum control, when attending an emergency situation
- Two tone colour keyway
- 127mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

• Drawers - 514 x 432 x 76 / 152 / 229 / 305mm (internal)

#### **Standard Colours**

- · Body Blue
- Drawers Light Taupe / Code Blue

#### **Options / Accessories**

• 8FL190 - Drawer Label Holders (x 10)

| Item     | Description                            | Overall Dimensions (w x d x h) | Weight |
|----------|--|--------------------------------|--------|
| 8LECCRP2 | 1 x 76mm, 3 x 152mm Drawers            | 965 x 575 x 984mm              | 84.0kg |
| 8LECCRP3 | 2 x 76mm, 1 x 152mm, 1 x 229mm Drawers | 965 x 575 x 984mm              | 88.0kg |
| 8LECCRP4 | 2 x 76mm, 1 x 152mm, 1 x 305mm Drawers | 965 x 575 x 1070mm             | 87.0kg |



Drawer Label Holders





|                                   | 8LECCRP2     | 8LECCRP3     | 8LECCRP4     |
|-----------------------------------|--------------|--------------|--------------|
| Top Cavity Tray                   |              | $\checkmark$ |              |
| Locking Side Bin (x 3)            |              | $\checkmark$ | $\checkmark$ |
| Waste Basket (28 Litre)           | $\checkmark$ |              |              |
| Lockable Sharps Container         | $\checkmark$ |              |              |
| Adjustable Defibrillator Platform | $\checkmark$ | $\checkmark$ | $\checkmark$ |
| Suction Pump Shelf                |              | $\checkmark$ |              |
| Backboard (Front Fitting)         |              | $\checkmark$ | $\checkmark$ |
| Backboard (Rear Fitting)          | $\checkmark$ |              |              |
| Oxygen Cylinder Holder            | √            | $\checkmark$ | $\checkmark$ |
| Two Hook Infusion Pole            | √            | $\checkmark$ | $\checkmark$ |
| Storage Bin                       | $\checkmark$ | $\checkmark$ | $\checkmark$ |
| Plastic Security Seals (x 100)    | $\checkmark$ | $\checkmark$ | $\checkmark$ |



# **Resuscitation Carts - Polymer, Metro Flexline™**

#### **Features**

- Lightweight polymer construction with smooth rounded corners to
   assist cleaning
- Incorporates Microban<sup>©</sup> offering antimicrobial protection
- Upper work surface provides a writing surface / space for preparation work
- Single seal mechanism secures the drawers (security seal compatible)
- Full extension lightweight polymer drawers, depths 76, 152 & 229mm (max load 7 20kg per drawer)
- Two tone colour keyway
- 127mm polymer castors, 2 x directional & 2 x total lock

## Dimensions (w x d x h)

- Overall Dimensions 819 x 568 x 914mm
- Drawers
  - 76mm 514 x 432 x 57mm (internal)
  - 152mm 514 x 432 x 133mm (internal)
  - 229mm 514 x 432 x 210mm (internal)

## **Standard Colours**

- Body Light Taupe
- Drawers Light Taupe / Red

## **Options / Accessories**

- 8FL141 76mm Egg Crate Style Drawer Divider Kit (10 x Vertical, 5 x Horizontal)
- 8FL142 152mm Egg Crate Style Drawer Divider Kit (10 x Vertical, 5 x Horizontal)
- 8FL190 Drawer Label Holders (x 10)
- 8FL234 100mm Diameter Oxygen Cylinder Holder
- 8FL315 Two Hook Infusion Pole
- 8LEC320 Plastic Security Seals (x 100)
- 8LEC304 Defibrillator Platform

## Airways / Difficult Intubation Trolley - Mild Steel

#### Features

- · Recessed liquid retaining top
- Available as non-locking or fitted with a tamper evident panel to secure the contents of the drawers (security seal compatible)
- 135mm deep drawers, fitted on full extension glides (max load 30kg per drawer)
- Carerails fitted to L/H side support a wide range of accessories
- Low level corner bumpers protect the fabric of the building
- 80mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall 615 x 500 x 915mm
- Drawers 480 x 375 x 135mm (internal)

#### Standard Colours

· Paint - White

#### **Options / Accessories**

Carerail Mounted - page 68



MW015 Tamper Evident Panel



| MW0 |  |
|-----|--|
|     |  |
|     |  |

| Item  | Description                             | Weight |
|-------|---|--------|
| MW010 | 5 x 125mm Drawers, Non Locking          | 55.0kg |
| MW015 | 5 x 125mm Drawers, Tamper Evident Panel | 57.0kg |

**Resuscitation Trolleys** 

ço



## **Caretray Trolley - Airways / Difficult Intubation**

#### **Features**

- Trolleys suitable for use with Airways / Difficult Intubation
- $\bullet\,$  Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Door fitted with a security seal lock (security seal compatible)
- Trolley pre-configured with the following accessories...
  - Alcohol Gel Holder
  - Airway Tube Holder
  - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

• Paint - Grey White

| Item                | Description                | Weight |
|---------------------|----------------------------|--------|
| CTAIRWM1/10/3S2D/SS | 3 x 100mm, 2 x 150mm Trays | 38.6kg |
| CTAIRWM1/10/4S1D/SS | 4 x 100mm, 1 x 150mm Trays | 38.3kg |
| CTAIRWM1/10/6S/SS   | 6 x 100mm                  | 39.4kg |

sales@bristolmaid.com

+44 (0)1258 484455

#### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
- Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77





# **Sackholders & Waste Bins**

- Removable Body Bins 320
- Free Standing Sackholders 325
- Fixed Body Plastic Sackholders Plastic 326
- Removable Body Plastic Sackholders Plastic 327



#### 20 Litre

| Item      | Description                            | Weight |
|-----------|--|--------|
| BR015/WHL | Hands Free                             | 8.5kg  |
| BR020/WHL | Hands Free & Silent Closing            | 8.5kg  |
| BR030/WHL | Hands Free, Rust Free & Silent Closing | 8.5kg  |

# Removable Body Bins - 20 / 50 Litre

#### **Features**

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- · Captive drip tray to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Choice of construction...
  - Mild steel tubular frame & sheet metal body, painted for extra durability
  - 'Rust Free' stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise when the lid closes

#### Dimensions (w x d x h)

- 20 Litre 267 x 346 x 552mm
- 50 Litre 444 x 397 x 551mm

#### **Standard Colours**

· Paint - Grey / White





'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)

www.bristolmaid.com



Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Removable Body Bin - 50ltr, Hands free with a green lid - BR045/GNL



#### 50 Litre

| Item      | Description                            | Weight |
|-----------|--|--------|
| BR045/WHL | Hands Free                             | 11.5kg |
| BR050/WHL | Hands Free & Silent Closing            | 11.5kg |
| BR060/WHL | Hands Free, Rust Free & Silent Closing | 11.5kg |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com



#### BR075/WHL

#### 75 Litre

| Item      | Description                            | Weight |
|-----------|--|--------|
| BR075/WHL | Hands Free                             | 15.0kg |
| BR080/WHL | Hands Free & Silent Closing            | 15.0kg |
| BR090/WHL | Hands Free, Rust Free & Silent Closing | 15.0kg |

## Removable Body Bins - 75 / 90 Litre

#### **Features**

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- · Captive drip tray to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- · Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Choice of construction...
  - Mild steel tubular frame & sheet metal body, painted for extra durability
  - 'Rust Free' stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise when the lid closes

#### Dimensions (w x d x h)

- 75 Litre 444 x 397 x 781mm
- 90 Litre 444 x 397 x 910mm

#### **Standard Colours**

· Paint - Grey / White





'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)



BR105/WHL

Sackholders & Waste Bins

Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Removable Body Bin - 75ltr, Hands free with a yellow lid - BR075/YEL



### 90 Litre

| Item      | Description                            | Weight |
|-----------|--|--------|
| BR105/WHL | Hands Free                             | 17.0kg |
| BR110/WHL | Hands Free & Silent Closing            | 17.0kg |
| BR120/WHL | Hands Free, Rust Free & Silent Closing | 17.0kg |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com






Positive bag retention through the use of a drop over clamp



Gas damper used on silent closing bins



Raised flange & recessed lid used on hands free bins

### **Removable Body Bins - 75 Litre, Front Opening**

### **Features**

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- Captive drip trav to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- · Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Choice of construction...
  - · Mild steel tubular frame & sheet metal body, painted for extra durability
  - 'Rust Free' stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise when the lid closes

### Dimensions (w x d x h)

75 Litre - 444 x 397 x 781mm

### Standard Colours

· Paint - Grev / White





'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)

Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Removable Body Bin - 75ltr. Hands Free & Silent Closing with a grey lid - BR140/GYL



BR135/WHL Hands Free BR140/WHL Hands Free & Silent Closing BR150/WHL Hands Free, Rust Free & Silent Closing

Description

| sales | @bristo | Imaid.com |  |
|-------|---------|-----------|--|

+44 (0)1258 484455

```
www.bristolmaid.com
```

Item

### Free Standing Sackholders - 50 / 75 Litre

### Features

- · Captive drip tray to collect spills
- · Positive bag retention through the use of an elasticated cord
- Stainless steel pedal & linkages
- · Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Mild steel tubular frame, painted for extra durability

### Dimensions (w x d x h)

- 50 Litre 442 x 397 x 551mm
- 75 Litre 442 x 397 x 780mm

### Standard Colours

· Paint - Grev



Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Free Standing Bag Holder - 50 Litre, Elastic Cord with brown lid - BH015/BNL



| Item      | Description                         | Weight |
|-----------|-------------------------------------|--------|
| BH015/WHL | Free Standing Sackholder - 50 Litre | 7.0kg  |
| BH030/WHL | Free Standing Sackholder - 75 Litre | 7.5kg  |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com



### Fixed Body Sackholders - Plastic, 20 Litre

### **Features**

- Meets NHS England 'Managing Healthcare Fire Safety' (HTM-05-01)
- Plastic construction with no metal parts
- Smooth surfaces to assist cleaning
- One-piece plastic lid, available in a range of colours to assist refuse identification (HTM07)
- Soft closing lid
- Bag positively secured using a drop over sack retention ring

### Dimensions (w x d x h)

Overall Dimensions - 350 x 220 x 660mm

### **Standard Colours**

- · Body White
- · Lid White

### Options

- Alternative Lid Colours Black / Orange / Yellow / Yellow & Black
- · Moulded In Graphics Available on request

Alternative lid colours available. Change suffix 'WH' on the item no. with any of the following, i.e. Fixed Body Sackholder - Plastic, Silent Closing, 20 Litre with a black lid - 5FSH20BL



```
www.bristolmaid.com
```

Item

# Sackholders & Waste Bins

# Removable Body Sackholders - Plastic, 42 & 70 Litre

### **Features**

- Meets NHS England 'Managing Healthcare Fire Safety' (HTM-05-03)
- One-piece removable body
- Plastic construction with no metal parts
- Smooth surfaces to assist cleaning
- One-piece plastic lid, available in a range of colours to assist refuse identification (HTM07)
- Soft operating lid on both opening & closing (moving parts enclosed)
- · Recessed drip tray to collect unwanted spills
- Bag positively secured using a drop over sack retention ring
- Non-marking rubber feet prevents the bin from sliding
- Wheels fitted to rear of bin assist moving when cleaning

### Dimensions (w x d x h)

• 42 Litre - 460 x 420 x 580mm

• 70 Litre - 460 x 450 x 790mm

### **Standard Colours**

- · Body White
- · Lid White

### Options

- Alternative Lid Colours Black/Orange/Yellow/Yellow & Black
- Moulded In Graphics Available on request



Alternative lid colours available. Change suffix 'WH' on the item no. with any of the following, i.e. Removable Body Sackholder - Plastic, Silent Closing, 42 Litre with a black lid - 5RSH42BL



| Item     | Description   | Weight |
|----------|---|--------|
| 5RSH42WH | Removable Body Sackholder - Plastic, Silent Closing, 42 Litre | 8.0kg  |
| 5RSH70WH | Removable Body Sackholder - Plastic, Silent Closing, 70 Litre | 9.0kg  |

sales@bristolmaid.com



# **Seating Solutions**

- Seat, Back Cushion & Headrest Options 356
  - Overnight Chair 358
    - Soft Seating 360
  - Visitor & Waiting Room Seating 362
    - Waiting Room Seating 364
      - Tech Chairs & Stools 370
    - PU Tech Chairs & Stools 372
  - Static Safe Tech Chairs & Stools 374
- Static Safe & Sterile Tech Chairs & Stools 376
- Medical Grade Clean Room Tech Chairs 378

- Office, Task & Operators Chairs 330
- 24 Hour Task & Operator Chairs 332
  - Bariatric Chairs 334
  - Saddle & Perch Stools 336
  - Low Back Patient Chairs 338
- Medium & High Back Patient Armchairs 339
  - Reclining Armchairs 340
  - Reclining Armchairs, Bariatric 350
  - Reclining Porters Armchairs 352
  - Rise & Recline Porters Armchairs 354





### **Office, Task & Operators Chairs - Century**

### **Features**

- Available as...
  - · Medium back
  - High back
  - Draughtsman
- $\bullet$  Medium  $\overset{-}{\&}$  high back available with fixed ring arms & non-braking castors
- Draughtsman has a medium back, foot ring & glides on a 27" base
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Two lever mechanism
- · Black polypropylene five star base
- · Pneumatic seat & back rest height adjustment
- Fully upholstered seat & back with deep CMHR foam
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Medium Back Chair in Aztec - 5160/FAC

### Vinyl - Ambla New Richmond in Colour



### Fabric - Camira Advantage





| ltem  | Description                     | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|-------|---------------------------------|--------------------------------|------------------------|--------|
| 5160  | Century, Medium Back            | 670 x 670 x 800-1040mm         | 440 - 570mm            | 10.0kg |
| 5161  | Century, Medium Back, Ring Arms | 670 x 670 x 800-1040mm         | 440 - 570mm            | 12.0kg |
| 5162  | Century, High Back              | 670 x 670 x 940-1070mm         | 440 - 570mm            | 11.0kg |
| 5163  | Century, High Back, Ring Arms   | 670 x 670 x 980-1110mm         | 440 - 570mm            | 13.0kg |
| 5160D | Century, Draughtsman            | 460 x 605 x 1010-1200mm        | 536 - 686mm            | 12.0kg |

FIRE



### 24 Hour Task & Operator Chairs - Eon

### **Features**

- BS 5459: Part 2: 2000. Specification for performance requirements & tests for office furniture. Office pedestal seating for use by persons weighing up to 150kg & for use up to 24 hours a day
   S1 level tested (severe) for heavy duty use
- Certified for use up to 150kg
- Available as...
  - Medium back
  - High back
- 24/7 Flax fire retardant heavy duty fabric upholstery
- Black polypropylene five-star base
- Easy lift back, ratchet height adjustable
- Independent seat & back mechanism with tension control
- Fully upholstered seat & back with deep CMHR foam
- Optional
  - Height adjustable, fold back or ring arms
  - 12" or 14" draughtsman's kit
- Maximum capacity 150kg

Replace suffix on the product code with any of the following to select your colour i.e. Eon, Medium Back in Zone - 5260/FYR09

### Fabric - Camira 24/7 Flax



| Item | Description      | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|------|------------------|--------------------------------|------------------------|--------|
| 5260 | Eon, Medium Back | 710 x 710 x 835-925mm          | 470 - 600mm            | 13.0kg |
| 5262 | Eon, High Back   | 710 x 710 x 925-1030mm         | 470 - 600mm            | 14.0kg |

### 24 Hour Task & Operator Chairs - Zenith

### **Features**

- BS 5459: Part 2: 2000. Specification for performance requirements & tests for office furniture. Office pedestal seating for use by persons weighing up to 150kg & for use up to 24 hours a day
- S1 level tested (severe) for heavy duty use
- Certified for use up to 150kg
- 24/7 Flax fire retardant heavy duty fabric upholstery
- Black polypropylene five-star base
- Easy lift back, ratchet height adjustable
- · Independent seat & back mechanism with tension control & seat slider
- Fully upholstered seat & back with deep CMHR foam
- Chrome swivel height adjustable arms & seat slider
- Maximum capacity 150kg

### **Options**

- 5/Z/HR Headrest
- 5/Z/PLS Pump-Up Lumbar Support

Replace suffix on the product code with any of the following to select your colour i.e. Zenith, High Back in Vernal - 59003/FYR02

### Fabric - Camira 24/7 Flax





59003/FYR02/5/Z/HR/5/Z/PLS shown with optional headrest & pump-up lumbar support

|       | Description       | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|-------|-------------------|--------------------------------|------------------------|--------|
| 59003 | Zenith, High Back | 710 x 710 x 1160-1240mm        | 450 - 530mm            | 21.0kg |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com







### **Bariatric High Back Chairs**

### **Features**

- Deep CMHR foam
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- · Black heavy duty steel five star base
- · Black heavy duty gas lift
- Black height adjustable soft touch armrests
- Easy lift ratchet height adjustable back
- · Independent seat & back mechanism with side tension control
- Maximum capacity 315kg

Replace suffix on the product code with any of the following to select your colour i.e. Bariatric High Back in Lime - 5660XL/VLI

### Vinyl - Ambla New Richmond in Colour



FBK

ItemDescriptionOverall Dimensions (w x d x h)Seat Height AdjustmentWeight5660XLHigh Back, Swivel750 x 750 x 1080-1190mm425 - 535mm29.0kg

FCT

FGE

### **Bariatric Chairs**

### **Features**

- Deep CMHR foam
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- · Fabric Fire retardant
- Beech, varnished four leg frame
- Maximum capacity 315kg

Alternative colours are available at extra cost. Replace suffix on the product code with any of the following to select your colour i.e. Bariatric, Beech Frame in Lime - 5671/VLI

### Vinyl - Ambla New Richmond in Colour







| 1 | FIRE<br>RETARDANT |
|---|-------------------|
|   | <u> </u>          |
|   |                   |
|   |                   |

| Item | Description            | Overall Dimensions (w x d x h) | Seat Height | Weight |
|------|------------------------|--------------------------------|-------------|--------|
| 5671 | Bariatric, Beech Frame | 844 x 580 x 890mm              | 490mm       | 16.0kg |



### Saddle Stools

### **Features**

- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Black polypropylene five star base
- Black gas lift with pneumatic seat height adjustment
- Free wheeling castors
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Saddle Stool in Lime - 5202/VLI

Vinyl - Ambla New Richmond in Colour



### Fabric - Camira Advantage



### **PU Saddle & Perch Stools**

### Features

- Soft touch black polyurethane seatBlack polypropylene five star base
- Black gas lift with pneumatic seat height adjustment
  Free wheeling castors
- Maximum capacity 120kg



|      | Description     | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|------|-----------------|--------------------------------|------------------------|--------|
| 5201 | PU Saddle Stool | 670 x 670 x 590-780mm          | 590 - 780mm            | 8.0kg  |
| 5205 | PU Perch Stool  | 670 x 670 x 650-730mm          | 510 - 690mm            | 8.0kg  |





### Low Back Patient Chairs - Alderbury

### **Features**

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- Ergonomic design to ensure correct sit & posture
- · Handgrips assist positioning from standing & seated positions
- Solid beech frame (PEFC accredited)
- · Mortice & tenon joint construction, glued & pinned
- Frame finished with anti-bacterial varnish
- · Plastic glides fitted to chair legs
- Contoured seat & back, seat section fixed
- Hygiene gap prevents the trapping of dirt & bacteria between the seat & back, provides access for cleaning
- Fire retardant foam (BS 5852: Part 2: 1982)
- Upholstered with fire retardant & antimicrobial vinyl
- Maximum capacity 160kg

### **Standard Colours**

Vinyl - Delft

### Dimensions (w x d x h)

- Overall 580 x 540 x 890mm
- Seat 470 x 430 at Height 470mm
- Width Between Arms 470mm
- Arm Height 655mm

Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Alderbury, Low Back Chair in Baltic - 5ND/LC/ALDE/VVTUR

Vinyl - Spradling Valencia C5



 Item
 Description
 Weight

 5ND/LC/ALDE/VVDEL
 Alderbury, Low Back Chair, Fixed Seat
 8.2kg

sales@bristolmaid.com

www.bristolmaid.com

# Medium & High Back Patient Armchairs - Alderbury

### **Features**

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- · Ergonomic design to ensure correct sit & posture
- Handgrips assist positioning from standing & seated positions
- · Solid beech frame (PEFC accredited)
- Mortice & tenon joint construction, glued & pinned
- · Frame finished with anti-bacterial varnish
- Plastic glides fitted to chair legs
- · Contoured seat & back
  - Fixed Seat is fixed to the frame
  - Removable Seat can be removed to assist cleaning (smooth vinyl upholstered underside) & is reversible to allow the seat section to be repositioned to provide additional clearance at the front of the chair
- Hygiene gap prevents the trapping of dirt & bacteria between the seat & back, provides access for cleaning
- Fire retardant foam (BS 5852: Part 2: 1982)
- · Upholstered with fire retardant & antimicrobial vinyl
- · Maximum capacity 160kg

### Dimensions (w x d x h)

- Medium Back 660 x 650 x 1000mm
- High Back 660 x 650 x 1195mm
- Seat 450 x 435/475\* at Height 500mm
- Width Between Arms 530mm
- Arm Height 655mm

\*Seat depth can be altered by repositioning the removable seat pad

### **Standard Colours**

• Vinyl - Delft

### Options

• HW - Housekeeping Wheels



Optional house keeping wheels assist movement of the chair during cleaning



5ND/HAR/ALDE/VVDEL

www.bristolmaid.com

Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Alderbury, High Back Armchair, Removable & Reversible Seat in Olive - 5ND/HAR/ALDE/VVTUR

Vinyl - Spradling Valencia C5



| Item               | Description                                     |        |
|--------------------|---|--------|
| Medium Back        |   |        |
| 5ND/MAF/ALDE/VVDEL | Alderbury, Medium Back Armchair, Fixed Seat     | 13.6kg |
| 5ND/MAR/ALDE/VVDEL | Alderbury, Medium Back Armchair, Removable Seat | 13.6kg |
| High Back          |   |        |
| 5ND/HAF/ALDE/VVDEL | Alderbury, High Back Armchair, Fixed Seat       | 14.4kg |
| 5ND/HAR/ALDE/VVDEL | Alderbury, High Back Armchair, Removable Seat   | 14.4kg |





### **Reclining Armchair - Shapwick**

### Features

- Available as either manual or electric operation
  - Manual Tilt in space recliner
  - Electric Single tilt in space reclining only or single tilt in space reclining with rise function
- · Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
  - Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting (fixed seat)
- Foam leg rest
- 4 x 50mm swivel braking castors
- Maximum capacity 127kg

### Dimensions (w x d x h)

- Overall Dimensions 737 x 902 x 1130mm
- Seat 510 x 510mm at Height 460mm

### **Standard Colours**

- Vinyl Delft
- VP Fabric Navy Blue

### **Options**

- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements

# ItemDescriptionWeight5RE/RAM/SHAP/VVDELReclining Armchair, Manual41.0kg5RE/RAE/SHAP/VVDELReclining Armchair, Electric55.0kg5RE/RRAE/SHAP/VVDELRise & Reclining Armchair, Electric57.0kg



Alternative colours available at additional cost.

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Shapwick, Reclining Armchair, Electric in Olive Vinyl & Green Vapour Permeable (VP) Multi-Way Stretch Fabric - 5RE/RAE/SHAP/VVOLI/VPGR



Fabric - VP Multi-Way Stretch



| 5RE/RRA | AE/SHA | P/VV | DEL |
|---------|--------|------|-----|
|         |        |      |     |

| Item                     | Description                         |
|--------------------------|-------------------------------------|
| 5RE/RAM/SHAP/VVDEL/VPNB  | Reclining Armchair, Manual          |
| 5RE/RAE/SHAP/VVDEL/VPNB  | Reclining Armchair, Electric        |
| 5RE/RRAE/SHAP/VVDEL/VPNB | Rise & Reclining Armchair, Electric |

www.bristolmaid.com

sales@bristolmaid.com

41.0kg 55.0kg 57.0kg





5RE/RAE/ALWE/VVDEL (fitted with optional pocket)

| Item                | Description                                   | Weight |
|---------------------|---|--------|
| 5RE/RAM/ALWE/VVDEL  | Reclining Armchair, Wide Arm, Manual          | 41.0kg |
| 5RE/RAE/ALWE/VVDEL  | Reclining Armchair, Wide Arm, Electric        | 55.0kg |
| 5RE/RRAE/ALWE/VVDEL | Rise & Reclining Armchair, Wide Arm, Electric | 57.0kg |

# **Reclining Armchair - Alweston, Wide Arm**

### **Features**

- Ideal for use in oncology, dialysis or haematology for intravenous cannulation, where a patient may be seated for a prolonged period of time, or alternatively in recovery following minor treatment
- Wide supportive cushioned arm (150mm)
- Available as either manual or electric operation
  - · Manual Tilt in space recliner
  - Electric Single tilt in space reclining only, or single tilt in space reclining with rise function
- · Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- · Upholstered with either
  - Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting (fixed seat)
- · Foam leg rest
- 4 x 50mm swivel braking castors
- Maximum capacity 127kg

### Dimensions (w x d x h)

### **Standard Colours**

- Overall Dimensions 737 x 902 x 1130mm
  Seat 510 x 510mm at Height 460mm
- Vinyl Delft
- VP Fabric Navy Blue

### **Options**

- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Table with Wireless Charging See page 327
- · Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements
- · USB Charging Socket Located on outside of arm
- · Pocket Positioned on the right-hand side of the chair



Delivery to Chosen Location
Unwrapping & Removal of Packaging
Product Positioning & Set-Up
Full Demonstration of Use to Staff

5RE/RAE/ALWE/VVLEI/VPBE (fitted with optional pocket)

Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VVNB' i.e. Ashington, Reclining Armchair, Waterfall Back, Wide Arm, Electric in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RAE/ASHI/VVOLI/VPGR

VINII - Spradling Valencia C5

Fabric - VP Multi-Way Stretch



|                          | 5RE/RAE/ALWE/VVBAL/VPGY                       |        |
|--------------------------|---|--------|
| Item                     | Description                                   | Weight |
| 5RE/RAM/ALWE/VVDEL/VPNB  | Reclining Armchair, Wide Arm, Manual          | 41.0kg |
| 5RE/RAE/ALWE/VVDEL/VPNB  | Reclining Armchair, Wide Arm, Electric        | 55.0kg |
| 5RE/RRAE/ALWE/VVDEL/VPNB | Rise & Reclining Armchair, Wide Arm, Electric | 57.0kg |





| Item                | Description                                       | Weight |
|---------------------|---|--------|
| 5RE/RAM/TWYF/VVDEL  | Reclining Armchair, Waterfall Back, Manual        | 43.0kg |
| 5RE/RAE/TWYF/VVDEL  | Reclining Armchair, Waterfall Back, Electric      | 57.0kg |
| 5RE/RRAE/TWYF/VVDEL | Rise & Recline Armchair, Waterfall Back, Electric | 59.0kg |

### **Reclining Armchair - Twyford, Waterfall Back**

### **Features**

- Available as either manual or electric operation
  - Manual Tilt in space recliner
  - Electric Single tilt in space reclining only or single tilt in space reclining with rise function
- Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V  $\sim,$  50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
  - · Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- · Easiform reflex foam provides comfort when sitting (fixed seat)
- · Foam leg rest
- Fixed deep padded waterfall back
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

### Dimensions (w x d x h)

- Overall Dimensions 765 x 865 x 1145mm
- Seat 520 x 510mm at Height 482mm

### **Options**

- Seat Pressure Management Coolform, Viscoform, Liquiform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements

### **Standard Colours**

- Vinyl Delft
- VP Fabric Navy Blue



- Delivery to Chosen Location
- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff



Alternative colours available at additional cost.

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Twyford, Reclining Armchair, Electric in Olive Vinyl & Green Vapour Permeable (VP) Multi-Way Stretch Fabric - 5RE/RAE/TWYF/VVOLI/VPGR





VPBE VPNB VPGY VPBL VPBU VPGR

Fabric - VP Multi-Way Stretch

| Item                     | Description                                       | Weight |
|--------------------------|---|--------|
| 5RE/RAM/TWYF/VVDEL/VPNB  | Reclining Armchair, Waterfall Back, Manual        | 43.0kg |
| 5RE/RAE/TWYF/VVDEL/VPNB  | Reclining Armchair, Waterfall Back, Electric      | 57.0kg |
| 5RE/RRAE/TWYF/VVDEL/VPNB | Rise & Recline Armchair, Waterfall Back, Electric | 59.0kg |

sales@bristolmaid.com

+44 (0)1258 484455

5RE/RRAE/TWYF/VVLIC/VPBU



- Ideal for use in oncology, dialysis or haematology for intravenous cannulation, where a patient may be seated for a prolonged period of time, or alternatively in recovery following minor treatment
- Wide supportive cushioned arm (150mm)
- · Available as either manual or electric operation
  - Manual Tilt in space recliner
  - Electric Single tilt in space reclining only, or single tilt in space reclining with rise function
- · Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- · Upholstered with either
  - · Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- · Easiform reflex foam provides comfort when sitting (fixed seat)
- · Foam leg rest
- Fixed deep padded waterfall back
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

### Dimensions (w x d x h)

- Overall Dimensions 765 x 865 x 1145mm
- Seat 520 x 510mm at Height 482mm

### **Standard Colours**

- Vinyl Delft
- VP Fabric Navy Blue

### **Options**

- Seat Pressure Management Coolform, Viscoform, Liquiform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Table with Wireless Charging See page 327
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements
- USB Charging Socket Located on outside of arm
- · Pocket Positioned on the right-hand side of the chair



5RE/RAM/ASHI/VVDELReclining Armchair, Waterfall Back, Wide Arm, Manual43.0kg5RE/RAE/ASHI/VVDELReclining Armchair, Waterfall Back, Wide Arm, Electric57.0kg5RE/RRAE/ASHI/VVDELRise & Recline Armchair, Waterfall Back, Wide Arm, Electric59.0kg

sales@bristolmaid.com



- Delivery to Chosen Location
  Unwrapping & Removal of Packaging Product Positioning & Set-Up
- Full Demonstration of Use to Staff

### 5RE/RAE/ASHI/VVSCH/VPGY (fitted with optional pocket)

**Seating Solutions** 

5RE/RAE/ASHI/VVROT/VPNB

Alternative colours available at additional cost.

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Alweston, Reclining Armchair, Wide Arm, Electric in Olive Vinyl & Green Vapour Permeable VP Multi-Way Stretch Fabric - 5RE/RAE/ALWE/VVOLI/VPGR



Fabric - VP Multi-Way Stretch



| Item                     | Description   | Weight |
|--------------------------|---|--------|
| 5RE/RAM/ASHI/VVDEL/VPNB  | Reclining Armchair, Waterfall Back, Wide Arm, Manual        | 43.0kg |
| 5RE/RAE/ASHI/VVDEL/VPNB  | Reclining Armchair, Waterfall Back, Wide Arm, Electric      | 57.0kg |
| 5RE/RRAE/ASHI/VVDEL/VPNB | Rise & Recline Armchair, Waterfall Back, Wide Arm, Electric | 59.0kg |

sales@bristolmaid.com





sales@bristolmaid.com

### **Reclining Armchair - Beckington, Waterfall Back**

### **Features**

- Electric reclining only, or dual tilt in space reclining with rise function
- Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Upholstered with flame retardant vinyl
- · Easiform reflex foam provides comfort when sitting
- · Foam leg rest
- · Fixed waterfall back
- · Side pocket & loop for handset
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

### Dimensions (w x d x h)

- Overall Dimensions 840 x 930 x 1160mm
- Seat 510 x 510mm at Height 483mm
- Back Height 710mm
- Arm Height 180mm

### Standard Colours

· Vinyl - Ultra



- Delivery to Chosen Location
- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- · Full Demonstration of Use to Staff

Alternative colours available at additional cost. Change suffix 'VMULT' on the product code with any of the following i.e. Beckington, Reclining Armchair, Waterfall Back, Electric in Garden - 5RE/RAE/BECK/VMGAR

### **Memphis Faux Leather**



| Item                | Description                                       | Weight |
|---------------------|---|--------|
| 5RE/RAE/BECK/VMULT  | Reclining Armchair, Waterfall Back, Electric      | 58.0kg |
| 5RE/RRAE/BECK/VMULT | Rise & Recline Armchair, Waterfall Back, Electric | 62.0kg |

+44 (0)1258 484455

www.bristolmaid.com





### Reclining Armchair - Hadleigh, Waterfall Back, Bariatric





### **Features**

- Specifically designed for comfort & pressure management control
- Dual standard rise & recliner, provides a full range of posture positive recline positions as well as support to a standing position
- Mains powered actuator (230V ~. 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Upholstered with a combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting
- · Seat & leg sections generously padded with memory foam & fitted with a removable cover
- · Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- Back blockers provided for additional lumbar support
- Flat wide arms for deep comfort & support
- · Loop for handset & deep side pocket
- Heavy duty glides
- Maximum capacity 254kg

### Dimensions (w x d x h)

- Chair 1000 / 1100 x 940 x 1145mm
- Seat 610 / 710mm x 510 x 460mm

### Standard Colours

- · Vinvl Delft
- VP Fabric Navy Blue

### **Options**

- Seat Pressure Management Coolform, Viscoform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- · USB Charging Socket Located on outside of arm
- · Full Length Arm Caps Positioned over the arm rests

5RE/RRAE/HADL610/VVDEL/VPNB Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 610mm 88.0kg 5RE/RRAE/HADL710/VVDEL/VPNB Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 710mm 88.0kg



- Delivery to Chosen LocationUnwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff

Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Hadleigh, Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 710mm in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RRAE/HADL710/VVOLI/VPGR



Fabric - VP Multi-Way Stretch



SRE/RAE/HADL610/VVROT/VPBU

### **Reclining Porters Armchair - Bathampton, Waterfall Back**





5RE/RPAH/BATH/VVDEL Reclining Porters Armchair, Waterfall Back, Electric 75.0kg

### **Features**

- Single tilt in space recliner, provides a full range of posture positive recline positions
- Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset (positioned at rear of chair/optional side of chair)
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
  - · Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multi-way stretch fabric. VP fabric applied to seat & back sections
- Upholstered seat & back sections are fitted with removable covers, allows seat section to be easily upgraded to aid pressure care
- · Removable internal Easiform reflex seat cushion
- Seat & leg sections generously padded with memory foam & fitted with a removable cover
- Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- Push handle at rear of the chair to assist when positioning/portering
- Slide out footplate to support feet
- 4 x 75mm swivel braking castors
- Maximum capacity 127kg

### Dimensions (w x d x h)

- Overall Dimensions 685 x 990 x 1220mm
- Seat 520 x 510mm at Height 470mm

### **Standard Colours**

- Vinyl Delft
- VP Fabric Navy Blue

### Options

- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform, Airform, Dynaform
- Back Cushions Comfort Lateral, Support Lateral, Large Profile with Waterfall, Large Profile Headrest with Comfort Lateral, Large Profile Headrest with Support Lateral
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements



Delivery to Chosen LocationUnwrapping & Removal of Packaging Product Positioning & Set-Up Full Demonstration of Use to Staff

5RE/RPAH/BATH/VVOLI/VPGR

Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Bathampton, Waterfall Back, Handset Operated in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RPAH/BATH/VVOLI/VPGR





|                          | Description  | Weight |
|--------------------------|--|--------|
| 5RE/RPAH/BATH/VVDEL/VPNB | Reclining Porters Armchair, Waterfall Back, Electric | 75.0kg |

VPGY

VPGR





### **Features**

- Dual tilt in space rise & recliner, provides a full range of posture positive recline positions as well as support to a standing position
- Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset (positioned at side of chair)
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- · Upholstered with either
  - Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinvl & vapour permeable (VP) multi-way stretch fabric. VP fabric applied to seat & back sections
- Upholstered seat & back sections are fitted with removable covers, allows seat section to be easily upgraded to aid pressure care
- Removable internal Easiform reflex seat cushion
- · Seat & leg sections generously padded with memory foam & fitted with a removable cover
- · Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- · Push handle at rear of the chair to assist when positioning/portering
- · Optional sliding footplate to support feet
- 4 x 100mm swivel braking castors
- Maximum capacity 127kg

### Dimensions (w x d x h)

- Overall Dimensions 775 x 1055 x 1245mm
- Seat 510 x 510mm at Height 480mm

### Standard Colours

- Vinyl Delft
- VP Fabric Navy Blue

### **Options**

- Sliding Footplate
- Seat Pressure Management Coolform, Viscoform, Liquiform. Celliform, Airform, Dynaform
- Back Cushions Comfort Lateral, Support Lateral, Large Profile with Waterfall, Large Profile Headrest with Comfort Lateral, Large Profile Headrest with Support Lateral
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile. Midline
- · Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements

5RE/RRPAH/DORC/VVDEL Rise & Recline Porters Armchair, Waterfall Back, Electric 76.5kg



Delivery to Chosen Location
Unwrapping & Removal of Packaging
Product Positioning & Set-Up
Full Demonstration of Use to Staff

5RE/RRPAH/DORC/VVBAL/VPNB (with optional footplate)

11

Alternative colours available at additional cost.

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Dorchester, Waterfall Back, Handset Operated in Leinen Vinyl & Grey VP Multi-Way Stretch Fabric - 5RE/RPAH/DORC/VVLEI/VPGY



Fabric - VP Multi-Way Stretch



5RE/RRPAH/DORC/VVLEI/VPNB (with optional sliding footplate)


### Description W

5RE/RRPAH/DORC/VVDEL/VPNB Rise & Recline Porters Armchair, Waterfall Back, Electric 76.5kg

### Seat Options For Pressure Management

Cushions are covered in multi-way stretch, vapour permeable fabric which is wipeable, removable, washable at 95°C, incontinence proof & breathable.

supported comfort.



for periods longer than one hour, require pressure relief &


High density foam provides support, whilst a top layer of

memory foam allows submersion into the cushion, providing

a large surface area to reduce pressure levels.


For those sitting for long periods with specific comfort issues, such as a high risk of pressure sores. Uses the latest medical grade silicone gel technology allowing pressure to be evenly distributed

The gel conforms to the shape of the body, equalises to body temperature, allows air to flow easily & is totally odourless & hypoallergenic.



For those sitting for long periods of time with specific comfort issues, such as a high risk of pressure sores & sheer forces.

Gelacell's unique grid design eliminates pressure points. contouring in response to the transfer of weight, evenly distributes pressure & reduces heat absorbtion, keeping the seat cooler & fresher



Airform



For those who have good posture, but are unable to alter their position easily & require high level pressure care.

This cushion provides the benefits of an air cushion, however without the sensation of alternating changes in air pressure.

Suitable for users 40 - 127kg.

**Dynaform** 

constantly adapting to movement.

Pressure Risk VERY HIGH

PUMP

For those at very high risk of developing pressure sores, or have a history of sores, where the priority is to prevent further pressure damage.

for long periods of time, require pressure relief & supported

High density foam provides support, topped with layers

of fluid liquid gel & memory foam to provide high levels of

pressure relief. The liquid gel helps to dissipate heat, whilst

comfort. Excellent for people with little fatty tissue.

Dynaform is a dynamic alternating air cushion, which constantly changes air levels within its cells to recreate the effect of pressure relieving & prevent pressure build up.

Air cell features

Alternating

air cushion

Reflex foam

- 6 dynamic cells 10cm wide
- 3 cell cycle repeated every 9 minutes

Suitable for users 40 - 121kg.

Please note the Dynaform cushion does need an electric supply at all times to operate.

Pressure risk levels detailed are for guidance only. Each individual user has different posture, medical needs & health issues which contribute to the performance of any pressure cushion. If in doubt, please consult a GP or Health Care Professional.

### **Back Cushion Options**



Waterfall Back Cushion Individual waterfall cushions, fully adjustable to accommodate maximum postural support



### Large Profile Headrest with Waterfall Back

For those with limited head control & requiring a soft adjustable lumbar support

### **Headrest Options**



**Head Pillow** 




**Comfort Lateral Back Cushion** Soft fibre filled removable lateral cushion for lower back & side lumber support, without feeling restricted



Large Profile Headrest with Comfort Lateral Back

For those with limited head control & requiring soft lateral support without feeling restricted

**Head Roll** 



Support Lateral Back Cushion A foam filled back providing firmer lateral support with a soft fibre filled headrest



Large profile headrest with Support Lateral Back For those with limited head control. requiring firm foam lateral lumbar & midline support

**Table with Wireless Charging** 360-degree rotation & angle adjustable table with a wipe clean surface. Supports QI wireless charging & fitted with a USB socket. Removable magnetic bar provides support & ensures correct positioning on the wireless receiver coil,

wipe clean surface (w 345 x h 245mm)





**Small Profile Headrest** 



Large Profile Headrest



**Midline Headrest** 

sales@bristolmaid.com

www.bristolmaid.com





### **Overnight Chair - Portland, Three Position**

### **Features**

- Multi-functional compact chair that provides comfortable seating by day, transforming into a bed by night
- Suitable for use in maternity units, children's wards, relative's rooms, hospices & care homes
- · Provides visitors with a comfortable place to sleep during longer visits
- · Chair is pulled out in three motions to a bed
  - Chair
  - Lounge Chair
  - Bed
- Upholstered using flame retardant antimicrobial vinyl
- Easiform reflex foam provides comfort when sitting or sleeping
- 50mm swivel braking castors
- Maximum capacity 160kg

### Dimensions (w x d x h)

- Overall Dimensions 870 x 840 x 870mm (when open overall bed length is 2060mm)
- Seat 690 x 540mm at Height 410mm

### **Standard Colours**

- Vinyl Delft
- · Vinyl & Patterned Fabric Vinyl Lavendel & Savanna Lavender



- Delivery to Chosen Location
- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff

Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Portland, Three Position in Olive - 5RE/OC/PORT/VVOLI

Vinyl - Spradling Valencia C5



ItemDescriptionWeight5RE/OC/PORT/VVDELThree Position, Vinyl Delft61.5kg5RE/OC/PORT/VVLAVSAVLThree Position, Vinyl Lavendel & Savanna Lavender61.5kg

sales@bristolmaid.com

+44 (0)1258 484455

www.bristolmaid.com







# Soft Seating - Zeta, Modular

#### **Features**

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- Zeta is a stylish, fully modular range of soft seating based around eight standard elements suitable for use in seating environments, ranging from reception areas, breakout spaces to waiting areas & cafeterias
- FSC® certified wood core, constructed from 18mm boards
- Traditional webbed seat for additional comfort & durability
- High density CMHR foam
- · Vinyl or fabric upholstery
  - · Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Chrome finish legs
- Arms available on single & double units with backs
- Backs are available with a raised feature to allow fluid run off (overall height 810mm)
- · Maximum capacity
  - Single 146kg

#### **Options**

- RA Right Arm
- LA Left Arm
- 2A Left & Right Arms

Replace suffix on the product code with any of the following to select your colour i.e. Single Unit in Cobalt - 5ZETA1/FCT

#### Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage



| Item        | Description                               | Overall Dimensions<br>(w x d x h) | Seat Height | Seat Depth | Weight |
|-------------|---|-----------------------------------|-------------|------------|--------|
| 5ZETA1      | Single Unit                               | 650 x 650 x 415mm                 | 415mm       | 650mm      | 12.0kg |
| 5ZETA2      | Double Unit                               | 1300 x 650 x 415mm                | 415mm       | 650mm      | 22.0kg |
| 5ZETA1/B    | Single Unit with Back                     | 650 x 650 x 765mm                 | 415mm       | 485mm      | 16.5kg |
| 5ZETA2/2B   | Double Unit with Back                     | 1300 x 650 x 765mm                | 415mm       | 485mm      | 31.0kg |
| 5ZETA1/CO   | Single Unit with Corner                   | 650 x 650 x 765mm                 | 415mm       | 485mm      | 16.5kg |
| 5ZETA2/B/RC | Double Unit with Right Corner & Left Back | 1300 x 650 x 765mm                | 415mm       | 485mm      | 35.5kg |
| 5ZETA2/B/LC | Double Unit with Left Corner & Right Back | 1300 x 650 x 765mm                | 415mm       | 485mm      | 35.5kg |
| 5ZETAQ      | Quarter Segment                           | 650 x 650 x 415mm                 | 415mm       | 650mm      | 12.0kg |
| 5ZETAH      | Half Segment                              | 1300 x 650 x 415mm                | 415mm       | 650mm      | 23.0kg |
| 5ZETAH/B    | Half Segment with Back                    | 1300 x 650 x 765mm                | 415mm       | 485mm      | 25.0kg |
| 5ZETAC45    | 45 Degree Curve                           | 950 x 650 x 415mm                 | 415mm       | 650mm      | 11.5kg |
| 5ZETAC45/B  | 45 Degree Curve with Back                 | 950 x 650 x 765mm                 | 415mm       | 485mm      | 16.5kg |
| 5ZETAC90    | 90 Degree Curve                           | 1625 x 650 x 415mm                | 415mm       | 650mm      | 22.0kg |
| 5ZETAC90/B  | 90 Degree Curve with Back                 | 1625 x 650 x 765mm                | 415mm       | 485mm      | 32.0kg |
| 5ZETA930C   | Single Corner Unit                        | 930 x 650 x 765mm                 | 415mm       | 650mm      | 28.0kg |

sales@bristolmaid.com

+44 (0)1258 484455

```
www.bristolmaid.com
```





# Visitor & Waiting Room Seating - Mata Chairs

#### **Features**

- BS EN 1729-1: Part 1 & 2: 2015 (functional dimensions & safety requirements)
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Available with...
  - Four legs
  - Sled base
- Moulded polypropylene shell available in seven colours
- Silver finish four leg
- Chrome finish sled base
- Four leg chair stackable five high
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Four Leg Chair in Olive - 5SM6L/SOE

Plastic Shell



sales@bristolmaid.com +44

+44 (0)1258 484455

# **Visitor & Waiting Room Seating - Prima Chairs**

#### **Features**

- BS EN 16139: 2013 Level 2 (strength, durability & safety)
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Certified for use up to 203kg
- Choice of vinyl, fabric upholstery or polypropylene seat & back
  - Vinyl Fire retardant & antimicrobial
  - · Fabric Fire retardant
- Stable four leg design
- · Black finish tubular frame
- Stacks up to ten high (transporter available)
- Built-in linking system (optional extension links)
- Maximum capacity 203kg

Replace suffix on the product code with any of the following to select your colour i.e. Four Leg Chair in Cobalt - 5301/FCT

#### Vinyl - Ambla New Richmond in Colour



Fabric - Camira Advantage



FKR SBK



FAC

FRD

5301/A/FRD

5301/VLI

| ltem   | Description         | Overall Dimensions (w x d x h) | Seat Height | Weight |
|--------|---------------------|--------------------------------|-------------|--------|
| 5301   | Four Leg            | 520 x 580 x 800mm              | 485mm       | 8.0kg  |
| 5301/A | Four Leg, Arm Rests | 620 x 580 x 800mm              | 485mm       | 10.0kg |





| ltem    | Description             | Overall Dimensions (w x d x h) | Seat Height | Weight |
|---------|-------------------------|--------------------------------|-------------|--------|
| 5NOVA/1 | Nova - Tub, Single Seat | 680 x 680 x 800mm              | 450mm       | 16.0kg |
| 5NOVA/2 | Nova - Tub, Double Seat | 1360 x 680 x 800mm             | 450mm       | 31.0kg |

# Waiting Room Seating - Nova, Tub

#### **Features**

- Nova single & double tub chairs are suitable for a wide range of seating environments, ranging from cafeterias to waiting areas
- FSC<sup>®</sup> certified wood core with moulded deep CMHR seat cushions
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Chrome finish legs
- Maximum capacity
  - Single 146kg
  - Double 292kg

Replace suffix on the product code with any of the following to select your colour i.e. Double Seat in Cobalt - 5NOVA/3/FCT

Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage



# Waiting Room Seating - Horizon, Beech Frame

#### **Features**

- Horizon chairs, sofas & tables are designed to suit a variety of waiting & occasional seating areas
- Beech frame with plywood seat shells
- Arms assist positioning from standing & seated positions
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - · Fabric Fire retardant
- Fixed seat, with a loose back cushion to assist cleaning
- Back cushion has a zipped cover, allowing for easy replacement if damaged
- · Maximum capacity
  - Single chair 146kg
  - Two seater sofa 120kg per seat
  - Three seater sofa 120kg per seat

Replace suffix on the product code with any of the following to select your colour i.e. Two Seater Sofa with Arms in Cobalt - 5HORI/2/A/FCT

Vinyl - Ambla New Richmond in Colour



Fabric - Camira Advantage







| Item      | Description                       | Overall Dimensions (w x d x h) | Seat Height | Weight |
|-----------|-----------------------------------|--------------------------------|-------------|--------|
| 5HORI/1/A | Horizon - Easy Chair, Arms        | 750 x 650 x 800mm              | 450mm       | 19.0kg |
| 5HORI/2/A | Horizon - Two Seater Sofa, Arms   | 1350 x 650 x 800mm             | 450mm       | 33.0kg |
| 5HORI/3/A | Horizon - Three Seater Sofa, Arms | 1950 x 650 x 800mm             | 450mm       | 47.0kg |

+44 (0)1258 484455 www





# Waiting Room Seating - Prima Beam, Fixed

#### **Features**

- Prima beam is suitable for all environments requiring durable & robust waiting room seating solutions
- Certified for use up to 32 stone (203kg) per seat
- Fixed seat
- Choice of vinyl, fabric upholstery or polypropylene seat & back
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Black finish beam & legs
- Drilled for floor fixing (fixings not supplied)
- Maximum capacity 203kg per seat

#### **Options**

• 510065 - Armrest (x1)

Replace suffix on the product code with any of the following to select your colour i.e. Fixed Beam Seat, Two Seat in Cobalt - 5PF332/FCT

#### Vinyl - Ambla New Richmond in Colour





#### **Fixed Seat**

| Item   | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|--------|-------------|--------------------------------|-------------|--------|
| 5PF332 | Two Seat    | 1200 x 560 x 790mm             | 480mm       | 20.0kg |
| 5PF333 | Three Seat  | 1800 x 560 x 790mm             | 480mm       | 25.0kg |
| 5PF334 | Four Seat   | 2400 x 560 x 790mm             | 480mm       | 30.0kg |
| 5PF335 | Five Seat   | 3000 x 560 x 790mm             | 480mm       | 35.0kg |

sales@bristolmaid.com

+44 (0)1258 484455



# **Options**

# Waiting Room Seating - Prima Beam, Tip-Up Seat

#### **Features**

- Prima beam is suitable for all environments requiring durable & robust waiting room seating solutions
- Certified for use up to 32 stone (203kg) per seat
- Tip-up seat (must be fixed to the floor)
- Choice of vinyl, fabric upholstery or polypropylene seat & back
  - · Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Black finish beam & leas
- Drilled for floor fixing (fixings not supplied)
- Maximum capacity 203kg per seat
- 510065 Armrest (x1)

Replace suffix on the product code with any of the following to select your colour i.e. Fixed Beam Seat. Two Seat in Cobalt - 5PF332/FCT

#### Vinyl - Ambla New Richmond in Colour





5PT338/FRD Shown with 510065 (optional arms sold individually) - contact us for details

#### **Tip-Up Seat**

| Item   | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|--------|-------------|--------------------------------|-------------|--------|
| 5PT336 | Two Seat    | 1200 x 490 x 790mm             | 490mm       | 20.0kg |
| 5PT337 | Three Seat  | 1800 x 490 x 790mm             | 490mm       | 25.0kg |
| 5PT338 | Four Seat   | 2400 x 490 x 790mm             | 490mm       | 30.0kg |
| 5PT339 | Five Seat   | 3000 x 490 x 790mm             | 490mm       | 35.0kg |

sales@bristolmaid.com

+44 (0)1258 484455



# **Tech Stools**

#### Features

- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
- Fabric Fire retardant
- Black polypropylene five star base
- Models
  - Medium/High Foot ring & glides
  - Low Non-braking castors
- Black gas lift with pneumatic seat height adjustment
- Maximum load 120kg

#### Options

• 5BLSTC - Brake loaded castors



| Item   | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|-------------|--------------------------------|------------------------|--------|
| 5TC110 | Low         | 670 x 670 x 440-570mm          | 440 - 570mm            | 8.0kg  |
| 5TC109 | Medium      | 670 x 670 x 570-690mm          | 570 - 690mm            | 10.0kg |
| 5TC108 | High        | 670 x 670 x 650-810mm          | 650 - 810mm            | 10.0kg |

Replace suffix on the product code with any of the following to select your colour i.e. Low Stool in Kingfisher - 5TC110/FKR

Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage



sales@bristolmaid.com

+44 (0)1258 484455 www.bris





| Item | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|------|----------------------------|--------------------------------|------------------------|--------|
| 5PU6 | Low, Castors               | 670 x 670 x 750-980mm          | 430 - 570mm            | 13.0kg |
| 5PU5 | Medium, Foot Ring & Glides | 710 x 710 x 830-1120mm         | 510 - 700mm            | 14.0kg |
| 5PU4 | High, Foot Ring & Glides   | 710 x 710 x 910-1270mm         | 580 - 840mm            | 14.0kg |

# **PU Tech Stools**

#### Features

- Soft touch black polyurethane seatBlack polypropylene five star base
- Models





| Item   | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5PU153 | Low, Castors               | 670 x 670 x 410-540mm          | 410 - 540mm            | 8.0kg  |
| 5PU152 | Medium, Foot Ring & Glides | 670 x 670 x 470-660mm          | 470 - 660mm            | 10.0kg |
| 5PU151 | High, Foot Ring & Glides   | 670 x 670 x 530-780mm          | 530 - 780mm            | 10.0kg |



|        | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC203 | Low, Castors               | 670 x 670 x 890-1090mm         | 440 - 580mm            | 13.0kg |
| 5TC202 | Medium, Foot Ring & Glides | 710 x 710 x 950-1210mm         | 520 - 710mm            | 14.0kg |
| 5TC201 | High, Foot Ring & Glides   | 710 x 710 x 980-1320mm         | 570 - 820mm            | 14.0kg |

Seating Solutions

# **Static Safe Tech Stools**

#### Features

- Hi-Tech ESD (anti-static) fire retardant black fabric upholstery
- Aluminium five star base
- Models
  - Medium/High Foot ring & glides
  - Low ESD free wheeling castors

5TC210

- Chrome gas lift with pneumatic seat height adjustment
- Supplied with ESD (anti-static) castors or glides
- Maximum capacity 120kg

#### Options

• Brake loaded ESD castors







670 x 670 x 540-790mm

www.bristolmaid.com

540 - 790mm

| 5 | TC209  |              | 5TC208                         |                        |        |
|---|--------|--------------|--------------------------------|------------------------|--------|
|   | Item   | Description  | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|   | 5TC210 | Low, Castors | 670 x 670 x 386-556mm          | 386 - 556mm            | 10.0kg |
|   |        |              |                                |                        |        |

5TC208 High, Foot Ring & Glides

11.0kg



## Static Safe & Sterile Tech Chairs

- ESD (anti-static) fire retardant & antimicrobial black vinyl upholstery
- Aluminium five star base
  - Medium/High Foot ring & glides
  - Low Castors
- Chrome gas lift with pneumatic seat & back rest height adjustment
- Supplied with ESD (anti-static) castors or glides
- Maximum capacity 120kg





**HEPA** filter

14.0kg

14.0kg

www.bristolmaid.com

520 - 710mm

570 - 8 20mm

5TC304 High, Foot Ring & Glides

5TC305 Medium, Foot Ring & Glides 710 x 710 x 950-1210mm

710 x 710 x 980-1320mm

# Static Safe & Sterile Tech Stools

#### Features

- ESD (anti-static) fire retardant & antimicrobial black vinyl upholstery
- Aluminium five star base
- Models
  - Medium/High Foot ring & glides
  - Low Castors







| Item   | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC310 | Low, Castors               | 670 x 670 x 386-556mm          | 386 - 556mm            | 10.0kg |
| 5TC309 | Medium, Foot Ring & Glides | 670 x 670 x 490-680mm          | 490 - 680mm            | 11.0kg |
| 5TC308 | High, Foot Ring & Glides   | 670 x 670 x 540-790mm          | 540 - 790mm            | 11.0kg |



|        | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC313 | Low, Castors               | 670 x 670 x 890-1090mm         | 440 - 580mm            | 13.0kg |
| 5TC312 | Medium, Foot Ring & Glides | 710 x 710 x 950-1210mm         | 520 - 710mm            | 14.0kg |
| 5TC311 | High, Foot Ring & Glides   | 710 x 710 x 980-1320mm         | 570 - 820mm            | 14.0kg |

# Medical Grade Clean Room Tech Chairs

- Medical grade fire retardant & antimicrobial black vinyl upholstery
- Aluminium five star base
  - Medium/High Foot ring & glides
  - Low Castors
- Chrome gas lift with pneumatic seat & back rest height adjustment
- Maximum capacity 120kg
- Independent seat & back mechanism
- Brake loaded castors



HEPA filter

sales@bristolmaid.com

# Other products that you may be interested in...

See Page 118

# Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep





# **Sepsis Trolleys**

- Sepsis Trolleys 382
  - Accessories 383
- Caretray Trolley Sepsis 386



## **Sepsis Trolleys**

#### **Features**

- Trollevs designed for use when treating Sepsis
- Supports delivery of 'The Sepsis Six' with clearly defined locations for equipment used in each of the stages
- · Door/s fitted with either:
  - Electronic (battery operated) Push Button Lock Customer can set their own separate user, master, submaster & technician access codes
  - High Security Bolt Lock Supplied with two keys
- Trolleys supplied with door shelves (useful storage for antibiotics etc) & an A4 document holder
  - Single Door 2 x door shelves / 1 x A4 document holder
- Double Door 5 x door shelves / 1 x A4 document holder
- Six labelled drawers (Oxygen / Blood Cultures / Intravenous Antibiotics / Intravenous Fluid /Haemoglobin & Lactates / Urine Measurement)
- Optional 'Sharps Box Holder' or 'Universal Sharps Box Bracket Mounting Plate' (fits to the I/h side of the trollev)
- Sharps Box Holder Accepts up to seven litre capacity sharps boxes (customer to fit - fixings supplied)
- Universal Sharps Box Bracket Mounting Plate Compatible with (customer to supply the plastic brackets & fit - fixings supplied)
  - Daniels Healthcare
    - POUDS<sup>®</sup> Wall-Bracket (up to 5L) SHARPSGUARD<sup>®</sup> mini, 1, 2.5. 3.75. 5 & com plus
    - POUDS<sup>®</sup> Wall-Bracket (7L & above) SHARPSGUARD<sup>®</sup> 7. 11.5 & 22
  - Frontier (Sharpsafe) 2, 3, 4 & 7
  - Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS
- · Upper work surface
- · Push handles incorporate integral buffers
- 125mm polymer castors
- . When not in use, trolleys storing medicines should be secured to the wall. See option MT/1/6F Security Clamp

#### Dimensions (w x d x h)

- Single Door Trolley 512 x 542 x 1072mm
- Double Door Trolley 712 x 542 x 1072mm
- Drawers
  - Single Door Trolley 338 x 297 x 115mm (internal)
  - Double Door Trolley 538 x 297 x 115mm (internal)

#### Standard Colours

- · Paint Grev
- Work Surface Orange

Sepsis Trolleys

Weight

53.4kg

69.1kg

# **Sharps Box Holder**

#### **Features**

SEP/SBH

- · Accepts up to seven litre capacity sharps boxes
- Sharps box secured using a rubber band fitted with a toggle
- Attaches to the left hand side of the unit (customer to fit fixings supplied)



#### **Features**

- Compatible with the following brackets (customer to supply the plastic brackets & fit fixings supplied)
   Daniels Healthcare
  - POUDS® Wall-Bracket (up to 5L) SHARPSGUARD® mini, 1, 2.5, 3.75, 5 & com plus
  - POUDS® Wall-Bracket (7L & above) SHARPSGUARD® 7, 11.5 & 22
  - Frontier (Sharpsafe) 2, 3, 4 & 7
  - Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS




sales@bristolmaid.com

+44 (0)1258 484455 wv

# **Security Clamp**

#### **Features**

- · Allows the trolley to be secured to the 'fabric' of the building
- Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door of the trolley when closed and locked
- · Recommended wall fitting height 1000mm



#### **Important Note - Security Clamps**

It is a UK requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area and broken into, but also encourages good nursing practice, giving the trolley a 'permanent home'.

'Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB- March 2005'


| Item    | Description    | Cable Length | Weight |
|---------|----------------|--------------|--------|
| MT/1/6F | Security Clamp | 950mm        | 0.8kg  |

sales@bristolmaid.com

# **Drawer Inserts**

#### **Features**

- Plastic drawer inserts, height 100mm
- Assists segregation of drawer contents
- · Suitable for single & double door trolleys
- Insert supplied with removable & repositionable dividers
  - SEP/DIS 2 x dividers providing 4 separate areas
  - SEP/DID 4 x dividers providing 6 separate areas
- Supplied for self-assembly
- Clip on label holders & plastic labels available to assist identification



| Item    | Description                                | Weight |
|---------|--|--------|
| SEP/DIS | Drawer Insert - Set - Single Door Trolleys | 0.5kg  |
| SEP/DID | Drawer Insert - Set - Double Door Trolleys | 0.6kg  |

# Label Holders & Labels

#### Features

- Plastic label holders & labels
- Assists identification of drawer contents
- Label holders snap onto drawer inserts/dividers
- Accepts ISO labels 40 x 60mm (w x h)
- Sold in packs of ten



| Item   | Description  |
|--------|--------------|
| SEP/LH | Label Holder |
| SEP/L  | Label        |


# Label - Sepsis Tiger

#### **Features**

- Sammy the sepsis tiger label developed for use on Paediatric trolleys
- Compatible with both single & double door trolleys



# Item Description SEP/TIG01 Label - Sammy Sepsis Tiger





# Item Description Weight CTSEPSM1/10/6S/HSB 6 x 100mm Trays, High Security Bolt Lock 36.8kg

sales@bristolmaid.com

# **Caretray Trolley - Sepsis**

#### **Features**

- · Trolleys designed for use in the treatment of Sepsis
- Supports delivery of 'The Sepsis Six' with clearly defined locations for equipment used in each of the stages
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Waste Bin & Bracket, 8 Litre
  - Sharps Box Holder
  - Labels, Sepsis Six
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (sepsis six labels supplied - oxygen / blood cultures / intravenous antibiotics / intravenous fluid / haemoglobin & lactates / urine measurement)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking
- When not in use, trolleys storing medicines should be secured to the wall using a
   Security Clamp (see locking & security accessories page 384)

#### Dimensions (w x d x h)

#### **Standard Colours**

Paint - Orange

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays

+44 (0)1258 484455

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

#### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67

- Carerail Mounted page 68
  Laptop Holders page 74
  Trays, Shelves, Dividers & Labels page 77





| Item               | Description                                  | Weight |
|--------------------|--|--------|
| CTSEPSM1/10/6S/EPB | 6 x 100mm Trays, Electronic Push Button Lock | 36.8kg |



# **Shelving & Racking Solutions**

Engineered for a Demanding Environment 390

- Polymer Racking MetroMax I<sup>®</sup> 392
- Polymer Racking MetroMax Q<sup>®</sup> 394
- Zinc Chrome Racking Olympic<sup>®</sup> 396
- High Density Storage Top-Track<sup>®</sup> 398
- High Density Storage qwikTRAK<sup>®</sup> 399

# **Engineered for a Demanding Environment**



# MetroMax I<sup>®</sup> Full Polymer System with Lifetime Warranty against corrosion





Microban<sup>®</sup> antimicrobial product protection helps keep shelves "cleaner between cleanings" by inhibiting the growth of mold, mildew, fungus and bacteria that cause odors and stains on the shelf surface.



MetroMax Q<sup>®</sup> Polymer / Epoxy Coated System with 15 Year Warranty against corrosion









#### **Stationary Units**

| 4 Shelves - 1590mm High |               | 5 Shelves - 1 | 5 Shelves - 1895mm High |              |
|-------------------------|---------------|---------------|-------------------------|--------------|
| Grid Shelves            | Solid Shelves | Grid Shelves  | Solid Shelves           |              |
| 8MX536SG                | 8MX536SS      | 85MX536SG     | 85MX536SS               | 914 x 610mm  |
| 8MX546SG                | 8MX546SS      | 85MX546SG     | 85MX546SS               | 1066 x 610mm |
| 8MX556SG                | 8MX556SS      | 85MX556SG     | 85MX556SS               | 1219 x 610mm |
| 8MX566SG                | 8MX566SS      | 85MX566SG     | 85MX566SS               | 1524 x 610mm |

#### **Mobile Units**

| 4 Shelves - 1 | 1730mm High   | 5 Shelves - 2 | 2035mm High   | Shelf Size (d x w) |
|---------------|---------------|---------------|---------------|--------------------|
| Grid Shelves  | Solid Shelves | Grid Shelves  | Solid Shelves |                    |
| 8MX536MG      | 8MX536MS      | 85MX536MG     | 85MX536MS     | 914 x 610mm        |
| 8MX546MG      | 8MX546MS      | 85MX546MG     | 85MX546MS     | 1066 x 610mm       |
| 8MX556MG      | 8MX556MS      | 85MX556MG     | 85MX556MS     | 1219 x 610mm       |
| 8MX566MG      | 8MX566MS      | 85MX566MG     | 85MX566MS     | 1524 x 610mm       |

sales@bristolmaid.com

# **Polymer Racking - MetroMax I®**

#### **Features**

- Corrosion proof polymer shelving & racking system
- Polymer shelves with stainless steel corners & a patented no tools assembly, adjustable in 25mm increments
- All polymer posts
- · Removable easy to clean polymer shelf mats (will fit most dishwashers) available as...
  - · Grid Shelves Open grid design allowing for the free circulation of air, improved visibility & light penetration
  - Solid Shelves Solid design offering a tear free surface, ideal for the storage of sterilised DIN/ISO travs within theatres
- Microban<sup>®</sup> product protection inhibits growth of bacteria & aids infection control
- · Available as stationary or mobile units
  - Stationary Fitted with adjustable feet to allow levelling on uneven floors
  - Mobile Fitted with 125mm swivel / braking castors
- Rack / unit capacity...
  - Stationary 907kg
  - Mobile 339kg
- Optional label holders, shelf dividers & ledges
- · Lifetime anti-corrosion warranty
- · Supplied for self-assembly



Solid Shelves

Contact us for more information

+44 (0)1258 484455



# **Label Holders**

#### **Features**


- · Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels

#### Dimensions (w x h)

• 102 x 27mm

#### **Standard Colours**

• Grey



# **Shelf Divider**

#### **Features**

• Corrosion proof polymer shelf divider (height 200mm)

+44 (0)1258 484455

- Use to create organised compartments or simply separate shelf contents
- Easily repositioned on grid & solid shelves
- Easy to clean

#### Dimensions (w x h)

• 610 x 200mm

#### **Standard Colours**

• Grey

sales@bristolmaid.com

| Item    | Description  | Weight | Item     | Description   | Weight |
|---------|--------------|--------|----------|---------------|--------|
| 89989PX | Label Holder | 0.01kg | 8MXD24-8 | Shelf Divider | 0.7kg  |



# **Stackable Shelf Ledges**

#### **Features**

- Corrosion proof stackable shelf ledge (height 51mm)
- Prevents larger items from protruding or falling
- Compatible with grid & solid shelves
- Easy to clean

#### **Standard Colours**

Chrome

| Item       | Description         | Weight |
|------------|---------------------|--------|
| Side Ledge |                     |        |
| 8MXLS24-2S | Side Ledge - 610mm  | 0.2kg  |
| Rear Ledge |                     |        |
| 8MXL36-2S  | Rear Ledge - 914mm  | 0.4kg  |
| 8MXL42-2S  | Rear Ledge - 1066mm | 0.6kg  |
| 8MXL48-2S  | Rear Ledge - 1219mm | 0.6kg  |
| 8MXL60-2S  | Rear Ledge - 1524mm | 0.8kg  |

| www.bristolmaid.com |
|---------------------|
|---------------------|





**Mobile Units** 

sales@bristolmaid.com

# Polymer Racking - MetroMax Q<sup>®</sup>

#### **Features**

- · Epoxy / polymer shelving & racking system
- · Epoxy coated shelves & posts with a patented no tools lever shelf release mechanism, adjustable in 25mm increments
- · Removable easy to clean polymer open grid shelf mats (will fit most dishwashers)
- Open grid design allowing for the free circulation of air, improved
- visibility & light penetration
- Microban® product protection inhibits growth of bacteria & aids infection control
- Available as stationary or mobile units
- Stationary Fitted with adjustable feet to allow levelling on uneven floors
- Mobile Fitted with 125mm swivel / braking castors
- Rack / unit capacity...
  - Stationary 907kg
  - Mobile 339kg
- Optional label holders, shelf dividers & ledges
- · Anti-corrosion warranty
  - Frame & Posts 15 years
  - Shelf Mats Lifetime
- · Supplied for self-assembly



Contact us for more information

## **Stationary Units**

| 4 Shelves             | Shelf Size (w x d)                | Height           | Weight           | 4 Shelves              | Shelf Size (w x d   |
|-----------------------|-----------------------------------|------------------|------------------|------------------------|---|
| 8Q536G2               | 914 x 610mm                       | 1590mm           | 30.0kg           | 8Q536G2M               | 914 x 610mm   |
| 8Q546G2               | 1066 x 610mm                      | 1590mm           | 33.6kg           | 8Q546G2M               | 1066 x 610mm  |
| 8Q556G2               | 1219 x 610mm                      | 1590mm           | 35.6kg           | 8Q556G2M               | 1219 x 610mm  |
| 8Q566G2               | 1524 x 610mm                      | 1590mm           | 44.4kg           | 8Q566G2M               | 1524 x 610mm  |
|                       |                                   |                  |                  |                        |   |
| 5 Shelves             | Shelf Size (w x d)                | Height           | Weight           | 5 Shelves              | Shelf Size (w x d   |
| 5 Shelves<br>85Q537G2 | Shelf Size (w x d)<br>914 x 610mm | Height<br>1895mm | Weight<br>36.7kg | 5 Shelves<br>85Q537G2N |   |
|                       |                                   |                  |                  |                        | 1 914 x 610mm   |
| 85Q537G2              | 914 x 610mm                       | 1895mm           | 36.7kg           | 85Q537G2N              | 1 914 x 610mm<br>1 1066 x 610mm   |
| 85Q537G2<br>85Q547G2  | 914 x 610mm<br>1066 x 610mm       | 1895mm<br>1895mm | 36.7kg<br>41.2kg | 85Q537G2N<br>85Q547G2N | <ol> <li>914 x 610mm</li> <li>1066 x 610mm</li> <li>1219 x 610mm</li> </ol> |

85Q566G2

| bile Unit | S *other siz       | *other sizes available on request |        |  |  |  |
|-----------|--------------------|-----------------------------------|--------|--|--|--|
| helves    | Shelf Size (w x d) | Height                            | Weight |  |  |  |
| 536G2M    | 914 x 610mm        | 1730mm                            | 34.6kg |  |  |  |
| 546G2M    | 1066 x 610mm       | 1730mm                            | 38.2kg |  |  |  |
| 556G2M    | 1219 x 610mm       | 1730mm                            | 40.2kg |  |  |  |
| 566G2M    | 1524 x 610mm       | 1730mm                            | 49.0kg |  |  |  |
| helves    | Shelf Size (w x d) | Height                            | Weight |  |  |  |
| Q537G2M   | 914 x 610mm        | 2035mm                            | 41.3kg |  |  |  |
| Q547G2M   | 1066 x 610mm       | 2035mm                            | 45.8kg |  |  |  |

Room Design & Installation Service Available

2035mm

48.3kg

59.3kg



# **Label Holders**

#### **Features**

- · Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels


#### Dimensions (w x h)

• 102 x 27mm

#### **Standard Colours**

• Grey

8Q04LH



# **Shelf Divider**

#### **Features**

- Corrosion proof polymer shelf divider (height 200mm)
- Use to create organised compartments or simply separate shelf contents
- Easily repositioned on grid & solid shelves
- · Easy to clean

#### Dimensions (w x h)

• 610 x 200mm

#### **Standard Colours**

• Grey

| Weight | Item     | Description   | Weight |
|--------|----------|---------------|--------|
| 0.01kg | 8MQD24-8 | Shelf Divider | 0.7kg  |



# **Stackable Shelf Ledges**

#### **Features**

- Corrosion proof stackable shelf ledge (height 51mm)
- Prevents larger items from protruding or falling
- Compatible with grid & solid shelves
- Easy to clean

#### **Standard Colours**

Chrome

| Item       | Description         | Weight |
|------------|---------------------|--------|
| Side Ledge |                     |        |
| 8MQLS24-2S | Side Ledge - 610mm  | 0.2kg  |
| Rear Ledge |                     |        |
| 8MQL24-2S  | Rear Ledge - 610mm  | 0.2kg  |
| 8MQL36-2S  | Rear Ledge - 914mm  | 0.4kg  |
| 8MQL42-2S  | Rear Ledge - 1066mm | 0.6kg  |
| 8MQL48-2S  | Rear Ledge - 1219mm | 0.6kg  |
| 8MQL60-2S  | Rear Ledge - 1524mm | 0.8kg  |

| Label Holder | 0.01kg |
|--------------|--------|
|              |        |

| sales@bristolmaid.com |  |
|-----------------------|--|
|-----------------------|--|

# DLYIVIPIC



#### **Stationary Units**

#### **Mobile Units**


| 4 Shelves | Shelf Size (w x d) | Height | Weight | 4 Shelves | Shelf Size (w x d) | Height | W |
|-----------|--------------------|--------|--------|-----------|--------------------|--------|---|
| 80536C    | 914 x 610mm        | 1590mm | 30.0kg | 80536CM   | 914 x 610mm        | 1730mm | 3 |
| 80546C    | 1066 x 610mm       | 1590mm | 33.6kg | 80546CM   | 1066 x 610mm       | 1730mm | 3 |
| 80556C    | 1219 x 610mm       | 1590mm | 35.6kg | 80556CM   | 1219 x 610mm       | 1730mm | 4 |
| 80966C    | 1524 x 610mm       | 1590mm | 44.4kg | 80966CM   | 1524 x 610mm       | 1730mm | 4 |
| 5 Shelves | Shelf Size (w x d) | Height | Weight | 5 Shelves | Shelf Size (w x d) | Height | N |
| 850537C   | 914 x 610mm        | 1895mm | 36.7kg | 850537CM  | 914 x 610mm        | 2035mm | 2 |
| 850547C   | 1066 x 610mm       | 1895mm | 41.2kg | 850547CM  | 1066 x 610mm       | 2035mm | 2 |
| 850557C   | 1219 x 610mm       | 1895mm | 43.7kg | 850557CM  | 1219 x 610mm       | 2035mm | 2 |
| 850967C   | 1524 x 610mm       | 1895mm | 54.7ka | 850967CM  | 1524 x 610mm       | 2035mm | Ę |

sales@bristolmaid.com

# Zinc Chrome Racking - Olympic<sup>®</sup>

#### **Features**

- Zinc chrome shelving & racking system
- · Chromate posts & shelves with a patented no tools tapered split sleeve system that snaps together on each post, adjustable in 25mm increments
- Available as stationary or mobile units
  - Stationary Fitted with adjustable feet to allow levelling on uneven floors
  - Mobile Fitted with 125mm swivel / braking castors
- Rack / unit capacity...
  - Stationary 907kg
  - Mobile 339kg
- Optional label holders, shelf dividers & ledges
- 1-year anti-corrosion warranty
- Supplied for self-assembly



Zinc Chrome Shelves

#### **Room Design & Installation Service Available** Contact us for more information



# **Label Holders**

#### **Features**

- Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels

#### Dimensions (w x h)

• 75 x 32mm

#### **Standard Colours**

Grey



# Shelf Divider

#### Features

• Keeps shelf contents orderly (height 203mm)

+44 (0)1258 484455


• Pressure fit to hold in place

#### **Dimensions** (w x h)

• 610 x 203mm

#### **Standard Colours**

Chrome



# **Shelf Ledges**

#### **Features**

- Prevents items from protruding or falling (height 100mm)
- Pressure fit to hold in place
- Ledges can be stacked

#### **Standard Colours**

Chrome

| Item     | Shelf Length | Weight |
|----------|--------------|--------|
| 8L24N-4C | 610mm        | 0.90kg |
| 8L36N-4C | 914mm        | 1.25kg |
| 8L42N-4C | 1066mm       | 1.48kg |
| 8L48N-4C | 1219mm       | 1.68kg |
| 8L60N-4C | 1524mm       | 1.81kg |

| ltem   | Description  | Weight | Item   | Description   | Weight |
|--------|--------------|--------|--------|---------------|--------|
| 89990P | Label Holder | 0.1kg  | 8DD24C | Shelf Divider | 0.7kg  |







# High Density Storage - Top-Track®

#### **Features**

- · Increases the storage capacity of a given area by up to 50%
- Maximises all available space
- Can eliminate the need to expand a given area through the use of optimised storage
- Creates moveable aisles
- · Eliminates hard to reach or 'dead' spaces
- The guide track is positioned above the shelving system
- Aluminium overhead track system acts as a guide for mobile units
- Floors are easy to clean
- Carts & trolleys can easily be easily manoeuvred into & out of the active aisle
- Mobile units are designed to accommodate medium-duty applications
- The maximum weight capacity of a Top-Track  $^{\!\otimes}$  mobile unit is 410kg



Overhead Guide Track

Room Design & Installation Service Available Contact us for more information

# High Density Storage - qwikTRAK®

#### **Features**

- Increases the storage capacity of a given area by up to 50%
- Maximises all available space
- Can eliminate the need to expand a given area through the use of optimised storage
- · Creates moveable aisles
- · Eliminates hard to reach or 'dead' spaces
- The guide track is positioned below the shelving system, & protects the floor from wear & tear (does not need to be mounted to the floor)
- · Compensates for uneven floor surfaces
- Track acts as a guide for mobile units
- · Mobile units are designed to accommodate heavy-duty applications
- Allows for double deep configurations, where two systems are joined together for even greater storage capacity
- The maximum weight capacity of a qwikTRAK® mobile unit is 970kg









Floor Mounted Guide Track



# **Sterile Services Equipment**

- Preparation Tables Variable Height 402
  - Preparation Tables Fixed Height 403
- Sterile Supplies Trolleys UN3291 Certified 404
- Sterile Supplies Trolleys UN3291 Certified, Trolley Washer Compatible 406

## **Preparation Tables - Stainless Steel, Variable Height**

#### **Features**

- Height adjustable preparation tables
- 304 grade stainless steel frame & work surface
- Uses a manual turn handle (folds away when not in use) • Static - 715 - 1015mm
  - Mobile 812 1112mm
- Available as static or mobile
  - Static Adjustable feet for levelling on uneven floors
  - Mobile 80mm pressed steel castors, 2 off braking
- Max capacity 70.0kg

#### Dimensions (w x d x h)

- Static 1500 x 750 x 715 1015mm
- Mobile 1500 x 750 x 812 1112mm

| Stand | dard | Col | lours |
|-------|------|-----|-------|
| Juan  | Jaiu | CO  | iou s |



| ltem  | Description  | Weight |
|-------|--|--------|
| PT050 | Preparation Table - Variable Height, Stainless Steel, Static | 49.0kg |
| PT110 | Preparation Table - Variable Height, Stainless Steel, Mobile | 53.0kg |

# **Preparation Tables - Stainless Steel, Fixed Height**

#### Features

- Fixed height preparation tables
- 304 grade stainless steel frame, lower shelf & work surface
- Worksurface height 950mm
- Available as static or mobile
  - Static Adjustable feet for levelling on uneven floors
  - Mobile 80mm pressed steel castors, 2 off braking
- Max capacity 70.0kg

#### Dimensions (w x d x h)

• Overall Dimensions - 1500 x 750 x 950mm

#### **Standard Colours**

Stainless Steel



| Item      | Description   | Weight |
|-----------|---|--------|
| PTS/115   | Preparation Table - Fixed Height, Stainless Steel, Static | 56.3kg |
| PTS/115/M | Preparation Table - Fixed Height, Stainless Steel, Mobile | 58.8kg |



AGV Compatible Trolleys & Carts For Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy distribution

# Sterile Supplies Trolleys - UN3291 Certified

#### Features

- Independently tested to UN3291 Certificate of Packaging Performance for the Transportation of Clinical Waste
- · Suitable for the transportation of DIN & ISO instrument trays
- · Mild steel frame
- · Lightweight aluminium cladding, shelves & door/s
- · Powder coated for durability
- Locking door/s…
  - Fitted with a flush fitting handle fitted with a CAM lock offering two-point locking (supplied with 2 keys)
  - When opened, doors fold back flat at either end of the trolley\*
- Three adjustable shelves, (max 40kg per shelf)
- Shelf size...
  - Small 534 x 614mm (w x d)
  - Large 1128 x 614mm (w x d)
- Removable stainless steel drip tray
- · Inset push handles at either end of trolley
- Low level bumper rail offers protection to both the trolley & fabric of the building
- Optional tow bar fitting...
  - Tow bar fitted to the front (swivel castor)
  - Antiluce fastener fitted to rear (fixed castor)
  - Tow bar position approx. 360-380mm above floor
  - Recommended max speed 7kph for towed trolleys
- 200mm heavy duty cushioned castors, 2 x fixed, 2 x swivel braking
  Maximum capacity...
  - Small 84k.0g per trolley
  - Large 109.0kg per trolley

\*Doors will not fold back flat when a tow bar is fitted

#### Dimensions (w x d x h) Standard Colours

- Small 805 x 774 x 1257mm Frame Grey
- Large 1399 x 774 x 1257mm Cladding / Shelves / Doors White

| Item        | Description  | Weight  |
|-------------|--|---------|
| Small       |  |         |
| ST/340/UN   | Sterile Supplies Trolley - UN3291 Certified, Small         | 74.5kg  |
| ST/340/UN/T | Sterile Supplies Trolley - UN3291 Certified, Towing, Small | 78.5kg  |
| Large       |  |         |
| ST/690/UN   | Sterile Supplies Trolley - UN3291 Certified, Large         | 110.0kg |
| ST/690/UN/T | Sterile Supplies Trolley - UN3291 Certified, Towing, Large | 114.0kg |



# Sterile Supplies Trolleys - UN3291 Certified, Trolley Washer Compatible



#### Features

- Independently tested to UN3291 Certificate of Packaging Performance for the Transportation of Clinical Waste
- · Compatible with trolley washers
- Suitable for the transportation of DIN & ISO instrument trays
- 304 grade stainless steel frame, cladding, shelves & doors
- Natural unpolished finish
- Locking door/s...
  - Fitted with a flush fitting handle fitted with a CAM lock offering two-point locking (supplied with 2 keys)
  - When opened, doors fold back flat at either end of the trolley
  - Spring loaded plungers hold the doors open, essential when the trolley is being processed in a trolley washer disinfector
- One fixed, two adjustable shelves (max 50kg per shelf)
- Shelf size 1227 x 632mm (w x d)
- Stainless steel drip tray...
  - Removable for cleaning
  - Can be located at the rear of the trolley during cleaning, held in position using spring loaded plungers
- · Inset push handles at either end of trolley
- High level corner bumpers offer protection to both the trolley & fabric of the building
- Optional tow bar fitting...
  - Tow bar fitted to the front (swivel castor)
  - Antiluce fastener fitted to rear (fixed castor)
  - Tow bar position approx. 360-380mm above floor
  - Recommended max speed 7kph for towed trolleys
- 200mm heavy duty sealed bearing stainless steel castors, 2 x fixed, 2 x swivel braking
- Maximum capacity 136.0kg per trolley


Doors will fold back flat with tow bar fitted

#### Dimensions (w x d x h) Standard Colours

- 1407 x 762 x 1288mm
- Stainless Steel Natural Unpolished

AGV Compatible Trolleys & Carts For Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy distribution

|         | Item          | Description   | Weight  |
|---------|---------------|---|---------|
| in      | STSS/690/UN   | Sterile Supplies Trolley - Stainless Steel, Trolley Washer Compatible         | 131.0kg |
| oplies, | STSS/690/UN/T | Sterile Supplies Trolley - Stainless Steel, Trolley Washer Compatible, Towing | 137.0kg |





High Level Corner Bumpers





# Storage Cupboards, Cabinets & Lockers

- Storage Cupboards 410
- Storage Cabinets Single & Double Door 412
- Storage Cabinets First Aid & Clinical Reagent Cabinets 413
  - Storage Cabinets Traka 21 414
    - Storage Cabinets Key 415
  - Storage Lockers Personal Effects 416











Shelves & Dividers

## Storage Cupboards - Single Door, 1870mm High

#### **Features**

- Single door 1870mm high cupboards
- Single-point locking using a flush fitting CAM lock fitted with a handle (supplied with 2 kevs)
- Available with either...
  - Plastic Travs
  - Shelves
  - Shelves & Dividers
- 0.9mm (20 Gauge) mild steel welded construction
- Modular design allows for easy interchange with HTM 63 storage systems

#### **Plastic Trays**

- Interchangeable with our popular Caretrav range
- Supplied with 7 x 100mm & 5 x 150mm travs
- · Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Travs feature an integral handle with a plastic clip-on label cover (coloured label strips) supplied - white, blue, green, vellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray

#### Shelves

- · Adjustable at 50mm increments
- Supplied with 8 x shelves

#### Shelves & Dividers

- · Adjustable at 50mm increments
- · Removable dividers 100mm high, adjustable at 25mm increments
- Supplied with 8 x shelves & 16 x dividers

#### Dimensions (w x d x h)

#### Cupboard

• 500 x 536 x 1870mm

Travs

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

#### Standard Colours

· Paint - White

| Options / | Accessories |
|-----------|-------------|
| optione,  |             |

- SC/JK/1870 Joining Kit
- SC/DIV Set of 16 Dividers
- CT/TRAY/S Shallow Trav

83.6kg

64.0kg

74.0kg

- CT/TRAY/D Deep Tray
- Weight **R/H Hinged** L/H Hinged Description SC/R/1870/7S5D SC/L/1870/7S5D Single Door, 1870mm High, Plastic Trays SC/R/1870/FS SC/L/1870/FS Single Door, 1870mm High, Flat Shelves SC/R/1870/FSD SC/L/1870/FSD Single Door, 1870mm High, Flat Shelves & Dividers

SC/L/1870/7S5D



**Joining Kits** 

#### **Features**

• Allows two units to be joined together

| Item       | Description |  |
|------------|-------------|--|
| SC/JK/1870 | Joining Kit |  |

# **Dividers**

#### **Features**

www.bristolmaid.com

- Spare / replacement dividers for SC cupboards supplied with shelves & dividers
- 100mm high dividers locate into slots on the shelves, adjustable at 25mm increments
- Compatible with 1870mm high units



| Item   | Description        |
|--------|--------------------|
| SC/DIV | Set of 16 Dividers |

SC/L/1870/FSD



# **Storage Cabinets - Single & Double Door**

#### Features

- Adjustable & removable shelves
- Flush fitting handle fitted with a CAM lock offering single-point locking (supplied with 2 keys)
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems

#### **Standard Colours**

• Paint - White





| Item        | Shelves | Overall Dimension (w x d x h) | Weight |
|-------------|---------|-------------------------------|--------|
| Single Door |         |                               |        |
| SC020       | 2       | 600 x 300 x 600mm             | 19.5kg |
| SC025       | 3       | 600 x 300 x 900mm             | 26.0kg |
| Double Door |         |                               |        |
| SC040       | 4       | 1000 x 300 x 600mm            | 32.5kg |
| SC045       | 6       | 1000 x 300 x 900mm            | 45.0kg |

# **Storage Cabinets - First Aid & Clinical Reagents**

#### **Features**

- Removable shelves
- Flush fitting handle fitted with a CAM lock offering single-point locking (supplied with 2 keys)
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes in back (fixings not included)
- Supplied with 'First Aid' & 'Clinical Reagents' self-adhesive labels

#### Dimensions (w x d x h)

- Small 300 x 200 x 400mm
- Medium 400 x 200 x 400mm
- Large 400 x 200 x 600mm

#### **Standard Colours**

• Paint - White








| Item  | Description                                    | Shelves | Weight |
|-------|--|---------|--------|
| GC005 | First Aid & Clinical Reagents Cabinet - Small  | 2       | 5.0kg  |
| GC010 | First Aid & Clinical Reagents Cabinet - Medium | 2       | 6.0kg  |
| GC015 | First Aid & Clinical Reagents Cabinet - Large  | 2       | 8.0kg  |







# How does it work?



Access User enters unique PIN code to gain access to the cabinet



Kev Removal Keys & keysets are individually locked within the cabinet. Green LED's indicate those keys which the user can access. The user cannot remove keys highlighted by the Red LED's which continue to remain locked in place

sales@bristolmaid.com



Key Replacement An amber light indicates the return position for the key

# Storage Cabinets - Keys, Traka21

#### Features

- User, key & access rights administration
- Screen based audit trail & reporting capability or export via USB port
- 21 robust, long-life iFobs with security seals
- 21 locking receptors with LEDs
- Robust & compact aluminium & ABS housing with no exposed hinges or obvious access points
- Manual override & door release functions in case of emergency
- Audible alarms

#### Operation

- Stand-alone 'Plug & Play' solution with advanced RFID technology
- Capacity for up to 21 keys
- Touchscreen interface
- · PIN access to designated keys or keysets which are individually locked in place
- · Keys are securely attached using special security seals
- Simple to configure through setup wizard
- No requirement for network connection or PC
- · Concealed robust fixing points allow secure wall fixing
- Mains operated (optional battery back-up)

#### Dimensions (w x d x h)

• 246 x 95 x 427mm

#### Standard Colours

· Aluminium & ABS - White

#### Options

- SKR/25 Security Key Rings (x 25)
- SKR/100 Security Key Rings (x 100)
- SKB/BAT Battery Back-Up



#### Traka21 iFob

Once attached to an iFob with a robust security seal, individual keys or keysets are assigned a specific receptor within the Traka cabinet & locked in place until released by an authorised user

| Item        | Description           | Weight |
|-------------|-----------------------|--------|
| SKC21/TRAKA | Key Cabinet - TRAKA21 | 4.0kg  |

# **Storage Cabinets - Keys**

#### **Features**

- · Cabinets for the control & tracking of keys
- Available with either a CAM or electronic push button lock
  - · CAM Supplied with two keys
  - · Electronic Push Button Battery operated
- Adjustable hook bars accommodate different key lengths
- Removable control index for out of hours peace of mind
- Colour coded hook bars & tags ensure guick identification. (both numbered)
- · Provided with security seals & loops to secure keys that require evidence of tampering or removal


#### Standard Colours

· Paint - Grev

#### **Electronic Push Button Lock**

- Uses a 10 button keypad and a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- · Four operating codes with differing levels of authority
  - Master code
  - Sub-master code
  - User code
  - Technician code
- · User codes are 4 digits
- · Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- · Low battery warning with failure override using 9v battery across LED's (+/-) & master code





SKC30

# 22 23 24 32 33 42 43 44 45 46

SKC100

Overall Dimensions (w x d x h)

185 x 110 x 260mm

215 x 110 x 305mm

300 x 110 x 355mm

300 x 110 x 450mm

380 x 110 x 550mm

Electronic Push Button Lock

#### **CAM Lock**

| Item   | Key Capacity | Overall Dimensions (w x d x h) | Weight | Item      | Key Capacity |
|--------|--------------|--------------------------------|--------|-----------|--------------|
| SKC20  | 20 Keys      | 185 x 80 x 260mm               | 2.0kg  | SKC20/CL  | 20 Keys      |
| SKC30  | 30 Keys      | 215 x 80 x 305mm               | 3.0kg  | SKC30/CL  | 30 Keys      |
| SKC48  | 48 Keys      | 300 x 80 x 355mm               | 4.0kg  | SKC48/CL  | 48 Keys      |
| SKC64  | 64 Keys      | 300 x 80 x 450mm               | 5.0kg  | SKC64/CL  | 64 Keys      |
| SKC100 | 100 Keys     | 380 x 80 x 550mm               | 8.0kg  | SKC100/CL | 100 Keys     |

Weight

2.0kg

3.0kg

4.0kg

5.0kg

8.0kg

Storage Cupboards,

Cabinets

Qo

Lockers

#### Alternative colours available

Replace suffix on the product code with any of the following to select your colour i.e. Single Bank with a Green Door - PL080/GD





GD Green Door



RD

Red Door

BD

Blue Door



YD Yellow Door

8 0



Hasp & Staple Lock

# Ø Ø 8 ø ø

## **Storage Lockers - Personal Effects**

#### Features

- Range of personal effects lockers available in single, double or triple banks in a number of compartment configurations
- · Choice of locking options
  - CAM lock (supplied with 2 keys)
  - Hasp & staple (accepts a padlock/not supplied)
- · Lock surround incorporates a 25mm round label holder
- Mounting holes provided in back (fixings not included)

0

0

8

· Optional stands & sloping tops

#### Dimensions (w x d x h)

## Internal Compartments (w x d x h)

- Single Bank 305 x 460 x 1780mm
- Double Bank 610 x 460 x 1780mm
- Triple Bank 915 x 460 x 1780mm

#### Standard Colours

- Body Grey
- · Doors Blue

0

0

0

0

0

0

A

- One 255 x 440 x 1680mm • Two - 255 x 440 x 815mm
- Three 255 x 440 x 525mm
- Four 255 x 440 x 385mm
- Five 255 x 440 x 295mm





Storage Cupboards, Cabinets & Lockers









PL145/YD

| Single Bank   |          |          |          | 2<br>2<br>2<br>2 |          |
|---------------|----------|----------|----------|------------------|----------|
| Compartments  | 1        | 2        | 3        | 4                | 5        |
| Coat Hook     | Yes      | Yes      | No       | No               | No       |
| Shelf         | Yes      | No       | No       | No               | No       |
| Hasp & Staple | PL005/BD | PL010/BD | PL015/BD | PL020/BD         | PL025/BD |
| CAM Lock      | PL080/BD | PL085/BD | PL090/BD | PL095/BD         | PL100/BD |





# Stand

#### **Features**

- Compatible with single, double & triple bank 1780mm high personal effects lockers
- Keeps lockers clear of the floor assisting cleaning & aiding air-circulation
- Fitted with adjustable feet
- Increases the overall height of the locker by 150mm

#### **Standard Colours**





PL915 - Sloping top

| Item  | Description         |  |
|-------|---------------------|--|
| PL900 | Stand - Single Bank |  |
| PL905 | Stand - Double Bank |  |
| PL910 | Stand - Triple Bank |  |

# **Sloping Top**

#### Features

- Compatible with single, double & triple bank 1780mm high personal effects lockers
- Prevents items being stored on the top of the locker & assists cleaning
- Increases the overall height of the locker by 150mm

#### **Standard Colours**

Sloping Top - Grey

| ltem  | Description               |
|-------|---------------------------|
| PL915 | Sloping Top - Single Bank |
| PL920 | Sloping Top - Double Bank |
| PL925 | Sloping Top - Triple Bank |



sales@bristolmaid.com



# **Theatre Equipment**

- Bowl Stands 420
- Bucket Stands 423
- Instrument Trolleys 424
  - Mayo Tables 426
- Operating Platforms 426
- Suture Storage Trolleys 427
- Sterile Container Storage 428



# **Bowl Stands - Fixed Height**

#### **Features**

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Bowl retaining rings fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 75mm pressed steel non-marking castors (1 x anti-static)

#### **Standard Colours**

Stainless Steel

#### **Options / Accessories**

- 2MECO079 360mm Diameter, 7 Litre Stainless Steel Bowl
- 2PLCO054 360mm Diameter, 7 Litre Polythene Bowl



2MECO079 7 Litre Stainless Steel Bowl



2PLCO054 7 Litre Polythene Bowl

| ltem   | Description  | Overall Dimensions (w x d x h) | Weight |
|--------|--|--------------------------------|--------|
| SBS/2  | Bowl Stand - Fixed Height, Double Bowl, Side by Side | 750 x 445 x 860mm              | 6.0kg  |
| SBST/2 | Bowl Stand - Fixed Height, Double Bowl, Tiered       | 445 x 445 x 860mm              | 5.0kg  |

www.bristolmaid.com

# **Bowl Stands - Variable Height**

#### **Features**

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Height adjustable upright using a handwheel (753 1100mm)
- Bowl retaining ring fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 50mm polypropylene/thermoplastic rubber twin-wheel castors (1 x anti-static)

#### **Standard Colours**

Stainless Steel

#### **Options / Accessories**

- 2MECO079 360mm Diameter, 7 Litre Stainless Steel Bowl
- 2PLCO054 360mm Diameter, 7 Litre Polythene Bowl



2MECO079 7 Litre Stainless Steel Bowl



2PLCO054 7 Litre Polythene Bowl

| Item   | Description                               | Overall Dimensions (w x d x h) | Weight |
|--------|---|--------------------------------|--------|
| SBSA/1 | Bowl Stand - Variable Height, Single Bowl | 650 x 650 x 753 - 1100mm       | 4.8kg  |
| SBSA/2 | Bowl Stand - Variable Height, Double Bowl | 746 x 650 x 753 - 1100mm       | 5.6kg  |



# Bowl Stands - Variable Height, Polyurethane, Easy Clean



| Item       | Description   | Weight |
|------------|---|--------|
| BO/P/VS/SP | Bowl Stand - Variable Height, Polyurethane, Easy Clean, Single Bowl | 9.0kg  |
| BO/P/VD/SP | Bowl Stand - Variable Height, Polyurethane, Easy Clean, Double Bowl | 9.8kg  |

#### Features

- · Five leg base with low centre of gravity to assist stability
- 25mm diameter scratch resistant stainless steel upright
- Polyurethane overmoulded base with smooth edges & contours provides...
  - Water-resistant finish, allowing for easy & effective cleaning (including steam cleaning)
  - Shield for the twin-wheel castors
  - Impact protection to patients & staff, reducing the risk of skin abrasions & bruising
  - · Protection to the fabric of the building
- Coloured indicator caps provide department identification & segregation, reducing the risk of cross contamination
- Height adjustable upright using a handwheel (820 1125mm)
- Bowl retaining rings fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 50mm polypropylene/thermoplastic rubber twin-wheel non-marking castors (1 x anti-static)

#### Dimensions (w x d x h)

- Single Bowl 620 x 620 x 820 1125mm
- Double Bowl 746 x 620 x 820 1125mm

#### **Standard Colours**

- Base Grey
- Coloured Indicator Cap Spearmint
- Bowl Retaining Rings Stainless Steel

#### **Options / Accessories**

- 2MECO079 Stainless Steel Bowl
- 2PLCO054 Polythene Bowl
- Coloured Indicator Cap Blue / Green / Red / Yellow



# **Bowl Stand - Kickabout, 360mm Diameter**

#### **Features**

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- · Bowl retaining ring fitted with removable anti-static spiral buffers to support bowls
- · Supplied with a 7 litre, 360mm diameter stainless steel bowl
- 75mm pressed steel non-marking castors (1 x anti-static, 2 x non-anti-static)

#### Dimensions (w x d x h)

#### Standard Colours

• 365 x 365 x 235mm

#### Stainless Steel

#### **Options / Accessories**

- 2MECO079 360mm Diameter, Stainless Steel Bowl
- 2PLCO054 360mm Diameter, Polythene Bowl



| Item   | Description                            | Weight |
|--------|--|--------|
| SBSL/1 | Bowl Stand - Kickabout, 360mm Diameter | 3.2kg  |

# **Bucket Stands - Kickabout, 15-Litre Bucket**

#### **Features**

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Bucket retaining ring fitted with removable anti-static spiral buffers to support the bucket
- Model available fitted with a handle
- Supplied with a 15-litre graduated stainless steel bucket
- 75mm pressed steel non-marking castors (1 x anti-static, 2 x non-anti-static)

#### Dimensions (w x d x h)

#### Standard Colours

- SBKSL/1 460 x 460 x 470mm Stainless Steel SBKSI H/1- 460 x 460 x 470mm
- Supplied with a 15-litre
  - stainless steel bucket

#### **Options / Accessories**

• 2MECO081 - 15-litre graduated stainless steel bucket



| Item     | Description                            | Weight |
|----------|--|--------|
| SBKSL/1  | Bucket Stand - 15 Litre Bucket         | 5.0kg  |
| SBKSLH/1 | Bucket Stand - 15 Litre Bucket, Handle | 5.5kg  |

Supplied with a

7 litre stainless

steel bowl

+44 (0)1258 484455



# Fixed Shelves, Sides Down, Stainless Steel, 450mm Deep

#### **Features**

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

#### **Options / Accessories**

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 75mm, Pressed Steel Castors (4 x braking)
- Drawer Units page 119
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

www.bristolmaid.com

# Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep

#### **Features**

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel anti-static non-marking castors

#### **Options / Accessories**

- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C015 100mm, Pressed Steel Castors (4 x braking)
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

| Item          | Size   | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|---------------|--------|--------------------|--------------------------------|--------|
| DTSF/66/2/SD  | Small  | 600 x 600mm        | 620 x 620 x 870mm              | 17.5kg |
| DTSF/96/2/SD  | Medium | 900 x 600mm        | 920 x 620 x 870mm              | 24.0kg |
| DTSF/126/2/SD | Large  | 1200 x 600mm       | 1220 x 620 x 870mm             | 29.6kg |



# Mayo Tables - Cantilever, Variable Height


#### **Features**

- · Seamlessly welded frame from 304 grade stainless steel for strength & ease of cleaning
- Available with two castors & skid feet or four castors
- Height adjustment...
  - SS/22/1 893 1480mm
  - SS/22/2 858 1450mm
- Supplied with a removable stainless steel tray, 512 x 387 x 20mm (w x d x h)
- 50mm polypropylene/thermoplastic rubber twin-wheel anti-static castors

#### Dimensions (w x d x h)

- SS/22/1 641 x 450 x 893-1480mm

Supplied with a stainless



# **Operating Platforms**

#### **Features**

- Seamlessly welded frame from 304 grade stainless steel for strength & ease of cleaning
- Aluminium embossed non-slip tread
- · Non-slip, anti-static rubber feet
- Max capacity 190kg

#### Dimensions (w x d x h)

- Small 400 x 300mm
- Large 600 x 300mm



Aluminium embossed non-slip tread



| ltem   | Description                            | Weight |
|--------|--|--------|
| Small  |  |        |
| OP/431 | Operating Platform - Small, 150mm High | 4.2kg  |
| OP/432 | Operating Platform - Small, 200mm High | 4.4kg  |
| Large  |  |        |
| OP/631 | Operating Platform - Large, 150mm High | 6.0kg  |
| OP/632 | Operating Platform - Large, 200mm High | 6.4kg  |



| Item    | Description   | Weight |
|---------|---|--------|
| SS/22/* | Mayo Table - Cantilever, Variable Height, Two Castors & Skid Feet | 5.4kg  |
| SS/22/2 | 2 Mayo Table - Cantilever, Variable Height, Four Castors          | 5.8kg  |

sales@bristolmaid.com

# **Suture Storage Trolleys**

#### **Features**

- Angled shelves present contents for instant identification & easy access
- Allows logical organisation, eliminates confusion, allowing for easier & faster restocking
- Open wire design minimises dust accumulation, allowing for free circulation of air, improved visibility & light penetration
- No tools shelf adjustment in 25mm increments
- Each trolley supplied with 12 x shelf retainers to assist with segregation of contents
- Delivered flat pack for self-assembly
- 125mm pressed steel castors, 2 off braking

#### **Options / Accessories**

- 89990P Label Holder, 75 x 32mm
- 8DCR17C Shelf Retainers



| Item                         | Description   | Overall Dimensions (w x d x h) | Weight |  |  |
|------------------------------|---|--------------------------------|--------|--|--|
| 1 x Flat, 4 x Angled Shelves |   |                                |        |  |  |
| 8DC15EC                      | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 610mm  | 610 x 457 x 1524mm             | 46.5kg |  |  |
| 8DC35EC                      | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 914mm  | 914 x 457 x 1524mm             | 49.0kg |  |  |
| 8DC55EC                      | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 1219mm | 1219 x 457 x 1524mm            | 51.0kg |  |  |
| 5 x Angled Shelves           |   |                                |        |  |  |
| 8DC16EC                      | Trolley - Suture Storage, 5 x Angled Shelves, 610mm             | 610 x 457 x 1524mm             | 50.0kg |  |  |
| 8DC36EC                      | Trolley - Suture Storage, 5 x Angled Shelves, 914mm             | 914 x 457 x 1778mm             | 51.5kg |  |  |
| 8DC56EC                      | Trolley - Suture Storage, 5 x Angled Shelves, 1219mm            | 1219 x 457 x 1778mm            | 52.5kg |  |  |

+44 (0)1258 484455



# **Sterile Container Storage**

#### **Features**

- Fully configurable freestanding modular racking suitable for the safe storage of sterile containers
- Compatible with both blue wrap & rigid containers, up to max 300mm x 600mm (width x depth)
- Containers stored on either...
  - Telescopic shelves
  - Wire grid shelves
- Anodised extruded aluminium uprights & cross members
- · Bolted construction using high quality moulded fixings provided rigidity & allows the system to be altered easily
- Available as either static or mobile units
  - Static Adjustable feet allow levelling on uneven floors
  - Mobile 50 / 80mm pressed steel castors
- Optional wall mounting brackets allow static units to be fixed securely to the wall, (when units are loaded over 100kg)
- Supplied for self-assembly

#### Dimensions (w x d x h)

Static

- Single Bay 394 x 642 x 1940mm
- Double Bay 737 x 642 x 1940mm

Mobile (50/80mm Castors)

- Single Bay 394 x 642 x 1996 / 2035mm
- Double Bay 737 x 642 x 1996 / 2035mm

#### **Standard Colours**

- Uprights & Cross Members Anodised Aluminium
- Wire Grids Polished Stainless Steel
- Telescopic Shelves Anodised Aluminium
- Plastic Mouldings Green / Grey

#### **Options / Accessories**

- 2HH-7532293 Wall Mounting Bracket (x1)
- 2HH-7532277 10mm Plastic Cover Cap (x1)

#### Static

| ltem             | Description   | Weight |
|------------------|---|--------|
| Single Bay       |   |        |
| 4HH-SCS/S1/2W5T  | Static, Single Bay, 2 x Wire Grid Shelves & 5 x Telescopic Shelves  | 30.2kg |
| 4HH-SCS/S1/8W    | Static, Single Bay, 8 x Wire Grid Shelves                           | 21.4kg |
| Double Bay       |   |        |
| 4HH-SCS/S2/4W10T | Static, Double Bay, 4 x Wire Grid Shelves & 10 x Telescopic Shelves | 56.9kg |
| 4HH-SCS/S2/16W   | Static, Double Bay, 16 x Wire Grid Shelves                          | 38.0kg |
|                  |   |        |

#### Mobile

4

Λ

| ltem              | Description   | Weight |
|-------------------|---|--------|
| Single Bay        |   |        |
| 4HH-SCS/M51/2W5T  | Mobile, Single Bay, 2 x Wire Grid Shelves & 5 x Telescopic Shelves  | 31.1kg |
| 4HH-SCS/M51/8W    | Mobile, Single Bay, 8 x Wire Grid Shelves                           | 22.3kg |
| Double Bay        |   |        |
| 4HH-SCS/M52/4W10T | Mobile, Double Bay, 4 x Wire Grid Shelves & 10 x Telescopic Shelves | 58.2kg |
| 4HH-SCS/M52/16W   | Mobile, Double Bay, 16 x Wire Grid Shelves                          | 39.3kg |

sales@bristolmaid.com +44 (0)1258 484455



2HH-7532255 - Wire Grid Shelf



2HH-7532335 - Telescopic Shelves



50mm Braking Pressed Steel Castor



Adjustable Feet



4HH-SCS/M51/8W Shown with BBraun Aesculap Sterile Containers (not supplied)

4HH-SCS/S2/16W Shown with BBraun Aesculap Sterile Containers (not supplied)

**Components** 

- 2HH-7532300-HH Base Unit (x 1)
- 2HH-7532310 Extension Unit (x 1)
- 2HH-7532335 Telescopic Shelf (x 1)
- 2HH-7532255 Wire Grid Shelf (x 1)
- 2HH-7532356 Runners (Pair)
- 2HH-7532320 Adjustable Foot Set (x 4)
- 2HH-7532321 Adjustable Foot Set (x 2) • 2HH-7532322 - 50mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532323 50mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)
- 2HH-7532324 80mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532325 80mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)

Room Design & Installation Service Available Contact us for more information

sales@bristolmaid.com



# Wireless Mobile Carts & Workstations

- Mobile Workstations Non-Powered 432
  - Options & Accessories 434
  - Caretray Trolley Laptop 436



# **Mobile Workstations - Non-Powered**

# Features

- Range of non-powered mobile workstations
- Models available for use with laptops & vesa display mounts & configurable with a keyboard tray
  - Laptop Mount Secures the laptop, (compatible with laptops up to 482mm in width)
  - Vesa Display Mount Height adjustable, compatible with 75 x 75mm / 100 x 100mm square hole patterns
  - Keyboard Tray Retractable tray with wrist support & left / right mouse pad (accepts full size & Medigenic keyboards with room for mouse storage)
- Height adjustable using a gas damper (780 1130mm)
- Extruded aluminium column provides a smooth linear height adjustment
- Safety feature prevents trapping of the user beneath the workstation
- · Cast aluminium four leg base
- · Acrylic work surface
- Configurable with a range of accessories
- 100mm twin-wheel castors (4 x braking)

#### Dimensions (w x d x h)

- Overall 630 x 630 x 780 1130mm
- Work Surface 500 x 465mm
- Keyboard Tray 500 x 255mm
- · Laptop Mount Accepts laptops up to 482mm

#### **Standard Colours**

- Base White
- Work Surface Blue
- Upright & Sheet White

| Item            | Description   | Weight |
|-----------------|---|--------|
| 5RC/MW/NP       | Mobile Workstation - Non-Powered                                    | 19.6kg |
| 5RC/MW/NP/LM    | Mobile Workstation - Non-Powered, Laptop Mount                      | 20.2kg |
| 5RC/MW/NP/VM    | Mobile Workstation - Non-Powered, Vesa Display Mount                | 22.4kg |
| 5RC/MW/NP/VM/KT | Mobile Workstation - Non-Powered, Vesa Display Mount, Keyboard Tray | 23.6kg |

# **Options / Accessories**

- 5RC/SB Storage Basket
- 5RC/DRK Side Din Rail Kit
- 5RC/PSH Power Supply Unit Holder
- GMBH/15 Glove / Mask Box Holder
- CT/ALC Alcohol Gel Holder
- SBH/15 Sharps Box Holder



# **Storage Basket**

# **Features**

- · Wire storage basket
- Mounts to the aluminium extrusion

# **Standard Colours**

• Paint - Grey White

# Dimensions (w x d x h)

• 380 x 180 x 206mm



| Item   | Description    | Weight |
|--------|----------------|--------|
| 5RC/SB | Storage Basket | 1.0kg  |

# **Power Supply Unit Holder**

# Features

- Provides a location to store the power supply unit
- Mounts to the underside of the top

# **Standard Colours**

White

# Side Din Rail Kit

# **Features**

- Aluminium DIN rail
- Provides a rail measuring 10 x 25mm to hang / position accessories
- Mounts to the right / left side of the unit

# **Standard Colours**

• Aluminium / White

#### Dimensions (w x d x h)

• 300 x 10 x 25mm





| Item    | Description       | Weight |
|---------|-------------------|--------|
| 5RC/DRK | Side Din Rail Kit | 0.6kg  |

Wireless Mobile Carts & Workstations

# **Glove / Mask Box Holder, One Box**

# Features

- Accommodates glove & mask boxes up to 150mm in width
- Box presented in a portrait orientation
- Fitted with a 15mm hook

# **Standard Colours**

Paint - Grey White

# Dimensions (w x d x h)

• 143 x 122 x 177mm

# Alcohol Gel Holder

# Features

- Accepts bottles with 90 x 90mm square base
- Fitted with a 15mm hook

# **Standard Colours**

Paint - Grey White

# Dimensions (w x d x h)

• 95 x 110 x 115mm

| _ |        |                    |        |
|---|--------|--------------------|--------|
| t | Item   | Description        | Weight |
|   | CT/ALC | Alcohol Gel Holder | 0.3kg  |

# **Sharps Box Holder**

# **Features**

- Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a toggle
- Fitted with a 15mm hook

# **Standard Colours**

· Paint - Grey White

# Dimensions (w x d x h)

• 127 x 118 x 230mm

Item

**SBH/15** 



Description

Sharps Box Holder

| <u> </u> |
|----------|
|          |
|          |
| (D)      |
| P        |
| S        |
| Wireless |
|          |
| Ξ        |
| 0        |
| Б        |
|          |
| P        |
| 0        |
| <b></b>  |
| -        |
| arts     |
| 0)       |
| 20       |
|          |
| $\geq$   |
| 8        |
| Ξ.       |
| <u> </u> |
| St       |
| 6        |
| T.       |
| 0        |
| 4        |
| 5        |
|          |



| Item    | Description                      | Weight |
|---------|----------------------------------|--------|
| GMBH/15 | Glove / Mask Box Holder, One Box | 0.6kg  |

| sales@bristolmaid.com   |  |
|-------------------------|--|
| sales@pristoirrialu.com |  |

Weight

0.5kg



| Item                    | Description   | Weight |
|-------------------------|---|--------|
| High Security Bolt Lock |   |        |
| CTLAPTM1/10/3S2D/HSB    | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 39.3kg |
| CTLAPTM1/10/4S1D/HSB    | 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 39.0kg |
| CTLAPTM1/10/6S/HSB      | 6 x 100mm Trays, High Security Bolt Lock            | 40.1kg |

# **Caretray Trolley - Laptop**

# Features

- Trolleys designed for use with laptops
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Laptop Holder, Removable, Swing Arm
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking
- When not in use, trolleys storing medicines should be secured to the wall using a
   Security Clamp (see locking & security accessories)

# Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

# **Standard Colours**

· Paint - Grey White

# **Options / Accessories**

- Defibrillator Platform page 67
- Carerail Mounted page 68
   Laptop Holders page 74

- Castors & Ring Buffers page 76
   Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTLAPTM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 39.3kg |
| CTLAPTM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 39.0kg |
| CTLAPTM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 40.1kg |





# Workstations & Speciality Storage

- Workstations 452
- Storage Units 453
- Catheter Storage Units 454





Other configurations available. Contact our customer services team to discuss your requirements

# Workstations - Mobile, Polymer, Metro Starsys™

#### **Features**

- · Adaptable workstations
- Available as single, double & triple width units
- Moulded body, drawers & doors from a lightweight polymer that will not rust, chip, dent, peel or corrode
- Laminate worktop (height 914mm)
- 125mm polymer castors, 2 off braking

# **Standard Colours**

- Body / Drawers / Doors Light Taupe
- Drawer / Door Handles Slate Blue





+44 (0)1258 484455

Item

# Storage Units - Mobile, Polymer, Metro Starsys™

# Features

- Full height storage units
- Available as single, double & triple width units
- Moulded body, drawers & hinged/tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode
- 125mm polymer castors, 2 off braking

# **Standard Colours**

# **Options / Accessories**

• 8SXREHAN - Push Handle

• 8SXLABKIT - Drawer / Door Label Holders (x 10)

- Body / Drawers / Doors Light Taupe
- Drawer / Door Handles Slate Blue





| Item             | Description                    | Overall Dimensions (w x d x h) | Weight  |
|------------------|--------------------------------|--------------------------------|---------|
| Drawers & Hinged | Doors                          |                                |         |
| 8SXRSGS2         | Single, Drawers                | 578 x 632 x 1990mm             | 187.7kg |
| 8SXRD76CMHD1     | Double, Drawers & Hinged Doors | 1080 x 709 x 1990mm            | 159.2kg |
| 8SXRTGS2         | Triple, Drawers & Hinged Doors | 1669 x 709 x 1990mm            | 187.8kg |
| Hinged Doors     |                                |                                |         |
| 8SXRS76CM5       | Single, Door                   | 578 x 709 x 1990mm             | 106.0kg |
| 8SXRD76CM3       | Double, Doors                  | 1080 x 709 x 1990mm            | 113.0kg |
| 8SXRTGS3         | Triple, Doors                  | 1669 x 709 x 1990mm            | 158.0kg |
| Tambour Doors    |                                |                                |         |
| 8SXRS83MGSTD     | Single, Tambour Door           | 578 x 709 x 2096mm             | 108.9kg |
| 8SXRD83MGSTD     | Double, Tambour Doors          | 1080 x 709 x 2096mm            | 152.9kg |





sales@bristolmaid.com



# Catheter Storage Units - Mobile, Polymer, Metro Starsys™





#### **Features**

- Full height catheter storage units available with either hinged or tambour door/s
- Catheters supported on retractable slides providing easy to access & safe storage
- Available as single & double width units
  - Single 180 Catheters
  - Double 360 Catheters
- Moulded body & hinged or tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode
- 125mm polymer castors, 2 off braking

#### **Standard Colours**

- Body / Drawers / Doors Light Taupe
- Drawer / Door Handles Slate Blue

#### **Options / Accessories**

- 8SXRCATHHK Catheter Hooks, (supports 12 catheters)
- 8SXRLABKIT Drawer / Door Label Holders, (x 10)
- 8SXREHAN Push Handle

Other configurations available. Contact our customer services team to discuss your requirements

| Item           | Description                                    | Overall Dimensions (w x d x h) | Weight  |
|----------------|--|--------------------------------|---------|
| Hinged Doors   |  |                                |         |
| 8SXRS76MCATH1  | Single, L/H Hinged Door, 180 x Catheters       | 578 x 632 x 1990mm             | 100.1kg |
| 8SXRD76MCATH1  | Double, R/H & L/H Hinged Door, 360 x Catheters | 1080 x 632 x 1990mm            | 136.3kg |
| Tambour Doors  |  |                                |         |
| 8SXRS83MCATHTD | Single, Tambour Door, 180 x Catheters          | 578 x 709 x 2096mm             | 109.3kg |
| 8SXRD83MCATHTD | Double, Tambour Doors, 360 x Catheters         | 1080 x 709 x 2096mm            | 146.1kg |

# Catheter Storage Units - Stationary, Polymer, Metro Starsys™

# **Features**

- · Full height catheter storage units available with either hinged or tambour door/s
- Catheters supported on retractable slides providing easy to access & safe storage
- · Available as single & double width units
  - Single 180 Catheters
  - Double 360 Catheters
- Moulded body & hinged or tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode

# **Standard Colours**

- Body / Drawers / Doors Light Taupe
- Drawer / Door Handles Slate Blue

# **Options / Accessories**

- 8SXRCATHHK Catheter Hooks, (supports 12 catheters)
- 8SXRLABKIT Drawer / Door Label Holders, (x 10)
- 8SXR-FLR Floor Mounting Kit
- 8SXR72UUFLR Unit to Unit Filler Kit\*\*
- 8SXR72BKFLR30 Back to Wall Kit (Half Depth Units)\*\*\*
- 8SXR72BKFLR Back to Wall Kit (Full Depth Units)\*\*\*

\*\*Covers the seam between two adjacent units

\*\*\*Fills the gap when a unit it is pulled away from the wall to be flush with the front of the worktop



| Item           | Description                                       | Overall Dimensions (w x d x h) | Weight  |
|----------------|---|--------------------------------|---------|
| Hinged Doors   |   |                                |         |
| 8SXRS72SCATH1  | Single, L/H Hinged Door, 180 x Catheters          | 540 x 594 x 1842mm             | 75.5kg  |
| 8SXRD72SCATH1  | Double, R/H & L/H Hinged Door, 360 x Catheters    | 1049 x 594 x 1842mm            | 143.6kg |
| Tambour Doors  |   |                                |         |
| 8SXRS80HCATHTD | Single, Half Depth, Tambour Door, 180 x Catheters | 540 x 457 x 2032mm             | 60.3kg  |
| 8SXRD80HCATHTD | Double, Half Depth, Tambour Door, 360 x Catheters | 360 1052 x 457 x 2032mm        | 113.9kg |
| 8SXRS80SCATHTD | Single, Full Depth, Tambour Door, 180 x Catheters | 180 540 x 620 x 2032mm         | 84.4kg  |
| 8SXRD80SCATHTD | Double, Full Depth, Tambour Door, 360 x Catheters | 1052 x 620 x 2032mm            | 140.2kg |







| 05400004                     | 07         | 4A1RS12            |     |                        |            | 4HH-HTM71/M51/600                    | 205        | 5RC/MW/NP                                     | 432        |
|------------------------------|------------|--------------------|-----|------------------------|------------|--------------------------------------|------------|---|------------|
| 2FACO024                     |            | 4A1RS15            |     | 4B1TD10/D              | 200        | 4HH-HTM71/M52/400                    | 205        | 5RC/MW/NP/LM                                  | 432        |
|                              | 209        | 4A1RS18<br>4A2DR09 |     | 4B1TD10/E<br>4B1TD10/F | 200<br>200 | 4HH-HTM71/M52/600                    | 205        | 5RC/MW/NP/VM<br>5RC/MW/NP/VM/KT               | 432<br>432 |
|                              | 209<br>207 | 4A2DR09<br>4A2DR10 |     | 4B11D10/F<br>4B11D10/G |            | 4HH-HTM71/S1/400<br>4HH-HTM71/S1/600 | 204<br>204 | 5RC/MW/NF/VM/KT                               | 432        |
| 2HH-240820<br>2HH-240820     |            | 4A2DR10<br>4A2DR12 |     | 4B1TD10/G<br>4B1TD12/H |            | 4HH-HTM71/S1/600<br>4HH-HTM71/S2/400 | 204        | 5RC/SB  | 434        |
|                              | 200        | 4A2DR12            |     | 4B1TD12/I              |            | 4HH-HTM71/S2/600                     | 204        | 5RE/OC/PORT/VVDEL                             | 358        |
|                              | 208        | 4A2DR18            |     | 4B1TD12/J              |            | 4HH-SCS/M51/2W5T                     | 428        | 5RE/OC/PORT/VVLAVSAVL                         | 358        |
|                              | 207        | 4A2DSTM            |     | 4B1VD09/A              | 199        | 4HH-SCS/M51/8W                       | 428        | 5RE/RAE/ALWE/VVDEL                            | 342        |
| 2HH-241301                   |            | 4A2DSTS            |     | 4B1VD09/B              | 199        | 4HH-SCS/M52/16W                      | 428        | 5RE/RAE/ALWE/VVDEL/VPNB                       | 343        |
| 2HH-241302                   |            | 4A2OF09            |     | 4B1VD09/C              | 199        | 4HH-SCS/M52/4W10T                    | 428        | 5RE/RAE/ASHI/VVDEL                            | 346        |
| 2HH-241302                   |            | 4A2OF10            |     | 4B1VD10/D              | 199        | 4HH-SCS/S1/2W5T                      | 428        | 5RE/RAE/ASHI/VVDEL/VPNB                       | 347        |
|                              | 206        | 4A2OF12            | 196 | 4B1VD10/E              | 199        | 4HH-SCS/S1/8W                        | 428        | 5RE/RAE/BECK/VMULT                            | 348        |
| 2HH-7011100                  | 206        | 4A2OF15            | 196 | 4B1VD10/F              | 199        | 4HH-SCS/S2/16W                       | 428        | 5RE/RAE/SHAP/VVDEL                            | 340        |
| 2HH-7012100                  | 206        | 4A2OF18            |     | 4B1VD12/G              | 199        | 4HH-SCS/S2/4W10T                     | 428        | 5RE/RAE/SHAP/VVDEL/VPNB                       | 341        |
| 2HH-7012200                  | 206        | 4A2OTTM            |     | 4B1VD12/H              |            | 5160                                 | 331        | 5RE/RAE/TWYF/VVDEL                            | 344        |
| 2HH-7014050                  | 206        | 4A2OTTS            |     | 4B1VD15/I              | 199        | 5161                                 | 331        | 5RE/RAE/TWYF/VVDEL/VPNB                       | 345        |
| 2HH-7014100                  |            | 4A2RS09            |     | 4B1VD15/J              | 199        | 5162                                 | 331        | 5RE/RAM/ALWE/VVDEL                            |            |
| 2HH-7014200                  |            | 4A2RS10            |     | 4B2DR09                | 198        | 5163                                 | 331        | 5RE/RAM/ALWE/VVDEL/VPNB                       | 343        |
| 2HH-7810102                  |            | 4A2RS12            |     | 4B2DR10                |            | 5201                                 | 337        | 5RE/RAM/ASHI/VVDEL                            |            |
| 2HH-7810103                  |            | 4A2RS15            |     |                        |            | 5202                                 | 336        | 5RE/RAM/ASHI/VVDEL/VPNB                       | 347        |
| 2HH-7810104                  |            | 4A2RS18            |     | 4B2DR15                | 198        | 5205                                 | 337        | 5RE/RAM/SHAP/VVDEL                            |            |
|                              |            | 4A3DR09            |     | 4B2DR18<br>4B2OF09     |            | 5260<br>5262                         | 332<br>332 | 5RE/RAM/SHAP/VVDEL/VPNB<br>5RE/RAM/TWYF/VVDEL | 341<br>344 |
| 2HH-7811204N                 | 208        | 4A3DR10<br>4A3DR12 |     | 4B2OF09<br>4B2OF10     |            | 5262                                 | 363        | 5RE/RAM/TWYF/VVDEL/VPNB                       | 344        |
| 2HH-7811205N<br>2HH-7811205N | 207<br>208 | 4A3DR12<br>4A3DR15 |     | 4B2OF10<br>4B2OF12     |            | 5671                                 | 14         | 5RE/RPAH/BATH/VVDEL                           |            |
| 2HH-7811205N                 |            | 4A3DR18            |     | 4B2OF12<br>4B2OF15     |            | 5671                                 | 335        | 5RE/RPAH/BATH/VVDEL/VPNB                      | 353        |
| 2HH-7811304                  | 207        | 4A30F09            |     | 4B2OF18                |            | 59003                                | 333        | 5RE/RRAE/ALWE/VVDEL                           | 342        |
|                              | 207        | 4A3OF10            |     | 4B2RS09                |            | 5160D                                | 331        | 5RE/RRAE/ALWE/VVDEL/VPNB                      | 343        |
|                              | 208        | 4A3OF12            |     |                        |            | 5301/A                               |            | 5RE/RRAE/ASHI/VVDEL                           |            |
| 2HH-7811306                  |            | 4A3OF15            |     | 4B2RS12                |            | 5660XL                               |            | 5RE/RRAE/ASHI/VVDEL/VPNB                      | 347        |
| 2HH-7910201N                 | 209        | 4A3OF18            | 196 | 4B2RS15                | 197        | 5660XL                               | 334        | 5RE/RRAE/BECK/VMULT                           | 348        |
| 2HH-7910204N                 | 209        | 4A3OTTM            | 203 | 4B2RS18                | 197        | 5FSH20WH                             | 326        | 5RE/RRAE/HADL610/VVDEL/VPNB                   | 350        |
| 2HH-K100187                  | 209        | 4A3OTTS            |     | 4B2TD09/AA *           | 200        | 5HORI/1/A                            | 365        | 5RE/RRAE/HADL710/VVDEL/VPNB                   | 350        |
|                              | 206        | 4A3RS09            |     | 4B2VD09/AA*            | 199        | 5HORI/2/A                            | 365        | 5RE/RRAE/SHAP/VVDEL                           | 340        |
| 2LOCK217                     |            | 4A3RS10            |     | 4B2VD10/DD *           | 199        | 5HORI/3/A                            | 365        | 5RE/RRAE/SHAP/VVDEL/VPNB                      | 341        |
| 2LOCK218                     | 273        | 4A3RS12            |     | 4B3DR09                | 198        | 5ND/HAF/ALDE/VVDEL                   | 339        | 5RE/RRAE/TWYF/VVDEL                           | 344        |
| 2LOKE011                     |            | 4A3RS15            |     | 4B3DR10                | 198        | 5ND/HAR/ALDE/VVDEL                   |            | 5RE/RRAE/TWYF/VVDEL/VPNB                      | 345        |
| 2MECO079                     |            | 4A3RS18            |     | 4B3DR12                |            | 5ND/LC/ALDE/VVDEL                    |            | 5RE/RRPAH/DORC/VVDEL                          | 354        |
| 2PLC0026                     | 65         | 4B1DR09            |     | 4B3DR15                | 198        | 5ND/MAF/ALDE/VVDEL                   |            | 5RE/RRPAH/DORC/VVDEL/VPNB                     | 355        |
|                              | 308        | 4B1DR10            |     | 4B3DR18                |            | 5ND/MAR/ALDE/VVDEL                   | 339        | 5RLDF0119                                     | 165        |
|                              | 420        | 4B1DR12<br>4B1DR15 |     | 4B3OF09<br>4B3OF10     | 196<br>196 | 5NOVA/1<br>5NOVA/2                   | 364<br>364 | 5RLDF0119DIGILOCK<br>5RLDF0219                | 165<br>165 |
| 4A1DR09<br>4A1DR10           |            | 4B1DR15<br>4B1DR18 |     | 4B3OF10<br>4B3OF12     |            | 5PF332                               | 367        | 5RLDF0219<br>5RLDF0219DIGILOCK                | 165        |
| 4A1DR10<br>4A1DR12           |            | 4B1DK18<br>4B1DSTM |     | 4B30F12<br>4B30F15     |            | 5PF333                               | 367        | 5RLDF0519                                     | 165        |
| 4A1DR12                      |            | 4B1DSTS            |     | 4B30F18                |            | 5PF334                               | 367        | 5RLDF0519DIGILOCK                             | 165        |
| 4A1DR18                      |            | 4B10F09            |     |                        |            | 5PF335                               | 367        | 5RLDG0519                                     | 165        |
| 4A1DSTM                      |            | 4B1OF10            |     |                        |            | 5PT336                               | 369        | 5RLDG0519DIGILOCK                             | 165        |
| 4A1DSTS                      |            | 4B1OF12            |     | 4B3RS12                |            | 5PT337                               | 369        | 5RSH42WH                                      | 327        |
|                              | 196        | 4B1OF15            | 196 | 4B3RS15                | 197        | 5PT338                               | 369        | 5RSH70WH                                      | 327        |
| 4A10F10                      | 196        | 4B1OF18            | 196 | 4B3RS18                | 197        | 5PT339                               | 369        | 5SM6L   | 362        |
| 4A10F12                      | 196        | 4B1RS09            |     | 4B3VD09/AAA **         | 199        | 5PU151                               | 373        | 5TC101  | 370        |
| 4A10F15                      |            | 4B1RS10            |     |                        |            | 5PU152                               |            | 5TC102  | 370        |
| 4A10F18                      |            | 4B1RS12            |     |                        |            | 5PU153                               |            |   | 370        |
| 4A1OTTM                      |            | 4B1RS15            |     | 4D1DSTCS               |            | 5PU4                                 | 372        | 5TC108  | 371        |
| 4A1OTTS                      |            | 4B1RS18            |     |                        |            | 5PU5                                 | 372        |   | 371        |
|                              | 197        | 4B1TD09/A          |     | 4D1TD09/B              | 201        | 5PU6                                 |            | 5TC110  | 371        |
| 4A1RS10                      | 197        | 4B1TD09/B          | 200 | 4HH-HTM71/M51/400      | 205        | 5RC/DRK                              | 434        | 5TC201  | 3/4        |
|                              |            |                    |     |                        |            |                                      |            |   |            |

0-9

|    | - |   |  |
|----|---|---|--|
|    |   |   |  |
| 1  |   |   |  |
|    |   | - |  |
| r. |   |   |  |
|    |   |   |  |
| c  |   |   |  |
|    |   |   |  |

|   | BC3CPDWL/BE | 29         |
|---|-------------|------------|
|   | BC3CPDWL/GW | 29         |
|   | BC3CPDWR/BE | 29         |
|   | BC3CPDWR/GW | 29         |
|   | BC3CWL/BE   | 29         |
|   | BC3CWL/GW   | 29         |
|   | BC3CWR/BE   | 29         |
|   | BC3CWR/GW   | 29         |
|   |             |            |
|   | BC3D/BE     | 22         |
|   | BC3D/GW     | 22         |
|   | BC3DPD/BE   | 22         |
|   | BC3DPD/GW   | 22         |
|   | BC3DPDWL/BE | 29         |
|   | BC3DPDWL/GW | 29         |
|   | BC3DPDWR/BE | 29         |
|   | BC3DPDWR/GW | 29         |
|   | BC3DWL/BE   | 29         |
|   | BC3DWL/GW   | 29         |
|   | BC3DWR/BE   | 29         |
|   | BC3DWR/GW   | 29         |
|   | BCC1CWL/BE  | 26         |
|   | BCC1CWL/GW  | 26         |
|   | BCC1CWR/BE  | 26         |
|   | BCC1CWR/GW  | 26         |
|   | BCC1DWL/BE  | 26         |
|   | BCC1DWL/GW  | 26         |
|   | BCC1DWR/BE  | 26         |
|   | BCC1DWR/GW  | 26         |
|   | BH015/WHL   | 325        |
|   | BH030/WHL   | 325        |
|   |             |            |
|   | BO/P/VD/SP  | 422        |
|   | BO/P/VS/SP  | 422        |
|   | BR015/WHL   | 320        |
|   | BR020/WHL   | 320        |
|   | BR030/WHL   | 320        |
|   | BR045/WHL   | 321        |
|   | BR050/WHL   | 321        |
|   | BR060/WHL   | 321        |
|   | BR075/WHL   | 322        |
|   | BR080/WHL   | 322        |
|   | BR090/WHL   | 322        |
|   | BR105/WHL   | 323        |
|   | BR110/WHL   | 323        |
|   | BR120/WHL   | 323        |
|   | BR135/WHL   | 324        |
|   | BR140/WHL   | 324        |
|   | BR150/WHL   | 324        |
| С |             |            |
|   | C001        | 121        |
|   | C002        | 121        |
|   | C002        | 121        |
|   | C003        | 76         |
|   | C005        | 76         |
|   | C005        | 76         |
|   | C008        | 309        |
|   |             |            |
|   | C013        | 121<br>121 |
|   | C015        |            |
|   | CB/100/A    | 86         |
|   |             |            |

|   | 8SXRSARTH   | 298  |
|---|-------------|------|
|   | 8SXRSCAST   | 298  |
|   |             |      |
|   | 8SXRSCCU    | 298  |
|   | 8SXRSDRS    | 298  |
|   | 8SXRSENT1   | 440  |
|   | 8SXRSGS2    | 441  |
|   | 8SXRSIV     |      |
|   |             | 298  |
|   | 8SXRSTRAMA  | 298  |
|   | 8SXRT38MW14 | 440  |
|   | 8SXRTGS2    | 441  |
|   | 8SXRTGS3    | 441  |
| В | 00/1110000  |      |
| В | DA /0/D     | 0.40 |
|   | BA/3/B      | 240  |
|   | BA/6        | 240  |
|   | BA/6/A      | 240  |
|   | BA/VH       | 238  |
|   | BA/VH/NC    | 238  |
|   | BA030       | 242  |
|   |             | 20   |
|   | BC1C/BE     |      |
|   | BC1C/GW     | 20   |
|   | BC1CPD/BE   | 20   |
|   | BC1CPD/GW   | 20   |
|   | BC1CPDWL/BE | 27   |
|   | BC1CPDWL/GW | 27   |
|   |             |      |
|   | BC1CPDWR/GW | 27   |
|   | BC1CWL/BE   | 27   |
|   | BC1CWL/GW   | 27   |
|   | BC1CWR/BE   | 27   |
|   | BC1CWR/GW   | 27   |
|   | BC1CWRPD/BE | 27   |
|   |             |      |
|   | BC1D/BE     | 20   |
|   | BC1D/GW     | 20   |
|   | BC1DPD/BE   | 20   |
|   | BC1DPD/GW   | 20   |
|   | BC1DPDWL/BE | 27   |
|   | BC1DPDWL/GW | 27   |
|   |             |      |
|   | BC1DPDWR/BE | 27   |
|   | BC1DPDWR/GW | 27   |
|   | BC1DWL/BE   | 27   |
|   | BC1DWL/GW   | 27   |
|   | BC1DWR/BE   | 27   |
|   | BC1DWR/GW   | 27   |
|   | BC2C/BE     | 21   |
|   |             |      |
|   | BC2C/GW     | 21   |
|   | BC2CWL/BE   | 28   |
|   | BC2CWL/GW   | 28   |
|   | BC2CWR/BE   | 28   |
|   | BC2CWR/GW   | 28   |
|   | BC2D/BE     | 21   |
|   |             |      |
|   | BC2D/GW     | 21   |
|   | BC2DWL/BE   | 28   |
|   | BC2DWL/GW   | 28   |
|   | BC2DWR/BE   | 28   |
|   | BC2DWR/GW   | 28   |
|   |             |      |
|   | BC3C/BE     | 22   |
|   | BC3C/GW     | 22   |
|   | BC3CPD/BE   | 22   |
|   | BC3CPD/GW   | 22   |
|   |             |      |
|   |             |      |

| 8MX546MS          | 392 |
|-------------------|-----|
| 8MX546SG          | 392 |
| 8MX546SS          | 392 |
| 8MX556MG          | 392 |
| 8MX556MS          | 392 |
| 8MX556SG          | 392 |
| 8MX556SS          | 392 |
| 8MX566MG          | 392 |
| 8MX566MS          | 392 |
|                   |     |
| 8MX566SG          | 392 |
| 8MX566SS          | 392 |
| 8MXD24-8          | 393 |
| 8MXL36-2S         | 393 |
| 8MXL42-2S         | 393 |
| 8MXL48-2S         | 393 |
| 8MXL60-2S         | 393 |
| 8MXLS24-2S        | 393 |
| 8MY1627-24BU      | 192 |
| 8MY1627-34BU      | 192 |
| 8MY2030-24BU      | 192 |
|                   |     |
| 8MY2030-34BU      | 192 |
| 8MY2636-25BU      | 192 |
| 8MY2636-35BU      | 192 |
| 80536C            | 396 |
| 80536CM           | 396 |
| 80546C            | 396 |
| 80546CM           | 396 |
| 80556C            | 396 |
| 80556CM           | 396 |
| 80966C            | 396 |
| 80966CM           | 396 |
| 80966CM<br>8Q04LH | 395 |
|                   |     |
| 8Q536G2           | 394 |
| 8Q536G2M          | 394 |
| 8Q546G2           | 394 |
| 8Q546G2M          | 394 |
| 8Q556G2           | 394 |
| 8Q556G2M          | 394 |
| 8Q566G2           | 394 |
| 8Q566G2M          | 394 |
| 8SXRD40MW10       | 440 |
| 8SXRD72SCATH1     | 443 |
| 8SXRD76CM3        | 441 |
| 8SXRD76CMHD1      | 441 |
| 8SXRD76MCATH1     | 441 |
|                   |     |
| 8SXRD80HCATHTD    | 443 |
| 8SXRD80SCATHTD    | 443 |
| 8SXRD83MCATHTD    | 442 |
| 8SXRD83MGSTD      | 441 |
| 8SXRDENT2         | 440 |
| 8SXRS72SCATH1     | 443 |
| 8SXRS76CM5        | 441 |
| 8SXRS76MCATH1     | 442 |
| 8SXRS80HCATHTD    | 443 |
| 8SXRS80SCATHTD    | 443 |
| 8SXRS80SCATHTD    | 443 |
|                   |     |
| 8SXRS83MGSTD      | 441 |
| 8SXRSANES         | 298 |
|                   |     |

|                      | 000 |  |
|----------------------|-----|--|
| 850537CM             | 396 |  |
| 850547C              | 396 |  |
| 850547CM             | 396 |  |
| 850557C              | 396 |  |
| 85O557CM             | 396 |  |
| 85O967C              | 396 |  |
| 85O967CM             | 396 |  |
| 85Q537G2             | 394 |  |
| 85Q537G2M            | 394 |  |
| 85Q547G2             | 394 |  |
| 85Q547G2M            | 394 |  |
| 85Q557G2             | 394 |  |
| 85Q557G2M            | 394 |  |
| 85Q567G2             | 394 |  |
| 85Q567G2M            | 394 |  |
| 89989PX              | 393 |  |
|                      |     |  |
| 89990P               | 397 |  |
| 8BC2030-2DMB         | 192 |  |
| 8BC2030-3DMB         | 192 |  |
| 8BC2636-2DMB         | 192 |  |
| 8BC2636-3DMB         | 192 |  |
| 8DC15EC              | 427 |  |
| 8DC16EC              | 427 |  |
| 8DC35EC              | 427 |  |
| 8DC36EC              | 427 |  |
| 8DC55EC              | 427 |  |
| 8DC56EC              | 427 |  |
| 8DD24C               | 397 |  |
| 8FL21P3369           | 314 |  |
| 8FLCCU1              | 296 |  |
| 8FLCLINE             | 296 |  |
| 8FLCOM-LT            | 296 |  |
| 8FLCRP4              | 296 |  |
| 8FLPED1              | 296 |  |
| 8FLTMENT2            | 296 |  |
|                      | 397 |  |
| 8L24N-4C             |     |  |
| 8L36N-4C             | 397 |  |
| 8L42N-4C             | 397 |  |
| 8L48N-4C             | 397 |  |
| 8L60N-4C             | 397 |  |
| 8LECCRP2             | 312 |  |
| 8LECCRP3             | 312 |  |
| 8LECCRP4             | 312 |  |
| 8MBC1210DB           | 300 |  |
| 8MBX1310ANES4        | 300 |  |
| 8MBX3210TRAUMA       | 300 |  |
| 8MQD24-8             | 395 |  |
| 8MQL24-2S            | 395 |  |
| 8MQL36-2S            | 395 |  |
| 8MQL42-2S            | 395 |  |
| 8MQL48-2S            | 395 |  |
| 8MQL60-2S            | 395 |  |
| 8MQLS24-2S           | 395 |  |
| 8MX536MG             | 392 |  |
| 8MX536MS             | 392 |  |
| 8MX536SG             | 392 |  |
| 8MX536SS             | 392 |  |
| 8MX536SS<br>8MX546MG |     |  |
| 0IVIA040IVIG         | 392 |  |
|                      |     |  |

| 5TC203         374           5TC208         375           5TC209         375           5TC204         376           5TC304         376           5TC305         376           5TC308         377           5TC303         377           5TC304         378           5TC305         377           5TC301         378           5TC311         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0112-071-000         11           5X-015-061-000         11           5X-015-071-000         11  | 5TC202          | 374 |
|---|-----------------|-----|
| 5TC208         375           STC209         375           STC210         376           STC210         376           STC305         376           STC308         377           STC309         377           STC310         377           STC311         378           STC312         378           STC313         378           STC314         378           STC313         378           STC313         378           STC314         378           STC313         378           STC3100-061-000/K         8           SX-0109-061-000/K         8           SX-0110-071-000/K         8           SX-0112-071-000/K         11           SX-0150-061-000         11           SX-0150-071-000         11 <tr< td=""><td></td><td></td></tr<>                                    |                 |     |
| 5TC209         375           5TC210         375           5TC304         376           5TC305         376           5TC306         377           5TC303         377           5TC311         378           5TC312         378           5TC313         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000         11           5X-015-071-000         11           5X-015-071-000         11           5X-015-071-000         11           5Z-TA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360<   |                 |     |
| 5TC210         375           5TC304         376           5TC305         376           5TC306         377           5TC303         377           5TC310         377           5TC311         378           5TC313         378           5TC313         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0100-071-000/K         8           5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-01142-061-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-114-071-000         11           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360  |                 |     |
| 5TC304         376           5TC305         376           5TC308         377           5TC303         377           5TC310         377           5TC311         378           5TC312         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0100-071-000/K         8           5X-0110-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         11           5X-0114-071-000         12           5X-0115-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-151-071-000         11           5X-52ETA1/B         360           5ZETA1/C         360           5ZETA2/B/L         360 <td></td> <td></td>                                |                 |     |
| 5TC305         376           5TC306         377           5TC308         377           5TC310         377           5TC311         378           5TC312         378           5TC313         378           5TC314         378           5TC313         378           5X-0100-061-000/K         8           5X-0110-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA4/B         360           5ZETA4/B  |                 |     |
| 5TC306         376           5TC308         377           5TC309         377           5TC310         377           5TC311         378           5TC312         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0109-071-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-217A/70         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B  |                 |     |
| 5TC308         377           5TC309         377           5TC310         378           5TC311         378           5TC312         378           5TC313         378           5TC311         378           5TC312         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0112-061-000         11           5X-0150-061-000         11           5X-0150-071-000         10           5X-0150-071-000         11   |                 |     |
| 5TC309         377           5TC310         377           5TC311         378           5TC311         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-061-000/K         8           5X-0109-071-000/K         9           5X-0109-071-000/K         9           5X-0109-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0150-061-000         11           5X-0150-061-000         11           5X-0150-071-000         11           5X-0150-071-070         10           5X-150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         10           5ZETA/D         360           5ZETA/B         360           5ZETA/B         360  |                 |     |
| 5TC310         377           5TC311         378           5TC312         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0109-071-000/K         9           5X-0110-071-000/K         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360   |                 |     |
| 5TC311         378           5TC312         378           5TC312         378           5X-0100-061-000/K         8           5X-0100-071-000/K         8           5X-0100-071-000/K         8           5X-0100-071-000/K         9           5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-01142-061-000         12           5X-0150-071-000         11           5X-120-87         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/RC         360           5ZETA4/S         360           5ZETA4/S         360           5ZETA4/S         360           5ZETA4/S         360   |                 |     |
| 5TC312         378           5TC312         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-061-000/K         8           5X-0109-071-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         11           5X-0150-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA1/C         360           5ZETA2/B/RC         360           5ZETA4/S   |                 |     |
| 5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0109-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5X-0151-071-000         16           5ZETA1/K         360           5ZETA2/B/C         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B/C         360           5ZETA4/B/C         360           5ZETA4/B/C         360  |                 |     |
| 5X-0100-061-000/K         8           5X-0100-071-000/K         8           5X-0109-071-000/K         8           5X-0109-071-000/K         9           5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-0110-071-0000         12           5X-0142-071-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0152-071-000         11           5X-0152-071-000         11           5X-0152-071-000         11           5X-0152-071-000         11           5ZETA1/B         360           5ZETA2/B         360           5ZETA2/B/C         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360 <tr< td=""><td></td><td></td></tr<> |                 |     |
| 5X-0100-071-000/K         8           5X-0109-061-000/K         8           5X-0109-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0142-071-000         12           5X-0150-061-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/RC         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA4/B         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360                       |                 |     |
| 5X-0109-061-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         12           5X-0142-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA2/B/RC         360           5ZETA2/2B         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360  |                 |     |
| 5X-0109-071-000/K         8           5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0142-071-000         12           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         12           5ZETA1/B         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC49/B         360           5ZETAC49/B         360           5ZETAC49/B         360           5ZETAC49/B         360           5ZETAC49/B         360 <t< td=""><td></td><td></td></t<>  |                 |     |
| 5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-0142-061-000         12           5X-0142-071-000         12           5X-0142-071-000         11           5X-0150-061-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA1/CO         360           5ZETA2/B/LC         360           5ZETAC90/B         360           5ZETA4/B         360   |                 |     |
| 5X-0110-071-000/K         9           5X-0142-061-000         12           5X-0142-071-000         11           5X-0150-061-000         11           5X-0150-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         10           5ZETA/IC         360           5ZETA/ZB         360           5ZETA/ZB/RC         360           5ZETA/SB/RC         360           5ZETA/SB/RC         360           5ZETA/SB/RC         360           5ZETA/SB/RC         392           85MX536MG         392           85  |                 |     |
| 5X-0142-061-000         12           5X-0142-071-000         12           5X-0150-061-000         11           5X-0150-061-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5Z-0151-071-000         11           5Z-0151-071-000         12           5Z-0151-071-000         11           5Z-0151-071-000         11           5Z-0151-071-000         11           5Z-0151-071-000         360           5Z-0151-071-000         360           5Z-0151-071-08         360           5Z-0151-071-08         360           5Z-072/2B         360  |                 |     |
| 5X-0142-071-000         12           5X-0150-061-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA1/CO         360           5ZETA2         360           5ZETA2         360           5ZETA2/B/LC         360           5ZETA2/B/RC         392           85MX536MG         392           85MX546MG         392           85MX546MG  |                 |     |
| 5X-0150-061-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-061-000         11           5Z-0151-071-000         11           5ZETA1/B         360           5ZETA1/C         360           5ZETA2         360           5ZETA2/2B         360           5ZETA2/2B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAH         360           5ZETAH         360           5ZETAK         360           5   |                 |     |
| 5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA1/B         360           5ZETA1/C         360           5ZETA2/B/C         360           5ZETA2/2B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETAC45/B         392           85MX566MG         392           85MX556MG         392           85MX566MG         392   |                 |     |
| 5X-0151-061-000         11           5X-0151-071-000         11           SZETA1/B         360           5ZETA1/CO         360           5ZETA1/CO         360           5ZETA1/CO         360           5ZETA2/B         360           5ZETA2/B         360           5ZETA2/B/AC         360           5ZETA2/B/AC         360           5ZETA2/B/AC         360           5ZETA2/B/AC         360           5ZETA/SJOC         392           85MX536MG         392           85MX546MG         392           85MX546MG         392           85MX556MG         392      <   | 5X-0150-061-000 | 11  |
| 5X-0151-071-000         11           SZETA1         360           SZETA1/B         360           SZETA1/C         360           SZETA1/C         360           SZETA2/B         360           SZETA2/B         360           SZETA2/B/C         360           SZETA2/B/C         360           SZETA2/B/C         360           SZETA2/B/C         360           SZETAC45/B         360           SZETAC45/B         360           SZETAC45/B         360           SZETAC45/B         360           SZETAH         360  | 5X-0150-071-000 | 11  |
| 52ETA1         360           5ZETA1/C         360           5ZETA2/Z         360           5ZETA2/ZB         360           5ZETA2/ZB         360           5ZETA2/ZB         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAH         360           5ZETAH/B         360           5ZETAH/B         360           5ZETAH/B         360           5ZETAK         360           5ZETAK         360           5ZETAH/B         360           5ZETAK         392  | 5X-0151-061-000 | 11  |
| 52ETA1         360           5ZETA1/C         360           5ZETA2/Z         360           5ZETA2/ZB         360           5ZETA2/ZB         360           5ZETA2/ZB         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAH         360           5ZETAH/B         360           5ZETAH/B         360           5ZETAH/B         360           5ZETAK         360           5ZETAK         360           5ZETAH/B         360           5ZETAK         392  | 5X-0151-071-000 | 11  |
| 5ZETA1/CO         360           5ZETA2/2         360           5ZETA2/2B         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         392           85MX546KG         392           85MX556KG         392           85MX556KG         392           85MX556KG         392           85MX556KM         392  |                 | 360 |
| 5ZETA1/CO         360           5ZETA2/2         360           5ZETA2/2B         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         392           85MX546KG         392           85MX556KG         392           85MX556KG         392           85MX556KG         392           85MX556KM         392  | 5ZETA1/B        | 360 |
| 52ETA2         360           52ETA2/2B         360           52ETA2/B/RC         360           52ETA2/B/RC         360           52ETA2/B/RC         360           52ETA2/B/RC         360           52ETA42/B/RC         360           52ETAC45/5         360           52ETAC45/8         360           52ETAC45/8         360           52ETAC45/8         360           52ETAH         392           85MX546MG         392           85MX556MG         392      8  |                 | 360 |
| 5ZETA2/28 360<br>5ZETA2/B/LC 360<br>5ZETA2/B/RC 360<br>5ZETAC45 360<br>5ZETAC45 360<br>5ZETAC45 360<br>5ZETAC90/B 360<br>5ZETAC90/B 360<br>5ZETAC90/B 360<br>5ZETAH/B 360<br>5ZETAH/B 360<br>5ZETA4/B 360<br>5ZETA4/B 360<br>85MX536MG 392<br>85MX536S 392<br>85MX546MG 392<br>85MX546MG 392<br>85MX546MG 392<br>85MX556MG 392<br>85MX566MG 392   |                 | 360 |
| 5ZETA2/B/LC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA42/B/RC         360           5ZETAC45/B         360           5ZETAC9/B         360           5ZETAC90/B         360           5ZETAQ         360           5ZETAQ         360           5ZETAG90/B         392           85MX536MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392  |                 |     |
| 5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/5         360           5ZETAC45/5         360           5ZETAC45/5         360           5ZETAC45/5         360           5ZETAC90/B         360           5ZETAM         360           5ZETAH         360           5ZETAK         392           85MX546MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX566MS         392           85MX566MS         392   |                 |     |
| 5ZETA930C         360           5ZETAC45         360           5ZETAC45         360           5ZETAC45         360           5ZETAC45         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC49/B         360           5ZETAH/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA/B         392           85MX546MG         392           85MX546MG         392           85MX556MS         392           85MX556MS         392           85MX556MG         392           85MX556MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392 </td <td></td> <td></td>   |                 |     |
| 5ZETAC45         360           5ZETAC45/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAH         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAQ         360           5ZETAQ         360           85MX536MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX566MG  |                 |     |
| 5ZETAC45/B         360           5ZETAC90         360           5ZETAC90/B         360           5ZETAH         360           5ZETAH         360           5ZETAM         360           5ZETAH         360           5ZETAM         360           5ZETAM         360           5ZETAM         360           5ZETAM         360           5ZETAM         360           5ZETAM         392           85MX536MS         392           85MX546MG         392           85MX546MS         392           85MX556MS         392           85MX566MS         392 </td <td></td> <td></td>   |                 |     |
| 5ZETAC90 360<br>5ZETAC90/B 360<br>5ZETAH/B 360<br>5ZETAH/B 360<br>5ZETAH/B 360<br>85MX536MG 392<br>85MX536MG 392<br>85MX536MS 392<br>85MX546MG 392<br>85MX546MS 392<br>85MX546MS 392<br>85MX546MS 392<br>85MX556MS 392<br>85MX556MS 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MS 392<br>85MX566MS 392<br>85MX566MS 392<br>85MX566MS 392<br>85MX566MS 392  |                 |     |
| 5ZETAC90/B 360<br>5ZETAH/B 360<br>5ZETAH/B 360<br>5ZETAQ 360<br>85MX536MG 392<br>85MX536S 392<br>85MX536S 392<br>85MX546MG 392<br>85MX546MG 392<br>85MX546MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX566MG 392<br>85MX566MG 392<br>85MX566MG 392<br>85MX566MG 392<br>85MX566MS 392<br>85MX566MS 392   |                 |     |
| 5ZETAH         360           5ZETAH/B         360           5ZETAQ         360           5ZETAQ         360           SSMX536MS         392           85MX536MS         392           85MX536MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX556MS         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566SMS         392           85MX566SMS         392           85MX566SMS         392           85MX566SMS         392           85MX566SS   |                 |     |
| 5ZETAH/B 360<br>5ZETAQ 360<br>85MX536MG 392<br>85MX536MS 392<br>85MX536SG 392<br>85MX546SG 392<br>85MX546MG 392<br>85MX546MS 392<br>85MX546SG 392<br>85MX556MS 392<br>85MX556MS 392<br>85MX556MS 392<br>85MX5566MG 392<br>85MX5566MG 392<br>85MX566MG 392<br>85MX566MS 392  |                 |     |
| 5ZETAQ         360           85MX536MG         392           85MX536MS         392           85MX536S         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392   |                 |     |
| 85MX536MG         392           85MX536MS         392           85MX536SG         392           85MX546MG         392           85MX546MG         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX556MG         392           85MX556MS         392           85MX566MS         392   |                 |     |
| 85MX536MS         392           85MX536SG         392           85MX536SG         392           85MX546MG         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX556MS         392           85MX556MS         392           85MX556SMS         392           85MX5566MS         392           85MX5566MS         392           85MX5566MS         392           85MX566MS         392           85MX566SS         392   |                 |     |
| 85MX536SG         392           85MX536SS         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX546MS         392           85MX556MS         392           85MX556MG         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MG         392           85MX566MG         392           85MX566MS         392   |                 |     |
| 85MX536SS         392           85MX546MG         392           85MX546MG         392           85MX546SG         392           85MX556MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566SG         392           85MX566SG         392   |                 |     |
| 85MX546MG         392           85MX546MS         392           85MX546SG         392           85MX556MG         392           85MX556MS         392           85MX556SMS         392           85MX556SMS         392           85MX556SMS         392           85MX556SMS         392           85MX556SMS         392           85MX566MS         392           85MX566SS         392  |                 |     |
| 85MX546MS         392           85MX546SG         392           85MX556MG         392           85MX556MG         392           85MX556MS         392           85MX556SG         392           85MX556MG         392           85MX556MS         392           85MX5566MG         392           85MX566MG         392           85MX566MG         392           85MX566MS         392           85MX566SS         392           85MX566SS         392  |                 |     |
| 85MX546SG         392           85MX546SS         392           85MX556MG         392           85MX556SG         392           85MX556SG         392           85MX556SG         392           85MX556SG         392           85MX556SG         392           85MX566MG         392           85MX566MG         392           85MX566SG         392           85MX566SG         392           85MX566SG         392           85MX566SG         392           85MX566SS         392           85MX566SS         392   |                 |     |
| 85MX546SS         392           85MX556MG         392           85MX556MS         392           85MX556SS         392           85MX556SG         392           85MX566MG         392           85MX566MS         392   |                 |     |
| 85MX556MG         392           85MX556MS         392           85MX556SS         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566SMS         392           85MX566SMS         392           85MX566SS         392   |                 |     |
| 85MX5566MS 392<br>85MX556SG 392<br>85MX556SS 392<br>85MX566MG 392<br>85MX566MS 392<br>85MX566GMS 392<br>85MX566SS 392   |                 |     |
| 85MX556SG         392           85MX556SS         392           85MX566MG         392           85MX566MG         392           85MX566SS         392           85MX566SS         392           85MX566SS         392   |                 |     |
| 85MX556SS         392           85MX566MG         392           85MX566MS         392           85MX566SG         392           85MX566SS         392           85MX566SS         392   |                 |     |
| 85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566SG         392           85MX566SS         392   |                 |     |
| 85MX566MS 392<br>85MX566SG 392<br>85MX566SS 392   |                 |     |
| 85MX566SG 392<br>85MX566SS 392  |                 |     |
| 85MX566SS 392   |                 |     |
|   |                 |     |
| 850537C 396   |                 |     |
|   | 850537C         | 396 |

|                    | 00010/00/000           |                         | 077 | 07010014/40/00/500   | 105 |                    | 000 |
|--------------------|------------------------|-------------------------|-----|----------------------|-----|--------------------|-----|
| CB/100/B 86        | CRSM2/18/16D 4         |                         | 277 | CTDISPM1/10/6S/EPB   | 125 |                    | 260 |
| CB005 88           | CRSM2/18/4S12D 4       |                         | 276 | CTDISPM1/10/6S/HSB   | 124 |                    | 260 |
| CB010 88           | CRSM2/18/6S12D 4       |                         | 292 | CTINFEM1/10/3S2D/EPB | 283 |                    | 260 |
| CB015 88           | CRSM2/18/8S10D 4       |                         | 292 | CTINFEM1/10/3S2D/HSB | 282 |                    | 261 |
| CB020 88           | CRSS1/16/2S6D 4        |                         |     | CTINFEM1/10/4S1D/HSB | 282 |                    | 261 |
| CD005 152          | CRSS1/16/3S6D 4        |                         | 279 | CTINFEM1/10/6S/EPB   | 283 |                    | 261 |
| CD010 152          | CRSS1/16/4S5D 4        |                         | 278 | CTINFEM1/10/6S/HSB   | 282 |                    | 288 |
| CD015 152          | CRSS1/16/5S5D 4        |                         | 279 | CTIVPRM1/10/3S2D/EPB | 285 |                    | 288 |
| CD020 152          | CRSS1/16/6S4D 4        |                         |     | CTIVPRM1/10/3S2D/HSB | 284 |                    | 310 |
| CD025 152          | CRSS1/16/8D 4          |                         | 279 | CTIVPRM1/10/4S1D/EPB | 285 |                    | 310 |
| CD030 153          | CRSS2/16/10S10D 4      |                         | 82  | CTIVPRM1/10/4S1D/HSB | 284 |                    | 310 |
| CD035 153          | CRSS2/16/12S8D 4       |                         | 82  | CTIVPRM1/10/6S/EPB   | 285 | CTS1/08/1S2D       | 40  |
| CD045 153          | CRSS2/16/16D 4         |                         | 83  | CTIVPRM1/10/6S/HSB   | 284 | CTS1/08/3D         | 40  |
| CD110 152          | CRSS2/16/4S12D 4       |                         | 83  | CTLAPTM1/10/3S2D/EPB | 437 | CTS1/08/3S1D       | 40  |
| CD115 152          | CRSS2/16/6S12D 4       |                         | 283 | CTLAPTM1/10/3S2D/HSB | 436 | CTS1/08/4S         | 40  |
| CD120 152          | CRSS2/16/8S10D 4       |                         | 280 | CTLAPTM1/10/4S1D/EPB | 437 | CTS1/10/1S3D       | 42  |
| CD125 152          |                        | 74 CTCRITM1/10/3S2D/EPB | 281 | CTLAPTM1/10/4S1D/HSB | 436 | CTS1/10/3S2D       | 42  |
| CD130 153          |                        | 73 CTCRITM1/10/3S2D/EPB | 281 | CTLAPTM1/10/6S/EPB   | 437 | CTS1/10/4D         | 42  |
| CD135 153          | CS010 1                |                         | 280 | CTLAPTM1/10/6S/HSB   | 436 | CTS1/10/4S1D       | 42  |
| CD145 153          | CS015 1                |                         |     | CTM1/08/1S2D         | 50  | CTS1/10/6S         | 42  |
| CD905 159          | CT/AIR 6               |                         | 281 | CTM1/08/3D           | 50  | CTS1/14/2S5D       | 44  |
| CD910 158          | CT/ALC 6               |                         | 61  | CTM1/08/3S1D         | 50  | CTS1/14/3S4D       | 44  |
| CD911 158          |                        | 35 CTD1/08EPB           | 62  | CTM1/08/4S           | 50  | CTS1/14/4S3D       | 44  |
| CD913 158          |                        | 92 CTD1/08HSB           | 63  | CTM1/10/1S3D         | 52  | CTS1/14/5S3D       | 44  |
| CF/F/PK10 252      | CT/BIN 6               | 8 CTD1/08SS             | 64  | CTM1/10/3S2D         | 52  | CTS1/14/6D         | 44  |
| CP 91              | CT/BOW 6               |                         | 61  | CTM1/10/4D           | 52  | CTS1/14/8S1D       | 44  |
| CRMM1/18/2S6D 56   | CT/BUC 6               |                         | 62  | CTM1/10/4S1D         | 52  | CTS2/08/3S4D       | 41  |
| CRMM1/18/3S6D 56   | CT/CR/2 7              |                         | 63  | CTM1/10/6S           | 52  | CTS2/08/4S3D       | 41  |
| CRMM1/18/4S5D 56   | CT/CR/5 7              | 0 CTD1/10SS             | 64  | CTM1/14/2S5D         | 54  | CTS2/08/6D         | 41  |
| CRMM1/18/5S5D 56   | CT/CYD 7               | 0 CTD1/14CA             | 61  | CTM1/14/3S4D         | 54  | CTS2/08/6S2D       | 41  |
| CRMM1/18/6S4D 56   | CT/DEFIB 6             | 7 CTD1/14EPB            | 62  | CTM1/14/4S3D         | 54  | CTS2/08/7S1D       | 41  |
| CRMM1/18/8D 56     | CT/DH/A4 6             |                         | 63  | CTM1/14/5S3D         | 54  | CTS2/08/8S         | 41  |
| CRMM2/18/10S10D 57 | CT/DIV/D 7             |                         | 64  | CTM1/14/6D           | 54  | CTS2/10/3S6D       | 43  |
| CRMM2/18/12S8D 57  | CT/DIV/S 7             |                         | 61  | CTM1/14/8S1D         | 54  | CTS2/10/4S5D       | 43  |
| CRMM2/18/16D 57    | CT/FS 7                |                         | 62  | CTM2/08/3S4D         | 51  | CTS2/10/5S4D       | 43  |
| CRMM2/18/4S12D 57  | CT/FWS 7               |                         | 63  | CTM2/08/4S3D         | 51  | CTS2/10/6S4D       | 43  |
| CRMM2/18/6S12D 57  | CT/IV 7                |                         | 64  | CTM2/08/6D           | 51  | CTS2/10/8D         | 43  |
| CRMM2/18/8S10D 57  | CT/LAB 7               |                         | 61  | CTM2/08/6S2D         | 51  | CTS2/10/8S2D       | 43  |
| CRMS1/16/2S6D 58   | CT/LAB/SEP 7           | 8 CTD2/10EPB            | 62  | CTM2/08/7S1D         | 51  | CTS2/14/10S6D      | 45  |
| CRMS1/16/3S6D 58   | CT/LHNR/S 7            | 5 CTD2/10HSB            | 63  | CTM2/08/8S           | 51  | CTS2/14/12D        | 45  |
| CRMS1/16/4S5D 58   | CT/LHR/S 7             | 4 CTD2/10SS             | 64  | CTM2/10/3S6D         | 53  | CTS2/14/16S2D      | 45  |
| CRMS1/16/5S5D 58   | CT/OB/MS 6             | 6 CTD2/14CA             | 61  | CTM2/10/4S5D         | 53  | CTS2/14/4S10D      | 45  |
| CRMS1/16/6S4D 58   | CT/OB/MS 3             | 11 CTD2/14EPB           | 62  | CTM2/10/5S4D         | 53  | CTS2/14/6S8D       | 45  |
| CRMS1/16/8D 58     | CT/OB/SS 6             | 6 CTD2/14HSB            | 63  | CTM2/10/6S4D         | 53  | CTS2/14/8S6D       | 45  |
| CRMS2/16/10S10D 59 | CT/PARM 7              | 1 CTD2/14SS             | 64  | CTM2/10/8D           | 53  | CTSEPSM1/10/6S/EPB | 387 |
| CRMS2/16/12S8D 59  | CT/SBB/15 7            | 1 CTDC1/08CA            | 61  | CTM2/10/8S2D         | 53  | CTSEPSM1/10/6S/HSB | 386 |
| CRMS2/16/16D 59    | CT/SPS 7               | 1 CTDC1/08EPB           | 62  | CTM2/14/10S6D        | 55  | CTVACIM1/10/3S2D   | 290 |
| CRMS2/16/4S12D 59  | CT/TRAY/D 7            | 7 CTDC1/08SS            | 64  | CTM2/14/12D          | 55  | CTVACIM1/10/4S1D   | 290 |
| CRMS2/16/6S12D 59  | CT/TRAY/S 7            | 7 CTDC1/10CA            | 61  | CTM2/14/16S2D        | 55  | CTVACIM1/10/6S     | 290 |
| CRMS2/16/8S10D 59  | CT/WB 7                | 2 CTDC1/10EPB           | 62  | CTM2/14/4S10D        | 55  | CWF/80             | 84  |
| CRSM1/18/2S6D 46   | CT/WF 7                | 2 CTDC1/10SS            | 64  | CTM2/14/6S8D         | 55  | CWF/80/SP          | 84  |
| CRSM1/18/3S6D 46   | CTAIRWM1/10/3S2D/SS 3  | 16 CTDC1/14CA           | 61  | CTM2/14/8S6D         | 55  | CYB/DE2            | 103 |
| CRSM1/18/4S5D 46   |                        | 16 CTDC1/14EPB          | 62  | CTMEDSM1/10/3S2D/CA  | 286 | CYB/F2             | 103 |
| CRSM1/18/5S5D 46   | CTAIRWM1/10/6S/SS 3    | 16 CTDC1/14SS           | 64  | CTMEDSM1/10/3S2D/EPB | 287 | CYB/G2             | 103 |
| CRSM1/18/6S4D 46   | CTANAEM1/10/3S2D/EPB 2 | 77 CTDISPM1/10/3S2D/EPB | 125 | CTMEDSM1/10/4S1D/CA  | 286 | CYB/J2             | 103 |
| CRSM1/18/8D 46     | CTANAEM1/10/3S2D/HSB 2 | 76 CTDISPM1/10/3S2D/HSB | 124 | CTMEDSM1/10/4S1D/EPB | 287 |                    | 103 |
| CRSM2/18/10S10D 47 |                        | 77 CTDISPM1/10/4S1D/EPB | 125 | CTMEDSM1/10/6S/CA    | 286 |                    | 103 |
| CRSM2/18/12S8D 47  | CTANAEM1/10/4S1D/HSB 2 | 76 CTDISPM1/10/4S1D/HSB | 124 | CTMEDSM1/10/6S/EPB   | 287 | CYHP/D1            | 272 |
|                    |                        |                         |     |                      |     |                    |     |

| F | 2 |  |
|---|---|--|
| L | 5 |  |
|   | D |  |
| s |   |  |

| G/RAIL-C                      | 308        | MRW/D4D/CAM              | 249      |
|-------------------------------|------------|--------------------------|----------|
| G/SB                          | 273        | MRW/D4D/CL               | 249      |
| GC005                         | 413        | MRW/LHNR/S               | 251      |
| GC010                         | 413        | MRW/LHR/S                | 250      |
| GC015                         | 413        | MRW/S1DOC/CAM            | 248      |
| GH005                         | 108        | MRW/S1DOC/CL             | 248      |
| GH015                         | 108        | MRW/S1DXP/CAM            | 248      |
| GMBH/15                       | 73         | MRW/S1DXP/CL             | 248      |
| GMBH/15                       | 107        | MRW/S2D/CAM              | 248      |
| GMBH/15                       | 435        | MRW/S2D/CL<br>MS/C4      | 248      |
| GMBH/WM                       | 107        | MS/C4<br>MS/S3           | 92<br>93 |
| GP005<br>GP015                | 193<br>193 | MS/S3<br>MS/S4           | 93<br>93 |
| GP015<br>GP025                | 193        | M3/34<br>MT/1/6F         | 93<br>65 |
| GP025<br>GP035                | 193        | MT/1/6F                  | 142      |
| 1                             | 195        | MT/1/6F                  | 384      |
| IS/P/2H/H/SP                  | 216        | MT/LHNR/S                | 130      |
| IS/P/21/1//SP                 | 216        | MT/LHR/S                 | 129      |
| IS/P/4H/H/SP                  | 216        | MT100                    | 125      |
| IS/P/4H/SP                    | 216        | MT100CL                  | 126      |
| IV001                         | 221        | MT100CL<br>MT110         | 120      |
| 1/001                         | 221        | MT110CL                  | 120      |
| M                             | 221        | MT120                    | 126      |
| MC/DIV                        | 159        | MT120CL                  | 126      |
| MC/JK/1870                    | 159        | MT130                    | 128      |
| MC/JK/870                     | 159        | MT140                    | 128      |
| MC/L/1870/7S5D                | 157        | MT150                    | 128      |
| MC/L/1870/FS                  | 157        | MT220                    | 127      |
| MC/L/1870/FSD                 | 157        | MT220CL                  | 127      |
| MC/L/870/3S2D                 | 156        | MT230                    | 127      |
| MC/L/870/FS                   | 156        | MT230CL                  | 127      |
| MC/R/1870/7S5D                | 157        | MW005                    | 294      |
| MC/R/1870/FS                  | 157        | MW010                    | 315      |
| MC/R/1870/FSD                 | 157        | MW015                    | 315      |
| MC/R/870/3S2D                 | 156        | 0                        |          |
| MC/R/870/FS                   | 156        | OBL/BE0/GW/S             | 32       |
| MC/ST                         | 158        | OBL/BE1/GW/S             | 32       |
| MR110                         | 257        | OBL/BE3/GW/S             | 32       |
| MR120                         | 257        | OBL/GW0/SP/S             | 32       |
| MR130                         | 256        | OBL/GW1/SP/S             | 32       |
| MR140                         | 256        | OBL/GW3/SP/S             | 32       |
| MR140CL                       | 256        | OBP/SP/S                 | 30       |
| MR210                         | 246        | OCL/BE0/GW/S             | 33       |
| MR210CL                       | 246        | OCL/BE1/GW/S             | 33       |
| MR310                         | 246        | OCL/BE3/GW/S             | 33       |
| MR310CL                       | 246        | OCL/GW0/SP/S             | 33       |
| MR410                         | 246        | OCL/GW1/SP/S             | 33       |
| MR410CL<br>MRF/L/LD/CAM       | 246<br>254 | OCL/GW3/SP/S<br>OCP/SP/S | 33<br>30 |
| MRF/L/LD/CAM<br>MRF/L/LD/CL   | 254<br>254 | OCP/SP/S<br>OP/431       | 426      |
| MRF/L/LD/CL<br>MRF/L/OC       | 254<br>254 | OP/431<br>OP/432         | 426      |
| MRF/L/OC<br>MRF/S/LD/CAM      | 254<br>254 | OP/432<br>OP/631         | 426      |
| MRF/S/LD/CAM<br>MRF/S/LD/CL   | 254<br>254 | OP/631<br>OP/632         | 426      |
| MRF/S/LD/CL<br>MRF/S/OC       | 254        | P                        | 420      |
| MRW/D2DOC/CAM                 | 249        | PC/1000/ST               | 158      |
| MRW/D2DOC/CAM<br>MRW/D2DOC/CL | 249        | PC/300/ST                | 158      |
| MRW/D2DXP/CAM                 | 249        | PC/500/ST                | 158      |
| MRW/D2DXP/CL                  | 249        | PC/600/ST                | 158      |
|                               |            |                          |          |

|   | EC1350/VBB     | 171 |
|---|----------------|-----|
|   | EC2205/VBB     |     |
|   |                | 176 |
|   | EC2210/VBB     | 176 |
|   | EC2215/VBB     | 176 |
|   | EC2220/VBB     | 176 |
|   | EC2225/VBB     | 177 |
|   |                |     |
|   | EC2230/VBB     | 177 |
|   | EC2235/VBB     | 177 |
|   | EC2240/VBB     | 177 |
|   | EC2305/VBB     | 178 |
|   |                |     |
|   | EC2310/VBB     | 178 |
|   | EC2315/VBB     | 178 |
|   | EC2320/VBB     | 178 |
|   | EC2325/VBB     | 179 |
|   | EC2330/VBB     | 179 |
|   | EC2335/VBB     | 179 |
|   |                |     |
|   | EC2340/VBB     | 179 |
|   | EC4305/VBB     | 180 |
|   | EC4310/VBB     | 180 |
|   | EC4315/VBB     | 180 |
|   | EC4320/VBB     | 180 |
|   | EC4325/VBB     | 180 |
|   |                |     |
|   | EC4330/VBB     | 181 |
|   | EC4335/VBB     | 181 |
|   | EC4340/VBB     | 181 |
|   | EC4345/VBB     | 181 |
|   | EC5000/VBB     | 184 |
|   |                |     |
|   | EC5000/VBB     | 264 |
|   | EC5305/VBB     | 182 |
|   | EC5305/VBB     | 262 |
|   | EC5310/VBB     | 182 |
|   | EC5310/VBB     | 262 |
|   | EC5315/VBB     | 182 |
|   |                |     |
|   | EC5315/VBB     | 262 |
|   | EC5320/VBB     | 182 |
|   | EC5320/VBB     | 262 |
|   | EC5325/VBB     | 182 |
|   | EC5325/VBB     | 262 |
|   |                |     |
|   | EC5330/VBB     | 183 |
|   | EC5330/VBB     | 263 |
|   | EC5335/VBB     | 183 |
|   | EC5335/VBB     | 263 |
|   | EC5340/VBB     | 183 |
|   | EC5340/VBB     | 263 |
|   |                |     |
|   | EC5345/VBB     | 183 |
|   | EC5345/VBB     | 263 |
| G |                |     |
|   | G/200/RS/VBB   | 268 |
|   | G/205/FS/VBB   | 270 |
|   | G/305/FS/VBB   | 269 |
|   |                |     |
|   | G/500/RS/CYL   | 272 |
|   | G/500/RS/VBB   | 10  |
|   | G/500/RS/VBB   | 271 |
|   | G/CARD/LOCK    | 273 |
|   | G/CARD/STATION | 273 |
|   | G/CLS          | 272 |
|   | G/CLS/500      |     |
|   |                | 272 |
|   | G/DP           | 272 |
|   |                |     |

| 2         DTSF/900/2/SD         424           2         DTSF/96/2/SD         118           2         DTSR/450/SR         114           0         DTSR/450/SR         114           0         DTSR/450/SR         114           0         DTSR/600/SR         114           1         DTSR/600/SR         114           3         DTSR/50/SR         114           0         DTSR/750/SR         114           0         DTSR/900/SR         114           0         DUS/R 119         DUL/R           0         DUS/F         119           0         E04         E04/U           185         E06         185           136         E07         185           14         E11         186           15         187         9           16         E11         186  |                                 |   |   |
|---|---------------------------------|---|---|
| E         E/SOCKET/4         308           7         E02         185           7         E04 - E04/U         185           9         E05         185           60         185         E06           60         185         E07           8         E07         185           9         E06         185           60         185         E07           8         E09         186           0         E11         186           0         E12         186           0         E13         186           0         E14         187           9         E21         187           9         E21         187           9         E21         187           9         E21         187           9         E8C020         18           8         EBC020         18           8         EBC020         18           8         EC1200/VBB         172           5         EC120/VBB         168           5         EC120/VBB         168           5         EC1220/VBB         168   | 2<br>2<br>0<br>0<br>1<br>1<br>3 | DTSF/96/2/SD<br>DTSF/96/2/SD<br>DTSR/450<br>DTSR/450/SR<br>DTSR/600/SR<br>DTSR/600/SR<br>DTSR/750<br>DTSR/900/SR<br>DTSR/900<br>DTSR/900/SR<br>DUL/F<br>DUL/R | 118<br>425<br>114<br>114<br>114<br>114<br>114<br>114<br>114<br>114<br>114<br>11 |
| E/SOCKET/4         308           F         E02         185           E03         185         E03         185           F         E04 - E04/U         185         E06         185           E07         185         E06         185           E07         185         E06         185           E07         185         E08         186           E09         186         E09         186           0         E11         186         13         186           0         E14         187         9         E13         187           9         E19         187         187         9         E21         187           9         E19         187         188           | 6                               |   | 119   |
| 9         E05         185           E06         185           E07         185           E08         186           E09         186           0         E10           0         E11           1         186           0         E12           1         E13           1         E13           1         E14           1         E13           1         E14           1         E15           1         E14           1         E15           1         E14           1         E17           1         E17           1         E18           1         E18           1         E13           18         E14           187         9           E14         187           9         E21           188         EBC015           18         EEC025           182         EC1205/VBB           188         EC120/VBB           188         EC1220/VBB           188         EC1220/VBB  |                                 | E/SOCKET/4<br>E02   | 185   |
| 0         E10         186           0         E11         186           0         E12         186           1         E13         186           0         E14         187           0         E14         187           0         E13         187           9         E13         187           9         E21         187           9         EBC015         18           8         EBC020         18           8         EBC030         18           8         EC010/VBB         188           8         EC010/VBB         188           5         EC120/VBB         168           5         EC120/VBB         168           5         EC120/VBB         168           5         EC1220/VBB         168           5         EC1230/VBB         169           6         EC1235/VBB         169           7         EC1305/VBB         169           6         EC1305/VBB         170           8         EC1310/VBB         170           8         EC1310/VBB         170           6 <td< td=""><td></td><td>E05<br/>E06<br/>E07<br/>E08</td><td>185<br/>185<br/>185<br/>185<br/>185<br/>186</td></td<> |                                 | E05<br>E06<br>E07<br>E08  | 185<br>185<br>185<br>185<br>185<br>186  |
| 0         E15         187           9         E18         187           9         E19         187           9         E21         187           9         E21         187           9         E8015         18           8         EBC020         18           8         EBC020         18           8         EBC020         18           8         EBC020         18           8         EBC030         18           8         EC010/VBB         172           5         EC120/VBB         168           5         EC120/VBB         168           5         EC120/VBB         168           5         EC125/VBB         168           5         EC1240/VBB         168           5         EC1240/VBB         169           6         EC1240/VBB         169           6         EC1240/VBB         169           6         EC1305/VBB         170           8         EC1310/VBB         170           8         EC1320/VBB         170           5         EC1320/VBB         170           6  | 0<br>0<br>1                     | E10<br>E11<br>E12<br>E13  | 186<br>186<br>186<br>186  |
| 9         EBC015         18           8         EBC020         18           8         EBC020         18           8         EBC030         18           8         EBC030         18           8         EBC030         18           8         EC010/VBB         172           5         EC120/VBB         168           5         EC1210/VBB         168           5         EC1225/VBB         168           5         EC1225/VBB         168           5         EC1225/VBB         169           6         EC125/VBB         169           6         EC125/VBB         169           6         EC1250/VBB         169           6         EC130/VBB         170           8         EC130/VBB         170           8         EC1310/VBB         170           8         EC1320/VBB         170           8         EC1320/VBB         170           8         EC1320/VBB         170           8         EC1320/VBB         170           8         EC1330/VBB         171           9         EC1330/VBB         171   | 0<br>9                          | E15<br>E18  | 187<br>187  |
| 8         EBC030         18           8         EC010/VB8         172           5         EC120/VB8         168           5         EC1210/VB8         168           5         EC125/VB8         168           5         EC1220/VB8         168           5         EC1225/VB8         168           5         EC1225/VB8         169           6         EC1230/VB8         169           6         EC1240/VB8         169           7         EC1305/VB8         170           8         EC1310/VB8         170           5         EC1320/VB8         170           5         EC1325/VB8         170           6         EC1330/VB8         171           7         EC1330/VB8         171           6         EC1340/VB8         171   | 9<br>8                          | EBC015<br>EBC020  | 18<br>18  |
| 5         EC1215/VBB         168           5         EC1220/VBB         168           6         EC1220/VBB         168           5         EC1220/VBB         168           5         EC1230/VBB         169           6         EC1240/VBB         169           7         EC1240/VBB         169           6         EC1240/VBB         169           6         EC1250/VBB         169           7         EC1305/VBB         170           8         EC1316/VBB         170           8         EC1320/VBB         170           6         EC1326/VBB         170           6         EC1330/VBB         171           6         EC1330/VBB         171           6         EC1330/VBB         171   | 8<br>8<br>5                     | EBC030<br>EC010/VBB   | 18<br>172   |
| 6         EC1235/VBB         169           7         EC1240/VBB         169           4         EC1240/VBB         169           5         EC1250/VBB         169           7         EC1305/VBB         170           4         EC1310/VBB         170           5         EC1320/VBB         170           6         EC1320/VBB         170           7         EC1330/VBB         170           6         EC1330/VBB         170           6         EC1330/VBB         171           6         EC1340/VBB         171   | 5<br>5                          | EC1215/VBB<br>EC1220/VBB  | 168<br>168  |
| 6         EC1250/VBB         169           7         EC1305/VBB         170           4         EC1310/VBB         170           5         EC1315/VBB         170           6         EC1315/VBB         170           7         EC1320/VBB         170           7         EC1320/VBB         170           7         EC1330/VBB         171           6         EC1336/VBB         171  | 6<br>7                          | EC1235/VBB<br>EC1240/VBB  | 169<br>169  |
| 5         EC1320/VBB         170           6         EC1325/VBB         170           7         EC1330/VBB         171           4         EC1336/VBB         171           6         EC1340/VBB         171  | 6<br>7<br>4                     | EC1250/VBB<br>EC1305/VBB  | 169<br>170  |
| 4 EC1335/VBB 171<br>6 EC1340/VBB 171  | 5<br>6                          | EC1320/VBB<br>EC1325/VBB  | 170<br>170  |
|   | 4<br>6                          | EC1335/VBB<br>EC1340/VBB  | 171<br>171  |

| CYHP/E1<br>CYT/CD1<br>CYT/D12<br>CYT/F1<br>CYT/FG1/ST<br>CYT/FG1/ST<br>CYT/J1/ST<br>CYWP/1  | 272<br>102<br>102<br>100<br>100<br>100<br>101<br>101   |
|---|--|
| D3<br>D5<br>D4F/1<br>DAR/1<br>DAR/1<br>DAR/1<br>DAR/1<br>D7<br>D7<br>D7<br>D7<br>D7<br>D7<br>D7<br>D7<br>D7<br>D7<br>D7<br>D7<br>D7 | 91<br>91<br>106<br>72<br>106<br>91<br>107<br>73<br>107<br>73<br>107<br>73<br>107<br>73<br>89<br>89<br>89<br>220<br>220<br>220<br>221<br>220<br>221<br>220<br>221<br>220<br>221<br>220<br>221<br>220<br>221<br>219<br>2219<br>2 |
|   |  |

D

www.bristolmaid.com

| SAC/WP  | 162  |
|---|--|
| SB/221  |  |
|   | 164  |
| SB/321  | 164  |
| SB/421  | 164  |
|   |  |
| SBH/15  | 73   |
| SBH/15  | 435  |
| SBH/WM  | 109  |
|   |  |
| SBKSL/1   | 423  |
| SBKSLH/1  | 423  |
| SBS/2   | 420  |
|   |  |
| SBSA/1  | 421  |
| SBSA/2  | 421  |
| SBSL/1  | 423  |
| SBST/2  |  |
|   | 420  |
| SC/DIV  | 411  |
| SC/JK/1870  | 411  |
| SC/L/1870/7S5D  | 410  |
|   |  |
| SC/L/1870/FS  | 410  |
| SC/L/1870/FSD   | 410  |
| SC/R/1870/7S5D  | 410  |
|   |  |
| SC/R/1870/FS  | 410  |
| SC/R/1870/FSD   | 410  |
| SC020   | 412  |
|   |  |
| SC025   | 412  |
| SC040   | 412  |
| SC045   | 412  |
| SEP/DID   | 384  |
|   |  |
| SEP/DIS   | 384  |
| SEP/L   | 385  |
| SEP/LH  | 385  |
|   |  |
| SEP/SBH   | 383  |
| SEP/TIG01   | 385  |
| SEP/USBBMP  | 383  |
| SEPD/6/EPB  |  |
|   | 382  |
| SEPD/6/HSB  | 382  |
| SEPS/6/EPB  | 382  |
| SEPS/6/HSB  | 382  |
|   |  |
|   |  |
| SKC100  | 415  |
| SKC100<br>SKC100/CL   |  |
| SKC100/CL   | 415<br>415   |
| SKC100/CL<br>SKC20  | 415<br>415<br>415  |
| SKC100/CL<br>SKC20<br>SKC20/CL  | 415<br>415<br>415<br>415   |
| SKC100/CL<br>SKC20  | 415<br>415<br>415  |
| SKC100/CL<br>SKC20<br>SKC20/CL  | 415<br>415<br>415<br>415   |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30  | 415<br>415<br>415<br>415<br>415<br>414<br>415                      |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL  | 415<br>415<br>415<br>415<br>414<br>415<br>415<br>415               |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48   | 415<br>415<br>415<br>415<br>414<br>415<br>415<br>415               |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL  | 415<br>415<br>415<br>415<br>414<br>415<br>415<br>415               |
| SKC100/CL<br>SKC20/CL<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48<br>SKC48/CL  | 415<br>415<br>415<br>415<br>414<br>415<br>415<br>415<br>415        |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48<br>SKC48/CL<br>SKC64  | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30/CL<br>SKC48<br>SKC48/CL<br>SKC48/CL<br>SKC64/CL  | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48<br>SKC48/CL<br>SKC64  | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48/CL<br>SKC48/CL<br>SKC64<br>SKC64/CL<br>SKC64/CL<br>SKC64/CL   | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48<br>SKC48/CL<br>SKC64<br>SKC64/CL<br>SR<br>SR<br>SS/22/1   | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48/CL<br>SKC64/CL<br>SKC64/CL<br>SR<br>SS/22/1<br>SS/22/2  | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48/CL<br>SKC48/CL<br>SKC48/CL<br>SKC48/CL<br>SKC4/CL<br>SKC4/CL<br>SKC4/CL<br>SKC4/CL<br>SK22/1<br>SS/22/2<br>SS/36                            | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48/CL<br>SKC64/CL<br>SKC64/CL<br>SR<br>SS/22/1<br>SS/22/2  | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30<br>SKC48/CL<br>SKC48/CL<br>SKC44/CL<br>SKC64<br>SKC64/CL<br>SKC64<br>SKC64/CL<br>SK22/2<br>SS/22/1<br>SS/22/2<br>SS/36<br>SS/40/MS                      | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48<br>SKC48/CL<br>SKC64<br>SKC64/CL<br>SKC64/CL<br>SS/22/1<br>SS/22/2<br>SS/22/2<br>SS/36<br>SS/40/MS<br>SS/48/3                               | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/TRAKA<br>SKC30/CL<br>SKC30/CL<br>SKC48<br>SKC48/CL<br>SKC48/CL<br>SKC48/CL<br>SKC48/CL<br>SKC4/CL<br>SKC4/CL<br>SK/22/1<br>SS/22/1<br>SS/22/2<br>SS/36<br>SS/40/MS<br>SS/40/MS<br>ST/340/UN | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/CL<br>SKC21/TRAKA<br>SKC30<br>SKC30/CL<br>SKC48<br>SKC48/CL<br>SKC64<br>SKC64/CL<br>SKC64/CL<br>SS/22/1<br>SS/22/2<br>SS/22/2<br>SS/36<br>SS/40/MS<br>SS/48/3                               | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |
| SKC100/CL<br>SKC20<br>SKC20/TRAKA<br>SKC30/CL<br>SKC30/CL<br>SKC48<br>SKC48/CL<br>SKC48/CL<br>SKC48/CL<br>SKC48/CL<br>SKC4/CL<br>SKC4/CL<br>SK/22/1<br>SS/22/1<br>SS/22/2<br>SS/36<br>SS/40/MS<br>SS/40/MS<br>ST/340/UN | 415<br>415<br>415<br>414<br>415<br>415<br>415<br>415<br>415<br>415 |

| PL125/BD       | 417 | R/C/T                     | 252 |
|----------------|-----|---------------------------|-----|
| PL130/BD       | 417 | R/C/TH                    | 252 |
| PL135/BD       | 417 | RC/2/15                   | 110 |
| PL140/BD       | 417 | RC/2/15                   | 110 |
| PL145/BD       | 417 | RC/5/15                   | 110 |
|                |     |                           |     |
| PL150/BD       | 417 | RC/5/15                   | 110 |
| PL900          | 417 | RD/B5/SS                  | 111 |
| PL905          | 417 | RD/U5/SS                  | 111 |
| PL910          | 417 | RD1                       | 91  |
| PL915          | 417 | RD2L                      | 91  |
| PL920          | 417 | RD2R                      | 91  |
| PL925          | 417 | RD3                       | 91  |
| PMH005         | 106 | RFID125/MC                | 25  |
|                |     |                           |     |
| PS/M/HC/4H     | 214 | RFID125/PC                | 25  |
| PS/M/SC/4H     | 214 | RFID125/SC                | 25  |
| PS/P/HC/4H     | 212 | RFID125/UC                | 25  |
| PS/P/SC/4H     | 212 | RFID1356/PC               | 25  |
| PT/LHNR/S      | 140 | RFID1356/SC               | 25  |
| PT/LHR/S       | 140 | RFID1356/SF               | 25  |
| PT050          | 402 | RFID1356/UC               | 25  |
| PT030          | 402 | RFID1356/UF               |     |
|                |     |                           | 25  |
| PTD/EPB/A18C12 | 135 | RPT                       | 23  |
| PTD/EPB/BD69   | 138 | S                         |     |
| PTD/EPB/BP9    | 136 | SA005                     | 143 |
| PTD/EPB/LP24   | 135 | SA025                     | 143 |
| PTD/EPB/MX60   | 139 | SA040                     | 163 |
| PTD/EPB/NC56   | 137 | SA040/CL                  | 163 |
| PTD/EPB/PA32   | 133 | SA045                     | 163 |
| PTD/HSB/A18C12 | 135 | SA045/CL                  | 163 |
|                |     |                           |     |
| PTD/HSB/BD69   | 138 | SAC/213/B                 | 161 |
| PTD/HSB/BP9    | 136 | SAC/213/B/CL              | 161 |
| PTD/HSB/LP24   | 135 | SAC/213/B/CL/VP           | 161 |
| PTD/HSB/MX60   | 139 | SAC/213/B/VP              | 161 |
| PTD/HSB/NC56   | 137 | SAC/213/S                 | 160 |
| PTD/HSB/PA32   | 133 | SAC/213/S/CL              | 160 |
| PTS/115        | 403 | SAC/213/S/CL/VP           | 160 |
| PTS/115/M      | 403 | SAC/213/S/VP              | 160 |
| PTS/EPB/BD46   | 138 | SAC/213/3/VF<br>SAC/313/B | 161 |
|                |     |                           |     |
| PTS/EPB/BP6    | 136 | SAC/313/B/CL              | 161 |
| PTS/EPB/DR6    | 134 | SAC/313/B/CL/VP           | 161 |
| PTS/EPB/LA16   | 135 | SAC/313/B/VP              | 161 |
| PTS/EPB/MX40   | 139 | SAC/313/S                 | 160 |
| PTS/EPB/PA24   | 132 | SAC/313/S/CL              | 160 |
| PTS/HSB/BD46   | 138 | SAC/313/S/CL/VP           | 160 |
| PTS/HSB/BP6    | 136 | SAC/313/S/VP              | 160 |
| PTS/HSB/DR6    | 134 | SAC/413/B                 | 161 |
|                |     |                           |     |
| PTS/HSB/LA16   | 135 | SAC/413/B/CL              | 161 |
| PTS/HSB/MX40   | 139 | SAC/413/B/CL/VP           | 161 |
| PTS/HSB/PA24   | 132 | SAC/413/B/VP              | 161 |
|                |     | SAC/413/S                 | 160 |
| R/BAG/BLUE     | 229 | SAC/413/S/CL              | 160 |
| R/BAG/GREEN    | 229 | SAC/413/S/CL/VP           | 160 |
| R/BAG/GREY     | 229 | SAC/413/S/VP              | 160 |
| R/BAG/ORANGE   | 229 | SAC/FS/2                  | 162 |
| R/BAG/ORANGE   | 229 | SAC/FS/2<br>SAC/FS/3      | 162 |
|                |     |                           |     |
| R/BAG/WHITE    | 229 | SAC/FS/4                  | 162 |
| R/BAG/YELLOW   | 229 | SAC/LT                    | 162 |
|                |     |                           |     |

+44 (0)1258 484455

R

| PC636/LH/LX          | 146 |  |
|----------------------|-----|--|
| PC636/LX             | 146 |  |
| PC639                | 146 |  |
| PC639/FS             | 148 |  |
| PC639/FS/LH          | 148 |  |
| PC639/FS/LH/LX       | 148 |  |
| PC639/FS/LX          | 148 |  |
| PC639/LH             | 146 |  |
| PC639/LH/LX          |     |  |
|                      | 146 |  |
| PC639/LH/MDS         | 150 |  |
| PC639/LH/MDS/LX      | 150 |  |
| PC639/LX             | 146 |  |
| PC639/MDS            | 150 |  |
| PC639/MDS/LX         | 150 |  |
| PC836                | 147 |  |
| PC836/FS             | 149 |  |
| PC836/FS/IND         | 149 |  |
| PC836/FS/IND/LX      | 149 |  |
| PC836/FS/LX          | 149 |  |
| PC836/IND            | 147 |  |
| PC836/IND/LX         | 147 |  |
| PC836/LX             | 147 |  |
| PC839                | 147 |  |
| PC839/FS             | 147 |  |
|                      |     |  |
| PC839/FS/IND         | 149 |  |
| PC839/FS/IND/LX      | 149 |  |
| PC839/FS/LX          | 149 |  |
| PC839/IND            | 147 |  |
| PC839/IND/LX         | 147 |  |
| PC839/LX             | 147 |  |
| PC905                | 159 |  |
| PCB                  | 159 |  |
| PE/30                | 164 |  |
| PL005/BD             | 417 |  |
| PL010/BD             | 417 |  |
| PL015/BD             | 417 |  |
| PL020/BD             | 417 |  |
| PL020/BD<br>PL025/BD | 417 |  |
| PL025/BD<br>PL030/BD | 417 |  |
|                      |     |  |
| PL035/BD             | 417 |  |
| PL040/BD             | 417 |  |
| PL045/BD             | 417 |  |
| PL050/BD             | 417 |  |
| PL055/BD             | 417 |  |
| PL060/BD             | 417 |  |
| PL065/BD             | 417 |  |
| PL070/BD             | 417 |  |
| PL075/BD             | 417 |  |
| PL080/BD             | 417 |  |
| PL085/BD             | 417 |  |
| PL090/BD             | 417 |  |
| PL095/BD             | 417 |  |
| PL095/BD<br>PL100/BD |     |  |
|                      | 417 |  |
| PL105/BD             | 417 |  |
| PL110/BD             | 417 |  |
| PL115/BD             | 417 |  |
| PL120/BD             | 417 |  |
|                      |     |  |

| PC/800/ST         | 158 |
|-------------------|-----|
| PC/RLBF           | 159 |
| PC/RLBS           | 159 |
| PC1036            | 147 |
| PC1036/FS         | 149 |
| PC1036/FS/IND     | 149 |
| PC1036/FS/IND/LX  | 149 |
| PC1036/FS/LX      | 149 |
| PC1030/1 3/LX     | 145 |
|                   |     |
| PC1036/IND/LX     | 147 |
| PC1036/LX         | 147 |
| PC1039            | 147 |
| PC1039/FS         | 149 |
| PC1039/FS/IND     | 149 |
| PC1039/FS/IND/LX  | 149 |
| PC1039/FS/LX      | 149 |
| PC1039/IND        | 147 |
| PC1039/IND/LX     | 147 |
| PC1039/IND/MDS    | 151 |
| PC1039/IND/MDS/LX |     |
|                   | 151 |
| PC1039/LX         | 147 |
| PC1039/MDS        | 151 |
| PC1039/MDS/LX     | 151 |
| PC336/FS          | 148 |
| PC336/FS/LH       | 148 |
| PC336/FS/LH/LX    | 148 |
| PC336/FS/LX       | 148 |
| PC536             | 146 |
| PC536/FS          | 148 |
| PC536/FS/LH       | 148 |
| PC536/FS/LH/LX    | 148 |
| PC536/FS/LX       | 148 |
| PC536/LH          | 146 |
|                   |     |
| PC536/LH/LX       | 146 |
| PC536/LH/MDS      | 150 |
| PC536/LH/MDS/LX   | 150 |
| PC536/LX          | 146 |
| PC536/MDS         | 150 |
| PC536/MDS/LX      | 150 |
| PC539             | 146 |
| PC539/FS          | 148 |
| PC539/FS/LH       | 148 |
| PC539/FS/LH/LX    | 148 |
| PC539/FS/LX       | 148 |
|                   |     |
| PC539/LH          | 146 |
| PC539/LH/LX       | 146 |
| PC539/LH/MDS      | 150 |
| PC539/LH/MDS/LX   | 150 |
| PC539/LX          | 146 |
| PC539/MDS         | 150 |
| PC539/MDS/LX      | 150 |
| PC636             | 146 |
| PC636/FS          | 148 |
| PC636/FS/LH       | 148 |
| PC636/FS/LH/LX    | 148 |
|                   |     |
| PC636/FS/LX       | 148 |
| PC636/LH          | 146 |
|                   |     |

|   | ŝ  |  |
|---|----|--|
|   |    |  |
| 6 | s' |  |
| 2 |    |  |
| 4 | Þ  |  |
|   | ł  |  |
|   | -  |  |
|   |    |  |

| ST/690/UN/T   | 404   |
|---|---|
| STSS/690/UN   | 406   |
| STSS/690/UN/T   | 406   |
| T<br>TBIN4/MS<br>TBIN4/MS<br>TBIN5/MS<br>TBIN5/MS<br>TBIN5/MS<br>TBIN5/MS<br>TBIN6/MS<br>TBIN6/MS<br>TBIN6/MS<br>TBIN6/SS<br>TBIN6/SS   | 66<br>311<br>66<br>66<br>311<br>66<br>66<br>311<br>66<br>23   |
| UB/9  | 109   |
| UB05  | 109   |
| UB010   | 109   |
| VH/950  | 96  |
| VH/950/LBT  | 96  |
| WT/11/MS<br>WT/13/MS<br>WT/33/SS/2<br>WT/33/SS/2<br>WT/33/SS/2<br>WT/33/SS/3/ZP005<br>WT/37/2/ZP010<br>WT/37/2/ZP010<br>WT/37/3/ZP010<br>WT/37/3/ZP010<br>WT/42/WT/COVER<br>WT/42/WT/COVER<br>WT/42/COVER/TBAR<br>WT/48/COVER/TBAR<br>WT/48/FBD/TBAR<br>WT/48/FBD/TBAR<br>WT/48/FBD/TBAR<br>WT/50/COVER/TBAR<br>WT/50/COVER<br>WT/50/COVER<br>WT/50/COVER<br>WT/50/FBD/TBAR<br>WT/50/FBD/TBAR<br>WT/50/FBD/TBAR<br>WT/50/FBD/TBAR<br>WT/50/FBD/TBAR<br>WT/50/FBD/TBAR<br>WT/50/FBD/TBAR<br>WT/50/FBD/TBAR<br>WT/50/FBD/TBAR<br>WT/6/2/MS<br>WT/6/2/PL/SS<br>WT/6/2/SS<br>WT/6/PL/SS<br>WT/6/PL/SS | 228<br>227<br>191<br>191<br>191<br>190<br>190<br>230<br>230<br>230<br>230<br>234<br>234<br>234<br>234<br>234<br>234<br>234<br>232<br>232<br>232 |

 WT/9/2/PL/MS
 226

 WT/9/PL/MS
 226

 WT/9/PL/MS
 226

 WT/SDPL/MS
 226

 WT/BOARD
 307

 WT/CBOARD
 307

 WT/CYL/DLF
 307

 WT/CPL/DF
 307

 WT/DEFIB
 306

 WT/SPS
 306

 WT200
 304

 WT250
 304

ZP005 121 ZP010 76

Ζ

# **PAINT & LAMINATE**

# Paint & Laminate

Painted surfaces should be cleaned using a soft damp cloth & polished dry. Do not use abrasive cloth. Persistent exposure to wet or humid conditions could result in rust marks, or cause the paint to lift. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.



# Laminate & Melamine

Laminates & melamine should be cleaned using a soft damp cloth & polished dry. Do not use abrasive cloth. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.




# VINYL

Vinyl - Ambla New Richmond in Colour



# Vinyl - Spradling Valencia<sup>™</sup> C5



# Faux Leather - Memphis



# **Care & Cleaning**

Vinyl's are resistant to most mild acids, alkalis, drinks & household stains. Clean with a damp, soapy cloth & rinse well with clean water. Dry the surface thoroughly. A soft brush with soapy water may be used to remove stubborn stains. Ensure surfaces are then rinsed & dried. Do not use solvents, bleaches\*, abrasives, synthetic detergents, wax polishes or aerosol sprays. Disinfectants may be used in dilution as specified by their manufacturers.


\* Chlorine Releasing Agents: Levels of 1-10,000 parts per million should pose no problems as long as the material does not endure long soakings at the higher concentrations. Rinse well with clean water. These include Haz-Tabs & other disinfectants commonly used in hospital infection control programmes. Surfaces are NaDCC compatible.

Swatches

# **FABRIC & PLASTIC**

# Fabric - Camira Advantage

Vacuum regularly. Wipe clean with a damp cloth or shampoo using proprietary upholstery shampoo.



# Fabric - VP Multi-Way Stretch

Wipe with mild detergent. Wash or dry clean using detergent up to 75°C. For disinfection use sodium hypo-chlorite at no more than 1,000 parts per million of active chlorine. Do not use biological or phenolic-based cleaners. Wash & dry completely. Some surface wrinkling may occur.



# FABRIC - Camira 24/7+

Vacuum regularly. Professionally dry clean or wipe clean with a damp cloth.



FYR02 - Vernal

FYR03 - Equinox

FYR21 - While

Swatches

# **Plastic Seat Shells**

Plastic should be cleaned using a soft damp cloth & polished dry. Do not use abrasive cloth. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.



# Panaz Patternistas for Childern

Vacuum regularly. Wipe clean with a damp cloth. Can be cleaned with water based proprietary cleaners & 5% sodium chlorite bleach solution, Haz-Tabs or Chlor-Clean.



SAVL - Savanna Lavender



# **TECHNICAL INFORMATION**

#### **Dimensions & Weight**

- Dimensions are stated in millimetres (width x depth x height)
- · Weights are stated in kilograms
- · Please ensure you are aware of the item size & weight prior to ordering

#### Castors

- All castors are fitted with non-marking tyres (with exception of some anti-static versions)
- Anti-static castors can be identified by either a yellow mark on the edge of the tyre, or where twin wheel castors are used, a different colouration to the tyre (the result of the addition of a conductive material)
- Castors & wheel bearings do not require specific maintenance other than regular cleaning
- The use of oil will breakdown the grease applied to the bearings, reducing the life of the castor
- We recommend that fittings are checked regularly & re-tightened as required
- Where anti-static castors are fitted, we recommend that the tyre is periodically cleaned to remove any dirt that may have built up. This will ensure a good contact between the tyre & the floor

# Locks & Keys

- Locks should operate smoothly. If they do not please arrange urgent attention, & where necessary, obtain a replacement
- Locks require minimal maintenance & should only be lubricated using a dry lubricant such as graphite which can be purchased from most locksmiths (NEVER USE OIL)
- When ordering replacement locks or keys, please ensure that the product code & any letters / numbers from the key & lock are available

#### Paint

- · We only use lead free epoxy polyester powder paint
- Where possible colours are applied in special booths to allow for overspray to be recycled
- Alternative colours are available at additional cost

#### **Stainless Steel**

- We use 304 & 430 grade stainless steel
  - · Grade 304 is non-magnetic & used on products which are not painted
  - · Grade 430 is magnetic & will be subsequently painted
- Please note that the term 'stainless' is not entirely accurate, as certain chemicals can 'eat' into the surface. In particular, 'iodine' which can badly corrode / stain, & as a result should not be used as a cleaning agent. Stainless steel products should always be dried after cleaning to reduce staining

#### Wood, Melamine & Laminate

- We only use FSC accredited materials. This provides assurances that the material is produced from well-managed forests & / or recycled materials
- We use a combination of melamine faced MDF / chipboard & high pressure laminates bonded MDF of chipboard substrates
- Edges are protected with 2mm ABS, T-Barb edging or sealed with varnish
- · Both melamine & laminate panels offer some resistance to water penetration under normal usage
- · We only use water based adhesives & varnishes

#### **Upholstered Panels & Fabrics**

- Vinyl upholstered panels are manufactured from 'Flame Retardant' materials, independently tested to BS 5852: 2006 Ignition Source 5, BS EN 1021-1: 2006, BS EN 1021-2: 2006 & BS 7176: 2007 Medium Hazard. Standard vinyl colours also contain an anti-microbial compound within its construction, which constantly migrates to the surface to replenish any material which is removed by cleaning or general usage
- Fabric upholstered panels are manufactured from 'Flame Retardant' materials, independently tested to BS 5852: 2006 Ignition Source 5, BS EN 1021-1:2006, BS EN 1021-2:2006 & BS 7176: 2007 Low Hazard
- Curtain screens are manufactured from 'Flame Retardant' materials, independently tested to BS 5438: Part 2

# **Caretray Trolley™ Range**

The Caretray Trolley was developed specifically to meet the needs of the Healthcare Professional across a wide variety of disciplines including Minor Treatment, A&E, Theatre, Phlebotomy, ITU, and Outpatients to name but a few. It has become a flexible, dependable partner, a trusted workhorse in hospitals throughout the UK and Overseas.

In its 25th year of production, the Caretray Trolley offers excellent value for money & is available as either pre-configured or can be specified to meet your specific requirements...

#### **Pre-Configured Trolleys**

- Chart Trolleys
- Dispensing Trolleys
- Phlebotomy Trolleys
- Procedure Trolleys
- Resuscitation Trolleys
- Sepsis Trolleys
- Laptop Trolleys

#### Configure Your Own Trolley (Mild & Stainless Steel)

- Low Level
- Standard Level
- High Level

#### Configure Your Own Rack (Mild & Stainless Steel)

- Mobile
- Static



#### **Electric Columns & Actuators**

 Electric columns & actuators comply with EMC Directive 2014/30/EU, Low Voltage Directive 2014/35/EU & RoHS2 Directive 2011/65/EU

#### **Medical Devices**

Where applicable, Hospital Metalcraft product is marked with the CE symbol to demonstrate compliance with the essential requirements of any European Directives (& transposing regulations) that may apply & for which a Declaration of Conformity is held. Hospital Metalcraft is registered with the MHRA (Medicines & Healthcare Regulatory Authority) in connection with product classified as 'Class 1 Medical Devices' for purposes of the Medical Devices Directive (93/42/EEC) with the registration number CA001358.

#### **WEEE Regulations**

Hospital Metalcraft meets its obligations as a business to business producer of Electrical and Electronic Equipment (EEE) through:

- Membership of the Veolia Environmental Services WEEE, (Waste Electrical & Electronic Equipment), Producer Compliance Scheme (Registration No. WEE/KE0677RX)
- Marking EEE placed on the market by Hospital Metalcraft with the mandatory crossed out wheeled bin symbol & Bristol Maid<sup>™</sup> logo
- Making information available, (for EEE placed on the market after April 2007), to assist those engaged in treatment, recovery & recycling activities
- Arranging for the collection, treatment, recovery & environmentally sound disposal of obligated WEEE

#### **Preventative Maintenance, Spares & Service Support**

Bristol Maid<sup>™</sup> products are designed & manufactured to cater for the demands of the healthcare environment, with the majority (under normal conditions) providing many years of trouble free service. Where annual maintenance or periodic inspections are required, we can offer comprehensive planned preventative maintenance through either a scheduled service contract or call out by our own trained Engineers. For your peace of mind & continued product reliability, please contact our customer service team for details on the following services.

- · Preventative maintenance contracts & packages
- Supply of spare parts
- · Site audit & equipment recommendations

#### **Cleaning Guidance & Instructions**

The maintenance of high standards of cleanliness on all surfaces & equipment is a crucial factor in the prevention & control of infection. Whilst all dust, dirt & moisture can harbour infection, the key risks are associated with contamination arising from contact with blood, body fluids, mucous membranes & damaged skin. In all cases surfaces, must be thoroughly cleaned & disinfected in order to destroy any pathogenic micro-organisms.

Bristol Maid products have been designed for use in demanding environments. Guidance & instructions can be found at - www.bristolmaid.com/cleaning-instructions





# Manufacturing Medical Furniture & Equipment Since 1953

With over 65 years of experience, Bristol Maid is a leading manufacturer and supplier of medical furniture & equipment.

From design and manufacture to supply and installation, Bristol Maid has a reputation for providing high product quality, excellent customer service and value for money.

We take pride in working in close partnership with our customers applying innovative design and technology to deliver effective solutions for today while anticipating the needs of tomorrow.





Hospital Metalcraft Limited Blandford Heights, Blandford Forum, Dorset. DT11 7TG, UK

Tel: +44 (0) 1258 451338 Fax: +44 (0) 1258 455056 sales@bristolmaid.com

www.bristolmaid.com







# **Family business**

Hospital Metalcraft Ltd is a family company. We own our factory, production lines and delivery lorries. We continually invest into our staff and the development of new product ranges.



# Healthcare focused

With an extensive portfolio that covers almost every discipline, we have solutions to cater for the demands of the healthcare environment.



# Manufactured in the UK

Hospital Metalcraft is proud to be a British manufacturing success story. Strict quality criteria are applied to both the Bristol Maid product range designed and produced within our own manufacturing facility, and complementary equipment from carefully selected partners.



# Delivered using our own lorries

Wherever possible we use our own fleet of lorries to ensure that your equipment arrives safe and secure.



# The best quality

With over 65 years' experience, we are focussed on the design and manufacture of high-quality products which will create value for our customers. With a worldwide customer base, the Bristol Maid brand is recognised by many thousands of satisfied customers.



# **Approved supplier**

We are an accredited NHS & public sector accredited supplier. This allows you to buy with confidence knowing your requirements for compliant and competitive tendering has already been met.

# Support you can rely on

Our external sales representatives, sales support engineers and our customer service team are here to help you at every stage, from advice with product selection, managing the delivery of your order through to full service and repair support.

# **HOW TO ORDER**

#### How to Order

- Orders are required in writing. Should you wish to order over the phone, payment is required in full (debit or credit card). Alternatively, we can issue a pro-forma invoice, payment is required in full prior to manufacture
- Orders should be submitted as an 'official purchase order' on company headed paper with all information clearly displayed
- Your order should contain the following information
  - Item No.
  - Description
  - Quantity
  - Delivery address
  - Invoice address
  - Full contact details, including name, email address, phone & fax numbers
- · Where a quotation has been issued, please reference this to ensure any agreed discounts are applied
- Item codes & descriptions in this catalogue supersede all previous editions
- · A full price list can be requested from our website
- Payments can be made by cheque, bank transfer, credit or debit card
- If you are a new customer, the first three orders will be processed on a pro-forma basis, paid in full prior to manufacture & delivery, after which a credit account may be applied for
- · For our terms and conditions refer to www.bristolmaid.com
- For further information and assistance, please contact our customer services team on +44 (0)1258 484455

# **Payment Options**

#### BACS

- · Payments to: HSBC Bank PLC
- Sort Code: 40-12-18
- Account Number: 11044052

#### Credit or Debit Card

• Major credit cards accepted (with exception of American Express)







7







# Bristol Maid Hospital Metalcraft<sup>™</sup> Limited Product Catalogue 2022









# **HOW TO ORDER**

#### **NHS & Public Sector Accredited Supplier**

Public sector organisations can take advantage of our nationally agreed discounts through the NHS Supply Chain and Shared Business Services. Buy with confidence knowing your requirement for compliant and competitive tendering has already been met.

#### **NHS Supply Chain**

Medical Healthcare Furniture - Ref: 2020/S 252-638069

#### **Shared Business Services**

Design, Furniture & Appliances - Ref: SBS/19/MA/TYX/9363/23

#### **National Services Scotland**

Medical Ward Equipment - Ref: NP846/21

#### **HSC Business Services**

Medical Ward Equipment - Ref: NP846/21

#### Your Partner for Healthcare Equipment & Furniture

Through many years' experience we can support you with the facilitation and project management of tenders and contracts.

- · Specifying Requirements We will analyse your needs to identify products that meet your requirements
- Product Demonstrations We will allow you the opportunity to gain 'hands on' experience of our products to ensure that it is suitable for your needs
- Project Planning We will work with you and your team to identify any constraints and agree delivery phasing to ensure that we are able to fully meet your expectations
- Scheduled Deliveries We will deliver your order using our own transport to an agreed delivery schedule and, if agreed place the product and remove any packaging
- Installation & Commissioning The majority of the products that we supply requires no installation or commissioning, however we
  can provide this where it is required
- Training We can supply training to ensure that the products are used and maintained in the correct manner
- Full Service & Repair Support Bristol Maid branded products are designed and manufactured to meet the demands of everyday use, which under normal conditions should provide many years of trouble-free use. However, there are times where attention is required due to wear and tear or misuse
- Maintenance Contracts Peace of mind service contracts to ensure that your product is periodically checked, and where necessary adjustments and repairs are made
- · Site Audit & Equipment Recommendations We will support you with site and equipment audits

#### **Custom Design & Manufacture**

Whether it is a trolley that is slightly lower, a work surface that is slightly larger, a cabinet that must fit into a specific space or a unique colour scheme on a bedside cabinet, we recognise

that our standard product range may not always to be able to provide a product that meets your specific requirements.

If you require a unique product, please contact us so we can work together to develop a solution that will exactly match your needs. Our design capability allows us to provide everything from minor product specification changes through to a full design & development service.

# MEDICAL FURNITURE & EQUIPMENT 2022

Founded in 1953, Bristol Maid has been manufacturing and supplying medical furniture & equipment for over 65 years. From design and manufacture to supply and installation, Bristol Maid has a reputation for providing high product quality, excellent customer service and value for money. We take pride in working in close partnership with our customers to deliver effective solutions based on innovative design and technology that will meet the demands of today while anticipating the solutions for tomorrow.

Our highly-trained staff are experienced in every aspect of customer support, from simple product selection to bespoke technical solutions, from on-site audits to complex project management arrangements, according to our customers' needs. Bristol Maid's commitment to quality is based on years of experience working with healthcare professionals in an exacting and regulated environment, and is supported by an ISO 9001:2015 quality management system that is registered with a UKAS accredited certification body.

If you require further assistance, please do not hesitate to contact our Customer Service team who are here to help; or alternatively can arrange for your Regional Sales Representative to visit and help find a solution that best suits your requirements.





# **Bariatric Equipment**

| Folding Wheelchairs            |  |
|--------------------------------|--|
| Portering Chairs               |  |
| Mobile Shower / Commode Chairs |  |
| Shower / Commode Chairs        |  |
| Chairs                         |  |

# **Bedside Cabinets & Overbed Tables**

| Polymer Bedside Cabinets                |
|---|
| Bedside Cabinets                        |
| Bedside Cabinet / Wardrobe Combinations |
| Polymer Overbed & Overchair Tables      |
| Laminate Overbed & Overchair Tables     |

# **Caretray Trolleys & Racks**

| Features                                 |
|--|
| Disassembly & Cleaning                   |
| Stainless Steel, Low Level Trolleys      |
| Stainless Steel, Standard Level Trolleys |
| Stainless Steel, High Level Trolleys     |
| Stainless Steel, Mobile Racks            |
| Stainless Steel, Static Racks            |
| Mild Steel, Low Level Trolleys           |
| Mild Steel, Standard Level Trolleys      |
| Mild Steel, High Level Trolleys          |
| Mild Steel, Mobile Racks                 |
| Mild Steel, Static Racks                 |
| Accessories                              |

# **Chart Trolleys, Workstations Room Dividers & Screens**

| Fixed Height Trolleys        |
|------------------------------|
| Fixed Height Workstations    |
| Variable Height Workstations |
| Chart Boards                 |
| Document Holders             |
| Room Dividers                |
| Curtain Screens              |
| Solid Screens                |
|                              |

# **Concealment Trolleys**

| Variable Height Trolleys |  |
|--------------------------|--|
| Body Trays & Covers      |  |

#### **Cylinder Storage & Transportation**

| Cylinder Trolleys - F/G Size |  |
|------------------------------|--|
| Cylinder Trolleys - J Size   |  |

# **Cylinder Storage & Transportation**

| Bulk Transportation Trolleys     |
|----------------------------------|
| Cylinder Trolleys - E, CD/D Size |
| Wall Brackets & Cylinder Holders |

# **Dispensers, Holders & Racks**

| Apron Dispensers, Pocket Mask Holders            |
|--|
| Glove / Mask Box Holders                         |
| Alcohol Gel Holders                              |
| Sharps Box Holder, Urine Bottle Holder, Catheter |
| Bag Holders                                      |
| Catheter Racks                                   |
| Drainage Racks                                   |
|  |

# **Dressing / Instrument Trolleys**

| Removable Shelves - Stainless Steel         |
|---|
| Removable Shelves - Mild Steel              |
| Fixed Shelves, Sides Up - Stainless Steel   |
| Fixed Shelves, Sides Down - Stainless Steel |
| Drawer Units                                |
| Options & Accessories                       |
|   |

# **Drug & Medicine Dispensing**

| Dispensing Trolleys                       | 124 |
|---|-----|
| Laptop Holders                            | 129 |
| CLIQ <sup>®</sup> & RFID Solutions        | 131 |
| Patient Administration Trolleys           | 132 |
| Drawer Trolleys                           | 134 |
| Tray Trolleys                             | 135 |
| Blister Pack (MDS) Trolleys               | 136 |
| Nomad <sup>®</sup> Cassette Trolleys      | 137 |
| Biodose <sup>®</sup> & MultiMeds Trolleys | 138 |
| Medinoxx Trolleys                         | 139 |
| Security Clamp                            | 142 |
| Distribution Trolleys                     | 143 |
|   |     |

# **Drug & Medicine Storage**

| Drug & Medicine Cabinets - Shelves & Door Shelves      | 146 |
|--|-----|
| Drug & Medicine Cabinets - Deep Shelves                | 148 |
| Drug & Medicine Cabinets - Blister Pack (MDS)          | 150 |
| Controlled Drug Cabinets                               | 152 |
| Auditable Access Control for Medicines Management      | 154 |
| Drug & Medicine Cupboards                              | 156 |
| Options  | 158 |
| Patients / Residents Own Medication Cabinets - Self    | 160 |
| Administration   |     |
| Patients / Residents Own Medication Cabinets - Blister | 163 |
| Pack (MDS)   |     |
|  |     |

# **Drug & Medicine Storage**

| Tablet Counter & Storage Boxes   | 164 |
|----------------------------------|-----|
| Pharmacy & Vaccine Refrigerators | 165 |

# **Examination / Treatment Couches & Chairs**

| Examination & Treatment Couches - Variable Height  | 168 |
|--|-----|
| Examination & Treatment Couches - Fixed Height     | 172 |
| Couch Steps  | 173 |
| Leg Rests  | 174 |
| Examination & Treatment Couches - Variable Height, | 176 |
| Bariatric  |     |
| Examination & Treatment Chairs - Variable Height   | 180 |
| Examination & Treatment Chairs - Variable Height,  | 181 |
| Bariatric  |     |
| Phlebotomy Chairs - Variable Height                | 182 |
| Phlebotomy Chairs - Variable Height Bariatric      | 183 |
| Phlebotomy Chairs - Fixed Height                   | 184 |
| Options  | 185 |
|  |     |

# **General Purpose Trolleys**

| Stainless Steel - Heavy Duty     | 190 |
|----------------------------------|-----|
| Stainless Steel - Light Duty     | 191 |
| Polymer - Shallow & Deep Shelves | 192 |
| Mild Steel - Laminate Shelves    | 193 |

# **HTM71 Storage Solutions**

| ScanCell <sup>®</sup> Carts            | 196 |
|--|-----|
| Distribution Trolleys                  | 201 |
| Open Transfer / Transport Trolleys     | 203 |
| HTM71 Racking                          | 204 |
| Trays, Baskets, Dividers & Accessories | 206 |
|  |     |

# Infusion / Pump Stands & Accessories

| Mobile Pump Stands            | 212 |
|-------------------------------|-----|
| Mobile Infusion Stands        | 216 |
| Infusion Poles                | 220 |
| Clamping Infusion Stands      | 221 |
| Wall Mounted Infusion Holders | 221 |

#### **Linen Trolleys & Bags**

| Soiled Linen Trolleys | 224 |
|-----------------------|-----|
| Linen Trolleys        | 227 |
| Hamper Trolleys       | 228 |
| Soiled Linen Bags     | 229 |
| Clean Linen Trolleys  | 230 |

# **Maternity Equipment**

| Baby Cribs - Variable Height, Battery Operation   |  |
|---|--|
| Baby Cribs - Fixed Height                         |  |
| Baby Cribs - Variable Height, Gas Strut Operation |  |

# **Medical Records**

| Medical Records Trolleys - Locking     |
|--|
| Medical Records Workstations - Locking |
| Laptop Holders                         |
| Carefile Filing Pockets                |
| Medical Records Trolleys - Filing      |
| Medical Records Transfer Trolleys      |

# Phlebotomy

| Caretray Trolley - Phlebotomy                  |
|--|
| Phlebotomy Chairs - Variable Height            |
| Phlebotomy Chairs - Variable Height, Bariatric |
| Phlebotomy Chairs - Fixed Height               |

# **Portering Chairs**

| Portering Chair - Rear Steer            |
|---|
| Portering Chair - Front Steer           |
| Portering Chair - Rear Steer, Bariatric |
| Options & Accessories                   |

# **Procedure Trolleys & Carts**

| Caretray Trolley - Anaesthesia             |
|--|
| Caretray Trolley - Cast                    |
| Caretray Trolley - Critical Care           |
| Caretray Trolley - Infection               |
| Caretray Trolley - IV Preparation          |
| Caretray Trolley - Medical / Surgical      |
| Caretray Trolley - Plaster                 |
| Caretray Trolley - Vaccination             |
| Caretray Trolley - Basics                  |
| Mobile Workstation - Mild Steel, Low Level |
| Procedure Carts - Polymer, Metro Flexline™ |
| Procedure Carts - Polymer, Metro Starsys™  |
| Procedure Carts - Polymer, Metro Basics™   |

#### **Resuscitation Trolleys & Carts**

| Resuscitation Trolleys - Mild Steel                 |
|---|
| Caretray Trolley - Resuscitation                    |
| Resuscitation Carts - Polymer, Metro Lifeline™      |
| Resuscitation Carts - Polymer, Metro Flexline™      |
| Airways / Difficult Intubation Trolley - Mild Steel |
| Caretray Trolley - Airways / Difficult Intubation   |

#### **Sackholders & Waste Bins**

| Removable Body Bins                  |  |
|--------------------------------------|--|
| Free Standing Sackholders            |  |
| Fixed Body Sackholders - Plastic     |  |
| Removable Body Sackholders - Plastic |  |

# **Seating Solutions**

| Office, Task & Operators Chairs            |
|--|
| 24 Hour Task & Operator Chairs             |
| Bariatric Chairs                           |
| Saddle & Perch Stools                      |
| Low Back Patient Chairs                    |
| Medium & High Back Patient Armchairs       |
| Reclining Armchairs                        |
| Reclining Armchairs, Bariatric             |
| Reclining Porters Armchairs                |
| Rise & Recline Porters Armchairs           |
| Seat, Back Cushion & Headrest Options      |
| Overnight Chair                            |
| Soft Seating                               |
| Visitor & Waiting Room Seating             |
| Waiting Room Seating                       |
| Tech Chairs & Stools                       |
| PU Tech Chairs & Stools                    |
| Static Safe Tech Chairs & Stools           |
| Static Safe & Sterile Tech Chairs & Stools |
| Medical Grade Clean Room Tech Chairs       |
|  |

# Sepsis Trolleys

# **Shelving & Racking Solutions**

| 390 |
|-----|
| 392 |
| 394 |
| 396 |
| 398 |
| 399 |
|     |

# **Sterile Services Equipment**

| Preparation Tables - Variable Height          | 402 |
|---|-----|
| Preparation Tables - Fixed Height             | 403 |
| Sterile Supplies Trolleys - UN3291 Certified  | 404 |
| Sterile Supplies Trolleys - UN3291 Certified, | 406 |
| Trolley Washer Compatible                     |     |

# Storage Cupboards, Cabinets & Lockers

| Storage Cupboards  | 410 |
|--|-----|
| Storage Cabinets - Single & Double Door                  | 412 |
| Storage Cabinets - First Aid & Clinical Reagents Storage | 413 |
| Cabinets - Traka 21                                      | 414 |
| Storage Cabinets - Keys                                  | 415 |
| Storage Lockers - Personal Effects                       | 416 |

# **Theatre Equipment**

| Bowl Stands               | 420 |
|---------------------------|-----|
| Bucket Stands             | 423 |
| Instrument Trolleys       | 424 |
| Mayo Tables               | 424 |
|                           |     |
| Operating Platforms       | 426 |
| Suture Storage Trolleys   | 427 |
| Sterile Container Storage | 428 |
|                           |     |

# Wireless Mobile Carts & Workstations

| Mobile Workstations - Non-Powered | 432 |
|-----------------------------------|-----|
| Options & Accessories             | 434 |
| Caretray Trolley - Laptop         | 436 |

# Workstations & Speciality Storage

| Workstations           | 440 |
|------------------------|-----|
| Storage Units          | 441 |
| Catheter Storage Units | 442 |



# **Bariatric Equipment**

- Folding Wheelchairs 8
  - Portering Chairs 10
- Mobile Shower / Commode Chairs 11
  - Shower / Commode Chairs 12
    - Chairs 13



## **Bariatric Folding Wheelchair - Minimaxx & Disc**

#### **Features**

- Maximum capacity 325kg
- · Powder coated chromium molybdenum steel frame
- Rigid folding mechanism
- Low seat height to assist access
- Viscoelastic foam cushion with two way stretch vapour permeable incontinence cover and non-slip base
- · Padded push handles to assist the attendant/porter
- Puncture proof wheels & castors
- Parking brakes as standard/optional disc brake for improved control
- Side/width adjustable polyurethane moulded footrests, 100kg capacity

#### Dimensions (w x d x h)

- Seat Width 610/710mm
- Seat Height 430mm
- Seat Depth 460mm
- Backrest Height 510mm

#### **Options / Accessories**

- 5X-0100-099-005 Butterfly Armrests, Pair
- 5X-0100-099-006 Angled Armrests, Pair
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0100-099-001 Universal Table
- 5X-0100-099-020 IV & O2 Holder, 100mm Dia

#### **Spares**

- 5X-0230-061-045 Viscoelastic Cushion with Incontinence Cover, 610mm
- 5X-0230-071-045 Viscoelastic Cushion with Incontinence Cover, 710mm





IV & O<sup>2</sup> Holder

Universal Table

Elevating Calf Support

sales@bristolmaid.com

| ltem              | Description          | Overall Dimensions (w x d x h) | Weight |
|-------------------|----------------------|--------------------------------|--------|
| 5X-0100-061-000/K | Minimaxx, 610mm      | 850 x 1100 x 990mm             | 27.0kg |
| 5X-0100-071-000/K | Minimaxx, 710mm      | 950 x 1100 x 990mm             | 27.5kg |
| 5X-0109-061-000/K | Minimaxx Disc, 610mm | 850 x 1100 x 990mm             | 28.0kg |
| 5X-0109-071-000/K | Minimaxx Disc, 710mm | 950 x 1100 x 990mm             | 28.5kg |

Bariatric Equipment

# **Bariatric Folding Wheelchair - Minimaxx Push Motor**

#### Features

- Maximum capacity 325kg
- · Powder coated chromium molybdenum steel frame
- · Rigid folding mechanism
- · Low seat height to assist access
- Viscoelastic foam cushion with two way stretch vapour permeable incontinence cover and non-slip base
- · Padded push handles to assist the attendant/porter
- · Puncture proof wheels & castors
- Parking brakes as standard/optional disc brake for improved control
- Side/width adjustable polyurethane moulded footrests, 100kg capacity
  Push motor to assist the attendant/porter, reduces the risk of sore
- shoulders & stress on wrists
- Compatible on flat surfaces up to 5 degrees
- Operated using a control panel fitted to the handle, offering drive & reverse modes with speed control from 0.6-3.1mph
- · One charge provides a range of 2.75miles
- Supplied with one removable battery and charger (battery has to be removed from the unit when charging)

#### Dimensions (w x d x h)

- Seat Width 610/710mm
- · Seat Height 430mm
- · Seat Depth 460mm
- Backrest Height 510mm

#### **Options / Accessories**

- 5X-0100-099-005 Butterfly Armrests, Pair
- 5X-0100-099-006 Angled Armrests, Pair
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0100-099-001 Universal Table
- 5X-0100-099-020 IV & O2 Holder, 100mm Dia

#### **Spares**

- 5X-0230-061-045 Viscoelastic Cushion with Incontinence Cover, 610mm
- 5X-0230-071-045 Viscoelastic Cushion with Incontinence Cover, 710mm

| Item              | Description                 | Overall Dimensions (w x d x h) | Weight |
|-------------------|-----------------------------|--------------------------------|--------|
| 5X-0110-061-000/K | Minimaxx, Push Motor, 610mm | 850 x 1100 x 990mm             | 50.0kg |
| 5X-0110-071-000/K | Minimaxx, Push Motor, 710mm | 950 x 1100 x 990mm             | 50.5kg |







Push Motor

Handle Mounted Control Panel

# Portering Chair - Bariatric, Rear Steer, Sliding Foot Rest





**Features** 

- · Maximum capacity 380kg
- Seamless & contoured upholstery using flame retardant vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles mounted to rear uprights
- Retractable diagonal seatbelt
- 200mm fixed wheels/200mm total stop braking castors (cushioned)
- We recommend that two members of staff are in attendance when used with a patient in excess of 190kg

#### Dimensions (w x d x h)

- Overall Dimensions 930 x 1000 x 1090mm
- Seat Width 740mm
- Seat Height 530mm
- Seat Depth 670mm
- Backrest Height 480mm

#### **Standard Colours**

- · Paint Grey
- Upholstery Bristol Blue

#### **Options / Accessories**

- G/500/CYL Cylinder Holder, CD/D
- G/CLS/500 Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/CARD/STATION Head of Station
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- 2LOKE011 Release Card

Alternative colours are available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Rear Steer Portering Chairs - Bariatric, Sliding Foot Rest in Aqua Marine - G/500/RS/VAQ

Vinyl - Ambla New Richmond in Colour



| Item         | Description                                 | Weight |
|--------------|---|--------|
| G/500/RS/VBB | Rear Steer, Sliding Foot Rest, Bristol Blue | 43.8kg |

+44 (0)1258 484455

# **Bariatric Mobile Shower / Commode Chairs**

#### **Features**

- Suitable for showering or toileting (supplied with a bucket & lid)
- Powder coated chromium molybdenum steel frame
- Available in two widths 610/710mm
- Rotating & removable straight or butterfly armrests to assist access
- Open back to allow easy sling access & assist positioning
- Polyurethane moulded seat & back
- Seat supplied with a removable lid
- Seat panel removable for cleaning
- Long open hole in seat (220mm at the widest point)
- Deep bucket with lid is removable from the front or back of the chair
- Side/width adjustable polyurethane moulded footrests, 100kg capacity
- 125mm swivel castors, 2 off braking
- Maximum capacity 325kg

#### Dimensions (w x d x h)

- Seat Width 610/710mm
- · Seat Height 550mm
- · Seat Depth 560mm
- Backrest Height 480mm

#### **Standard Colours**

- Paint White
- Polyurethane Blue

#### **Options / Accessories**

- 5X-0150-999-022 Seat Panel 610mm
- 5X-0150-999-023 Seat Panel 660/710mm
- 5X-0140-999-001 Commode Bucket & Lid
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0156-099-010 Butterfly Armrests, Pair
- 5X-0150-099-010/K Straight Armrests, Pair

| Item               | Description   | Overall Dimensions (w x d x h) | Weight |
|--------------------|---|--------------------------------|--------|
| Standard Armrest   |   |                                |        |
| 5X-0150-061-000    | Bariatric Mobile Shower / Commode Chair - Standard Armrests, 610mm  | 670 x 1120 x 990mm             | 27.0kg |
| 5X-0150-071-000    | Bariatric Mobile Shower / Commode Chair - Standard Armrests, 710mm  | 770 x 1120 x 990mm             | 27.5kg |
| Butterfly Armrests |   |                                |        |
| 5X-0151-061-000    | Bariatric Mobile Shower / Commode Chair - Butterfly Armrests, 610mm | 840 x 1120 x 990mm             | 27.5kg |
| 5X-0151-071-000    | Bariatric Mobile Shower / Commode Chair - Butterfly Armrests, 710mm | 940 x 1120 x 990mm             | 28.0kg |



+44 (0)1258 484455

CE



#### 5X-0142-061-000

# **Bariatric Shower / Commode Chair**

#### **Features**

- Suitable for showering or toileting (supplied with a bucket & lid) Powder coated chromium molybdenum steel frame
- Available in two widths 610/710mm
- Seat height adjustable 420-550mm
- · Long armrests to assist access
- Open back to allow easy sling access & assist positioning
- Polyurethane moulded seat & back
- Seat supplied with a removable lid
- Seat panel removable for cleaning
- Long open hole in seat (220mm at the widest point)
- Optional deep bucket is removable from the front or back of the chair
- Maximum capacity 325kg

#### Dimensions (w x d x h) Standard Colours

- Seat Width 610/710mmSeat Height 420-550mm
- Paint WhitePolyurethane Blue
- Seat Depth 400-520mm

#### **Options / Accessories**

- 5X-0150-999-022 Seat Panel 610mm
- 5X-0150-999-023 Seat Panel 660/710mm
- 5X-0140-999-001 Commode Bucket & Lid



| Item            | Description                              | Overall Dimensions (w x d x h) | Weight |
|-----------------|--|--------------------------------|--------|
| 5X-0142-061-000 | Bariatric Shower / Commode Chair - 610mm | 730 x 1010 x 790 - 940mm       | 15.0kg |
| 5X-0142-071-000 | Bariatric Shower / Commode Chair - 710mm | 830 x 1010 x 790 - 940mm       | 16.0kg |

# **Bariatric High Back Chairs**

#### **Features**

- Deep CMHR foam
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- · Fabric Fire retardant
- Black heavy duty steel five star base
- Black heavy duty gas lift
- · Black height adjustable soft touch armrests
- Easy lift ratchet height adjustable back
- Independent seat & back mechanism with side tension control
- Maximum capacity 315kg

Replace suffix on the product code with any of the following to select your colour i.e. Bariatric High Back in Lime - 5660XL/VLI

#### Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage









# **Bariatric Chairs**

#### Features

- Deep CMHR foam
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- · Beech, varnished four leg frame
- Maximum capacity 315kg

Alternative colours are available at extra cost. Replace suffix on the product code with any of the following to select your colour i.e. Bariatric, Beech Frame in Lime - 5671/VLI

Vinyl - Ambla New Richmond in Colour



Fabric - Camira Advantage



| ltem | Description            | Overall Dimensions (w x d x h) | Seat Height | Weight |
|------|------------------------|--------------------------------|-------------|--------|
| 5671 | Bariatric, Beech Frame | 844 x 580 x 890mm              | 490mm       | 16.0kg |

# Other Products you may be interested in...

See Pages 176 & 178

Examination & Treatment Couch Two Section, Bariatric Examination & Treatment Couch Three Section, Bariatric







# **Bedside Cabinets** & Overbed Tables

- Polymer Bedside Cabinets 18
  - Bedside Cabinets 20
- Bedside Cabinet / Wardrobe Combinations 26
  - Polymer Overbed & Overchair Tables 30
  - Laminate Overbed & Overchair Tables 32

# eClean™



# **Polymer Bedside Cabinet - Easy Clean**

#### **Features**

- BS 3962 Part 6:1980 Assessment of resistance to mechanical damage
- BS EN 12720:2009 Furniture assessment of surface resistance to cold liquids
- BS EN 12721:2009 Furniture assessment of surface resistance to wet heat
- BS EN 12722:2009 Furniture assessment of surface resistance to dry heat
- BS 5852:2006 Methods of the test for assessment of the ignitability
- BS 4875-7:2006 Domestic & contract storage furniture
- Developed as part of Design Bugs Out & evaluated by the Department of Health
- Manufactured from recyclable polyethylene
- Open backless design ensures that all areas are accessible, whilst surfaces allow for easy & effective cleaning
- Distinct upper, middle & lower zones
  - Upper Nurse zone with integral rim to contain any spillages. Ideally positioned for staff to use while standing
  - Middle Patient zone provides storage area for patient use
  - Lower Patient zone provides shelf for washbowl, footwear etc.
- Drawers can be inserted from either side, allowing the cabinet to be positioned either side of the bed
- Upper & middle drawers can be fitted with an L15 125KHz RFID lock, providing secure storage for patient medicines / property\*
- 100mm polymer twin-wheel braking castors

\*cabinets with locks are supplied with one patient card - a programming card will be required per ward / directorate as per hospital policy

| Item   | Description                        | Weight |
|--------|------------------------------------|--------|
| EBC015 | No Locks                           | 29.0kg |
| EBC020 | Upper Drawer with Lock*            | 29.1kg |
| EBC025 | Middle Drawer with Lock*           | 29.1kg |
| EBC030 | Upper & Middle Drawers with Locks* | 29.2kg |

#### Dimensions (w x d x h)

- Overall Dimensions 480 x 665 x 1237mm
- Drawers (internal)
  - Upper 149 x 222 x 217mm
    Middle 338 x 491 x 197mm

  - Lower 338 x 491 x 220mm

### **Standard Colours**

- Body Grey White
- Drawers Spearmint

#### **Options / Accessories**

• RFID Lock Cards - page 25





BE

CC

ΒA

# **Bedside Cabinet - Single Upper Drawer**

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper & personal effects drawer
- Lower storage offered through...
  - Cupboard with double doors (adjustable shelf model dependant)
  - Drawer (recessed shelf model dependant)
- Available as either two-tone / single colour
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions: 485 x 485 x 905mm
- Internal Dimensions...
  - Upper Drawer 413 x 381 x 182mm
  - Personal Drawer 413 x 381 x 102mm
  - Lower Drawer 413 x 381 x 505 / 297mm
  - Cupboard 445 x 443 x 537 / 394mm

#### **Standard Colours**

- Two-Tone
  - · Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour Beech

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White | Beech     | Description  | Weight |
|------------|-----------|--|--------|
| BC1C/GW    | BC1C/BE   | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 34.0kg |
| BC1CPD/GW  | BC1CPD/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors  | 35.8kg |
| BC1D/GW    | BC1D/BE   | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 40.4kg |
| BC1DPD/GW  | BC1DPD/BE | Single Upper Drawer, Personal Drawer & Lower Drawer                      | 42.1kg |

NO

NC

sales@bristolmaid.com +

# **Bedside Cabinet - Double Upper Drawer**

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawers suitable for the storage of patient's medication or personal effects (1 x fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for either drawer
- Lower storage offered through...
  - Cupboard with double doors (adjustable shelf)
  - Drawer (recessed shelf)
- Available as either two-tone / single colour
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions 485 x 485 x 905mm
- Internal Dimensions...
  - Upper Drawer 180 x 381 x 182mm
  - Lower Drawer 413 x 381 x 505mm
  - Cupboard 445 x 443 x 537mm

#### **Standard Colours**

- Two-Tone
  - Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- · Single Colour Beech

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23



| <b>Grey White</b> | Beech   | Description  | Weight |
|-------------------|---------|--|--------|
| BC2C/GW           | BC2C/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 34.2kg |
| BC2D/GW           | BC2D/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf                    | 40.5kg |



Alternative Colours Available



+44 (0)1258 484455







BC3C/GW

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper section suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- · Rear side doors allow access to contents
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper section & personal effects drawer
- Lower storage offered through...
  - Cupboard with double doors (adjustable shelf model dependant)
  - Drawer (recessed shelf model dependant)
- Available as either two-tone / single colour
- 50mm twin-wheel polypropylene castors

#### **Dimensions** (w x d x h)

- Overall Dimensions 485 x 485 x 905mm
- Internal Dimensions...
  - Upper Section 445 x 135 x 249mm
  - Rear Side Section 290 x 445 x 249mm
  - Personal Drawer 413 x 381 x 102mm
  - Lower Drawer 413 x 381 x 505 / 297mm
  - Cupboard 445 x 443 x 537 / 394mm

#### **Standard Colours**

Two-Tone

www.bristolmaid.com

- · Body Grey White
- Top / Drawer / Door Fronts Spearmint
- · Single Colour Beech

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White | Beech     | Description  | Weight |
|------------|-----------|--|--------|
| BC3C/GW    | BC3C/BE   | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 31.8kg |
| BC3CPD/GW  | BC3CPD/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors  | 33.5kg |
| BC3D/GW    | BC3D/BE   | Upper Section, Lower Drawer with Recessed Shelf                    | 38.1kg |
| BC3DPD/GW  | BC3DPD/BE | Upper Section, Personal Drawer & Lower Drawer                      | 39.8kg |



**Overbed Tables** 

**Bedside Cabinets &** 

# **Towel Holder**

#### Features

- · Self-adhesive hook
- Supplied for customer to fit



| Item | Description  | Weight |
|------|--------------|--------|
| TH   | Towel Holder | 0.1kg  |

# **Removable Plastic Top**

#### **Features**

- Raised edges on three sides reduces misplacement of objects
- Clips to the top of the bedside cabinet, & is removable to assist cleaning
- Not compatible with bedside cabinets fitted with wardrobes



RPT shown with BC1C/GW

| Item | Description           | Weight |
|------|-----------------------|--------|
| RPT  | Removable Plastic Top | 0.8kg  |

# **Electronic Push Button & RFID Locks Solutions**

Through either electronic push button or contactless RFID, we can offer a secure locking solution to suit your specific requirements.

#### **RFID Locks**

- Available as 125KHz or 13.56Mhz frequencies
- Easy programming using a dedicated card
- Three / four user levels offer differing levels of authority Administration / Master / Sub-Master / User
- Available with clamshell cards or fobs
- Compatible with Abloy CLIQ RFID chip inserted into the key







EBC Lock (L15)

PS Solo RFID Lock (PL)

ABLOY RFID Lock (AL)

|                        | EBC Lock (L15)          | PS Solo Lock (PL)       | ABLOY RFID Lock (AL)   |
|------------------------|-------------------------|-------------------------|--|
| Frequency              | 125KHz                  | 125KHz                  | 13,56Mhz   |
| Programming            | Programming Card        | Programming Card        | Programming Card   |
| Modes                  | Public & Private        | Public & Private        | Public & Private   |
| Audit Trail            | No                      | No                      | No   |
| User Levels            | 3 - Master, Staff, User | 3 - Master, Staff, User | 4 - Administration, Master, Staff, User                        |
| No. of Cards / Fobs    | 50 x Cards              | 50 x Cards              | 50 x Cards / Fobs  |
| Batteries              | 1 x CR123               | 1 x CR123               | 3 x AAA  |
| Openings (Battery Set) | 20,000                  | 20,000                  | 4,000  |
| Low Battery Warning    | Audible Alarm           | Audible Alarm           | Audible Alarm  |
| Low Battery Override   | Fails Unlocked          | Fails Unlocked          | Micro USB - External Power Source (power bank / phone charger) |
| Guarantee              | 1 Year                  | 2 Years                 | 1 Year   |
| Abloy CLIQ Compatible* | Yes                     | Yes                     | Yes  |

\*RFID chip inserted into the CLIQ key

#### **RFID Cards & Fobs**

- Compatible with 125KHz & 13,56Mhz frequencies
- Available as clamshell cards or fobs
- Can be fitted to a key ring
- A programming card is required for each hospital or ward as per hospital procedure

#### Dimensions (w x d x h)

- Clamshell Card 86 x 54 x 2mm
- Key Fob 27 x 35 x 5mm



#### Frequency - 13.56 mhz

| Item        | Description                 | Lock Compatibility |
|-------------|-----------------------------|--------------------|
| RFID1356/PC | 13,56KHz Programming Card   | AL                 |
| RFID1356/SC | 13,56KHz Staff Card (Green) | AL                 |
| RFID1356/UC | 13,56KHz User Card (Blue)   | AL                 |
| RFID1356/SF | 13,56KHz Staff Fob (Green)  | AL                 |
| RFID1356/UF | 13,56KHz User Fob (Blue)    | AL                 |

#### Frequency - 125KHz

| Item       | Description                                     | Lock Compatibility |
|------------|---|--------------------|
| RFID125/PC | 125KHz Programming Card - Staff<br>Pair (Black) | PL, L15            |
| RFID125/SC | 125KHz Staff Card (Green)                       | PL, L15            |
| RFID125/MC | 125KHz Master Card (Red)                        | PL, L15            |
| RFID125/UC | 125KHz User Card (Blue)                         | PL, L15            |

#### **Electronic Push Button Lock**

- Uses a 10-button keypad & # button
- Easy programming & code change
- Four user levels offer differing levels of authority...
  - Master / Technician 8 digits
    Staff / User 4 digits
- Available with a key override, allowing the lock to be opened immediately in the event of a low battery





CodeLock Key Overide (CLOR)

Electronic Push Button Lock (CL)

|                        | CodeLock Key Overide (CLOR)         | Electronic Push Button Lock (CL)    |
|------------------------|-------------------------------------|-------------------------------------|
| Programming            | 10-Button Keypad & # Button         | 10-Button Keypad & # Button         |
| Modes                  | Public & Private                    | Public & Private                    |
| Audit Trail            | No                                  | No                                  |
| User Levels            | 4 - Master, Technician, Staff, User | 4 - Master, Technician, Staff, User |
| User / Staff Codes     | Twenty                              | Single                              |
| Batteries              | 2 x AAA                             | 2 x AAA                             |
| Openings (Battery Set) | 100,000                             | 15,000                              |
| Low Battery Warning    | Red Flashing LED                    | Red Flashing LED                    |
| User / Staff Codes     | 9V PP3 Battery + Code               | 9V PP3 Battery + Code               |
| Guarantee              | 2 Years                             | 2 Years                             |
| Abloy CLIQ Compatible  | N/A                                 | N/A                                 |

## **Compact Bedside Cabinet / Wardrobe Combination - Single Upper Drawer**





Alternative Colours Available





BCC1CWL/GW

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper drawer
- Lower storage offered through...
  - · Cupboard with double doors (adjustable shelf)
  - Drawer (recessed shelf)
- · Available as either two-tone / single colour
- · Wardrobe supplied with a hanging rail, available right / left hand opening
- 50mm twin-wheel polypropylene castors

#### **Dimensions** (w x d x h)

- Overall Dimensions 633 x 489 x 1420mm
- Internal Dimensions...
  - Upper Drawer 328 x 381 x 182mm
  - Lower Drawer 328 x 381 x 505mm
  - Cupboard 360 x 443 x 537mm
  - Wardrobe 215 x 443 x 1420mm

#### **Standard Colours**

- Two Tone
  - Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour Beech

#### **Options / Accessories**

www.bristolmaid.com

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White | Beech      | Description  | Weight |
|------------|------------|--|--------|
| Right Hand |            |  |        |
| BCC1CWR/GW | BCC1CWR/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 50.6kg |
| BCC1DWR/GW | BCC1DWR/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 56.9kg |
| Left Hand  |            |  |        |
| BCC1CWL/GW | BCC1CWL/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 50.6kg |
| BCC1DWL/GW | BCC1DWL/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 56.9kg |

sales@bristolmaid.com

# **Bedside Cabinet / Wardrobe Combination - Single Upper Drawer**

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper drawer & personal effects drawer
- Lower storage offered through...
  - Cupboard with double doors (adjustable shelf)
- Drawer (recessed shelf model dependant)
- Available as either two-tone / single colour
- · Wardrobe supplied with a hanging rail, available right / left hand opening
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions: 718 x 489 x 1420mm
- Internal Dimensions...
  - Upper Drawer: 413 x 381 x 182mm
  - Personal Drawer: 413 x 381 x 102mm
  - Lower Drawer: 413 x 381 x 505 / 297mm
  - Cupboard: 445 x 443 x 537 / 394mm
  - Wardrobe: 215 x 443 x 1420mm

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White  | Beech       | Description  | Weight |
|-------------|-------------|--|--------|
| Right Hand  |             |  |        |
| BC1CWR/GW   | BC1CWR/BE   | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 54.5kg |
| BC1CPDWR/GW | BC1CPDWR/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors  | 56.2kg |
| BC1DWR/GW   | BC1DWR/BE   | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 60.8kg |
| BC1DPDWR/GW | BC1DPDWR/BE | Single Upper Drawer, Personal Drawer & Lower Drawer                      | 62.5kg |
| Left Hand   |             |  |        |
| BC1CWL/GW   | BC1CWL/BE   | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 54.5kg |
| BC1CPDWL/GW | BC1CPDWL/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors  | 56.2kg |
| BC1DWL/GW   | BC1DWL/BE   | Single Upper Drawer, Lower Drawer with Recessed Shelf                    | 60.8kg |
| BC1DPDWL/GW | BC1DPDWL/BE | Single Upper Drawer, Personal Drawer & Lower Drawer                      | 62.5kg |

#### **Standard Colours**

- Two Tone
  - Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour- Beech



Alternative Colours Available CP GE RB LA VI BI GW GP BF CA NO NC BA CC BE





BC2CWL/GW

Alternative Colours Available



## **Bedside Cabinet / Wardrobe Combination - Double Upper Drawer**

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper drawers suitable for the storage of patient's medication or personal effects (1 x fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for either drawer
- Lower storage offered through...
  - · Cupboard with double doors (adjustable shelf)
  - Drawer (recessed shelf)
- Available as either two-tone / single colour
- Wardrobe supplied with a hanging rail, available right / left hand opening
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions: 718 x 489 x 1420mm
- Internal Dimensions...
  - Upper Drawer: 180 x 381 x 182mm
  - Lower Drawer: 413 x 381 x 505mm
  - Cupboard: 445 x 443 x 537mm
  - Wardrobe: 215 x 443 x 1420mm

#### **Standard Colours**

- Two Tone
  - Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour Beech

#### **Options / Accesssories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White | Beech     | Description  | Weight |
|------------|-----------|--|--------|
| Right Hand |           |  |        |
| BC2CWR/GW  | BC2CWR/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 55.0kg |
| BC2DWR/GW  | BC2DWR/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf                    | 61.3kg |
| Left Hand  |           |  |        |
| BC2CWL/GW  | BC2CWL/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 55.0kg |
| BC2DWL/GW  | BC2DWL/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf                    | 61.3kg |

sales@bristolmaid.com +44

# Bedside Cabinet / Wardrobe Combination - Upper Section & Rear Side Doors

#### **Features**

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients
- Melamine faced chipboard / laminate construction
- Upper section suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- · Rear side doors allow access to contents
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button or 125KHz / 13.56Mhz RFID locks available for the upper section & personal effects drawer
- Lower storage offered through...
  - · Cupboard with double doors (adjustable shelf model dependant)
  - Drawer (recessed shelf model dependant)
- · Available as either two-tone / single colour
- Wardrobe supplied with a hanging rail, available right / left hand opening
- 50mm twin-wheel polypropylene castors

#### Dimensions (w x d x h)

- Overall Dimensions: 718 x 489 x 1420mm
- Internal Dimensions...
  - Upper Section: 445 x 135 x 249mm
  - Rear Side Section: 290 x 445 x 249mm
  - Personal Drawer: 413 x 381 x 102mm
  - Lower Drawer: 413 x 381 x 505 / 297mm
  - Cupboard: 445 x 443 x 537 / 394mm
  - Wardrobe: 215 x 443 x 1420mm

#### **Standard Colours**

- Two Tone
  - · Body Grey White
  - Top / Drawer / Door Fronts Spearmint
- Single Colour Beech

#### **Options / Accessories**

- Electronic Push Button & RFID Locks page 24/25
- RFID Lock Cards & Fobs page 25
- Towel Holder & Removable Plastic Top page 23

| Grey White  | Beech       | Description  | Weight |
|-------------|-------------|--|--------|
| Right Hand  |             |  |        |
| BC3CWR/GW   | BC3CWR/BE   | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 52.1kg |
| BC3CPDWR/GW | BC3CPDWR/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors  | 53.8kg |
| BC3DWR/GW   | BC3DWR/BE   | Upper Section, Lower Drawer with Recessed Shelf                    | 58.4kg |
| BC3DPDWR/GW | BC3DPDWR/BE | Upper Section, Personal Drawer & Lower Drawer                      | 60.1kg |
| Left Hand   |             |  |        |
| BC3CWL/GW   | BC3CWL/BE   | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 52.1kg |
| BC3CPDWL/GW | BC3CPDWL/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors  | 53.8kg |
| BC3DWL/GW   | BC3DWL/BE   | Upper Section, Lower Drawer with Recessed Shelf                    | 58.4kg |
| BC3DPDWL/GW | BC3DPDWL/BE | Upper Section, Personal Drawer & Lower Drawer                      | 60.1kg |



BC3CPDWR/GW

#### Alternative Colours Available





# **Overbed & Overchair Table - Polymer, Easy Clean**

#### **Features**

- BS 2483: 1977 Specification for overbed tables (selected tests)
- Finish Performance Tests (FIRA Standard 6250 2005)
  - BS 3962-6: 1980 Methods of test for finishes for wooden furniture
  - BS EN 12720: 2009 + A1: 2013 Assessment of surface resistance to cold liquids
  - BS EN 12721: 2009 + A1: 2013 Assessment of surface resistance to wet heat
  - BS EN 12722: 2009 + A1: 2013 Assessment of surface resistance to dry heat
- Height adjustable using a gas damper (760 1070mm)
- Safety feature prevents trapping of the patient between the bed / chair & table
- · Extruded aluminium column provides a smooth linear height adjustment
- Cantilever tubular base available as...
  - C shaped suitable for use with beds & chairs
  - · Straight leg suitable for use with beds
- Recyclable polypropylene one-piece top, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
- · Raised edges reduces misplacement of objects & contains spillages
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors
- Maximum load 25kg

#### Dimensions (w x d x h)

• Table Top - 880 x 440mm

#### **Standard Colours**

Top - Grey White

- Upper Extrusion Spearmint
- Lower Extrusion / Base Grey White

#### **Options / Accessories**

• B - Braking Castors (x 2)

| Item     | Description  | Overall Dimensions (w x d x h) | Weight |
|----------|--|--------------------------------|--------|
| OBP/SP/S | Overbed Table - Height Adjustable, Polymer, Easy Clean   | 895 x 440 x 760 - 1070mm       | 14.0kg |
| OCP/SP/S | Overchair Table - Height Adjustable, Polymer, Easy Clean | 918 x 440 x 760 - 1070mm       | 13.5kg |

sales@bristolmaid.com

+44 (0)1258 484455



# **Overbed & Overchair Table - Laminate**

#### **Features**

• BS 2483: 1977 Specification for overbed tables (selected tests)

 BS EN 438-2: 2016 High-pressure decorative laminates (HPL). Sheets based on thermosetting resins (usually called laminates). Determination of properties

- Height adjustable using a gas damper (735 1045mm)
- Safety feature prevents trapping of the patient between the bed / chair & table
- Extruded aluminium column provides a smooth linear height adjustment
- Cantilever tubular base available as...
  - · C shaped suitable for use with beds & chairs
  - · Straight leg suitable for use with beds
- · Laminate top with varnished edges available in three configurations...
  - Flat top
  - One varnished beech lip
  - Three varnished beech lips
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors
- · Maximum load 25kg

#### Dimensions (w x d x h)

- Table Top 880 x 395mm
- Overbed Table 895 x 440 x 735 -1045mm
- Overchair Table 918 x 440 x 735 -1045mm

#### **Standard Colours**

- Top Grey White or Beech
- Upper Extrusion Spearmint or Grey White
- · Lower Extrusion / Base Grey White

#### **Options / Accessories**

• B - Braking Castors (x 2)

#### **Overbed Tables**

| Grey White   | Beech        | Description | Weight |
|--------------|--------------|-------------|--------|
| OBL/GW0/SP/S | OBL/BE0/GW/S | Flat Top    | 14.5kg |
| OBL/GW1/SP/S | OBL/BE1/GW/S | One Lip     | 15.4kg |
| OBL/GW3/SP/S | OBL/BE3/GW/S | Three Lips  | 15.8kg |



sales@bristolmaid.com

OBL/GW1/SP/S

+44 (0)1258 484455

HEIGHT



www.bristolmaid.com



# **Caretray Trolleys & Racks**

- Mild Steel, Low Level Trolleys 50
- Mild Steel, Standard Level Trolleys 52
  - Mild Steel, High Level Trolleys 54
    - Mild Steel, Mobile Racks 56
    - Mild Steel, Static Racks 58
      - Accessories 60

- Features 36
- Disassembly & Cleaning 38
- Stainless Steel, Low Level Trolleys 40
- Stainless Steel, Standard Level Trolleys 42
  - Stainless Steel, High Level Trolleys 44
    - Stainless Steel, Mobile Racks 46
    - Stainless Steel, Static Racks 48

# **Features**

# Carerails

 Fitted to both sides of the trolley, compatible with a wide range of accessories

# Frame

- · Fully welded eliminating gaps
- 304 grade stainless steel or mild steel that is powder coated in a durable epoxy polyester paint

# **Side Panels**

- One-piece moulding
- Smooth surfaces & rounded edges assist cleaning (compatible with steam cleaning)
- Runners with positive stops prevent trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)

# **Easy Clean**

- One-piece mouldings
- Smooth surfaces & rounded edges assist cleaning (compatible with steam cleaning)
- Quick & easy tool-less disassembly assists thorough or deep cleaning programmes



# Тор

- One-piece moulding
- Smooth surfaces & rounded edges assist cleaning (compatible with steam cleaning)
- Raised edges on three sides provides defined grab points to act as push handles, whilst an integral rim at the front contains spillages
- Maximum load 30kg

# **Trays & Dividers**

- Translucent & impact-resistant (100 / 150mm depths)
- Smooth surfaces & rounded edges assist cleaning (compatible with steam cleaning & able to withstand high temperature machine washing)
- Trays feature an integral handle with plastic clip on label cover, (suitable for use with our perforated coloured label strips)
- Trays can be subdivided in both length & width using simple slot in dividers
- Maximum load 5kg

# **Castors & Ring Buffers**

- 100mm pressed steel or polymer castors (model dependant)
- Two braking castors fitted as standard
- Ring buffers offer protection to both the trolley & fabric of the building
- Smooth surfaces assist cleaning

Racks

øð

# Accessories

The Caretray Trolley has been developed specifically to meet the needs of Healthcare Professionals across a wide variety of disciplines including Minor Treatment, A&E. Theatre. Phlebotomy, ITU & Outpatients. Through a wide range of accessories, you can configure your trolley to suit your specific needs.

# Locking & Security - Page 60

- · Solid grade laminate / clear acrylic doors
- · Range of locking options
- Security clamp & seals to safeguard contents

# Laptop Holders - Page 74

- · Suitable for use with Electronic Prescribing & Medicines Administration (EPMA) & Electronic Patient Records (EPR)
- Offers height, rotation & angle adjustment
- Laptop can be secured using a security clamp



## **Overbridge & Tilt Bins - Page 66**

- Provides high level visible storage & assists organisation
- Tilt-bins available with 4, 5 or 6 clear compartments
- Compartments tilt to 45 degrees, to assist access of contents



### **Defibrillator Platform - Page 67**

- · Compatible with all popular defibrillators
- · Adjustable, allowing compatibility with all popular models
- Platform rotates 360 degrees, whilst the arm swings 90 degrees from the worksurface

# Castors & Ring Buffers - Page 76

- Range of polymer castors
- Available in various braking & non-braking configurations
- Ring buffers offer protection to both the trolley & fabric of the building



# **Carerail Mounted - Page 68**

- · An extensive range of accessories can be attached/suspended from the Carerails
- Allows the trolley to be configured to meet specific user requirements
- Accessories range from sharps & glove box holders, catheter racks & waste bins, through to work flaps & suction pump shelves



# Trays, Shelves, Dividers & Labels - Page 77

- · Additional dividers assist segregation
- · Coloured labels ensure that trays can be easily identified
- · Flat shelves provide an alternative to trays



sales@bristolmaid.com

# **Disassembly & Cleaning**

The design of the trolley allows for quick & easy tool-less disassembly providing access for all areas when conducting thorough or deep cleaning programmes. The simple guide below illustrates the disassembly & re-assembly process...






Remove the trays from the unit



# **Remove Top**

The top is held in place with two thumb screws, located on the underside of the top. Once removed, the top can be removed with a firm tug



## **Remove Back Panel** (if fitted)

Once the top has been removed, the back panel can be lifted clear of the cross members on the rear of the frame



# **Remove Side Panels**

The side panels are removed by lifting off the positioning inserts. Double height trolleys require thumb screws at the top of the lower panel to be removed.



# Clean

The frame, side panels & top can now be cleaned



# Reassembly

The trolley is re-assembled by repeating the steps in reverse

# Caretray Trolleys - Stainless Steel, Low Level, Single & Double Column





CTS1/08/3S1D - shown with optional CT/WF, CT/ALC & GMBH/15



Easy disassembly & cleaning - page 38

#### **Features**

- · Low level, single & double column easy clean procedure trolleys
- Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- · 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 807mm
  - Double Column 1005 x 550 x 807mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 782mm
- Double Column 915 x 455 (internal) at height 782mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

 Item
 Description
 Weight

 CTS1/08/1S2D
 1 x 100mm, 2 x 150mm Trays
 21.2kg

 CTS1/08/3D
 3 x 150mm Trays
 21.5kg

 CTS1/08/3D1
 3 x 100mm, 1 x 150mm Trays
 22.3kg

 CTS1/08/4S
 4 x 100mm Trays
 21.9kg

# **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67

- Carerail Mounted page 68
  Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77



CTS2/08/4S3D

| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTS2/08/3S4D | 3 x 100mm, 4 x 150mm Trays | 37.3kg |
| CTS2/08/4S3D | 4 x 100mm, 3 x 150mm Trays | 36.9kg |
| CTS2/08/6D   | 6 x 150mm Trays            | 36.5kg |
| CTS2/08/6S2D | 6 x 100mm, 2 x 150mm Trays | 38.1kg |
| CTS2/08/7S1D | 7 x 100mm, 1 x 150mm Trays | 37.7kg |
| CTS2/08/8S   | 8 x 100mm Trays            | 37.3kg |

CTS2/08/4S3D - shown with optional CT/IV, GMBH/15 & CT/ALC

### Caretray Trolleys - Stainless Steel, Standard Level, Single & Double Column



| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTS1/10/1S3D | 1 x 100mm, 3 x 150mm Trays | 27.9kg |
| CTS1/10/3S2D | 3 x 100mm, 2 x 150mm Trays | 29.0kg |
| CTS1/10/4D   | 4 x 150mm Trays            | 28.3kg |
| CTS1/10/4S1D | 4 x 100mm, 1 x 150mm Trays | 28.7kg |
| CTS1/10/6S   | 6 x 100mm                  | 29.8kg |



CTS1/10/3S2D - shown with optional CT/IV, CT/CR2, GMBH/15 & CT/ALC



Easy disassembly & cleaning - page 38

#### **Features**

- Standard level, single & double column easy clean procedure trolleys
- Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 993mm
  - Double Column 1005 x 550 x 993mm
- Plastic Top
  - Single Column 455 x 455 (internal) at height 968mm
  - Double Column 915 x 455 (internal) at height 968mm

#### Trays

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

aid.com +44 (0)1258 484455

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
- Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77



CTS2/10/6S4D

| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTS2/10/3S6D | 3 x 100mm, 6 x 150mm Trays | 46.6kg |
| CTS2/10/4S5D | 4 x 100mm, 5 x 150mm Trays | 46.2kg |
| CTS2/10/5S4D | 5 x 100mm, 4 x 150mm Trays | 45.9kg |
| CTS2/10/6S4D | 6 x 100mm, 4 x 150mm Trays | 47.3kg |
| CTS2/10/8D   | 8 x 150mm Trays            | 45.8kg |
| CTS2/10/8S2D | 8 x 100mm, 2 x 150mm Trays | 46.6kg |



CTS2/10/6S4D - shown with optional CT/IV, CT/CR2, GMBH/15 & CT/ALC



#### CTS1/14/4S3D

| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTS1/14/2S5D | 2 x 100mm, 5 x 150mm Trays | 36.2kg |
| CTS1/14/3S4D | 3 x 100mm, 4 x 150mm Trays | 35.8kg |
| CTS1/14/4S3D | 4 x 100mm, 3 x 150mm Trays | 35.4kg |
| CTS1/14/5S3D | 5 x 100mm, 3 x 150mm Trays | 37.8kg |
| CTS1/14/6D   | 6 x 150mm Trays            | 35.0kg |
| CTS1/14/8S1D | 8 x 100mm, 1 x 150mm Trays | 37.8kg |



CTS1/14/4S3D - shown with optional CT/DH/A4



Easy disassembly & cleaning - page 38

#### Features

Caretray Trolleys - Stainless Steel, High Level, Single & Double Column

- High level, single & double column easy clean procedure trolleys
- Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 1410mm
  - Double Column 1005 x 550 x 1410mm
- Plastic Top
  - Single Column 455 x 455 (internal) at height 1385mm
- Double Column 915 x 455 (internal) at height 1385mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

- Locking & Security page 60
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



CTS2/14/10S6D - shown with optional CT/DH/A4

| Item          | Description                 | Weight |
|---------------|-----------------------------|--------|
| CTS2/14/4S10D | 2 x 100mm, 10 x 150mm Trays | 60.6kg |
| CTS2/14/6S8D  | 6 x 100mm, 8 x 150mm Trays  | 62.8kg |
| CTS2/14/8S6D  | 8 x 100mm, 6 x 150mm Trays  | 63.0kg |
| CTS2/14/10S6D | 10 x 100mm, 6 x 150mm Trays | 65.0kg |
| CTS2/14/12D   | 12 x 150mm Trays            | 61.3kg |
| CTS2/14/16S2D | 16 x 100mm, 2 x 150mm Trays | 66.6kg |



CTS2/14/10S6D





#### CRSM1/18/6S4D

| Item          | Description                | Weight |
|---------------|----------------------------|--------|
| CRSM1/18/2S6D | 2 x 100mm, 6 x 150mm Trays | 41.7kg |
| CRSM1/18/3S6D | 3 x 100mm, 6 x 150mm Trays | 44.1kg |
| CRSM1/18/4S5D | 4 x 100mm, 5 x 150mm Trays | 42.8kg |
| CRSM1/18/5S5D | 5 x 100mm, 5 x 150mm Trays | 44.3kg |
| CRSM1/18/6S4D | 6 x 100mm, 4 x 150mm Trays | 43.9kg |
| CRSM1/18/8D   | 8 x 150mm Trays            | 43.3kg |

Easy disassembly & cleaning - page 38

# Caretray Racks - Stainless Steel, Mobile, Single & Double Column

#### **Features**

- Mobile, single & double column easy clean racks
- · Quick & easy disassembly, assists thorough & deep cleaning
- · Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recvclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each trav
- · Polyurethane ring buffers offer protection to both the rack & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 1782mm
  - Double Column 1005 x 550 x 1782mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 1757mm
- Double Column 915 x 455 (internal) at height 1757mm
- Travs
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

sales@bristolmaid.com

- Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



CRSM2/18/8S10D -shown with optional C004

| Item            | Description                  | Weight |
|-----------------|------------------------------|--------|
| CRSM2/18/4S12D  | 4 x 100mm, 12 x 150mm Trays  | 72.1kg |
| CRSM2/18/6S12D  | 6 x 100mm, 12 x 150mm Trays  | 77.1kg |
| CRSM2/18/8S10D  | 8 x 100mm, 10 x 150mm Trays  | 75.4kg |
| CRSM2/18/10S10D | 10 x 100mm, 10 x 150mm Trays | 78.4kg |
| CRSM2/18/12S8D  | 12 x 100mm, 8 x 150mm Trays  | 77.7kg |
| CRSM2/18/16D    | 16 x 150mm Trays             | 74.6kg |







### Features

Caretray Racks - Stainless Steel, Static, Single & Double Column

- Static, single & double column easy clean racks
- Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- · Adjustable feet allow levelling on uneven floors

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 536 x 536 x 1685mm
  - Double Column 995 x 536 x 1685mm
- Plastic Top
  - Single Column 455 x 455 (internal) at height 1645mm
- Double Column 915 x 455 (internal) at height 1645mm

#### Trays

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray -355 x 405 x 150mm (internal)

CRSS1/16/6S4D

| Item          | Description                | Weight |
|---------------|----------------------------|--------|
| CRSS1/16/2S6D | 2 x 100mm, 6 x 150mm Trays | 38.4kg |
| CRSS1/16/3S6D | 3 x 100mm, 6 x 150mm Trays | 40.8kg |
| CRSS1/16/4S5D | 4 x 100mm, 5 x 150mm Trays | 39.5kg |
| CRSS1/16/5S5D | 5 x 100mm, 5 x 150mm Trays | 41.0kg |
| CRSS1/16/6S4D | 6 x 100mm, 4 x 150mm Trays | 40.6kg |
| CRSS1/16/8D   | 8 x 150mm Trays            | 40.0kg |

Easy disassembly & cleaning - page 38

www.bristolmaid.com

• Trays, Shelves, Dividers & Labels - page 77

CRSS2/16/8S10D

| Item            | Description                  | Weight |
|-----------------|------------------------------|--------|
| CRSS2/16/4S12D  | 4 x 100mm, 12 x 150mm Trays  | 71.1kg |
| CRSS2/16/6S12D  | 6 x 100mm, 12 x 150mm Trays  | 74.0kg |
| CRSS2/16/8S10D  | 8 x 100mm, 10 x 150mm Trays  | 72.3kg |
| CRSS2/16/10S10D | 10 x 100mm, 10 x 150mm Trays | 75.3kg |
| CRSS2/16/12S8D  | 12 x 100mm, 8 x 150mm Trays  | 74.6kg |
| CRSS2/16/16D    | 16 x 150mm Trays             | 71.5kg |



# Caretray Trolleys - Mild Steel, Low Level, Single & Double Column





CTM1/08/3S1D - shown with optional CT/WF, SBH/15 & CT/ALC



Easy disassembly & cleaning - page 38

#### **Features**

- Low level, single & double column easy clean procedure trollevs
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows trav insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant travs & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- 100 / 150mm deep travs, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each trav
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 807mm
  - Double Column 1005 x 550 x 807mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 782mm
- Double Column 915 x 455 (internal) at height 782mm
- Travs
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

### Standard Colours

Paint - Spearmint

Item Description Weight CTM1/08/1S2D 1 x 100mm, 2 x 150mm Trays 20.9kg CTM1/08/3D 3 x 150mm Trays 21.2kg CTM1/08/3S1D 3 x 100mm, 1 x 150mm Trays 22.0kg CTM1/08/4S 4 x 100mm Trays 21.6kg

sales@bristolmaid.com

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67

- Carerail Mounted page 68
  Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77

CTM2/08/4S3D

| Item         | Description                | Weight |
|--------------|----------------------------|--------|
| CTM2/08/3S4D | 3 x 100mm, 4 x 150mm Trays | 36.8kg |
| CTM2/08/4S3D | 4 x 100mm, 3 x 150mm Trays | 36.4kg |
| CTM2/08/6D   | 6 x 150mm Trays            | 36.0kg |
| CTM2/08/6S2D | 6 x 100mm, 2 x 150mm Trays | 37.6kg |
| CTM2/08/7S1D | 7 x 100mm, 1 x 150mm Trays | 37.2kg |
| CTM2/08/8S   | 8 x 100mm Trays            | 36.8kg |

CTM2/08/4S3D - shown with optional CT/WF, SBH/15, CT/ALC

# Caretray Trolleys - Mild Steel, Standard Level, Single & Double Column





#### Easy disassembly & cleaning - page 38

#### **Features**

- Standard level, single & double column easy clean procedure trolleys
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 993mm
  - Double Column 1005 x 550 x 993mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 968mm
- Double Column 915 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

### **Standard Colours**

Paint - Spearmint

Description Weight Item CTM1/10/1S3D 1 x 100mm, 3 x 150mm Trays 27.5kg CTM1/10/3S2D 3 x 100mm, 2 x 150mm Trays 28.6kg CTM1/10/4D 4 x 150mm Travs 27.9kg CTM1/10/4S1D 4 x 100mm, 1 x 150mm Trays 28.3kg CTM1/10/6S 6 x 100mm 29.4kg

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
- Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77

CTM2/10/6S4D

| Item   | Description   | Weight                           |
|--|---|----------------------------------|
| CTM2/10/3S6D   | 3 x 100mm, 6 x 150mm Trays  | 46.0kg                           |
| CTM2/10/4S5D   | 4 x 100mm, 5 x 150mm Trays  | 45.6kg                           |
| CTM2/10/5S4D   | 5 x 100mm, 4 x 150mm Trays  | 45.3kg                           |
| CTM2/10/6S4D   | 6 x 100mm, 4 x 150mm Trays  | 46.7kg                           |
| CTM2/10/8D   | 8 x 150mm Trays   | 45.2kg                           |
| CTM2/10/8S2D   | 8 x 100mm, 2 x 150mm Trays  | 46.0kg                           |
| CTM2/10/4S5D<br>CTM2/10/5S4D<br>CTM2/10/6S4D<br>CTM2/10/8D | 4 x 100mm, 5 x 150mm Trays<br>5 x 100mm, 4 x 150mm Trays<br>6 x 100mm, 4 x 150mm Trays<br>8 x 150mm Trays | 45.6k<br>45.3k<br>46.7k<br>45.2k |

and the second

CTM2/10/6S4D - shown with optional CT/IV, CT/CR2, GMBH/15 & CT/ALC



sales@bristolmaid.com

# High level, single & double column easy clean procedure trolleys Quick & easy disassembly, assists thorough & deep cleaning

- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly) Recvclable polypropylene one-piece top & side panels, smooth
- surfaces & radiused edges assists cleaning (compatible with steam cleaning)
- Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
- Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant travs & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- 100 / 150mm deep travs, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each trav
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

### Dimensions (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 1410mm
  - Double Column 1005 x 550 x 1410mm
- Plastic Top

**Features** 

- Single Column 455 x 455 (internal) at height 1385mm
- Double Column 915 x 455 (internal) at height 1385mm
- Travs
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

## Standard Colours

Paint - Spearmint

Item Description CTM1/14/2S5D 2 x 100mm, 5 x 150mm Trays 36.2kg CTM1/14/3S4D 3 x 100mm, 4 x 150mm Trays 35.8kg 4 x 100mm, 3 x 150mm Trays CTM1/14/4S3D 35.4kg CTM1/14/5S3D 5 x 100mm. 3 x 150mm Travs 37.3kg CTM1/14/6D 6 x 150mm Trays 35.0kg CTM1/14/8S1D 8 x 100mm, 1 x 150mm Trays 37.2kg Easy disassembly & cleaning - page 38

- Locking & Security page 60

- Carerail Mounted Accessories page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



CTM2/14/10S6D - shown with optional CTD2/14SS (Double door fitted with Security Seal Lock)

| Item          | Description                 | Weight |
|---------------|-----------------------------|--------|
| CTM2/14/4S10D | 2 x 100mm, 10 x 150mm Trays | 60.6kg |
| CTM2/14/6S8D  | 6 x 100mm, 8 x 150mm Trays  | 62.8kg |
| CTM2/14/8S6D  | 8 x 100mm, 6 x 150mm Trays  | 62.3kg |
| CTM2/14/10S6D | 10 x 100mm, 6 x 150mm Trays | 65.0kg |
| CTM2/14/12D   | 12 x 150mm Trays            | 61.3kg |
| CTM2/14/16S2D | 16 x 100mm, 2 x 150mm Trays | 66.6kg |







CRMM1/18/6S4D

Easy disassembly & cleaning - page 38

| Item          | Description                | Weight |
|---------------|----------------------------|--------|
| CRMM1/18/2S6D | 2 x 100mm, 6 x 150mm Trays | 41.2kg |
| CRMM1/18/3S6D | 3 x 100mm, 6 x 150mm Trays | 43.6kg |
| CRMM1/18/4S5D | 4 x 100mm, 5 x 150mm Trays | 42.3kg |
| CRMM1/18/5S5D | 5 x 100mm, 5 x 150mm Trays | 43.8kg |
| CRMM1/18/6S4D | 6 x 100mm, 4 x 150mm Trays | 43.4kg |
| CRMM1/18/8D   | 8 x 150mm Trays            | 42.8kg |

#### **Features**

- Mobile, single & double column easy clean racks
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- · 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the rack & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 x 1782mm
  - Double Column 1005 x 550 x 1782mm
- · Plastic Top
  - Single Column 455 x 455 (internal) at height 1757mm
- Double Column 915 x 455 (internal) at height 1757mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

### **Standard Colours**

· Paint - Spearmint

sales@bristolmaid.com +

- Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



CRMM2/18/8S10D

| Item            | Description                  | Weight |
|-----------------|------------------------------|--------|
| CRMM2/18/4S12D  | 4 x 100mm, 12 x 150mm Trays  | 73.3kg |
| CRMM2/18/6S12D  | 6 x 100mm, 12 x 150mm Trays  | 76.3kg |
| CRMM2/18/8S10D  | 8 x 100mm, 10 x 150mm Trays  | 74.6kg |
| CRMM2/18/10S10D | 10 x 100mm, 10 x 150mm Trays | 77.6kg |
| CRMM2/18/12S8D  | 12 x 100mm, 8 x 150mm Trays  | 76.9kg |
| CRMM2/18/16D    | 16 x 150mm Trays             | 73.8kg |





| Item          | Description                | Weight |
|---------------|----------------------------|--------|
| CRMS1/16/2S6D | 2 x 100mm, 6 x 150mm Trays | 38.4kg |
| CRMS1/16/3S6D | 3 x 100mm, 6 x 150mm Trays | 39.5kg |
| CRMS1/16/4S5D | 4 x 100mm, 5 x 150mm Trays | 39.5kg |
| CRMS1/16/5S5D | 5 x 100mm, 5 x 150mm Trays | 41.0kg |
| CRMS1/16/6S4D | 6 x 100mm, 4 x 150mm Trays | 40.6kg |
| CRMS1/16/8D   | 8 x 150mm Trays            | 39.6kg |
|               |                            |        |



Easy disassembly & cleaning - page 38

# Caretray Racks - Mild Steel, Static, Single & Double Column

#### **Features**

- Mobile, single & double column easy clean racks
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- · Adjustable feet allow levelling on uneven floors

#### **Dimensions**

- Overall Dimensions
  - Single Column 536 x 536 x 1670-1685mm
  - Double Column 995 x 536 x 1670-1685mm
- Plastic Top
  - Single Column 455 x 455 (internal) at height 1645mm
  - Double Column 915 x 455 (internal) at height 1645mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

### **Standard Colours**

· Paint - Spearmint

sales@bristolmaid.com

www.bristolmaid.com

• Trays, Shelves, Dividers & Labels - page 77



| Item            | Description                  | Weight |
|-----------------|------------------------------|--------|
| CRMS2/16/4S12D  | 4 x 100mm, 12 x 150mm Trays  | 70.2kg |
| CRMS2/16/6S12D  | 6 x 100mm, 12 x 150mm Trays  | 73.2kg |
| CRMS2/16/8S10D  | 8 x 100mm, 10 x 150mm Trays  | 72.3kg |
| CRMS2/16/10S10D | 10 x 100mm, 10 x 150mm Trays | 75.3kg |
| CRMS2/16/12S8D  | 12 x 100mm, 8 x 150mm Trays  | 74.6kg |
| CRMS2/16/16D    | 16 x 150mm Trays             | 71.5kg |
|                 |                              |        |





# Locking

# **Features**

- · Doors available as either
  - Solid Grade Laminate Single / double column units
  - Clear Acrylic Single column units only
- Back panels are solid grade laminate, secures the contents of the trolley (removable to assist cleaning)
- Available for both single & double column units (mild & stainless steel)
  - Single Column Units Single door
  - Double Column Units Double door, features a return to ensure that only one lock is required to secure both doors
- 270-degree opening hinges allow the doors to fold flat against the sides of the trolley
- Four locking options
  - CAM Lock with Wing Handle Solid grade laminate / clear acrylic
  - Electronic Push Button Lock Solid grade laminate / clear acrylic
  - High Security Bolt Lock Solid grade laminate only
  - Security Seal Lock Solid grade laminate / clear acrylic
- Factory fitted



Door Open - CTM1/10/3S2D/CTD1/10EPB

SS

EPB

# **CAM Lock with Wing Handle**

# Features

- Supplied with two keys
  Wing handle can operate independently of the key
  Locks keyed to differ, alike on request



CTM1/08/3S1D/CTD1/08CA

| Description                   | Solid Grade Laminate | Clear Acrylic | Weight |
|-------------------------------|----------------------|---------------|--------|
| Low Level, Single Column      | CTD1/08CA            | CTDC1/08CA    | 2.7kg  |
| Low Level, Double Column      | CTD2/08CA            | -             | 5.3kg  |
| Standard Level, Single Column | CTD1/10CA            | CTDC1/10CA    | 3.5kg  |
| Standard Level, Double Column | CTD2/10CA            | -             | 6.8kg  |
| High Level, Single Column     | CTD1/14CA            | CTDC1/14CA    | 5.3kg  |
| High Level, Double Column     | CTD2/14CA            | -             | 10.3kg |



CTM1/10/3S2D/CTDC1/10CA



# **Electronic Push Button Lock**

# **Features**

- Uses a 10-button keypad
- Easy programming & code change
- Four operating codes with differing levels of authority Master / Sub-Master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code



CTM1/08/3S1D/CTD1/08EPB

| Description                   | Solid Grade Laminate | Clear Acrylic | Weight |
|-------------------------------|----------------------|---------------|--------|
| Low Level, Single Column      | CTD1/08EPB           | CTDC1/08EPB   | 2.9kg  |
| Low Level, Double Column      | CTD2/08EPB           | -             | 5.5kg  |
| Standard Level, Single Column | CTD1/10EPB           | CTDC1/10EPB   | 3.7kg  |
| Standard Level, Double Column | CTD2/10EPB           | -             | 7.0kg  |
| High Level, Single Column     | CTD1/14EPB           | CTDC1/14EPB   | 5.5kg  |
| High Level, Double Column     | CTD2/14EPB           | -             | 10.5kg |



www.bristolmaid.com

KitLock

# **High Security Bolt-Lock**

# Features

Supplied with two keys



CTD1/10HSB shown on CTM1/10/3S2D





| Description                   | Solid Grade Laminate | Weight |
|-------------------------------|----------------------|--------|
| Low Level, Single Column      | CTD1/08HSB           | 2.9kg  |
| Low Level, Double Column      | CTD2/08HSB           | 5.6kg  |
| Standard Level, Single Column | CTD1/10HSB           | 3.6kg  |
| Standard Level, Double Column | CTD2/10HSB           | 7.1kg  |
| High Level, Single Column     | CTD1/14HSB           | 5.3kg  |
| High Level, Double Column     | CTD2/14HSB           | 10.6kg |

www.bristolmaid.com

# **Security Seal Lock**

# **Features**

• Accepts standard plastic security seals (order separately)



| Description                   | Solid Grade Laminate | Clear Acrylic | Weight |
|-------------------------------|----------------------|---------------|--------|
| Low Level, Single Column      | CTD1/08SS            | CTDC1/08SS    | 2.6kg  |
| Low Level, Double Column      | CTD2/08SS            |               | 5.2kg  |
| Standard Level, Single Column | CTD1/10SS            | CTDC1/10SS    | 3.7kg  |
| Standard Level, Double Column | CTD2/10SS            |               | 6.7kg  |
| High Level, Single Column     | CTD1/14SS            | CTDC1/14SS    | 5.2kg  |
| High Level, Double Column     | CTD2/14SS            |               | 10.2kg |

CTM1/10/3S2D/CTDC1/10SS

# **Security Clamp**


### **Features**

- Allows the trolley to be secured to the 'fabric of the building'
- Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door of the trolley when closed & locked
- Recommended wall fitting height approx. 1000mm

# **Security Seals**

### **Features**

- Snap on plastic security seals
- Pack of 100





#### Important Note - Security Clamps

It is a requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area & broken into, but also encourages good nursing practice, giving the trolley a "permanent home".

'Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB-March 2005"

| Item    | Description    | Cable Length | Weight |
|---------|----------------|--------------|--------|
| MT/1/6F | Security Clamp | 950mm        | 0.75kg |

| ltem     | Description    |
|----------|----------------|
| 2PLCO026 | Security Seals |

www.bristolmaid.com

sales@bristolmaid.com

+44 (0)1258 484455

# Overbridge

### **Features**

- Attaches to the rear uprights of the trolley (increases overall height by 625mm)
- Available as either painted mild steel or self-colour stainless steel
- Provides a fixing location for tilt bins (order separately)
- Compatible with single column trolleys only
- Factory fitted only

### **Standard Colours**

- Paint Grey White
- Stainless Steel Self Colour



| Item     | Description                  | Weight |
|----------|------------------------------|--------|
| CT/OB/MS | Overbridge - Mild Steel      | 2.0kg  |
| CT/OB/SS | Overbridge - Stainless Steel | 2.0kg  |

# Tilt Bins

- Available with 4, 5 or 6 compartments
- Clear bins provide visible storage & assists organisation
- Bins tilt to 45° allowing easy access of contents, & can be removed to assist refiling / cleaning
- Label holder for easy identification
- Retaining bar secures bins in place, preventing from tipping forward during transit
- Supplied with a mild / stainless steel mounting

# Dimensions (w x d x h)

- TBIN4 143 x 127 x 190mm
- TBIN5 108 x 95 x 146mm
- TBIN6 89 x 67 x 95mm

### **Standard Colours**

- Plastic Cream
- Paint Grey White
- Stainless Steel Self Colour



Label Holder

Tilts to 45°

| Mild Steel | Stainless Steel | Description      | Dimensions (w x d x h) | Weight |
|------------|-----------------|------------------|------------------------|--------|
| TBIN4/MS   | TBIN4/SS        | Four Compartment | 600 x 171 x 207mm      | 2.3kg  |
| TBIN5/MS   | TBIN5/SS        | Five Compartment | 600 x 135 x 164mm      | 1.6kg  |
| TBIN6/MS   | TBIN6/SS        | Six Compartment  | 600 x 93 x 113mm       | 0.9kg  |

Racks

+44 (0)1258 484455 www.br

# **Adjustable Defibrillator Platform**

### **Features**

- Compatible with all popular defibrillators
- Defibrillator supported on a 365mm wide platform & is held securely in place using a clamp (adjusts using a handwheel 205 380mm)
- Platform mounted on an arm that swings 90 degrees from the worksurface & can be locked in position using a tensioning knob (rotates 360 degrees)
- Base of platform sits 130mm above the top of the trolley
- Fitted to the rear left hand upright of the trolley
- Maximum capacity 6kg
- Factory fitted

# **Standard Colours**

• Paint - Grey White

# Dimensions (w x d x h)

• 477 x 365 x 642mm (increases overall height by 170mm)

| Item     | Description                       | Weight |
|----------|-----------------------------------|--------|
| CT/DEFIB | Adjustable Defibrillator Platform | 5.5kg  |

MAXIMUM SAFE

# **Carerail Mounted**

# **Alcohol Gel Holder**

#### **Features**

- Accepts bottles with 90 x 90mm square base
- Fitted with a 15mm hook

### **Standard Colours**

• Paint - Grey White

# Dimensions (w x d x h)

• 95 x 110 x 115mm



| Item   | Description        | Weight |
|--------|--------------------|--------|
| CT/ALC | Alcohol Gel Holder | 0.3kg  |

# **Airway Tube Holder**

### **Features**

• Internal dimensions of holder 95 x 50 x 510mm

ol Maid 🙎

- Outer sleeve removable for cleaning
- Fitted with a 15mm hook

# **Standard Colours**

Paint - Grey White

### **Dimensions** (w x d x h)

• 140 x 70 x 555mm

# **Debris Bin & Bracket**

### **Features**

- Plastic debris bin & bracket
- Fitted with a 15mm hook

# **Standard Colours**

• Paint - Grey White

# Dimensions (w x d x h)

- Overall Dimensions 210 x 164 x 109mm
- Internal Dimensions 195 x 130 x 95mm





| Item   | Description          | Weight |
|--------|----------------------|--------|
| CT/BIN | Debris Bin & Bracket | 0.4kg  |

sales@bristolmaid.com +44

Item CT/AIR

+44 (0)1258 484455

Airway Tube Holder

Weight

1.4kg

# **Bowl Support & Bracket**

### **Features**

- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- Bowl retaining ring fitted with removable anti-static spiral buffers to support bowls

### **Standard Colours**

- Paint Grey White
- Stainless Steel

### Dimensions (w x d x h)

• 510 x 360 x 325mm

#### Accessories

- 2MECO079 360mm Diameter, Stainless Steel Bowl
- 2PLCO054 360mm Diameter, Polythene Bowl



| Item   | Description            | Weight |
|--------|------------------------|--------|
| CT/BOW | Bowl Support & Bracket | 2.5kg  |

# **Bucket Support & Bracket**

### **Features**

- Suitable for use with 15-litre graduated stainless steel buckets (order separately)
- Bucket retaining ring fitted with removable anti-static spiral buffers to support the bucket

### **Standard Colours**

- Paint Grey White
- Stainless Steel

### Dimensions (w x d x h)

• 510 x 360 x 335mm

#### Accessories

 2MECO081 - 15-litre Graduated Stainless Steel Bucket



| Item   | Description              | Weight |
|--------|--------------------------|--------|
| CT/BUC | Bucket Support & Bracket | 3.9kg  |

# **Document Holder**

# Features

- Suitable for A4 ring binders or foolscap notes (pocket size 274 x 63mm)
- Fitted with a 15mm hook

# **Standard Colours**

· Paint - Grey White

### **Dimensions** (w x d x h)

• 275 x 82 x 402mm



| Item     | Description     | Weight |
|----------|-----------------|--------|
| CT/DH/A4 | Document Holder | 1.7kg  |

# **Catheter Racks**

#### **Features**

- Available with 2 or 5 sections
- Clear perspex panel provides visibility of contents, & can be removed to assist cleaning
- Fitted with a 15mm hook

### Standard Colours

· Paint - Grey White

### **Dimensions** (w x d x h)

• Sections - 64 x 75 x 400mm (internal)



| Item   | Description  | Dimensions (w x d x h ) | Weight |
|--------|--------------|-------------------------|--------|
| CT/CR2 | Two Section  | 130 x 108 x 405mm       | 1.4kg  |
| CT/CR5 | Five Section | 325 x 108 x 405mm       | 3.0kg  |

# Cylinder Holder, C/CD/D Size

### **Features**

- Compatible with C/CD/D size cylinders (max diameter 100mm)
- Cylinder held in position using a chain, whilst a rubber pad protects wear points & reduces noise

### **Standard Colours**

· Paint - Grey White

#### Dimensions (w x d x h)

• 145 x 135 x 400mm



#### Oxygen cylinder not supplied

| Item   | Description                  | Weight |
|--------|------------------------------|--------|
| CT/CYD | Cylinder Holder, C/CD/D Size | 1.4kg  |

# **Infusion Pole & Bracket**

## **Features**

- Two hook infusion pole & bracket (pole length 1,160mm)
- Attaches to upper carerail
- Maximum capacity on hooks 2kg

### **Standard Colours**

- Paint Grey White
- Stainless Steel

### **Dimensions** (w x d x h)

• 238 x 104 x 1160mm



| Item  | Description             | Weight |
|-------|-------------------------|--------|
| CT/IV | Infusion Pole & Bracket | 1.6kg  |

Racks

ංඊ

**Caretray Trolleys** 

# **Phlebotomy Armrest**

#### **Features**

- Fully adjustable easy clean moulded polyurethane armrest
- Attaches to lower carerail

# Standard Colours

- · Paint Grey White
- Stainless Steel

### Dimensions (w x d x h)

Armrest Pad - 350 x 120 x 40mm



| Item    | Description        | Weight |
|---------|--------------------|--------|
| CT/PARM | Phlebotomy Armrest | 2.1kg  |

# **Sharps Box Bracket, Universal**

### **Features**

- · Compatible with the following brackets (customer to supply & fit the plastic brackets - fixings provided)
- Daniels Healthcare
  - POUDS® Wall-Bracket (up to 5L) SHARPSGUARD® mini. 1. 2.5. 3.75. 5 & com plus
  - POUDS® Wall-Bracket (7L & above) SHARPSGUARD® 7.11.5 & 22
- Frontier (Sharpsafe) 2. 3. 4 & 7
- Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS
- Fitted with a 15mm hook

## Standard Colours

· Paint - Grey White

# **Dimensions** (w x d x h)

• 75 x 32 x 152mm





Plastic brackets not supplied

Sharps Box Bracket, Universal

Description



Weiaht

0.2kg





| Item   | Description        | Weight |
|--------|--------------------|--------|
| CT/SPS | Suction Pump Shelf | 5.0kg  |

# **Suction Pump Shelf**

# **Features**

- · Compatible with all popular suction pumps
- Shelf provides an area 415 x 300mm (w x d)
- · 30mm upstand at front of shelf
- Corner bumpers offer protection to both the trollev & fabric of the building
- Floor to base of shelf at height 180mm
- Maximum capacity 15kg

# Standard Colours

· Paint - Grey White

# **Dimensions** (w x d x h)

• 460 x 245 x 390mm

Item CT/SBB/15

# Waste Bin & Bracket

#### **Features**

- 8 litre plastic swing bin & metal support frame
- Fitted with a 15mm hook

### Standard Colours

- · Paint Grey White
- Bin Cream

#### Dimensions (w x d x h)

• 8 Litre Plastic Swing Bin - 220 x 190 x 370mm



| Item  | Description         | Weight |
|-------|---------------------|--------|
| CT/WB | Waste Bin & Bracket | 0.6kg  |

# **Work Flap**

### **Features**

- Fold down flap (provides work area 300 x 455mm)
- Machined from solid grade laminate allowing for easy cleaning
- Maximum capacity 2kg

# **Standard Colours**

- Paint Grey White
- Solid Grade Laminate Grey White

## **Dimensions** (w x d x h)

• 325 x 465 x 70mm



### **Features**

- Suitable for use with roll type aprons, maximum width 415mm
- Fitted with a 15mm hook

# **Standard Colours**

• Paint - Grey White

#### Dimensions (w x d x h)

• 434 x 187 x 222mm





| Item     | Description           | Weight |
|----------|-----------------------|--------|
| DAR/1/15 | Apron Dispenser, Roll | 2.2kg  |

Item CT/WF

Work Flap

Weight

3.2kg

# Glove, Mask Box Holder, One Box Glove, Mask Box Holder, Three Box

· Plastic coated wire frame

Standard Colours

Dimensions (w x d x h)

• DG/3/15 - 290 x 141 x 405mm

· Paint - White

• Fitted with a 15mm square hook

· Boxes presented in landscape orientation

• Accommodates glove & mask boxes up to 150mm in width

4. 5. - 2 4

**Features** 

#### **Features**

- Accommodates glove & mask boxes up to 150mm in width
- · Box presented in portrait orientation
- Fitted with a 15mm square hook

#### **Standard Colours**

· Paint - Grey White

#### **Dimensions** (w x d x h)

• GMBH/15 - 143 x 122 x 177mm



| Item    | Description               | Weight |
|---------|---------------------------|--------|
| GMBH/15 | One Box, 15mm Square Hook | 0.6kg  |

| Item    | Description                 | Weight |
|---------|-----------------------------|--------|
| DG/3/15 | Three Box, 15mm Square Hook | 0.8kg  |

# **Sharps Box Holder**

### **Features**

- · Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a toggle
- Fitted with a 15mm hook

### **Standard Colours**

· Paint - Grey White

### Dimensions (w x d x h)

• 127 x 118 x 230mm



| Item   | Description       | Weight |
|--------|-------------------|--------|
| SBH/15 | Sharps Box Holder | 0.5kg  |



# Laptop Holder, Removeable, Swing Arm

### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with all standard height trolleys
- Fitted to the rear right upright of the trolley
- Laptop supported in a cradle, fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

\*Not compatible with Dell laptops using the N17 wedge lock system

# **Standard Colours**

• Paint - Grey White



| Item     | Description          | Weight |
|----------|----------------------|--------|
| CT/LHR/S | Removable, Swing Arm | 5.4kg  |

# Laptop Holder, Non-Removeable, Swing Arm

### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with all standard height trolleys
- Fitted to the rear right upright of the trolley
- Laptop supported in a cradle fitted with a steel mounting bracket, fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

\*Not compatible with Dell laptops using the N17 wedge lock system

# **Standard Colours**

· Paint - Grey White



www.bristolmaid.com

| Item      | Description              | Weight |
|-----------|--------------------------|--------|
| CT/LHNR/S | Non-Removable, Swing Arm | 5.7kg  |

Caretray Trolleys



C004



# **Polymer Castors**

# **Features**

- Offers reduced rolling resistance making the trolley easier to manoeuvre & lowers noise levels
- Smooth surfaces assist cleaning
- Supplied as a set of 4 castors
  - 2 x braking / 2 x non-braking
  - 4 x braking
  - 4 x non-braking
- Factory fitted only

| Item | Description   | Weight |
|------|---|--------|
| C004 | 100mm, Grey White Polymer, 2 x Braking / 2 x Non-Braking (set of 4) | 2.0kg  |
| C005 | 100mm, Grey White Polymer, 4 x Braking (set of 4)                   | 2.0kg  |
| C006 | 100mm, Grey White Polymer, 4 x Non-Braking (set of 4)               | 2.0kg  |

# **Ring Buffers**

## **Features**

- Set of four ring buffers
- Offers protection to both the trolley & fabric of the building
- Increases the length & width of trolley by approx. 80mm

www.bristolmaid.com

Factory fitted



| Item  | Description                | Weight |
|-------|----------------------------|--------|
| ZP010 | Ring Buffers (set of four) | 0.3kg  |

# Trays

# **Features**

- 100/150mm deep translucent high impact plastic trays with integral handle (max load per tray 5kg)
- Each tray is supplied with plastic clip on label cover & four removable dividers (2 x long & 2 x short)
- Trays are supplied without coloured label strips (see CT/LAB)

### Internal Dimensions (w x d x h)

- CT/TRAY/S 355 x 405 x 100mm
- CT/TRAY/D 355 x 405 x 150mm





| Item      | Description        | Weight |
|-----------|--------------------|--------|
| CT/TRAY/S | 100mm Shallow Tray | 1.5kg  |
| CT/TRAY/D | 150mm Deep Tray    | 1.8kg  |

# **Dividers**

# **Features**

- 100/150mm deep translucent dividers
- Supplied as a set of 2 x long & 2 x short



www.bristolmaid.com



| Item     | Description                              | Weight |
|----------|--|--------|
| CT/DIV/S | Set of 4 dividers for 100mm Shallow Tray | 0.6kg  |
| CT/DIV/D | Set of 4 dividers for 150mm Deep Tray    | 0.4kg  |

# Labels, Coloured

### **Features**


- Coloured labels Red / Yellow / Green / Blue / White
- Labels position beneath the clip-on label cover on the plastic trays
- Supplied as a set of 6 sheets



# Labels, Sepsis Six

### Features

- · Labels to support the delivery of 'The Sepsis Six'
- Printed Oxygen / Blood Cultures / Intravenous Antibiotics / Intravenous Fluid / Haemoglobin & Lactates / Urine Measurement
- Labels position beneath the clip-on label cover on the plastic trays
- Supplied as a set of 6 labels



| ltem   | Description       |
|--------|-------------------|
| CT/LAB | Labels - Coloured |

| Item       | Description         |
|------------|---------------------|
| CT/LAB/SEP | Labels - Sepsis Six |

sales@bristolmaid.com +44 (0)1258 484455

## Flat Shelf, Removable

#### Features

- Solid grade laminate (max load per shelf 5kg)
- Supported on the runner guides (non-pull out)

## **Standard Colours**

• Solid Grade Laminate - Grey White

## Dimensions (w x d x h)

• 413 x 457 x 4mm



CT/FS Flat Shelf - Removable 1.1kg

## Full Width Shelf, Removable

## **Features**

- Solid grade laminate (max load per shelf 5kg)
- Supported on the trolley frame

## **Standard Colours**

• Solid Grade Laminate - Grey White

## Dimensions (w x d x h)

• Overall Dimensions - 480 x 480 x 4mm





#### www.bristolmaid.com



# Chart Trolleys Workstations Room Dividers & Screens

- Fixed Height Trolleys 82
- Fixed Height Workstations 84
- Variable Height Workstations 86
  - Chart Boards 88
  - Document Holders 89
    - Room Dividers 90
    - Curtain Screens 92
      - Solid Screens 93





## **Caretray Trolley - Chart**

#### Features

- Standard level, single & double column easy clean fixed height chart trolleys
- Angle adjustable chartboard fitted with a spring clip (adjustable 0 50 degrees)
- Usable chart area...
  - Single Column 500 x 530mm (w x h)
  - Double Column 800 x 530mm (w x h)
- · Lower lip on chartboard to support pens etc.
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side (2 x L/H, 1 x R/H) support a wide range of accessories
- Drop down work flap fitted to R/H side
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- 100mm pressed steel castors, 2 off braking

## Dimensions (w x d x h)

- Overall Dimensions
  - Single Column 550 x 550 620 x 1040 1450mm
  - Double Column 1010 x 550 620 x 1040 1450mm
- Usable Chart Area
  - Single Column 500 x 530mm
  - Double Column 800 x 530mm
- Work Flap 420 x 290mm (w x d)
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

## **Standard Colours**

- Paint White
- Chartboard / Work Flap White

| Item             | Description                | Weight |
|------------------|----------------------------|--------|
| CTCHARM1/10/3S2D | 3 x 100mm, 2 x 150mm Trays | 33.2kg |
| CTCHARM1/10/4S1D | 4 x 100mm, 1 x 150mm Trays | 32.9kg |

## **Options / Accessories**

- Carerail Mounted page 68
- Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77



Angle adjustable chartboard (adjustable 0 - 50 degrees)

| Item             | Description                | Weight |
|------------------|----------------------------|--------|
| CTCHARM2/10/6S4D | 6 x 100mm, 4 x 150mm Trays | 53.2kg |
| CTCHARM2/10/8S2D | 8 x 100mm, 2 x 150mm Trays | 51.2kg |





## **Chart Workstations - Fixed Height**

## **Features**

- Fixed height chart workstation
- Angle adjustable chartboard fitted with a spring clip (0 50 degrees)
- Usable chart area 800 x 530mm (w x h)
- Lower lip on chartboard to support pens etc.
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- · Models available with
  - Chartboard only
  - Storage pocket
- Carerails fitted to each side of the trolley (2 x L/H, 1 x R/H) support a wide range of accessories
- Drop down work flap fitted to R/H side
- 100mm pressed steel castors, 2 off braking

## Dimensions (w x d x h)

- Overall Dimensions 1010 x 550 620 x 1040 1450mm
- Usable Chart Area 800 x 530mm
- Storage Pocket 825 x 450 x 150mm (internal)
- Work Flap 420 x 290mm (w x d)

## **Standard Colours**

- Paint White
- · Chartboard / Work Flap White

| Item      | Description                                      | Weight |
|-----------|--|--------|
| CWF/80    | Chart Workstation - Fixed Height                 | 23.0kg |
| CWF/80/SP | Chart Workstation - Fixed Height, Storage Pocket | 33.0kg |

## **Options / Accessories**

- Carerail Mounted page 68
- Castors & Ring Buffers page 76

Angle adjustable chartboard (0 - 50 degrees)





## **Chart Workstations - Variable Height**

#### **Features**

- Height adjustable using a gas damper (795 1075mm)
- Angle adjustable chartboard fitted with spring clips (adjustable 0 - 90 degrees)
- Chartboard available in two sizes (w x h)
  - Small 790 x 700mm
  - Large 790 x 800mm
- Lower lip on chartboard to support pens etc
- · Slide out shelf at front
- Document storage pocket on side to accept patients notes, ring binders, etc
- 80mm swivel castors, 2 off braking

## Dimensions (w x d x h)

- Overall Dimensions 800 x 700 / 800 x 1095mm
- Document Storage Pocket 75 x 255 x 360mm (internal)

#### **Standard Colours**

- Paint White
- · Chartboard White

## **Options / Accessories**

• ZP010 - 100mm Ring Buffers

| Item     | Description  | Weight |
|----------|--|--------|
| CB/100/A | Chart Workstation - Variable Height, Small, 700mm High Board | 40.0kg |
| CB/100/B | Chart Workstation - Variable Height, Large, 800mm High Board | 41.0kg |

## Other products that you may be interested in...

See Pages 374 & 375

## **Static Safe Tech Chairs**



+44 (0)1258 484455

www.bristolmaid.com

sales@bristolmaid.com

87



## Chart Boards - A4 & A3

## Features

- A4 one spring clip (holds one chart 220 x 280mm)
  A3 two spring clips (holds one chart 450 x 280mm or two charts 220 x 280mm)
- Available with a 30 or 50mm mild steel square hook

## **Standard Colours**

White

| Item    | Description   | Overall Dimensions (w x d x h) | Weight |
|---------|---------------|--------------------------------|--------|
| A4 Char | t Board       |                                |        |
| CB005   | A4, 30mm Hook | 230 x 80 x 370mm               | 0.5kg  |
| CB010   | A4, 50mm Hook | 230 x 100 x 370mm              | 0.5kg  |
| A3 Char | t Board       |                                |        |
| CB015   | A3, 30mm Hook | 460 x 80 x 370mm               | 0.8kg  |
| CB020   | A3, 50mm Hook | 460 x 100 x 370mm              | 0.8kg  |

sales@bristolmaid.com



| Item      | Description      | Overall Dimensions (w x d x h) | Weight |
|-----------|------------------|--------------------------------|--------|
| A4 Docume | ent Holder       |                                |        |
| DH/A4/WF  | A4, Wall Fitting | 280 x 70 x 400mm               | 1.4kg  |
| DH/A4/50  | A4, 50mm Hook    | 280 x 125 x 405mm              | 1.4kg  |
| A3 Docume | ent Holder       |                                |        |
| DH/A3/WF  | A3, Wall Fitting | 405 x 105 x 400mm              | 2.9kg  |
| DH/A3/50  | A3, 50mm Hook    | 405 x 160 x 405mm              | 3.0kg  |

## **Document Holders - A4 & A3**

## Features

- A4 size suitable for A4 ring binders or foolscap notes (pocket size 260 x 65mm)
- A3 size suitable for A3 notes & x-rays (pocket size 398 x 100mm)
- Open front design allows for easy patient identification
- Available for wall fitting or with a 50mm square hook

## **Standard Colours**

• Paint - Grey White





## **Room Dividers**

#### Features

- · Fixed height mobile room divider/workstation
- Ideal for use in critical care areas A&E. ICU. ICT. HDU
- Modular design allows the unit to be configured to suit your specific requirements
- Frame, worktop & upstand/panel available as...
  - Full width upstand
  - Full width upstand & right hand mounted panel
  - Full width upstand & left hand mounted panel
- Full width panel
- · Four configurations of lower units...
  - Three drawer unit, upper drawer fitted with an electronic push button lock (2 x 170, 1 x 285mm drawers)
  - Five drawer unit, upper drawer fitted with an electronic push button lock (5 x 125mm drawers)
  - · Cupboard with double doors & removable shelf
  - Cover panel (creates a space within the unit)
- Easy clean solid grade laminate work surface & up-stand/panel
- · Upper panel fitted with
  - Triple glove box holder
  - Alcohol gel holder (compatible with 500ml bottles, 70 x 70 x 225mm)
  - Apron roll dispenser
- 100mm ring buffers to protect the fabric of the building
- 100mm swivel castors, 2 off braking

## Dimensions (w x d x h)

- RD1 1768 x 698 x 1000mm
- RD2R/RD2L 1768 x 698 x 1588mm
- RD3 1768 x 698 x 1350mm


## **Standard Colours**

- · Paint White
- · Worksurface / Up-Stand / Panel White

## **Options / Accessories**

- CC Colour coded drawers available on request
- DTSF/RD/450/2 Dressing Trolley
  - · Seamlessly welded 304 grade stainless steel frame
  - · Fixed shelves with down turned sides
  - Shelf Heights 375 & 840mm
  - Max load per shelf 40kg

## Configure Your Workstation in Four Easy Steps...







MS/C4

+44 (0)1258 484455

sales@bristolmaid.com

## Screen - Curtain, Four Section, Folding

## Features

- PVC fire retardant curtain (machine washable)
- Available as four section
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors

## Dimensions (w x d x h)

- Overall Screen Height 1708mm
- Lowest Screen Height 288mm
- Folded 650 x 420 x 1790mm

## **Standard Colours**

- Paint White
- Curtain Grey

## **Spares**

• WS/SPARES/001 - Set of 4 replacement curtains

www.bristolmaid.com

Overall Dimensions (w x d x h)

Weight

14.6kg

## Mobile Screen - Solid, Three & Four Section, Folding

## Features

- Fire retardant folding compact laminate solid screens
- Available as either three or four section
- Smooth surface assists cleaning
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors

## Dimensions (w x d x h)

- Overall Screen Height 1680mm
- Lowest Screen Height 180mm
- Folded 600 x 435 x 1680mm

## **Standard Colours**

- · Paint Spearmint
- · Screen White



| Item  | Description                                   | Overall Dimensions (w x d x h) | Weight |
|-------|---|--------------------------------|--------|
| MS/S3 | Mobile Screen - Solid, Three Section, Folding | 1802 x 435 x 1680mm            | 20.8kg |
| MS/S4 | Mobile Screen - Solid, Four Section, Folding  | 2403 x 435 x 1680mm            | 27.8kg |





# **Concealment Trolleys**

- Variable Height Trolleys 96
  - Body Trays & Covers 97



## **Variable Height**

#### **Features**

- · Height adjustable
  - Fixed Tray 580 1280mm
  - Loose Tray 610 1310mm
- Scissor action provides vertical lift, ideal for use in confined areas where space is limited
- Battery powered actuator
  - Rechargeable battery is charged through the handset using a UK 3 pin plug
  - Audible warning when battery power is at 50%
  - Controlled using a handset
- Liquid retaining tray
  - Fixed Tray 1900 x 540 x 25mm
  - Loose Tray 2100 x 560 x 70mm
- Drop down side-rails, distance between rails 660mm
- Bumper bar / push handle at either end
- · Heavy duty reinforced PVC cover (locates over side rails)
- 200mm heavy duty castors, 2 x braking, 2 x directional locking
- Maximum capacity 190kg

## **Standard Colours**

• Paint - White



**Concealment Trolleys** 

2150 x 770 x 580 - 1280mm

2270 x 770 x 610 - 1310mm

www.bristolmaid.com

Variable Height, Fixed Body Tray

Variable Height, Loose Body Tray

VH/950

VH/950/LBT

## **Body Tray**

#### Features

- · Heavy duty body tray, mild steel frame & stainless steel tray
- Integral runners
- Fits most standard mortuary racks (check prior to ordering)
- Maximum capacity 190kg

## Dimensions (w x d x h)

• Overall Dimensions - 2100 x 550 x 70mm

## Standard Colours

· Paint - White

Item



## **Features**

· Heavy duty reinforced PVC cover (locates over side rails)





| Item     | Description                       | Weight |
|----------|-----------------------------------|--------|
| 2FACO024 | Cover - PVC Heavy Duty Reinforced | 5.0kg  |

www.bristolmaid.com

Shown with VH/950



# Cylinder Storage & Transportation

- Cylinder Trolleys F/G Size 100
  - Cylinder Trolleys J Size 101
- Bulk Transportation Trolleys 102
- Cylinder Trolleys E, CD/D, Size 102
- Wall Brackets & Cylinder Holders 103

## Cylinder Trolley F/G Size

#### **Features**

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinder
- 160mm anti-static wheels

## Dimensions (w x d x h)

• 510 x 405 x 1145mm

#### **Standard Colours**

Paint - Grey



| Item    | Description | Weight |
|---------|-------------|--------|
| CYT/FG1 | F/G Size    | 10.0kg |

## Cylinder Trolley with Stabilisers F/G Size

#### **Features**

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinder
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 160mm anti-static wheels, 100mm swivel castors

## Dimensions (w x d x h)

• 520 x 700 x 1145mm

#### **Standard Colours**

Paint - Grey



| Item       | Description | Weight |
|------------|-------------|--------|
| CYT/FG1/ST | F/G Size    | 14.0kg |

www.bristolmaid.com

## **Cylinder Trolley J** Size

#### **Features**

- Compatible with J (240 cu ft), AK, L & HL cylinders, max diameter 245mm, height 1550mm
- One piece construction
- Drop over chains secure the cylinder
- 160mm anti-static wheels

## Dimensions (w x d x h)

• 510 x 475 x 1290mm

#### Standard Colours

Paint - Grey



| ľ | tem    | Description |
|---|--------|-------------|
| ( | CYT/J1 | J Size      |

## **J** Size

- **Features**
- Compatible with J (240 cu ft), AK, L & HL cylinders, max diameter, max diameter 245mm, height 1550mm

**Cylinder Trolley with Stabilisers** 

- One piece construction
- Drop over chains secure the cylinder
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 160mm anti-static wheels, 100mm swivel castors

## Dimensions (w x d x h)

• 510 x 740 x 1290mm

#### **Standard Colours**

· Paint - Grey



| Item      | Description | Weight |
|-----------|-------------|--------|
| CYT/J1/ST | J Size      | 14.8kg |

+44 (0)1258 484455

Weight 10.5kg

## **Bulk Transportation Trolley** D & E Cylinders, 12 x Capacity

#### **Features**

- Compatible with D, AD, CD, DD, PD, RD, ZD & E (24 cu ft) cylinders (12 x capacity) max diameter 105mm, height 870mm
- One piece construction
- Cylinders locate into individual sections, supported on a non slip rubber base
- 125mm swivel castors, 2 off braking

#### Dimensions (w x d x h)

• 460 x 590 x 1005mm

#### Standard Colours

· Paint - Grey



| ltem     | Description                    | Weight |
|----------|--------------------------------|--------|
| CYT/DE12 | D & E Cylinders, 12 x Capacity | 23kg   |

## **Cylinder Trolley** F Size

#### **Features**

- · Compatible with E (24 cu ft) cylinders, max diameter 115mm, height 900mm
- One piece construction
- · Recessed base plate and drop over chain to secure the cylinder
- 40mm hook allows trolley to hang at bed end, a useful feature when transporting patients
- 100mm anti-static wheels

## Dimensions (w x d x h)

• 320 x 235 x 980mm

#### Standard Colours

· Paint - Grey



| Item   | Description | Weight |
|--------|-------------|--------|
| CYT/E1 | E Size      | 4.3kg  |

## **Cylinder Trolley CD/D Size**

## **Features**

- Compatible with CD & D cylinders, max diameter 115mm, height 540mm
- One piece construction
- Recessed base plate and drop over chain to secure the cvlinder
- 40mm hook allows trolley to hang at bed end, a useful feature when transporting patients
- 100mm anti-static wheels

## Dimensions (w x d x h)

• 320 x 235 x 980mm

## Standard Colours

· Paint - Grey



| Item    | Description | Weight |
|---------|-------------|--------|
| CYT/CD1 | CD/D Size   | 4.3kg  |

sales@bristolmaid.com


+44 (0)1258 484455

## Wall Brackets Cylinder

#### Features

- Max cylinder diameters
- D/E 112mm, F 148mm, G 185mm, J 235mm
- One piece construction
- Drop over chains to secure the cylinder
- · Holes for wall mounting (fixings not provided)
- Should be positioned on the wall at approx half cylinder height

## **Standard Colours**



| Item    | Description       | Dimensions (w x d x h) | Weight |
|---------|-------------------|------------------------|--------|
| CYB/DE2 | 2 x D/E Cylinders | 359 x 107 x 88mm       | 1.1kg  |
| CYB/F2  | 2 x F Cylinders   | 435 x 126 x 88mm       | 1.4kg  |
| CYB/G2  | 2 x G Cylinders   | 511 x 145 x 88mm       | 1.7kg  |
| CYB/J2  | 2 x J Cylinders   | 613 x 170 x 88mm       | 2.1kg  |

## Cylinder Holder D & E Size, Vertical

#### Features

- Compatible with...
  - D Size D, AD, CD, DD, PD, RD & ZD cylinders (max 112mm diameter)
  - E Size E (max 112mm diameter)
- One piece construction
- Carry handle and extended hook
- Suitable for beds with 25mm wide side/end members
- Holder can be hung and/or bolted to 25mm section bed-frame

#### Dimensions (w x d x h)

- CYH/D1 200 x 160 x 290mm
- CYH/E1 200 x 160 x 595mm

CYH/D1

+44 (0)1258 484455

Description

D Size, Vertical

E Size, Vertical

#### **Standard Colours**

• Paint - Grey

Item

CYH/D1

CYH/E1

## Cylinder Holder Wall Plate

## **Features**

- Painted stainless steel wall plate to support cylinder holders
- Compatible with CYH/D1, CYH/E1, CYHP/D1 & CYHP/E1

## Dimensions (w x d x h)

• 240 x 47 x 254mm

## **Standard Colours**

• Paint - Grey





CYH/F1

Weight

1.4kg

2.5kg



## **Dispensers, Holders & Racks**

- Apron Dispensers, Pocket Mask Holders 106
  - Glove / Mask Box Holders 107
    - Alcohol Gel Holders 108
- Sharps Box Holders, Urine Bottle Holder, Catheter Bag Holders 109
  - Catheter Racks 110
  - Drainage Racks 111

## **Apron Dispensers - Roll**

#### **Features**

- Apron roll holders, maximum width 415mm
- Holes for wall mounting (fixings not provided)

## Dimensions (w x d x h)

- One Roll 420 x 168 x 268mm
- Three Roll 420 x 168 x 615mm

#### Standard Colours

· Paint - White



| Item  | Description | Weight |
|-------|-------------|--------|
| DAR/1 | One Roll    | 2.0kg  |
| DAR/3 | Three Rolls | 8.2kg  |

## **Apron Dispenser - Flat Pack**

## Features

- Holder to accept flat packed aprons, maximum width 400mm
- Packs of aprons fold over the centrally positioned divider
- Holes for wall mounting (fixings not provided)

## Dimensions (w x d x h)

• 410 x 105 x 290mm

## Standard Colours

· Paint - White



| Item  | Description | Weight |
|-------|-------------|--------|
| DAF/1 | Flat Pack   | 2.0kg  |

## Pocket Mask Holder

### **Features**

- Holder to accept Laerdal & Ambu Pocket Masks (hard cases)
- Allows the case to be guickly removed from the holder
- Supplied with self adhesive tape/holes for wall mounting (fixings not provided)

## Dimensions (w x d x h)

• 110 x 40 x 110mm

#### Standard Colours

· Paint - White



Pocket mask not supplied



| Item   | Description        | Weight |
|--------|--------------------|--------|
| PMH005 | Pocket Mask Holder | 0.1kg  |

& Racks

**Dispensers, Holders** 

sales@bristolmaid.com

## Glove / Mask Box Holders Three Box

#### **Features**

- Plastic coated wire frame
- Accommodates glove & mask boxes up to 150mm in width
- Boxes presented in landscape orientation
- Available for wall mounting or fitted with a 15mm square hook
   (fixings not provided)

## Dimensions (w x d x h)

- DG/3 290 x 115 x 405mm
- DG/3/15 290 x 141 x 405mm

## **Standard Colours**

· Paint - White



| Item    | Description                 | Weight |
|---------|-----------------------------|--------|
| DG/3    | Three Box, Wall Mounting    | 0.7kg  |
| DG/3/15 | Three Box, 15mm Square Hook | 0.8kg  |

## Glove / Mask Box Holders One Box

## Features

- Accommodates glove & mask boxes up to 150mm in width
- · Box presented in portrait orientation
- Available for wall mounting or fitted with a 15mm square hook (fixings not provided)

## Dimensions (w x d x h)

- GMBH/WM 143 x 100 x 152mm
- GMBH/15 143 x 122 x 177mm

## **Standard Colours**

· Paint - Grey White



## **Features**

- Packs of sterile gloves separated by two
   adjustable dividers
- Full width card holder for identification
- Clear front panel assists with stock level inspection
- Hinged lid assists replenishment, whilst top loading ensures stock rotation
- Holes for wall mounting (fixings not provided)

## Dimensions (w x d x h)

• 565 x 260 x 315mm

## **Standard Colours**

• Paint - White



| Item    | Description              | Weight |
|---------|--------------------------|--------|
| SS/48/3 | Three Column, Adjustable | 7.0kg  |

Item

GMBH/WM

GMBH/15

Weight

0.5ka

Description

One Box. Wall Mounted

One Box, 15mm Square Hook

## **Alcohol Gel Holders - Wire Frame**

#### **Features**

- Compatible with popular 500ml bottles, 70mm x 70mm x 225mm (not supplied)
- Stainless steel construction
- · Plastic drip tray
- Available with holes for wall mounting (fixings not provided) or with a 30mm square hook

### **Standard Colours**

Stainless Steel



## Item Description Overall Dimensions (w x d x h) Weight GH005 Lever Operated, Wall Mounted, Angled 77 x 300 x 278mm 0.4kg GH015 30mm Hook 77 x 117 x 180mm 0.3kg

sales@bristolmaid.com

## Syringe Dispenser - Five Column, Adjustable

#### **Features**

- Syringes separated by four adjustable dividers
- Full width card holder for identification
- Clear front panel assists with stock level inspection
- Hinged lid assists replenishment, whilst top loading ensures stock rotation
- Holes for wall mounting (fixings not provided)

## Dimensions (w x d x h)

• 565 x 317 x 320mm

#### **Standard Colours**

• Paint - White



| Item  | Description             | Weight |
|-------|-------------------------|--------|
| SS/36 | Five Column, Adjustable | 9.2kg  |

## **Sharps Box Holder One Box, Wall Mounted**

## Features

- · Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a togale
- Holes for wall mounting (fixings not provided

## Dimensions (w x d x h)

• 127 x 100 x 190mm

#### Standard Colours

· Paint - White



## **Urine Bottle Holder** Hanging

### **Features**

- · Compatible with bottles up to 115mm Diameter
- Plastic coated wire frame
- 30mm square hook

## Dimensions (w x d x h)

• 150 x 182 x 282mm

## Standard Colours

· Plastic - White



| Item | Description  | Weight |
|------|--------------|--------|
| UB/9 | Urine Bottle | 0.3kg  |

## **Catheter Bag Holders Free Standing/Hanging**

## **Features**

- · Plastic coated wire frame
- · Accepts most widely used bags
- Available freestanding or hanging (30mm square hook)
- · Upper clip to hold the inlet tube
- Lower clip angles the drainage tube upwards to prevent drips should the tap leak

## Standard Colours

· Plastic - White

www.bristolmaid.com



| Item  | Description  | Overall Dimensions (w x d x h) | Weight |
|-------|--------------|--------------------------------|--------|
| UB005 | Freestanding | 250 x 250 x 550mm              | 0.6kg  |
| UB010 | Hanging      | 250 x 77 x 410mm               | 0.4kg  |

## **Catheter Racks - Two & Five Section**

#### Features

• Internal dimensions of each section 65 x 75 x 405mm

RC/2

- Removable clear perspex panel to aid cleaning
- Available with holes for wall mounting (fixings not provided) or with a 15mm square hook

#### **Standard Colours**

• Paint - White



| Item Description |                                     | Overall Dimensions (w x d x h) | Weight |
|------------------|-------------------------------------|--------------------------------|--------|
| Two Section      |                                     |                                |        |
| RC/2             | Wall Fitting                        | 135 x 85 x 410mm               | 1.3kg  |
| RC/2/15          | 15mm Hook                           | 135 x 115 x 410mm              | 1.4kg  |
| Five Section     |                                     |                                |        |
| RC/5             | Wall Fitting                        | 320 x 85 x 410mm               | 3.3kg  |
| RC/5/15          | RC/5/15 15mm Hook 320 x 115 x 410mm |                                | 3.4kg  |

sales@bristolmaid.com

## **Catheter Racks - One Section**

## Features

- Outer sleeve removable for cleaning
- 20mm square hook

## Dimensions (w x d x h)

• 185 x 55 x 555mm

#### Standard Colours

· Paint - White



| Item     | Description | Weight |
|----------|-------------|--------|
| SS/40/MS | One Section | 2.0kg  |

www.bristolmaid.com


## **Drainage Racks - Urinal**

#### Features

- Holds five urine bottles
- Stainless steel construction
- Supplied with a stainless steel drip tray
- Holes for wall mounting (fixings not provided)

## Dimensions (w x d x h)

• 535 x 205 x 415mm



## Drainage Racks - Bedpan

## Features

- Holds five bedpans
- Stainless steel construction
- Supplied with a stainless steel drip tray
- Holes for wall mounting (fixings not provided)

## Dimensions (w x d x h)

• 530 x 295 x 375mm



| Item     | Description             | Weight |
|----------|-------------------------|--------|
| RD/U5/SS | Urinal, Stainless Steel | 1.7kg  |

| Description             | Weight |
|-------------------------|--------|
| Bedpan, Stainless Steel | 3.3kg  |

Item RD/B5/SS



# **Dressing / Instrument Trolleys**

- Removable Shelves Stainless Steel 114
  - Removable Shelves Mild Steel 115
- Fixed Shelves, Sides Up Stainless Steel 116
- Fixed Shelves, Sides Down Stainless Steel 117
  - Drawer Units 119
  - Options & Accessories 120



## Removable Shelves, Stainless Steel, 450mm Deep

## Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Flat or reversible 304 grade stainless steel shelves, cushioned on anti-static spiral buffers (removable for cleaning)
- Reversible shelves allow the end user to choose between a flat shelf or one with three upturned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

## **Options / Accessories**

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C013 75mm, Pressed Steel Castors (4 x braking) refer to page 121
- Drawer Units page 119 Factory Fitted Only
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only

| Flat Shelves | Reversible Shelves | Size    | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|--------------|--------------------|---------|--------------------|--------------------------------|--------|
| DTSR/450     | DTSR/450/SR        | Small   | 450 x 450mm        | 470 x 470 x 890mm              | 11.5kg |
| DTSR/600     | DTSR/600/SR        | Medium  | 600 x 450mm        | 620 x 470 x 890mm              | 16.0kg |
| DTSR/750     | DTSR/750/SR        | Large   | 750 x 450mm        | 770 x 470 x 890mm              | 18.0kg |
| DTSR/900     | DTSR/900/SR        | X-Large | 900 x 450mm        | 920 x 470 x 890mm              | 20.0kg |

sales@bristolmaid.com

```
www.bristolmaid.com
```



## Removable Shelves, Mild Steel, 450mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- · Seamlessly welded mild steel frame
- Flat or reversible 304 grade stainless steel shelves, cushioned on anti-static spiral buffers (removable for cleaning)
- Reversible shelves allow the end user to choose between a flat shelf or a shelf with three upturned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

#### **Standard Colours**

Paint - White

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C013 75mm, Pressed Steel Castors (4 x braking) refer to page 121
- Drawer Units page 119 Factory Fitted Only
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only

| Flat Shelves | Reversible Shelves | Size   | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|--------------|--------------------|--------|--------------------|--------------------------------|--------|
| DTMR/450     | DTMR/450/SR        | Small  | 450 x 450mm        | 470 x 470 x 890mm              | 11.5kg |
| DTMR/600     | DTMR/600/SR        | Medium | 600 x 450mm        | 620 x 470 x 890mm              | 16.0kg |



## Fixed Shelves, Sides Up, Stainless Steel, 450mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with three sides up & one down
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

#### **Options / Accessories**

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C013 75mm, Pressed Steel Castors (4 x braking) refer to page 121
- Drawer Units page 119 Factory Fitted Only
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only

| Item     | Size    | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|----------|---------|--------------------|--------------------------------|--------|
| DTSF/450 | Small   | 450 x 450mm        | 470 x 470 x 890mm              | 11.5kg |
| DTSF/600 | Medium  | 600 x 450mm        | 620 x 470 x 890mm              | 16.0kg |
| DTSF/750 | Large   | 750 x 450mm        | 770 x 470 x 890mm              | 18.0kg |
| DTSF/900 | X-Large | 900 x 450mm        | 920 x 470 x 890mm              | 20.0kg |

sales@bristolmaid.com

DTSF/450



## Fixed Shelves, Sides Down, Stainless Steel, 450mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trollevs for instruments. dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- · Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

#### **Options / Accessories**

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C013 75mm, Pressed Steel Castors (4 x braking) refer to page 121
- Drawer Units page 119 Factory Fitted Only
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only

| sales@bristolmaid.com |
|-----------------------|
|                       |

18.0kg

21.0kg

770 x 470 x 870mm

920 x 470 x 870mm

750 x 450mm

DTSF/750/2/SD

Large DTSF/900/2/SD X-Large 900 x 450mm



## Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep

#### Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel anti-static non-marking castors
- \* Drawers not suitable for this product

- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking) Factory Fitted Only - refer to page 121
- C015 100mm, Pressed Steel Castors (4 x braking) Factory Fitted
   Only refer to page 121
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers Factory Fitted Only refer to page 121

| Item          | Size   | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|---------------|--------|--------------------|--------------------------------|--------|
| DTSF/66/2/SD  | Small  | 600 x 600mm        | 620 x 620 x 870mm              | 17.5kg |
| DTSF/96/2/SD  | Medium | 900 x 600mm        | 920 x 620 x 870mm              | 24.0kg |
| DTSF/126/2/SD | Large  | 1200 x 600mm       | 1220 x 620 x 870mm             | 29.6kg |





2HH-7011100

#### DTSR/450/DUS/R

| Item     | Description                                   | Trolley Compatibility | Weight |
|----------|---|-----------------------|--------|
| Removal  | ble Shelf Trolleys                            |                       |        |
| DUS/R    | Drawer Unit - Small, Removable Shelf Trolleys | 450mm & 600mm         | 5.1kg  |
| DUL/R    | Drawer Unit - Large, Removable Shelf Trolleys | 750mm & 900mm         | 6.9kg  |
| Fixed Sh | elf Trolleys                                  |                       |        |
| DUS/F    | Drawer Unit - Small, Fixed Shelf Trolleys     | 450mm & 600mm         | 4.7kg  |
| DUL/F    | Drawer Unit - Large, Fixed Shelf Trolleys     | 450mm & 600mm         | 6.5kg  |

## **Drawer Units**

#### Features

- Compatible with removable & fixed shelf ranges, with exception of 600mm deep trolleys
- · Available in two sizes
  - Small 300 x 400mm
  - Large 600 x 400mm
- Each drawer unit is supplied with a 100mm deep HTM71 modular basket
- Stainless steel basket support & drawer carriage
- Precision telescopic runners
- Solid grade laminate drawer front
- Maximum load 5kg
- Trolleys accept only one tier of drawers
- Drawers are not lockable
- Factory fitted & must be ordered at the same time as the trolley

#### Internal Dimensions (w x d x h)

- Small 235 x 335mm (internal)
- Large 535 x 335mm (internal)

#### **Options / Accessories**

- 2HH-7014100 Basket, 300 x 400 x 100mm
- 2HH-7012100 Basket, 600 x 400 x 100mm
- 2HH-7811206N Divider, 300 x 100mm
- 2HH-7811205N Divider, 400 x 100mm
- 2HH-7811204N Divider, 600 x 100mm
- Additional Trays, Baskets, Dividers & Accessories page 206

#### sales@bristolmaid.com

## **Reversible Shelves**

#### Features

- Allows the end user to choose between a flat shelf or one with three upturned sides
- Removable for cleaning
- Compatible with DTMR & DTSR ranges only



Description

**Reversible Shelves** 

## **Clamping Infusion Stands - Trolley**

#### **Features**

- Clamp secures to 25mm vertical round tube (designed to fit our range of dressing trolleys)
- When fitting to other trolleys, we strongly recommend that a risk assessment is carried out to assess stability
- When fitted protrudes 45mm
- Supplied with a two pigtail hook infusion pole
- Maximum capacity on hooks 2kg

#### Dimensions (w x d x h)

• 55 x 220 x 1200mm





| Item | Description                        | Weight |
|------|------------------------------------|--------|
| DP/C | Clamping Infusion Stands - Trolley | 1.8kg  |

Item

SR

Weight

5-10kg

## Castors

#### Features

- Polymer & pressed steel castors
- Anti-static & non-marking
- Supplied as a set of 4 castors
- 80mm polymer castors increase the height of the trolley by approx. 5mm
- 100mm polymer castors increase the height of the trolley by approx. 30mm
- 100mm steel castors increase the height of the trolley by approx. 30mm
- · Factory fitted



#### Item Description

- C001 Castor 80mm, Grey Polymer, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C002 Castor 100mm, Grey Polymer, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C003 Castor 100mm, Pressed Steel, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C013 Castor 75mm, Pressed Steel (4 x braking) Factory Fitted Only
- C015 Castor 100mm, Pressed Steel (4 x braking) Factory Fitted Only

| <b>Ring B</b> | uffers |
|---------------|--------|
|---------------|--------|

#### Features

- · Set of four ring buffers
- Offers protection to both the trolley & fabric of the building
- Increases the length & width of trolley by approx. 80mm
- Factory fitted





| Item  | Description                                      | Weight |
|-------|--|--------|
| ZP005 | Ring Buffers (set of four) - Factory Fitted Only | 0.3kg  |



# **Drug & Medicine Dispensing**

- Blister Pack (MDS) Trolleys 136
- Nomad<sup>®</sup> Cassette Trolleys 137
- Biodose<sup>®</sup> & MultiMeds Trolleys 138
  - Medinoxx Trolleys 139
    - Security Clamp 142
  - Distribution Trolleys 143

- Dispensing Trolleys 124
  - Laptop Holders 129
- CLIQ<sup>®</sup> & RFID Solutions 131
- Patient Administration Trolleys 132
  - Drawer Trolleys 134
    - Tray Trolleys 135





CTDISPM1/10/3S2D/HSB

| Description   | Weight   |
|---|--|
|   |  |
| 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 35.5kg   |
| 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 35.2kg   |
| 6 x 100mm Trays, High Security Bolt Lock            | 36.3kg   |
|   | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock<br>4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock |

## **Caretray Trolley - Dispensing**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines
- Suitable for use with Electronic Prescribing & Medicines Administration (EPMA) & Electronic Patient Records (EPR)
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Waste Bin & Bracket, 8 Litre
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking
- When not in use, trolleys storing medicines should be secured to the wall using a Security Clamp (see locking & security accessories see page 142)

#### Dimensions (w x d x h)

#### **Standard Colours**

• Paint - Spearmint

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Carerail Mounted page 68

- Laptop Holders page 74
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77





| CTDISPM1/ | /10/3S2D/EI | B |
|-----------|-------------|---|
|           |             |   |

| Item                     | Description   | Weight |
|--------------------------|---|--------|
| Electronic Push Button L | ock   |        |
| CTDISPM1/10/3S2D/EPB     | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 35.5kg |
| CTDISPM1/10/4S1D/EPB     | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 35.2kg |
| CTDISPM1/10/6S/EPB       | 6 x 100mm Trays, Electronic Push Button Lock            | 36.3kg |

sales@bristolmaid.com



| Item                 | Description                         | Weight |
|----------------------|-------------------------------------|--------|
| <b>High Security</b> | Bolt Lock                           |        |
| MT100                | Small, High Security Bolt Lock      | 26.0kg |
| MT110                | Medium, High Security Bolt Lock     | 35.0kg |
| MT120                | Large, High Security Bolt Lock      | 42.5kg |
| Electronic Pus       | sh Button Lock                      |        |
| MT100CL              | Small, Electronic Push Button Lock  | 26.0kg |
| MT110CL              | Medium, Electronic Push Button Lock | 35.0kg |
| MT120CL              | Large, Electronic Push Button Lock  | 42.5kg |

## **Dispensing Trolleys - Wooden, Standard Capacity**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines
  Lockable hinged lids...
  - Upper Assisted using a gas strut, providing document support
  - Lower Provides a work area, at height 980mm (max load 4.0kg)
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Hinged work flap provides a work area at height 930mm (max load 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Two slide out storage / waste bins
- 100mm pressed steel castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

- Overall...
  - Small 750 x 490 x 1130mm
  - Medium 890 x 490 x 1130mm
  - Large 1030 x 525 x 1130mm
- Internal...
  - Small 530 x 427 x 263mm sloping to 112mm at front
  - Medium 683 x 427 x 263mm sloping to 112mm at front
  - Large 805 x 461 x 263mm sloping to 112mm at front

#### **Standard Colours**

- Frame Grey
- Shelves / Dividers White
- Cabinet Spearmint
- Shelf / Flap Grey

- MT/1/6F Security Wall Clamp
- ZP005 Ring Buffers
- MT/LHR/S Laptop Holder Removable, Swing Arm
- MT/LHNR/S Laptop Holder Non-Removable, Swing Arm



| Item                        | Description                         | Weight |  |
|-----------------------------|-------------------------------------|--------|--|
| High Secur                  | ity Bolt Lock                       |        |  |
| MT220                       | Medium, High Security Bolt Lock     | 42.0kg |  |
| MT230                       | Large, High Security Bolt Lock      | 50.0kg |  |
| Electronic Push Button Lock |                                     |        |  |
| MT220CL                     | Medium, Electronic Push Button Lock | 42.0kg |  |
| MT230CL                     | Large, Electronic Push Button Lock  | 50.0kg |  |

## **Dispensing Trolleys - Wooden, High Capacity**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines
  Lockable hinged lids...
  - Upper Assisted using a gas strut, providing document support
  - Lower Provides a work area, at height 921mm (max load 4.0kg)
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Hinged work flap provides a work area at height 850mm (max load 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Two slide out storage / waste bins
- 125mm pressed steel castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

- Overall...
  - Medium 760 x 560 x 1290mm
  - Large 905 x 560 x 1290mm
- Internal...
  - Medium 530 x 461 x 513mm sloping to 137mm at the front
  - Large 683 x 461 x 513mm sloping to 137mm at the front

#### **Standard Colours**

- Frame Grey
- · Shelves / Dividers White
- Cabinet Spearmint
- Shelf / Flap Grey

- MT/1/6F Security Wall Clamp
- MT/LHR/S Laptop Holder Removable, Swing Arm
- MT/LHNR/S Laptop Holder Non-Removable, Swing Arm



## Features

- Trolleys suitable for the storage & dispensing of drugs & medicines
- Lockable hinged lid provides document support
- Fitted with a high security bolt lock (supplied with two keys)

**Dispensing Trolleys - Metal, Standard Capacity** 

- Hinged writing flap provides a work area at height 907mm (max capacity 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Two slide out storage / waste bins
- 100mm pressed steel castors
- When not in use, trolley should be secured to the building see optional  $\ensuremath{\mathsf{MT}}\xspace{-1.5mu}$

#### Dimensions (w x d x h)

- Overall...
  - Small 750 x 459 x 1112mm
  - Medium 891 x 459 x 1112mm
  - Large 1030 x 492 x 1112mm
- Internal...
  - Small 489 x 457 x 292mm sloping to 127mm at the front
  - Medium 706 x 457 x 292mm sloping to 127mm at the front
  - Large 832 x 490 x 292mm sloping to 127mm at the front

#### **Standard Colours**

- Frame Grey
- Cabinet White
- Shelf / Flap Spearmint

- MT/1/6F Security Wall Clamp
- ZP005 Ring Buffers
- MT/LHR/S Laptop Holder Removable, Swing Arm
- MT/LHNR/S Laptop Holder Non-Removable, Swing Arm

| Item  | Description                     | Weight |
|-------|---------------------------------|--------|
| MT130 | Small, High Security Bolt Lock  | 22.0kg |
| MT140 | Medium, High Security Bolt Lock | 25.0kg |
| MT150 | Large, High Security Bolt Lock  | 31.0kg |

## Laptop Holder - Removable, Swing Arm

#### Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MT range of dispensing trolleys
- Fitted to the rear right hand upright of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle, fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports
- \*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**

Paint - Grey



| Item     | Description                          | Weight |
|----------|--------------------------------------|--------|
| MT/LHR/S | Laptop Holder - Removable, Swing Arm | 3.7kg  |

## Laptop Holder - Non-Removable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MT range of dispensing trolleys
- Fitted to the rear right hand upright of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle fitted with a steel mounting bracket, fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports
- \*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**


Paint - Grey

| Item      | Description                              | Weight |
|-----------|--|--------|
| MT/LHNR/S | Laptop Holder - Non-Removable, Swing Arm | 5.0kg  |

fixed using either nylock or wing nuts

Laptop supported in a cradle fitted with a steel mounting bracket that can be

**Drug & Medicine Dispensing** 





#### **Single Door**

## High Security Bolt LockElectronic Push Button LockDescriptionWeightPTS/HSB/PA24PTS/EPB/PA24Patient Administration, 24 Compartments51.4kg

## **Pharmacy Trolleys - Patient Administration**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using individual compartments dedicated for each patient
- Two locking options...
  - · High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Supplied with 4 / 8 trays subdivided into equal compartments...
  - 24 Compartment 4 trays subdivided into 6 equal compartments, each compartment - 190mm x 132mm x 220mm (5.5 Litres). Label holders provided for each compartment
  - 32 Compartment 8 trays subdivided into 4 equal compartments, each compartment - 192mm x 143mm x 220mm (6.1 Litres). Label holders provided for each compartment
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional  $\rm MT/1/6F$

#### Dimensions (w x d x h)

#### **Standard Colours**

• Work Surface - Spearmint

· Paint - Grev

- Single Door 627 x 579 x 1165mm
- Double Door 820 x 579 x 1165mm

#### **Options / Accessories**

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm



PTS/EPB/PA24



#### **Double Door**

| High Security Bolt Lock | Electronic Push Button Lock | Description                             | Weight |
|-------------------------|-----------------------------|---|--------|
| PTD/HSB/PA32            | PTD/EPB/PA32                | Patient Administration, 32 Compartments | 68.8kg |



Trays tilt to assist access to contents



PTD/HSB/PA32



| Item        | Description                              | Weight |
|-------------|--|--------|
| PTS/HSB/DR6 | High Security Bolt Lock, Six Drawers     | 53.0kg |
| PTS/EPB/DR6 | Electronic Push Button Lock, Six Drawers | 53.0kg |

## **Pharmacy Trolleys - Drawers**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using drawers
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- 3 adjustable door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

- Overall Dimensions 511 x 542 x 1072mm
- Drawer 335 x 295 x 113mm (internal)

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint

#### **Options / Accessories**

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- SEP/DIS Drawer Inserts, single door trolley refer to page 384



PTS/EPB/DR6



| ecurity Bolt Lock Electro | onic Push Button Lock | Description                 | Weight |
|---------------------------|-----------------------|-----------------------------|--------|
| Door                      |                       |                             |        |
| SB/LA16 PTS/E             | PB/LA16               | Sixteen LA Trays            | 39.0kg |
| e Door                    |                       |                             |        |
| SB/A18C12 PTD/E           | PB/A18C12             | Eighteen A & Twelve C Trays | 57.0kg |
| SB/LP24 PTD/E             | PB/LP24               | Twenty Four LP Trays        | 57.0kg |
| SB/A18C12 PTD/E           |                       | б ў                         |        |

## **Pharmacy Trolleys - Trays**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines
   using individual patient trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Trays supported on adjustable shelves
- Tray Sizes (w x d x h)...
  - LA 86 x 268 x 57mm (aluminium)
  - A 85 x 181 x 54mm (plastic)
  - C 129 x 181 x 54mm (plastic)
  - LP 115 x 266 x 95mm (plastic)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- Adjustable door shelves...
  - Single Door 3 door shelves
  - Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional  $\ensuremath{\mathsf{MT}}\xspace{-1.5mu}$

#### Dimensions (w x d x h)

- Single Door 511 x 542 x 1072mm
- Double Door 711 x 568 x 1072mm

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm





PTS/HSB/BP6 (frames not included)

## Pharmacy Trolleys - Blister Pack (MDS)

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using the blister packed monitored dosage system (MDS)
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Frames supported on three adjustable slide out shelves (frames not supplied)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- Adjustable door shelves...
  - Single Door 3 door shelves
  - Double Door 6 door shelves
- Fold over MDS dispensing surface
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

#### **Options / Accessories**

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- Single Door 511 x 542 x 1072mm
  Double Door 711 x 568 x 1072mm

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint



Tab prevents the shelf from being removed

| High Security Bolt Lock | Electronic Push Button Lock | Description               | Weight |
|-------------------------|-----------------------------|---------------------------|--------|
| Single Door             |                             |                           |        |
| PTS/HSB/BP6             | PTS/EPB/BP6                 | Blister Pack, Six Frames  | 38.0kg |
| Double Door             |                             |                           |        |
| PTD/HSB/BP9             | PTD/EPB/BP9                 | Blister Pack, Nine Frames | 42.0kg |

## Pharmacy Trolleys - Nomad<sup>®</sup> Cassettes

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using the Nomad® cassette monitored dosage system
- Two locking options...
  - · High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Cassettes supported on three adjustable shelves (cassettes not supplied)
- · Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- · Six adjustable door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

• Overall Dimensions 711 x 568 x 1072mm



- · Paint Grey
- Work Surface Spearmint

#### **Options / Accessories**

High Security Bolt Lock

PTD/HSB/NC56

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm

**Electronic Push Button Lock** 

PTD/EPB/NC56



sales@bristolmaid.com +44 (0)1258 484455

PTD/EPB/NC56

Weight

48.0kg

(-REER

16

Description

Nomad<sup>®</sup>, Fifty Six Cassettes

## Pharmacy Trolleys - Biodose® & MultiMeds



| High Security Bolt Lock | Electronic Push Button Lock | Description                                | Weight |
|-------------------------|-----------------------------|--|--------|
| Single Door             |                             |  |        |
| PTS/HSB/BD46            | PTS/EPB/BD46                | Biodose <sup>®</sup> & MultiMeds, 46 Trays | 37.3kg |
| Double Door             |                             |  |        |
| PTD/HSB/BD69            | PTD/EPB/BD69                | Biodose <sup>®</sup> & MultiMeds, 69 Trays | 52.0kg |

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using the Biodose<sup>®</sup> & MultiMeds monitored dosage system
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Trays supported on fixed runners (trays not supplied)
- Upper work surface
- Push handles incorporate integral buffers to offer protection to both the trolley & fabric of the building
- · Adjustable door shelves...
  - Single Door 3 door shelves
  - Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors
- When not in use, trolley should be secured to the building see optional  $\ensuremath{\mathsf{MT}}\xspace{1/6}$

#### Dimensions (w x d x h)

- Single Door 581 x 568 x 1072mm
- Double Door 812 x 568 x 1072mm

#### **Standard Colours**

- Paint Grey
- Work Surface Spearmint

#### **Options / Accessories**

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm

sales@bristolmaid.com +44 (0)1258 484455



| High Security Bolt Lock | Electronic Push Button Lock | Description        | Weight |
|-------------------------|-----------------------------|--------------------|--------|
| Single Door             |                             |                    |        |
| PTS/HSB/MX40            | PTS/EPB/MX40                | Medinoxx, 40 Trays | 41.4kg |
| Double Door             |                             |                    |        |
| PTD/HSB/MX60            | PTD/EPB/MX60                | Medinoxx, 60 Trays | 55.2kg |

## **Pharmacy Trolleys - Medinoxx**

#### **Features**

- Trolleys suitable for the storage & dispensing of drugs & medicines using the Medinoxx monitored dosage system (compatible with 7 x 4 medi-trays)
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Programmable four-digit code
- Trays supported on fixed runners (trays not supplied)
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the building
- · Adjustable door shelves
  - Single Door 3 door shelves
  - Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building see optional MT/1/6F

#### Dimensions (w x d x h)

- Single Door 581 x 568 x 1072mm
- Double Door 812 x 568 x 1072mm

### **Standard Colours**

- Paint Grey
- Work Surface Spearmint

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm



PTS/EPB/MX40

## Laptop Holder - Removable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the PT range of single & double door pharmacy trolleys
- Fitted to an upright bolted to the rear left-hand corner of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle, fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

#### Standard Colours

• Paint - Grey



| Item     | Description          | Weight |
|----------|----------------------|--------|
| PT/LHR/S | Removable, Swing Arm | 4.7kg  |

## Laptop Holder - Non-Removable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the PT range of single & double door pharmacy trolleys
- Fitted to an upright bolted to the rear left-hand corner of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle fitted with a steel mounting bracket, fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg

**Standard Colours** 

· Paint - Grey

• Design allows easy access to laptop ports



Laptop supported in a cradle fitted with a steel mounting bracket that can be fixed using either nylock or wing nuts

| ço         |
|------------|
| S          |
| Medicine   |
| <u>0</u>   |
| ne         |
|            |
| Dispensing |
| ĕ          |
| <u></u>    |
| J.         |
|            |
|            |


Drug

| Item      | Description              | Weight |
|-----------|--------------------------|--------|
| PT/LHNR/S | Non-Removable, Swing Arm | 5.0kg  |

## **Fitting Your Security Clamp**



Securely fasten one end of the clamp to the wall



The 'hook shaped' end is trapped by the door when closed & locked



sales@bristolmaid.com

## **Security Clamp**

#### **Features**

- · Allows the trolley to be secured to the 'fabric of the building'
- · Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door/lid of the trolley when closed & locked
- Recommended wall fitting height approx 1000mm



**Important Note - Security Clamps** 

It is a requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area & broken into, but also encourages good nursing practice, giving the trolley a "permanent home". "Safe & Secure Handling of Medicines - Royal Pharmaceutical Society,

GB-March 2005"

| Item    | Description    | Cable Length | Weight |
|---------|----------------|--------------|--------|
| MT/1/6F | Security Clamp | 950mm        | 0.75kg |



| Item  | Description                       | Weight |
|-------|-----------------------------------|--------|
| SA005 | Distribution Trolley - 15 Drawers | 49.0kg |

Distribution Trolley - 18 Drawers

SA025

| rawers | 56.0kg |  |
|--------|--------|--|
|        |        |  |
|        |        |  |

#### www.bristolmaid.com

## **Distribution Trolleys - Lockable Drawers**

#### **Features**

- Trolleys suitable for the distribution of drugs & medicines from pharmacy to the ward & also with self-administration / patient own drug schemes
- Provides secure drug delivery using drawers which when locked are securely fastened to the trolley
- Drawers are fitted with a CAM lock (supplied with two keys)
- Max load per drawer 3kg
- Supplied as standard with locks to differ, allowing each ward to have its own key & be responsible for removing the drawer from the trolley (remaining key should be kept at the pharmacy)
- Master keys available (additional cost)
- Upper work surface
- · Push handles at either end
- · Document pocket fitted to end of trolley
- Low level corner bumpers offer protection to both the trolley & fabric of the building
- 125mm pressed steel castors. 2 off braking

#### Dimensions (w x d x h)

- SA005 1026 x 523 x 1002mm
- SA025 1180 x 523 x 1002mm
- Drawer 145 x 423 x 230mm (Internal)

#### **Standard Colours**

- Paint Grev
- · Work Surface Spearmint

#### **Options / Accessories**

· Master Keys - contact us for details



# **Drug & Medicine Storage**

- Drug & Medicine Cabinets Shelves & Door Shelves 146
  - Drug & Medicine Cabinets Deep Shelves 148
  - Drug & Medicine Cabinets Blister Pack (MDS) 150
    - Controlled Drug Cabinets 152
- Auditable Access Control for Medicines Management 154
  - Drug & Medicine Cupboards 156
    - Options 158
- Patients / Residents Own Medication Cabinets Self-Administration 160
- Patients / Residents Own Medication Cabinets Blister Pack (MDS) 163
  - Tablet Counter & Storage Boxes 164
  - Pharmacy & Vaccine Refrigerators 165



## Drug & Medicine Cabinets - Shelves & Door Shelves

#### **Features**

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Adjustable & removable shelves (190mm deep) & door shelves (75mm deep)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems
- Compatible with our range of CD cabinets for the storage of Controlled Drugs (see page 152)

#### **Options / Accessories**

- DGL Digital Push Button Lock (Not available on independent locking cabinets)
- SL Suited Locks (contact us for more information)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/800/ST Sloping Top (800mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

#### **Standard Colours**

• Paint - White

| Right Hand Opening |          | Left Hand Opening |             | * Transformer casing adds 105mm to over |                        |                        | overall height |
|--------------------|----------|-------------------|-------------|---|------------------------|------------------------|----------------|
| Light*             | No Light | Light*            | No Light    | Overall Dimensions (w x d x h)          | Shelves / Door Shelves | Compatible CD Cabinets | Weight         |
| PC536              | PC536/LX | PC536/LH          | PC536/LH/LX | 500 x 300 x 600mm                       | 3/3                    | CD005, CD010           | 29.9kg         |
| PC539              | PC539/LX | PC539/LH          | PC539/LH/LX | 500 x 300 x 900mm                       | 4 / 4                  | CD005, CD010, CD015    | 34.0kg         |
| PC636              | PC636/LX | PC636/LH          | PC636/LH/LX | 600 x 300 x 600mm                       | 3/3                    | CD005, CD010           | 33.9kg         |
| PC639              | PC639/LX | PC639/LH          | PC639/LH/LX | 600 x 300 x 900mm                       | 4 / 4                  | CD005, CD010, CD015    | 38.0kg         |

#### Single Door Cabinets

+44 (0)1258 484455

sales@bristolmaid.com



PC836/LX

#### \*Transformer casing adds 105mm to overall height \*\*CD010 not compatible with PC836/IND & PC836/IND/LX **Double Door Cabinets** Single Lock (Secures Both Doors) Each Door Fitted With A Lock \*\*\*CD015 not compatible with PC839/IND & PC839/IND/LX Overall Dimensions (w x d x h) Shelves / Door Shelves Compatible CD Cabinets Weight PC836/LX PC836 PC836/IND PC836/IND/LX 800 x 300 x 600mm 6/5 CD005, CD010\*\* 40.9kg PC839 PC839/LX PC839/IND PC839/IND/LX 800 x 300 x 900mm 8/8 CD005, CD010, CD015\*\*\* 45.9kg PC1036 PC1036/LX PC1036/IND PC1036/IND/LX 1000 x 300 x 600mm 6/6 CD005, CD010 48.9kg PC1039 PC1039/LX PC1039/IND PC1039/IND/LX 1000 x 300 x 900mm 8/8 CD005, CD010, CD015 53.9kg

+44 (0)1258 484455







sales@bristolmaid.com

#### **Single Door Cabinets**

## **Drug & Medicine Cabinets - Deep Shelves**

#### **Features**

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Adjustable & removable deep shelves (250mm deep)
- Key profile unique to Hospital Metalcraft Ltd. offering 250,000 differs
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- · Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems
- · Compatible with our range of CD cabinets for the storage of Controlled Drugs (see page 152)

#### **Options / Accessories**

- DGL Digital Push Button Lock (Not available on independent locking cabinets)
- SL Suited Locks (contact us for more information)
- PC/300/ST Sloping Top (300mm Wide Cabinets)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/800/ST Sloping Top (800mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

#### Standard Colours

· Paint - White

| Right Ha | Right Hand Opening |             | Left Hand Opening |                                | * Transformer casing adds 105mm to |                        |                | o overall height |
|----------|--------------------|-------------|-------------------|--------------------------------|------------------------------------|------------------------|----------------|------------------|
| Light*   | No Light           | Light*      | No Light          | Overall Dimensions (w x d x h) | Shelves                            | Compatible CD Cabinets | Nomad Capacity | Weight           |
| PC336/FS | PC336/FS/LX        | PC336/FS/LH | PC336/FS/LH/LX    | 300 x 300 x 600mm              | 3                                  | CD005                  | 13 x Cassettes | 16.9kg           |
| PC536/FS | PC536/FS/LX        | PC536/FS/LH | PC536/FS/LH/LX    | 500 x 300 x 600mm              | 3                                  | CD005, CD010           | 24 x Cassettes | 29.9kg           |
| PC636/FS | PC636/FS/LX        | PC636/FS/LH | PC636/FS/LH/LX    | 600 x 300 x 600mm              | 3                                  | CD005, CD010           | 32 x Cassettes | 33.9kg           |
| PC539/FS | PC539/FS/LX        | PC539/FS/LH | PC539/FS/LH/LX    | 500 x 300 x 900mm              | 4                                  | CD005, CD010, CD015    | 35 x Cassettes | 34.9kg           |
| PC639/FS | PC639/FS/LX        | PC639/FS/LH | PC639/FS/LH/LX    | 600 x 300 x 900mm              | 4                                  | CD005, CD010, CD015    | 44 x Cassettes | 38.9kg           |



#### **Double Door Cabinets**

| Doub   | le Door Cabin       | ets             |                       | _                              |         | *Transformer<br>**CD010 not compatible with | casing adds 105mm to or<br>PC836/FS/IND & PC836 |            |
|--------|---------------------|-----------------|-----------------------|--------------------------------|---------|---|---|------------|
| Single | e Lock (Secures Doo | ors) Each Doc   | or Fitted With A Lock |                                |         | *** CD015 not compatible with               | n PC839/FS/IND & PC839                          | /FS/IND/LX |
| Light* | No Light            | Light*          | No Light              | Overall Dimensions (w x d x h) | Shelves | Compatible CD Cabinets                      | Nomad Capacity                                  | Weight     |
| PC836  | S/FS PC836/FS/L     | X PC836/FS/IND  | PC836/FS/IND/LX**     | 800 x 300 x 600mm              | 6       | CD005, CD010**                              | -   | 40.9kg     |
| PC839  | /FS PC839/FS/L      | X PC839/FS/IND  | PC839/FS/IND/LX***    | 800 x 300 x 900mm              | 8       | CD005, CD010, CD015***                      | -   | 45.9kg     |
| PC103  | 86/FS PC1036/FS/    | LX PC1036/FS/IN | D PC1036/FS/IND/LX    | 1000 x 300 x 600mm             | 6       | CD005, CD010                                | 58 Cassettes                                    | 48.9kg     |
| PC103  | 39/FS PC1039/FS/    | LX PC1039/FS/IN | D PC1039/FS/IND/LX    | 1000 x 300 x 900mm             | 8       | CD005, CD010, CD015                         | 80 Cassettes                                    | 53.9kg     |

sales@bristolmaid.com



PC536/MDS (Frames not included)

## **Drug & Medicine Cabinets - Blister Pack (MDS)**

#### **Features**

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Hooks at rear of cabinet support frames for the popular blister pack Monitored Dosage System (MDS)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/ EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems

#### **Standard Colours**

· Paint - White

| Distant.                             |              | 1.0011       | 10              | * Transform                    | er casing adds 105mm t | o overall height |
|--------------------------------------|--------------|--------------|-----------------|--------------------------------|------------------------|------------------|
| Right Hand Opening Left Hand Opening |              |              |                 |                                |                        |                  |
| Light*                               | No Light     | Light*       | No Light        | Overall Dimensions (w x d x h) | MDS Capacity           | Weight           |
| PC536/MDS                            | PC536/MDS/LX | PC536/LH/MDS | PC536/LH/MDS/LX | 500 x 300 x 600mm              | 2 x Frames             | 26.9kg           |
| PC539/MDS                            | PC539/MDS/LX | PC539/LH/MDS | PC539/LH/MDS/LX | 500 x 300 x 900mm              | 4 x Frames             | 31.4kg           |
| PC639/MDS                            | PC639/MDS/LX | PC639/LH/MDS | PC639/LH/MDS/LX | 600 x 300 x 900mm              | 6 x Frames             | 35.2kg           |

#### **Single Door Cabinets**

### **Options / Accessories**

Double Door Cabinets
Single Lock (Secures Both Doors)

• DGL - Digital Push Button Lock (Not available on independent locking cabinets)

PC1039/MDS PC1039/MDS/LX PC1039/IND/MDS PC1039/IND/MDS/LX 1000 x 300 x 900mm

- SL Suited Locks (contact us for more information)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)



52.1kg

8 x Frames



# Features Exceeds the statutory requires of Drugs (Safe Custody) Reg CD140 & CD145) BS 3621 anti-pick & drill, key identification card) Single-point locking Key profile unique to Hospita Compatible with CLIQ®, integ cylinders & keys, combining Can be supplied with a microsystem by others (CD905), ow with volt-free contact to allow -240y ~ 50Hz -60Hz 0.3A

# **Controlled Drug Cabinets - Single Point Locking**

- Exceeds the statutory requirement for the storage of controlled drugs "The Misuse of Drugs (Safe Custody) Regulations 1973" (with the exception of CD040, CD045, CD140 & CD145)
- BS 3621 anti-pick & drill, key retaining six pin cylinder lock (supplied with 2 keys & identification card)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs per profile
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Can be supplied with a micro-switch to allow connection to remote lighting/nurse call system by others (CD905), or supplied with a 12v DC LED lighting/warning system with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer casing adds 105mm to the overall height
- 1.5mm (16 gauge) mild steel, seamless welded construction
- Mounting holes provided in back (fixings not included)
- Adjustable & removable shelves
- Security of models CD005, CD010 & CD015 can be enhanced by concealing inside one of our 'Drug & Medicine Cabinets'
- Modular design allows for easy interchange with HTM 63 storage systems
- Supplied right hand hinging as standard, left hand available on request

### **Standard Colours**

• Paint - White

### **Options / Accessories**

- CD905 Micro-switch to connect to remote light box / alarm
- CD910 / CD911 / CD913 12v DC LED Lighting/Warning System (for compatibility see tables below)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

R/H Door L/H Door Overall Dimensions (w x d x h) Shelves Lighting / Warning System (Optional) Weight CD005 210 x 270 x 300mm CD910 -1 8.0kg CD010 CD110 CD910 335 x 270 x 300mm 1 9.0kg CD015 CD115 335 x 270 x 550mm 2 CD910 17.0kg CD020 CD120 500 x 300 x 550mm 2 CD910 20.0kg CD025 CD125 500 x 300 x 850mm 3 CD910 37.0kg

CD125







CD135

| R/H Door | L/H Door | Overall Dimensions (w x d x h) | Shelves | Lighting / Warning System (Optional) | Weight |
|----------|----------|--------------------------------|---------|--------------------------------------|--------|
| CD030    | CD130    | 500 x 450 x 850mm              | 3       | CD911                                | 42.0kg |
| CD035    | CD135    | 500 x 450 x 1250mm             | 4       | CD911                                | 59.0kg |
| CD045    | CD145    | 500 x 465 x 1900mm             | 6       | CD913                                | 91.0kg |



sales@bristolmaid.com

# Auditable Access Control for Medicines Management



### What Is CLIQ<sup>®</sup>?

CLIQ<sup>®</sup> is the very latest in user-friendly, web-based security management. Part of the Protec<sup>2</sup> CLIQ<sup>®</sup> from Abloy range, CLIQ<sup>®</sup> is designed for remote sites, enabling the commissioning, programming and full administration of security management centrally via web-based software.

### How Does It Work?

CLIQ<sup>®</sup> integrates electronics into mechanical cylinders and keys, providing a security solution that can be managed remotely. In short, this works like an electronic handshake: when the key is inserted it engages with the cylinder, sending a random number back to the key. The key replies with an encrypted response which is decrypted by the cylinder and if valid the cylinder then allows access. It is important to remember that the mechanical combination must match as well as the digital combination. CLIQ<sup>®</sup> allows a single key to be used with multiple cabinets or trolleys.

### **System Features**

- Secure Provides dual security through the patented Protec<sup>2</sup> CLIQ<sup>®</sup> from Abloy mechanical disc cylinders and keys along with DES encryption, a recognised electronic standard
- · Easily Programmable Keys can be added quickly and programmed to specific time and date schedules
- · Fully Auditable- Audit trails can be accessed from both the cylinder and key
- · Reduces Keys One key can be programmed to open multiple cabinets or trolleys
- The CLIQ® environment works in an Amazon cloud hosted server
- · maintained and supported by Abloy allowing remote user access via a secure
- · connection from anywhere in the world
- Cost Effective Low installation & maintenance costs along with a long lifecycle (requires no wires or system updates, power for the cylinder is provided by the key)
- Easily Retrofitable CLIQ<sup>®</sup> can be used with a variety of locking solutions including euro lock profile cases, padlocks, camlocks & furniture locks

### **System Benefits**

- · Save valuable nursing time which can be used to deliver patient care
- · Reduces delay in administering medication
- · Improves drug & medicine security
- · No more unlocked drug & medicine cabinets
- · Eliminates the need to replace locks as a result of lost keys
- · One key can be programmed to open multiple cabinets or trolleys
- The key is made personal to that member of staff, resulting in increased ownership and a reduction in the number of keys lost
- Easy to apply remote management, saves money and has the ability to cover several hospitals
- Time and date based access rights, with a full audit trail allow for full visibility of when visits were made
- New keys can be added to the system quickly and easily, removal of lost ones is hassle-free
- · Reduction in time spent looking for the person holding the key

# **Process Flow**

A JA



### **During The Shift**

### Operation

The key will only open cabinets where access permissions have been granted. Each time the key is used in a cabinet fitted with an PROTEC<sup>2</sup> cylinder an audit trail is created on both the key and cylinder



### **Start of Shift**

### Key Validation & Authorisation

Recurring validation forces users to place their key into the Wall PD. User authorisations are uploaded to the key



### End of Shift

### Automated Data Transfer

Key is placed into the Wall PD. Previous key transactions are automatically downloaded from the key and stored in the audit trail archive

# **System Components**



### **CLIQ®** Kev

- · Provides dual security through the patented PROTEC<sup>2</sup> mechanical disc cylinders and keys along with DES encryption, a recognised electronic standard
- Thin, light but incredibly strong key

### **PROTEC<sup>2</sup>** Lock Barrel

A bump proof mechanism

as well as padlocks

· Exceptionally smooth key entry and operation




· Can be master keyed and keyed alike with cylinders, cam and furniture locks



### Wall Programming Device

- · Updates keys remotely (linking the server and users keys)
- Offers high level data encryption
- Allows the use of recurring validation, forcing users to validate their keys at the start of their shift whilst automatically downloading audit trails which are relayed to the server

### Traka21



- · Keys are accessed using PIN number technology
- Holds up to 21 keys or keysets



### **Questions & Answers** Answers to common questions about CLIQ®

### O How secure is it?

- A CLIQ<sup>®</sup> uses a DES encryption (Data Encryption Standard) so it is just as safe as online banking
- Q Can I use this with an existing Ablov master key system?
- A Yes: CLIQ<sup>®</sup> can be integrated with an existing PROTEC<sup>2</sup> master key system
- Q What is the benefit?
- A CLIQ<sup>®</sup> will give you some of the benefits of access control, such as adding, deleting and amending access rights, take audit trails from keys and cylinders as well as setting the key to expire on a specific date at a specific time. Also, lost or stolen key can be easily be deleted electronically

### Q Is the system expensive to install?

A The system has a low installation & maintenance cost along with a long life cycle. CLIQ<sup>®</sup> cylinders are designed to fit into Euro Profile mechanical lock cases. No wiring is required to the cabinet. The wall programming device requires an Ethernet connection and power through a 12-24V adapter or alternatively power over Ethernet

### Q Who maintains the system?

- A System maintenance is normally assigned to a member of the hospital staff who will be responsible for programming/issuing of keys along with interrogation of audit trails
- Q How many audit trails can be collected from a single key?
- A The programmable user key stores the last 1.800 transactions and the last 20 foreign transactions
- Q How large is the cylinder memory?
- A 200 key groups; each key group holding 65,000 individual unique key combinations. This provides up to 13 million unique possible key combinations in one cylinder and can retain a black list of 3.000 keys
- Q Can I take an audit trail from a cylinder?
- A Yes, the cylinder stores an audit trail of the last 1,000 key insertions

### Q What happens if a key is lost?

A We recommend that time schedules are set to each key, so if the key is not returned it will expire. Keys can be removed from the cylinder with the programming key or directly from the system. If subsequently found keys can be added back and re-issued. In the event that an audit trail is required for a lost key. lock cylinders can be interrogated using a programming key (each cylinder stores an audit trail of the last 1,000 key insertions)

### Q How long do the batteries last?

- A The battery is located in the key, not the cylinder. Under normal usage a battery will operate up to 20.000 cvcles or 2 vears
- Q How is the information sent?
- A Via the internet to a wall mounted PD
- Q Can the information be taken off the Wall mounted PD by anyone?
- A No, there is no information stored on the PD, it is stored on a server until the correct key is inserted and only then is the information securely transferred
- Q Is it possible to access the system remotely?
- A The CLIQ<sup>®</sup> environment works in an Amazon cloud hosted server environment ISO27001 Accredited. maintained and supported by Abloy allowing remote user access via a secure connection from anywhere in the world.
- Q Where is the server?
  - A Amazon Republic Of Ireland
- Q Can I host this myself?
- A Yes, but there would be additional costs
- Q What if there is no internet connection on a particular site?
- A mobile PD is available which will work on the GSM network
- Q How can I better control the issuing of keys?
- A Through the use of Traka21 (key safe) keys can be automatically assigned to specific users. Keys are issued using either the staff members own security card or via pin number technology.
- Q Is it possible to integrate the CLIQ<sup>®</sup> key with other systems?
- A Unique RFID tag identifiers can be inserted into the key allowing it to be used with other MIFARE card based RFID systems
  - Self-Administration Schemes
  - Medicines Management
  - Patient Records



# **Drug & Medicine Cupboards - Single Door**

### **Features**

- Independently tested to BS2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Single door cupboards available as either 870 / 1870mm high
- Available with either...
  - Plastic Trays
  - Shelves
  - Shelves & Dividers
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ<sup>®</sup>, integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- Worksurface supplied with 870mm cupboards (removable to assist with positioning under existing worktops)
- 0.9mm (20 Gauge) mild steel welded construction
- Modular design allows for easy interchange with HTM 63 storage systems

### Dimensions (w x d x h)

- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal

### **Standard Colours**

· Paint - White

### **Options / Accessories**

- SL Suited Locks (contact us for more information)
- DGL Digital Push Button Lock
- MC/JK/870 Joining Kit & Full Width Worktop (required when connecting 870mm cupboards)
- MC/JK/1870 Joining Kit (required when connecting 1870mm cupboards)
- MC/ST Sloping Top (1870mm cupboards)
- MC/DIV Set of 16 Dividers

| R/H Hinged    | L/H Hinged    | Overall Dimensions (w x d x h) | Plastic Trays / Shelves | Weight |
|---------------|---------------|--------------------------------|-------------------------|--------|
| MC/R/870/3S2D | MC/L/870/3S2D | 500 x 536 x 870mm              | Plastic Trays           | 45.6kg |
| MC/R/870/FS   | MC/L/870/FS   | 500 x 536 x 870mm              | Shelves                 | 37.0kg |



MC/R/1870/7S5D MC/L/1870/7S5D

MC/L/1870/FS

MC/L/1870/FSD

MC/R/1870/FS

MC/R/1870/FSD





### **Plastic Trays**

- Interchangeable with our popular Caretray range
  - 870mm Cupboards 3 x 100mm, 2 x 150mm trays
  - 1870mm Cupboards 7 x 100mm, 5 x 150mm trays
  - Supplied with 7 x 100mm & 5 x 150mm trays
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray



- Adjustable at 50mm increments
  - 870mm Cupboards 3 x shelves
  - 1870mm Cupboards 8 x shelves
- Suitable for the storage of MDS frames (if required)
  - 870mm Cupboards 6 x frames
  - 1870mm Cupboards 16 x frames



### **Shelves & Dividers**

- Adjustable at 50mm increments
- Removable dividers 100mm high adjustable at 25mm increments
  - 870mm Cupboards N/A
  - 1870mm Cupboards 8 x shelves & 16 x dividers

sales@bristolmaid.com

**Plastic Trays** 

Shelves & Dividers

Shelves

Overall Dimensions (w x d x h)

500 x 536 x 1870mm

500 x 536 x 1870mm

500 x 536 x 1870mm

88.6kg

69.0kg

79.0kg

Plastic Trays / Shelves / Shelves & Dividers

# **Sloping Top**

### Features

- Compatible with 300/500/600/800 & 1000mm wide PC cabinets & 500mm wide MC cupboards
- · Prevents items being stored on the top of the cabinet & assists cleaning
- Increases the overall height of the cabinet by 193mm
- Void within the top provides a space to locate the transformer for the lighting/warning system
- · Access to the transformer provided via a removable plate

### **Standard Colours**

· Paint - White


# Warning Light System

### **Features**

- Suitable for use with Bristol Maid CD Controlled Drug cabinets
- · Provides a visual indication when the cabinet is open
- LED illuminates...
  - Green when the door is closed
  - Red when the door is open
  - Flashes red when open for 2 1/2 minutes
- Operates using a door activated micro-switch connected to a 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height



| Item     | Description  |
|----------|--|
| CD Range | Compatible   |
| CD910    | 12v DC LED Lighting/Warning System - Refer To Page 152 For Compatibility |
| CD911    | 12v DC LED Lighting/Warning System - Refer To Page 153 For Compatibility |
| CD913    | 12v DC LED Lighting/Warning System - Refer To Page 153 For Compatibility |



| Item       | Description             | Weight |
|------------|-------------------------|--------|
| ntem       | Description             | weight |
| PC/300/ST  | 300mm Wide PC Cabinets  | 1.3kg  |
| PC/500/ST  | 500mm Wide PC Cabinets  | 1.5kg  |
| PC/600/ST  | 600mm Wide PC Cabinets  | 1.6kg  |
| PC/800/ST  | 800mm Wide PC Cabinets  | 1.7kg  |
| PC/1000/ST | 1000mm Wide PC Cabinets | 1.8kg  |
| MC/ST      | 500mm Wide MC Cupboards | 2.9kg  |

# **Joining Kits - MC Range of Cabinets**

### **Features**

- · Allows two units to be joined together
- 870mm joining kit supplied with a full width worktop (2 units)



| Item       | Description                      | Compatibility     |
|------------|----------------------------------|-------------------|
| MC/JK/870  | Joining Kit & Full Width Worktop | 870mm High Units  |
| MC/JK/1870 | Joining Kit                      | 1870mm High Units |

# **Remote Light Box**

### **Features**

- Connects to all Bristol Maid PC, MS & CD cabinets fitted with a light)
- Provides a remote visual indication when the cabinet is open
- · Models available for surface or flush mounting
- Units fitted without a light require a separate micro-switch & a 12v DC supply



| Item    | Description                        | <b>Overall Dimensions</b> |
|---------|------------------------------------|---------------------------|
| PC/RLBS | Remote Light Box - Surface Mounted | 100 x 50 x 150mm          |
| PC/RLBF | Remote Light Plate - Flush Mounted | 100 x 50 x 150mm          |

# **Digital Push Button Lock**

### **Features**

- Requires a five digit code to operate (set by customer)
- Cabinets fitted with this lock will not conform to BS2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises



Digital Push Button Lock

# Dividers - MC Range of Cabinets

### **Features**

DGL

- Spare / replacement dividers for MC cupboards supplied with shelves & dividers
- 100mm high dividers locate into slots on the shelves, adjustable at 25mm increments
- · Compatible with 870mm & 1870mm high units



### Shelf not supplied

| Item   | Descrip  |
|--------|----------|
| MC/DIV | Set of 1 |

| Description        |  |
|--------------------|--|
| Set of 16 Dividers |  |

# **Micro Switch**

### **Features**

• Micro-switch to allow connection to



CD905

### Item Description

CD905Micro-switch (CD Range Compatible)PC905Micro-switch (PC Range Compatible)

# **Buzzer Alarm**

### Features

- · Replacement circuit board with buzzer
- Compatible for PC, MS & CD cabinets fitted with a light unit
- Buzzer pulses on/off after the door has been open for 2 ½ minutes
- Factory fit only to cabinets fitted with lights

| Item | Description  |
|------|--------------|
| PCB  | Buzzer Alarm |









SAC/313/S/VP

### Features

- Suitable for use with...
  - Drugs retained in original packs
  - Individual MDS packs with dividers & name cards
- Nomad cassettes
- Roll type blister packs
- Available as side or bottom hinged fitted with either a CAM or electronic push button lock
  - CAM Supplied with two keys & uses 40 different master suites (for security we recommend that each ward/wing operates a different series, master keys available at additional cost)
  - Electronic Push Button Battery operated
- Available with a solid door or with an inset perspex visibility panel
- 0.9mm (20 gauge) mild steel, spot welded construction
- Mounting holes in back/bottom (fixings not included)
- Range of mounting plates (optional)
- Shelves not included (optional)
- Supplied with a self-adhesive name card holder (length 60mm)

### **Standard Colours**

• Paint - White

### **Electronic Push Button Lock**

- Uses a 10-button keypad
- Easy programming & code change
- Four operating codes with differing levels of authority - Master / Sub-Master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code



| CAM Lock |            | Electric Push Button Lock |              |                          |                                |                     |        |  |  |
|----------|------------|---------------------------|--------------|--------------------------|--------------------------------|---------------------|--------|--|--|
|          | Solid Door | Perspex Visibility Panel  | Solid Door   | Perspex Visibility Panel | Overall Dimensions (w x d x h) | Shelf Compatibility | Weight |  |  |
|          | SAC/213/S  | SAC/213/S/VP              | SAC/213/S/CL | SAC/213/S/CL/VP          | 210 x 155 x 315mm              | SAC/FS/2            | 2.6kg  |  |  |
|          | SAC/313/S  | SAC/313/S/VP              | SAC/313/S/CL | SAC/313/S/CL/VP          | 310 x 155 x 315mm              | SAC/FS/3            | 3.4kg  |  |  |
|          | SAC/413/S  | SAC/413/S/VP              | SAC/413/S/CL | SAC/413/S/CL/VP          | 410 x 155 x 315mm              | SAC/FS/4            | 4.1kg  |  |  |

### **Side Hinged**

| sales@bristolmaid.com | +44 (0)1258 484455 |
|-----------------------|--------------------|

4455 www.bristolmaid.com

### Options

- SAC/FS/2 Shelf, 200mm Wide
- SAC/FS/3 Shelf, 300mm Wide
- SAC/FS/4 Shelf, 400mm Wide
- SAC/WP Mounting Plate, Wall
- SAC/LT Mounting Plate, Locker Top



### **Bottom Hinged**

CAM Lock **Electric Push Button Lock** SAC/FS/2 2.6kg SAC/213/B SAC/213/B/VP SAC/213/B/CL SAC/213/B/CL/VP 210 x 155 x 315mm SAC/313/B SAC/313/B/VP SAC/313/B/CL SAC/313/B/CL/VP 310 x 155 x 315mm SAC/FS/3 3.4kg SAC/413/B SAC/413/B/VP SAC/413/B/CL SAC/413/B/CL/VP 410 x 155 x 315mm SAC/FS/4 4.1kg SAC/213/B/VP

sales@bristolmaid.com

+44 (0)1258 484455

www.bristolmaid.com

# **Shelves**

### **Features**

- Flat shelves (100mm deep)
- Removable to assist cleaning, allow the storage of larger items

### **Standard Colours**

Paint - White

# **Mounting Plates**

### **Features**


- Compatible with all SAC cabinets using existing holes (fixings not provided)
- SAC/WP Cabinet secured to wall plate using a small thumbscrew. Allows the cabinet to be repositioned
- SAC/LT Mounts to upper surface of traditional bedside cabinets

### Standard Colours



| Item     | Description                                 | Weight |
|----------|---|--------|
| SAC/FS/2 | Shelf - 200mm Wide Cabinets (SAC/213 Range) | 0.2kg  |
| SAC/FS/3 | Shelf - 300mm Wide Cabinets (SAC/313 Range) | 0.2kg  |
| SAC/FS/4 | Shelf - 400mm Wide Cabinets (SAC/413 Range) | 0.3kg  |

| Item   | Description | Overall Dimensions (w x d x h) | Weight |
|--------|-------------|--------------------------------|--------|
| SAC/WP | Wall        | 190 x 12 x 162mm               | 0.5kg  |
| SAC/LT | Locker Top  | 190 x 150 x 150mm              | 0.9kg  |



SAC/FS/4

sales@bristolmaid.com

```
www.bristolmaid.com
```

# Patients / Residents Own Medication Cabinets - Monitored Dosage, MDS Frames



### **Features**

- Suitable for the storage of MDS frames with dividers & name cards (max 100mm depth)
- · Models available to hold one or two frames (not included)
  - Small One frame
  - Large Two frames
- Available with either a CAM or electronic push button lock
  - CAM Supplied with two keys, master keys available at additional cost
  - · Electronic Push Button Battery operated
- Mounting holes provided in back (fixings not included)
- Supplied with an adjustable shelf (SA045 & SA045/CL only)
- Constructed from melamine faced chipboard offering an extremely strong & durable, yet aesthetically pleasing appearance
- MDS frames not supplied

### **Standard Colours**

• Body / Door - Beech Woodgrain

### **Electronic Push Button Lock**

- Uses a 10-button keypad
- Easy programming & code change
- Four operating codes with differing levels of authority - Master / Sub-Master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

02

34

56

28

00

SA040 - frames not included

| CAM Lock | Digital Push Button Lock | Overall Dimensions (w x d x h) | MDS Capacity | Weight |
|----------|--------------------------|--------------------------------|--------------|--------|
| SA040    | SA040/CL                 | 325 x 245 x 420mm              | 1 x Frame    | 7.5kg  |
| SA045    | SA045/CL                 | 325 x 245 x 775mm              | 2 x Frames   | 12.1kg |

+44 (0)1258 484455

# **Tablet Counter**

### **Features**

- Triangular stainless steel tablet counter
- Etched legend to interpret rows (max 22) into number of tablets

# **Storage Boxes - Lightweight**


### **Features**

- Lightweight aluminium riveted construction
- Fitted with handle
- Single-point locking using a CAM lock (supplied with 2 keys)
- Cannot be wall mounted

### **Standard Colours**

• Paint - White







| Item  | Description    | Dimensions (w x d x h) |
|-------|----------------|------------------------|
| PE/30 | Tablet Counter | 205 x 195 x 10mm       |

+44 (0)1258 484455

sales@bristolmaid.com

www.bristolmaid.com





Powered by

5RLDF0219

| CAM Lock          | Digital Push Button | Overall Dimensions (w x d x h) | Capacity   | Description          | Shelves | Weight |
|-------------------|---------------------|--------------------------------|------------|----------------------|---------|--------|
| Solid Door        |                     |                                |            |                      |         |        |
| 5RLDF0119         | 5RLDF0119DIGILOCK   | 450 x 510 x 538mm              | 36 Litres  | Bench/Wall Mounting* | 2       | 30.0kg |
| 5RLDF0219         | 5RLDF0219DIGILOCK   | 450 x 510 x 740mm              | 66 Litres  | Bench/Wall Mounting* | 3       | 34.0kg |
| 5RLDF0519         | 5RLDF0519DIGILOCK   | 600 x 600 x 835mm              | 150 Litres | Under Bench          | 3       | 51.0kg |
| Glass Door        |                     |                                |            |                      |         |        |
| 5RLDG0519         | 5RLDG0519DIGILOCK   | 600 x 600 x 835mm              | 150 Litres | Under Bench          | 3       | 54.0kg |
| *Integral Wall Mo | ounting Brackets    |                                |            |                      |         |        |
|                   |                     |                                |            | 44 (0) 4050 404455   |         |        |

### Features

- Compliant with the Department of Health Green Book
- Recommendations in the General Pharmaceutical Council Standards
- Royal College of Veterinary Surgeons Guidelines
- Designed specifically for the safe storage of vaccines and pharmaceuticals
- Available as bench top, under bench, wall fitting or freestanding
- Temperature range +2°C to +8°C
- Second wandering temperature probe (can be placed anywhere within the refrigerator, ideal for pharmacists needing to comply with MHRA regulations)
- IntelliCold<sup>®</sup> touch screen controller
  - Automatic 24/7 temperature & fridge operation data recording (micro SD card required)
  - Min/max temperature control, display & recording
  - Displays temperature from both probes
  - High/low temperature alarm
  - Door open alarm
  - Mains failure alarm (requires batteries)
  - Alarm battery back-up
- Internal LED light
- · Access port & remote alarm contacts
- · Fan assisted cooling
- · Off cycle auto defrost
- Low GWP environmentally friendly refrigerant
- · Self-closing doors with door open alarm
- Available with either a CAM or digital push button lock
  - CAM Supplied with two keys
  - Digital Push Button Mains operated (uses a 4-digit pin)
- · Adjustable shelves
- Micro SD card (optional)
- UKAS traceable calibration service (optional)
- · Repeat alarm (optional)

### **Standard Colours**

- Exterior White / Blue
- · Interior White



# Examination / Treatment Couches & Chairs

- Examination & Treatment Couches Variable Height 168
  - Examination & Treatment Couches Fixed Height 172
    - Couch Steps 173
      - Leg Rests 174
- Examination & Treatment Couches Variable Height, Bariatric 176
  - Examination & Treatment Chairs Variable Height 180
  - Examination & Treatment Chairs Variable Height, Bariatric 181
    - Bariatric Examination & Treatment Couches 192
      - Phlebotomy Chairs Variable Height 182
    - Phlebotomy Chairs Variable Height, Bariatric 183
      - Phlebotomy Chairs Fixed Height 184
        - Options 185



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Hydraulic in Aquamarine - EC1205/VAQ



|            |   |        | 1170mm  | 615mi  |
|------------|---|--------|---------|--------|
| Item       | Description                                   | Weight | TT7 Omm | 015111 |
| EC1205/VBB | Hydraulic                                     | 68.8kg |         |        |
| EC1210/VBB | Electric, Hand Switch                         | 68.8kg |         |        |
| EC1215/VBB | Electric, Foot Switch                         | 68.8kg |         |        |
| EC1220/VBB | Electric with Electric Backrest - Hand Switch | 71.1kg |         |        |
| EC1225/VBB | Electric with Electric Backrest - Foot Switch | 71.1kg |         |        |
|            |   |        |         |        |

# Examination & Treatment Couch Two Section

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

### Dimensions (w x d x h)

• 630 x 1790 x 490 - 950mm

### **Standard Colours**

- · Paint White
- Vinyl Bristol Blue

### **Options**

- E02 Paper Roll Holder
- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E11 Lithotomy Stirrups (Pair)
- E12 Leg Supports (Pair)

630mm

# **Examination & Treatment Couch Two Section, Mobile**



### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- · Height adjustment from 490 950mm
- · Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- · Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinvl
- · Swing down side support rails fitted to either side
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

### Dimensions (w x d x h)

• 630 x 1790 x 490 - 950mm

### Standard Colours

- · Paint White
- · Vinyl Bristol Blue

### Options

- · E02 Paper Roll Holder
- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E11 Lithotomy Stirrups (Pair)
- E12 Leg Supports (Pair)

| Foot Switch | Alternative | E Colours avai | ilable at ex<br>gi.e. Electr | tra cost. Changic, Hand Switc                | ge suffix 'V   | EC1235/V                                      | oduct code |             | HADI      | EIGHT                      |
|-------------|-------------|----------------|------------------------------|--|----------------|---|------------|-------------|-----------|----------------------------|
| Hand Switch | VBB         | VAQ            | VFR                          | VGH  | VSD            | VOB   | VLI        | VCH         | VCB       | VPE                        |
| 1170mr      |             | 615mm          |                              | Item<br>EC1230/V                             | D              | escription                                    |            |             |           | Weight<br>75.0kg           |
|             |             |                | 630mm                        | EC1235/V<br>EC1235/V<br>EC1240/V<br>EC1245/V | BB EI<br>BB EI | ectric, Hanc<br>ectric, Foot<br>ectric with E | Switch     | krest - Har | nd Switch | 75.0kg<br>75.0kg<br>75.3kg |

EC1250/VBB

Electric with Electric Backrest - Foot Switch

77.3kg



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric, Foot Switch in Onyx Black - EC1315/VOB



| Item       | Description                                   | Weight | 690mm | 470mm |
|------------|---|--------|-------|-------|
| EC1305/VBB | Hydraulic                                     | 70.2kg |       |       |
| EC1310/VBB | Electric, Hand Switch                         | 70.2kg |       |       |
| EC1315/VBB | Electric, Foot Switch                         | 70.2kg |       |       |
| EC1320/VBB | Electric with Electric Backrest - Hand Switch | 72.5kg |       |       |
| EC1325/VBB | Electric with Electric Backrest - Foot Switch | 72.5kg |       |       |

# Examination & Treatment Couch Three Section

### Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch

Planned

Preventative

Maintenance Contact us for

details

- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

### Dimensions (w x d x h)

• 630 x 1790 x 490 - 950mm

### **Standard Colours**

- · Paint White
- Vinyl Bristol Blue

### Options

- E02 Paper Roll Holder
- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole

VCB

VPE

630mm

Hand Switch

615mm

# Examination & Treatment Couch Three Section, Mobile

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- · Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- 125mm swivel castors featuring a central locking & steering facility
- · Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless & contoured upholstery using flame retardant & antimicrobial vinyl
- · Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

### Dimensions (w x d x h)

• 630 x 1790 x 490 - 950mm

### **Standard Colours**

- · Paint White
- Vinyl Bristol Blue

### **Options**

- E02 Paper Roll Holder
- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E05 Breather Hole & Plug
- · E06 Battery Back Up



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Coastal Blue - EC1345/VCB



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC1330/VBB | Hydraulic                                     | 76.4kg |
| EC1335/VBB | Electric, Hand Switch                         | 76.4kg |
| EC1340/VBB | Electric, Foot Switch                         | 76.4kg |
| EC1345/VBB | Electric with Electric Backrest - Hand Switch | 78.7kg |
| EC1350/VBB | Electric with Electric Backrest - Foot Switch | 78.7kg |
|            |   |        |

Hand Switch

Foot Switch

630mm



### **Features**

- Fixed height approx. 786mm
- Adjustable angle backrest 0 to +90 degrees using a scissor jack
- · Adjustable feet allow levelling on uneven floors
- Seamless upholstery using flame retardant anti-microbial vinyl
- Paper roll holder fitted to foot end (suitable for rolls max 510mm wide x 200mm diameter)
- Supplied boxed for self-assembly
- Maximum capacity 125kg

### Dimensions (w x d x h)

• 610 x 1860 x 795mm

### **Standard Colours**

- Paint White
- · Legs Chrome
- Vinyl Bristol Blue





sales@bristolmaid.com

Chairs లర **Examination / Treatment Couches** 

CE

PETAPDAN

www.bristolmaid.com

# Couch Steps One & Two Tread

### Features

- Non-slip rubber feet & tread
- Tread dimension 460 x 250mm
- Tread Heights
  - One Tread 180mm
  - Two Tread 185 & 365mm
- One Step available with a handle for additional safety
- Handle height max. 1005mm
- Maximum capacity 125kg

### **Standard Colours**

• Paint - White





| ltem  | Description                    | Overall Dimensions (w x d x h) | Weight |
|-------|--------------------------------|--------------------------------|--------|
| CS005 | Couch Step - One Tread         | 519 x 284 x 228mm              | 4.4kg  |
| CS010 | Couch Step - One Tread, Handle | 519 x 284 x 1005mm             | 6.0kg  |
| CS015 | Couch Step - Two Tread         | 519 x 599 x 405mm              | 9.6kg  |

www.bristolmaid.com





# Leg Rest - Adjustable

### **Features**

- Independently adjustable legs
- Height & angle adjustment from 355 510mm
- Seamless upholstery using flame retardant anti-microbial vinyl
- Maximum capacity 25kg

### Dimensions (w x d x h)

- Overall Dimensions 460 x 305 x 350 510mm
- Upholstered Pad 460 x 305mm

### **Standard Colours**

- Paint White
- Vinyl Bristol Blue



| Item      | Description           | Weight |
|-----------|-----------------------|--------|
| CS/28/VBB | Leg Rest - Adjustable | 3.5kg  |

VOB

VLI

VCH

VCB

VSD

VPE

VBB

VAQ

VFR

VGH

# Other products that you may be interested in...

See Pages 52 & 53

Caretray Trolleys - Mild Steel, Standard Level, Single & Double Column







Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Peach - EC2215/VPE



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC2205/VBB | Electric, Hand Switch                         | 79.0kg |
| EC2210/VBB | Electric, Foot Switch                         | 79.0kg |
| EC2215/VBB | Electric with Electric Backrest - Hand Switch | 81.3kg |
| EC2220/VBB | Electric with Electric Backrest - Foot Switch | 81.3kg |

# **Examination & Treatment Couch Two Section, Bariatric**

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- · Height adjustment from 520 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

### Dimensions (w x d x h)

• 850 x 1790 x 520 - 950mm

### Standard Colours

- · Paint White
- · Vinyl Bristol Blue

### **Options**

- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E18 Battery Back Up
- · E19 Paper Roll Holder
- E21 Breather Hole & Plug

### Planned Preventative Maintenance Contact us for details

1170mm

850mm

Hand Switch

615mm

VPF

Chairs

ංඊ

**Examination / Treatment Couches** 

# Examination & Treatment Couch Two Section, Mobile, Bariatric

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- · Height adjustment from 520 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- 230V  $\sim,$  50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinyl
- · Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

### Dimensions (w x d x h)

• 850 x 1790 x 520 - 950mm

### **Standard Colours**

- Paint White
- · Vinyl Bristol Blue

### Options

- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

|   |  | CE                   |
|---|--|----------------------|
|   |  |                      |
|   |  | HEIGHT<br>ADJUSTABLE |
|   | EC2225/VBB   |                      |
| · | Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Foot Switch in Bristol Blue - EC2240/VBB |                      |
|   |  |                      |



| 1170mm |
|--------|
|        |

| 615mm |       |            |   |        |
|-------|-------|------------|---|--------|
|       |       | ltem       | Description                                   | Weight |
|       | 8     | EC2225/VBB | Electric, Hand Switch                         | 85.2kg |
|       | 850mm | EC2230/VBB | Electric, Foot Switch                         | 85.2kg |
|       | З     | EC2235/VBB | Electric with Electric Backrest - Hand Switch | 87.5kg |
|       |       | EC2240/VBB | Electric with Electric Backrest - Foot Switch | 87.5kg |

Foot Switch

Hand Switch



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC2305/VBB | Electric, Hand Switch                         | 81.8kg |
| EC2310/VBB | Electric, Foot Switch                         | 81.8kg |
| EC2315/VBB | Electric with Electric Backrest - Hand Switch | 84.1kg |
| EC2320/VBB | Electric with Electric Backrest - Foot Switch | 84.1kg |
|            |   |        |

# **Examination & Treatment Couch Three Section, Bariatric**

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- · Height adjustment from 520 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinvl
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

### Dimensions (w x d x h)

850 x 1790 x 520 - 950mm

### Standard Colours

- Paint White
- Vinvl Bristol Blue

### Options

- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- · E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

Planned Preventative Maintenance Contact us for details

600mm

470mm

850mm

615mm

Chairs ංඊ **Examination / Treatment Couches** 

# **Examination & Treatment Couch Three Section, Mobile, Bariatric**

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- · Height adjustment from 520 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections (optional electric backrest)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- 230V  $\sim,$  50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinyl
- · Swing down side support rails fitted to either side
- · Frames prepared to accept optional accessories
- Maximum capacity 325kg

### Dimensions (w x d x h)

• 850 x 1790 x 520 - 950mm

### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

### Options

- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug



- Foot Switch

| Hand | Switch |
|------|--------|
|------|--------|

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric, Foot Switch in Flame Red - EC2330/VFR

www.bristolmaid.com



| 470mm | 615mm |       |            |   |        |
|-------|-------|-------|------------|---|--------|
|       |       |       | Item       | Description                                   | Weight |
|       |       | 8     | EC2325/VBB | Electric, Hand Switch                         | 88.0kg |
|       |       | 850mm | EC2330/VBB | Electric, Foot Switch                         | 88.0kg |
|       |       | З     | EC2335/VBB | Electric with Electric Backrest - Hand Switch | 90.3kg |
|       |       |       | EC2340/VBB | Electric with Electric Backrest - Foot Switch | 90.3kg |

EC2340/VBB



HEIGHT

690mm



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric - Foot Switch in Peach - EC4315/VPE

| VBB     | VAQ        | VFR           | VGH         | VSD          | VOB        | VLI    | VCH  | VCB    | VPE |       |
|---------|------------|---------------|-------------|--------------|------------|--------|------|--------|-----|-------|
| Item    |            | Description   | 1           |              |            | Weight |      |        |     |       |
| EC4305/ | <b>VBB</b> | Hydraulic     |             |              |            | 66.5kg | 360m | m 470m | m   | 820mm |
| EC4310/ | VBB        | Electric - Ha | nd Switch   |              |            | 66.5kg |      |        |     |       |
| EC4315/ | <b>VBB</b> | Electric - Fo | ot Switch   |              |            | 66.5kg |      |        |     |       |
| EC4320/ | VBB        | Electric with | Electric Ba | ackrest - Ha | and Switch | 68.8kg |      |        |     |       |
| EC4325/ | VBB        | Electric with | Electric Ba | ackrest - Fo | oot Switch | 68.8kg |      |        |     |       |
|         |            |               |             |              |            |        |      |        |     |       |

# Examination & Treatment Chair Drop End, Three Section

### Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

### Dimensions (w x d x h)

• 630 x 1710 x 490 - 950mm

### **Standard Colours**

- Paint White
- Vinyl Bristol Blue

### Options

- E02 Paper Roll Holder
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E09 Bowl

www.bristolmaid.com

630mm

- E10 Bowl Support
- E13 Phlebotomy Armrests (Pair)
- E14 Side Support Loops (Pair)



Foot Switch



Hand Switch

# Examination & Treatment Chair Drop End, Three Section, Bariatric

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- · Height adjustment from 520 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- · Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant antimicrobial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

### Dimensions (w x d x h)

• 850 x 1710 x 520 - 950mm

### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

### Options

- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E09 Bowl
- E10 Bowl Support
- E13 Phlebotomy Armrests (Pair)
- · E14 Side Support Loops (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- · E21 Breather Hole & Plug



Foot Switch

Hand Switch







| 820mm |       |            |   |        |
|-------|-------|------------|---|--------|
|       |       | Item       | Description                                   | Weight |
|       | 85    | EC4330/VBB | Electric - Hand Switch                        | 76.7kg |
|       | 850mm | EC4335/VBB | Electric - Foot Switch                        | 76.7kg |
|       | 2     | EC4340/VBB | Electric with Electric Backrest - Hand Switch | 79.0kg |
|       |       | EC4345/VBB | Electric with Electric Backrest - Foot Switch | 79.0kg |
|       |       |            |   |        |

sales@bristolmaid.com +44

470mm

HEIGHT



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric - Foot Switch in Graphite - EC5315/VGH



| Item       | Description  | Weight   | 360mm  |
|------------|--|--|--|
| EC5305/VBB | Hydraulic  | 67.0kg   |  |
| EC5310/VBB | Electric - Hand Switch                               | 67.0kg   |  |
| EC5315/VBB | Electric - Foot Switch                               | 67.0kg   |  |
| EC5320/VBB | Electric with Electric Backrest - Hand Switch        | 69.3kg   |  |
| EC5325/VBB | Electric with Electric Backrest - Foot Switch        | 69.3kg   |  |
|            | EC5305/VBB<br>EC5310/VBB<br>EC5315/VBB<br>EC5320/VBB | EC5305/VBBHydraulicEC5310/VBBElectric - Hand SwitchEC5315/VBBElectric - Foot SwitchEC5320/VBBElectric with Electric Backrest - Hand Switch | EC5305/VBBHydraulic67.0kgEC5310/VBBElectric - Hand Switch67.0kgEC5315/VBBElectric - Foot Switch67.0kgEC5320/VBBElectric with Electric Backrest - Hand Switch69.3kg |

# Phlebotomy Chair Variable Height

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (selfassembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

## Dimensions (w x d x h)

• 770 x 1710 x 490 - 950mm

### **Standard Colours**

- Paint White
- Vinyl Bristol Blue

### Options

E06 - Battery Back Up

470mm

820mm

630mm

# Phlebotomy Chair Variable Height, Bariatric

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- Height adjustment from 490 950mm
- Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- · Clearance beneath frame allows access for mobile hoists
- 230V  $\sim,$  50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (selfassembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- · Frames prepared to accept optional accessories
- Maximum capacity 325kg

### Dimensions (w x d x h)

• 990 x 1710 x 490 - 950mm

### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

### Options


· E18 - Battery Back Up



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Onyx Black - EC5340/VOB



360mm 470mm 820mm



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC5330/VBB | Electric - Hand Switch                        | 77.2kg |
| EC5335/VBB | Electric - Foot Switch                        | 77.2kg |
| EC5340/VBB | Electric with Electric Backrest - Hand Switch | 79.5kg |
| EC5345/VBB | Electric with Electric Backrest - Foot Switch | 79.5kg |

Planned Preventative

Maintenance

Contact us for

details



| Item       | Description  | Seat Height | Weight |
|------------|--------------|-------------|--------|
| EC5000/VBB | Fixed Height | 485mm       | 15.0kg |

# **Phlebotomy Chair Fixed Height**

### **Features**

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Certified for use up to 32 stone (203kg)
- Stable four leg design
- · Black finish tubular frame
- · Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrest
- · Armrest can be positioned on either side of the chair (additional armrest available if required)
- Maximum capacity 203kg

### Dimensions (w x d x h)

• 620 x 580 x 800mm

### **Standard Colours**

- · Paint Black
- · Vinyl Bristol Blue

### **Options**

• E15 - Phlebotomy Armrest (Single)

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Fixed Height Phlebotomy Chair in Aquamarine - EC5000/VAQ



www.bristolmaid.com

Chairs

Examination / Treatment Couches &



# Paper Roll Holder

· Fitted to the head section

E02

 Suitable for rolls max 510mm wide x 200mm diameter



# Side Support Rails (Pair)

### **E03**

- Swing down-side support rails
- Can be added to either one or both sides of the couch to provide additional safety for patients



# Cushioned Side Support Rails (Pair)

### E04 - E04/U

- Vinyl upholstered swing down-side support rails
- Can be added to either one or both sides of the couch to provide additional safety for patients
- Upgrade option for factory fitted side support rails (E04/U)



# Breather Hole & Plug E05

- Breather hole located in backrest section to improve patient comfort when in the prone position
- Supplied with a plug / insert



# **Battery Back Up**

### **E06**

- Allows the couch to be lowered in a controlled manner in the event of a power failure
- Provides enough charge for 1-3 cycles



# Height Adjustable Armrests (Pair)

### **E07**

 Height adjustable easy clean moulded polyurethane armrests





## **Infusion Pole**

### E08

- Height adjustable two hook infusion poleCannot be fitted in conjunction with side
- support rails E03, E04, E04/U



# **Bowl**

 360mm diameter 304 Grade stainless steel bowl



# **Bowl Support**

### E10

- Support to hold 360mm diameter 304 Grade stainless steel bowl
- Order bowl separately



# Lithotomy Stirrups (Pair) E11

• Traditional lithotomy stirrups



# Leg Supports (Pair) E12

• Fully adjustable easy clean moulded polyurethane leg supports



# Phlebotomy **Armrests (Pair)** E13

· Fully adjustable easy clean moulded polyurethane armrests



# Side Support Loops (Pair)

 Side support loops to assist mount / dismount & provide patient security



# Phlebotomy Armrest (Single)

# E15

- Fully adjustable easy clean moulded polyurethane armrest
- Compatible with EC5000



# **Battery Back Up**

# E18 (Bariatric Couch)

- Allows the couch to be lowered in a controlled manner in the event of a power failure
- Provides enough charge for 1-3 cycles



# Paper Roll Holder

# E19 (Bariatric Couch)

- Fitted to the head section
- Suitable for rolls max 510mm wide x 200mm diameter



# Breather Hole & Plug

# E21 (Bariatric Couch)

- Breather hole located in backrest section to improve patient comfort when in the prone position
- Supplied with a plug / insert



# **General Purpose Trolleys**

- Stainless Steel Trolleys Heavy Duty 190
  - Stainless Steel Light Duty 191
  - Polymer Shallow & Deep Shelves 192
    - Mild Steel Laminate Shelves 193





# General Purpose Trolley Stainless Steel, Heavy Duty

#### Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with three sides up & one down (shelf size 920 x 535mm)
- Available with two or three shelves, shelf heights
  - Two shelf 280 & 860mm
  - Three shelf 220, 540 & 860mm
- Maximum capacity per shelf 75kg
- 125mm pressed steel castors

# Dimensions (w x d x h)

- No Ring Buffers 1000 x 570 x 1020mm
- Ring Buffers 1030 x 600 x 1020mm

## **Standard Colours**

Stainless Steel

#### Options

• C008 - 125mm Castors, Heavy Duty Pressed Steel, 2 x Braking, 2 x Non-Braking (set of 4) - Factory fitted

| Item          | Description   | Weight |
|---------------|---|--------|
| WT/37/2       | Stainless Steel - Heavy Duty, Two Shelves                 | 28.0kg |
| WT/37/3       | Stainless Steel - Heavy Duty, Three Shelves               | 37.5kg |
| WT/37/2/ZP010 | Stainless Steel - Heavy Duty, Two Shelves, Ring Buffers   | 28.2kg |
| WT/37/3/ZP010 | Stainless Steel - Heavy Duty, Three Shelves, Ring Buffers | 37.7kg |

# General Purpose Trolley Stainless Steel, Light Duty

### Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed liquid retaining sunken shelves (shelf size 825 x 430mm)
- Available with two or three shelves, shelf heights
  - Two shelf 250 & 870mm
  - Three shelf 250, 560 & 870mm
- Maximum capacity per shelf 50kg
- 100mm pressed steel castors

# Dimensions (w x d x h)

- No Ring Buffers 1100 x 540 x 1040mm
- Ring Buffers 1100 x 580 x 1040mm

# **Standard Colours**

Stainless Steel



www.bristolmaid.com

| Item             | Description   | Weight |
|------------------|---|--------|
| WT/33/SS/2       | Stainless Steel - Light Duty, Two Shelves                 | 16.0kg |
| WT/33/SS/3       | Stainless Steel - Light Duty, Three Shelves               | 22.0kg |
| WT/33/SS/2/ZP005 | Stainless Steel - Light Duty, Two Shelves, Ring Buffers   | 16.2kg |
| WT/33/SS/3/ZP005 | Stainless Steel - Light Duty, Three Shelves, Ring Buffers | 22.2kg |

# **General Purpose Trollev Polymer, Shallow Shelves**

#### **Features**

- Easy to clean corrosion resistant high-density polyethylene shelves (not suitable for use with automated washers)
- 11mm raised edge retains spills & prevents items from sliding off during transport
- Microban® product protection inhibits growth of bacteria & aids infection control
- Ergonomic handle & recessed shelves provide additional legroom allowing you to push the cart without hitting your shins & feet
- Chrome plated posts offers shelf adjustment in 25mm increments
- Available with two / three shelves in three sizes
- Small 620 x 415mm (internal) at heights 185, 515 & 845mm (lower shelves) adiustable)
  - Easy grip handle & utility trav
  - Maximum capacity per shelf 68kg / per cart 136kg
  - 100mm pressed steel castors
- Medium 693 x 541mm (internal) at heights 185, 515 & 845mm (lower shelves) adjustable)
  - Handle has built in recessed support channel to manage odd shaped items
  - Maximum capacity per shelf 68kg / per cart 136kg
  - 100mm pressed steel castors
- Large 833 x 652mm (internal) at heights & 218, 548 & 878mm (lower shelves) adiustable)
  - Handle has built in slot to hold scanners & labellers
  - Maximum capacity per shelf 90kg / per cart 180kg
  - 125mm pressed steel castors
- Self-assembly (supplied boxed)

# **General Purpose Trolley Polymer, Deep Shelves**

# Features

- Easy to clean corrosion resistant high-density polyethylene shelves (not suitable for use with automated washers)
- 70mm deep shelves retains spills & prevents items from sliding off during transport
- Microban® product protection inhibits growth of bacteria & aids infection control
- · Chrome plated posts offers shelf adjustment in 25mm increments
- Available with two / three shelves in two sizes
  - Small 730 x 469mm (internal) at heights 140, 470 & 800mm (lower shelves) adiustable)
  - Large 877 x 616mm (internal) at heights 140, 470 & 800mm (lower shelves adjustable)
- Max capacity per shelf 68kg / per cart 136kg
- 100mm pressed steel castors
- Self-assembly (supplied boxed)



# Standard Colours

Standard Colours

Small - 850 x 547 x 1067mm

Large - 997 x 691 x 1067mm

· Uprights - Chrome Shelves - Blue

8BC2030-2DMB

8BC2030-3DMB

8BC2636-2DMB

8BC2636-3DMB

Item

- · Uprights Chrome
- · Shelves Blue

| ltem              | Shelves     | Weight |
|-------------------|-------------|--------|
| Small - 805 x 470 | x 902mm     |        |
| 8MY1627-24BU      | 2           | 12.9kg |
| 8MY1627-34BU      | 3           | 15.4kg |
| Medium - 877 x 6  | 600 x 902mi | m      |
| 8MY2030-24BU      | 2           | 14.9kg |
| 8MY2030-34BU      | 3           | 18.8kg |
| Large - 1023 x 70 | 04 x 937mm  | 1      |
| 8MY2636-25BU      | 2           | 18.1kg |
| 8MY2636-35BU      | 3           | 23.1kg |





sales@bristolmaid.com +44 (0)1258 484455 Shelves

2

3

2

3

www.bristolmaid.com

Weight

13.0ka

17.0kg

15.9kg

21.1kg

# General Purpose Trolley Mild Steel, Laminate Shelves

### Features

- Laminate shelves
- Available with two or three shelves in two sizes
- Shelf heights
  - Two shelf 210 & 845mm
  - Three shelf 210, 525 & 845mm
- Shelf size
  - Small 705 x 445mm
  - Large 855 x 445mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel castors

# **Standard Colours**

- Paint Grey
- Shelf Spearmint

# Options

• ZP005 - Ring Buffers





| Item         | Description                                | Weight |
|--------------|--|--------|
| Small - 972  | x 545 x 1035mm                             |        |
| GP005        | Mild Steel - Two Laminate Shelves, Small   | 18.0kg |
| GP015        | Mild Steel - Three Laminate Shelves, Small | 23.0kg |
| Large - 1122 | 2 x 545 x 1035mm                           |        |
| GP025        | Mild Steel - Two Laminate Shelves, Large   | 20.0kg |
| GP035        | Mild Steel - Three Laminate Shelves, Large | 26.5kg |



# **HTM71 Storage Solutions**

- ScanCell<sup>®</sup> Carts 196
- Distribution Trolleys 201
- Open Transfer / Transport Trolleys 203
  - HTM71 Racking 204
- Trays, Baskets, Dividers & Accessories 206





# ScanCell<sup>®</sup> Carts - Open Fronts



#### **Features**

- · Mobile open front storage carts for use with 600 x 400mm trays & baskets
- Available as single, double & treble columns in a range of heights
- Lightweight easyclean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides. & back
- Fitted with angled / horizontal side panels to suit 400 or 600mm insertions
  - Type A 600mm Insertions
  - Type B 400mm Insertions
- 125mm polymer castors, 2 off braking

# Standard Colours

Trays, Baskets, Dividers & Accessories Page 206 - 209

· Body - Grey

# Options

- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
- 4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)
- 4PR1A Push Rail, Type A (Cart Heights 1503mm, 1803mm)
- 4PR1B Push Rail, Type B (Cart Heights 1503mm, 1803mm)



(trays & baskets not included)

(trays & baskets not included)

# ScanCell<sup>®</sup> Carts - Roller Shutters



### **Features**

- Mobile storage carts fitted with roller shutters for use with 600 x 400mm trays & baskets
- · Available as single, double & treble columns in a range of heights
- Lightweight easyclean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back, polypropylene roller shutter
- · Fitted with angled / horizontal side panels to suit 400 or 600mm insertions
  - Type A 600mm Insertions
  - Type B 400mm Insertions
- 125mm polymer castors, 2 off braking

Trays, Baskets, Dividers & Accessories Page 206 - 209

# Body - Grey Options

Standard Colours

- 411P1A Litility Pail Type A (Cart Heights 005mm 1055
- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
   4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)









# **ScanCell® Carts - Doors**



Trays, Baskets, Dividers

& Accessories

Page 206 - 209

#### **Features**

- Mobile storage carts fitted with doors for use with 600 x 400mm trays & baskets
- Available as single, double & treble columns in a range of heights
- Lightweight easy clean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back, post form / plastic finish doors
- Fitted with angled / horizontal side panels to suit 400 or 600mm insertions
- Type A 600mm Insertion Single door for each column
  - Single Column R/H Hinged
  - Double Column R/H & L/H Hinged
  - Treble Column 1 x R/H & 2 x L/H Hinged
- Type B 400mm Insertion Double doors for each column
- 125mm polymer castors, 2 off braking

# **Standard Colours**

Body - Grey

## Options

- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
- 4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)
- 4PR1A Push Rail, Type A (Cart Heights 1503mm, 1803mm)
- 4PR1B Push Rail, Type B (Cart Heights 1503mm, 1803mm)



HTM71 Storage Solutions

sales@bristolmaid.com +44 (0)1258 484455

www.bristolmaid.com

# ScanCell<sup>®</sup> Carts - Vario Drawers



#### **Features**

- Mobile storage carts fitted with vario drawers & supplied with 600 x 400mm trays
- Available as single, double & treble columns in a range of heights
- Lightweight easy clean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back
- Type B 400mm Insertions
- 125mm polymer castors, 2 off braking

# Dimensions (w x d x h)

## **Standard Colours**

- Single Column 662 x 494 x Height
- Body / Drawers Grey
  Drawer Handle Blue
- Double Column 1324 x 494 x Height
- Treble Column 1986 x 494 x Height

# Options

- 4UR1B Utility Rail (Cart Height 905mm, 1055mm, 1205mm)
- 4PR1B Push Rail (Cart Height 1503mm)

| Item                    | 50mm Tray | 100mm Tray | 200mm Tray | 300mm Tray |
|-------------------------|-----------|------------|------------|------------|
| Height - 905mm          |           |            |            |            |
| 4B1VD09/A               | 1         | 2          | 1          |            |
| 4B1VD09/B               |           | 4          |            |            |
| 4B1VD09/C               | 2         | 3          |            |            |
| 4B2VD09/AA*             | 2         | 4          | 2          |            |
| 4B3VD09/AAA**           | 3         | 6          | 3          |            |
| Height - 1055mm         |           |            |            |            |
| 4B1VD10/D               | 1         | 3          | 1          |            |
| 4B1VD10/E               |           | 5          |            |            |
| 4B1VD10/F               | 2         | 4          |            |            |
| 4B2VD10/DD*             | 2         | 6          | 2          |            |
| 4B3VD10/DDD**           | 3         | 9          | 3          |            |
| Height - 1205mm         |           |            |            |            |
| 4B1VD12/G               | 1         | 4          | 1          |            |
| 4B1VD12/H               | 2         | 2          | 2          |            |
| Height - 1503mm         |           |            |            |            |
| 4B1VD15/I               | 1         | 6          | 1          |            |
| 4B1VD15/J               | 1         | 4          | 1          | 1          |
| *double column **treble | e column  |            |            |            |



**TYPE B - 400mm Insertion** 







#### sales@bristolmaid.com



# **ScanCell® Carts - Telescopic Drawers**

#### **Features**

- Mobile storage carts with telescopic drawers supplied with 600 x 400mm trays & baskets
- Available as single, double & treble columns in a range of heights
- Lightweight easy clean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back
- Type B 400mm Insertions
- Fitted with horizontal side panels to suit 400 insertions
- 125mm polymer castors, 2 off braking

# Dimensions (w x d x h)

- Single Column 662 x 494 x Height
- Double Column 1324 x 494 x Height

# **Standard Colours**

• Body / Drawers - Grey

# Options

- Trays, Baskets, Dividers & Accessories Page 206 - 209
- lowe
- 4UR1B Utility Rail (Cart Height 905mm, 1055mm, 1205mm)

| Item            | 50mm Tray | 100mm Tray | 200mm Basket |
|-----------------|-----------|------------|--------------|
| Height - 905mm  |           |            |              |
| 4B1TD09/A       | 1         | 2          | 1            |
| 4B1TD09/B       |           | 4          |              |
| 4B1TD09/C       | 2         | 3          |              |
| 4B2TD09/AA *    | 2         | 4          | 2            |
| Height - 1055mm |           |            |              |
| 4B1TD10/D       | 1         | 3          | 1            |
| 4B1TD10/E       |           | 5          |              |
| 4B1TD10/F       | 2         | 4          |              |
| 4B1TD10/G       | 10        |            |              |
| Height - 1205mm |           |            |              |
| 4B1TD12/H       | 1         | 4          | 1            |
| 4B1TD12/I       | 2         | 5          |              |
| 4B1TD12/J       | 6         | 3          |              |

\*double column

# **ScanCell® Carts - Compact**



#### **Features**

- · Mobile storage carts with telescopic drawers supplied with 300 x 400mm trays & baskets
- Lightweight easyclean design
- Injection moulded ABS top / chassis with aluminium corner profiles, laminate sides & back
- Type D 400mm Insertions
- 125mm polymer castors, 2 off braking

# Dimensions (w x d x h) Options

- 462 x 494 x 905mm
- 4UR1B Utility Rail (Type B)

#### Standard Colours

· Body / Drawers - Grev



4D1TD09/A shown with optional utility rail

| Item      | 50mm Tray | 100mm Tray | 200mm Basket |
|-----------|-----------|------------|--------------|
| 4D1TD09/A | 1         | 2          | 1            |
| 4D1TD09/B |           | 4          |              |

# **Distribution Trolleys - Compact**

# **Features**

- Mobile distribution trollevs for use with 300 x 400mm travs & baskets
- Available in mild & 304 grade stainless steel
  - · Mild Steel Welded & powder coated
  - · Stainless Steel Seamlessly welded for ease of cleaning
- Runners allow travs to be accessed from either side (4 stops required per trav / basket)
  - Type D 400mm insertion
- Trays & baskets can be supported on the upper & lower frame (optional solid grade laminate worktop)
- 75mm pressed steel castors

# Dimensions (w x d x h)

• 360 x 460 x 850mm

# Standard Colours

· Paint - Grey

# **Options**

- ZP010 Ring Buffers
- 4D1SLWT- Solid Grade Laminate Worktop







4D1DSTCM (trays & baskets not included)

| Item     | Description                             |
|----------|---|
| 4D1DSTCM | Compact, Single Column, Mild Steel      |
| 4D1DSTCS | Compact, Single Column, Stainless Steel |



# **Distribution Trolleys**

### **Features**

- Mobile distribution trolleys for use with 600 x 400mm trays & baskets
- · Available in mild & 304 grade stainless steel
  - Mild Steel Welded & powder coated
  - · Stainless Steel Seamlessly welded for ease of cleaning
- Runners allow trays to be accessed from either side (4 stops required per tray / basket)
  - Type A 600mm insertion
  - Type B 400mm insertion
- Trays & baskets can be supported on the upper & lower frame (optional solid grade laminate worktop)
- 75mm pressed steel castors

# **Standard Colours**

• Paint - Grey

# Options

- ZP010 Ring Buffers
- 4A1SLWT Solid Grade Laminate Worktop, Single Column (Type A)
- 4A2SLWT Solid Grade Laminate Worktop, Double Column (Type A)
- 4B1SLWT Solid Grade Laminate Worktop, Single Column (Type B)

#### TYPE A - 600mm Insertion

#### **TYPE B - 400mm Insertion**

Trays, Baskets, Dividers

& Accessories

Page 206 - 209



| ltem            | Description                            | Overall Dimensions (w x d x h) |
|-----------------|--|--------------------------------|
| Mild Steel      |  |                                |
| 4A1DSTM         | Single Column, Type A, Mild Steel      | 460 x 655 x 850mm              |
| 4A2DSTM         | Double Column, Type A, Mild Steel      | 895 x 655 x 850mm              |
| 4B1DSTM         | Single Column, Type B, Mild Steel      | 660 x 455 x 850mm              |
| Stainless Steel |  |                                |
| 4A1DSTS         | Single Column, Type A, Stainless Steel | 460 x 655 x 850mm              |
| 4A2DSTS         | Double Column, Type A, Stainless Steel | 895 x 655 x 850mm              |
| 4B1DSTS         | Single Column, Type B, Stainless Steel | 660 x 455 x 850mm              |
|                 |  |                                |

sales@bristolmaid.com

665 x 725 x 1555mm

1100 x 725 x 1555mm

1535 x 725 x 1635mm

665 x 725 x 1555mm

1100 x 725 x 1555mm

1535 x 725 x 1635mm

Solutions

# **Open Transfer / Transport Trolleys**

# Features

- Mobile transfer / transport trolleys for use with 600 x 400mm trays & baskets
- Available in mild & 304 grade stainless steel
  - Mild Steel Welded & powder coated
  - Stainless Steel Seamlessly welded for ease of cleaning
- Runners allow trays to be accessed from either side (4 stops required per tray / basket)
  - Type A 600mm insertion
- Push handles at either end
- · Ring buffers to protect the trolley / fabric of the building
- Castors
  - Single / Double Column 125mm pressed steel castors, 2 off braking
  - Treble Column 200mm heavy duty pressed steel castors, 2
     off fixed, 2 off swivel braking

# **Standard Colours**

· Paint - Grey

# Options

Item

Mild Steel 4A10TTM

4A2OTTM

4A3OTTM

4A2OTTS

4A3OTTS

Stainless Steel 4A10TTS

- 4A1OTTNC Nylon Cover (Single Column)
- 4A2OTTNC Nylon Cover (Double Column)

Description

Single Column, Type A, Mild Steel

Double Column, Type A, Mild Steel

Treble Column, Type A, Mild Steel

Single Column, Type A, Stainless Steel

Double Column, Type A, Stainless Steel

Treble Column, Type A, Stainless Steel

• 4A3OTTNC - Nylon Cover (Treble Column)













2HH-7532293 Wall Mounting Bracket



2HH-7532277 10mm Plastic Cover Cap

sales@bristolmaid.com

4HH-HTM71/S1/400 Runners, trays, baskets, dividers & angle brackets not included (order separately)

# Components

- 2HH-7532301-HH Base Unit, Type A, 600mm Insertion (x 1)
- 2HH-7532311 Extension Unit, Type A, 600mm Insertion (x 1)
- 2HH-7532356 Runners, Type A, 600mm Insertion (pair)
- 2HH-7532329 Angle Adaptor Plate, Type A, 600mm Insertion (pair)
- 2HH-7532304-HH Base Unit, Type B, 400mm Insertion (x 1)
- 2HH-7532314 Extension Unit, Type B, 400mm Insertion (x 1)
- 2HH-7532355 Runners, Type B, 400mm Insertion (pair)
- 2HH-7532328 Angle Adaptor Plate, Type B, 400mm Insertion (pair)
- 2HH-7532320 Adjustable Foot Set (x 4)
- 2HH-7532321 Adjustable Foot Set (x 2)
- 2HH-7532322 50mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532323 50mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)
- 2HH-7532324 80mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532325 80mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)

# HTM71 Racking

#### • Fully configurable freestanding modular racking for use with 600 x 400mm modular trays & baskets

- Available to suit 600 or 400mm insertions
  - Type A 600mm insertion
  - Type B 400mm insertion
- · Trays & baskets supported on adjustable moulded plastic/glass fibre runners
- Runners can be positioned horizontal or angled (approx. 25 degrees requires angle adaptor brackets)
- · Integral stops on the runners allow trays/baskets to be accessible in two positions
  - Straight (approx. 5 degrees)
  - Angled (approx. 25 degrees)
- · Runners offer pass through facility allowing access on either side
- · Anodised extruded aluminium uprights & cross members
- Bolted construction using high quality moulded fixings provides rigidity & allows the system to be altered easily
- · Available as either static or mobile units
  - · Static Adjustable feet allow levelling on uneven floors
  - Mobile 50/80mm pressed steel castors
- Optional wall mounting brackets allow static units to be fixed securely to the wall, (when units are loaded over 100kg
- Supplied for self-assembly

# Dimensions (w x d x h)

#### Static

Type A - 600mm Insertion

- Single Bay 494 x 642 x 1940mm
- Double Bay 937 x 642 x 1940mm

Type B - 400mm Insertion

- Single Bay 694 x 442 x 1940mm
- Double Bay 1337 x 442 x 1940mm

#### **Standard Colours**

- Uprights & Cross Members Anodised Aluminium
- Plastic Mouldings Green/Grey

#### **Static - Adjustable Feet**

Item Description Weight Type A - 600mm Insertion 4HH-HTM71/S1/600 Static, Single Bay, 600mm Insertion 10.3kg 4HH-HTM71/S2/600 Static, Double Bay, 600mm Insertion 16.1kg Type B - 400mm Insertion 4HH-HTM71/S1/400 Static, Single Bay, 400mm Insertion 10.2kg 4HH-HTM71/S2/400 Static, Double Bay, 400mm Insertion 16.3kg

- Mobile
  - Type A 600mm Insertion
    - Single Bay 494 x 642 x 1996mm
    - Double Bay 937 x 642 x 1996mm

#### Type B - 400mm Insertion

- Single Bay 694 x 442 x 1996mm
- Double Bay 1337 x 442 x 1996mm

#### Options

- 2HH-7532293 Wall Mounting Bracket
- 2HH-7532277 10mm Plastic Cover Cap

| +44 ( | (0)1258 484455 |
|-------|----------------|
|-------|----------------|



#### Mobile - 50mm Castors

| Item                     | Description                         | Weight |
|--------------------------|-------------------------------------|--------|
| Type A - 600mm Insertion |                                     |        |
| 4HH-HTM71/M51/600        | Mobile, Single Bay, 600mm Insertion | 11.1kg |
| 4HH-HTM71/M52/600        | Mobile, Double Bay, 600mm Insertion | 17.3kg |
| Type B - 400mm Insertion |                                     |        |
| 4HH-HTM71/M51/400        | Mobile, Single Bay, 400mm Insertion | 11.0kg |
| 4HH-HTM71/M52/400        | Mobile, Double Bay, 400mm Insertion | 17.5kg |



Runners, travs, baskets, dividers not included (order separately)

V LITTLE CONTRACTOR 

utation and the second

THE OWNER WATCHING THE PARTY OF THE PARTY OF

NUMBER OF STREET, STRE

manumunum

Runners, trays, baskets, dividers

not included (order separately)

TANK AND REAL PROPERTY AND REAL PROPERTY.

annin mark 

4HH-HTM71/M52/400

NUMBER OF STREET, STRE

......

111111

OTT DATA

TTLLIT

ITING

HUHB TZ

TETHON

....

1111111111111111111

# **Trays & Baskets** Type D - Size 300 x 400mm

Item







2HH-7014200

Basket - ABS, 200mm Dimensions (w x d x h): 300 x 400 x 200mm Material: ABS Grey

Description



2HH-7012100



2HH-7012200

Basket - ABS, 100mm Dimensions (w x d x h): 600 x 400 x 100mm Material: ABS Grey

Description

Tray - ABS, 50mm

Material: ABS Grev

Tray - ABS, 100mm

Material: ABS Grey

Basket - ABS, 200mm Dimensions (w x d x h): 600 x 400 x 200mm Material: ABS Grey



Tray/Lid - ABS, 15mm Dimensions (w x d x h): 600 x 400 x 15mm Material: Plastic Grey

2HH-K561039D


# **Dividers**

# Type D - Size 300 x 400mm

| Item         | Description  | Tray/Basket Compatibility  | Item        | Description  | Tray/Basket Compatibility  |
|--------------|--|----------------------------|-------------|--|----------------------------|
| 2HH-7810103  | Divider - PC<br>Dimensions (w x h): 300 x 50mm<br>Material: Polycarbonate, Transparent | 2HH-7014050                |             | Divider - ABS<br>Dimensions (w x h): 400 x 200mm<br>Material: ABS Grey   | 2HH-7014200<br>2HH-7012200 |
| 2HH-7810102  | Divider - PC<br>Dimensions (w x h): 400 x 50mm<br>Material: Polycarbonate, Transparent | 2HH-7014050<br>2HH-7011050 | 2HH-7811305 | Sliding Divider - SAN<br>Dimensions (w x h): 170 x 45mm<br>Material: Plastic, Transparent                            | 2HH-7014050<br>2HH-7011050 |
| 2HH-7811206N | Divider - ABS<br>Dimensions (w x h): 300 x 100mm<br>Material: ABS Grey                 | 2HH-7014100                | 2HH-241302  | Sliding Divider - SAN<br>(including cardholder)<br>Dimensions (w x h): 170 x 45mm<br>Material: Plastic, Transparent  | 2HH-7014050<br>2HH-7011050 |
| 2HH-7811205N | Divider - ABS<br>Dimensions (w x h): 400 x 100mm<br>Material: ABS Grey                 | 2HH-7014100<br>2HH-7012100 | 2HH-241301  | Sliding Divider - SAN<br>Dimensions (w x h): 170 x 100mm<br>Material: Plastic, Transparent                           | 2HH-7014100                |
|              | Divider - ABS<br>Dimensions (w x h): 300 x 200mm<br>Material: ABS Grey                 | 2HH-7014200                | 2HH-240820  | Sliding Divider - SAN<br>(including cardholder)<br>Dimensions (w x h): 170 x 100mm<br>Material: Plastic, Transparent | 2HH-7014100                |
| 2HH-7811306  |  |                            | 2HH-240821  |  |                            |

# **Dividers**


# Type A & B - Size - 600 x 400mm



# Accessories

| Item         | Description   | Tray/Basket Compatibility |
|--------------|---|---------------------------|
| 2HH-240501   | T/Clip  | All Trays & Baskets       |
| 2HH-7910204N | Safety Stop - Long<br>(Scan Modul Compatible)<br>Material: Plastic  | All Trays & Baskets       |
| 2HH-7910201N | Safety Stop - Short<br>(H+H System Compatible)<br>Material: Plastic | All Trays & Baskets       |
| 2HH-100072   | Label Holder<br>Credit Card Size: (85 x 54mm)<br>Material: Plastic  | All Trays & Baskets       |
|              | Label<br>Credit Card Size: (85 x 54mm)<br>Material: Plastic, White  | All Trays & Baskets       |

2HH-K100187



# Infusion / Pump Stands & Accessories

- Mobile Pump Stands 212
- Mobile Infusion Stands 216
- Clamping Infusion Stands 221
- Wall Mounted Infusion Holders 221



# Mobile Pump Stands - Polyurethane, Easy Clean

#### **Features**

- Range of stands designed specifically to hold blood warmers, pumps & syringe drivers
- · Five leg base with low centre of gravity to assist stability
- 32mm diameter 304 grade stainless steel upright & 25mm diameter pump support bows
- Polyurethane overmoulded base with smooth edges & contours provides...
  - Water-resistant finish, allowing for easy & effective cleaning (including steam cleaning)
  - · Shield for the castors
  - Impact protection to patients & staff, reducing the risk of skin abrasions & bruising
  - Protection to the fabric of the building
- Height adjustable upright using a handwheel (1,300 2,000mm)
- · Four U shaped hooks
- Maximum recommended load...
  - Hooks 12kg per stand plus
  - Fixed Upright/s 22kg per stand
- 75mm polypropylene / thermoplastic rubber twin-wheel castors (1 x anti-static)

# Dimensions (w x d x h)

- Overall 620 x 620 x 1,300 2,000mm
- Vertical Clamping Area
  - PS/P/SC/4H 975mm
  - PS/P/HC/4H 771mm (upright) 2 x 646mm (support bows)

# **Standard Colours**

- Base Grey
- Upright & Support Bows Stainless Steel
- Vertical Clamping Area
- PS/P/SC/4H 975mm
- PS/P/HC/4H 771mm (upright) 2 x 646mm (support bows)

ItemDescriptionWeightPS/P/SC/4HMobile Pump Stand - Polyurethane, Easy Clean, Standard Capacity, Four Hook11.0kgPS/P/HC/4HMobile Pump Stand - Polyurethane, Easy Clean, High Capacity, Four Hook13.0kg

sales@bristolmaid.com +44 (0)1258 484455





# Mobile Pump Stands - Mild Steel, Four Leg

#### **Features**

- Range of stands designed specifically to hold blood warmers, pumps & syringe drivers
- · Fully welded four leg base
- 32mm diameter 304 grade stainless steel upright & 25mm diameter pump support bows
- Height adjustable upright using a handwheel (1,270 2,000mm)
- Four U-shaped hooks
- Maximum recommended load...
  - Hooks 12kg per stand plus
  - Fixed Upright/s 22kg per stand
- 75mm polypropylene / thermoplastic rubber twin-wheel castors (1 x anti-static)

# Dimensions (w x d x h)

- Overall 560 x 621 x 1,270 2,000mm
- Vertical Clamping Area
  - PS/M/SC/4H 975mm
  - PS/M/HC/4H 771mm (upright) 2 x 646mm (support bows)

# **Standard Colours**

- Base Grey
- Upright & Support Bows Stainless Steel

| Item       | Description  | Weight |
|------------|--|--------|
| PS/M/SC/4H | Mobile Pump Stand - Mild Steel, Four Leg, Standard Capacity, Four Hook | 14.7kg |
| PS/M/HC/4H | Mobile Pump Stand - Mild Steel, Four Leg, High Capacity, Four Hook     | 16.7kg |



www.bristolmaid.com



IS/P/2H/H/SP IS/P/2H/H/BL IS/P/2H/H/GR IS/P/2H/H/RD IS/P/2H/H/YE Two Hook, Handle

IS/P/4H/H/SP IS/P/4H/H/BL IS/P/4H/H/GR IS/P/4H/H/RD IS/P/4H/H/YE Four Hook, Handle

sales@bristolmaid.com

# Mobile Infusion Stands - Polyurethane, Easy Clean

#### **Features**

- Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile Infusion Stands
- · Five leg base with low centre of gravity to assist stability
- · 25mm diameter 304 grade stainless steel upright
- Polyurethane overmoulded base with smooth edges & contours provides...
  - Water-resistant finish, allowing for easy & effective cleaning (including steam cleaning)
  - Shield for the castors
  - Impact protection to patients & staff, reducing the risk of skin abrasions & bruising
  - Protection to the fabric of the building
- Coloured indicator caps provide department identification & segregation, reducing the risk of cross contamination
- Height adjustable upright using a handwheel (1,260 2,100mm)
- Available with two or four pigtail hooks
- Maximum recommended load...
  - Hooks 4kg per stand plus
  - Fixed Upright 4kg per stand
- Optional handle keeps feet clear of the base, assists manoeuvrability
   & improves patient safety
- 50mm polypropylene / thermoplastic rubber twin-wheel non-marking castors (1 x anti-static)

# Dimensions (w x d x h)

- Overall 620 x 620 x 1,260 2,100mm
- Vertical Clamping Area
  - IS/P/2H & IS/P/4H 970mm
  - IS/P/2H/H & IS/P/4H/H 840mm

# **Standard Colours**

- · Base Grey
- · Upright & Handle Stainless Steel
- · Coloured Indicator Cap Spearmint

# Options

• Coloured Indicator Cap - Blue / Green / Red / Yellow



8.2kg

8.3kg





# **Mobile Infusion Stands - Stainless Steel**

#### **Features**

- Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile
   Infusion Stands
- · Seamless welded 304 grade stainless steel construction
- Height adjustable upright using a handwheel (1,270 2,000mm)
- Available with two or four pigtail hooks
- 9kg plastic coated weighted base available for additional stability
- Maximum recommended load...
- · Hooks 4kg per stand plus
- Fixed Upright
  - Non-Weighted Base 2kg per stand
  - Weighted Base 4kg per stand
- 50mm polypropylene / thermoplastic rubber twin wheel castors (1 x anti-static)

# Dimensions (w x d x h)

• 650 x 650 x 1,270 - 2,000mm

### **Standard Colours**


- Base / Upright Stainless Steel
- Weighted Base Grey



sales@bristolmaid.com



www.bristolmaid.com



| Item      | Description                                  | Weight |
|-----------|--|--------|
| DPKS/1    | Mobile Infusion Stand - Knockdown, Two Hook  | 3.8kg  |
| DPKS/1/4H | Mobile Infusion Stand - Knockdown, Four Hook | 3.9kg  |

sales@bristolmaid.com

| Item  | Description             | Weight |
|-------|-------------------------|--------|
| DP/2H | Two Hook Infusion Pole  | 0.6kg  |
| DP/4H | Four Hook Infusion Pole | 0.8kg  |

+44 (0)1258 484455

www.bristolmaid.com

# **Clamping Infusion Stands - Trolley**

#### **Features**

- Clamp secures to 25mm vertical round tube (designed to fit our range of dressing trolleys)
- When fitting to other trolleys, we strongly recommend that a risk assessment is carried out to assess stability
- When fitted protrudes 45mm
- Supplied with a two pigtail hook infusion pole
- Maximum load on hooks 2kg

# Dimensions (w x d x h)

• 55 x 220 x 1200mm





| Item | Description                        | Weight |
|------|------------------------------------|--------|
| DP/C | Clamping Infusion Stands - Trolley | 1.8kg  |

# Wall Mounted Infusion Holders

#### Features

CE

- Wall mounted infusion bag holder
- Available with one/two U shaped hooks
- 25mm diameter "D" shaped tube to accept pumps (max load 10kg)
- Maximum load per hook 2kg
- Fixings not included

# Dimensions (w x d x h)

- IV001 344 x 75 x 500mm
- IV002 344 x 75 x 570mm



| Item  | Description                            | Weight |
|-------|--|--------|
| IV001 | One Hook, Wall Mounted Infusion Holder | 1.2kg  |
| IV002 | Two Hook, Wall Mounted Infusion Holder | 1.6kg  |

Infusion &

CE



## **Linen Trolleys & Bags**

- Soiled Linen Trolleys 224
  - Linen Trolleys 227
  - Hamper Trolleys 228
  - Soiled Linen Bags 229
- Clean Linen Trolleys 230

#### Soiled Linen Trolleys - Stainless Steel, Cantilever Frame



#### **Features**

- · Seamlessly welded from from 304 grade stainless steel for strength & ease of cleaning
- Suitable for use with fabric linen bags using a tubular C shaped ring\*
- · Bags supported on an easy clean solid base
- · Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors

\*Recommended minimum bag circumference 1360mm, height 750mm

#### Dimensions (w x d x h)

- Single 440 x 460 x 820mm
- Double 835 x 460 x 820mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

#### **Options / Accessories**

Linen Bags:

- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grev
- R/BAG/ORANGE Orange

Linen Bags page 229



Pedal operated lid

WT/6/2/PL/SS

Weight

6.5kg

10.0kg

8.5kg

```
www.bristolmaid.com
```

WT/6/PL/SS

Description

Stainless Steel Cantilever, Single

Stainless Steel Cantilever, Double

Stainless Steel Cantilever, Single, Pedal Operated

Item

WT/6/SS

WT/6/2/SS

WT/6/PL/SS

#### Soiled Linen Trolleys - Mild Steel, Cantilever Frame

#### Features

- · Welded mild steel frame
- · Suitable for use with fabric linen bags using a tubular C shaped ring\*
- Bags supported on an easy clean solid base
- Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors

\*Recommended minimum bag circumference 1360mm, height 750mm

#### Dimensions (w x d x h)

- Single 440 x 460 x 820mm
- Double 835 x 460 x 820mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

#### **Standard Colours**

· Paint - Spearmint

#### **Options / Accessories**

Linen Bags:

- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grey
- R/BAG/ORANGE Orange

Linen Bags page 229



Mild Steel Cantilever, Single, Pedal Operated

Mild Steel Cantilever, Double, Pedal Operated

WT/6/PL/MS

WT/6/2/PL/MS

www.bristolmaid.com

8.5ka

14.5kg



#### Soiled Linen Trolleys - Mild Steel, Sprung Frame

#### **Features**

- Welded mild steel frame
- Suitable for use with either fabric or disposable polythene bags using a sprung rod frame\*
- · Bags supported on an easy clean solid base
- Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors

\*Factory set for standard 1360mm circumference bags. Can be manually adjusted to suit other bags. Recommended minimum bag height 750mm

#### Dimensions (w x d x h)

- Single 440 x 460 x 825mm
- Double 835 x 460 x 825mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

#### **Standard Colours**

Paint - Spearmint

#### **Options / Accessories**

- Linen Bags:
  - R/BAG/WHITE White
  - R/BAG/BLUE Blue
  - R/BAG/RED Red
  - R/BAG/GREEN Green
  - R/BAG/YELLOW Yellow
  - R/BAG/GREY Grey
  - R/BAG/ORANGE Orange

Linen Bags page 229

www.bristolmaid.com

10.0kg

14.5kg

WT/9/2/MS

Mild Steel, Sprung Frame, Double

WT/9/2/PL/MS Mild Steel, Sprung Frame, Double, Pedal Operated

#### Soiled Linen Trolley - Mild Steel, Folding

#### **Features**

- · Folds to minimize space
- Suitable for use with fabric linen bags
- 50mm polypropylene / thermoplastic rubber twin-wheel castors
- \*Recommended minimum bag circumference 1360mm, height 750mm

#### Dimensions (w x d x h)

- Overall Dimensions 430 x 430 x 1030mm
- Folded 430 x 200 x 1040mm

#### **Standard Colours**

· Paint - Grey

#### **Options / Accessories**

- Linen Bags:
  - R/BAG/WHITE White
  - R/BAG/BLUE Blue
  - R/BAG/RED Red
  - R/BAG/GREEN Green
  - R/BAG/YELLOW Yellow
  - R/BAG/GREY Grey
  - R/BAG/ORANGE Orange

#### Linen Bags page 229

#### Linen Trolley - Bed Changing

#### **Features**

- Welded mild steel frame
- Three removable, sunken, liquid retaining trays
- Tray heights 230, 590 & 960mm
- Dividing panel separates clean & soiled linen
- Suitable for use with fabric linen bags using a tubular C shaped ring\*
- Bag supported on a base frame
- 100mm pressed steel castors

\*Recommended minimum bag circumference 1360mm, height 750mm

#### Dimensions (w x d x h)

- Overall Dimensions 830 x 490 x 1040mm
- Tray 390 x 400 x 25mm

#### **Standard Colours**

· Paint - Grey

#### **Options / Accessories**

- ZP005 100mm Ring Buffers
- Linen Bags:
  - R/BAG/WHITE White
  - R/BAG/BLUE Blue
  - R/BAG/RED Red
  - R/BAG/GREEN Green
  - R/BAG/YELLOW Yellow
  - R/BAG/GREY Grey
  - R/BAG/ORANGE Orange



| Item    | Description         | Weight |
|---------|---------------------|--------|
| WT/7/MS | Folding, Mild Steel | 3.25kg |

| Item     | Description  | Weight |
|----------|--------------|--------|
| WT/13/MS | Bed Changing | 15.5kg |

WT/7/MS



#### Hamper Trolleys - Mild Steel

#### Features

- Welded mild steel construction
- Suitable for use with fabric hamper bags supported on 4 hooks
- Bag supported on an easy clean solid base
- 100mm pressed steel castors

#### Dimensions (w x d x h)

• Overall Dimensions - 815 x 514 x 1000mm

#### Standard Colours

Paint - Spearmint

#### **Options / Accessories**

- ZP005 100mm Ring Buffers 2FACO013 Hamper Bag (190 litre)

| Item     | Description                 | Weight |
|----------|-----------------------------|--------|
| WT/11/MS | Hamper Trolley - Mild Steel | 12.5kg |

www.bristolmaid.com

#### **Soiled Linen Bags**

#### Features

- Washable 7oz polyester soiled linen bag
  Circumference 1360mm, height 750mm, 120 litres
- Fitted with eyelets and draw cord with plastic toggle
  Available in a range of colours to allow compliance with local, national & international guidelines





| Item         | Description      | Weight |
|--------------|------------------|--------|
| R/BAG/WHITE  | White Linen Bag  | 0.3kg  |
| R/BAG/BLUE   | Blue Linen Bag   | 0.3kg  |
| R/BAG/RED    | Red Linen Bag    | 0.3kg  |
| R/BAG/GREEN  | Green Linen Bag  | 0.3kg  |
| R/BAG/YELLOW | Yellow Linen Bag | 0.3kg  |
| R/BAG/GREY   | Grey Linen Bag   | 0.3kg  |
| R/BAG/ORANGE | Orange Linen Bag | 0.3kg  |



R/BAG/BLUE



R/BAG/YELLOW



R/BAG/GREY



R/BAG/RED



R/BAG/GREEN



R/BAG/ORANGE



#### Clean Linen Trolleys - Small (390 Litre)

#### Features

- Suitable for ward clean linen storage
- Mild steel frame
- Lightweight aluminium cladding & shelves
- Powder coated for durability
- Two adjustable shelves plus base (max 45kg per shelf)
- Shelf size 1000 x 450mm (w x d)
- · Push handles at either end
- Available with a flexible cover or fold back doors, preventing dust ingress
- 125mm pressed steel castors
- Maximum capacity 135kg per trolley

#### Dimensions (w x d x h)

- Overall Dimensions 1175 x 530 x 1065mm
- Internal Dimensions 1000 x 500 x 810mm

#### **Standard Colours**

• Paint - Grey

#### Options

- C008 125mm Castors, Heavy Duty Pressed Steel, 2 x Braking, 2 x Non-Braking (set of 4) Factory fitted
- ZP010 100mm Ring Buffers

| Item        | Description                                  | Weight |
|-------------|--|--------|
| WT/42       | Clean Linen Trolley - Small                  | 39.7kg |
| WT/42/COVER | Clean Linen Trolley - Small, Flexible Cover  | 42.4kg |
| WT/42/FBD   | Clean Linen Trolley - Small, Fold Back Doors | 46.0kg |



#### **Flexible Cover**

- · Fire retardant cover
- Secured to the top of the trolley & drops down to enclose the front
- Weighted bottom hem with velcro style side fastenings



#### **Fold Back Doors**

- Fold back doors encloses the front
- When opened, doors fold against the push handles at either end of the trolley
- Drop down bolt secures base of doors
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)

# -inen Trolleys & Bags

#### AGV Compatible Trolleys & Carts Bristol Maid<sup>™</sup> is your ideal partner for transport trolleys & carts for Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy

distribution



WT/50

sales@bristolmaid.com

| Item             | Description  | Weight  |
|------------------|--|---------|
| Non-Towing       |  |         |
| WT/50            | Clean Linen Trolley - Medium                           | 87.0kg  |
| WT/50/COVER      | Clean Linen Trolley - Medium, Flexible Cover           | 90.0kg  |
| WT/50/FBD        | Clean Linen Trolley - Medium, Fold Back Doors          | 97.3kg  |
| Towing           |  |         |
| WT/50/TBAR       | Clean Linen Trolley - Medium, Tow Bar                  | 89.5kg  |
| WT/50/COVER/TBAR | Clean Linen Trolley - Medium, Flexible Cover, Tow Bar  | 92.5kg  |
| WT/50/FBD/TBAR   | Clean Linen Trolley - Medium, Fold Back Doors, Tow Bar | 100.8kg |

#### Clean Linen Trolleys - Medium (690 Litre)

#### **Features**

- Widely used for bulk linen delivery, part of a linen exchange system, or as an alternative to dedicated ward linen rooms
- Mild steel frame
- Lightweight aluminium cladding & shelves
- Powder coated for durability
- Three adjustable shelves plus base (max 45kg per shelf)
- Shelf size 1180 x 460mm (w x d)
- · Inset push handles at either end
- Low level bumper bar offers protection to both the trolley & fabric of the building
- Available with a flexible cover or fold back doors, preventing dust ingress
- Optional tow bar fitting
- 200mm heavy duty castors, 2 x fixed, 2 x swivel
- Maximum capacity 180kg per trolley

#### Dimensions (w x d x h)

- Overall Dimensions 1430 x 600 x 1495 mm
- Internal Dimensions 1180 x 500 x 1180 mm

#### **Standard Colours**

• Paint - Grey



#### **Tow Bar Fitting**

- Tow bar fitted to the front (swivel castor)
- Antiluce fastener fitted to rear (fixed castor)
- Tow bar position approx. 360 380mm above floor
- Recommended max speed 7kph for towed trolleys



#### **Flexible Cover**

- · Fire retardant cover
- Secured to the top of the trolley & drops down to enclose the front
- Weighted bottom hem with velcro style side fastenings



WT/50/FBD

#### **Fold Back Doors**

- · Fold back doors, encloses the front
- When opened, doors fold back flat at either end of the trolley\*
- Drop down bolt secures base of doors
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)

\*Doors will not fold back flat when a tow bar is fitted

#### AGV Compatible Trolleys & Carts

For Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy distribution



| Item             | Description   | Weight  |
|------------------|---|---------|
| Non-Towing       |   |         |
| WT/48            | Clean Linen Trolley - Large                           | 99.6kg  |
| WT/48/COVER      | Clean Linen Trolley - Large, Flexible Cover           | 102.9kg |
| WT/48/FBD        | Clean Linen Trolley - Large, Fold Back Doors          | 111.5kg |
| Towing           |   |         |
| WT/48/TBAR       | Clean Linen Trolley - Large, Tow Bar                  | 103.1kg |
| WT/48/COVER/TBAR | Clean Linen Trolley - Large, Flexible Cover, Tow Bar  | 106.4kg |
| WT/48/FBD/TBAR   | Clean Linen Trolley - Large, Fold Back Doors, Tow Bar | 115.0kg |

#### Clean Linen Trolleys - Large (1,100 Litre)

#### **Features**

- Widely used for bulk linen delivery, part of a linen exchange system, or as an alternative to dedicated ward linen rooms
- Mild steel frame
- Lightweight aluminium cladding & shelves
- Powder coated for durability
- Three adjustable shelves plus base (max 45kg per shelf)
- Shelf size 1430 x 555mm (w x d)
- · Inset push handles at either end
- Low level bumper bar offers protection to both the trolley & fabric of the building
- Available with a flexible cover or fold back doors, preventing dust ingress
- Optional tow bar fitting
- 200mm heavy duty castors, 2 x fixed, 2 x swivel
- Maximum capacity 180kg per trolley

#### Dimensions (w x d x h)

- Overall Dimensions 1710 x 700 x1620mm
- Internal Dimensions 1420 x 610 x 1295mm

#### **Standard Colours**

· Paint - Grey



#### **Tow Bar Fitting**

- Tow bar fitted to the front (swivel castor)
- Antiluce fastener fitted to rear (fixed castor)
- Tow bar position approx. 360 380mm above floor
- Recommended max speed 7kph for towed trolleys

www.bristolmaid.com



# 0 e 6 o

#### **Flexible Cover**

- · Fire retardant cover
- Secured to the top of the trolley & drops down to enclose the front
- Weighted bottom hem with velcro style side fastenings



#### **Fold Back Doors**

- Fold back doors, encloses the front
- When opened, doors fold back flat at either end of the trolley\*
- Drop down bolt secures base of doors
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)

\*Doors will not fold back flat when a tow bar is fitted



# **Maternity Equipment**

- Baby Cribs Variable Height, Battery Operated 238
  - Baby Cribs Fixed Height 240
- Baby Crib Variable Height, Gas Strut Operation 242



#### **Baby Crib - Variable Height, Battery Operation**

#### **Features**

- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Rechargeable battery powered height adjustment (700 910mm)
- Crib can be tilted at either end (horizontal to approx. 6 degrees)
- Push handles at either end
- Lower cupboard with fixed shelf
- Two lift up work flaps (max capacity 2.0kg)
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- Supplied with one rechargeable battery & handset\*
- 100mm pressed steel braking castors

\*Battery must be removed from the unit for charging (see options / accessories for additional batteries)

#### Dimensions (w x d x h)

- Overall Dimensions 940 x 510 x 700 910mm
- Mattress 600 x 330 x 60mm
- Cupboard 210 x 265 x 315mm (internal, either side of the column)
- Work Flap 415 x 305mm

#### **Standard Colours**

- · Paint Spearmint
- Cabinet / Lower Shelf / Work Flaps Spearmint / Grey

#### **Accessories / Spares**

- 2LIDE033 Additional Rechargeable Battery
- 2LIDE037 Charger for Rechargable Battery
- 2PLCO014 Clear Plastic Crib, Solid Side (Spare / Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare / Replacement)

| Item     | Description                       | Weight |
|----------|-----------------------------------|--------|
| BA/VH    | Battery Operated, With Charger    | 49.5kg |
| BA/VH/NC | Battery Operated, Without Charger | 49.5kg |

```
184455 v
```

```
www.bristolmaid.com
```





Suppled with a handset



Rechargeable battery located at rear of crib



#### **Baby Cribs - Fixed Height**

#### Features

- Available with...
  - Lower shelf
  - · Lower shelf & small cupboard (with adjustable shelf)
  - Large cupboard (with adjustable shelf)
- Crib can be tilted at either end (horizontal to approx. 6 degrees)
- Two lift up work flaps (max capacity 2.0kg)
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- 100mm pressed steel castors, 2 off braking

#### **Dimensions** (w x d x h)

- Overall Dimensions 952 x 500 x 929mm
- Mattress 600 x 330 x 60mm
- Large Shelf 810 x 400mm at height 270mm (BA/3/B)
- Small Shelf 460 x 355mm at height 270mm (BA/6)
- Small Cupboard 355 x 420 x 345mm (internal) at height 270mm (BA/6)
- Large Cupboard 710 x 420 x 345mm (internal) at height 270mm (BA/6/A)
- Work Flap 415 x 305mm

#### **Standard Colours**

- · Paint Grey White
- Cabinet / Lower Shelf / Work Flaps Spearmint / Grey White

#### **Accessories / Spares**

- 2PLCO014 Clear Plastic Crib, Solid Side (Spare / Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare / Replacement)

| Item   | Description                  | Weight |
|--------|------------------------------|--------|
| BA/3/B | Fixed Height, Lower Shelf    | 23.5kg |
| BA/6   | Fixed Height, Small Cupboard | 34.0kg |
| BA/6/A | Fixed Height, Large Cupboard | 42.0kg |





#### **Baby Crib - Variable Height, Gas Strut Operation**

#### **Features**

- Gas damper powered height adjustment (710 1025mm)
- Crib can be tilted at either end (horizontal to approx. 6 degrees)
- · Push handles at either end
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- 75mm polymer twin-wheel braking castors

#### **Dimensions** (w x d x h)

- Overall Dimensions 930 x 476 x 775 1025mm
- Mattress 600 x 330 x 60mm

#### **Standard Colours**

- · Paint Spearmint / Grey
- · Work Flap Spearmint

#### **Accessories / Spares**

- BA902 Slide Out Debris Bin
- BA903 Wire Basket Internal Dimensions 420 x 170 x 220mm
- BA904 Lift Up Flap Dimensions 400 x 300mm
- 2PLCO014 Clear Plastic Crib (Spare / Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare / Replacement)





#### Other products that you may be interested in...



#### **Overnight Chair - Portland, Three Position**





### **Medical Records**

- Medical Records Trolleys Locking 246
- Medical Records Workstations Locking 248
  - Laptop Holders 250
  - Carefile Filing Pockets 252
  - Medical Records Trolleys Filing 254
  - Medical Records Transfer Trolleys 256





#### **Medical Records Trolleys - Locking**

#### **Features**

- Suitable for use with our popular narrow filing pockets (order separately)
- Filing capacity
  - Small 470mm
  - Medium 695mm
  - Large 940mm
- Locking lid/s act as a work surface when open & is fitted with either a...
  - CAM Lock Supplied with 2 keys
  - Electronic Push Button Lock Programmable four-digit code
- Maximum load on locking lid/s 2.0kg
- Document pocket/s fitted to end of trolley
- Optional ring buffers offer protection to both the trolley & fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

#### Dimensions (w x d x h)

- Small 671 x 522 x 1056mm (width open 1190mm)
- Medium 920 x 522 x 1056mm (width open 1420mm)
- Large 1240 x 522 x 1056mm (width open 2085mm)

#### **Standard Colours**

- Paint Spearmint / Grey
- Lid / Flap Spearmint

#### **Options / Accessories**

- CF/F/PK10 Foolscap Filing Pockets (x 10)
- R/C/T Paper Tabs (x 500)
- R/C/TH Tab Holders (x 25)
- ZP010 Ring Buffers Factory Fitted Only
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C005 100mm Polymer Castor, 4 x Braking (set of 4) Factory Fitted Only

| High Security Bolt Lock | Electronic Push Button Lock | Description | Weight |
|-------------------------|-----------------------------|-------------|--------|
| MR210                   | MR210CL                     | Small       | 27.2kg |
| MR310                   | MR310CL                     | Medium      | 34.4kg |
| MR410                   | MR410CL                     | Large       | 46.4kg |



02

56 78

00

\*\*\*\*\*\*\*\*\*\*\*\*

#### www.bristolmaid.com



#### **Electronic Push Button Lock**

Uses a 10-button keypad

11-12

- Easy programming & code change
- Four operating codes with differing levels of authority -Master / Sub-Master / User / Technician

MRW/S1DXP/CAM

- User codes are 4 digits
- · Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

#### **Medical Records Workstations - Locking**

#### Features

- Suitable for use with our filing pockets (order separately)
- Provides 375mm filing capacity per drawer
- Available as single & double column units with...
  - One / two drawers & x-ray pigeon hole (5 / 10 sections)
  - One / two drawers & open cupboard
  - Two / four drawers
- Two locking options...
  - · CAM Lock Supplied with two keys
- · Electronic Push Button Lock Programmable four-digit code
- Where relevant one lock secures all drawers
- Double drawer units fitted with an anti-tilt mechanism preventing two drawers from being opened at the same time
- Max load per drawer 16.0kg
- Push handle & document pocket fitted to end of trolley
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Single Column 604 x 538 x 1028mm
- Double Column 1056 x 538 x 1034mm

#### **Standard Colours**

- Paint Grey White
- Cabinet Spearmint / Grey White

MRW/S1DOC/CAM

#### Single Column

MRW/S2D/CAM

fitted with MRW/LHR/S

(laptop not included)

| High Security Bolt Lock | Electronic Push Button Lock | Description                   | Weight |
|-------------------------|-----------------------------|-------------------------------|--------|
| MRW/S1DXP/CAM           | MRW/S1DXP/CL                | One Drawer, X-Ray Pigeon Hole | 34.2kg |
| MRW/S1DOC/CAM           | MRW/S1DOC/CL                | One Drawer, Open Cupboard     | 34.0kg |
| MRW/S2D/CAM             | MRW/S2D/CL                  | Two Drawers                   | 40.6kg |



#### **Options / Accessories**

MRW/D2DXP/CAM

- MRW/LHR/S Laptop Holder, Removable, Swing Arm
   MRW/LHNR/S Laptop Holder, Non-Removable, Swing Arm
- CF/F/PK10 Foolscap Filing Pockets (x 10)
- R/C/T Paper Tabs (x 500)
- R/C/TH Tab Holders (x 25)
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
   C005 100mm Polymer Castor, 4 x Braking (set of 4) Factory Fitted Only

#### **Double Column**

| High Security Bolt Lock | Electronic Push Button Lock | Description                    | Weight |
|-------------------------|-----------------------------|--------------------------------|--------|
| MRW/D2DXP/CAM           | MRW/D2DXP/CL                | Two Drawers, X-Ray Pigeon Hole | 62.0kg |
| MRW/D2DOC/CAM           | MRW/D2DOC/CL                | Two Drawers, Open Cupboard     | 58.0kg |
| MRW/D4D/CAM             | MRW/D4D/CL                  | Four Drawers                   | 72.0kg |

sales@bristolmaid.com +44 (0)1258 484455

MRW/D2DOC/CL

MRW/D4D/CAM

-2

#### **Laptop Holders**



#### Removeable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MRW range of medical record workstations
- Fitted to an upright bolted to the rear left-hand corner of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle, fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

\*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**

• Paint - Grey White



| Item      | Description          | Weight |
|-----------|----------------------|--------|
| MRW/LHR/S | Removable, Swing Arm | 3.7kg  |



#### Non-Removeable, Swing Arm

#### **Features**

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MRW range of medical record workstations
- Fitted to an upright bolted to the rear left-hand corner of the trolley (upright & holder supplied for customer to fit)
- Laptop supported in a cradle fitted with a steel mounting bracket, fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)\*
- Adjustable to accept laptops 280 435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Design allows easy access to laptop ports

\*Not compatible with Dell laptops using the N17 wedge lock system

#### **Standard Colours**

• Paint - Grey White



| Item       | Description              | Weight |
|------------|--------------------------|--------|
| MRW/LHNR/S | Non-Removable, Swing Arm | 5.0kg  |




#### **Carefile Filing Pockets**

#### **Features**

- Filing pockets suitable with Bristol Maid trolleys / workstations & most standard filing suspension systems
- Provides a durable alternative to paper or card types
- Suitable for runners with 380mm spacings / centres
- Provides an internal pocket (w x h) 360 x 260mm
- Fire retardant PVC coated nylon pockets
- Unique 'stay-put' tab holders, with pink / blue reversible card tabs
- Supplied in packs of ten

#### **Standard Colours**

• Nylon Pockets - Grey



| Item      | Description                    |
|-----------|--------------------------------|
| CF/F/PK10 | Foolscap Filing Pockets (x 10) |
| R/C/T     | Paper Tabs (x 500)             |
| R/C/TH    | Tab Holders (x 25)             |

#### Other products that you may be interested in...

See Page 339

#### Medium & High Back Patient Armchairs - Alderbury

5ND/MAR/ALDE/VVDEL





#### **Medical Records Filing Trolleys - Tiered**

#### **Features**

- Suitable for use with ring binders, case note & lever arch files (max file height 350mm)
- Available in two sizes with either an open front or locking door/s
- Small 900mm filing capacity
- Large 1800mm filing capacity Adjustable dividers (small 4 / large 8)
- Doors fold against the sides of the trolley (not flat)
- Large trolleys use one lock to secure both doors
- Two locking options...
  - · CAM Lock Supplied with two keys
  - Electronic Push Button Lock Programmable four-digit code
- Optional ring buffers offer protection to both the trolley & fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

#### Dimensions (w x d x h)

- Small 546 x 440 x 1016mm
- Large 998 x 440 x 1016mm

#### **Standard Colours**

- · Paint Spearmint / Grey
- Top Spearmint

| Item         | Description  | Weight |
|--------------|--|--------|
| Small        |  |        |
| MRF/S/OC     | Small - Open Cupboard                              | 21.5kg |
| MRF/S/LD/CAM | Small - Locking Door, CAM Lock                     | 26.5kg |
| MRF/S/LD/CL  | Small - Locking Door, Electronic Push Button Lock  | 26.5kg |
| Large        |  |        |
| MRF/L/OC     | Large - Open Cupboard                              | 30.0kg |
| MRF/L/LD/CAM | Large - Locking Doors, CAM Lock                    | 41.0kg |
| MRF/L/LD/CL  | Large - Locking Doors, Electronic Push Button Lock | 41.0kg |
|              |  |        |

- MRF/LHR/S Laptop Holder Removable, Swing Arm
- MRF/LHNR/S Laptop Holder Non-Removable, Swing Arm
- ZP010 Ring Buffers Factory Fitted Only
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4) Factory Fitted Only
- C005 100mm Polymer Castor, 4 x Braking (set of 4) Factory Fitted Only
- C006 100mm Polymer Castor, 4 x Non-Braking (set of 4) Factory Fitted Only

\$

• Low battery warning with failure override using 9v battery

across LED's (+/-) & master code



#### Capacity

|                | Case Note Files  | Ring Binders     | Lever Arch Files |
|----------------|------------------|------------------|------------------|
|                | 235 x 50 x 320mm | 250 x 40 x 315mm | 290 x 75 x 315mm |
| Small - 900mm  | 18               | 22               | 12               |
| Large - 1800mm | 36               | 44               | 24               |

#### Medical Records Transfer Trolleys Single Box



#### **Features**

- Suitable for use the transportation of bulk case notes & x-rays
- Records stored in an upper box
- Lower tray
- Two locking options...
  - CAM lock supplied with 2 keys
  - Electronic Push Button lock Programmable four-digit code
- Trolley max capacity 50kg
- Optional ring buffers offer protection to both the trolley & fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

#### Dimensions (w x d x h)

- Overall Dimensions 817 x 523 x 975mm
- Box (internal) 600 x 450 x 250mm
- Top of Box (to floor) 805mm

#### **Standard Colours**

• Paint - Spearmint / Grey

#### **Options**

• ZP005 - Ring Buffers

| Item    | Description                              | Weight |
|---------|--|--------|
| Open    |  |        |
| MR130   | Open Top                                 | 16.0kg |
| MR140   | Locking Lid, CAM Lock                    | 19.0kg |
| MR140CL | Locking Lid, Electronic Push Button Lock | 19.0kg |

#### Medical Records Transfer Trolley Double Box, Extending

#### **Features**

- Trolley for the transportation of bulk case notes & x-rays
- · Records stored in two removable boxes
- Clearance between upper & lower boxes, 235mm
- Max capacity 25kg (per box)
- · 80mm ring buffers fitted to front of unit
- 125mm pressed steel castors & 200mm wheels (braking)

#### Dimensions (w x d x h)

- Overall Dimensions 875 x 598 x 1040mm
- Box (internal) 568 x 368 x 316mm
- Top of Box (to floor)
  - Lower 440mm
  - Upper 975mm

#### **Standard Colours**

· Paint - Spearmint

Alternative

Colours

Page 450

Description

Double Box, Extending

Item

MR110

#### **Options / Accessories**

• 2PLCO271 - Removable Box

#### Medical Records Transfer Trolley Double Box, Duo

#### **Features**

- · Unique 'two' trolley system for the bulk transportation of case notes & x-rays
- Records stored in fixed boxes, each supplied with an adjustable divider
- · Clearance between upper & lower boxes, 87mm
- Trolleys are loaded independently, before being linked together & transported
- Trolleys link together using a simple clevis pin
- Max capacity 25kg (per box)
- 100mm pressed steel castors (2 off braking) & 200mm wheels (braking)

#### Dimensions (w x d x h)

- Overall Dimensions 866 x 522 x 970mm
- Box (internal) 600 x 420 x 255mm
- Top of Box (to floor)
- Lower 527mm
- Upper 875mm

#### **Standard Colours**

Paint - Spearmint / Grey



www.bristolmaid.com

| em   | Description     | Weight |
|------|-----------------|--------|
| R120 | Double Box, Duo | 34.6kg |

Weight

31.6kg

Ite

MF



# Phlebotomy

- Caretray Trolley Phlebotomy 260
- Phlebotomy Chairs Variable Height 262
- Phlebotomy Chairs Variable Height, Bariatric 263
  - Phlebotomy Chairs Fixed Height 264



## **Caretray Trolley - Phlebotomy**

### **Features**

- Trollevs designed for use with Phlebotomy
- Quick & easy disassembly, assists thorough & deep cleaning
- Available in mild or stainless steel (requires no on-site assembly)
  - Fully welded mild steel frame, powder coated in a durable epoxy paint
  - · Seamlessly welded, one piece 304 grade stainless steel frame
- Ergonomic push handle fitted to one side of the trolley ensures correct posture & keeps feet away from the frame, ideal for mobile phlebotomy or in cases where the trolley is wheeled for long periods of time
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Available with a solid grade laminate door & back panel providing security for the contents of the travs
- Two locking options...
  - · CAM Lock with Wing Handle Supplied with two keys
  - · Electronic Push Button Lock Programmable four-digit code
- Trollev pre-configured with the following accessories...
  - Alcohol Gel Holder
  - Debris Bin & Bracket
  - Glove / Mask Box Holder, One Box
  - · Sharps Box Holder
- Translucent, impact-resistant travs & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each trav
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### CTPHLEM1/10/3S2D

### Mild Steel

| Description                | No Door          | Door with CAM Lock  | Door with Electronic<br>Push Button Lock | Weight |
|----------------------------|------------------|---------------------|--|--------|
| 3 x 100mm, 2 x 150mm Trays | CTPHLEM1/10/3S2D | CTPHLEM1/10/3S2D/CA | CTPHLEM1/10/3S2D/EPB                     | 31.1kg |
| 4 x 100mm, 1 x 150mm Trays | CTPHLEM1/10/4S1D | CTPHLEM1/10/4S1D/CA | CTPHLEM1/10/4S1D/EPB                     | 30.8kg |
| 6 x 100mm Trays            | CTPHLEM1/10/6S   | CTPHLEM1/10/6S/CA   | CTPHLEM1/10/6S/EPB                       | 31.9kg |

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com

### Dimensions (w x d x h)

- Overall Dimensions 610 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

### **Standard Colours**

- Mild Steel Paint Spearmint
- Stainless Steel Natural Finish

### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
- Defibrillator Platform page 67
- Carerail Mounted Accessories page 68
- Laptop Holders page 74
- Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77



### **Stainless Steel**

| Description                | No Door          | Door with CAM Lock  | Door with Electronic<br>Push Button Lock | Weight |
|----------------------------|------------------|---------------------|--|--------|
| 3 x 100mm, 2 x 150mm Trays | CTPHLES1/10/3S2D | CTPHLES1/10/3S2D/CA | CTPHLES1/10/3S2D/EPB                     | 31.1kg |
| 4 x 100mm, 1 x 150mm Trays | CTPHLES1/10/4S1D | CTPHLES1/10/4S1D/CA | CTPHLES1/10/4S1D/EPB                     | 30.8kg |
| 6 x 100mm Trays            | CTPHLES1/10/6S   | CTPHLES1/10/6S/CA   | CTPHLES1/10/6S/EPB                       | 31.9kg |

CTPHLES1/10/3S2D/CA



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric - Foot Switch in Graphite - EC5315/VGH



| Item       | Description                                   | Weight |
|------------|---|--------|
| EC5305/VBB | Hydraulic                                     | 67.0kg |
| EC5310/VBB | Electric - Hand Switch                        | 67.0kg |
| EC5315/VBB | Electric - Foot Switch                        | 67.0kg |
| EC5320/VBB | Electric with Electric Backrest - Hand Switch | 69.3kg |
| EC5325/VBB | Electric with Electric Backrest - Foot Switch | 69.3kg |
|            |   |        |

## Phlebotomy Chair Three Section

### Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (selfassembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

### Dimensions (w x d x h)

• 770 x 1710 x 490 - 950mm

### **Options**

· E06 - Battery Back Up

### **Standard Colours**

- Paint White
- Vinyl Bristol Blue

sales@bristolmaid.com +44 (0)1258 484455

360mm

470mm

630mm

820mm

## Phlebotomy Chair Three Section, Bariatric

### **Features**

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- · Height adjustment from 490 950mm
- · Adjustable angle backrest
  - Manual -25 to +90 degrees
  - Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- · One adjustable foot fitted to allow levelling on uneven floors
- · Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant anti-microbial vinyl
- · Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (selfassembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

### Dimensions (w x d x h)

• 990 x 1710 x 520 - 950mm

### **Standard Colours**

- · Paint White
- · Vinyl Bristol Blue

### Options

· E18 - Battery Back Up



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Onyx Black - EC5340/VOB

www.bristolmaid.com





| Item       | Description                                   | Weight |
|------------|---|--------|
| EC5330/VBB | Electric - Hand Switch                        | 77.2kg |
| EC5335/VBB | Electric - Foot Switch                        | 77.2kg |
| EC5340/VBB | Electric with Electric Backrest - Hand Switch | 79.5kg |
| EC5345/VBB | Electric with Electric Backrest - Foot Switch | 79.5kg |



## Phlebotomy Chair Fixed Height

### Features

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- BS 4875: Part1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Certified for use up to 32 stone (203kg)
- Stable four leg design
- Black finish tubular frame
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrest
- Armrest can be positioned on either side of the chair (additional armrest available if required)
- Maximum capacity 203kg

### Dimensions (w x d x h)

• 620 x 580 x 800mm

### **Standard Colours**

- Paint Black
- Vinyl Bristol Blue

### Options

• E15 - Phlebotomy Armrest (Single)

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Fixed Height Phlebotomy Chair in Aquamarine - EC5000/VAQ





# **Portering Chairs**

- Portering Chair Rear Steer 268
- Portering Chair Front Steer 269
- Portering Chair Bariatric, Rear Steer 271
  - Options & Accessories 272



## Portering Chair Rear Steer, Nesting, Sliding Foot Rest

### Features

- Nests to occupy reduced floor space & encourage tidy storage
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- Design ensures chair will not tip should the patient stand on the step
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- · Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 200mm fixed wheels/150mm total stop braking castors (cushioned)
- Maximum capacity 190.0kg

### Dimensions (w x d x h)

- Overall Dimensions 660 x 840 x 1040mm
- Seat 540 x 420mm at Height 530mm
- Back 490 x 455mm

### **Standard Colours**

- Paint Grey
- Upholstery Bristol Blue

### **Options / Accessories**

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cylinder Holder, E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/SB Lap Seatbelt, Factory Fitted
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card

Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Rear Steer, Nesting, Sliding Foot Rest in Bristol Blue - G/200/RS/VBB



| ltem         | Description                            | Weight |
|--------------|--|--------|
| G/200/RS/VBB | Rear Steer, Nesting, Sliding Foot Rest | 29.0kg |

sales@bristolmaid.com

+44 (0)1258 484455

X-Ray/Notes Pocket

## Portering Chair Front Steer, Nesting, Sliding Foot Rest

### **Features**

- · Nests to occupy reduced floor space & encourage tidy storage
- · Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- · Design ensures chair will not tip should patient stand on the step
- X-ray/notes pocket on rear of seat
- · Two hooks for patients belongings
- · Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 125mm swivel castors/200mm braking fixed wheels (cushioned)
- Maximum capacity 190.0kg
- Ergonomically designed over swept push handle gives greater control making chair easier to push/use

### Dimensions (w x d x h)

- Overall Dimensions 660 x 1005 x 1090mm
- Seat 540 x 420mm at Height 600mm
- Back 490 x 455mm

### **Standard Colours**

- · Paint Grey
- Upholstery Bristol Blue

### **Options / Accessories**

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cylinder Holder, E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/SB Lap Seatbelt, Factory Fitted
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card



X-Ray/Notes Pocket



Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Front Steer, Nesting, Sliding Foot Rest in Lime - G/305/FS/VLI

|  | VGH |  |  |  |
|--|-----|--|--|--|

| Item         | Description                             | Weight |
|--------------|---|--------|
| G/305/FS/VBB | Front Steer, Nesting, Sliding Foot Rest | 29.0kg |

sales@bristolmaid.com

www.bristolmaid.com

+44 (0)1258 484455



## **Portering Chair** Front Steer, Nesting, Hinged Foot Rests

### **Features**

- Nests to occupy reduced floor space & encourage tidy storage
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & helpers
- · Folding/hinging foot rests aids access as well as supporting feet during transit
- · X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 125mm swivel castors/200mm braking fixed wheels (cushioned)
- Maximum capacity 190.0kg

### Dimensions (w x d x h)

- Overall Dimensions 660 x 710 x 1055mm
- Seat 540 x 420mm at Height 530mm
- Back 490 x 455mm

### **Standard Colours**

- · Paint Grev
- Upholstery Bristol Blue

### **Options / Accessories**

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cvlinder Holder, E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- · G/SB Lap Seatbelt, Factory Fitted
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card



| Item         | Description                             | Weight |
|--------------|---|--------|
| G/205/FS/VBB | Front Steer, Nesting, Hinged Foot Rests | 27.5kg |

+44 (0)1258 484455

X-Rav/Notes Pocket

## Portering Chair Bariatric, Rear Steer, Sliding Foot Rest

### **Features**

- · Heavy duty bariatric portering chair
- · Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients & helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- · Design ensures chair will not tip should the patient stand on the step
- · X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- · Fittings to accept IV poles mounted to rear uprights
- · Retractable diagonal seatbelt
- 200mm fixed wheels/200mm total stop braking castors (cushioned)
- We recommend that two members of staff are in attendance when used with a patient in excess of 190kg
- Maximum capacity 380.0kg

### Dimensions (w x d x h)

- Overall Dimensions 930 x 1000 x 1090mm
- Seat 740 x 670mm at Height 530mm
- Back 700 x 480mm

### **Standard Colours**

- · Paint Grey
- Upholstery Bristol Blue

### **Options / Accessories**

- G/500/RS/CYL Cylinder Holder, CD/D
- G/CLS/500 Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card



X-Ray/Notes Pocket





Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Front Steer, Nesting, Sliding Foot Rest in Lime - G/305/FS/VLI

|  | VGH |  |  |  |
|--|-----|--|--|--|

| Item         | Description                              | Weight |
|--------------|--|--------|
| G/500/RS/VBB | Bariatric, Rear Steer, Sliding Foot Rest | 43.8kg |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com

**Portering Chairs** 

## **Cylinder Holders**

### **Features**

- Compatible with 100mm diameter cylinders (C/CD/D/E)
- . Mounts to IV support tube at rear of chair using a handwheel
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models

### **Standard Colours**

Paint - Grey



| Item         | Description             | Weight |
|--------------|-------------------------|--------|
| CYHP/D1      | C/CD/D Size             | 1.0kg  |
| CYHP/E1      | E Size                  | 1.8kg  |
| G/500/RS/CYL | Bariatric - C/CD/D Size | 1.0kg  |

## Calf / Leg Support - Pair

### **Features**

- Moulded calf/leg support on captive retractable support tube
- Retracts beneath the seat when not in use
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models
- Supplied as a pair

### **Standard Colours**

Stainless Steel/Black



| Item      | Description                      | Weight |
|-----------|----------------------------------|--------|
| G/CLS     | Calf/Leg Support, Pair           | 2.3kg  |
| G/CLS/500 | Bariatric Calf/Leg Support, Pair | 2.3kg  |

## **Infusion Pole**

### **Features**

- Stainless steel infusion pole with a U shaped hook (length of pole 980mm)
- Mounts to IV support tube at rear of chair using a handwheel
- · Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models



## Lap Seatbelt

### **Features**

- Lap seat belt with snap buckle
- Compatible with G/200/RS, G/205/FS & G/305/FS models
- · Factory fitted

### **Standard Colours**

Black





| Item | Description  | Weight |
|------|--------------|--------|
| G/SB | Lap Seatbelt | 0.4kg  |

## **Card/Coin Lock**

### Features

- Compatible with old & new one pound coin (euro compatible)
- Coin/card remains in the lock until the chair is connected to another chair or the head of station
- Plastic release card available for use by porters
- Factory fitted
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models
- Replacement card/coin lock cassette available (requires cassette removal tool)



| Item        | Description                                     | Weight |
|-------------|---|--------|
| G/CARD/LOCK | Card/Coin Lock, Factory Fitted                  | 1.0kg  |
| 2LOKE011    | Plastic Release Card                            | 0.01kg |
| 2LOCK217    | Card/Coin Lock Cassette, New £1 Coin Compatible | 0.2kg  |
| 2LOCK218    | Cassette Removal Tool                           | 0.01kg |

## ItemDescriptionWeightG/CARD/STATIONHead of Station1.6kg

### **Features**

- Head of station for chairs fitted with card/coin locks
- Requires fitting to the fabric of the building (wall fixings not supplied)
- Chain length 1000mm
- One required for each row of chairs, (G/200/RS, G/205/FS & G/305/FS)
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models





## **Procedure Trolleys & Carts**

- Caretray Trolley Anaesthesia 276
  - Caretray Trolley Cast 278
- Caretray Trolley Critical Care 280
- Caretray Trolley Infection 282
- Caretray Trolley IV Preparation 284
- Caretray Trolley Medical / Surgical 286
  - Caretray Trolley Plaster 288
- Caretray Trolley Vaccination290Caretray Trolley Basics292Mobile Workstation Mild Steel, Low Level294Procedure Carts Polymer, Metro Flexline™296Procedure Carts Polymer, Metro Starsys™298Procedure Carts Polymer, Metro Basics™300

#### CTANAEM1/10/3S2D/HSB

| Item                    | Description   | Weight |
|-------------------------|---|--------|
| High Security Bolt Lock |   |        |
| CTANAEM1/10/3S2D/HSB    | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 41.4kg |
| CTANAEM1/10/4S1D/HSB    | 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 41.1kg |
| CTANAEM1/10/6S/HSB      | 6 x 100mm Trays, High Security Bolt Lock            | 42.2kg |

## **Caretrav Trollev - Anaesthesia**

### **Features**

- Trollevs suitable for use with Anaesthesia
- Quick & easy disassembly, assists thorough & deep cleaning
- · Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trollev pre-configured with the following accessories...
  - Overbridge, Mild Steel
  - · Four Compartment Tilt Bin, Mild Steel
  - · Five Compartment Tilt Bin, Mild Steel
  - Waste Bin & Bracket, 8 Litre
- · Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- · Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- · 100mm polymer castors, 2 off braking

### Standard Colours · Paint - Grey White

### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Travs
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

- Locking & Security page 60
  Overbridge & Tilt Bins page 66
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTANAEM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 41.4kg |
| CTANAEM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 41.1kg |
| CTANAEM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 42.2kg |







| Item                      | Description   | Weight |
|---------------------------|---|--------|
| CAM Lock with Wing Handle |   |        |
| CTCASTM1/10/3S2D/CA       | 3 x 100mm, 2 x 150mm Trays, CAM Lock with Wing Handle | 43.5kg |
| CTCASTM1/10/4S1D/CA       | 4 x 100mm, 1 x 150mm Trays, CAM Lock with Wing Handle | 43.2kg |
| CTCASTM1/10/6S/CA         | 6 x 100mm Trays, CAM Lock with Wing Handle            | 44.3kg |

## **Caretray Trolley - Cast**

### **Features**

- Trolleys suitable for use within A&E / Plaster Room
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · CAM Lock with Wing Handle Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Overbridge, Mild Steel
  - Four Compartment Tilt Bin, Mild Steel
  - Bowl Support & Bracket
  - Stainless Steel Bowl, 7 Litre
  - Waste Bin & Bracket, 8 Litre
  - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

### **Standard Colours**

· Paint - Grey White

### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

sales@bristolmaid.com +44 (0)1258 484455

- Locking & Security page 60
  Overbridge & Tilt Bins page 66
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTCASTM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 43.7kg |
| CTCASTM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 43.4kg |
| CTCASTM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 44.5kg |



CTCASTM1/10/3S2D/EPB



| Item                      | Description   | Weight |
|---------------------------|---|--------|
| CAM Lock with Wing Handle |   |        |
| CTCRITM1/10/3S2D/CA       | 3 x 100mm, 2 x 150mm Trays, CAM Lock with Wing Handle | 44.6kg |
| CTCRITM1/10/4S1D/CA       | 4 x 100mm, 1 x 150mm Trays, CAM Lock with Wing Handle | 44.3kg |
| CTCRITM1/10/6S/CA         | 6 x 100mm Trays, CAM Lock with Wing Handle            | 45.4kg |

## **Caretray Trolley - Critical Care**

### **Features**

- Trolleys suitable for use with Critical Care
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · CAM Lock with Winged Handle Supplied with two keys
  - Electric Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Overbridge, Mild Steel
  - Five Compartment Tilt Bin, Mild Steel
  - Six Compartment Tilt Bin, Mild Steel
  - Alcohol Gel Holder
  - Infusion Pole & Bracket
  - Waste Bin & Bracket, 8 Litre
  - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

### **Standard Colours**

· Paint - Grey White

### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

- Locking & Security page 60
  Overbridge & Tilt Bins page 66
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTCRITM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 44.8kg |
| CTCRITM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 44.5kg |
| CTCRITM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 45.6kg |



CTCRITM1/10/3S2D/EPB

σ



CTINFEM1/10/3S2D/HSB

| Description   | Weight   |
|---|--|
|   |  |
| 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 35.5kg   |
| 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 35.2kg   |
| 6 x 100mm Trays, High Security Bolt Lock            | 36.3kg   |
|   | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock<br>4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock |

sales@bristolmaid.com

## **Caretray Trolley - Infection**

### **Features**

- Trolleys designed for use when treating patients that have been infected with multi-resistant pathogens
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Waste Bin & Bracket, 8 Litre
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

### **Standard Colours**

- Dimensions (w x d x h)
- Paint Grey White
   Overall Dimensions 550 x 550 x 993mm
  - Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
  - Plastic Top 455 x 455 (internal) at height 968mm
  - Trays
    - Shallow Tray 355 x 405 x 100mm (internal)
    - Deep Tray 355 x 405 x 150mm (internal)

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
- Laptop Holders page 74
- Castors & Ring Buffers page 76
   Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTINFEM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 35.5kg |
| CTINFEM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 35.2kg |
| CTINFEM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 36.3kg |



CTINFEM1/10/3S2D/EPB



| Item                    | Description   | Weight |
|-------------------------|---|--------|
| High Security Bolt Lock |   |        |
| CTIVPRM1/10/3S2D/HSB    | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 39.8kg |
| CTIVPRM1/10/4S1D/HSB    | 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 39.5kg |
| CTIVPRM1/10/6S/HSB      | 6 x 100mm Trays, High Security Bolt Lock            | 40.6kg |

sales@bristolmaid.com

## **Caretray Trolley - IV Preparation**

### **Features**

- Trolleys designed for use when preparing & administering intravenous infusions
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Infusion Pole & Bracket
  - Waste Bin & Bracket, 8 Litre
  - Glove / Mask Box Holder, One Box
  - · Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking
- When not in use, trolleys storing medicines should be secured to the wall using a Security Clamp (see locking & security accessories)

## • Paint - Grey White

### **Dimensions** (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

+44 (0)1258 484455

- Locking & Security page 60
   Overbridge & Tilt Bins page 66
   Defibrillator Platform page 67
   Carerail Mounted page 68

- Laptop Holders page 74
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTIVPRM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 39.8kg |
| CTIVPRM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 39.5kg |
| CTIVPRM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 40.6kg |



+44 (0)1258 484455



| Item                      | Description   | Weight |
|---------------------------|---|--------|
| CAM Lock with Wing Handle |   |        |
| CTMEDSM1/10/3S2D/CA       | 3 x 100mm, 2 x 150mm Trays, CAM Lock with Wing Handle | 43.2kg |
| CTMEDSM1/10/4S1D/CA       | 4 x 100mm, 1 x 150mm Trays, CAM Lock with Wing Handle | 42.9kg |
| CTMEDSM1/10/6S/CA         | 6 x 100mm Trays, CAM Lock with Wing Handle            | 44.0kg |

CTMEDSM1/10/3S2D/CA

## **Caretray Trolley - Medical / Surgical**

### **Features**

- Trolleys designed for use within the Medical / Surgical Environment
- · Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- · Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - · CAM Lock with Wing Handle Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Overbridge, Mild Steel
  - Four Compartment Tilt Bin, Mild Steel
  - Infusion Pole & Bracket
  - Waste Bin & Bracket, 8 Litre
  - Glove / Mask Box Holder, One Box
  - · Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

### **Standard Colours**

· Paint - Grey White

### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with Overbridge) 600 x 550 x 1620mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

sales@bristolmaid.com +44 (0)1258 484455

- Locking & Security page 60
  Overbridge & Tilt Bins page 66
  Carerail Mounted page 68
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTMEDSM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 43.4kg |
| CTMEDSM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 43.1kg |
| CTMEDSM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 44.0kg |



CTMEDSM1/10/3S2D/EPB

O



## **Caretray Trolley - Plaster**

### **Features**

- Trolleys suitable for use within A&E / Plaster Room
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Trolley pre-configured with the following accessories...
  - Bucket Support & Bracket
  - Stainless Steel Graduated Bucket, 15 Litre
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

### Dimensions (w x d x h)

- Overall Dimensions 1005 x 550 x 993mm
- Plastic Top 915 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

### **Standard Colours**

• Paint - Grey White

| Item             | Description                | Weight |
|------------------|----------------------------|--------|
| CTPLASM2/10/6S4D | 6 x 100mm, 4 x 150mm Trays | 50.6kg |
| CTPLASM2/10/8S2D | 8 x 100mm, 2 x 150mm Trays | 49.9kg |

sales@bristolmaid.com +44 (0)1258 484455

- Carerail Mounted page 68
  Laptop Holders page 74
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77







## **Caretray Trolley - Vaccination**

### **Features**

- Trolleys suitable for use with vaccination
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Trolley pre-configured with the following accessories...
  - Alcohol Gel Holder
  - Debris Bin & Bracket
  - Glove, Mask Box Holder, One Box
  - Sharps Box Holder
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured label strips supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

### **Standard Colours**

### Dimensions (w x d x h)

- Paint Grey White
- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

| Item             | Description                | Weight |
|------------------|----------------------------|--------|
| CTVACIM1/10/3S2D | 3 x 100mm, 2 x 150mm Trays | 32.0kg |
| CTVACIM1/10/4S1D | 4 x 100mm, 1 x 150mm Trays | 31.7kg |
| CTVACIM1/10/6S   | 6 x 100mm Trays            | 32.8kg |

- Locking & Security page 60
   Overbridge & Tilt Bins page 66
   Defibrillator Platform page 67
   Carerail Mounted page 68

- Laptop Holders page 74
  Castors & Ring Buffers page 76
  Trays, Shelves, Dividers & Labels page 77







|          | C |  |
|----------|---|--|
| <b>S</b> |   |  |

| Item             | Description            | Weight |
|------------------|------------------------|--------|
| CTBASIM1/10      | Top Panel              | 14.2kg |
| CTBASIM1/10/1FWS | 1 x Full Width Shelf   | 15.5kg |
| CTBASIM1/10/2FWS | 2 x Full Width Shelves | 16.8kg |

## **Caretray Trolley - Basics**

### **Features**

- Standard level single column easy clean procedure trolleys
- Design allows for the trolley to be upgraded with side panels & travs
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top...
  - Smooth surfaces & radiused edges to assist cleaning (compatible with steam cleaning)
- Raised edges on three sides, providing defined grab points to act as push handles whilst an integral rim at the front contains spillages (max load 30kg)
- · Solid grade laminate full width shelves, removable & supported on the trolley frame (max load per shelf 5kg)
- · Carerails fitted to either side support a wide range of accessories
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Full Width Shelves 420 x 420 (internal) at heights 193 & 554mm

### **Standard Colours**

· Paint - Spearmint

- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
  Laptop Holders page 74

- Castors & Ring Buffers page 76
   Trays, Shelves, Dividers & Labels page 77





CTBASIM1/10/2FWS



### Mobile Workstation - Mild Steel, Low Level

### Features

- Recessed liquid retaining top
- Single Point Locking Using a CAM lock Supplied with 2 keys
- Anti-tilt mechanism, allows only one drawer to be opened at any one time
- 85, 170 & 255mm deep drawers, fitted on full extension glides (max load 30kg per drawer)
- Two tone colour keyway, drawer fronts available in a range of colours
- Low level corner bumpers protect the fabric of the building
- 80mm pressed steel castors, 2 off braking

### Dimensions (w x d x h)

- Overall 615 x 500 x 810mm
- Drawers 480 x 375 x 85 / 170 / 255mm (internal)

### **Standard Colours**

- Body Light Grey
- Drawers Dark Grey

| Item  | Description                            | Weight |
|-------|--|--------|
| MW005 | 2 x 85mm, 1 x 170mm, 1 x 255mm Drawers | 49.0kg |

```
www.bristolmaid.com
```





Single Point Locking Using a CAM lock - supplied with 2 keys)

Low level corner bumpers to protect the fabric of the building



## Procedure Carts - Polymer, Metro Flexline™

#### **Features**

- · Lightweight polymer construction with smooth rounded corners to assist cleaning
- Incorporates Microban<sup>®</sup> offering antimicrobial protection
- Upper work surface provides a writing surface / space for preparation work
- Model specific locking
  - CAM Lock Supplied with 2 keys
  - Keyless Entry Touch Pad Lock
  - Security Seal Lock Security seal compatible
- Full extension lightweight polymer drawers, depths 76, 152 & 229mm (max load 7 - 20kg per drawer)
- Two tone colour keyway, drawer handles available in a range of colours
- 127mm polymer castors, 2 off directional & 2 off total lock

#### **Dimensions** (w x d x h)

- Drawers
  - 76mm 514 x 432 x 57mm (internal)
  - 152mm 514 x 432 x 133mm (internal)
  - 229mm 514 x 432 x 210mm (internal)

#### Standard Colours

- Body Light Taupe
- Drawers Model Dependant

#### **Options / Accessories**

• 8FL190 - Drawer Label Holders (x 10)

| Item      | Description  | Overall Dimensions (w x d x h) | Weight |
|-----------|--|--------------------------------|--------|
| 8FLCRP4   | Emergency - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Code Blue Handles                  | 1056 x 607 x 1273mm            | 91.0kg |
| 8FLPED1   | Paediatric Emergency - 8 x 76mm, 1 x 152mm Drawers, Coloured Handles                   | 1056 x 607 x 1273mm            | 91.0kg |
| 8FLTMENT2 | Treatment - 2 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Slate Blue Handles                 | 819 x 568 x 1064mm             | 91.0kg |
| 8FLCOM-LT | Laptop - 1 x Keyboard Tray, 1 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Slate Blue Handles | 819 x 568 x 1064mm             | 83.0kg |
| 8FLXCLINE | Central Line - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles              | 819 x 568 x 1064mm             | 91.0kg |
| 8FLXCCU1  | Critical Care - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Red Handles                    | 840 x 569 x 1807mm             | 90.0kg |



8FLCRP4



|   | Emergency<br>8FLCRP4 | Paediatric Emergency<br>8FLPED1 | Treatment<br>8FLTMENT2 | Laptop<br>8FLCOM-LT | Central Line<br>8FLCLINE | Critical Care<br>8FLCCU1 |
|---|----------------------|---------------------------------|------------------------|---------------------|--------------------------|--------------------------|
| Locking & Security                      |                      |                                 |                        |                     |                          |                          |
| CAM Lock                                |                      |                                 | $\checkmark$           |                     |                          |                          |
| Security Seal Lock                      | $\checkmark$         | √                               |                        |                     |                          |                          |
| Keyless Entry Touchpad Lock             |                      |                                 |                        | $\checkmark$        | $\checkmark$             | $\checkmark$             |
| Drawer Seal Lock Bar                    |                      | √                               |                        |                     |                          |                          |
| Plastic Security Seals (x 100)          | $\checkmark$         | $\checkmark$                    |                        |                     |                          |                          |
| Overbridge & Tilt Bins                  |                      |                                 |                        |                     |                          |                          |
| Overbridge with Top Shelf & Hanger Rail |                      |                                 |                        |                     |                          | $\checkmark$             |
| Six Compartment Tilt Bin for Overbridge |                      |                                 |                        |                     |                          | $\checkmark$             |
| Defibrillator Platform                  |                      |                                 |                        |                     |                          |                          |
| Adjustable Defibrillator Platform       | $\checkmark$         | $\checkmark$                    |                        |                     |                          |                          |
| Drawer Trays & Dividers                 |                      |                                 |                        |                     |                          |                          |
| 76mm Drawer Divider Kit                 |                      |                                 | 2                      |                     | 2                        |                          |
| 152 / 229mm Drawer Divider Kit          |                      |                                 | 2                      |                     | 1                        |                          |
| 76mm Drawer Divider Tray & Dividers     | 3                    |                                 |                        |                     |                          |                          |
| 152mm Drawer Divider Tray & Dividers    | 2                    |                                 |                        |                     |                          |                          |
| Side Bins                               |                      |                                 |                        |                     |                          |                          |
| Side Bin, Locking                       | 4                    |                                 |                        |                     | 2                        | 4                        |
| Side Bin, Non-Locking                   |                      |                                 |                        | 3                   | 2                        | 4                        |
| Laptop Holders                          |                      |                                 |                        |                     |                          |                          |
| Laptop Arm                              |                      |                                 |                        | $\checkmark$        |                          |                          |
| Accessories                             |                      |                                 |                        |                     |                          |                          |
| Waste Basket (28 Litre)                 |                      |                                 |                        |                     | $\checkmark$             |                          |
| Oxygen Cylinder Holder                  | $\checkmark$         | √                               |                        |                     |                          |                          |
| Lockable Sharps Container               |                      |                                 |                        |                     | $\checkmark$             |                          |
| Cord Manager                            | $\checkmark$         |                                 |                        |                     |                          |                          |
| Hospital Grade 6, Outlet Strip & Holder | $\checkmark$         |                                 |                        |                     |                          |                          |
| Suction Pump Shelf                      | $\checkmark$         | √                               |                        |                     |                          |                          |
| Rear Fitting Backboard                  | $\checkmark$         | √                               |                        |                     |                          |                          |
| Pull Out-Side Shelf                     | $\checkmark$         | √                               |                        |                     |                          |                          |
| Two Hook Infusion Pole                  | $\checkmark$         | $\checkmark$                    |                        |                     |                          |                          |



## Procedure Carts - Polymer, Metro Starsys<sup>™</sup>

#### **Features**

- · Lightweight polymer construction with smooth rounded corners to assist cleaning
- Incorporates Microban® offering antimicrobial protection
- Upper work surface provides a writing surface / space for preparation work
- CAM lock Supplied with 2 keys\*
- Full extension lightweight polymer drawers, depths 76, 152 & 229mm (max load 7 - 20kg per drawer)
- Two tone colour keyway, drawer handles available in a range of colours
- 127mm polymer castors, 2 off braking
- \*model specific

#### Dimensions (w x d x h)

- Drawers
  - 76mm 412 x 400 x 60mm (internal)
  - 152mm 412 x 400 x 136mm (internal)
  - 229mm 412 x 400 x 225mm (internal)

#### Standard Colours

- Body Light Taupe
- Drawers Handles Model Dependant

#### **Options / Accessories**

• 8FL190 - Drawer Label Holders (x 10)



| Item       | Description   | Overall Dimensions (w x d x h) | Weight  |
|------------|---|--------------------------------|---------|
| 8SXRSANES  | Anaesthesia - 4 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock        | 882 x 633 x 1136mm             | 104.0kg |
| 8SXRSDRS   | Dressing - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock | 577 x 633 x 1060mm             | 83.0kg  |
| 8SXRSTRAMA | Trauma - 3 x 76mm, 4 x 152mm Drawers, Slate Blue Handles, CAM Lock              | 882 x 633 x 1136mm             | 116.0kg |
| 8SXRSCCU   | Critical Care - 1 x 76mm, 1 x 152mm, 3 x 229mm Drawers, Slate Blue Handles      | 882 x 633 x 1212mm             | 107.0kg |
| 8SXRSARTH  | Arthroscopic - 6 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock       | 882 x 633 x 1136mm             | 110.0kg |
| 8SXRSIV    | IV Therapy - 5 x 76mm, 2 x 229mm Drawers, Green Handles, CAM Lock               | 882 x 633 x 1136mm             | 106.0kg |
| 8SXRSCAST  | Cast - 2 x 76mm, 2 x 152mm, 2 x 229mm Drawers, Slate Blue Handles               | 882 x 633 x 1212mm             | 109.0kg |



|   | Anaesthesia<br>8SXRSANES | Dressing<br>8SXRSDRS | Trauma<br>8SXRSTRAMA | Critical Care<br>8SXRSCCU | Arthroscopic<br>8SXRSARTH | IV Therapy<br>8SXRSIV | Cast<br>8SXRSCAST |
|---|--------------------------|----------------------|----------------------|---------------------------|---------------------------|-----------------------|-------------------|
| Locking & Security  |                          |                      |                      |                           |                           |                       |                   |
| CAM Lock  |                          |                      |                      | $\checkmark$              |                           |                       | $\checkmark$      |
| Overbridge & Accessory Rails                              |                          |                      |                      |                           |                           |                       |                   |
| Overbridge with One Shelf & Accessory Rail                |                          |                      | ~                    |                           | $\checkmark$              |                       |                   |
| Overbridge with Two Accessory Rails                       |                          |                      | ~                    | $\checkmark$              | $\checkmark$              | $\checkmark$          | $\checkmark$      |
| Drawer Trays & Dividers                                   |                          |                      |                      |                           |                           |                       | $\checkmark$      |
| 76mm Drawer Divider Kit                                   |                          |                      |                      | $\checkmark$              |                           |                       | $\checkmark$      |
| Pods & Tilt-Bins  |                          |                      |                      |                           |                           |                       |                   |
| Left Swing Out, 3 x Tilt-Bins                             |                          |                      |                      | $\checkmark$              |                           |                       | $\checkmark$      |
| Left Swing Out, 3 x Locking Tilt-Bins (CAM Lock)          |                          |                      |                      |                           |                           | $\checkmark$          |                   |
| Left Swing Out Locking, 3 x Locking Tilt-Bins (CAM Lock)  | ~                        |                      |                      |                           |                           |                       |                   |
| Left Swing Out, Top Locking Tilt-Bin (CAM Lock)           |                          |                      | √                    |                           | √                         | $\checkmark$          |                   |
| Right Fixed Pod, 1 x Tilt-Bin                             |                          |                      |                      | $\checkmark$              |                           |                       | √                 |
| Right Swing Out   |                          |                      |                      |                           | √                         |                       |                   |
| Right Swing Out Locking, 3 x Locking Tilt-Bins (CAM Lock) | $\checkmark$             |                      |                      |                           |                           |                       |                   |
| Right Swing Out, Top Locking Tilt-Bin (CAM Lock)          |                          |                      | √                    |                           |                           |                       |                   |
| Middle Non-Locking Tilt-Bin                               |                          |                      | √                    |                           |                           |                       |                   |
| Pod Tray & Dividers                                       |                          |                      | ✓                    | $\checkmark$              | $\checkmark$              | $\checkmark$          | $\checkmark$      |
| Tilt-Bin  |                          |                      |                      |                           | √                         |                       |                   |
| Accessories   |                          |                      |                      |                           |                           |                       |                   |
| Waste Basket & Holder                                     |                          |                      | $\checkmark$         | $\checkmark$              | ✓                         |                       | $\checkmark$      |
| Universal Clamp   | ~                        |                      |                      |                           |                           |                       |                   |
| Laryngoscope Blade / Specimen Holder                      | ~                        |                      |                      |                           |                           |                       |                   |
| Half Size Utility Bin (Metal)                             | $\checkmark$             |                      |                      |                           |                           |                       |                   |
| Full Size Utility Bin (Metal)                             | $\checkmark$             |                      |                      |                           |                           |                       |                   |
| Label / Tape Dispenser (x 3)                              | $\checkmark$             |                      |                      |                           |                           |                       |                   |
| Plastic Utility Bins with Clear Covers (x 6)              | $\checkmark$             |                      |                      | $\checkmark$              | $\checkmark$              |                       |                   |
| Sharps Container Bracket                                  |                          |                      | $\checkmark$         |                           | $\checkmark$              |                       |                   |
| Two Hook Infusion Pole                                    |                          |                      | $\checkmark$         |                           |                           | $\checkmark$          |                   |
| Wire Basket   |                          |                      | $\checkmark$         | $\checkmark$              | $\checkmark$              |                       | $\checkmark$      |



 $\overline{\mathbf{O}}$ 

# Procedure Carts - Metro Basix™

#### **Features**

σ

- · Recessed liquid retaining plastic work surface
- Model specific locking
  - CAM Lock Supplied with 2 keys
  - Push Button Lock
  - Security Seal Lock Security seal compatible
- Full extension drawers, depths 76, 152 & 230mm (max load 7 24kg per drawer)
- Two tone colour keyway, drawer handles available in a range of colours
- Low level bumpers fitted to Trauma & Anaesthesia model
- 127mm polymer castors, 2 off directional & 2 off total lock

#### Dimensions (w x d x h)

- Compact
  - Overall Dimensions 489 x 470 x 876mm
  - Drawers 381 x 432 x 76 / 152 / 230mm (internal)
- Trauma / Anaesthesia
  - Overall Dimensions 762 x 613 x 1440mm
  - Drawers 432 x 565 x 76 / 152 / 230mm (internal)

#### **Standard Colours**

- Body Light Taupe
- Drawers Dark Blue (see below for alternative colours)

#### **Options / Accessories**

• 8MBA172 - Label Holders (x 12)

Alternative drawer colours are available. Change suffix 'DB' on the product code with any of the following i.e. - Compact Basix Cart with red drawers 8MBC1210RE



| $\checkmark$ | Item             | Description  |        |
|--------------|------------------|--|--------|
| $\checkmark$ | 8MBC1210DB       | Compact - 1 x 76mm, 2 x 152mm, 1 x 230mm Drawers     | 63.0kg |
| $\checkmark$ | 8MBX3210TRAUMADB | Trauma - 3 x 76mm, 2 x 152mm, 1 x 230mm Drawers      | 78.0kg |
|              | 8MBX1310ANES4DB  | Anaesthesia - 1 x 76mm, 3 x 152mm, 1 x 230mm Drawers | 81.0kg |





# **Resuscitation Trolleys & Carts**

- Resuscitation Trolleys Mild Steel 316
  - Caretray Trolley Resuscitation 322
- Resuscitation Carts Polymer, Metro Lifeline<sup>™</sup> 324
- Resuscitation Carts Polymer, Metro Flexline<sup>™</sup> 326
- Airways / Difficult Intubation Trolley Mild Steel 327
- Caretray Trolley Airways / Difficult Intubation 328



Drop down cupboard door prevents drawers from being opened. Can be secured using security seals.



# **Resuscitation Trolleys - Mild Steel**

#### **Features**

- Popular resuscitation trolley
- Drawer contents are secured when the lower cupboard is closed (security seal compatible)
- Upper drawer supplied with a plastic drawer liner, with dividers
- Recessed liquid retaining stainless steel work surface (usable 630 x 445)
- Lift up work flap (usable 305 x 455mm)
- Two hook infusion pole mounted to rear of trolley
- Low level bumper bar
- · Push handles one end
- 125mm pressed steel castors

#### Dimensions (w x d x h)

- Drawers 600 x 430 x 70 / 96 / 155mm (internal)
- Lower Cupboard 625 x 475 x 295mm (internal)

#### **Standard Colours**

- · Paint Red
- Lift Up Flap White

| Item  | Description                                       | Overall Dimensions (w x d x h) | Weight |
|-------|---|--------------------------------|--------|
| WT200 | 2 x 70mm & 1 x 155mm Deep Drawers, Lower Cupboard | 825 x 575 x 999mm              | 74.0kg |
| WT210 | 3 x 96mm Deep Drawers, Lower Cupboard             | 825 x 575 x 999mm              | 75.1kg |
| WT250 | 4 x 70mm & 1 x 155mm Deep Drawers, Lower Cupboard | 825 x 575 x 1165mm             | 86.0kg |

www.bristolmaid.com



WT200



# **Adjustable Defibrillator Platform**

#### **Features**

- Compatible with all popular defibrillators
- Defibrillator supported on a 365mm wide platform & is held securely in place using a clamp (adjusts using a handwheel 205 380mm)
- Platform mounted on an arm that swings 90 degrees from the worksurface & can be locked in position using a tensioning knob (rotates 360 degrees)
- Base of platform sits 150mm above the top of the trolley
- Fitted to the rear left-hand side of the trolley
- Maximum capacity 6kg
- Factory fitted

#### Dimensions (w x d x h)

 477 x 365 x 642mm (increases overall height by 185mm)

#### **Standard Colours**

· Paint - Grey White



#### WT/DEFIB

 Item
 Description

 WT/DEFIB
 Adjustable Defibrillator Platform

# **Suction Pump Shelf**

#### **Features**

- Compatible with all popular suction pumps
- Shelf provides an area 415 x 200mm (w x d)
- 30mm upstand at front of shelf
- Floor to base of shelf at height 260mm
- Corner bumpers offer protection to both the trolley & fabric of the building
- Maximum capacity 15kg
- Factory fitted

#### Dimensions (w x d x h)

• Overall Dimensions - 420 x 205 x 380mm

#### **Standard Colours**

Paint - Red



WT/SPS shown on WT/200

| Item   | Description        | Weight |
|--------|--------------------|--------|
| WT/SPS | Suction Pump Shelf | 5.0kg  |

Weight

5.5kg

# **Cardiac Board**

#### Features

- Laminate cardiac board approx. 500 x 400 x 17mm
- Positioned at the rear of the trolley for easy access
- Factory fitted



#### Features

- WT/CYL/D/E compatible with AD, CD, D, DD, PD, RD, ZD & E cylinders, max diameter 100mm
- WT/CYL/F compatible with F, LF, VF, HX & ZX cylinders, max diameter 145mm
- Maximum two-cylinder holders per trolley
- Factory fitted



- Plastic drawer liner with removable dividers
- One tidy supplied as standard with each trolley





| • |  |
|---|--|
|   |  |
|   |  |

| Item      | Description   | Weight |
|-----------|---------------|--------|
| WT/CBOARD | Cardiac Board | 1.5kg  |

| Item       | Description                  | Weight |
|------------|------------------------------|--------|
| WT/CYL/D/E | D/E Size, Max Diameter 100mm | 3.0kg  |
| WT/CYL/F   | F Size, Max Diameter 145mm   | 3.4kg  |

| Item  | Description | Weight |
|-------|-------------|--------|
| WT/DT | Drawer Tidy | 1.6kg  |

| em | Description | Weight |
|----|-------------|--------|



# Four Gang Socket

#### **Features**

- Switched & fused four gang socket
- Supplied with 3 metres of beige cable & fitted with a UK 3 pin plug
- Fitted to the rear of the trolley, complete with a cleat to secure the cable when not in use
- Factory fitted

# **Guard Rail**

#### Features

- Stainless steel guard rail
- Increases the height of the trolley by 60mm
- Factory fitted

# **Security Seals**

#### **Features**

- Snap on plastic securityPack of 100







| Item       | Description      | Weight |
|------------|------------------|--------|
| E/SOCKET/4 | Four Gang Socket | 1.4kg  |

| Item     | Description | Weight |
|----------|-------------|--------|
| G/RAIL-C | Guard Rail  | 1.5kg  |

| Item     | Description    |
|----------|----------------|
| 2PLCO026 | Security Seals |

# **Heavy Duty Castors**

#### Features

- 125mm heavy duty pressed steel castorsOffers additional cushioning, lowering noise levels
- Factory fitted



6

www.bristolmaid.com



# **Caretray Trolley - Resuscitation**

#### **Features**

- Trolleys suitable for use with Resuscitation
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Door fitted with a security seal lock (security seal compatible)
- Trolley pre-configured with the following accessories...
  - Cylinder Holder, C/CD/D
  - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

#### **Options / Accessories**

• Paint - Red

Defibrillator Platform - page 67

| Item                 | Description                | Weight |
|----------------------|----------------------------|--------|
| CTRESUSM1/10/3S2D/SS | 3 x 100mm, 2 x 150mm Trays | 38.3kg |
| CTRESUSM1/10/4S1D/SS | 4 x 100mm, 1 x 150mm Trays | 38.0kg |
| CTRESUSM1/10/6S/SS   | 6 x 100mm Trays            | 39.1kg |

# **Overbridge**

#### Features

- · Fixes onto the rear uprights of the trolley
- Provides a fixing location for tilt bins (order separately)
- Factory fitted

#### Standard Colours

· Paint - Grey White



# **Tilt Bins**

### **Features**

- · Available with 4.5 or 6 compartments
- · Clear bins provide visible storage & assists organisation
- Bins tilt to 45° allowing easy access of contents, & can be removed to assist refilling / cleaning
- · Label holder for easy identification
- · Retaining bar secures bins in place, preventing from tipping forward during transit
- Supplied with a mild steel mounting

#### Dimensions (w x d x h)

TBIN6/MS - 89 x 67 x 95mm

#### Standard Colours

**Compartment Dimensions** 

45° Angle Tilt

Label Holder

+44 (0)1258 484455

- · Plastic Cream
- TBIN4/MS 143 x 127 x 190mm • TBIN5/MS - 108 x 95 x 146mm
- · Paint Grey White



Middle: TBIN5/MS Bottom: TBIN4/MS

| Item     | Description      | Overall Dimensions (w x d x h) | Weight |
|----------|------------------|--------------------------------|--------|
| TBIN4/MS | Four Compartment | 600 x 171 x 207mm              | 2.3kg  |
| TBIN5/MS | Five Compartment | 600 x 135 x 164mm              | 1.6kg  |
| TBIN6/MS | Six Compartment  | 600 x 93 x 113mm               | 0.9kg  |

Weight

2.0kg

www.bristolmaid.com



## Resuscitation Carts - Polymer, Metro Lifeline™

#### **Features**

- Range of pre-configured emergency carts
- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Incorporates Microban<sup>®</sup> offering antimicrobial protection
- Single seal mechanism secures the top compartment, drawers & tiltout side bins (security seal compatible)
- Top storage compartment features a clear removable cover, & provides instant access to first-line equipment
- Strapless defibrillator platform on swing arm adjusts to accommodate various models
- Full extension lightweight polymer drawers, depths 76, 152, 229 & 305mm, (max load 7 20kg per drawer)
- Ergonomic handle with fifth wheel steering provides maximum control, when attending an emergency situation
- Two tone colour keyway
- 127mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

• Drawers - 514 x 432 x 76 / 152 / 229 / 305mm (internal)

#### **Standard Colours**

- · Body Blue
- Drawers Light Taupe / Code Blue

#### **Options / Accessories**

• 8FL190 - Drawer Label Holders (x 10)

| Item     | Description                            | Overall Dimensions (w x d x h) | Weight |
|----------|--|--------------------------------|--------|
| 8LECCRP2 | 1 x 76mm, 3 x 152mm Drawers            | 965 x 575 x 984mm              | 84.0kg |
| 8LECCRP3 | 2 x 76mm, 1 x 152mm, 1 x 229mm Drawers | 965 x 575 x 984mm              | 88.0kg |
| 8LECCRP4 | 2 x 76mm, 1 x 152mm, 1 x 305mm Drawers | 965 x 575 x 1070mm             | 87.0kg |



Drawer Label Holders





|                                   | 8LECCRP2     | 8LECCRP3     | 8LECCRP4     |
|-----------------------------------|--------------|--------------|--------------|
| Top Cavity Tray                   |              | $\checkmark$ |              |
| Locking Side Bin (x 3)            |              | $\checkmark$ | $\checkmark$ |
| Waste Basket (28 Litre)           | $\checkmark$ |              |              |
| Lockable Sharps Container         | $\checkmark$ |              |              |
| Adjustable Defibrillator Platform | $\checkmark$ | $\checkmark$ | $\checkmark$ |
| Suction Pump Shelf                |              | $\checkmark$ |              |
| Backboard (Front Fitting)         |              | $\checkmark$ | $\checkmark$ |
| Backboard (Rear Fitting)          | $\checkmark$ |              |              |
| Oxygen Cylinder Holder            | ✓            | $\checkmark$ | $\checkmark$ |
| Two Hook Infusion Pole            | $\checkmark$ | $\checkmark$ | $\checkmark$ |
| Storage Bin                       | $\checkmark$ | $\checkmark$ | $\checkmark$ |
| Plastic Security Seals (x 100)    | $\checkmark$ | $\checkmark$ | $\checkmark$ |



# **Resuscitation Carts - Polymer, Metro Flexline™**

#### **Features**

- Lightweight polymer construction with smooth rounded corners to
   assist cleaning
- Incorporates Microban<sup>©</sup> offering antimicrobial protection
- Upper work surface provides a writing surface / space for preparation work
- Single seal mechanism secures the drawers (security seal compatible)
- Full extension lightweight polymer drawers, depths 76, 152 & 229mm (max load 7 20kg per drawer)
- Two tone colour keyway
- 127mm polymer castors, 2 x directional & 2 x total lock

#### Dimensions (w x d x h)

- Overall Dimensions 819 x 568 x 914mm
- Drawers
  - 76mm 514 x 432 x 57mm (internal)
  - 152mm 514 x 432 x 133mm (internal)
  - 229mm 514 x 432 x 210mm (internal)

#### **Standard Colours**

- Body Light Taupe
- Drawers Light Taupe / Red

#### **Options / Accessories**

- 8FL141 76mm Egg Crate Style Drawer Divider Kit (10 x Vertical, 5 x Horizontal)
- 8FL142 152mm Egg Crate Style Drawer Divider Kit (10 x Vertical, 5 x Horizontal)
- 8FL190 Drawer Label Holders (x 10)
- 8FL234 100mm Diameter Oxygen Cylinder Holder
- 8FL315 Two Hook Infusion Pole
- 8LEC320 Plastic Security Seals (x 100)
- 8LEC304 Defibrillator Platform

# Airways / Difficult Intubation Trolley - Mild Steel

#### Features

- · Recessed liquid retaining top
- Available as non-locking or fitted with a tamper evident panel to secure the contents of the drawers (security seal compatible)
- 135mm deep drawers, fitted on full extension glides (max load 30kg per drawer)
- Carerails fitted to L/H side support a wide range of accessories
- Low level corner bumpers protect the fabric of the building
- 80mm pressed steel castors, 2 off braking

#### Dimensions (w x d x h)

- Overall 615 x 500 x 915mm
- Drawers 480 x 375 x 135mm (internal)

#### Standard Colours

· Paint - White

#### **Options / Accessories**

Carerail Mounted - page 68



MW015 Tamper Evident Panel



| MW010 |  |
|-------|--|
|       |  |
|       |  |

| Item  | Description                             | Weight |
|-------|---|--------|
| MW010 | 5 x 125mm Drawers, Non Locking          | 55.0kg |
| MW015 | 5 x 125mm Drawers, Tamper Evident Panel | 57.0kg |

**Resuscitation Trolleys** 

ço



# **Caretray Trolley - Airways / Difficult Intubation**

#### **Features**

- Trolleys suitable for use with Airways / Difficult Intubation
- $\bullet\,$  Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Door fitted with a security seal lock (security seal compatible)
- Trolley pre-configured with the following accessories...
  - Alcohol Gel Holder
  - Airway Tube Holder
  - Glove / Mask Box Holder, One Box
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking

#### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

• Paint - Grey White

| Item                | Description                | Weight |
|---------------------|----------------------------|--------|
| CTAIRWM1/10/3S2D/SS | 3 x 100mm, 2 x 150mm Trays | 38.6kg |
| CTAIRWM1/10/4S1D/SS | 4 x 100mm, 1 x 150mm Trays | 38.3kg |
| CTAIRWM1/10/6S/SS   | 6 x 100mm                  | 39.4kg |

sales@bristolmaid.com

+44 (0)1258 484455

#### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67
- Carerail Mounted page 68
- Laptop Holders page 74
  Castors & Ring Buffers page 76
- Trays, Shelves, Dividers & Labels page 77





# **Sackholders & Waste Bins**

- Removable Body Bins 320
- Free Standing Sackholders 325
- Fixed Body Plastic Sackholders Plastic 326
- Removable Body Plastic Sackholders Plastic 327



#### 20 Litre

| Item      | Description                            | Weight |
|-----------|--|--------|
| BR015/WHL | Hands Free                             | 8.5kg  |
| BR020/WHL | Hands Free & Silent Closing            | 8.5kg  |
| BR030/WHL | Hands Free, Rust Free & Silent Closing | 8.5kg  |

# Removable Body Bins - 20 / 50 Litre

#### **Features**

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- · Captive drip tray to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Choice of construction...
  - Mild steel tubular frame & sheet metal body, painted for extra durability
  - 'Rust Free' stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise when the lid closes

#### Dimensions (w x d x h)

- 20 Litre 267 x 346 x 552mm
- 50 Litre 444 x 397 x 551mm

#### **Standard Colours**

· Paint - Grey / White





'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)

www.bristolmaid.com



Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Removable Body Bin - 50ltr, Hands free with a green lid - BR045/GNL



#### 50 Litre

| Item      | Description                            | Weight |
|-----------|--|--------|
| BR045/WHL | Hands Free                             | 11.5kg |
| BR050/WHL | Hands Free & Silent Closing            | 11.5kg |
| BR060/WHL | Hands Free, Rust Free & Silent Closing | 11.5kg |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com



#### BR075/WHL

#### 75 Litre

| Item      | Description                            | Weight |
|-----------|--|--------|
| BR075/WHL | Hands Free                             | 15.0kg |
| BR080/WHL | Hands Free & Silent Closing            | 15.0kg |
| BR090/WHL | Hands Free, Rust Free & Silent Closing | 15.0kg |

# Removable Body Bins - 75 / 90 Litre

#### **Features**

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- · Captive drip tray to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- · Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Choice of construction...
  - Mild steel tubular frame & sheet metal body, painted for extra durability
  - 'Rust Free' stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise when the lid closes

#### Dimensions (w x d x h)

- 75 Litre 444 x 397 x 781mm
- 90 Litre 444 x 397 x 910mm

#### **Standard Colours**

· Paint - Grey / White





'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)



BR105/WHL

Sackholders & Waste Bins

Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Removable Body Bin - 75ltr, Hands free with a yellow lid - BR075/YEL



#### 90 Litre

| Item      | Description                            | Weight |
|-----------|--|--------|
| BR105/WHL | Hands Free                             | 17.0kg |
| BR110/WHL | Hands Free & Silent Closing            | 17.0kg |
| BR120/WHL | Hands Free, Rust Free & Silent Closing | 17.0kg |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com






Positive bag retention through the use of a drop over clamp



Gas damper used on silent closing bins



Raised flange & recessed lid used on hands free bins

# **Removable Body Bins - 75 Litre, Front Opening**

#### **Features**

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- Captive drip trav to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- · Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Choice of construction...
  - · Mild steel tubular frame & sheet metal body, painted for extra durability
  - 'Rust Free' stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise when the lid closes

#### Dimensions (w x d x h)

75 Litre - 444 x 397 x 781mm

#### Standard Colours

· Paint - Grev / White





'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)

Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Removable Body Bin - 75ltr. Hands Free & Silent Closing with a grey lid - BR140/GYL



BR135/WHL Hands Free BR140/WHL Hands Free & Silent Closing BR150/WHL Hands Free, Rust Free & Silent Closing

Description

| sales | @bristo | Imaid.com |  |
|-------|---------|-----------|--|

+44 (0)1258 484455

```
www.bristolmaid.com
```

Item

# Free Standing Sackholders - 50 / 75 Litre

#### Features

- · Captive drip tray to collect spills
- · Positive bag retention through the use of an elasticated cord
- Stainless steel pedal & linkages
- · Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Mild steel tubular frame, painted for extra durability

#### Dimensions (w x d x h)

- 50 Litre 442 x 397 x 551mm
- 75 Litre 442 x 397 x 780mm

#### Standard Colours

· Paint - Grev



Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Free Standing Bag Holder - 50 Litre, Elastic Cord with brown lid - BH015/BNL



| Item      | Description                         | Weight |
|-----------|-------------------------------------|--------|
| BH015/WHL | Free Standing Sackholder - 50 Litre | 7.0kg  |
| BH030/WHL | Free Standing Sackholder - 75 Litre | 7.5kg  |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com



# Fixed Body Sackholders - Plastic, 20 Litre

#### **Features**

- Meets NHS England 'Managing Healthcare Fire Safety' (HTM-05-01)
- Plastic construction with no metal parts
- Smooth surfaces to assist cleaning
- One-piece plastic lid, available in a range of colours to assist refuse identification (HTM07)
- Soft closing lid
- Bag positively secured using a drop over sack retention ring

#### Dimensions (w x d x h)

Overall Dimensions - 350 x 220 x 660mm

#### Standard Colours

- · Body White
- · Lid White

#### Options

- Alternative Lid Colours Black / Orange / Yellow / Yellow & Black
- · Moulded In Graphics Available on request

Alternative lid colours available. Change suffix 'WH' on the item no. with any of the following, i.e. Fixed Body Sackholder - Plastic, Silent Closing, 20 Litre with a black lid - 5FSH20BL



```
www.bristolmaid.com
```

Item

# Sackholders & Waste Bins

# Removable Body Sackholders - Plastic, 42 & 70 Litre

#### **Features**

- Meets NHS England 'Managing Healthcare Fire Safety' (HTM-05-03)
- One-piece removable body
- Plastic construction with no metal parts
- Smooth surfaces to assist cleaning
- One-piece plastic lid, available in a range of colours to assist refuse identification (HTM07)
- Soft operating lid on both opening & closing (moving parts enclosed)
- · Recessed drip tray to collect unwanted spills
- Bag positively secured using a drop over sack retention ring
- Non-marking rubber feet prevents the bin from sliding
- Wheels fitted to rear of bin assist moving when cleaning

#### Dimensions (w x d x h)

• 42 Litre - 460 x 420 x 580mm

• 70 Litre - 460 x 450 x 790mm

#### **Standard Colours**

- · Body White
- · Lid White

#### Options

- Alternative Lid Colours Black/Orange/Yellow/Yellow & Black
- Moulded In Graphics Available on request



Alternative lid colours available. Change suffix 'WH' on the item no. with any of the following, i.e. Removable Body Sackholder - Plastic, Silent Closing, 42 Litre with a black lid - 5RSH42BL



| Item     | Description   | Weight |
|----------|---|--------|
| 5RSH42WH | Removable Body Sackholder - Plastic, Silent Closing, 42 Litre | 8.0kg  |
| 5RSH70WH | Removable Body Sackholder - Plastic, Silent Closing, 70 Litre | 9.0kg  |

sales@bristolmaid.com



# **Seating Solutions**

- Seat, Back Cushion & Headrest Options 356
  - Overnight Chair 358
    - Soft Seating 360
  - Visitor & Waiting Room Seating 362
    - Waiting Room Seating 364
      - Tech Chairs & Stools 370
    - PU Tech Chairs & Stools 372
  - Static Safe Tech Chairs & Stools 374
- Static Safe & Sterile Tech Chairs & Stools 376
- Medical Grade Clean Room Tech Chairs 378

- Office, Task & Operators Chairs 330
- 24 Hour Task & Operator Chairs 332
  - Bariatric Chairs 334
  - Saddle & Perch Stools 336
  - Low Back Patient Chairs 338
- Medium & High Back Patient Armchairs 339
  - Reclining Armchairs 340
  - Reclining Armchairs, Bariatric 350
  - Reclining Porters Armchairs 352
  - Rise & Recline Porters Armchairs 354





## **Office, Task & Operators Chairs - Century**

#### **Features**

- Available as...
  - · Medium back
  - High back
  - Draughtsman
- $\bullet$  Medium  $\overset{-}{\&}$  high back available with fixed ring arms & non-braking castors
- Draughtsman has a medium back, foot ring & glides on a 27" base
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Two lever mechanism
- · Black polypropylene five star base
- · Pneumatic seat & back rest height adjustment
- Fully upholstered seat & back with deep CMHR foam
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Medium Back Chair in Aztec - 5160/FAC

#### Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage





| ltem  | Description                     | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|-------|---------------------------------|--------------------------------|------------------------|--------|
| 5160  | Century, Medium Back            | 670 x 670 x 800-1040mm         | 440 - 570mm            | 10.0kg |
| 5161  | Century, Medium Back, Ring Arms | 670 x 670 x 800-1040mm         | 440 - 570mm            | 12.0kg |
| 5162  | Century, High Back              | 670 x 670 x 940-1070mm         | 440 - 570mm            | 11.0kg |
| 5163  | Century, High Back, Ring Arms   | 670 x 670 x 980-1110mm         | 440 - 570mm            | 13.0kg |
| 5160D | Century, Draughtsman            | 460 x 605 x 1010-1200mm        | 536 - 686mm            | 12.0kg |

FIRE



## 24 Hour Task & Operator Chairs - Eon

#### **Features**

- BS 5459: Part 2: 2000. Specification for performance requirements & tests for office furniture. Office pedestal seating for use by persons weighing up to 150kg & for use up to 24 hours a day
   S1 level tested (severe) for heavy duty use
- Certified for use up to 150kg
- Available as...
  - Medium back
  - High back
- 24/7 Flax fire retardant heavy duty fabric upholstery
- Black polypropylene five-star base
- Easy lift back, ratchet height adjustable
- Independent seat & back mechanism with tension control
- Fully upholstered seat & back with deep CMHR foam
- Optional
  - Height adjustable, fold back or ring arms
  - 12" or 14" draughtsman's kit
- Maximum capacity 150kg

Replace suffix on the product code with any of the following to select your colour i.e. Eon, Medium Back in Zone - 5260/FYR09

#### Fabric - Camira 24/7 Flax



| Item | Description      | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|------|------------------|--------------------------------|------------------------|--------|
| 5260 | Eon, Medium Back | 710 x 710 x 835-925mm          | 470 - 600mm            | 13.0kg |
| 5262 | Eon, High Back   | 710 x 710 x 925-1030mm         | 470 - 600mm            | 14.0kg |

## 24 Hour Task & Operator Chairs - Zenith

#### **Features**

- BS 5459: Part 2: 2000. Specification for performance requirements & tests for office furniture. Office pedestal seating for use by persons weighing up to 150kg & for use up to 24 hours a day
- S1 level tested (severe) for heavy duty use
- Certified for use up to 150kg
- 24/7 Flax fire retardant heavy duty fabric upholstery
- Black polypropylene five-star base
- Easy lift back, ratchet height adjustable
- · Independent seat & back mechanism with tension control & seat slider
- Fully upholstered seat & back with deep CMHR foam
- Chrome swivel height adjustable arms & seat slider
- Maximum capacity 150kg

#### **Options**

- 5/Z/HR Headrest
- 5/Z/PLS Pump-Up Lumbar Support

Replace suffix on the product code with any of the following to select your colour i.e. Zenith, High Back in Vernal - 59003/FYR02

#### Fabric - Camira 24/7 Flax





59003/FYR02/5/Z/HR/5/Z/PLS shown with optional headrest & pump-up lumbar support

|       | Description       | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|-------|-------------------|--------------------------------|------------------------|--------|
| 59003 | Zenith, High Back | 710 x 710 x 1160-1240mm        | 450 - 530mm            | 21.0kg |

sales@bristolmaid.com

+44 (0)1258 484455 www.bristolmaid.com







## **Bariatric High Back Chairs**

#### **Features**

- Deep CMHR foam
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- · Black heavy duty steel five star base
- · Black heavy duty gas lift
- Black height adjustable soft touch armrests
- Easy lift ratchet height adjustable back
- · Independent seat & back mechanism with side tension control
- Maximum capacity 315kg

Replace suffix on the product code with any of the following to select your colour i.e. Bariatric High Back in Lime - 5660XL/VLI

#### Vinyl - Ambla New Richmond in Colour



FBK

ItemDescriptionOverall Dimensions (w x d x h)Seat Height AdjustmentWeight5660XLHigh Back, Swivel750 x 750 x 1080-1190mm425 - 535mm29.0kg

FCT

FGE

## **Bariatric Chairs**

#### **Features**

- Deep CMHR foam
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- · Fabric Fire retardant
- Beech, varnished four leg frame
- Maximum capacity 315kg

Alternative colours are available at extra cost. Replace suffix on the product code with any of the following to select your colour i.e. Bariatric, Beech Frame in Lime - 5671/VLI

#### Vinyl - Ambla New Richmond in Colour







| 1 | FIRE<br>RETARDANT |
|---|-------------------|
|   | <u> </u>          |
|   |                   |
|   |                   |

| Item | Description            | Overall Dimensions (w x d x h) | Seat Height | Weight |
|------|------------------------|--------------------------------|-------------|--------|
| 5671 | Bariatric, Beech Frame | 844 x 580 x 890mm              | 490mm       | 16.0kg |



## Saddle Stools

#### **Features**

- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Black polypropylene five star base
- Black gas lift with pneumatic seat height adjustment
- Free wheeling castors
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Saddle Stool in Lime - 5202/VLI

Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage



## **PU Saddle & Perch Stools**

#### Features

- Soft touch black polyurethane seatBlack polypropylene five star base
- Black gas lift with pneumatic seat height adjustment
  Free wheeling castors
- Maximum capacity 120kg



|      | Description     | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|------|-----------------|--------------------------------|------------------------|--------|
| 5201 | PU Saddle Stool | 670 x 670 x 590-780mm          | 590 - 780mm            | 8.0kg  |
| 5205 | PU Perch Stool  | 670 x 670 x 650-730mm          | 510 - 690mm            | 8.0kg  |





## Low Back Patient Chairs - Alderbury

#### **Features**

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- Ergonomic design to ensure correct sit & posture
- · Handgrips assist positioning from standing & seated positions
- Solid beech frame (PEFC accredited)
- · Mortice & tenon joint construction, glued & pinned
- Frame finished with anti-bacterial varnish
- · Plastic glides fitted to chair legs
- Contoured seat & back, seat section fixed
- Hygiene gap prevents the trapping of dirt & bacteria between the seat & back, provides access for cleaning
- Fire retardant foam (BS 5852: Part 2: 1982)
- Upholstered with fire retardant & antimicrobial vinyl
- Maximum capacity 160kg

#### **Standard Colours**

Vinyl - Delft

#### Dimensions (w x d x h)

- Overall 580 x 540 x 890mm
- Seat 470 x 430 at Height 470mm
- Width Between Arms 470mm
- Arm Height 655mm

Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Alderbury, Low Back Chair in Baltic - 5ND/LC/ALDE/VVTUR

Vinyl - Spradling Valencia C5



 Item
 Description
 Weight

 5ND/LC/ALDE/VVDEL
 Alderbury, Low Back Chair, Fixed Seat
 8.2kg

sales@bristolmaid.com

# Medium & High Back Patient Armchairs - Alderbury

#### **Features**

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- · Ergonomic design to ensure correct sit & posture
- Handgrips assist positioning from standing & seated positions
- · Solid beech frame (PEFC accredited)
- Mortice & tenon joint construction, glued & pinned
- · Frame finished with anti-bacterial varnish
- Plastic glides fitted to chair legs
- · Contoured seat & back
  - Fixed Seat is fixed to the frame
  - Removable Seat can be removed to assist cleaning (smooth vinyl upholstered underside) & is reversible to allow the seat section to be repositioned to provide additional clearance at the front of the chair
- Hygiene gap prevents the trapping of dirt & bacteria between the seat & back, provides access for cleaning
- Fire retardant foam (BS 5852: Part 2: 1982)
- · Upholstered with fire retardant & antimicrobial vinyl
- · Maximum capacity 160kg

#### Dimensions (w x d x h)

- Medium Back 660 x 650 x 1000mm
- High Back 660 x 650 x 1195mm
- Seat 450 x 435/475\* at Height 500mm
- Width Between Arms 530mm
- Arm Height 655mm

\*Seat depth can be altered by repositioning the removable seat pad

#### **Standard Colours**

• Vinyl - Delft

#### Options

• HW - Housekeeping Wheels



Optional house keeping wheels assist movement of the chair during cleaning



5ND/HAR/ALDE/VVDEL

www.bristolmaid.com

Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Alderbury, High Back Armchair, Removable & Reversible Seat in Olive - 5ND/HAR/ALDE/VVTUR

Vinyl - Spradling Valencia C5



| Item               | Description                                     |        |
|--------------------|---|--------|
| Medium Back        |   |        |
| 5ND/MAF/ALDE/VVDEL | Alderbury, Medium Back Armchair, Fixed Seat     | 13.6kg |
| 5ND/MAR/ALDE/VVDEL | Alderbury, Medium Back Armchair, Removable Seat | 13.6kg |
| High Back          |   |        |
| 5ND/HAF/ALDE/VVDEL | Alderbury, High Back Armchair, Fixed Seat       | 14.4kg |
| 5ND/HAR/ALDE/VVDEL | Alderbury, High Back Armchair, Removable Seat   | 14.4kg |





## **Reclining Armchair - Shapwick**

#### Features

- Available as either manual or electric operation
  - Manual Tilt in space recliner
  - Electric Single tilt in space reclining only or single tilt in space reclining with rise function
- · Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
  - Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting (fixed seat)
- Foam leg rest
- 4 x 50mm swivel braking castors
- Maximum capacity 127kg

## Dimensions (w x d x h)

- Overall Dimensions 737 x 902 x 1130mm
- Seat 510 x 510mm at Height 460mm

## **Standard Colours**

- Vinyl Delft
- VP Fabric Navy Blue

## **Options**

- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements

# ItemDescriptionWeight5RE/RAM/SHAP/VVDELReclining Armchair, Manual41.0kg5RE/RAE/SHAP/VVDELReclining Armchair, Electric55.0kg5RE/RRAE/SHAP/VVDELRise & Reclining Armchair, Electric57.0kg



Alternative colours available at additional cost.

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Shapwick, Reclining Armchair, Electric in Olive Vinyl & Green Vapour Permeable (VP) Multi-Way Stretch Fabric - 5RE/RAE/SHAP/VVOLI/VPGR



Fabric - VP Multi-Way Stretch



| 5RE/RRA | AE/SHA | P/VVE | DEL |
|---------|--------|-------|-----|
|         |        |       |     |

| Item                     | Description                         |
|--------------------------|-------------------------------------|
| 5RE/RAM/SHAP/VVDEL/VPNB  | Reclining Armchair, Manual          |
| 5RE/RAE/SHAP/VVDEL/VPNB  | Reclining Armchair, Electric        |
| 5RE/RRAE/SHAP/VVDEL/VPNB | Rise & Reclining Armchair, Electric |

www.bristolmaid.com

sales@bristolmaid.com

41.0kg 55.0kg 57.0kg





5RE/RAE/ALWE/VVDEL (fitted with optional pocket)

| Item                | Description                                   | Weight |
|---------------------|---|--------|
| 5RE/RAM/ALWE/VVDEL  | Reclining Armchair, Wide Arm, Manual          | 41.0kg |
| 5RE/RAE/ALWE/VVDEL  | Reclining Armchair, Wide Arm, Electric        | 55.0kg |
| 5RE/RRAE/ALWE/VVDEL | Rise & Reclining Armchair, Wide Arm, Electric | 57.0kg |

# **Reclining Armchair - Alweston, Wide Arm**

#### **Features**

- Ideal for use in oncology, dialysis or haematology for intravenous cannulation, where a patient may be seated for a prolonged period of time, or alternatively in recovery following minor treatment
- Wide supportive cushioned arm (150mm)
- Available as either manual or electric operation
  - · Manual Tilt in space recliner
  - Electric Single tilt in space reclining only, or single tilt in space reclining with rise function
- · Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- · Upholstered with either
  - Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting (fixed seat)
- · Foam leg rest
- 4 x 50mm swivel braking castors
- Maximum capacity 127kg

#### Dimensions (w x d x h)

#### **Standard Colours**

- Overall Dimensions 737 x 902 x 1130mm
  Seat 510 x 510mm at Height 460mm
- Vinyl Delft
- VP Fabric Navy Blue

#### **Options**

- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Table with Wireless Charging See page 327
- · Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements
- · USB Charging Socket Located on outside of arm
- · Pocket Positioned on the right-hand side of the chair



Delivery to Chosen Location
Unwrapping & Removal of Packaging
Product Positioning & Set-Up
Full Demonstration of Use to Staff

5RE/RAE/ALWE/VVLEI/VPBE (fitted with optional pocket)

Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VVNB' i.e. Ashington, Reclining Armchair, Waterfall Back, Wide Arm, Electric in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RAE/ASHI/VVOLI/VPGR

VINII - Spradling Valencia C5

Fabric - VP Multi-Way Stretch



| 5RE/RAE/ALWE/VVBAL/VPGY  |   |        |
|--------------------------|---|--------|
| Item                     | Description                                   | Weight |
| 5RE/RAM/ALWE/VVDEL/VPNB  | Reclining Armchair, Wide Arm, Manual          | 41.0kg |
| 5RE/RAE/ALWE/VVDEL/VPNB  | Reclining Armchair, Wide Arm, Electric        | 55.0kg |
| 5RE/RRAE/ALWE/VVDEL/VPNB | Rise & Reclining Armchair, Wide Arm, Electric | 57.0kg |





| Item                | Description                                       | Weight |
|---------------------|---|--------|
| 5RE/RAM/TWYF/VVDEL  | Reclining Armchair, Waterfall Back, Manual        | 43.0kg |
| 5RE/RAE/TWYF/VVDEL  | Reclining Armchair, Waterfall Back, Electric      | 57.0kg |
| 5RE/RRAE/TWYF/VVDEL | Rise & Recline Armchair, Waterfall Back, Electric | 59.0kg |

## **Reclining Armchair - Twyford, Waterfall Back**

#### **Features**

- Available as either manual or electric operation
  - Manual Tilt in space recliner
  - Electric Single tilt in space reclining only or single tilt in space reclining with rise function
- Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V  $\sim,$  50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
  - · Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- · Easiform reflex foam provides comfort when sitting (fixed seat)
- · Foam leg rest
- Fixed deep padded waterfall back
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

#### Dimensions (w x d x h)

- Overall Dimensions 765 x 865 x 1145mm
- Seat 520 x 510mm at Height 482mm

#### **Options**

- Seat Pressure Management Coolform, Viscoform, Liquiform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements

## **Standard Colours**

- Vinyl Delft
- VP Fabric Navy Blue



- Delivery to Chosen Location
- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff



Alternative colours available at additional cost.

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Twyford, Reclining Armchair, Electric in Olive Vinyl & Green Vapour Permeable (VP) Multi-Way Stretch Fabric - 5RE/RAE/TWYF/VVOLI/VPGR





VPBE VPNB VPGY VPBL VPBU VPGR

Fabric - VP Multi-Way Stretch

| Item                     | Description                                       | Weight |
|--------------------------|---|--------|
| 5RE/RAM/TWYF/VVDEL/VPNB  | Reclining Armchair, Waterfall Back, Manual        | 43.0kg |
| 5RE/RAE/TWYF/VVDEL/VPNB  | Reclining Armchair, Waterfall Back, Electric      | 57.0kg |
| 5RE/RRAE/TWYF/VVDEL/VPNB | Rise & Recline Armchair, Waterfall Back, Electric | 59.0kg |

sales@bristolmaid.com

+44 (0)1258 484455

5RE/RRAE/TWYF/VVLIC/VPBU



- Ideal for use in oncology, dialysis or haematology for intravenous cannulation, where a patient may be seated for a prolonged period of time, or alternatively in recovery following minor treatment
- Wide supportive cushioned arm (150mm)
- · Available as either manual or electric operation
  - Manual Tilt in space recliner
  - Electric Single tilt in space reclining only, or single tilt in space reclining with rise function
- · Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- · Upholstered with either
  - · Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- · Easiform reflex foam provides comfort when sitting (fixed seat)
- · Foam leg rest
- Fixed deep padded waterfall back
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

#### Dimensions (w x d x h)

- Overall Dimensions 765 x 865 x 1145mm
- Seat 520 x 510mm at Height 482mm

#### **Standard Colours**

- Vinyl Delft
- VP Fabric Navy Blue

#### **Options**

- Seat Pressure Management Coolform, Viscoform, Liquiform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- · Table with Wireless Charging See page 327
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements
- USB Charging Socket Located on outside of arm
- · Pocket Positioned on the right-hand side of the chair



5RE/RAM/ASHI/VVDELReclining Armchair, Waterfall Back, Wide Arm, Manual43.0kg5RE/RAE/ASHI/VVDELReclining Armchair, Waterfall Back, Wide Arm, Electric57.0kg5RE/RRAE/ASHI/VVDELRise & Recline Armchair, Waterfall Back, Wide Arm, Electric59.0kg

sales@bristolmaid.com



- Delivery to Chosen Location
  Unwrapping & Removal of Packaging Product Positioning & Set-Up
- Full Demonstration of Use to Staff

#### 5RE/RAE/ASHI/VVSCH/VPGY (fitted with optional pocket)

**Seating Solutions** 

5RE/RAE/ASHI/VVROT/VPNB

Alternative colours available at additional cost.

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Alweston, Reclining Armchair, Wide Arm, Electric in Olive Vinyl & Green Vapour Permeable VP Multi-Way Stretch Fabric - 5RE/RAE/ALWE/VVOLI/VPGR



Fabric - VP Multi-Way Stretch



| Item                     | Description   | Weight |
|--------------------------|---|--------|
| 5RE/RAM/ASHI/VVDEL/VPNB  | Reclining Armchair, Waterfall Back, Wide Arm, Manual        | 43.0kg |
| 5RE/RAE/ASHI/VVDEL/VPNB  | Reclining Armchair, Waterfall Back, Wide Arm, Electric      | 57.0kg |
| 5RE/RRAE/ASHI/VVDEL/VPNB | Rise & Recline Armchair, Waterfall Back, Wide Arm, Electric | 59.0kg |

sales@bristolmaid.com





sales@bristolmaid.com

## **Reclining Armchair - Beckington, Waterfall Back**

## **Features**

- Electric reclining only, or dual tilt in space reclining with rise function
- Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Upholstered with flame retardant vinyl
- · Easiform reflex foam provides comfort when sitting
- · Foam leg rest
- · Fixed waterfall back
- · Side pocket & loop for handset
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

## Dimensions (w x d x h)

- Overall Dimensions 840 x 930 x 1160mm
- Seat 510 x 510mm at Height 483mm
- Back Height 710mm
- Arm Height 180mm

#### Standard Colours

· Vinyl - Ultra



- Delivery to Chosen Location
- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- · Full Demonstration of Use to Staff

Alternative colours available at additional cost. Change suffix 'VMULT' on the product code with any of the following i.e. Beckington, Reclining Armchair, Waterfall Back, Electric in Garden - 5RE/RAE/BECK/VMGAR

#### **Memphis Faux Leather**



| Item                | Description                                       | Weight |
|---------------------|---|--------|
| 5RE/RAE/BECK/VMULT  | Reclining Armchair, Waterfall Back, Electric      | 58.0kg |
| 5RE/RRAE/BECK/VMULT | Rise & Recline Armchair, Waterfall Back, Electric | 62.0kg |

+44 (0)1258 484455





## Reclining Armchair - Hadleigh, Waterfall Back, Bariatric





#### **Features**

- Specifically designed for comfort & pressure management control
- Dual standard rise & recliner, provides a full range of posture positive recline positions as well as support to a standing position
- Mains powered actuator (230V ~. 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Upholstered with a combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting
- · Seat & leg sections generously padded with memory foam & fitted with a removable cover
- · Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- Back blockers provided for additional lumbar support
- Flat wide arms for deep comfort & support
- · Loop for handset & deep side pocket
- Heavy duty glides
- Maximum capacity 254kg

#### Dimensions (w x d x h)

- Chair 1000 / 1100 x 940 x 1145mm
- Seat 610 / 710mm x 510 x 460mm

#### Standard Colours

- · Vinvl Delft
- VP Fabric Navy Blue

#### **Options**

- Seat Pressure Management Coolform, Viscoform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- · USB Charging Socket Located on outside of arm
- · Full Length Arm Caps Positioned over the arm rests

5RE/RRAE/HADL610/VVDEL/VPNB Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 610mm 88.0kg 5RE/RRAE/HADL710/VVDEL/VPNB Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 710mm 88.0kg



- Delivery to Chosen LocationUnwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff

Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Hadleigh, Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 710mm in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RRAE/HADL710/VVOLI/VPGR



Fabric - VP Multi-Way Stretch



SRE/RAE/HADL610/VVROT/VPBU

## **Reclining Porters Armchair - Bathampton, Waterfall Back**





5RE/RPAH/BATH/VVDEL Reclining Porters Armchair, Waterfall Back, Electric 75.0kg

#### **Features**

- Single tilt in space recliner, provides a full range of posture positive recline positions
- Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset (positioned at rear of chair/optional side of chair)
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
  - · Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multi-way stretch fabric. VP fabric applied to seat & back sections
- Upholstered seat & back sections are fitted with removable covers, allows seat section to be easily upgraded to aid pressure care
- · Removable internal Easiform reflex seat cushion
- Seat & leg sections generously padded with memory foam & fitted with a removable cover
- Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- Push handle at rear of the chair to assist when positioning/portering
- Slide out footplate to support feet
- 4 x 75mm swivel braking castors
- Maximum capacity 127kg

#### Dimensions (w x d x h)

- Overall Dimensions 685 x 990 x 1220mm
- Seat 520 x 510mm at Height 470mm

#### **Standard Colours**

- Vinyl Delft
- VP Fabric Navy Blue

#### Options

- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform, Airform, Dynaform
- Back Cushions Comfort Lateral, Support Lateral, Large Profile with Waterfall, Large Profile Headrest with Comfort Lateral, Large Profile Headrest with Support Lateral
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements



Delivery to Chosen LocationUnwrapping & Removal of Packaging Product Positioning & Set-Up Full Demonstration of Use to Staff

5RE/RPAH/BATH/VVOLI/VPGR

Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Bathampton, Waterfall Back, Handset Operated in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RPAH/BATH/VVOLI/VPGR





|                          | Description  | Weight |
|--------------------------|--|--------|
| 5RE/RPAH/BATH/VVDEL/VPNB | Reclining Porters Armchair, Waterfall Back, Electric | 75.0kg |

VPGY

VPGR





#### **Features**

- Dual tilt in space rise & recliner, provides a full range of posture positive recline positions as well as support to a standing position
- Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset (positioned at side of chair)
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- · Upholstered with either
  - Flame retardant antimicrobial vinyl throughout
  - A combination of flame retardant antimicrobial vinvl & vapour permeable (VP) multi-way stretch fabric. VP fabric applied to seat & back sections
- Upholstered seat & back sections are fitted with removable covers, allows seat section to be easily upgraded to aid pressure care
- Removable internal Easiform reflex seat cushion
- · Seat & leg sections generously padded with memory foam & fitted with a removable cover
- · Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- · Push handle at rear of the chair to assist when positioning/portering
- · Optional sliding footplate to support feet
- 4 x 100mm swivel braking castors
- Maximum capacity 127kg

#### Dimensions (w x d x h)

- Overall Dimensions 775 x 1055 x 1245mm
- Seat 510 x 510mm at Height 480mm

#### Standard Colours

- Vinyl Delft
- VP Fabric Navy Blue

#### **Options**

- Sliding Footplate
- Seat Pressure Management Coolform, Viscoform, Liquiform. Celliform, Airform, Dynaform
- Back Cushions Comfort Lateral, Support Lateral, Large Profile with Waterfall, Large Profile Headrest with Comfort Lateral, Large Profile Headrest with Support Lateral
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile. Midline
- · Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements

5RE/RRPAH/DORC/VVDEL Rise & Recline Porters Armchair, Waterfall Back, Electric 76.5kg



Delivery to Chosen Location
Unwrapping & Removal of Packaging
Product Positioning & Set-Up
Full Demonstration of Use to Staff

5RE/RRPAH/DORC/VVBAL/VPNB (with optional footplate)

11

Alternative colours available at additional cost.

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Dorchester, Waterfall Back, Handset Operated in Leinen Vinyl & Grey VP Multi-Way Stretch Fabric - 5RE/RPAH/DORC/VVLEI/VPGY



Fabric - VP Multi-Way Stretch



5RE/RRPAH/DORC/VVLEI/VPNB (with optional sliding footplate)


## Description W

5RE/RRPAH/DORC/VVDEL/VPNB Rise & Recline Porters Armchair, Waterfall Back, Electric 76.5kg

## Seat Options For Pressure Management

Cushions are covered in multi-way stretch, vapour permeable fabric which is wipeable, removable, washable at 95°C, incontinence proof & breathable.

supported comfort.



for periods longer than one hour, require pressure relief &


High density foam provides support, whilst a top layer of

memory foam allows submersion into the cushion, providing

a large surface area to reduce pressure levels.


For those sitting for long periods with specific comfort issues, such as a high risk of pressure sores. Uses the latest medical grade silicone gel technology allowing pressure to be evenly distributed

The gel conforms to the shape of the body, equalises to body temperature, allows air to flow easily & is totally odourless & hypoallergenic.



For those sitting for long periods of time with specific comfort issues, such as a high risk of pressure sores & sheer forces.

Gelacell's unique grid design eliminates pressure points. contouring in response to the transfer of weight, evenly distributes pressure & reduces heat absorbtion, keeping the seat cooler & fresher



Airform



For those who have good posture, but are unable to alter their position easily & require high level pressure care.

This cushion provides the benefits of an air cushion, however without the sensation of alternating changes in air pressure.

Suitable for users 40 - 127kg.

**Dynaform** 

constantly adapting to movement.

Pressure Risk VERY HIGH

PUMP

For those at very high risk of developing pressure sores, or have a history of sores, where the priority is to prevent further pressure damage.

for long periods of time, require pressure relief & supported

High density foam provides support, topped with layers

of fluid liquid gel & memory foam to provide high levels of

pressure relief. The liquid gel helps to dissipate heat, whilst

comfort. Excellent for people with little fatty tissue.

Dynaform is a dynamic alternating air cushion, which constantly changes air levels within its cells to recreate the effect of pressure relieving & prevent pressure build up.

Air cell features

Alternating

air cushion

Reflex foam

- 6 dynamic cells 10cm wide
- 3 cell cycle repeated every 9 minutes

Suitable for users 40 - 121kg.

Please note the Dynaform cushion does need an electric supply at all times to operate.

Pressure risk levels detailed are for guidance only. Each individual user has different posture, medical needs & health issues which contribute to the performance of any pressure cushion. If in doubt, please consult a GP or Health Care Professional.

## **Back Cushion Options**



Waterfall Back Cushion Individual waterfall cushions, fully adjustable to accommodate maximum postural support



#### Large Profile Headrest with Waterfall Back

For those with limited head control & requiring a soft adjustable lumbar support

## **Headrest Options**



**Head Pillow** 




**Comfort Lateral Back Cushion** Soft fibre filled removable lateral cushion for lower back & side lumber support, without feeling restricted



Large Profile Headrest with Comfort Lateral Back

For those with limited head control & requiring soft lateral support without feeling restricted

**Head Roll** 



Support Lateral Back Cushion A foam filled back providing firmer lateral support with a soft fibre filled headrest



Large profile headrest with Support Lateral Back For those with limited head control. requiring firm foam lateral lumbar & midline support

**Table with Wireless Charging** 360-degree rotation & angle adjustable table with a wipe clean surface. Supports QI wireless charging & fitted with a USB socket. Removable magnetic bar provides support & ensures correct positioning on the wireless receiver coil,

wipe clean surface (w 345 x h 245mm)





**Small Profile Headrest** 



Large Profile Headrest



**Midline Headrest** 

sales@bristolmaid.com





## **Overnight Chair - Portland, Three Position**

#### **Features**

- Multi-functional compact chair that provides comfortable seating by day, transforming into a bed by night
- Suitable for use in maternity units, children's wards, relative's rooms, hospices & care homes
- · Provides visitors with a comfortable place to sleep during longer visits
- · Chair is pulled out in three motions to a bed
  - Chair
  - Lounge Chair
  - Bed
- Upholstered using flame retardant antimicrobial vinyl
- Easiform reflex foam provides comfort when sitting or sleeping
- 50mm swivel braking castors
- Maximum capacity 160kg

#### Dimensions (w x d x h)

- Overall Dimensions 870 x 840 x 870mm (when open overall bed length is 2060mm)
- Seat 690 x 540mm at Height 410mm

#### **Standard Colours**

- Vinyl Delft
- · Vinyl & Patterned Fabric Vinyl Lavendel & Savanna Lavender



- Delivery to Chosen Location
- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff

Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Portland, Three Position in Olive - 5RE/OC/PORT/VVOLI

Vinyl - Spradling Valencia C5



ItemDescriptionWeight5RE/OC/PORT/VVDELThree Position, Vinyl Delft61.5kg5RE/OC/PORT/VVLAVSAVLThree Position, Vinyl Lavendel & Savanna Lavender61.5kg

sales@bristolmaid.com

+44 (0)1258 484455







## Soft Seating - Zeta, Modular

#### **Features**

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- Zeta is a stylish, fully modular range of soft seating based around eight standard elements suitable for use in seating environments, ranging from reception areas, breakout spaces to waiting areas & cafeterias
- FSC® certified wood core, constructed from 18mm boards
- Traditional webbed seat for additional comfort & durability
- High density CMHR foam
- · Vinyl or fabric upholstery
  - · Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Chrome finish legs
- Arms available on single & double units with backs
- Backs are available with a raised feature to allow fluid run off (overall height 810mm)
- · Maximum capacity
  - Single 146kg

#### **Options**

- RA Right Arm
- LA Left Arm
- 2A Left & Right Arms

Replace suffix on the product code with any of the following to select your colour i.e. Single Unit in Cobalt - 5ZETA1/FCT

#### Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage



| Item        | Description                               | Overall Dimensions<br>(w x d x h) | Seat Height | Seat Depth | Weight |
|-------------|---|-----------------------------------|-------------|------------|--------|
| 5ZETA1      | Single Unit                               | 650 x 650 x 415mm                 | 415mm       | 650mm      | 12.0kg |
| 5ZETA2      | Double Unit                               | 1300 x 650 x 415mm                | 415mm       | 650mm      | 22.0kg |
| 5ZETA1/B    | Single Unit with Back                     | 650 x 650 x 765mm                 | 415mm       | 485mm      | 16.5kg |
| 5ZETA2/2B   | Double Unit with Back                     | 1300 x 650 x 765mm                | 415mm       | 485mm      | 31.0kg |
| 5ZETA1/CO   | Single Unit with Corner                   | 650 x 650 x 765mm                 | 415mm       | 485mm      | 16.5kg |
| 5ZETA2/B/RC | Double Unit with Right Corner & Left Back | 1300 x 650 x 765mm                | 415mm       | 485mm      | 35.5kg |
| 5ZETA2/B/LC | Double Unit with Left Corner & Right Back | 1300 x 650 x 765mm                | 415mm       | 485mm      | 35.5kg |
| 5ZETAQ      | Quarter Segment                           | 650 x 650 x 415mm                 | 415mm       | 650mm      | 12.0kg |
| 5ZETAH      | Half Segment                              | 1300 x 650 x 415mm                | 415mm       | 650mm      | 23.0kg |
| 5ZETAH/B    | Half Segment with Back                    | 1300 x 650 x 765mm                | 415mm       | 485mm      | 25.0kg |
| 5ZETAC45    | 45 Degree Curve                           | 950 x 650 x 415mm                 | 415mm       | 650mm      | 11.5kg |
| 5ZETAC45/B  | 45 Degree Curve with Back                 | 950 x 650 x 765mm                 | 415mm       | 485mm      | 16.5kg |
| 5ZETAC90    | 90 Degree Curve                           | 1625 x 650 x 415mm                | 415mm       | 650mm      | 22.0kg |
| 5ZETAC90/B  | 90 Degree Curve with Back                 | 1625 x 650 x 765mm                | 415mm       | 485mm      | 32.0kg |
| 5ZETA930C   | Single Corner Unit                        | 930 x 650 x 765mm                 | 415mm       | 650mm      | 28.0kg |

sales@bristolmaid.com

+44 (0)1258 484455

```
www.bristolmaid.com
```





## Visitor & Waiting Room Seating - Mata Chairs

#### **Features**

- BS EN 1729-1: Part 1 & 2: 2015 (functional dimensions & safety requirements)
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Available with...
  - Four legs
  - Sled base
- Moulded polypropylene shell available in seven colours
- Silver finish four leg
- Chrome finish sled base
- Four leg chair stackable five high
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Four Leg Chair in Olive - 5SM6L/SOE

Plastic Shell



sales@bristolmaid.com +44

+44 (0)1258 484455

## **Visitor & Waiting Room Seating - Prima Chairs**

#### **Features**

- BS EN 16139: 2013 Level 2 (strength, durability & safety)
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Certified for use up to 203kg
- Choice of vinyl, fabric upholstery or polypropylene seat & back
  - Vinyl Fire retardant & antimicrobial
  - · Fabric Fire retardant
- Stable four leg design
- · Black finish tubular frame
- Stacks up to ten high (transporter available)
- Built-in linking system (optional extension links)
- Maximum capacity 203kg

Replace suffix on the product code with any of the following to select your colour i.e. Four Leg Chair in Cobalt - 5301/FCT

#### Vinyl - Ambla New Richmond in Colour



Fabric - Camira Advantage



FKR SBK



FAC

FRD

5301/A/FRD

5301/VLI

| ltem   | Description         | Overall Dimensions (w x d x h) | Seat Height | Weight |
|--------|---------------------|--------------------------------|-------------|--------|
| 5301   | Four Leg            | 520 x 580 x 800mm              | 485mm       | 8.0kg  |
| 5301/A | Four Leg, Arm Rests | 620 x 580 x 800mm              | 485mm       | 10.0kg |





| ltem    | Description             | Overall Dimensions (w x d x h) | Seat Height | Weight |
|---------|-------------------------|--------------------------------|-------------|--------|
| 5NOVA/1 | Nova - Tub, Single Seat | 680 x 680 x 800mm              | 450mm       | 16.0kg |
| 5NOVA/2 | Nova - Tub, Double Seat | 1360 x 680 x 800mm             | 450mm       | 31.0kg |

## Waiting Room Seating - Nova, Tub

#### **Features**

- Nova single & double tub chairs are suitable for a wide range of seating environments, ranging from cafeterias to waiting areas
- FSC<sup>®</sup> certified wood core with moulded deep CMHR seat cushions
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Chrome finish legs
- Maximum capacity
  - Single 146kg
  - Double 292kg

Replace suffix on the product code with any of the following to select your colour i.e. Double Seat in Cobalt - 5NOVA/3/FCT

Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage



## Waiting Room Seating - Horizon, Beech Frame

#### **Features**

- Horizon chairs, sofas & tables are designed to suit a variety of waiting & occasional seating areas
- Beech frame with plywood seat shells
- Arms assist positioning from standing & seated positions
- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
  - · Fabric Fire retardant
- Fixed seat, with a loose back cushion to assist cleaning
- Back cushion has a zipped cover, allowing for easy replacement if damaged
- · Maximum capacity
  - Single chair 146kg
  - Two seater sofa 120kg per seat
  - Three seater sofa 120kg per seat

Replace suffix on the product code with any of the following to select your colour i.e. Two Seater Sofa with Arms in Cobalt - 5HORI/2/A/FCT

Vinyl - Ambla New Richmond in Colour



Fabric - Camira Advantage







| Item      | Description                       | Overall Dimensions (w x d x h) | Seat Height | Weight |
|-----------|-----------------------------------|--------------------------------|-------------|--------|
| 5HORI/1/A | Horizon - Easy Chair, Arms        | 750 x 650 x 800mm              | 450mm       | 19.0kg |
| 5HORI/2/A | Horizon - Two Seater Sofa, Arms   | 1350 x 650 x 800mm             | 450mm       | 33.0kg |
| 5HORI/3/A | Horizon - Three Seater Sofa, Arms | 1950 x 650 x 800mm             | 450mm       | 47.0kg |

+44 (0)1258 484455 www





## Waiting Room Seating - Prima Beam, Fixed

#### **Features**

- Prima beam is suitable for all environments requiring durable & robust waiting room seating solutions
- Certified for use up to 32 stone (203kg) per seat
- Fixed seat
- Choice of vinyl, fabric upholstery or polypropylene seat & back
  - Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Black finish beam & legs
- Drilled for floor fixing (fixings not supplied)
- Maximum capacity 203kg per seat

#### **Options**

• 510065 - Armrest (x1)

Replace suffix on the product code with any of the following to select your colour i.e. Fixed Beam Seat, Two Seat in Cobalt - 5PF332/FCT

#### Vinyl - Ambla New Richmond in Colour





#### **Fixed Seat**

| Item   | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|--------|-------------|--------------------------------|-------------|--------|
| 5PF332 | Two Seat    | 1200 x 560 x 790mm             | 480mm       | 20.0kg |
| 5PF333 | Three Seat  | 1800 x 560 x 790mm             | 480mm       | 25.0kg |
| 5PF334 | Four Seat   | 2400 x 560 x 790mm             | 480mm       | 30.0kg |
| 5PF335 | Five Seat   | 3000 x 560 x 790mm             | 480mm       | 35.0kg |

sales@bristolmaid.com

+44 (0)1258 484455

www.bristolmaid.com



# **Options**

# Waiting Room Seating - Prima Beam, Tip-Up Seat

#### **Features**

- Prima beam is suitable for all environments requiring durable & robust waiting room seating solutions
- · Certified for use up to 32 stone (203kg) per seat
- Tip-up seat (must be fixed to the floor)
- Choice of vinyl, fabric upholstery or polypropylene seat & back
  - · Vinyl Fire retardant & antimicrobial
  - Fabric Fire retardant
- Black finish beam & leas
- Drilled for floor fixing (fixings not supplied)
- Maximum capacity 203kg per seat
- 510065 Armrest (x1)

Replace suffix on the product code with any of the following to select your colour i.e. Fixed Beam Seat. Two Seat in Cobalt - 5PF332/FCT

#### Vinyl - Ambla New Richmond in Colour





5PT338/FRD Shown with 510065 (optional arms sold individually) - contact us for details

#### **Tip-Up Seat**

| Item   | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|--------|-------------|--------------------------------|-------------|--------|
| 5PT336 | Two Seat    | 1200 x 490 x 790mm             | 490mm       | 20.0kg |
| 5PT337 | Three Seat  | 1800 x 490 x 790mm             | 490mm       | 25.0kg |
| 5PT338 | Four Seat   | 2400 x 490 x 790mm             | 490mm       | 30.0kg |
| 5PT339 | Five Seat   | 3000 x 490 x 790mm             | 490mm       | 35.0kg |

sales@bristolmaid.com

+44 (0)1258 484455



# **Tech Stools**

#### Features

- Vinyl or fabric upholstery
  - Vinyl Fire retardant & antimicrobial
- Fabric Fire retardant
- Black polypropylene five star base
- Models
  - Medium/High Foot ring & glides
  - Low Non-braking castors
- Black gas lift with pneumatic seat height adjustment
- Maximum load 120kg

#### Options

• 5BLSTC - Brake loaded castors



| Item   | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|-------------|--------------------------------|------------------------|--------|
| 5TC110 | Low         | 670 x 670 x 440-570mm          | 440 - 570mm            | 8.0kg  |
| 5TC109 | Medium      | 670 x 670 x 570-690mm          | 570 - 690mm            | 10.0kg |
| 5TC108 | High        | 670 x 670 x 650-810mm          | 650 - 810mm            | 10.0kg |

Replace suffix on the product code with any of the following to select your colour i.e. Low Stool in Kingfisher - 5TC110/FKR

Vinyl - Ambla New Richmond in Colour



#### Fabric - Camira Advantage



sales@bristolmaid.com

+44 (0)1258 484455 www.bris





| Item | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|------|----------------------------|--------------------------------|------------------------|--------|
| 5PU6 | Low, Castors               | 670 x 670 x 750-980mm          | 430 - 570mm            | 13.0kg |
| 5PU5 | Medium, Foot Ring & Glides | 710 x 710 x 830-1120mm         | 510 - 700mm            | 14.0kg |
| 5PU4 | High, Foot Ring & Glides   | 710 x 710 x 910-1270mm         | 580 - 840mm            | 14.0kg |

www.bristolmaid.com

# **PU Tech Stools**

#### Features

- Soft touch black polyurethane seatBlack polypropylene five star base
- Models





| Item   | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5PU153 | Low, Castors               | 670 x 670 x 410-540mm          | 410 - 540mm            | 8.0kg  |
| 5PU152 | Medium, Foot Ring & Glides | 670 x 670 x 470-660mm          | 470 - 660mm            | 10.0kg |
| 5PU151 | High, Foot Ring & Glides   | 670 x 670 x 530-780mm          | 530 - 780mm            | 10.0kg |



|        | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC203 | Low, Castors               | 670 x 670 x 890-1090mm         | 440 - 580mm            | 13.0kg |
| 5TC202 | Medium, Foot Ring & Glides | 710 x 710 x 950-1210mm         | 520 - 710mm            | 14.0kg |
| 5TC201 | High, Foot Ring & Glides   | 710 x 710 x 980-1320mm         | 570 - 820mm            | 14.0kg |

Seating Solutions

# **Static Safe Tech Stools**

#### Features

- Hi-Tech ESD (anti-static) fire retardant black fabric upholstery
- Aluminium five star base
- Models
  - Medium/High Foot ring & glides
  - Low ESD free wheeling castors

5TC210

- Chrome gas lift with pneumatic seat height adjustment
- Supplied with ESD (anti-static) castors or glides
- Maximum capacity 120kg

#### Options

• Brake loaded ESD castors







670 x 670 x 540-790mm

www.bristolmaid.com

540 - 790mm

| 5 | TC209  |              | 5TC208                         |                        |        |
|---|--------|--------------|--------------------------------|------------------------|--------|
|   | Item   | Description  | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|   | 5TC210 | Low, Castors | 670 x 670 x 386-556mm          | 386 - 556mm            | 10.0kg |
|   |        |              |                                |                        |        |

5TC208 High, Foot Ring & Glides

11.0kg



#### Static Safe & Sterile Tech Chairs

- ESD (anti-static) fire retardant & antimicrobial black vinyl upholstery
- Aluminium five star base
  - Medium/High Foot ring & glides
  - Low Castors
- Chrome gas lift with pneumatic seat & back rest height adjustment
- Supplied with ESD (anti-static) castors or glides
- Maximum capacity 120kg





**HEPA** filter

14.0kg

14.0kg

www.bristolmaid.com

520 - 710mm

570 - 8 20mm

5TC304 High, Foot Ring & Glides

5TC305 Medium, Foot Ring & Glides 710 x 710 x 950-1210mm

710 x 710 x 980-1320mm

# Static Safe & Sterile Tech Stools

#### Features

- ESD (anti-static) fire retardant & antimicrobial black vinyl upholstery
- Aluminium five star base
- Models
  - Medium/High Foot ring & glides
  - Low Castors







| Item   | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC310 | Low, Castors               | 670 x 670 x 386-556mm          | 386 - 556mm            | 10.0kg |
| 5TC309 | Medium, Foot Ring & Glides | 670 x 670 x 490-680mm          | 490 - 680mm            | 11.0kg |
| 5TC308 | High, Foot Ring & Glides   | 670 x 670 x 540-790mm          | 540 - 790mm            | 11.0kg |



|        | Description                | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC313 | Low, Castors               | 670 x 670 x 890-1090mm         | 440 - 580mm            | 13.0kg |
| 5TC312 | Medium, Foot Ring & Glides | 710 x 710 x 950-1210mm         | 520 - 710mm            | 14.0kg |
| 5TC311 | High, Foot Ring & Glides   | 710 x 710 x 980-1320mm         | 570 - 820mm            | 14.0kg |

## Medical Grade Clean Room Tech Chairs

- Medical grade fire retardant & antimicrobial black vinyl upholstery
- Aluminium five star base
  - Medium/High Foot ring & glides
  - Low Castors
- Chrome gas lift with pneumatic seat & back rest height adjustment
- Maximum capacity 120kg
- Independent seat & back mechanism
- Brake loaded castors



HEPA filter

sales@bristolmaid.com

www.bristolmaid.com

# Other products that you may be interested in...

See Page 118

### Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep





# **Sepsis Trolleys**

- Sepsis Trolleys 382
  - Accessories 383
- Caretray Trolley Sepsis 386



#### **Sepsis Trolleys**

#### **Features**

- Trollevs designed for use when treating Sepsis
- Supports delivery of 'The Sepsis Six' with clearly defined locations for equipment used in each of the stages
- · Door/s fitted with either:
  - Electronic (battery operated) Push Button Lock Customer can set their own separate user, master, submaster & technician access codes
  - High Security Bolt Lock Supplied with two keys
- Trolleys supplied with door shelves (useful storage for antibiotics etc) & an A4 document holder
  - Single Door 2 x door shelves / 1 x A4 document holder
- Double Door 5 x door shelves / 1 x A4 document holder
- Six labelled drawers (Oxygen / Blood Cultures / Intravenous Antibiotics / Intravenous Fluid /Haemoglobin & Lactates / Urine Measurement)
- Optional 'Sharps Box Holder' or 'Universal Sharps Box Bracket Mounting Plate' (fits to the I/h side of the trollev)
- Sharps Box Holder Accepts up to seven litre capacity sharps boxes (customer to fit - fixings supplied)
- Universal Sharps Box Bracket Mounting Plate Compatible with (customer to supply the plastic brackets & fit - fixings supplied)
  - Daniels Healthcare
    - POUDS<sup>®</sup> Wall-Bracket (up to 5L) SHARPSGUARD<sup>®</sup> mini, 1, 2.5. 3.75. 5 & com plus
    - POUDS<sup>®</sup> Wall-Bracket (7L & above) SHARPSGUARD<sup>®</sup> 7. 11.5 & 22
  - Frontier (Sharpsafe) 2, 3, 4 & 7
  - Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS
- · Upper work surface
- · Push handles incorporate integral buffers
- 125mm polymer castors
- . When not in use, trolleys storing medicines should be secured to the wall. See option MT/1/6F Security Clamp

#### Dimensions (w x d x h)

- Single Door Trolley 512 x 542 x 1072mm
- Double Door Trolley 712 x 542 x 1072mm
- Drawers
  - Single Door Trolley 338 x 297 x 115mm (internal)
  - Double Door Trolley 538 x 297 x 115mm (internal)

#### Standard Colours

- · Paint Grev
- Work Surface Orange

Sepsis Trolleys

Weight

53.4kg

69.1kg

# **Sharps Box Holder**

#### **Features**

SEP/SBH

- · Accepts up to seven litre capacity sharps boxes
- Sharps box secured using a rubber band fitted with a toggle
- Attaches to the left hand side of the unit (customer to fit fixings supplied)



#### **Features**

- Compatible with the following brackets (customer to supply the plastic brackets & fit fixings supplied)
   Daniels Healthcare
  - POUDS® Wall-Bracket (up to 5L) SHARPSGUARD® mini, 1, 2.5, 3.75, 5 & com plus
  - POUDS® Wall-Bracket (7L & above) SHARPSGUARD® 7, 11.5 & 22
  - Frontier (Sharpsafe) 2, 3, 4 & 7
  - Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS




sales@bristolmaid.com

+44 (0)1258 484455 wv

# **Security Clamp**

#### **Features**

- · Allows the trolley to be secured to the 'fabric' of the building
- Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door of the trolley when closed and locked
- · Recommended wall fitting height 1000mm



#### **Important Note - Security Clamps**

It is a UK requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area and broken into, but also encourages good nursing practice, giving the trolley a 'permanent home'.

'Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB- March 2005'


| Item    | Description    | Cable Length | Weight |
|---------|----------------|--------------|--------|
| MT/1/6F | Security Clamp | 950mm        | 0.8kg  |

sales@bristolmaid.com

# **Drawer Inserts**

#### **Features**

- Plastic drawer inserts, height 100mm
- Assists segregation of drawer contents
- · Suitable for single & double door trolleys
- Insert supplied with removable & repositionable dividers
  - SEP/DIS 2 x dividers providing 4 separate areas
  - SEP/DID 4 x dividers providing 6 separate areas
- Supplied for self-assembly
- Clip on label holders & plastic labels available to assist identification



| Item    | Description                                | Weight |
|---------|--|--------|
| SEP/DIS | Drawer Insert - Set - Single Door Trolleys | 0.5kg  |
| SEP/DID | Drawer Insert - Set - Double Door Trolleys | 0.6kg  |

# Label Holders & Labels

#### Features

- Plastic label holders & labels
- Assists identification of drawer contents
- Label holders snap onto drawer inserts/dividers
- Accepts ISO labels 40 x 60mm (w x h)
- Sold in packs of ten



| Item   | Description  |
|--------|--------------|
| SEP/LH | Label Holder |
| SEP/L  | Label        |


# Label - Sepsis Tiger

#### **Features**

- Sammy the sepsis tiger label developed for use on Paediatric trolleys
- Compatible with both single & double door trolleys



# Item Description SEP/TIG01 Label - Sammy Sepsis Tiger





# Item Description Weight CTSEPSM1/10/6S/HSB 6 x 100mm Trays, High Security Bolt Lock 36.8kg

sales@bristolmaid.com

# **Caretray Trolley - Sepsis**

#### **Features**

- · Trolleys designed for use in the treatment of Sepsis
- Supports delivery of 'The Sepsis Six' with clearly defined locations for equipment used in each of the stages
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Waste Bin & Bracket, 8 Litre
  - Sharps Box Holder
  - Labels, Sepsis Six
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (sepsis six labels supplied - oxygen / blood cultures / intravenous antibiotics / intravenous fluid / haemoglobin & lactates / urine measurement)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking
- When not in use, trolleys storing medicines should be secured to the wall using a
   Security Clamp (see locking & security accessories page 384)

#### Dimensions (w x d x h)

#### **Standard Colours**

Paint - Orange

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays

+44 (0)1258 484455

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

#### **Options / Accessories**

- Locking & Security page 60
- Overbridge & Tilt Bins page 66
  Defibrillator Platform page 67

- Carerail Mounted page 68
  Laptop Holders page 74
  Trays, Shelves, Dividers & Labels page 77





| Item               | Description                                  | Weight |
|--------------------|--|--------|
| CTSEPSM1/10/6S/EPB | 6 x 100mm Trays, Electronic Push Button Lock | 36.8kg |



# **Shelving & Racking Solutions**

Engineered for a Demanding Environment 390

- Polymer Racking MetroMax I<sup>®</sup> 392
- Polymer Racking MetroMax Q<sup>®</sup> 394
- Zinc Chrome Racking Olympic<sup>®</sup> 396
- High Density Storage Top-Track<sup>®</sup> 398
- High Density Storage qwikTRAK<sup>®</sup> 399

# **Engineered for a Demanding Environment**



# MetroMax I<sup>®</sup> Full Polymer System with Lifetime Warranty against corrosion





Microban<sup>®</sup> antimicrobial product protection helps keep shelves "cleaner between cleanings" by inhibiting the growth of mold, mildew, fungus and bacteria that cause odors and stains on the shelf surface.



MetroMax Q<sup>®</sup> Polymer / Epoxy Coated System with 15 Year Warranty against corrosion









#### **Stationary Units**

| 4 Shelves -  | 1590mm High   | 5 Shelves - 1 | 895mm High    | Shelf Size (w x d) |
|--------------|---------------|---------------|---------------|--------------------|
| Grid Shelves | Solid Shelves | Grid Shelves  | Solid Shelves |                    |
| 8MX536SG     | 8MX536SS      | 85MX536SG     | 85MX536SS     | 914 x 610mm        |
| 8MX546SG     | 8MX546SS      | 85MX546SG     | 85MX546SS     | 1066 x 610mm       |
| 8MX556SG     | 8MX556SS      | 85MX556SG     | 85MX556SS     | 1219 x 610mm       |
| 8MX566SG     | 8MX566SS      | 85MX566SG     | 85MX566SS     | 1524 x 610mm       |

#### **Mobile Units**

| 4 Shelves - 1 | 1730mm High   | 5 Shelves - 2 | 2035mm High   | Shelf Size (d x w) |
|---------------|---------------|---------------|---------------|--------------------|
| Grid Shelves  | Solid Shelves | Grid Shelves  | Solid Shelves |                    |
| 8MX536MG      | 8MX536MS      | 85MX536MG     | 85MX536MS     | 914 x 610mm        |
| 8MX546MG      | 8MX546MS      | 85MX546MG     | 85MX546MS     | 1066 x 610mm       |
| 8MX556MG      | 8MX556MS      | 85MX556MG     | 85MX556MS     | 1219 x 610mm       |
| 8MX566MG      | 8MX566MS      | 85MX566MG     | 85MX566MS     | 1524 x 610mm       |

sales@bristolmaid.com

# **Polymer Racking - MetroMax I®**

#### **Features**

- Corrosion proof polymer shelving & racking system
- Polymer shelves with stainless steel corners & a patented no tools assembly, adjustable in 25mm increments
- All polymer posts
- · Removable easy to clean polymer shelf mats (will fit most dishwashers) available as...
  - · Grid Shelves Open grid design allowing for the free circulation of air, improved visibility & light penetration
  - Solid Shelves Solid design offering a tear free surface, ideal for the storage of sterilised DIN/ISO travs within theatres
- Microban<sup>®</sup> product protection inhibits growth of bacteria & aids infection control
- · Available as stationary or mobile units
  - Stationary Fitted with adjustable feet to allow levelling on uneven floors
  - Mobile Fitted with 125mm swivel / braking castors
- Rack / unit capacity...
  - Stationary 907kg
  - Mobile 339kg
- Optional label holders, shelf dividers & ledges
- · Lifetime anti-corrosion warranty
- · Supplied for self-assembly



Solid Shelves

Contact us for more information

+44 (0)1258 484455

www.bristolmaid.com



## **Label Holders**

#### **Features**


- · Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels

#### Dimensions (w x h)

• 102 x 27mm

#### **Standard Colours**

• Grey



# **Shelf Divider**

#### **Features**

• Corrosion proof polymer shelf divider (height 200mm)

+44 (0)1258 484455

- Use to create organised compartments or simply separate shelf contents
- Easily repositioned on grid & solid shelves
- Easy to clean

#### Dimensions (w x h)

• 610 x 200mm

#### **Standard Colours**

• Grey

sales@bristolmaid.com

| Item    | Description  | Weight | Item     | Description   | Weight |
|---------|--------------|--------|----------|---------------|--------|
| 89989PX | Label Holder | 0.01kg | 8MXD24-8 | Shelf Divider | 0.7kg  |



# **Stackable Shelf Ledges**

#### **Features**

- Corrosion proof stackable shelf ledge (height 51mm)
- Prevents larger items from protruding or falling
- Compatible with grid & solid shelves
- Easy to clean

#### **Standard Colours**

Chrome

| Item       | Description         | Weight |
|------------|---------------------|--------|
| Side Ledge |                     |        |
| 8MXLS24-2S | Side Ledge - 610mm  | 0.2kg  |
| Rear Ledge |                     |        |
| 8MXL36-2S  | Rear Ledge - 914mm  | 0.4kg  |
| 8MXL42-2S  | Rear Ledge - 1066mm | 0.6kg  |
| 8MXL48-2S  | Rear Ledge - 1219mm | 0.6kg  |
| 8MXL60-2S  | Rear Ledge - 1524mm | 0.8kg  |

| www.bristolmaid.com |
|---------------------|
|---------------------|





**Mobile Units** 

sales@bristolmaid.com

# Polymer Racking - MetroMax Q<sup>®</sup>

#### **Features**

- · Epoxy / polymer shelving & racking system
- · Epoxy coated shelves & posts with a patented no tools lever shelf release mechanism, adjustable in 25mm increments
- · Removable easy to clean polymer open grid shelf mats (will fit most dishwashers)
- Open grid design allowing for the free circulation of air, improved
- visibility & light penetration
- Microban® product protection inhibits growth of bacteria & aids infection control
- Available as stationary or mobile units
- Stationary Fitted with adjustable feet to allow levelling on uneven floors
- Mobile Fitted with 125mm swivel / braking castors
- Rack / unit capacity...
  - Stationary 907kg
  - Mobile 339kg
- Optional label holders, shelf dividers & ledges
- · Anti-corrosion warranty
  - Frame & Posts 15 years
  - Shelf Mats Lifetime
- · Supplied for self-assembly



Contact us for more information

#### **Stationary Units**

| 4 Shelves             | Shelf Size (w x d)                | Height           | Weight           | 4 Shelves              | Shelf Size (w x d   |
|-----------------------|-----------------------------------|------------------|------------------|------------------------|---|
| 8Q536G2               | 914 x 610mm                       | 1590mm           | 30.0kg           | 8Q536G2M               | 914 x 610mm   |
| 8Q546G2               | 1066 x 610mm                      | 1590mm           | 33.6kg           | 8Q546G2M               | 1066 x 610mm  |
| 8Q556G2               | 1219 x 610mm                      | 1590mm           | 35.6kg           | 8Q556G2M               | 1219 x 610mm  |
| 8Q566G2               | 1524 x 610mm                      | 1590mm           | 44.4kg           | 8Q566G2M               | 1524 x 610mm  |
|                       |                                   |                  |                  |                        |   |
| 5 Shelves             | Shelf Size (w x d)                | Height           | Weight           | 5 Shelves              | Shelf Size (w x d   |
| 5 Shelves<br>85Q537G2 | Shelf Size (w x d)<br>914 x 610mm | Height<br>1895mm | Weight<br>36.7kg | 5 Shelves<br>85Q537G2N |   |
|                       |                                   |                  |                  |                        | 1 914 x 610mm   |
| 85Q537G2              | 914 x 610mm                       | 1895mm           | 36.7kg           | 85Q537G2N              | 1 914 x 610mm<br>1 1066 x 610mm   |
| 85Q537G2<br>85Q547G2  | 914 x 610mm<br>1066 x 610mm       | 1895mm<br>1895mm | 36.7kg<br>41.2kg | 85Q537G2N<br>85Q547G2N | <ol> <li>914 x 610mm</li> <li>1066 x 610mm</li> <li>1219 x 610mm</li> </ol> |

85Q566G2

| bile Unit | S *other siz       | zes available | on request |
|-----------|--------------------|---------------|------------|
| helves    | Shelf Size (w x d) | Height        | Weight     |
| 536G2M    | 914 x 610mm        | 1730mm        | 34.6kg     |
| 546G2M    | 1066 x 610mm       | 1730mm        | 38.2kg     |
| 556G2M    | 1219 x 610mm       | 1730mm        | 40.2kg     |
| 566G2M    | 1524 x 610mm       | 1730mm        | 49.0kg     |
| helves    | Shelf Size (w x d) | Height        | Weight     |
| Q537G2M   | 914 x 610mm        | 2035mm        | 41.3kg     |
| Q547G2M   | 1066 x 610mm       | 2035mm        | 45.8kg     |

Room Design & Installation Service Available

2035mm

48.3kg

59.3kg

www.bristolmaid.com



# **Label Holders**

#### **Features**

- · Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels


#### Dimensions (w x h)

• 102 x 27mm

#### **Standard Colours**

• Grey

8Q04LH



# **Shelf Divider**

#### **Features**

- Corrosion proof polymer shelf divider (height 200mm)
- Use to create organised compartments or simply separate shelf contents
- Easily repositioned on grid & solid shelves
- · Easy to clean

#### Dimensions (w x h)

• 610 x 200mm

#### **Standard Colours**

• Grey

| Weight | Item     | Description   | Weight |
|--------|----------|---------------|--------|
| 0.01kg | 8MQD24-8 | Shelf Divider | 0.7kg  |



# **Stackable Shelf Ledges**

#### **Features**

- Corrosion proof stackable shelf ledge (height 51mm)
- Prevents larger items from protruding or falling
- Compatible with grid & solid shelves
- Easy to clean

#### **Standard Colours**

Chrome

| Item       | Description         | Weight |
|------------|---------------------|--------|
| Side Ledge |                     |        |
| 8MQLS24-2S | Side Ledge - 610mm  | 0.2kg  |
| Rear Ledge |                     |        |
| 8MQL24-2S  | Rear Ledge - 610mm  | 0.2kg  |
| 8MQL36-2S  | Rear Ledge - 914mm  | 0.4kg  |
| 8MQL42-2S  | Rear Ledge - 1066mm | 0.6kg  |
| 8MQL48-2S  | Rear Ledge - 1219mm | 0.6kg  |
| 8MQL60-2S  | Rear Ledge - 1524mm | 0.8kg  |

| Label Holder | 0.01kg |
|--------------|--------|
|              |        |

| sales@bristolmaid.com |  |
|-----------------------|--|
|-----------------------|--|

# DLYIVIPIC



#### **Stationary Units**

#### **Mobile Units**


| 4 Shelves | Shelf Size (w x d) | Height | Weight | 4 Shelves | Shelf Size (w x d) | Height | W |
|-----------|--------------------|--------|--------|-----------|--------------------|--------|---|
| 80536C    | 914 x 610mm        | 1590mm | 30.0kg | 80536CM   | 914 x 610mm        | 1730mm | 3 |
| 8O546C    | 1066 x 610mm       | 1590mm | 33.6kg | 80546CM   | 1066 x 610mm       | 1730mm | 3 |
| 80556C    | 1219 x 610mm       | 1590mm | 35.6kg | 80556CM   | 1219 x 610mm       | 1730mm | 4 |
| 80966C    | 1524 x 610mm       | 1590mm | 44.4kg | 80966CM   | 1524 x 610mm       | 1730mm | 4 |
| 5 Shelves | Shelf Size (w x d) | Height | Weight | 5 Shelves | Shelf Size (w x d) | Height | N |
| 850537C   | 914 x 610mm        | 1895mm | 36.7kg | 850537CM  | 914 x 610mm        | 2035mm | 2 |
| 850547C   | 1066 x 610mm       | 1895mm | 41.2kg | 850547CM  | 1066 x 610mm       | 2035mm | 2 |
| 850557C   | 1219 x 610mm       | 1895mm | 43.7kg | 850557CM  | 1219 x 610mm       | 2035mm | 2 |
| 85O967C   | 1524 x 610mm       | 1895mm | 54.7kg | 85O967CM  | 1524 x 610mm       | 2035mm | Ę |

sales@bristolmaid.com

# Zinc Chrome Racking - Olympic<sup>®</sup>

#### **Features**

- Zinc chrome shelving & racking system
- · Chromate posts & shelves with a patented no tools tapered split sleeve system that snaps together on each post, adjustable in 25mm increments
- Available as stationary or mobile units
  - Stationary Fitted with adjustable feet to allow levelling on uneven floors
  - Mobile Fitted with 125mm swivel / braking castors
- Rack / unit capacity...
  - Stationary 907kg
  - Mobile 339kg
- Optional label holders, shelf dividers & ledges
- 1-year anti-corrosion warranty
- Supplied for self-assembly



Zinc Chrome Shelves

#### **Room Design & Installation Service Available** Contact us for more information



# **Label Holders**

#### **Features**

- Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels

#### Dimensions (w x h)

• 75 x 32mm

#### **Standard Colours**

Grey



# Shelf Divider

#### Features

• Keeps shelf contents orderly (height 203mm)

+44 (0)1258 484455


• Pressure fit to hold in place

#### **Dimensions** (w x h)

• 610 x 203mm

#### **Standard Colours**

Chrome



# **Shelf Ledges**

#### **Features**

- Prevents items from protruding or falling (height 100mm)
- Pressure fit to hold in place
- Ledges can be stacked

#### **Standard Colours**

Chrome

| Item     | Shelf Length | Weight |
|----------|--------------|--------|
| 8L24N-4C | 610mm        | 0.90kg |
| 8L36N-4C | 914mm        | 1.25kg |
| 8L42N-4C | 1066mm       | 1.48kg |
| 8L48N-4C | 1219mm       | 1.68kg |
| 8L60N-4C | 1524mm       | 1.81kg |

| ltem   | Description  | Weight | Item   | Description   | Weight |
|--------|--------------|--------|--------|---------------|--------|
| 89990P | Label Holder | 0.1kg  | 8DD24C | Shelf Divider | 0.7kg  |







# High Density Storage - Top-Track®

#### **Features**

- · Increases the storage capacity of a given area by up to 50%
- Maximises all available space
- Can eliminate the need to expand a given area through the use of optimised storage
- Creates moveable aisles
- · Eliminates hard to reach or 'dead' spaces
- The guide track is positioned above the shelving system
- Aluminium overhead track system acts as a guide for mobile units
- Floors are easy to clean
- Carts & trolleys can easily be easily manoeuvred into & out of the active aisle
- Mobile units are designed to accommodate medium-duty applications
- The maximum weight capacity of a Top-Track  $^{\!\otimes}$  mobile unit is 410kg



Overhead Guide Track

Room Design & Installation Service Available Contact us for more information

# High Density Storage - qwikTRAK®

#### **Features**

- Increases the storage capacity of a given area by up to 50%
- Maximises all available space
- Can eliminate the need to expand a given area through the use of optimised storage
- · Creates moveable aisles
- · Eliminates hard to reach or 'dead' spaces
- The guide track is positioned below the shelving system, & protects the floor from wear & tear (does not need to be mounted to the floor)
- · Compensates for uneven floor surfaces
- Track acts as a guide for mobile units
- · Mobile units are designed to accommodate heavy-duty applications
- Allows for double deep configurations, where two systems are joined together for even greater storage capacity
- The maximum weight capacity of a qwikTRAK® mobile unit is 970kg









Floor Mounted Guide Track



# **Sterile Services Equipment**

- Preparation Tables Variable Height 402
  - Preparation Tables Fixed Height 403
- Sterile Supplies Trolleys UN3291 Certified 404
- Sterile Supplies Trolleys UN3291 Certified, Trolley Washer Compatible 406

#### **Preparation Tables - Stainless Steel, Variable Height**

#### **Features**

- Height adjustable preparation tables
- 304 grade stainless steel frame & work surface
- Uses a manual turn handle (folds away when not in use) • Static - 715 - 1015mm
  - Mobile 812 1112mm
- Available as static or mobile
  - Static Adjustable feet for levelling on uneven floors
  - Mobile 80mm pressed steel castors, 2 off braking
- Max capacity 70.0kg

#### Dimensions (w x d x h)

- Static 1500 x 750 x 715 1015mm
- Mobile 1500 x 750 x 812 1112mm

| Stand | dard | Col | ours |
|-------|------|-----|------|
| Juan  | Jaiu | CU  | ouis |



| Item  | Description  | Weight |
|-------|--|--------|
| PT050 | Preparation Table - Variable Height, Stainless Steel, Static | 49.0kg |
| PT110 | Preparation Table - Variable Height, Stainless Steel, Mobile | 53.0kg |

# **Preparation Tables - Stainless Steel, Fixed Height**

#### Features

- Fixed height preparation tables
- 304 grade stainless steel frame, lower shelf & work surface
- Worksurface height 950mm
- Available as static or mobile
  - Static Adjustable feet for levelling on uneven floors
  - Mobile 80mm pressed steel castors, 2 off braking
- Max capacity 70.0kg

#### Dimensions (w x d x h)

• Overall Dimensions - 1500 x 750 x 950mm

#### **Standard Colours**

Stainless Steel



| Item      | Description   | Weight |
|-----------|---|--------|
| PTS/115   | Preparation Table - Fixed Height, Stainless Steel, Static | 56.3kg |
| PTS/115/M | Preparation Table - Fixed Height, Stainless Steel, Mobile | 58.8kg |



AGV Compatible Trolleys & Carts For Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy distribution

# Sterile Supplies Trolleys - UN3291 Certified

#### Features

- Independently tested to UN3291 Certificate of Packaging Performance for the Transportation of Clinical Waste
- · Suitable for the transportation of DIN & ISO instrument trays
- · Mild steel frame
- · Lightweight aluminium cladding, shelves & door/s
- · Powder coated for durability
- Locking door/s…
  - Fitted with a flush fitting handle fitted with a CAM lock offering two-point locking (supplied with 2 keys)
  - When opened, doors fold back flat at either end of the trolley\*
- Three adjustable shelves, (max 40kg per shelf)
- Shelf size...
  - Small 534 x 614mm (w x d)
  - Large 1128 x 614mm (w x d)
- Removable stainless steel drip tray
- · Inset push handles at either end of trolley
- Low level bumper rail offers protection to both the trolley & fabric of the building
- Optional tow bar fitting...
  - Tow bar fitted to the front (swivel castor)
  - Antiluce fastener fitted to rear (fixed castor)
  - Tow bar position approx. 360-380mm above floor
  - Recommended max speed 7kph for towed trolleys
- 200mm heavy duty cushioned castors, 2 x fixed, 2 x swivel braking
  Maximum capacity...
  - Small 84k.0g per trolley
  - Large 109.0kg per trolley

\*Doors will not fold back flat when a tow bar is fitted

#### Dimensions (w x d x h) Standard Colours

- Small 805 x 774 x 1257mm Frame Grey
- Large 1399 x 774 x 1257mm Cladding / Shelves / Doors White

| Item Description |  |         |  |
|------------------|--|---------|--|
| Small            |  |         |  |
| ST/340/UN        | Sterile Supplies Trolley - UN3291 Certified, Small         | 74.5kg  |  |
| ST/340/UN/T      | Sterile Supplies Trolley - UN3291 Certified, Towing, Small | 78.5kg  |  |
| Large            |  |         |  |
| ST/690/UN        | Sterile Supplies Trolley - UN3291 Certified, Large         | 110.0kg |  |
| ST/690/UN/T      | Sterile Supplies Trolley - UN3291 Certified, Towing, Large | 114.0kg |  |



# Sterile Supplies Trolleys - UN3291 Certified, Trolley Washer Compatible



#### Features

- Independently tested to UN3291 Certificate of Packaging Performance for the Transportation of Clinical Waste
- · Compatible with trolley washers
- Suitable for the transportation of DIN & ISO instrument trays
- 304 grade stainless steel frame, cladding, shelves & doors
- Natural unpolished finish
- Locking door/s...
  - Fitted with a flush fitting handle fitted with a CAM lock offering two-point locking (supplied with 2 keys)
  - When opened, doors fold back flat at either end of the trolley
  - Spring loaded plungers hold the doors open, essential when the trolley is being processed in a trolley washer disinfector
- One fixed, two adjustable shelves (max 50kg per shelf)
- Shelf size 1227 x 632mm (w x d)
- Stainless steel drip tray...
  - Removable for cleaning
  - Can be located at the rear of the trolley during cleaning, held in position using spring loaded plungers
- · Inset push handles at either end of trolley
- High level corner bumpers offer protection to both the trolley & fabric of the building
- Optional tow bar fitting...
  - Tow bar fitted to the front (swivel castor)
  - Antiluce fastener fitted to rear (fixed castor)
  - Tow bar position approx. 360-380mm above floor
  - Recommended max speed 7kph for towed trolleys
- 200mm heavy duty sealed bearing stainless steel castors, 2 x fixed, 2 x swivel braking
- Maximum capacity 136.0kg per trolley


Doors will fold back flat with tow bar fitted

#### Dimensions (w x d x h) Standard Colours

- 1407 x 762 x 1288mm
- Stainless Steel Natural Unpolished

AGV Compatible Trolleys & Carts For Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy distribution

|         | Item          | Description   | Weight  |
|---------|---------------|---|---------|
| in      | STSS/690/UN   | Sterile Supplies Trolley - Stainless Steel, Trolley Washer Compatible         | 131.0kg |
| oplies, | STSS/690/UN/T | Sterile Supplies Trolley - Stainless Steel, Trolley Washer Compatible, Towing | 137.0kg |





High Level Corner Bumpers





# Storage Cupboards, Cabinets & Lockers

- Storage Cupboards 410
- Storage Cabinets Single & Double Door 412
- Storage Cabinets First Aid & Clinical Reagent Cabinets 413
  - Storage Cabinets Traka 21 414
    - Storage Cabinets Key 415
  - Storage Lockers Personal Effects 416











Shelves & Dividers

# Storage Cupboards - Single Door, 1870mm High

#### **Features**

- Single door 1870mm high cupboards
- Single-point locking using a flush fitting CAM lock fitted with a handle (supplied with 2 kevs)
- Available with either...
  - Plastic Travs
  - Shelves
  - Shelves & Dividers
- 0.9mm (20 Gauge) mild steel welded construction
- Modular design allows for easy interchange with HTM 63 storage systems

#### **Plastic Trays**

- Interchangeable with our popular Caretrav range
- Supplied with 7 x 100mm & 5 x 150mm travs
- · Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Travs feature an integral handle with a plastic clip-on label cover (coloured label strips) supplied - white, blue, green, vellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray

#### Shelves

- · Adjustable at 50mm increments
- Supplied with 8 x shelves

#### Shelves & Dividers

- · Adjustable at 50mm increments
- · Removable dividers 100mm high, adjustable at 25mm increments
- Supplied with 8 x shelves & 16 x dividers

#### Dimensions (w x d x h)

#### Cupboard

• 500 x 536 x 1870mm

Travs

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

#### Standard Colours

· Paint - White

| Options / | Accessories |
|-----------|-------------|
| optione,  |             |

- SC/JK/1870 Joining Kit
- SC/DIV Set of 16 Dividers
- CT/TRAY/S Shallow Trav

83.6kg

64.0kg

74.0kg

- CT/TRAY/D Deep Tray
- Weight **R/H Hinged** L/H Hinged Description SC/R/1870/7S5D SC/L/1870/7S5D Single Door, 1870mm High, Plastic Trays SC/R/1870/FS SC/L/1870/FS Single Door, 1870mm High, Flat Shelves SC/R/1870/FSD SC/L/1870/FSD Single Door, 1870mm High, Flat Shelves & Dividers

SC/L/1870/7S5D



**Joining Kits** 

#### **Features**

• Allows two units to be joined together

| Item       | Description |  |
|------------|-------------|--|
| SC/JK/1870 | Joining Kit |  |

# **Dividers**

#### **Features**

www.bristolmaid.com

- Spare / replacement dividers for SC cupboards supplied with shelves & dividers
- 100mm high dividers locate into slots on the shelves, adjustable at 25mm increments
- Compatible with 1870mm high units



| Item   | Description        |
|--------|--------------------|
| SC/DIV | Set of 16 Dividers |

SC/L/1870/FSD



# **Storage Cabinets - Single & Double Door**

#### Features

- Adjustable & removable shelves
- Flush fitting handle fitted with a CAM lock offering single-point locking (supplied with 2 keys)
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems

#### **Standard Colours**

• Paint - White





| Item Shelves Overall Dimension (w x d x h) |   | Weight             |        |
|--|---|--------------------|--------|
| Single Door                                |   |                    |        |
| SC020                                      | 2 | 600 x 300 x 600mm  | 19.5kg |
| SC025                                      | 3 | 600 x 300 x 900mm  | 26.0kg |
| Double Door                                |   |                    |        |
| SC040                                      | 4 | 1000 x 300 x 600mm | 32.5kg |
| SC045                                      | 6 | 1000 x 300 x 900mm | 45.0kg |

# **Storage Cabinets - First Aid & Clinical Reagents**

#### **Features**

- Removable shelves
- Flush fitting handle fitted with a CAM lock offering single-point locking (supplied with 2 keys)
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes in back (fixings not included)
- Supplied with 'First Aid' & 'Clinical Reagents' self-adhesive labels

#### Dimensions (w x d x h)

- Small 300 x 200 x 400mm
- Medium 400 x 200 x 400mm
- Large 400 x 200 x 600mm

#### **Standard Colours**

• Paint - White








| Item  | Description                                    | Shelves | Weight |
|-------|--|---------|--------|
| GC005 | First Aid & Clinical Reagents Cabinet - Small  | 2       | 5.0kg  |
| GC010 | First Aid & Clinical Reagents Cabinet - Medium | 2       | 6.0kg  |
| GC015 | First Aid & Clinical Reagents Cabinet - Large  | 2       | 8.0kg  |







# How does it work?



Access User enters unique PIN code to gain access to the cabinet



Kev Removal Keys & keysets are individually locked within the cabinet. Green LED's indicate those keys which the user can access. The user cannot remove keys highlighted by the Red LED's which continue to remain locked in place

sales@bristolmaid.com



Key Replacement An amber light indicates the return position for the key

# Storage Cabinets - Keys, Traka21

#### Features

- User, key & access rights administration
- Screen based audit trail & reporting capability or export via USB port
- 21 robust, long-life iFobs with security seals
- 21 locking receptors with LEDs
- Robust & compact aluminium & ABS housing with no exposed hinges or obvious access points
- Manual override & door release functions in case of emergency
- Audible alarms

#### Operation

- Stand-alone 'Plug & Play' solution with advanced RFID technology
- · Capacity for up to 21 keys
- Touchscreen interface
- · PIN access to designated keys or keysets which are individually locked in place
- · Keys are securely attached using special security seals
- Simple to configure through setup wizard
- No requirement for network connection or PC
- · Concealed robust fixing points allow secure wall fixing
- Mains operated (optional battery back-up)

#### Dimensions (w x d x h)

• 246 x 95 x 427mm

#### Standard Colours

· Aluminium & ABS - White

#### Options

- SKR/25 Security Key Rings (x 25)
- SKR/100 Security Key Rings (x 100)
- SKB/BAT Battery Back-Up



#### Traka21 iFob

Once attached to an iFob with a robust security seal, individual keys or keysets are assigned a specific receptor within the Traka cabinet & locked in place until released by an authorised user

| Item        | Description           | Weight |
|-------------|-----------------------|--------|
| SKC21/TRAKA | Key Cabinet - TRAKA21 | 4.0kg  |

# **Storage Cabinets - Keys**

#### Features

- · Cabinets for the control & tracking of keys
- Available with either a CAM or electronic push button lock
  - · CAM Supplied with two keys
  - · Electronic Push Button Battery operated
- Adjustable hook bars accommodate different key lengths
- Removable control index for out of hours peace of mind
- Colour coded hook bars & tags ensure guick identification. (both numbered)
- · Provided with security seals & loops to secure keys that require evidence of tampering or removal


#### Standard Colours

· Paint - Grev

#### **Electronic Push Button Lock**

- Uses a 10 button keypad and a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- · Four operating codes with differing levels of authority
  - Master code
  - Sub-master code
  - User code
  - Technician code
- · User codes are 4 digits
- · Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- · Low battery warning with failure override using 9v battery across LED's (+/-) & master code





SKC30

# 22 23 24 32 33 42 43 44 45 46

SKC100

Overall Dimensions (w x d x h)

185 x 110 x 260mm

215 x 110 x 305mm

300 x 110 x 355mm

300 x 110 x 450mm

380 x 110 x 550mm

Electronic Push Button Lock

#### **CAM Lock**

| Item   | Key Capacity | Overall Dimensions (w x d x h) | Weight | Item      | Key Capacity |
|--------|--------------|--------------------------------|--------|-----------|--------------|
| SKC20  | 20 Keys      | 185 x 80 x 260mm               | 2.0kg  | SKC20/CL  | 20 Keys      |
| SKC30  | 30 Keys      | 215 x 80 x 305mm               | 3.0kg  | SKC30/CL  | 30 Keys      |
| SKC48  | 48 Keys      | 300 x 80 x 355mm               | 4.0kg  | SKC48/CL  | 48 Keys      |
| SKC64  | 64 Keys      | 300 x 80 x 450mm               | 5.0kg  | SKC64/CL  | 64 Keys      |
| SKC100 | 100 Keys     | 380 x 80 x 550mm               | 8.0kg  | SKC100/CL | 100 Keys     |

Weight

2.0kg

3.0kg

4.0kg

5.0kg

8.0kg

Storage Cupboards,

Cabinets

Qo

Lockers

#### Alternative colours available

Replace suffix on the product code with any of the following to select your colour i.e. Single Bank with a Green Door - PL080/GD





GD Green Door



RD

Red Door

BD

Blue Door



YD Yellow Door

8 0



Hasp & Staple Lock

# Ø Ø 8 ø ø

# **Storage Lockers - Personal Effects**

#### **Features**

- Range of personal effects lockers available in single, double or triple banks in a number of compartment configurations
- · Choice of locking options
  - CAM lock (supplied with 2 keys)
  - Hasp & staple (accepts a padlock/not supplied)
- · Lock surround incorporates a 25mm round label holder
- Mounting holes provided in back (fixings not included)

0

0

8

· Optional stands & sloping tops

#### Dimensions (w x d x h)

# Internal Compartments (w x d x h)

- Single Bank 305 x 460 x 1780mm
- Double Bank 610 x 460 x 1780mm
- Triple Bank 915 x 460 x 1780mm

#### Standard Colours

- Body Grey
- · Doors Blue

0

0

0

0

0

0

A

- One 255 x 440 x 1680mm • Two - 255 x 440 x 815mm
- Three 255 x 440 x 525mm
- Four 255 x 440 x 385mm
- Five 255 x 440 x 295mm





Storage Cupboards, Cabinets & Lockers









PL145/YD

| Single Bank   |          |          |          | 2<br>2<br>2<br>2 |          |
|---------------|----------|----------|----------|------------------|----------|
| Compartments  | 1        | 2        | 3        | 4                | 5        |
| Coat Hook     | Yes      | Yes      | No       | No               | No       |
| Shelf         | Yes      | No       | No       | No               | No       |
| Hasp & Staple | PL005/BD | PL010/BD | PL015/BD | PL020/BD         | PL025/BD |
| CAM Lock      | PL080/BD | PL085/BD | PL090/BD | PL095/BD         | PL100/BD |





# Stand

#### **Features**

- Compatible with single, double & triple bank 1780mm high personal effects lockers
- Keeps lockers clear of the floor assisting cleaning & aiding air-circulation
- Fitted with adjustable feet
- Increases the overall height of the locker by 150mm

#### **Standard Colours**





PL915 - Sloping top

| Item  | Description         |  |
|-------|---------------------|--|
| PL900 | Stand - Single Bank |  |
| PL905 | Stand - Double Bank |  |
| PL910 | Stand - Triple Bank |  |

# **Sloping Top**

### Features

- Compatible with single, double & triple bank 1780mm high personal effects lockers
- Prevents items being stored on the top of the locker & assists cleaning
- Increases the overall height of the locker by 150mm

#### **Standard Colours**

Sloping Top - Grey

| ltem  | Description               |
|-------|---------------------------|
| PL915 | Sloping Top - Single Bank |
| PL920 | Sloping Top - Double Bank |
| PL925 | Sloping Top - Triple Bank |



sales@bristolmaid.com



# **Theatre Equipment**

- Bowl Stands 420
- Bucket Stands 423
- Instrument Trolleys 424
  - Mayo Tables 426
- Operating Platforms 426
- Suture Storage Trolleys 427
- Sterile Container Storage 428



# **Bowl Stands - Fixed Height**

#### **Features**

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Bowl retaining rings fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 75mm pressed steel non-marking castors (1 x anti-static)

#### **Standard Colours**

Stainless Steel

#### **Options / Accessories**

- 2MECO079 360mm Diameter, 7 Litre Stainless Steel Bowl
- 2PLCO054 360mm Diameter, 7 Litre Polythene Bowl



2MECO079 7 Litre Stainless Steel Bowl



2PLCO054 7 Litre Polythene Bowl

| ltem   | Description  | Overall Dimensions (w x d x h) | Weight |
|--------|--|--------------------------------|--------|
| SBS/2  | Bowl Stand - Fixed Height, Double Bowl, Side by Side | 750 x 445 x 860mm              | 6.0kg  |
| SBST/2 | Bowl Stand - Fixed Height, Double Bowl, Tiered       | 445 x 445 x 860mm              | 5.0kg  |

www.bristolmaid.com

# **Bowl Stands - Variable Height**

#### **Features**

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Height adjustable upright using a handwheel (753 1100mm)
- Bowl retaining ring fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 50mm polypropylene/thermoplastic rubber twin-wheel castors (1 x anti-static)

#### **Standard Colours**

Stainless Steel

#### **Options / Accessories**

- 2MECO079 360mm Diameter, 7 Litre Stainless Steel Bowl
- 2PLCO054 360mm Diameter, 7 Litre Polythene Bowl



2MECO079 7 Litre Stainless Steel Bowl



2PLCO054 7 Litre Polythene Bowl

| Item   | Description                               | Overall Dimensions (w x d x h) | Weight |
|--------|---|--------------------------------|--------|
| SBSA/1 | Bowl Stand - Variable Height, Single Bowl | 650 x 650 x 753 - 1100mm       | 4.8kg  |
| SBSA/2 | Bowl Stand - Variable Height, Double Bowl | 746 x 650 x 753 - 1100mm       | 5.6kg  |



# Bowl Stands - Variable Height, Polyurethane, Easy Clean



| Item       | Description   | Weight |
|------------|---|--------|
| BO/P/VS/SP | Bowl Stand - Variable Height, Polyurethane, Easy Clean, Single Bowl | 9.0kg  |
| BO/P/VD/SP | Bowl Stand - Variable Height, Polyurethane, Easy Clean, Double Bowl | 9.8kg  |

#### Features

- · Five leg base with low centre of gravity to assist stability
- 25mm diameter scratch resistant stainless steel upright
- Polyurethane overmoulded base with smooth edges & contours provides...
  - Water-resistant finish, allowing for easy & effective cleaning (including steam cleaning)
  - Shield for the twin-wheel castors
  - Impact protection to patients & staff, reducing the risk of skin abrasions & bruising
  - · Protection to the fabric of the building
- Coloured indicator caps provide department identification & segregation, reducing the risk of cross contamination
- Height adjustable upright using a handwheel (820 1125mm)
- Bowl retaining rings fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 50mm polypropylene/thermoplastic rubber twin-wheel non-marking castors (1 x anti-static)

#### Dimensions (w x d x h)

- Single Bowl 620 x 620 x 820 1125mm
- Double Bowl 746 x 620 x 820 1125mm

#### **Standard Colours**

- Base Grey
- Coloured Indicator Cap Spearmint
- Bowl Retaining Rings Stainless Steel

#### **Options / Accessories**

- 2MECO079 Stainless Steel Bowl
- 2PLCO054 Polythene Bowl
- Coloured Indicator Cap Blue / Green / Red / Yellow



# **Bowl Stand - Kickabout, 360mm Diameter**

#### **Features**

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- · Bowl retaining ring fitted with removable anti-static spiral buffers to support bowls
- · Supplied with a 7 litre, 360mm diameter stainless steel bowl
- 75mm pressed steel non-marking castors (1 x anti-static, 2 x non-anti-static)

#### Dimensions (w x d x h)

#### Standard Colours

• 365 x 365 x 235mm

#### Stainless Steel

#### **Options / Accessories**

- 2MECO079 360mm Diameter, Stainless Steel Bowl
- 2PLCO054 360mm Diameter, Polythene Bowl



| Item   | Description                            | Weight |
|--------|--|--------|
| SBSL/1 | Bowl Stand - Kickabout, 360mm Diameter | 3.2kg  |

# **Bucket Stands - Kickabout, 15-Litre Bucket**

#### **Features**

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Bucket retaining ring fitted with removable anti-static spiral buffers to support the bucket
- Model available fitted with a handle
- Supplied with a 15-litre graduated stainless steel bucket
- 75mm pressed steel non-marking castors (1 x anti-static, 2 x non-anti-static)

#### Dimensions (w x d x h)

#### Standard Colours

- SBKSL/1 460 x 460 x 470mm Stainless Steel SBKSI H/1- 460 x 460 x 470mm
- Supplied with a 15-litre
  - stainless steel bucket

#### **Options / Accessories**

• 2MECO081 - 15-litre graduated stainless steel bucket



| Item     | Description                            | Weight |
|----------|--|--------|
| SBKSL/1  | Bucket Stand - 15 Litre Bucket         | 5.0kg  |
| SBKSLH/1 | Bucket Stand - 15 Litre Bucket, Handle | 5.5kg  |

Supplied with a

7 litre stainless

steel bowl

+44 (0)1258 484455



# Fixed Shelves, Sides Down, Stainless Steel, 450mm Deep

#### **Features**

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

#### **Options / Accessories**

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 75mm, Pressed Steel Castors (4 x braking)
- Drawer Units page 119
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

www.bristolmaid.com

# Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep

#### **Features**

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel anti-static non-marking castors

#### **Options / Accessories**

- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C015 100mm, Pressed Steel Castors (4 x braking)
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

| Item          | Size   | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|---------------|--------|--------------------|--------------------------------|--------|
| DTSF/66/2/SD  | Small  | 600 x 600mm        | 620 x 620 x 870mm              | 17.5kg |
| DTSF/96/2/SD  | Medium | 900 x 600mm        | 920 x 620 x 870mm              | 24.0kg |
| DTSF/126/2/SD | Large  | 1200 x 600mm       | 1220 x 620 x 870mm             | 29.6kg |



# Mayo Tables - Cantilever, Variable Height


#### **Features**

- · Seamlessly welded frame from 304 grade stainless steel for strength & ease of cleaning
- · Available with two castors & skid feet or four castors
- Height adjustment...
  - SS/22/1 893 1480mm
  - SS/22/2 858 1450mm
- Supplied with a removable stainless steel tray, 512 x 387 x 20mm (w x d x h)
- 50mm polypropylene/thermoplastic rubber twin-wheel anti-static castors

#### Dimensions (w x d x h)

- SS/22/1 641 x 450 x 893-1480mm

Supplied with a stainless



# **Operating Platforms**

#### **Features**

- Seamlessly welded frame from 304 grade stainless steel for strength & ease of cleaning
- Aluminium embossed non-slip tread
- · Non-slip, anti-static rubber feet
- Max capacity 190kg

#### Dimensions (w x d x h)

- Small 400 x 300mm
- Large 600 x 300mm



Aluminium embossed non-slip tread



| ltem   | Description                            | Weight |
|--------|--|--------|
| Small  |  |        |
| OP/431 | Operating Platform - Small, 150mm High | 4.2kg  |
| OP/432 | Operating Platform - Small, 200mm High | 4.4kg  |
| Large  |  |        |
| OP/631 | Operating Platform - Large, 150mm High | 6.0kg  |
| OP/632 | Operating Platform - Large, 200mm High | 6.4kg  |



| Item    | Description   | Weight |
|---------|---|--------|
| SS/22/  | Mayo Table - Cantilever, Variable Height, Two Castors & Skid Feet | 5.4kg  |
| SS/22/2 | 2 Mayo Table - Cantilever, Variable Height, Four Castors          | 5.8kg  |

sales@bristolmaid.com

# **Suture Storage Trolleys**

#### **Features**

- Angled shelves present contents for instant identification & easy access
- Allows logical organisation, eliminates confusion, allowing for easier & faster restocking
- Open wire design minimises dust accumulation, allowing for free circulation of air, improved visibility & light penetration
- No tools shelf adjustment in 25mm increments
- Each trolley supplied with 12 x shelf retainers to assist with segregation of contents
- Delivered flat pack for self-assembly
- 125mm pressed steel castors, 2 off braking

#### **Options / Accessories**

- 89990P Label Holder, 75 x 32mm
- 8DCR17C Shelf Retainers



| Item               | Description   | Overall Dimensions (w x d x h) | Weight |  |
|--------------------|---|--------------------------------|--------|--|
| 1 x Flat, 4        | x Angled Shelves  |                                |        |  |
| 8DC15EC            | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 610mm  | 610 x 457 x 1524mm             | 46.5kg |  |
| 8DC35EC            | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 914mm  | 914 x 457 x 1524mm             | 49.0kg |  |
| 8DC55EC            | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 1219mm | 1219 x 457 x 1524mm            | 51.0kg |  |
| 5 x Angled Shelves |   |                                |        |  |
| 8DC16EC            | Trolley - Suture Storage, 5 x Angled Shelves, 610mm             | 610 x 457 x 1524mm             | 50.0kg |  |
| 8DC36EC            | Trolley - Suture Storage, 5 x Angled Shelves, 914mm             | 914 x 457 x 1778mm             | 51.5kg |  |
| 8DC56EC            | Trolley - Suture Storage, 5 x Angled Shelves, 1219mm            | 1219 x 457 x 1778mm            | 52.5kg |  |

+44 (0)1258 484455



# **Sterile Container Storage**

#### **Features**

- Fully configurable freestanding modular racking suitable for the safe storage of sterile containers
- Compatible with both blue wrap & rigid containers, up to max 300mm x 600mm (width x depth)
- Containers stored on either...
  - Telescopic shelves
  - Wire grid shelves
- Anodised extruded aluminium uprights & cross members
- · Bolted construction using high quality moulded fixings provided rigidity & allows the system to be altered easily
- Available as either static or mobile units
  - Static Adjustable feet allow levelling on uneven floors
  - Mobile 50 / 80mm pressed steel castors
- Optional wall mounting brackets allow static units to be fixed securely to the wall, (when units are loaded over 100kg)
- Supplied for self-assembly

#### Dimensions (w x d x h)

Static

- Single Bay 394 x 642 x 1940mm
- Double Bay 737 x 642 x 1940mm

Mobile (50/80mm Castors)

- Single Bay 394 x 642 x 1996 / 2035mm
- Double Bay 737 x 642 x 1996 / 2035mm

#### **Standard Colours**

- Uprights & Cross Members Anodised Aluminium
- Wire Grids Polished Stainless Steel
- Telescopic Shelves Anodised Aluminium
- Plastic Mouldings Green / Grey

#### **Options / Accessories**

- 2HH-7532293 Wall Mounting Bracket (x1)
- 2HH-7532277 10mm Plastic Cover Cap (x1)

#### Static

| ltem             | Description   | Weight |
|------------------|---|--------|
| Single Bay       |   |        |
| 4HH-SCS/S1/2W5T  | Static, Single Bay, 2 x Wire Grid Shelves & 5 x Telescopic Shelves  | 30.2kg |
| 4HH-SCS/S1/8W    | Static, Single Bay, 8 x Wire Grid Shelves                           | 21.4kg |
| Double Bay       |   |        |
| 4HH-SCS/S2/4W10T | Static, Double Bay, 4 x Wire Grid Shelves & 10 x Telescopic Shelves | 56.9kg |
| 4HH-SCS/S2/16W   | Static, Double Bay, 16 x Wire Grid Shelves                          | 38.0kg |
|                  |   |        |

#### Mobile

4

Λ

| ltem              | Description   | Weight |
|-------------------|---|--------|
| Single Bay        |   |        |
| 4HH-SCS/M51/2W5T  | Mobile, Single Bay, 2 x Wire Grid Shelves & 5 x Telescopic Shelves  | 31.1kg |
| 4HH-SCS/M51/8W    | Mobile, Single Bay, 8 x Wire Grid Shelves                           | 22.3kg |
| Double Bay        |   |        |
| 4HH-SCS/M52/4W10T | Mobile, Double Bay, 4 x Wire Grid Shelves & 10 x Telescopic Shelves | 58.2kg |
| 4HH-SCS/M52/16W   | Mobile, Double Bay, 16 x Wire Grid Shelves                          | 39.3kg |

sales@bristolmaid.com +44 (0)1258 484455



2HH-7532255 - Wire Grid Shelf



2HH-7532335 - Telescopic Shelves



50mm Braking Pressed Steel Castor



Adjustable Feet



4HH-SCS/M51/8W Shown with BBraun Aesculap Sterile Containers (not supplied)

4HH-SCS/S2/16W Shown with BBraun Aesculap Sterile Containers (not supplied)

**Components** 

- 2HH-7532300-HH Base Unit (x 1)
- 2HH-7532310 Extension Unit (x 1)
- 2HH-7532335 Telescopic Shelf (x 1)
- 2HH-7532255 Wire Grid Shelf (x 1)
- 2HH-7532356 Runners (Pair)
- 2HH-7532320 Adjustable Foot Set (x 4)
- 2HH-7532321 Adjustable Foot Set (x 2) • 2HH-7532322 - 50mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532323 50mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)
- 2HH-7532324 80mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532325 80mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)

Room Design & Installation Service Available Contact us for more information

sales@bristolmaid.com



# Wireless Mobile Carts & Workstations

- Mobile Workstations Non-Powered 432
  - Options & Accessories 434
  - Caretray Trolley Laptop 436



# **Mobile Workstations - Non-Powered**

#### Features

- Range of non-powered mobile workstations
- Models available for use with laptops & vesa display mounts & configurable with a keyboard tray
  - Laptop Mount Secures the laptop, (compatible with laptops up to 482mm in width)
  - Vesa Display Mount Height adjustable, compatible with 75 x 75mm / 100 x 100mm square hole patterns
  - Keyboard Tray Retractable tray with wrist support & left / right mouse pad (accepts full size & Medigenic keyboards with room for mouse storage)
- Height adjustable using a gas damper (780 1130mm)
- Extruded aluminium column provides a smooth linear height adjustment
- Safety feature prevents trapping of the user beneath the workstation
- · Cast aluminium four leg base
- · Acrylic work surface
- Configurable with a range of accessories
- 100mm twin-wheel castors (4 x braking)

#### Dimensions (w x d x h)

- Overall 630 x 630 x 780 1130mm
- Work Surface 500 x 465mm
- Keyboard Tray 500 x 255mm
- · Laptop Mount Accepts laptops up to 482mm

#### **Standard Colours**

- Base White
- Work Surface Blue
- Upright & Sheet White

| Item            | Description   | Weight |
|-----------------|---|--------|
| 5RC/MW/NP       | Mobile Workstation - Non-Powered                                    | 19.6kg |
| 5RC/MW/NP/LM    | Mobile Workstation - Non-Powered, Laptop Mount                      | 20.2kg |
| 5RC/MW/NP/VM    | Mobile Workstation - Non-Powered, Vesa Display Mount                | 22.4kg |
| 5RC/MW/NP/VM/KT | Mobile Workstation - Non-Powered, Vesa Display Mount, Keyboard Tray | 23.6kg |

#### **Options / Accessories**

- 5RC/SB Storage Basket
- 5RC/DRK Side Din Rail Kit
- 5RC/PSH Power Supply Unit Holder
- GMBH/15 Glove / Mask Box Holder
- CT/ALC Alcohol Gel Holder
- SBH/15 Sharps Box Holder



# **Storage Basket**

#### **Features**

- · Wire storage basket
- Mounts to the aluminium extrusion

#### **Standard Colours**

• Paint - Grey White

#### Dimensions (w x d x h)

• 380 x 180 x 206mm



| Item   | Description    | Weight |
|--------|----------------|--------|
| 5RC/SB | Storage Basket | 1.0kg  |

# **Power Supply Unit Holder**

#### Features

- Provides a location to store the power supply unit
- Mounts to the underside of the top

#### **Standard Colours**

White

# Side Din Rail Kit

#### **Features**

- Aluminium DIN rail
- Provides a rail measuring 10 x 25mm to hang / position accessories
- Mounts to the right / left side of the unit

#### **Standard Colours**

• Aluminium / White

#### Dimensions (w x d x h)

• 300 x 10 x 25mm





| Item    | Description       | Weight |
|---------|-------------------|--------|
| 5RC/DRK | Side Din Rail Kit | 0.6kg  |

Wireless Mobile Carts & Workstations

# **Glove / Mask Box Holder, One Box**

#### Features

- Accommodates glove & mask boxes up to 150mm in width
- Box presented in a portrait orientation
- Fitted with a 15mm hook

#### **Standard Colours**

Paint - Grey White

#### Dimensions (w x d x h)

• 143 x 122 x 177mm

# Alcohol Gel Holder

#### Features

- Accepts bottles with 90 x 90mm square base
- Fitted with a 15mm hook

#### **Standard Colours**

Paint - Grey White

#### Dimensions (w x d x h)

• 95 x 110 x 115mm

| _ |        |                    |        |
|---|--------|--------------------|--------|
| t | Item   | Description        | Weight |
|   | CT/ALC | Alcohol Gel Holder | 0.3kg  |

# **Sharps Box Holder**

### **Features**

- Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a toggle
- Fitted with a 15mm hook

# **Standard Colours**

· Paint - Grey White

### Dimensions (w x d x h)

• 127 x 118 x 230mm

Item

**SBH/15** 



Description

Sharps Box Holder

| _          |
|------------|
| <u> </u>   |
|            |
| 2          |
| <u> </u>   |
| Wireless   |
| Ś          |
| S          |
| $\leq$     |
|            |
| <u> </u>   |
| <u>o</u>   |
|            |
| (P         |
| C          |
| ۵Ú         |
| 3          |
| arts       |
|            |
| <b>Š</b> 0 |
| <          |
| 8          |
| <u> </u>   |
| <b>_</b>   |
| St         |
| at .       |
| 2          |
| 5          |
| 2          |
| 2          |
|            |



| Item    | Description                      | Weight |
|---------|----------------------------------|--------|
| GMBH/15 | Glove / Mask Box Holder, One Box | 0.6kg  |

| sales@bristolmaid.com   |  |
|-------------------------|--|
| sales@pristoirrialu.com |  |

Weight

0.5kg



| Item                    | Description   | Weight |
|-------------------------|---|--------|
| High Security Bolt Lock |   |        |
| CTLAPTM1/10/3S2D/HSB    | 3 x 100mm, 2 x 150mm Trays, High Security Bolt Lock | 39.3kg |
| CTLAPTM1/10/4S1D/HSB    | 4 x 100mm, 1 x 150mm Trays, High Security Bolt Lock | 39.0kg |
| CTLAPTM1/10/6S/HSB      | 6 x 100mm Trays, High Security Bolt Lock            | 40.1kg |

# **Caretray Trolley - Laptop**

#### Features

- Trolleys designed for use with laptops
- Quick & easy disassembly, assists thorough & deep cleaning
- Fully welded mild steel frame, powder coated in a durable epoxy paint (requires no on-site assembly)
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
  - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
  - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerails fitted to either side support a wide range of accessories
- Solid grade laminate door & back panel providing security for the contents of the trays
- Two locking options...
  - High Security Bolt Lock Supplied with two keys
  - Electronic Push Button Lock Programable four-digit code
- Trolley pre-configured with the following accessories...
  - Laptop Holder, Removable, Swing Arm
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- 100 / 150mm deep trays, can be subdivided in both length & width using removable dividers (max load 5kg)
- 2 x long & 2 x short dividers supplied with each tray
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- 100mm polymer castors, 2 off braking
- When not in use, trolleys storing medicines should be secured to the wall using a
   Security Clamp (see locking & security accessories)

### Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Plastic Top 455 x 455 (internal) at height 968mm
- Trays
  - Shallow Tray 355 x 405 x 100mm (internal)
  - Deep Tray 355 x 405 x 150mm (internal)

#### **Standard Colours**

· Paint - Grey White

#### **Options / Accessories**

- Defibrillator Platform page 67
- Carerail Mounted page 68
   Laptop Holders page 74

- Castors & Ring Buffers page 76
   Trays, Shelves, Dividers & Labels page 77



| Item                        | Description   | Weight |
|-----------------------------|---|--------|
| Electronic Push Button Lock |   |        |
| CTLAPTM1/10/3S2D/EPB        | 3 x 100mm, 2 x 150mm Trays, Electronic Push Button Lock | 39.3kg |
| CTLAPTM1/10/4S1D/EPB        | 4 x 100mm, 1 x 150mm Trays, Electronic Push Button Lock | 39.0kg |
| CTLAPTM1/10/6S/EPB          | 6 x 100mm Trays, Electronic Push Button Lock            | 40.1kg |





# Workstations & Speciality Storage

- Workstations 452
- Storage Units 453
- Catheter Storage Units 454





Other configurations available. Contact our customer services team to discuss your requirements

# Workstations - Mobile, Polymer, Metro Starsys™

# **Features**

- · Adaptable workstations
- Available as single, double & triple width units
- Moulded body, drawers & doors from a lightweight polymer that will not rust, chip, dent, peel or corrode
- Laminate worktop (height 914mm)
- 125mm polymer castors, 2 off braking

# **Standard Colours**

- Body / Drawers / Doors Light Taupe
- Drawer / Door Handles Slate Blue





+44 (0)1258 484455

Item

# Storage Units - Mobile, Polymer, Metro Starsys™

# Features

- Full height storage units
- Available as single, double & triple width units
- Moulded body, drawers & hinged/tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode
- 125mm polymer castors, 2 off braking

# **Standard Colours**

# **Options / Accessories**

• 8SXREHAN - Push Handle

• 8SXLABKIT - Drawer / Door Label Holders (x 10)

- Body / Drawers / Doors Light Taupe
- Drawer / Door Handles Slate Blue





| Item                   | Description                    | Overall Dimensions (w x d x h) | Weight  |  |  |  |
|------------------------|--------------------------------|--------------------------------|---------|--|--|--|
| Drawers & Hinged Doors |                                |                                |         |  |  |  |
| 8SXRSGS2               | Single, Drawers                | 578 x 632 x 1990mm             | 187.7kg |  |  |  |
| 8SXRD76CMHD1           | Double, Drawers & Hinged Doors | 1080 x 709 x 1990mm            | 159.2kg |  |  |  |
| 8SXRTGS2               | Triple, Drawers & Hinged Doors | 1669 x 709 x 1990mm            | 187.8kg |  |  |  |
| Hinged Doors           |                                |                                |         |  |  |  |
| 8SXRS76CM5             | Single, Door                   | 578 x 709 x 1990mm             | 106.0kg |  |  |  |
| 8SXRD76CM3             | Double, Doors                  | 1080 x 709 x 1990mm            | 113.0kg |  |  |  |
| 8SXRTGS3               | Triple, Doors                  | 1669 x 709 x 1990mm            | 158.0kg |  |  |  |
| Tambour Doors          |                                |                                |         |  |  |  |
| 8SXRS83MGSTD           | Single, Tambour Door           | 578 x 709 x 2096mm             | 108.9kg |  |  |  |
| 8SXRD83MGSTD           | Double, Tambour Doors          | 1080 x 709 x 2096mm            | 152.9kg |  |  |  |





sales@bristolmaid.com



# Catheter Storage Units - Mobile, Polymer, Metro Starsys™





# **Features**

- Full height catheter storage units available with either hinged or tambour door/s
- Catheters supported on retractable slides providing easy to access & safe storage
- Available as single & double width units
  - Single 180 Catheters
  - Double 360 Catheters
- Moulded body & hinged or tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode
- 125mm polymer castors, 2 off braking

# **Standard Colours**

- Body / Drawers / Doors Light Taupe
- Drawer / Door Handles Slate Blue

## **Options / Accessories**

- 8SXRCATHHK Catheter Hooks, (supports 12 catheters)
- 8SXRLABKIT Drawer / Door Label Holders, (x 10)
- 8SXREHAN Push Handle

Other configurations available. Contact our customer services team to discuss your requirements

| Item           | Description                                    | Overall Dimensions (w x d x h) | Weight  |
|----------------|--|--------------------------------|---------|
| Hinged Doors   |  |                                |         |
| 8SXRS76MCATH1  | Single, L/H Hinged Door, 180 x Catheters       | 578 x 632 x 1990mm             | 100.1kg |
| 8SXRD76MCATH1  | Double, R/H & L/H Hinged Door, 360 x Catheters | 1080 x 632 x 1990mm            | 136.3kg |
| Tambour Doors  |  |                                |         |
| 8SXRS83MCATHTD | Single, Tambour Door, 180 x Catheters          | 578 x 709 x 2096mm             | 109.3kg |
| 8SXRD83MCATHTD | Double, Tambour Doors, 360 x Catheters         | 1080 x 709 x 2096mm            | 146.1kg |

# Catheter Storage Units - Stationary, Polymer, Metro Starsys™

# **Features**

- · Full height catheter storage units available with either hinged or tambour door/s
- Catheters supported on retractable slides providing easy to access & safe storage
- · Available as single & double width units
  - Single 180 Catheters
  - Double 360 Catheters
- Moulded body & hinged or tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode

# **Standard Colours**

- Body / Drawers / Doors Light Taupe
- Drawer / Door Handles Slate Blue

# **Options / Accessories**

- 8SXRCATHHK Catheter Hooks, (supports 12 catheters)
- 8SXRLABKIT Drawer / Door Label Holders, (x 10)
- 8SXR-FLR Floor Mounting Kit
- 8SXR72UUFLR Unit to Unit Filler Kit\*\*
- 8SXR72BKFLR30 Back to Wall Kit (Half Depth Units)\*\*\*
- 8SXR72BKFLR Back to Wall Kit (Full Depth Units)\*\*\*

\*\*Covers the seam between two adjacent units

\*\*\*Fills the gap when a unit it is pulled away from the wall to be flush with the front of the worktop



| Item           | Description                                       | Overall Dimensions (w x d x h) | Weight  |
|----------------|---|--------------------------------|---------|
| Hinged Doors   |   |                                |         |
| 8SXRS72SCATH1  | Single, L/H Hinged Door, 180 x Catheters          | 540 x 594 x 1842mm             | 75.5kg  |
| 8SXRD72SCATH1  | Double, R/H & L/H Hinged Door, 360 x Catheters    | 1049 x 594 x 1842mm            | 143.6kg |
| Tambour Doors  |   |                                |         |
| 8SXRS80HCATHTD | Single, Half Depth, Tambour Door, 180 x Catheters | 540 x 457 x 2032mm             | 60.3kg  |
| 8SXRD80HCATHTD | Double, Half Depth, Tambour Door, 360 x Catheters | 360 1052 x 457 x 2032mm        | 113.9kg |
| 8SXRS80SCATHTD | Single, Full Depth, Tambour Door, 180 x Catheters | 180 540 x 620 x 2032mm         | 84.4kg  |
| 8SXRD80SCATHTD | Double, Full Depth, Tambour Door, 360 x Catheters | 1052 x 620 x 2032mm            | 140.2kg |







| 05400004                     | 07         | 4A1RS12            |     |                        |            | 4HH-HTM71/M51/600                    | 205        | 5RC/MW/NP                                     | 432        |
|------------------------------|------------|--------------------|-----|------------------------|------------|--------------------------------------|------------|---|------------|
| 2FACO024                     |            | 4A1RS15            |     | 4B1TD10/D              | 200        | 4HH-HTM71/M52/400                    | 205        | 5RC/MW/NP/LM                                  | 432        |
|                              | 209        | 4A1RS18<br>4A2DR09 |     | 4B1TD10/E<br>4B1TD10/F | 200<br>200 | 4HH-HTM71/M52/600                    | 205        | 5RC/MW/NP/VM<br>5RC/MW/NP/VM/KT               | 432<br>432 |
|                              | 209<br>207 | 4A2DR09<br>4A2DR10 |     | 4B11D10/F<br>4B11D10/G |            | 4HH-HTM71/S1/400<br>4HH-HTM71/S1/600 | 204<br>204 | 5RC/MW/NF/VM/KT                               | 432        |
| 2HH-240820<br>2HH-240820     |            | 4A2DR10<br>4A2DR12 |     | 4B1TD10/G<br>4B1TD12/H |            | 4HH-HTM71/S1/600<br>4HH-HTM71/S2/400 | 204        | 5RC/SB  | 434        |
|                              | 200        | 4A2DR12            |     | 4B1TD12/I              |            | 4HH-HTM71/S2/600                     | 204        | 5RE/OC/PORT/VVDEL                             | 358        |
|                              | 208        | 4A2DR18            |     | 4B1TD12/J              |            | 4HH-SCS/M51/2W5T                     | 428        | 5RE/OC/PORT/VVLAVSAVL                         | 358        |
|                              | 207        | 4A2DSTM            |     | 4B1VD09/A              | 199        | 4HH-SCS/M51/8W                       | 428        | 5RE/RAE/ALWE/VVDEL                            | 342        |
| 2HH-241301                   |            | 4A2DSTS            |     | 4B1VD09/B              | 199        | 4HH-SCS/M52/16W                      | 428        | 5RE/RAE/ALWE/VVDEL/VPNB                       | 343        |
| 2HH-241302                   |            | 4A2OF09            |     | 4B1VD09/C              | 199        | 4HH-SCS/M52/4W10T                    | 428        | 5RE/RAE/ASHI/VVDEL                            | 346        |
| 2HH-241302                   |            | 4A2OF10            |     | 4B1VD10/D              | 199        | 4HH-SCS/S1/2W5T                      | 428        | 5RE/RAE/ASHI/VVDEL/VPNB                       | 347        |
|                              | 206        | 4A2OF12            | 196 | 4B1VD10/E              | 199        | 4HH-SCS/S1/8W                        | 428        | 5RE/RAE/BECK/VMULT                            | 348        |
| 2HH-7011100                  | 206        | 4A2OF15            | 196 | 4B1VD10/F              | 199        | 4HH-SCS/S2/16W                       | 428        | 5RE/RAE/SHAP/VVDEL                            | 340        |
| 2HH-7012100                  | 206        | 4A2OF18            |     | 4B1VD12/G              | 199        | 4HH-SCS/S2/4W10T                     | 428        | 5RE/RAE/SHAP/VVDEL/VPNB                       | 341        |
| 2HH-7012200                  | 206        | 4A2OTTM            |     | 4B1VD12/H              |            | 5160                                 | 331        | 5RE/RAE/TWYF/VVDEL                            | 344        |
| 2HH-7014050                  | 206        | 4A2OTTS            |     | 4B1VD15/I              | 199        | 5161                                 | 331        | 5RE/RAE/TWYF/VVDEL/VPNB                       | 345        |
| 2HH-7014100                  |            | 4A2RS09            |     | 4B1VD15/J              | 199        | 5162                                 | 331        | 5RE/RAM/ALWE/VVDEL                            |            |
| 2HH-7014200                  |            | 4A2RS10            |     | 4B2DR09                | 198        | 5163                                 | 331        | 5RE/RAM/ALWE/VVDEL/VPNB                       | 343        |
| 2HH-7810102                  |            | 4A2RS12            |     | 4B2DR10                |            | 5201                                 | 337        | 5RE/RAM/ASHI/VVDEL                            |            |
| 2HH-7810103                  |            | 4A2RS15            |     |                        |            | 5202                                 | 336        | 5RE/RAM/ASHI/VVDEL/VPNB                       | 347        |
| 2HH-7810104                  |            | 4A2RS18            |     | 4B2DR15                | 198        | 5205                                 | 337        | 5RE/RAM/SHAP/VVDEL                            |            |
|                              |            | 4A3DR09            |     | 4B2DR18<br>4B2OF09     |            | 5260<br>5262                         | 332<br>332 | 5RE/RAM/SHAP/VVDEL/VPNB<br>5RE/RAM/TWYF/VVDEL | 341<br>344 |
| 2HH-7811204N                 | 208        | 4A3DR10<br>4A3DR12 |     | 4B2OF09<br>4B2OF10     |            | 5262                                 | 363        | 5RE/RAM/TWYF/VVDEL/VPNB                       | 344        |
| 2HH-7811205N<br>2HH-7811205N | 207<br>208 | 4A3DR12<br>4A3DR15 |     | 4B2OF10<br>4B2OF12     |            | 5671                                 | 14         | 5RE/RPAH/BATH/VVDEL                           |            |
| 2HH-7811205N                 |            | 4A3DR18            |     | 4B2OF12<br>4B2OF15     |            | 5671                                 | 335        | 5RE/RPAH/BATH/VVDEL/VPNB                      | 353        |
| 2HH-7811304                  | 207        | 4A30F09            |     | 4B2OF18                |            | 59003                                | 333        | 5RE/RRAE/ALWE/VVDEL                           | 342        |
|                              | 207        | 4A3OF10            |     | 4B2RS09                |            | 5160D                                | 331        | 5RE/RRAE/ALWE/VVDEL/VPNB                      | 343        |
|                              | 208        | 4A3OF12            |     |                        |            | 5301/A                               |            | 5RE/RRAE/ASHI/VVDEL                           |            |
| 2HH-7811306                  |            | 4A3OF15            |     | 4B2RS12                |            | 5660XL                               |            | 5RE/RRAE/ASHI/VVDEL/VPNB                      | 347        |
| 2HH-7910201N                 | 209        | 4A3OF18            | 196 | 4B2RS15                | 197        | 5660XL                               | 334        | 5RE/RRAE/BECK/VMULT                           | 348        |
| 2HH-7910204N                 | 209        | 4A3OTTM            | 203 | 4B2RS18                | 197        | 5FSH20WH                             | 326        | 5RE/RRAE/HADL610/VVDEL/VPNB                   | 350        |
| 2HH-K100187                  | 209        | 4A3OTTS            |     | 4B2TD09/AA *           | 200        | 5HORI/1/A                            | 365        | 5RE/RRAE/HADL710/VVDEL/VPNB                   | 350        |
|                              | 206        | 4A3RS09            |     | 4B2VD09/AA*            | 199        | 5HORI/2/A                            | 365        | 5RE/RRAE/SHAP/VVDEL                           | 340        |
| 2LOCK217                     |            | 4A3RS10            |     | 4B2VD10/DD *           | 199        | 5HORI/3/A                            | 365        | 5RE/RRAE/SHAP/VVDEL/VPNB                      | 341        |
| 2LOCK218                     | 273        | 4A3RS12            |     | 4B3DR09                | 198        | 5ND/HAF/ALDE/VVDEL                   | 339        | 5RE/RRAE/TWYF/VVDEL                           | 344        |
| 2LOKE011                     |            | 4A3RS15            |     | 4B3DR10                | 198        | 5ND/HAR/ALDE/VVDEL                   |            | 5RE/RRAE/TWYF/VVDEL/VPNB                      | 345        |
| 2MECO079                     |            | 4A3RS18            |     | 4B3DR12                |            | 5ND/LC/ALDE/VVDEL                    |            | 5RE/RRPAH/DORC/VVDEL                          | 354        |
| 2PLC0026                     | 65         | 4B1DR09            |     | 4B3DR15                | 198        | 5ND/MAF/ALDE/VVDEL                   |            | 5RE/RRPAH/DORC/VVDEL/VPNB                     | 355        |
|                              | 308        | 4B1DR10            |     | 4B3DR18                |            | 5ND/MAR/ALDE/VVDEL                   | 339        | 5RLDF0119                                     | 165        |
|                              | 420        | 4B1DR12<br>4B1DR15 |     | 4B3OF09<br>4B3OF10     | 196<br>196 | 5NOVA/1<br>5NOVA/2                   | 364<br>364 | 5RLDF0119DIGILOCK<br>5RLDF0219                | 165<br>165 |
| 4A1DR09<br>4A1DR10           |            | 4B1DR15<br>4B1DR18 |     | 4B3OF10<br>4B3OF12     |            | 5PF332                               | 367        | 5RLDF0219<br>5RLDF0219DIGILOCK                | 165        |
| 4A1DR10<br>4A1DR12           |            | 4B1DK18<br>4B1DSTM |     | 4B30F12<br>4B30F15     |            | 5PF333                               | 367        | 5RLDF0519                                     | 165        |
| 4A1DR12                      |            | 4B1DSTS            |     | 4B30F18                |            | 5PF334                               | 367        | 5RLDF0519DIGILOCK                             | 165        |
| 4A1DR18                      |            | 4B10F09            |     |                        |            | 5PF335                               | 367        | 5RLDG0519                                     | 165        |
| 4A1DSTM                      |            | 4B1OF10            |     |                        |            | 5PT336                               | 369        | 5RLDG0519DIGILOCK                             | 165        |
| 4A1DSTS                      |            | 4B1OF12            |     | 4B3RS12                |            | 5PT337                               | 369        | 5RSH42WH                                      | 327        |
|                              | 196        | 4B1OF15            | 196 | 4B3RS15                | 197        | 5PT338                               | 369        | 5RSH70WH                                      | 327        |
| 4A10F10                      | 196        | 4B1OF18            | 196 | 4B3RS18                | 197        | 5PT339                               | 369        | 5SM6L   | 362        |
| 4A10F12                      | 196        | 4B1RS09            |     | 4B3VD09/AAA **         | 199        | 5PU151                               | 373        | 5TC101  | 370        |
| 4A10F15                      |            | 4B1RS10            |     |                        |            | 5PU152                               |            | 5TC102  | 370        |
| 4A10F18                      |            | 4B1RS12            |     |                        |            | 5PU153                               |            |   | 370        |
| 4A1OTTM                      |            | 4B1RS15            |     | 4D1DSTCS               |            | 5PU4                                 | 372        | 5TC108  | 371        |
| 4A1OTTS                      |            | 4B1RS18            |     |                        |            | 5PU5                                 | 372        |   | 371        |
|                              | 197        | 4B1TD09/A          |     | 4D1TD09/B              | 201        | 5PU6                                 |            | 5TC110  | 371        |
| 4A1RS10                      | 197        | 4B1TD09/B          | 200 | 4HH-HTM71/M51/400      | 205        | 5RC/DRK                              | 434        | 5TC201  | 3/4        |
|                              |            |                    |     |                        |            |                                      |            |   |            |

0-9

|    | - |   |  |
|----|---|---|--|
|    |   |   |  |
| 1  |   |   |  |
|    |   | - |  |
| r. |   |   |  |
|    |   |   |  |
| c  |   |   |  |
|    |   |   |  |

|   | BC3CPDWL/BE | 29         |
|---|-------------|------------|
|   | BC3CPDWL/GW | 29         |
|   | BC3CPDWR/BE | 29         |
|   | BC3CPDWR/GW | 29         |
|   | BC3CWL/BE   | 29         |
|   | BC3CWL/GW   | 29         |
|   | BC3CWR/BE   | 29         |
|   | BC3CWR/GW   | 29         |
|   |             |            |
|   | BC3D/BE     | 22         |
|   | BC3D/GW     | 22         |
|   | BC3DPD/BE   | 22         |
|   | BC3DPD/GW   | 22         |
|   | BC3DPDWL/BE | 29         |
|   | BC3DPDWL/GW | 29         |
|   | BC3DPDWR/BE | 29         |
|   | BC3DPDWR/GW | 29         |
|   | BC3DWL/BE   | 29         |
|   | BC3DWL/GW   | 29         |
|   | BC3DWR/BE   | 29         |
|   | BC3DWR/GW   | 29         |
|   | BCC1CWL/BE  | 26         |
|   | BCC1CWL/GW  | 26         |
|   | BCC1CWR/BE  | 26         |
|   | BCC1CWR/GW  | 26         |
|   | BCC1DWL/BE  | 26         |
|   | BCC1DWL/GW  | 26         |
|   | BCC1DWR/BE  | 26         |
|   | BCC1DWR/GW  | 26         |
|   | BH015/WHL   | 325        |
|   | BH030/WHL   | 325        |
|   |             |            |
|   | BO/P/VD/SP  | 422        |
|   | BO/P/VS/SP  | 422        |
|   | BR015/WHL   | 320        |
|   | BR020/WHL   | 320        |
|   | BR030/WHL   | 320        |
|   | BR045/WHL   | 321        |
|   | BR050/WHL   | 321        |
|   | BR060/WHL   | 321        |
|   | BR075/WHL   | 322        |
|   | BR080/WHL   | 322        |
|   | BR090/WHL   | 322        |
|   | BR105/WHL   | 323        |
|   | BR110/WHL   | 323        |
|   | BR120/WHL   | 323        |
|   | BR135/WHL   | 324        |
|   | BR140/WHL   | 324        |
|   | BR150/WHL   | 324        |
| С |             |            |
|   | C001        | 121        |
|   | C002        | 121        |
|   | C002        | 121        |
|   | C003        | 76         |
|   | C005        | 76         |
|   | C005        | 76         |
|   | C008        | 309        |
|   |             |            |
|   | C013        | 121<br>121 |
|   | C015        |            |
|   | CB/100/A    | 86         |
|   |             |            |

|   | 8SXRSARTH   | 298  |
|---|-------------|------|
|   | 8SXRSCAST   | 298  |
|   |             |      |
|   | 8SXRSCCU    | 298  |
|   | 8SXRSDRS    | 298  |
|   | 8SXRSENT1   | 440  |
|   | 8SXRSGS2    | 441  |
|   | 8SXRSIV     |      |
|   |             | 298  |
|   | 8SXRSTRAMA  | 298  |
|   | 8SXRT38MW14 | 440  |
|   | 8SXRTGS2    | 441  |
|   | 8SXRTGS3    | 441  |
| В | 00/1110000  |      |
| в | DA /0/D     | 0.40 |
|   | BA/3/B      | 240  |
|   | BA/6        | 240  |
|   | BA/6/A      | 240  |
|   | BA/VH       | 238  |
|   | BA/VH/NC    | 238  |
|   | BA030       | 242  |
|   |             | 20   |
|   | BC1C/BE     |      |
|   | BC1C/GW     | 20   |
|   | BC1CPD/BE   | 20   |
|   | BC1CPD/GW   | 20   |
|   | BC1CPDWL/BE | 27   |
|   | BC1CPDWL/GW | 27   |
|   |             |      |
|   | BC1CPDWR/GW | 27   |
|   | BC1CWL/BE   | 27   |
|   | BC1CWL/GW   | 27   |
|   | BC1CWR/BE   | 27   |
|   | BC1CWR/GW   | 27   |
|   | BC1CWRPD/BE | 27   |
|   |             |      |
|   | BC1D/BE     | 20   |
|   | BC1D/GW     | 20   |
|   | BC1DPD/BE   | 20   |
|   | BC1DPD/GW   | 20   |
|   | BC1DPDWL/BE | 27   |
|   | BC1DPDWL/GW | 27   |
|   |             |      |
|   | BC1DPDWR/BE | 27   |
|   | BC1DPDWR/GW | 27   |
|   | BC1DWL/BE   | 27   |
|   | BC1DWL/GW   | 27   |
|   | BC1DWR/BE   | 27   |
|   | BC1DWR/GW   | 27   |
|   | BC2C/BE     | 21   |
|   |             |      |
|   | BC2C/GW     | 21   |
|   | BC2CWL/BE   | 28   |
|   | BC2CWL/GW   | 28   |
|   | BC2CWR/BE   | 28   |
|   | BC2CWR/GW   | 28   |
|   | BC2D/BE     | 21   |
|   |             |      |
|   | BC2D/GW     | 21   |
|   | BC2DWL/BE   | 28   |
|   | BC2DWL/GW   | 28   |
|   | BC2DWR/BE   | 28   |
|   | BC2DWR/GW   | 28   |
|   |             |      |
|   | BC3C/BE     | 22   |
|   | BC3C/GW     | 22   |
|   | BC3CPD/BE   | 22   |
|   | BC3CPD/GW   | 22   |
|   |             |      |
|   |             |      |

| 8MX546MS          | 392 |
|-------------------|-----|
| 8MX546SG          | 392 |
| 8MX546SS          | 392 |
| 8MX556MG          | 392 |
| 8MX556MS          | 392 |
| 8MX556SG          | 392 |
| 8MX556SS          | 392 |
| 8MX566MG          | 392 |
| 8MX566MS          | 392 |
|                   |     |
| 8MX566SG          | 392 |
| 8MX566SS          | 392 |
| 8MXD24-8          | 393 |
| 8MXL36-2S         | 393 |
| 8MXL42-2S         | 393 |
| 8MXL48-2S         | 393 |
| 8MXL60-2S         | 393 |
| 8MXLS24-2S        | 393 |
| 8MY1627-24BU      | 192 |
| 8MY1627-34BU      | 192 |
| 8MY2030-24BU      | 192 |
|                   |     |
| 8MY2030-34BU      | 192 |
| 8MY2636-25BU      | 192 |
| 8MY2636-35BU      | 192 |
| 80536C            | 396 |
| 80536CM           | 396 |
| 80546C            | 396 |
| 80546CM           | 396 |
| 80556C            | 396 |
| 80556CM           | 396 |
| 80966C            | 396 |
| 80966CM           | 396 |
| 80966CM<br>8Q04LH | 395 |
|                   |     |
| 8Q536G2           | 394 |
| 8Q536G2M          | 394 |
| 8Q546G2           | 394 |
| 8Q546G2M          | 394 |
| 8Q556G2           | 394 |
| 8Q556G2M          | 394 |
| 8Q566G2           | 394 |
| 8Q566G2M          | 394 |
| 8SXRD40MW10       | 440 |
| 8SXRD72SCATH1     | 443 |
| 8SXRD76CM3        | 441 |
| 8SXRD76CMHD1      | 441 |
| 8SXRD76MCATH1     | 441 |
|                   |     |
| 8SXRD80HCATHTD    | 443 |
| 8SXRD80SCATHTD    | 443 |
| 8SXRD83MCATHTD    | 442 |
| 8SXRD83MGSTD      | 441 |
| 8SXRDENT2         | 440 |
| 8SXRS72SCATH1     | 443 |
| 8SXRS76CM5        | 441 |
| 8SXRS76MCATH1     | 442 |
| 8SXRS80HCATHTD    | 443 |
| 8SXRS80SCATHTD    | 443 |
| 8SXRS80SCATHTD    | 443 |
|                   |     |
| 8SXRS83MGSTD      | 441 |
| 8SXRSANES         | 298 |
|                   |     |

|                      | 000 |  |
|----------------------|-----|--|
| 850537CM             | 396 |  |
| 850547C              | 396 |  |
| 850547CM             | 396 |  |
| 850557C              | 396 |  |
| 85O557CM             | 396 |  |
| 85O967C              | 396 |  |
| 85O967CM             | 396 |  |
| 85Q537G2             | 394 |  |
| 85Q537G2M            | 394 |  |
| 85Q547G2             | 394 |  |
| 85Q547G2M            | 394 |  |
| 85Q557G2             | 394 |  |
| 85Q557G2M            | 394 |  |
| 85Q567G2             | 394 |  |
| 85Q567G2M            | 394 |  |
| 89989PX              | 393 |  |
|                      |     |  |
| 89990P               | 397 |  |
| 8BC2030-2DMB         | 192 |  |
| 8BC2030-3DMB         | 192 |  |
| 8BC2636-2DMB         | 192 |  |
| 8BC2636-3DMB         | 192 |  |
| 8DC15EC              | 427 |  |
| 8DC16EC              | 427 |  |
| 8DC35EC              | 427 |  |
| 8DC36EC              | 427 |  |
| 8DC55EC              | 427 |  |
| 8DC56EC              | 427 |  |
| 8DD24C               | 397 |  |
| 8FL21P3369           | 314 |  |
| 8FLCCU1              | 296 |  |
| 8FLCLINE             | 296 |  |
| 8FLCOM-LT            | 296 |  |
| 8FLCRP4              | 296 |  |
| 8FLPED1              | 296 |  |
| 8FLTMENT2            | 296 |  |
|                      | 397 |  |
| 8L24N-4C             |     |  |
| 8L36N-4C             | 397 |  |
| 8L42N-4C             | 397 |  |
| 8L48N-4C             | 397 |  |
| 8L60N-4C             | 397 |  |
| 8LECCRP2             | 312 |  |
| 8LECCRP3             | 312 |  |
| 8LECCRP4             | 312 |  |
| 8MBC1210DB           | 300 |  |
| 8MBX1310ANES4        | 300 |  |
| 8MBX3210TRAUMA       | 300 |  |
| 8MQD24-8             | 395 |  |
| 8MQL24-2S            | 395 |  |
| 8MQL36-2S            | 395 |  |
| 8MQL42-2S            | 395 |  |
| 8MQL48-2S            | 395 |  |
| 8MQL60-2S            | 395 |  |
| 8MQLS24-2S           | 395 |  |
| 8MX536MG             | 392 |  |
| 8MX536MS             | 392 |  |
| 8MX536SG             | 392 |  |
| 8MX536SS             | 392 |  |
| 8MX536SS<br>8MX546MG |     |  |
| 0IVIA040IVIG         | 392 |  |
|                      |     |  |

| 5TC203         374           5TC208         375           5TC209         375           5TC204         376           5TC304         376           5TC305         376           5TC308         377           5TC303         377           5TC304         378           5TC305         377           5TC301         378           5TC311         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0112-071-000         11           5X-015-061-000         11           5X-015-071-000         11  | 5TC202          | 374 |
|---|-----------------|-----|
| 5TC208         375           STC209         375           STC210         376           STC304         376           STC305         376           STC308         377           STC309         377           STC310         377           STC311         378           STC312         378           STC313         378           STC314         378           STC313         378           STC313         378           STC313         378           STC314         378           STC315         377           STC314         378           STC313         378           STC314         378           STC317         378           STC3100-01-000/K         8           SX-0110-071-000/K         9           SX-0112-071-000         11           SX-015-071-000         11 </td <td></td> <td></td>   |                 |     |
| 5TC209         375           5TC210         375           5TC304         376           5TC305         376           5TC306         377           5TC303         377           5TC311         378           5TC312         378           5TC313         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000         11           5X-015-071-000         11           5X-015-071-000         11           5X-015-071-000         11           5Z-TA2/B/C         360           5Z-ETA2/B/C         360           5Z-ETA2/B/C         360           5Z-ETA2/B/C         360           5Z-ETA2/B/C  |                 |     |
| 5TC210         375           5TC304         376           5TC305         376           5TC306         377           5TC303         377           5TC310         377           5TC311         378           5TC313         378           5TC313         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0100-071-000/K         8           5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-01142-061-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-114-071-000         11           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360  |                 |     |
| 5TC304         376           5TC305         376           5TC308         377           5TC303         377           5TC310         377           5TC311         378           5TC312         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0100-071-000/K         8           5X-0110-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         11           5X-0110-071-000/K         11           5X-0110-071-000/K         11           5X-0110-071-000/K         11           5X-0115-071-000         11           5X-0115-071-000         11           5X-015-071-000         11           5X-151-071-000         11           5X-217A/2B         360           5ZETA1/B         360           5ZETA2/B/L         360           5ZETA2/B/L         360           5ZETA2/B/L         360           5ZETA2/B/L         360           5ZETA2/B/L         360           5ZETA2/B/L <t< td=""><td></td><td></td></t<>                      |                 |     |
| 5TC305         376           5TC306         377           5TC308         377           5TC310         377           5TC311         378           5TC312         378           5TC313         378           5TC314         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA4/B         360           5ZETA4/B  |                 |     |
| 5TC306         376           5TC308         377           5TC309         377           5TC310         377           5TC311         378           5TC312         378           5TC313         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0109-071-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-217A/70         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B  |                 |     |
| 5TC308         377           5TC309         377           5TC310         378           5TC311         378           5TC312         378           5TC313         378           5TC311         378           5TC312         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0112-061-000         11           5X-0150-061-000         11           5X-0150-071-000         10           5X-0150-071-000         11   |                 |     |
| 5TC309         377           5TC310         377           5TC311         378           5TC311         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-061-000/K         8           5X-0109-071-000/K         9           5X-0109-071-000/K         9           5X-0109-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0150-061-000         11           5X-0150-061-000         11           5X-0150-071-000         11           5X-0150-071-070         10           5X-150-071-000         11           5X-0150-071-000         11           5X-0150-071-000         10           5ZETA/D         360           5ZETA/B         360           5ZETA/B         360  |                 |     |
| 5TC310         377           5TC311         378           5TC312         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0109-071-000/K         9           5X-0110-071-000/K         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360   |                 |     |
| 5TC311         378           5TC312         378           5TC312         378           5X-0100-061-000/K         8           5X-0100-071-000/K         8           5X-0100-071-000/K         8           5X-0100-071-000/K         9           5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-01142-061-000         11           5X-0150-071-000         11           5X-120-87         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/RC         360           5ZETA4/S         360           5ZETA4/S         360           5ZETA4/S         360           5ZETA4/S         360   |                 |     |
| 5TC312         378           5TC312         378           5TC313         378           5X-0100-061-000/K         8           5X-0109-061-000/K         8           5X-0109-071-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         9           5X-0110-061-000/K         11           5X-0150-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA1/C         360           5ZETA2/B/RC         360           5ZETA4/S   |                 |     |
| 5TC313         378           5X-0100-061-000/K         8           5X-0109-071-000/K         8           5X-0109-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5X-0151-071-000         16           5ZETA1/K         360           5ZETA2/B/C         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B/C         360           5ZETA4/B/C         360           5ZETA4/B/C         360  |                 |     |
| 5X-0100-061-000/K         8           5X-0100-071-000/K         8           5X-0109-071-000/K         8           5X-0109-071-000/K         9           5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-0110-071-0000         12           5X-0142-071-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0152-071-000         11           5X-0152-071-000         11           5X-0152-071-000         11           5X-0152-071-000         11           5ZETA1/B         360           5ZETA2/B         360           5ZETA2/B/C         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360 <tr< td=""><td></td><td></td></tr<> |                 |     |
| 5X-0100-071-000/K         8           5X-0109-061-000/K         8           5X-0109-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0142-071-000         12           5X-0150-061-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/RC         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA4/B         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360                       |                 |     |
| 5X-0109-061-000/K         8           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         12           5X-0142-071-000         11           5X-0150-071-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA2/B/RC         360           5ZETA2/2B         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360  |                 |     |
| 5X-0109-071-000/K         8           5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-0110-071-000/K         9           5X-0142-071-000         12           5X-0150-071-000         11           5Z-141/B         360           5ZETA1/CO         360           5ZETA2/2/B         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC49/B         360           5ZETAC49/B         360           5ZETAC49/B         360           5ZETAC49/B         360           5ZETAC49/B         360           5ZETAC49/B         360           <                                    |                 |     |
| 5X-0110-061-000/K         9           5X-0110-071-000/K         9           5X-0142-061-000         12           5X-0142-071-000         12           5X-0142-071-000         11           5X-0150-061-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA1/CO         360           5ZETA2/B/LC         360           5ZETAC90/B         360           5ZETA4/B         360   |                 |     |
| 5X-0110-071-000/K         9           5X-0142-061-000         12           5X-0142-071-000         11           5X-0150-061-000         11           5X-0150-061-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         10           5ZETA/IC         360           5ZETA/Z2         360           5ZETA/Z2B         360           5ZETA/ZB/RC         360           5ZETA/S36MG         392           85MX536MG         392           85MX546MG         392           85MX546MS         392           85MX546MS         392           85MX546MS<  |                 |     |
| 5X-0142-061-000         12           5X-0142-071-000         12           5X-0150-061-000         11           5X-0150-061-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5Z-0151-071-000         11           5Z-0151-071-000         12           5Z-0151-071-000         11           5Z-0151-071-000         11           5Z-0151-071-000         11           5Z-0151-071-000         360           5Z-0151-071-000         360           5Z-0151-071-08         360           5Z-072         360           5Z-072/2B         360           5  |                 |     |
| 5X-0142-071-000         12           5X-0150-061-000         11           5X-0150-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA1/CO         360           5ZETA2         360           5ZETA2         360           5ZETA2/B/LC         360           5ZETA4/B         360           5ZETA4/B  |                 |     |
| 5X-0150-061-000         11           5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-061-000         11           5Z-0151-071-000         11           5ZETA1/B         360           5ZETA1/C         360           5ZETA2         360           5ZETA2/2B         360           5ZETA2/2B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAH         360           5ZETAH         360           5ZETAK         360           5   |                 |     |
| 5X-0150-071-000         11           5X-0151-061-000         11           5X-0151-071-000         11           5ZETA1/B         360           5ZETA1/B         360           5ZETA1/C         360           5ZETA2/B/C         360           5ZETA2/2B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETA2/B/C         360           5ZETAC45/B         392           85MX566MG         392           85MX556MG         392           85MX556MG         392   |                 |     |
| 5X-0151-061-000         11           5X-0151-071-000         11           SZETA1/B         360           5ZETA1/CO         360           5ZETA1/CO         360           5ZETA1/CO         360           5ZETA2/2B         360           5ZETA2/2B         360           5ZETA2/2B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA/2/B/RC         360           5ZETA/30C         360           5ZETA/4         360           5ZETA/   | 5X-0150-061-000 | 11  |
| 5X-0151-071-000         11           SZETA1         360           SZETA1/B         360           SZETA1/C         360           SZETA1/C         360           SZETA2/B         360           SZETA2/B         360           SZETA2/B/C         360           SZETA2/B/C         360           SZETA2/B/RC         360           SZETA2/B/RC         360           SZETAC45/B         360           SZETAC45/B         360           SZETAC45/B         360           SZETAH         360  | 5X-0150-071-000 | 11  |
| 52ETA1         360           5ZETA1/C         360           5ZETA2/2B         360           5ZETA2/2B         360           5ZETA2/2B         360           5ZETA2/2B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAH/B         360           5ZETAH/B         360           5ZETAK         360           5ZETAK         360           5ZETAH/B         360           5ZETAK         392  | 5X-0151-061-000 | 11  |
| 52ETA1         360           5ZETA1/C         360           5ZETA2/2B         360           5ZETA2/2B         360           5ZETA2/2B         360           5ZETA2/2B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAH/B         360           5ZETAH/B         360           5ZETAK         360           5ZETAK         360           5ZETAH/B         360           5ZETAK         392  | 5X-0151-071-000 | 11  |
| 5ZETA1/CO         360           5ZETA2/2         360           5ZETA2/2B         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         392           85MX546KG         392           85MX556KG         392           85MX556KG         392           85MX556KG         392           85MX556KM         392  |                 | 360 |
| 5ZETA1/CO         360           5ZETA2/2         360           5ZETA2/2B         360           5ZETA2/B/LC         360           5ZETA2/B/LC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC90         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         360           5ZETAC45/B         392           85MX546KG         392           85MX556KG         392           85MX556KG         392           85MX556KG         392           85MX556KM         392  | 5ZETA1/B        | 360 |
| 52ETA2         360           52ETA2/2B         360           52ETA2/B/RC         360           52ETA2/B/RC         360           52ETA2/B/RC         360           52ETA2/B/RC         360           52ETA42/B/RC         360           52ETAC45/5         360           52ETAC45/8         360           52ETAC45/8         360           52ETAC45/8         360           52ETAH         392           85MX546MG         392           85MX556MG         392      8  |                 | 360 |
| 5ZETA2/28 360<br>5ZETA2/B/LC 360<br>5ZETA2/B/RC 360<br>5ZETAC45 360<br>5ZETAC45 360<br>5ZETAC45 360<br>5ZETAC90/B 360<br>5ZETAC90/B 360<br>5ZETAC90/B 360<br>5ZETAH/B 360<br>5ZETAH/B 360<br>5ZETA4/B 360<br>5ZETA4/B 360<br>85MX536MG 392<br>85MX536S 392<br>85MX546MG 392<br>85MX546MG 392<br>85MX546MG 392<br>85MX556MG 392<br>85MX556MS 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX566MG 392<br>85MX566MG 392   |                 | 360 |
| 5ZETA2/B/LC         360           5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETA42/B/RC         360           5ZETAC45/B         360           5ZETAC9/B         360           5ZETAC90/B         360           5ZETAQ         360           5ZETAQ         360           5ZETAG90/B         392           85MX536MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392  |                 |     |
| 5ZETA2/B/RC         360           5ZETA2/B/RC         360           5ZETAC45/5         360           5ZETAC45/5         360           5ZETAC45/5         360           5ZETAC45/5         360           5ZETAC90/B         360           5ZETAM         360           5ZETAH         360           5ZETAH         360           5ZETAK         360           5ZETAH         360           5ZETAH         360           5ZETAH         360           5ZETAK         392           85MX546MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX566MS         392           85MX566MS         392   |                 |     |
| 5ZETA930C         360           5ZETAC45         360           5ZETAC45         360           5ZETAC45         360           5ZETAC45         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC49/B         360           5ZETAH/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA4/B         360           5ZETA/B         392           85MX546MG         392           85MX546MG         392           85MX556MS         392           85MX556MS         392           85MX556MG         392           85MX556MG         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392 </td <td></td> <td></td>   |                 |     |
| 5ZETAC45         360           5ZETAC45/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAH         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAC90/B         360           5ZETAQ         360           5ZETAQ         360           85MX536MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX566MG  |                 |     |
| 5ZETAC45/B         360           5ZETAC90         360           5ZETAC90/B         360           5ZETAH         360           5ZETAH         360           5ZETAM         360           5ZETAH         360           5ZETAM         360           5ZETAM         360           5ZETAM         360           5ZETAM         360           5ZETAM         360           5ZETAM         392           85MX536MS         392           85MX546MG         392           85MX546MS         392           85MX556MS         392           85MX566MS         392 </td <td></td> <td></td>   |                 |     |
| 5ZETAC90 360<br>5ZETAC90/B 360<br>5ZETAH/B 360<br>5ZETAH/B 360<br>5ZETAH/B 360<br>85MX536MG 392<br>85MX536MG 392<br>85MX536MS 392<br>85MX546MG 392<br>85MX546MS 392<br>85MX546MS 392<br>85MX546MS 392<br>85MX556MS 392<br>85MX556MS 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MS 392<br>85MX566MS 392<br>85MX566MS 392<br>85MX566MS 392<br>85MX566MS 392  |                 |     |
| 5ZETAC90/B 360<br>5ZETAH/B 360<br>5ZETAH/B 360<br>5ZETAQ 360<br>85MX536MG 392<br>85MX536S 392<br>85MX536S 392<br>85MX546MG 392<br>85MX546MG 392<br>85MX546MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX556MG 392<br>85MX566MG 392<br>85MX566MG 392<br>85MX566MG 392<br>85MX566MG 392<br>85MX566MS 392<br>85MX566MS 392   |                 |     |
| 5ZETAH         360           5ZETAH/B         360           5ZETAQ         360           5ZETAQ         360           SSMX536MS         392           85MX536MS         392           85MX536MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX556MS         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566SMS         392           85MX566SMS         392           85MX566SMS         392           85MX566SMS         392           85MX566SS   |                 |     |
| 5ZETAH/B 360<br>5ZETAQ 360<br>85MX536MG 392<br>85MX536MS 392<br>85MX536SG 392<br>85MX546SG 392<br>85MX546MG 392<br>85MX546MS 392<br>85MX546SG 392<br>85MX556MS 392<br>85MX556MS 392<br>85MX556SM 392<br>85MX5566MG 392<br>85MX5566MG 392<br>85MX566MG 392<br>85MX566MS 392  |                 |     |
| 5ZETAQ         360           85MX536MG         392           85MX536MS         392           85MX536S         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX556MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566SS         392   |                 |     |
| 85MX536MG         392           85MX536MS         392           85MX536SG         392           85MX546MG         392           85MX546MG         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX556MG         392           85MX556MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566SS         392   |                 |     |
| 85MX536MS         392           85MX536SG         392           85MX536SG         392           85MX546MG         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX546MS         392           85MX556MS         392           85MX556MS         392           85MX556SMS         392           85MX5566MS         392           85MX5566MS         392           85MX5566MS         392           85MX566MS         392           85MX566SS         392   |                 |     |
| 85MX536SG         392           85MX536SS         392           85MX546MG         392           85MX546MG         392           85MX546MG         392           85MX546MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX556MS         392           85MX566MS         392   |                 |     |
| 85MX536SS         392           85MX546MG         392           85MX546MS         392           85MX546SS         392           85MX556MG         392           85MX566MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566SS         392   |                 |     |
| 85MX546MG         392           85MX546MS         392           85MX546SG         392           85MX556MG         392           85MX556MG         392           85MX556SMS         392           85MX556SMS         392           85MX5566MG         392           85MX5566MG         392           85MX5566MG         392           85MX566MS         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566SS         392  |                 |     |
| 85MX546MS         392           85MX546SG         392           85MX556MG         392           85MX556MG         392           85MX556MS         392           85MX556SG         392           85MX556MG         392           85MX556MS         392           85MX556MS         392           85MX566MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566SS         392   |                 |     |
| 85MX546SG         392           85MX546SS         392           85MX556MG         392           85MX556SG         392           85MX556SG         392           85MX556SG         392           85MX556SG         392           85MX556SG         392           85MX566MG         392           85MX566MG         392           85MX566SG         392           85MX566SG         392           85MX566SG         392           85MX566SG         392           85MX566SS         392           85MX566SS         392   |                 |     |
| 85MX546SS         392           85MX556MG         392           85MX556MS         392           85MX556SS         392           85MX556SG         392           85MX566MG         392           85MX566MS         392   |                 |     |
| 85MX556MG         392           85MX556MS         392           85MX556SS         392           85MX566MG         392           85MX566MG         392           85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566MS         392           85MX566SS         392   |                 |     |
| 85MX5566MS 392<br>85MX556SG 392<br>85MX556SS 392<br>85MX566MG 392<br>85MX566MS 392<br>85MX566GMS 392<br>85MX566SS 392   |                 |     |
| 85MX556SG         392           85MX556SS         392           85MX566MG         392           85MX566MG         392           85MX566SS         392           85MX566SS         392           85MX566SS         392   |                 |     |
| 85MX556SS         392           85MX566MG         392           85MX566MS         392           85MX566SG         392           85MX566SS         392           85MX566SS         392   |                 |     |
| 85MX566MG         392           85MX566MS         392           85MX566MS         392           85MX566SG         392           85MX566SS         392   |                 |     |
| 85MX566MS 392<br>85MX566SG 392<br>85MX566SS 392   |                 |     |
| 85MX566SG 392<br>85MX566SS 392  |                 |     |
| 85MX566SS 392   |                 |     |
|   |                 |     |
| 850537C 396   |                 |     |
|   | 850537C         | 396 |

|                    | 00010/00/000           |                         | 077 | 07010014/40/00/500   | 105 |                    | 000 |
|--------------------|------------------------|-------------------------|-----|----------------------|-----|--------------------|-----|
| CB/100/B 86        | CRSM2/18/16D 4         |                         | 277 | CTDISPM1/10/6S/EPB   | 125 |                    | 260 |
| CB005 88           | CRSM2/18/4S12D 4       |                         | 276 | CTDISPM1/10/6S/HSB   | 124 |                    | 260 |
| CB010 88           | CRSM2/18/6S12D 4       |                         | 292 | CTINFEM1/10/3S2D/EPB | 283 |                    | 260 |
| CB015 88           | CRSM2/18/8S10D 4       |                         | 292 | CTINFEM1/10/3S2D/HSB | 282 |                    | 261 |
| CB020 88           | CRSS1/16/2S6D 4        |                         |     | CTINFEM1/10/4S1D/HSB | 282 |                    | 261 |
| CD005 152          | CRSS1/16/3S6D 4        |                         | 279 | CTINFEM1/10/6S/EPB   | 283 |                    | 261 |
| CD010 152          | CRSS1/16/4S5D 4        |                         | 278 | CTINFEM1/10/6S/HSB   | 282 |                    | 288 |
| CD015 152          | CRSS1/16/5S5D 4        |                         | 279 | CTIVPRM1/10/3S2D/EPB | 285 |                    | 288 |
| CD020 152          | CRSS1/16/6S4D 4        |                         |     | CTIVPRM1/10/3S2D/HSB | 284 |                    | 310 |
| CD025 152          | CRSS1/16/8D 4          |                         | 279 | CTIVPRM1/10/4S1D/EPB | 285 |                    | 310 |
| CD030 153          | CRSS2/16/10S10D 4      |                         | 82  | CTIVPRM1/10/4S1D/HSB | 284 |                    | 310 |
| CD035 153          | CRSS2/16/12S8D 4       |                         | 82  | CTIVPRM1/10/6S/EPB   | 285 | CTS1/08/1S2D       | 40  |
| CD045 153          | CRSS2/16/16D 4         |                         | 83  | CTIVPRM1/10/6S/HSB   | 284 | CTS1/08/3D         | 40  |
| CD110 152          | CRSS2/16/4S12D 4       |                         | 83  | CTLAPTM1/10/3S2D/EPB | 437 | CTS1/08/3S1D       | 40  |
| CD115 152          | CRSS2/16/6S12D 4       |                         | 283 | CTLAPTM1/10/3S2D/HSB | 436 | CTS1/08/4S         | 40  |
| CD120 152          | CRSS2/16/8S10D 4       |                         | 280 | CTLAPTM1/10/4S1D/EPB | 437 | CTS1/10/1S3D       | 42  |
| CD125 152          |                        | 74 CTCRITM1/10/3S2D/EPB | 281 | CTLAPTM1/10/4S1D/HSB | 436 | CTS1/10/3S2D       | 42  |
| CD130 153          |                        | 73 CTCRITM1/10/3S2D/EPB | 281 | CTLAPTM1/10/6S/EPB   | 437 | CTS1/10/4D         | 42  |
| CD135 153          | CS010 1                |                         | 280 | CTLAPTM1/10/6S/HSB   | 436 | CTS1/10/4S1D       | 42  |
| CD145 153          | CS015 1                |                         |     | CTM1/08/1S2D         | 50  | CTS1/10/6S         | 42  |
| CD905 159          | CT/AIR 6               |                         | 281 | CTM1/08/3D           | 50  | CTS1/14/2S5D       | 44  |
| CD910 158          | CT/ALC 6               |                         | 61  | CTM1/08/3S1D         | 50  | CTS1/14/3S4D       | 44  |
| CD911 158          |                        | 35 CTD1/08EPB           | 62  | CTM1/08/4S           | 50  | CTS1/14/4S3D       | 44  |
| CD913 158          |                        | 92 CTD1/08HSB           | 63  | CTM1/10/1S3D         | 52  | CTS1/14/5S3D       | 44  |
| CF/F/PK10 252      | CT/BIN 6               | 8 CTD1/08SS             | 64  | CTM1/10/3S2D         | 52  | CTS1/14/6D         | 44  |
| CP 91              | CT/BOW 6               |                         | 61  | CTM1/10/4D           | 52  | CTS1/14/8S1D       | 44  |
| CRMM1/18/2S6D 56   | CT/BUC 6               |                         | 62  | CTM1/10/4S1D         | 52  | CTS2/08/3S4D       | 41  |
| CRMM1/18/3S6D 56   | CT/CR/2 7              |                         | 63  | CTM1/10/6S           | 52  | CTS2/08/4S3D       | 41  |
| CRMM1/18/4S5D 56   | CT/CR/5 7              | 0 CTD1/10SS             | 64  | CTM1/14/2S5D         | 54  | CTS2/08/6D         | 41  |
| CRMM1/18/5S5D 56   | CT/CYD 7               | 0 CTD1/14CA             | 61  | CTM1/14/3S4D         | 54  | CTS2/08/6S2D       | 41  |
| CRMM1/18/6S4D 56   | CT/DEFIB 6             | 7 CTD1/14EPB            | 62  | CTM1/14/4S3D         | 54  | CTS2/08/7S1D       | 41  |
| CRMM1/18/8D 56     | CT/DH/A4 6             |                         | 63  | CTM1/14/5S3D         | 54  | CTS2/08/8S         | 41  |
| CRMM2/18/10S10D 57 | CT/DIV/D 7             |                         | 64  | CTM1/14/6D           | 54  | CTS2/10/3S6D       | 43  |
| CRMM2/18/12S8D 57  | CT/DIV/S 7             |                         | 61  | CTM1/14/8S1D         | 54  | CTS2/10/4S5D       | 43  |
| CRMM2/18/16D 57    | CT/FS 7                |                         | 62  | CTM2/08/3S4D         | 51  | CTS2/10/5S4D       | 43  |
| CRMM2/18/4S12D 57  | CT/FWS 7               |                         | 63  | CTM2/08/4S3D         | 51  | CTS2/10/6S4D       | 43  |
| CRMM2/18/6S12D 57  | CT/IV 7                |                         | 64  | CTM2/08/6D           | 51  | CTS2/10/8D         | 43  |
| CRMM2/18/8S10D 57  | CT/LAB 7               |                         | 61  | CTM2/08/6S2D         | 51  | CTS2/10/8S2D       | 43  |
| CRMS1/16/2S6D 58   | CT/LAB/SEP 7           | 8 CTD2/10EPB            | 62  | CTM2/08/7S1D         | 51  | CTS2/14/10S6D      | 45  |
| CRMS1/16/3S6D 58   | CT/LHNR/S 7            | 5 CTD2/10HSB            | 63  | CTM2/08/8S           | 51  | CTS2/14/12D        | 45  |
| CRMS1/16/4S5D 58   | CT/LHR/S 7             | 4 CTD2/10SS             | 64  | CTM2/10/3S6D         | 53  | CTS2/14/16S2D      | 45  |
| CRMS1/16/5S5D 58   | CT/OB/MS 6             | 6 CTD2/14CA             | 61  | CTM2/10/4S5D         | 53  | CTS2/14/4S10D      | 45  |
| CRMS1/16/6S4D 58   | CT/OB/MS 3             | 11 CTD2/14EPB           | 62  | CTM2/10/5S4D         | 53  | CTS2/14/6S8D       | 45  |
| CRMS1/16/8D 58     | CT/OB/SS 6             | 6 CTD2/14HSB            | 63  | CTM2/10/6S4D         | 53  | CTS2/14/8S6D       | 45  |
| CRMS2/16/10S10D 59 | CT/PARM 7              | 1 CTD2/14SS             | 64  | CTM2/10/8D           | 53  | CTSEPSM1/10/6S/EPB | 387 |
| CRMS2/16/12S8D 59  | CT/SBB/15 7            | 1 CTDC1/08CA            | 61  | CTM2/10/8S2D         | 53  | CTSEPSM1/10/6S/HSB | 386 |
| CRMS2/16/16D 59    | CT/SPS 7               | 1 CTDC1/08EPB           | 62  | CTM2/14/10S6D        | 55  | CTVACIM1/10/3S2D   | 290 |
| CRMS2/16/4S12D 59  | CT/TRAY/D 7            | 7 CTDC1/08SS            | 64  | CTM2/14/12D          | 55  | CTVACIM1/10/4S1D   | 290 |
| CRMS2/16/6S12D 59  | CT/TRAY/S 7            | 7 CTDC1/10CA            | 61  | CTM2/14/16S2D        | 55  | CTVACIM1/10/6S     | 290 |
| CRMS2/16/8S10D 59  | CT/WB 7                | 2 CTDC1/10EPB           | 62  | CTM2/14/4S10D        | 55  | CWF/80             | 84  |
| CRSM1/18/2S6D 46   | CT/WF 7                | 2 CTDC1/10SS            | 64  | CTM2/14/6S8D         | 55  | CWF/80/SP          | 84  |
| CRSM1/18/3S6D 46   | CTAIRWM1/10/3S2D/SS 3  | 16 CTDC1/14CA           | 61  | CTM2/14/8S6D         | 55  | CYB/DE2            | 103 |
| CRSM1/18/4S5D 46   |                        | 16 CTDC1/14EPB          | 62  | CTMEDSM1/10/3S2D/CA  | 286 | CYB/F2             | 103 |
| CRSM1/18/5S5D 46   | CTAIRWM1/10/6S/SS 3    | 16 CTDC1/14SS           | 64  | CTMEDSM1/10/3S2D/EPB | 287 | CYB/G2             | 103 |
| CRSM1/18/6S4D 46   | CTANAEM1/10/3S2D/EPB 2 | 77 CTDISPM1/10/3S2D/EPB | 125 | CTMEDSM1/10/4S1D/CA  | 286 | CYB/J2             | 103 |
| CRSM1/18/8D 46     | CTANAEM1/10/3S2D/HSB 2 | 76 CTDISPM1/10/3S2D/HSB | 124 | CTMEDSM1/10/4S1D/EPB | 287 |                    | 103 |
| CRSM2/18/10S10D 47 |                        | 77 CTDISPM1/10/4S1D/EPB | 125 | CTMEDSM1/10/6S/CA    | 286 |                    | 103 |
| CRSM2/18/12S8D 47  | CTANAEM1/10/4S1D/HSB 2 | 76 CTDISPM1/10/4S1D/HSB | 124 | CTMEDSM1/10/6S/EPB   | 287 | CYHP/D1            | 272 |
|                    |                        |                         |     |                      |     |                    |     |

| F | 2 |  |
|---|---|--|
| L | 5 |  |
|   | D |  |
| s |   |  |

| G/RAIL-C           | 308        | MRW/D4D/CAM                  | 249        |
|--------------------|------------|------------------------------|------------|
| G/SB               | 273        | MRW/D4D/CL                   | 249        |
| GC005              | 413        | MRW/LHNR/S                   | 251        |
| GC010              | 413        | MRW/LHR/S                    | 250        |
| GC015              | 413        | MRW/S1DOC/CAM                | 248        |
| GH005              | 108        | MRW/S1DOC/CL                 | 248        |
| GH015              | 108        | MRW/S1DXP/CAM                | 248<br>248 |
| GMBH/15<br>GMBH/15 | 73<br>107  | MRW/S1DXP/CL<br>MRW/S2D/CAM  | 248        |
| GMBH/15<br>GMBH/15 | 435        | MRW/S2D/CAM<br>MRW/S2D/CL    | 248        |
| GMBH/WM            | 435        | MS/C4                        | 248<br>92  |
| GP005              | 193        | MS/C4<br>MS/S3               | 92         |
| GP003<br>GP015     | 193        | MS/S3<br>MS/S4               | 93         |
| GP013<br>GP025     | 193        | M3/34<br>MT/1/6F             | 65         |
| GP035              | 193        | MT/1/6F                      | 142        |
| 1                  | 155        | MT/1/6F                      | 384        |
| IS/P/2H/H/SP       | 216        | MT/LHNR/S                    | 130        |
| IS/P/2H/SP         | 216        | MT/LHR/S                     | 129        |
| IS/P/4H/H/SP       | 216        | MT100                        | 126        |
| IS/P/4H/SP         | 216        | MT100CL                      | 126        |
| IV001              | 221        | MT10002                      | 126        |
| IV002              | 221        | MT110CL                      | 126        |
| М                  |            | MT120                        | 126        |
| MC/DIV             | 159        | MT120CL                      | 126        |
| MC/JK/1870         | 159        | MT130                        | 128        |
| MC/JK/870          | 159        | MT140                        | 128        |
| MC/L/1870/7S5D     | 157        | MT150                        | 128        |
| MC/L/1870/FS       | 157        | MT220                        | 127        |
| MC/L/1870/FSD      | 157        | MT220CL                      | 127        |
| MC/L/870/3S2D      | 156        | MT230                        | 127        |
| MC/L/870/FS        | 156        | MT230CL                      | 127        |
| MC/R/1870/7S5D     | 157        | MW005                        | 294        |
| MC/R/1870/FS       | 157        | MW010                        | 315        |
| MC/R/1870/FSD      | 157        | MW015                        | 315        |
| MC/R/870/3S2D      | 156        | 0                            |            |
| MC/R/870/FS        | 156        | OBL/BE0/GW/S                 | 32         |
| MC/ST              | 158        | OBL/BE1/GW/S                 | 32         |
| MR110<br>MR120     | 257        | OBL/BE3/GW/S<br>OBL/GW0/SP/S | 32<br>32   |
|                    | 257        | OBL/GW0/SP/S<br>OBL/GW1/SP/S | 32         |
| MR130<br>MR140     | 256<br>256 | OBL/GW1/SP/S<br>OBL/GW3/SP/S | 32<br>32   |
| MR140<br>MR140CL   | 256        | OBL/GW3/3F/3<br>OBP/SP/S     | 30         |
| MR1400L<br>MR210   | 246        | OCL/BE0/GW/S                 | 33         |
| MR210CL            | 246        | OCL/BE1/GW/S                 | 33         |
| MR310              | 246        | OCL/BE3/GW/S                 | 33         |
| MR310CL            | 246        | OCL/GW0/SP/S                 | 33         |
| MR410              | 246        | OCL/GW1/SP/S                 | 33         |
| MR410CL            | 246        | OCL/GW3/SP/S                 | 33         |
| MRF/L/LD/CAM       | 254        | OCP/SP/S                     | 30         |
| MRF/L/LD/CL        | 254        | OP/431                       | 426        |
| MRF/L/OC           | 254        | OP/432                       | 426        |
| MRF/S/LD/CAM       | 254        | OP/631                       | 426        |
| MRF/S/LD/CL        | 254        | OP/632                       | 426        |
| MRF/S/OC           | 254        | P                            |            |
| MRW/D2DOC/CAM      | 249        | PC/1000/ST                   | 158        |
| MRW/D2DOC/CL       | 249        | PC/300/ST                    | 158        |
| MRW/D2DXP/CAM      | 249        | PC/500/ST                    | 158        |
| MRW/D2DXP/CL       | 249        | PC/600/ST                    | 158        |
|                    |            |                              |            |

|   | EC1350/VBB     | 171 |
|---|----------------|-----|
|   | EC2205/VBB     |     |
|   |                | 176 |
|   | EC2210/VBB     | 176 |
|   | EC2215/VBB     | 176 |
|   | EC2220/VBB     | 176 |
|   | EC2225/VBB     | 177 |
|   |                |     |
|   | EC2230/VBB     | 177 |
|   | EC2235/VBB     | 177 |
|   | EC2240/VBB     | 177 |
|   | EC2305/VBB     | 178 |
|   |                |     |
|   | EC2310/VBB     | 178 |
|   | EC2315/VBB     | 178 |
|   | EC2320/VBB     | 178 |
|   | EC2325/VBB     | 179 |
|   | EC2330/VBB     | 179 |
|   | EC2335/VBB     | 179 |
|   |                |     |
|   | EC2340/VBB     | 179 |
|   | EC4305/VBB     | 180 |
|   | EC4310/VBB     | 180 |
|   | EC4315/VBB     | 180 |
|   | EC4320/VBB     | 180 |
|   | EC4325/VBB     | 180 |
|   |                |     |
|   | EC4330/VBB     | 181 |
|   | EC4335/VBB     | 181 |
|   | EC4340/VBB     | 181 |
|   | EC4345/VBB     | 181 |
|   | EC5000/VBB     | 184 |
|   |                |     |
|   | EC5000/VBB     | 264 |
|   | EC5305/VBB     | 182 |
|   | EC5305/VBB     | 262 |
|   | EC5310/VBB     | 182 |
|   | EC5310/VBB     | 262 |
|   | EC5315/VBB     | 182 |
|   |                |     |
|   | EC5315/VBB     | 262 |
|   | EC5320/VBB     | 182 |
|   | EC5320/VBB     | 262 |
|   | EC5325/VBB     | 182 |
|   | EC5325/VBB     | 262 |
|   |                |     |
|   | EC5330/VBB     | 183 |
|   | EC5330/VBB     | 263 |
|   | EC5335/VBB     | 183 |
|   | EC5335/VBB     | 263 |
|   | EC5340/VBB     | 183 |
|   | EC5340/VBB     | 263 |
|   |                |     |
|   | EC5345/VBB     | 183 |
|   | EC5345/VBB     | 263 |
| G |                |     |
|   | G/200/RS/VBB   | 268 |
|   | G/205/FS/VBB   | 270 |
|   | G/305/FS/VBB   | 269 |
|   | G/500/RS/CYL   | 209 |
|   |                |     |
|   | G/500/RS/VBB   | 10  |
|   | G/500/RS/VBB   | 271 |
|   | G/CARD/LOCK    | 273 |
|   | G/CARD/STATION | 273 |
|   | G/CLS          | 272 |
|   | G/CLS/500      |     |
|   |                | 272 |
|   | G/DP           | 272 |
|   |                |     |

| 2         DTSF/900/2/SD         424           2         DTSF/96/2/SD         118           2         DTSR/450/SR         114           0         DTSR/60/SR         114           0         DTSR/60/SR         114           0         DTSR/60/SR         114           1         DTSR/50/SR         114           3         DTSR/50/SR         114           0         DUS/F         119           0         DUJ/R         119           0         DUS/F         119           0         DUS/F         119           0         DUS/F         119           0         DUS/F         119           0         E04         E04/U           185         E06         185           136         E07         185           14         E11         186           15         187         187           16         E11         186  |                                 |   |   |
|---|---------------------------------|---|---|
| E         E/SOCKET/4         308           7         E02         185           7         E04 - E04/U         185           9         E05         185           60         185         E06           60         185         E07           8         E07         185           9         E06         185           60         185         E07           8         E09         186           0         E11         186           0         E12         186           0         E13         186           0         E14         187           9         E21         187           8         EBC020         18           8         EBC020         18           8         EC010/VBB         172           5         EC120/VBB         168           5         EC1225/VBB         168  | 2<br>2<br>0<br>0<br>1<br>1<br>3 | DTSF/96/2/SD<br>DTSF/96/2/SD<br>DTSR/450<br>DTSR/450/SR<br>DTSR/600/SR<br>DTSR/600/SR<br>DTSR/750<br>DTSR/900/SR<br>DTSR/900<br>DTSR/900/SR<br>DUL/F<br>DUL/R | 118<br>425<br>114<br>114<br>114<br>114<br>114<br>114<br>114<br>114<br>114<br>11 |
| E/SOCKET/4         308           F         E02         185           E03         185         E03         185           F         E04 - E04/U         185         E06         185           E07         185         E06         185           E07         185         E06         185           E07         185         E08         186           E09         186         E09         186           0         E11         186         13         186           0         E14         187         9         E13         187           9         E19         187         187         9         E21         187           9         E19         187         188           | 6                               |   | 119   |
| 9         E05         185           E06         185           E07         185           E08         186           E09         186           0         E10           0         E11           1         186           0         E12           1         E13           1         E13           1         E14           1         E13           1         E14           1         E15           1         E14           1         E15           1         E14           1         E14           1         E15           187         9           E14         187           9         E21           188         EBC015           18         EBC025           182         188           EC010/VBB         172           17         EC1205/VBB           168         EC1225/VBB           168         EC1225/VBB           168         EC1220/VBB           169         EC1240/VBB           160         EC1240/VBB   |                                 | E/SOCKET/4<br>E02   | 185   |
| 0         E10         186           0         E11         186           0         E12         186           1         E13         186           0         E14         187           0         E14         187           0         E13         187           9         E13         187           9         E21         187           9         EBC015         18           8         EBC020         18           8         EBC030         18           8         EC010/VBB         188           8         EC010/VBB         188           5         EC120/VBB         168           5         EC120/VBB         168           5         EC120/VBB         168           5         EC1220/VBB         168           5         EC1220/VBB         168           5         EC1230/VBB         169           6         EC1235/VBB         169           7         EC1305/VBB         170           8         EC1315/VBB         170           8         EC1320/VBB         170           6 <td< td=""><td></td><td>E05<br/>E06<br/>E07<br/>E08</td><td>185<br/>185<br/>185<br/>185<br/>185<br/>186</td></td<> |                                 | E05<br>E06<br>E07<br>E08  | 185<br>185<br>185<br>185<br>185<br>186  |
| 0         E15         187           9         E18         187           9         E19         187           9         E21         187           9         E21         187           9         E8015         18           8         EBC020         18           8         EBC020         18           8         EBC020         18           8         EBC020         18           8         EBC030         18           8         EC010/VBB         172           5         EC120/VBB         168           5         EC120/VBB         168           5         EC125/VBB         168           5         EC125/VBB         168           5         EC1240/VBB         169           6         EC1240/VBB         169           6         EC130/VBB         170           6         EC130/VBB         170           8         EC130/VBB         170           8         EC130/VBB         170           8         EC1320/VBB         170           5         EC1320/VBB         170           6  | 0<br>0<br>1                     | E10<br>E11<br>E12<br>E13  | 186<br>186<br>186<br>186  |
| 9         EBC015         18           8         EBC020         18           8         EBC020         18           8         EBC030         18           8         EBC030         18           8         EBC030         18           8         EC010/VBB         172           5         EC1210/VBB         168           5         EC1210/VBB         168           5         EC1225/VBB         168           5         EC1225/VBB         168           5         EC1225/VBB         169           6         EC125/VBB         169           6         EC125/VBB         169           6         EC1250/VBB         169           6         EC130/VBB         170           8         EC130/VBB         170           8         EC1310/VBB         170           8         EC1320/VBB         170           8         EC1320/VBB         170           8         EC1320/VBB         170           8         EC1320/VBB         170           8         EC1330/VBB         171           9         EC1330/VBB         171  | 0<br>9                          | E15<br>E18  | 187<br>187  |
| 8         EBC030         18           8         EC010/VBB         172           5         EC120/VBB         168           5         EC1210/VBB         168           5         EC125/VBB         168           5         EC1220/VBB         168           5         EC1225/VBB         168           5         EC1225/VBB         168           5         EC1225/VBB         169           6         EC1230/VBB         169           6         EC1240/VBB         169           6         EC1240/VBB         169           6         EC1240/VBB         169           6         EC1240/VBB         169           6         EC1305/VBB         170           8         EC1310/VBB         170           8         EC1310/VBB         170           5         EC1325/VBB         170           6         EC1330/VBB         171           7         EC1330/VBB         171           6         EC1340/VBB         171  | 9<br>8                          | EBC015<br>EBC020  | 18<br>18  |
| 5         EC1215/VBB         168           5         EC1220/VBB         168           6         EC1220/VBB         168           5         EC1220/VBB         168           5         EC1230/VBB         169           6         EC1240/VBB         169           7         EC1240/VBB         169           6         EC1240/VBB         169           6         EC1250/VBB         169           7         EC1305/VBB         170           8         EC1316/VBB         170           8         EC1320/VBB         170           6         EC1326/VBB         170           6         EC1330/VBB         171           6         EC1330/VBB         171           6         EC1330/VBB         171   | 8<br>8<br>5                     | EBC030<br>EC010/VBB   | 18<br>172   |
| 6         EC1235/VBB         169           7         EC1240/VBB         169           4         EC1240/VBB         169           5         EC1250/VBB         169           7         EC1305/VBB         170           4         EC1310/VBB         170           5         EC1320/VBB         170           6         EC1320/VBB         170           6         EC1335/VBB         170           7         EC1330/VBB         170           6         EC1335/VBB         170           6         EC1330/VBB         171           6         EC1340/VBB         171  | 5<br>5                          | EC1215/VBB<br>EC1220/VBB  | 168<br>168  |
| 6         EC1250/VBB         169           7         EC1305/VBB         170           4         EC1310/VBB         170           5         EC1315/VBB         170           6         EC1315/VBB         170           7         EC1320/VBB         170           7         EC1320/VBB         170           7         EC1330/VBB         171           6         EC1336/VBB         171  | 6<br>7                          | EC1235/VBB<br>EC1240/VBB  | 169<br>169  |
| 5         EC1320/VBB         170           6         EC1325/VBB         170           7         EC1330/VBB         171           4         EC1336/VBB         171           6         EC1340/VBB         171  | 6<br>7<br>4                     | EC1250/VBB<br>EC1305/VBB  | 169<br>170  |
| 4 EC1335/VBB 171<br>6 EC1340/VBB 171  | 5<br>6                          | EC1320/VBB<br>EC1325/VBB  | 170<br>170  |
|   | 4<br>6                          | EC1335/VBB<br>EC1340/VBB  | 171<br>171  |

| CYT/J1/ST 10<br>CYWP/1 10   |   |
|---|---|
| D3         91           D5         91           DAF/1         10           DAR/1         10           DAR/115         72           DAR/3         10           DG/3         10           DH/A3/50         85           DH/A4/50         82           DP/2H         22           DP/2H         22           DP/2H         22           DPKS/14H         22           DPWK/14H         21           DPWM/14H         21           DPWS/1/4H         21           DTMR/450/SR         11 <td< td=""><td>0606       06       073       0</td></td<> | 0606       06       073       0 |

D

www.bristolmaid.com

| SAC/WP                           | 162        |
|----------------------------------|------------|
|                                  |            |
| SB/221                           | 164        |
| SB/321                           | 164        |
| SB/421                           | 164        |
|                                  |            |
| SBH/15                           | 73         |
| SBH/15                           | 435        |
| SBH/WM                           | 109        |
|                                  |            |
| SBKSL/1                          | 423        |
| SBKSLH/1                         | 423        |
| SBS/2                            | 420        |
|                                  |            |
| SBSA/1                           | 421        |
| SBSA/2                           | 421        |
| SBSL/1                           | 423        |
|                                  |            |
| SBST/2                           | 420        |
| SC/DIV                           | 411        |
| SC/JK/1870                       | 411        |
|                                  |            |
| SC/L/1870/7S5D                   | 410        |
| SC/L/1870/FS                     | 410        |
| SC/L/1870/FSD                    | 410        |
|                                  |            |
| SC/R/1870/7S5D                   | 410        |
| SC/R/1870/FS                     | 410        |
| SC/R/1870/FSD                    | 410        |
|                                  |            |
| SC020                            | 412        |
| SC025                            | 412        |
|                                  |            |
| SC040                            | 412        |
| SC045                            | 412        |
| SEP/DID                          | 384        |
|                                  |            |
| SEP/DIS                          | 384        |
| SEP/L                            | 385        |
| SEP/LH                           | 385        |
|                                  |            |
| SEP/SBH                          | 383        |
| SEP/TIG01                        | 385        |
| SEP/USBBMP                       | 383        |
|                                  |            |
| SEPD/6/EPB                       | 382        |
| SEPD/6/HSB                       | 382        |
| SEPS/6/EPB                       | 382        |
|                                  |            |
| SEPS/6/HSB                       | 382        |
| SKC100                           | 415        |
| SKC100/CL                        | 415        |
|                                  |            |
| SKC20                            | 415        |
| SKC20/CL                         | 415        |
| SKC21/TRAKA                      | 414        |
| SKC30                            |            |
|                                  | 415        |
| SKC30/CL                         | 415        |
| SKC48                            | 415        |
|                                  |            |
| SKC48/CL                         | 415        |
| SKC64                            | 415        |
| SKC64/CL                         | 415        |
| SR                               |            |
|                                  | 120        |
| SS/22/1                          | 426        |
| SS/22/2                          | 426        |
|                                  |            |
|                                  | 108        |
| SS/36                            |            |
| SS/40/MS                         | 110        |
| SS/40/MS                         |            |
| SS/40/MS<br>SS/48/3              | 107        |
| SS/40/MS<br>SS/48/3<br>ST/340/UN | 107<br>404 |
| SS/40/MS<br>SS/48/3              | 107        |
| SS/40/MS<br>SS/48/3<br>ST/340/UN | 107<br>404 |

| PL125/BD       | 417 | R/C/T                     | 252 |
|----------------|-----|---------------------------|-----|
| PL130/BD       | 417 | R/C/TH                    | 252 |
| PL135/BD       | 417 | RC/2/15                   | 110 |
| PL140/BD       | 417 | RC/2/15                   | 110 |
| PL145/BD       | 417 | RC/5/15                   | 110 |
|                |     |                           |     |
| PL150/BD       | 417 | RC/5/15                   | 110 |
| PL900          | 417 | RD/B5/SS                  | 111 |
| PL905          | 417 | RD/U5/SS                  | 111 |
| PL910          | 417 | RD1                       | 91  |
| PL915          | 417 | RD2L                      | 91  |
| PL920          | 417 | RD2R                      | 91  |
| PL925          | 417 | RD3                       | 91  |
| PMH005         | 106 | RFID125/MC                | 25  |
|                |     |                           |     |
| PS/M/HC/4H     | 214 | RFID125/PC                | 25  |
| PS/M/SC/4H     | 214 | RFID125/SC                | 25  |
| PS/P/HC/4H     | 212 | RFID125/UC                | 25  |
| PS/P/SC/4H     | 212 | RFID1356/PC               | 25  |
| PT/LHNR/S      | 140 | RFID1356/SC               | 25  |
| PT/LHR/S       | 140 | RFID1356/SF               | 25  |
| PT050          | 402 | RFID1356/UC               | 25  |
| PT030          | 402 | RFID1356/UF               |     |
|                |     |                           | 25  |
| PTD/EPB/A18C12 | 135 | RPT                       | 23  |
| PTD/EPB/BD69   | 138 | S                         |     |
| PTD/EPB/BP9    | 136 | SA005                     | 143 |
| PTD/EPB/LP24   | 135 | SA025                     | 143 |
| PTD/EPB/MX60   | 139 | SA040                     | 163 |
| PTD/EPB/NC56   | 137 | SA040/CL                  | 163 |
| PTD/EPB/PA32   | 133 | SA045                     | 163 |
| PTD/HSB/A18C12 | 135 | SA045/CL                  | 163 |
|                |     |                           |     |
| PTD/HSB/BD69   | 138 | SAC/213/B                 | 161 |
| PTD/HSB/BP9    | 136 | SAC/213/B/CL              | 161 |
| PTD/HSB/LP24   | 135 | SAC/213/B/CL/VP           | 161 |
| PTD/HSB/MX60   | 139 | SAC/213/B/VP              | 161 |
| PTD/HSB/NC56   | 137 | SAC/213/S                 | 160 |
| PTD/HSB/PA32   | 133 | SAC/213/S/CL              | 160 |
| PTS/115        | 403 | SAC/213/S/CL/VP           | 160 |
| PTS/115/M      | 403 | SAC/213/S/VP              | 160 |
| PTS/EPB/BD46   | 138 | SAC/213/3/VF<br>SAC/313/B | 161 |
|                |     |                           |     |
| PTS/EPB/BP6    | 136 | SAC/313/B/CL              | 161 |
| PTS/EPB/DR6    | 134 | SAC/313/B/CL/VP           | 161 |
| PTS/EPB/LA16   | 135 | SAC/313/B/VP              | 161 |
| PTS/EPB/MX40   | 139 | SAC/313/S                 | 160 |
| PTS/EPB/PA24   | 132 | SAC/313/S/CL              | 160 |
| PTS/HSB/BD46   | 138 | SAC/313/S/CL/VP           | 160 |
| PTS/HSB/BP6    | 136 | SAC/313/S/VP              | 160 |
| PTS/HSB/DR6    | 134 | SAC/413/B                 | 161 |
|                |     |                           |     |
| PTS/HSB/LA16   | 135 | SAC/413/B/CL              | 161 |
| PTS/HSB/MX40   | 139 | SAC/413/B/CL/VP           | 161 |
| PTS/HSB/PA24   | 132 | SAC/413/B/VP              | 161 |
|                |     | SAC/413/S                 | 160 |
| R/BAG/BLUE     | 229 | SAC/413/S/CL              | 160 |
| R/BAG/GREEN    | 229 | SAC/413/S/CL/VP           | 160 |
| R/BAG/GREY     | 229 | SAC/413/S/VP              | 160 |
| R/BAG/ORANGE   | 229 | SAC/FS/2                  | 162 |
| R/BAG/ORANGE   | 229 | SAC/FS/2<br>SAC/FS/3      | 162 |
|                |     |                           |     |
| R/BAG/WHITE    | 229 | SAC/FS/4                  | 162 |
| R/BAG/YELLOW   | 229 | SAC/LT                    | 162 |
|                |     |                           |     |

+44 (0)1258 484455

R

| PC636/LH/LX          | 146 |  |
|----------------------|-----|--|
| PC636/LX             | 146 |  |
| PC639                | 146 |  |
| PC639/FS             | 148 |  |
| PC639/FS/LH          | 148 |  |
| PC639/FS/LH/LX       | 148 |  |
| PC639/FS/LX          | 148 |  |
| PC639/LH             | 146 |  |
| PC639/LH/LX          |     |  |
|                      | 146 |  |
| PC639/LH/MDS         | 150 |  |
| PC639/LH/MDS/LX      | 150 |  |
| PC639/LX             | 146 |  |
| PC639/MDS            | 150 |  |
| PC639/MDS/LX         | 150 |  |
| PC836                | 147 |  |
| PC836/FS             | 149 |  |
| PC836/FS/IND         | 149 |  |
| PC836/FS/IND/LX      | 149 |  |
| PC836/FS/LX          | 149 |  |
| PC836/IND            | 147 |  |
| PC836/IND/LX         | 147 |  |
| PC836/LX             | 147 |  |
| PC836/LX<br>PC839    | 147 |  |
|                      |     |  |
| PC839/FS             | 149 |  |
| PC839/FS/IND         | 149 |  |
| PC839/FS/IND/LX      | 149 |  |
| PC839/FS/LX          | 149 |  |
| PC839/IND            | 147 |  |
| PC839/IND/LX         | 147 |  |
| PC839/LX             | 147 |  |
| PC905                | 159 |  |
| PCB                  | 159 |  |
| PE/30                | 164 |  |
| PL005/BD             | 417 |  |
| PL010/BD             | 417 |  |
| PL015/BD             | 417 |  |
| PL015/BD<br>PL020/BD | 417 |  |
| PL020/BD<br>PL025/BD |     |  |
|                      | 417 |  |
| PL030/BD             | 417 |  |
| PL035/BD             | 417 |  |
| PL040/BD             | 417 |  |
| PL045/BD             | 417 |  |
| PL050/BD             | 417 |  |
| PL055/BD             | 417 |  |
| PL060/BD             | 417 |  |
| PL065/BD             | 417 |  |
| PL070/BD             | 417 |  |
| PL075/BD             | 417 |  |
| PL080/BD             | 417 |  |
| PL085/BD             | 417 |  |
| PL000/BD             | 417 |  |
| PL090/BD<br>PL095/BD |     |  |
|                      | 417 |  |
| PL100/BD             | 417 |  |
| PL105/BD             | 417 |  |
| PL110/BD             | 417 |  |
| PL115/BD             | 417 |  |
| PL120/BD             | 417 |  |
|                      |     |  |

| PC/800/ST             | 158 |
|-----------------------|-----|
| PC/RLBF               | 159 |
| PC/RLBS               | 159 |
| PC1036                | 147 |
| PC1036/FS             | 149 |
| PC1036/FS/IND         | 149 |
| PC1036/FS/IND/LX      | 149 |
| PC1036/FS/LX          | 149 |
| PC1030/1 3/LX         | 145 |
|                       |     |
| PC1036/IND/LX         | 147 |
| PC1036/LX             | 147 |
| PC1039                | 147 |
| PC1039/FS             | 149 |
| PC1039/FS/IND         | 149 |
| PC1039/FS/IND/LX      | 149 |
| PC1039/FS/LX          | 149 |
| PC1039/IND            | 147 |
| PC1039/IND/LX         | 147 |
| PC1039/IND/MDS        | 151 |
| PC1039/IND/MDS/LX     | 151 |
| PC1039/LX             | 147 |
| PC1039/MDS            | 151 |
| PC1039/MDS/LX         | 151 |
| PC336/FS              | 148 |
|                       |     |
| PC336/FS/LH           | 148 |
| PC336/FS/LH/LX        | 148 |
| PC336/FS/LX           | 148 |
| PC536                 | 146 |
| PC536/FS              | 148 |
| PC536/FS/LH           | 148 |
| PC536/FS/LH/LX        | 148 |
| PC536/FS/LX           | 148 |
| PC536/LH              | 146 |
| PC536/LH/LX           | 146 |
| PC536/LH/MDS          | 150 |
| PC536/LH/MDS/LX       | 150 |
| PC536/LX              | 146 |
| PC536/MDS             | 150 |
| PC536/MDS/LX          | 150 |
| PC530/MD3/LX<br>PC539 | 146 |
|                       |     |
| PC539/FS              | 148 |
| PC539/FS/LH           | 148 |
| PC539/FS/LH/LX        | 148 |
| PC539/FS/LX           | 148 |
| PC539/LH              | 146 |
| PC539/LH/LX           | 146 |
| PC539/LH/MDS          | 150 |
| PC539/LH/MDS/LX       | 150 |
| PC539/LX              | 146 |
| PC539/MDS             | 150 |
| PC539/MDS/LX          | 150 |
| PC636                 | 146 |
| PC636/FS              | 140 |
|                       |     |
| PC636/FS/LH           | 148 |
| PC636/FS/LH/LX        | 148 |
| PC636/FS/LX           | 148 |
| PC636/LH              | 146 |
|                       |     |

|   | ŝ  |  |
|---|----|--|
|   |    |  |
| 6 | s' |  |
| ÷ |    |  |
| 4 | Þ  |  |
|   | ł  |  |
|   |    |  |
|   |    |  |

| ST/690/UN/T<br>STSS/690/UN<br>STSS/690/UN/T  | 404<br>406<br>406   |
|--|---|
| TBIN4/MS<br>TBIN4/MS<br>TBIN4/SS<br>TBIN5/MS<br>TBIN5/MS<br>TBIN5/SS<br>TBIN5/MS<br>TBIN6/MS<br>TBIN6/MS<br>TBIN6/MS<br>TBIN6/SS   | 66<br>311<br>66<br>66<br>311<br>66<br>66<br>311<br>66<br>23   |
| U UB/9<br>UB005<br>UB010   | 109<br>109<br>109   |
| V<br>VH/950<br>VH/950/LBT  | 96<br>96  |
| WT/11/MS<br>WT/13/MS<br>WT/33/SS/2<br>WT/33/SS/2<br>WT/33/SS/2<br>WT/33/SS/2<br>WT/33/SS/3/ZP005<br>WT/37/2<br>WT/37/2<br>WT/37/2<br>WT/37/2<br>WT/37/3/ZP010<br>WT/37/3/ZP010<br>WT/42/WT/COVER<br>WT/42/WT/COVER<br>WT/48/COVER/TBAR<br>WT/48/COVER/TBAR<br>WT/48/FBD/TBAR<br>WT/48/FBD/TBAR<br>WT/48/FBD/TBAR<br>WT/50/COVER/TBAR<br>WT/50/COVER<br>WT/50/COVER<br>WT/50/COVER<br>WT/50/CD/TBAR<br>WT/50/CD/TBAR<br>WT/50/CD/TBAR<br>WT/50/CD/TBAR<br>WT/50/CD/TBAR<br>WT/6/2/MS<br>WT/6/2/PL/SS<br>WT/6/2/SS<br>WT/6/PL/SS<br>WT/6/PL/SS<br>WT/6/PL/SS | 228<br>227<br>191<br>191<br>191<br>190<br>190<br>230<br>230<br>230<br>230<br>234<br>234<br>234<br>234<br>234<br>234<br>234<br>232<br>232<br>232 |

 WT/9/2/PL/MS
 226

 WT/9/PL/MS
 226

 WT/9/PL/MS
 226

 WT/SDPL/MS
 226

 WT/BOARD
 307

 WT/CBOARD
 307

 WT/CYL/DLF
 307

 WT/CPL/DLF
 307

 WT/DEFIB
 306

 WT/SPS
 306

 WT200
 304

 WT250
 304

ZP005 121 ZP010 76

Ζ

# **PAINT & LAMINATE**

# Paint & Laminate

Painted surfaces should be cleaned using a soft damp cloth & polished dry. Do not use abrasive cloth. Persistent exposure to wet or humid conditions could result in rust marks, or cause the paint to lift. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.



# Laminate & Melamine

Laminates & melamine should be cleaned using a soft damp cloth & polished dry. Do not use abrasive cloth. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.




# VINYL

Vinyl - Ambla New Richmond in Colour



# Vinyl - Spradling Valencia<sup>™</sup> C5



# Faux Leather - Memphis



# **Care & Cleaning**

Vinyl's are resistant to most mild acids, alkalis, drinks & household stains. Clean with a damp, soapy cloth & rinse well with clean water. Dry the surface thoroughly. A soft brush with soapy water may be used to remove stubborn stains. Ensure surfaces are then rinsed & dried. Do not use solvents, bleaches\*, abrasives, synthetic detergents, wax polishes or aerosol sprays. Disinfectants may be used in dilution as specified by their manufacturers.


\* Chlorine Releasing Agents: Levels of 1-10,000 parts per million should pose no problems as long as the material does not endure long soakings at the higher concentrations. Rinse well with clean water. These include Haz-Tabs & other disinfectants commonly used in hospital infection control programmes. Surfaces are NaDCC compatible.

Swatches

# **FABRIC & PLASTIC**

# Fabric - Camira Advantage

Vacuum regularly. Wipe clean with a damp cloth or shampoo using proprietary upholstery shampoo.



# Fabric - VP Multi-Way Stretch

Wipe with mild detergent. Wash or dry clean using detergent up to 75°C. For disinfection use sodium hypo-chlorite at no more than 1,000 parts per million of active chlorine. Do not use biological or phenolic-based cleaners. Wash & dry completely. Some surface wrinkling may occur.



# FABRIC - Camira 24/7+

Vacuum regularly. Professionally dry clean or wipe clean with a damp cloth.



FYR03 - Equinox

FYR21 - While

Swatches

# **Plastic Seat Shells**

Plastic should be cleaned using a soft damp cloth & polished dry. Do not use abrasive cloth. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.



# Panaz Patternistas for Childern

Vacuum regularly. Wipe clean with a damp cloth. Can be cleaned with water based proprietary cleaners & 5% sodium chlorite bleach solution, Haz-Tabs or Chlor-Clean.



SAVL - Savanna Lavender

# **TECHNICAL INFORMATION**

## **Dimensions & Weight**

- Dimensions are stated in millimetres (width x depth x height)
- · Weights are stated in kilograms
- · Please ensure you are aware of the item size & weight prior to ordering

## Castors

- All castors are fitted with non-marking tyres (with exception of some anti-static versions)
- Anti-static castors can be identified by either a yellow mark on the edge of the tyre, or where twin wheel castors are used, a different colouration to the tyre (the result of the addition of a conductive material)
- Castors & wheel bearings do not require specific maintenance other than regular cleaning
- The use of oil will breakdown the grease applied to the bearings, reducing the life of the castor
- We recommend that fittings are checked regularly & re-tightened as required
- Where anti-static castors are fitted, we recommend that the tyre is periodically cleaned to remove any dirt that may have built up. This will ensure a good contact between the tyre & the floor

# Locks & Keys

- Locks should operate smoothly. If they do not please arrange urgent attention, & where necessary, obtain a replacement
- Locks require minimal maintenance & should only be lubricated using a dry lubricant such as graphite which can be purchased from most locksmiths (NEVER USE OIL)
- When ordering replacement locks or keys, please ensure that the product code & any letters / numbers from the key & lock are available

#### Paint

- · We only use lead free epoxy polyester powder paint
- Where possible colours are applied in special booths to allow for overspray to be recycled
- Alternative colours are available at additional cost

#### **Stainless Steel**

- We use 304 & 430 grade stainless steel
  - · Grade 304 is non-magnetic & used on products which are not painted
  - · Grade 430 is magnetic & will be subsequently painted
- Please note that the term 'stainless' is not entirely accurate, as certain chemicals can 'eat' into the surface. In particular, 'iodine' which can badly corrode / stain, & as a result should not be used as a cleaning agent. Stainless steel products should always be dried after cleaning to reduce staining

### Wood, Melamine & Laminate

- We only use FSC accredited materials. This provides assurances that the material is produced from well-managed forests & / or recycled materials
- We use a combination of melamine faced MDF / chipboard & high pressure laminates bonded MDF of chipboard substrates
- Edges are protected with 2mm ABS, T-Barb edging or sealed with varnish
- · Both melamine & laminate panels offer some resistance to water penetration under normal usage
- · We only use water based adhesives & varnishes

# **Upholstered Panels & Fabrics**

- Vinyl upholstered panels are manufactured from 'Flame Retardant' materials, independently tested to BS 5852: 2006 Ignition Source 5, BS EN 1021-1: 2006, BS EN 1021-2: 2006 & BS 7176: 2007 Medium Hazard. Standard vinyl colours also contain an anti-microbial compound within its construction, which constantly migrates to the surface to replenish any material which is removed by cleaning or general usage
- Fabric upholstered panels are manufactured from 'Flame Retardant' materials, independently tested to BS 5852: 2006 Ignition Source 5, BS EN 1021-1:2006, BS EN 1021-2:2006 & BS 7176: 2007 Low Hazard
- Curtain screens are manufactured from 'Flame Retardant' materials, independently tested to BS 5438: Part 2

**Fechnical Information** 

#### **Electric Columns & Actuators**

 Electric columns & actuators comply with EMC Directive 2014/30/EU, Low Voltage Directive 2014/35/EU & RoHS2 Directive 2011/65/EU

## **Medical Devices**

Where applicable, Hospital Metalcraft product is marked with the CE symbol to demonstrate compliance with the essential requirements of any European Directives (& transposing regulations) that may apply & for which a Declaration of Conformity is held. Hospital Metalcraft is registered with the MHRA (Medicines & Healthcare Regulatory Authority) in connection with product classified as 'Class 1 Medical Devices' for purposes of the Medical Devices Directive (93/42/EEC) with the registration number CA001358.

# **WEEE Regulations**

Hospital Metalcraft meets its obligations as a business to business producer of Electrical and Electronic Equipment (EEE) through:

- Membership of the Veolia Environmental Services WEEE, (Waste Electrical & Electronic Equipment), Producer Compliance Scheme (Registration No. WEE/KE0677RX)
- Marking EEE placed on the market by Hospital Metalcraft with the mandatory crossed out wheeled bin symbol & Bristol Maid<sup>™</sup> logo
- Making information available, (for EEE placed on the market after April 2007), to assist those engaged in treatment, recovery & recycling activities
- Arranging for the collection, treatment, recovery & environmentally sound disposal of obligated WEEE

## **Preventative Maintenance, Spares & Service Support**

Bristol Maid<sup>™</sup> products are designed & manufactured to cater for the demands of the healthcare environment, with the majority (under normal conditions) providing many years of trouble free service. Where annual maintenance or periodic inspections are required, we can offer comprehensive planned preventative maintenance through either a scheduled service contract or call out by our own trained Engineers. For your peace of mind & continued product reliability, please contact our customer service team for details on the following services.

- · Preventative maintenance contracts & packages
- Supply of spare parts
- · Site audit & equipment recommendations

### **Cleaning Guidance & Instructions**

The maintenance of high standards of cleanliness on all surfaces & equipment is a crucial factor in the prevention & control of infection. Whilst all dust, dirt & moisture can harbour infection, the key risks are associated with contamination arising from contact with blood, body fluids, mucous membranes & damaged skin. In all cases surfaces, must be thoroughly cleaned & disinfected in order to destroy any pathogenic micro-organisms.

Bristol Maid products have been designed for use in demanding environments. Guidance & instructions can be found at - www.bristolmaid.com/cleaning-instructions





# **Caretray Trolley™ Range**

The Caretray Trolley was developed specifically to meet the needs of the Healthcare Professional across a wide variety of disciplines including Minor Treatment, A&E, Theatre, Phlebotomy, ITU, and Outpatients to name but a few. It has become a flexible, dependable partner, a trusted workhorse in hospitals throughout the UK and Overseas.

In its 25th year of production, the Caretray Trolley offers excellent value for money & is available as either pre-configured or can be specified to meet your specific requirements...

#### **Pre-Configured Trolleys**

- Chart Trolleys
- Dispensing Trolleys
- Phlebotomy Trolleys
- Procedure Trolleys
- Resuscitation Trolleys
- Sepsis Trolleys
- Laptop Trolleys

#### Configure Your Own Trolley (Mild & Stainless Steel)

- Low Level
- Standard Level
- High Level

#### Configure Your Own Rack (Mild & Stainless Steel)

- Mobile
- Static





# Manufacturing Medical Furniture & Equipment Since 1953

With over 65 years of experience, Bristol Maid is a leading manufacturer and supplier of medical furniture & equipment.

From design and manufacture to supply and installation, Bristol Maid has a reputation for providing high product quality, excellent customer service and value for money.

We take pride in working in close partnership with our customers applying innovative design and technology to deliver effective solutions for today while anticipating the needs of tomorrow.





Hospital Metalcraft Limited Blandford Heights, Blandford Forum, Dorset. DT11 7TG, UK

Tel: +44 (0) 1258 451338 Fax: +44 (0) 1258 455056 sales@bristolmaid.com

www.bristolmaid.com

