

Medisinsk utstyr og helseartikler

Caretray™ Trolleys & Racks

Disassembly & Cleaning	44
Stainless Steel Trolleys, Low Level	46
Stainless Steel Trolleys, Standard Level	48
Stainless Steel Trolleys, High Level	50
Stainless Steel Racks, Mobile	52
Stainless Steel Racks, Static	54
Mild Steel Trolleys, Low Level	56
Mild Steel Trolleys, Standard Level	58
Mild Steel Trolleys, High Level	60
Mild Steel Racks. Mobile	62

Mild Steel Racks, Static 64

Features 42

Options 66

Caretray™ Trolley, Rack & Cabinet Range

The Caretray™ Trolley was developed specifically to meet the needs of the Healthcare Professional across a wide variety of disciplines including Minor Treatment, A&E, Theatre, Phlebotomy, ITU, and Outpatients to name but a few. The Caretray™ trolley is your flexible, dependable partner and has become a trusted workhorse in hospitals throughout the UK and Overseas.

The range now comprises:

- Caretray™ Trolleys Low, standard & high level single or double column mild or stainless steel
- Caretray™ Racks Mobile & static racks single or double column mild or stainless steel
- Drug & Medicine Cupboards 870mm or 1870mm high
- Storage Cupboards 1870mm high
- Chart Trolleys Standard level single column mild steel
- Phlebotomy Trolleys Standard level single column mild or stainless steel
- · Resuscitation Trolleys Standard level single column mild steel



- www.gymo.no - post@gymo.no - (+47) 4000 3522



Distributør i Norge: Gymo AS - www.gymo.no - post@gymo.no - (+47) 4000 3522

Caretray[™] Trolley Features

Carerails

• Fitted to both sides of the trolley, compatible with a wide range of accessories

Side Panels

- · One piece moulded side panels
- Smooth surfaces, rounded edges & contours assists cleaning (compatible with steam cleaning)
- Runners with positive stops prevent trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)

Easy Clean

- · One piece mouldings used throughout
- Smooth surfaces, rounded edges & contours assists cleaning
- Steam cleaning compatible
- Quick & easy tool-less disassembly assists thorough or deep cleaning programmes



-Top

- One piece moulded top prevents ingress of dirt
- Smooth surfaces, rounded edges & contours assists cleaning (compatible with steam cleaning)
- Raised edges on three sides provides defined grab points to act as push handles, whilst an integral rim at the front contains spillages
- Maximum load of 30kg

Frame

- · Frames are fully welded eliminating any gaps
- Available in either 304 grade stainless steel or mild steel that is blast finished & powder coated in a durable epoxy polyester paint



Ring Buffers

- Ring buffers offer protection to both the trolley & fabric of the building
- · Smooth surfaces assist cleaning

Trays / Dividers / Plastic Clip On Label Cover / Labels

- · Translucent, impact-resistant trays & dividers
- Smooth surfaces, rounded edges & contours assists cleaning (compatible with steam cleaning & able to withstand high temperature machine washing)
- Trays feature an integral handle with plastic clip on label cover, (suitable for use with our perforated coloured label strips)
- Trays are available in 100 or 150mm depths & can be subdivided in both length & width using simple slot in dividers
- Each tray is supplied with 4 removable dividers (2 x long & 2 x short)
- Maximum load of 5kg



Castors

- Trolleys are fitted with 100mm pressed steel castors
- · Two braking castors are fitted as standard
- Optional 100mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lower noise levels & smooth surfaces to assist cleaning



Optional Polymer Castors

Doors & Back Panel

- · Solid grade laminate doors & back panel
- Smooth surfaces, rounded edges & contours assists cleaning (compatible with steam cleaning)
- 270 degree opening hinges allow the doors to fold flat against the sides of the trolley
- Three locking options...











Caretray[™] Trolley Disassembly & Cleaning

The design of the trolley allows for quick & easy tool-less disassembly providing access for all areas when conducting thorough or deep cleaning programmes. The simple guide below illustrates the disassembly & re-assembly process...

Remove Trays

Remove the trays from the unit

2 Remove Top

The top is held in place with two thumb screws. These are located on the underside of the top. Once these are removed the top can be removed with a firm tug

Remove Back Panel (if fitted)

Once the top has been removed the back panel can be lifted clear of the cross members on the rear of the frame









Remove Side Panels

The side panels can now be removed by lifting off of the positioning inserts. Double height trolleys require thumb screws at the top of the lower panel to be removed.

5 Clean

The frame, side panels & top can now be cleaned

6 Reassemble

The trolley can now be re-assembled by repeating the steps in reverse









Caretray[™] Trolleys - Stainless Steel, Low Level, Single & Double Column





CTS1/08/3S1D - Shown with optional CT/WF, CT/AGH & GMBH/15



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Features

- Low level single & double column easy clean procedure trollevs. Height makes it suitable to be stored under bench
- Designed for guick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing travs from being accidentally withdrawn
- · Carerail fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover. (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white. blue, green, vellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building: smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 782mm
 - Double Column 915 x 455 (internal) at height 782mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Description Weight

Single Column - w 550 x d 550 x h 807mm CTS1/08/1S2D 1 x 100mm, 2 x 150mm Trays 21.2kg CTS1/08/3D 3 x 150mm Trays 21.5kg CTS1/08/3S1D 3 x 100mm, 1 x 150mm Trays 22.3kg CTS1/08/4S 4 x 100mm Trays 21.9kg

Options

- · CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- . CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR/2 & CT/CR/5 Catheter Racks
- CT/CYD Cylinder Holder C/CD/D Size
- Single & Double Doors See page 72 . CT/FS - Flat Shelf, Removable
- . CT/IV Infusion Pole & Bracket
- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box. 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- Polymer Castors See page 73

Item

Optional Polymer Castors

- Easier to manoeuvre
- Reduced noise levels
- Easy to clean

See page 73







Item	Description	Weight
Double Column - w 1005 x d 550 x h 807r		
CTS2/08/3S4D	3 x 100mm, 4 x 150mm Trays	37.3kg
CTS2/08/4S3D	4 x 100mm, 3 x 150mm Trays	36.9kg
CTS2/08/6D	6 x 150mm Trays	36.5kg
CTS2/08/6S2D	6 x 100mm, 2 x 150mm Trays	38.1kg
CTS2/08/7S1D	7 x 100mm, 1 x 150mm Trays	37.7kg
CTS2/08/8S	8 x 100mm Trays	37.3kg

Distributør i Norge: Gymo AS - www.gymo.no - post@gymo.no - (+47) 4000 3522

Caretray[™] Trolleys - Stainless Steel, Standard Level, Single & Double Column





CTS1/10/3S2D - Shown with optional CT/IV, CT/CR2, GMBH/15 & CT/AGH



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Features

- Standard level single & double column easy clean procedure trolleys
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides and integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerails fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 968mm
 - Double Column 915 x 455 (internal) at height 968mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Item	Description	Weight
Single Column - w 550	c d 550 x h 993mm	
CTS1/10/1S3D	1 x 100mm, 3 x 150mm Trays	27.9kg
CTS1/10/3S2D	3 x 100mm, 2 x 150mm Trays	29.0kg
CTS1/10/4D	4 x 150mm Trays	28.3kg
CTS1/10/4S1D	4 x 100mm, 1 x 150mm Trays	28.7kg
CTS1/10/6S	6 x 100mm	29.8kg

Options

- CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR/2 & CT/CR/5 Catheter Racks
- CT/CYD Cylinder Holder C/CD/D Size
- Single & Double Doors See page 72
- CT/FS Flat Shelf, Removable
- CT/IV Infusion Pole & Bracket

- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box. 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- Laptop Holders See page 70
- CT/LH/C5 Motion C5 Holder, Swing Arm
- CT/OB/SS Overbridge
- Tilt Bins See page 71
- Polymer Castors See page 73

Optional Polymer Castors

- Easier to manoeuvre
- Reduced noise levels
- Easy to clean

See page 73







Item	Description	Weight
Double Column - w 1005 x d 550 x h	993mm	
CTS2/10/3S6D	3 x 100mm, 6 x 150mm Trays	46.6kg
CTS2/10/4S5D	4 x 100mm, 5 x 150mm Trays	46.2kg
CTS2/10/5S4D	5 x 100mm, 4 x 150mm Trays	45.9kg
CTS2/10/6S4D	6 x 100mm, 4 x 150mm Trays	47.3kg
CTS2/10/8D	8 x 150mm Trays	45.8kg
CTS2/10/8S2D	8 x 100mm, 2 x 150mm Trays	46.6kg

Distributør i Norge: Gymo AS

CTS1/14/4S3D

Item Description Weight Single Column - w 550 x d 550 x h 1410mm CTS1/14/2S5D 2 x 100mm, 5 x 150mm Trays 36.2kg CTS1/14/3S4D 3 x 100mm, 4 x 150mm Trays 35.8kg CTS1/14/4S3D 4 x 100mm, 3 x 150mm Trays 35.4kg CTS1/14/5S3D 37.8ka 5 x 100mm, 3 x 150mm Travs CTS1/14/6D 6 x 150mm Trays 35.0kg CTS1/14/8S1D 8 x 100mm, 1 x 150mm Trays 37.8kg

Caretray™ Trolleys - Stainless Steel, High Level, Single & Double Column



CTS1/14/4S3D - Shown with optional CT/D4/A4



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Features

- High level single & double column easy clean procedure trolleys
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerail fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 1385mm
 - Double Column 915 x 455 (internal) at height 1385mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Options

- CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- OT/BIT//CT BoodinonCTToldol, /CT
- CT/CR/2 & CT/CR/5 Catheter Racks
- CT/CYD Cylinder Holder C/CD/D Size
- CT/FS Flat Shelf, Removable
- Single & Double Doors See page 72

- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box. 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)
- C006 100mm Polymer Castor, 4 x Non-Braking (set of 4)

Optional Polymer Castors

- Easier to manoeuvre
- Reduced noise levels
- Easy to clean

See page 73



Item	Description	Weight
Double Column - w 100	5 x d 550 x h 1410mm	
CTS2/14/4S10D	4 x 100mm, 10 x 150mm Trays	60.6kg
CTS2/14/6S8D	6 x 100mm, 8 x 150mm Trays	62.8kg
CTS2/14/8S6D	8 x 100mm, 6 x 150mm Trays	63.0kg
CTS2/14/10S6D	10 x 100mm, 6 x 150mm Trays	65.0kg
CTS2/14/12D	12 x 150mm Trays	61.3kg
CTS2/14/16S2D	16 x 100mm, 2 x 150mm Trays	66.6kg

CTS2/14/10S6D

Caretray™ Racks - Stainless Steel, Mobile, Single & Double Column



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Features

- Single & double column easy clean easy clean rack
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 1757mm
 - Double Column 915 x 455 (internal) at height 1757mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Options

- CT/FS Flat Shelf, Removable
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)
- C006 100mm Polymer Castor, 4 x Non-Braking (set of 4)

Item	Description	Weight
Single Column - w 550 >	c d 550 x h 1782mm	
CRSM1/18/2S6D	2 x 100mm, 6 x 150mm Trays	41.7kg
CRSM1/18/3S6D	3 x 100mm, 6 x 150mm Trays	44.1kg
CRSM1/18/4S5D	4 x 100mm, 5 x 150mm Trays	42.8kg
CRSM1/18/5S5D	5 x 100mm, 5 x 150mm Trays	44.3kg
CRSM1/18/6S4D	6 x 100mm, 4 x 150mm Trays	43.9kg
CRSM1/18/8D	8 x 150mm Trays	43.3kg

CRSM1/18/6S4D

Optional Polymer Castors

- Easier to manoeuvre
- · Reduced noise levels
- Easy to clean

See page 73



Quick & easy disassembly when conducting thorough or deep cleaning programmes



CRSM2/18/8S10D fitted with C004 (polymer castors)





Distributør i Norge: Gymo AS - www.gymo.no - post@gymo.no - (+47) 4000 3522

CRSS1/16/6S4D

Caretray™ Racks - Stainless Steel, Static, Single & Double Column



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Features

- Single & double column easy clean easy clean rack
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Adjustable feet allow levelling on uneven floors

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 1645mm
 - Double Column 915 x 455 (internal) at height 1645mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray -355 x 405 x 150mm (internal)

Options

• CT/FS - Flat Shelf, Removable

Item	Description	Weight	
Single Column - w 536 x d 5	536 x h 1670 - 1685mm		
CRSS1/16/2S6D	2 x 100mm, 6 x 150mm Trays	38.4kg	
CRSS1/16/3S6D	3 x 100mm, 6 x 150mm Trays	40.8kg	
CRSS1/16/4S5D	4 x 100mm, 5 x 150mm Trays	39.5kg	
CRSS1/16/5S5D	5 x 100mm, 5 x 150mm Trays	41.0kg	
CRSS1/16/6S4D	6 x 100mm, 4 x 150mm Trays	40.6kg	
CRSS1/16/8D	8 x 150mm Trays	40.0kg	
CRSS1/16/8D	8 x 150mm Trays	40.0kg	



Quick & easy disassembly when conducting thorough or deep cleaning programmes



Item	Description	Weight
Double Column - w 995 x	d 536 x w 1670 - 1685mm	
CRSS2/16/4S12D	4 x 100mm, 12 x 150mm Trays	71.1kg
CRSS2/16/6S12D	6 x 100mm, 12 x 150mm Trays	74.0kg
CRSS2/16/8S10D	8 x 100mm, 10 x 150mm Trays	72.3kg
CRSS2/16/10S10D	10 x 100mm, 10 x 150mm Trays	75.3kg
CRSS2/16/12S8D	12 x 100mm, 8 x 150mm Trays	74.6kg
CRSS2/16/16D	16 x 150mm Travs	71.5kg



Distributør i Norge: Gymo AS - www.gymo.no - post@gymo.no - (+47) 4000 3522

Caretray™ Trolleys - Mild Steel, Low Level, Single & Double Column





CTM1/08/3S1D - Shown with optional CT/WF, SBH/15 & CT/AGH



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Features

- Low level single & double column easy clean procedure trollevs. Height makes it suitable to be stored underbench
- Designed for guick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing travs from being accidentally withdrawn
- Carerail fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic travs feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white. blue, green, vellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building: smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 782mm
 - Double Column 915 x 455 (internal) at height 782mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Spearmint



Item	Description	Weight
Single Column - w	7 550 x d 550 x h 807mm	
CTM1/08/1S2D	1 x 100mm, 2 x 150mm Trays	20.9kg
CTM1/08/3D	3 x 150mm Trays	21.2kg
CTM1/08/3S1D	3 x 100mm, 1 x 150mm Trays	22.0kg
CTM1/08/4S	4 x 100mm Trays	21.6kg

Options

- . CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- . CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR/2 & CT/CR/5 Catheter Racks
- CT/CYD Cylinder Holder C/CD/D Size
- Single & Double Doors See page 72
- CT/FS Flat Shelf, Removable

- . CT/IV Infusion Pole & Bracket
- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box. 15mm Hook
- SBH/15 Sharps Box Holder One Box. 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)
- C006 100mm Polymer Castor, 4 x Non-Braking (set of 4)

Optional Polymer Castors

- Easier to manoeuvre
- Reduced noise levels
- Easy to clean

See page 73

Quick & easy disassembly when conducting thorough or deep cleaning programmes







Distributør i Norge: Gymo AS

CTM2/08/4S3D

- www.gymo.no - post@gymo.no - (+47) 4000 3522

Caretray™ Trolleys - Mild Steel, Standard Level, Single & Double Column



Item	Description	Weight
Single Column - w	/ 550 x d 550 x h 993mm	
CTM1/10/1S3D	1 x 100mm, 3 x 150mm Trays	27.5kg
CTM1/10/3S2D	3 x 100mm, 2 x 150mm Trays	28.6kg
CTM1/10/4D	4 x 150mm Trays	27.9kg
CTM1/10/4S1D	4 x 100mm, 1 x 150mm Trays	28.3kg
CTM1/10/6S	6 x 100mm	29.4kg



CTM1/10/3S2D - Shown with optional CT/IV, CT/CR2, GMBH/15 & CT/AGH



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Features

- Standard level single & double column easy clean procedure trolleys
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides and integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerails fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 968mm
 - Double Column 915 x 455 (internal) at height 968mm
- Travs
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

Paint - Spearmint



Options

- CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR/2 & CT/CR/5 Catheter Racks
- CT/CYD Cylinder Holder C/CD/D Size
- Single & Double Doors See page 72
- CT/FS Flat Shelf, Removable
- CT/IV Infusion Pole & Bracket

- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- Laptop Holders See page 70
- CT/LH/C5 Motion C5 Holder, Swing Arm
- CT/OB/MS Overbridge
- Tilt Bins See page 71
- Polymer Castors See page 73









Item	Description	Weight
Double Column - w 1005	c d 550 x h 993mm	
CTM2/10/3S6D	3 x 100mm, 6 x 150mm Trays	46.0kg
CTM2/10/4S5D	4 x 100mm, 5 x 150mm Trays	45.6kg
CTM2/10/5S4D	5 x 100mm, 4 x 150mm Trays	45.3kg
CTM2/10/6S4D	6 x 100mm, 4 x 150mm Trays	46.7kg
CTM2/10/8D	8 x 150mm Trays	45.2kg
CTM2/10/8S2D	8 x 100mm, 2 x 150mm Trays	46.0kg

Distributør i Norge: Gymo AS



Item	Description	Weight
Single Column - w 5	50 x d 550 x h 1410mm	
CTM1/14/2S5D	2 x 100mm, 5 x 150mm Trays	36.2kg
CTM1/14/3S4D	3 x 100mm, 4 x 150mm Trays	35.8kg
CTM1/14/4S3D	4 x 100mm, 3 x 150mm Trays	35.4kg
CTM1/14/5S3D	5 x 100mm, 3 x 150mm Trays	37.3kg
CTM1/14/6D	6 x 150mm Trays	35.0kg
CTM1/14/8S1D	8 x 100mm, 1 x 150mm Trays	37.2kg

Caretray[™] Trolleys - Mild Steel, High Level, Single & Double Column



CTM1/14/4S3D - Shown with optional CT/DH/A4



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Options

- . CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- . CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR/2 & CT/CR/5 Catheter Racks
- CT/CYD Cylinder Holder C/CD/D Size
- Single & Double Doors See page 72
- . CT/FS Flat Shelf, Removable

Features

- High level single & double column easy clean procedure trolleys
- Designed for guick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerail' fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white. blue, green, vellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 1385mm
 - Double Column 915 x 455 (internal) at height 1385mm
- Travs
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Spearmint



- . CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box. 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)
- C006 100mm Polymer Castor, 4 x Non-Braking (set of 4)

Optional Polymer Castors

- Easier to manoeuvre
- Reduced noise levels
- Easy to clean

See page 73



CTM2/14/8S6D - Shown with optional CTD2/14SS (double door with security seal)

Item	Description	Weight
Double Column - w 1	005 x d 550 x h 1410mm	
CTM2/14/4S10D	4 x 100mm, 10 x 150mm Trays	60.6kg
CTM2/14/6S8D	6 x 100mm, 8 x 150mm Trays	62.8kg
CTM2/14/8S6D	8 x 100mm, 6 x 150mm Trays	62.3.kg
CTM2/14/10S6D	10 x 100mm, 6 x 150mm Trays	65.0kg
CTM2/14/12D	12 x 150mm Trays	61.3kg
CTM2/14/16S2D	16 x 100mm, 2 x 150mm Trays	66.6kg

Distributør i Norge: Gymo AS - www.gymo.no - post@gymo.no - (+47) 4000 3522 CTM2/14/10S6D



Features

- Single & double column easy clean rack
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- 100/150mm deep translucent high impact plastic travs feature an integral handle with plastic clip on label cover. (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white, blue, green, vellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 1757mm
 - Double Column 915 x 455 (internal) at height 1757mm
- Travs Shallow Tray - 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Spearmint

Options

- CT/FS Flat Shelf, Removable
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)
- C006 100mm Polymer Castor, 4 x Non-Braking (set of 4)

Item	Description	Weight	
Single Column - w 550	x d 550 x h 1782mm		
CRMM1/18/2S6D	2 x 100mm, 6 x 150mm Trays	41.2kg	
CRMM1/18/3S6D	3 x 100mm, 6 x 150mm Trays	43.6kg	
CRMM1/18/4S5D	4 x 100mm, 5 x 150mm Trays	42.3kg	
CRMM1/18/5S5D	5 x 100mm, 5 x 150mm Trays	43.8kg	
CRMM1/18/6S4D	6 x 100mm, 4 x 150mm Trays	43.4kg	
CRMM1/18/8D	8 x 150mm Trays	42.8kg	



CRMS1/16/6S4D

Caretray™ Racks - Mild Steel, Static, Single & Double Column



Quick & easy disassembly when conducting thorough or deep cleaning programmes

Features

- Single & double column easy clean rack
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Adjustable feet allow levelling on uneven floors

Dimensions

- Plastic Top
 - Single Column 455 x 455 (internal) at height 1645mm
 - Double Column 915 x 455 (internal) at height 1645mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Spearmint

Options

• CT/FS - Flat Shelf, Removable



Item	Description	Weight
Single Column - w 536	x d 536 x h 1670 - 1685mm	
CRMS1/16/2S6D	2 x 100mm, 6 x 150mm Trays	38.4kg
CRMS1/16/3S6D	3 x 100mm, 6 x 150mm Trays	39.5kg
CRMS1/16/4S5D	4 x 100mm, 5 x 150mm Trays	39.5kg
CRMS1/16/5S5D	5 x 100mm, 5 x 150mm Trays	41.0kg
CRMS1/16/6S4D	6 x 100mm, 4 x 150mm Trays	40.6kg
CRMS1/16/8D	8 x 150mm Trays	39.6kg



Item	Description	Weight
Double Column - w 995 x	d 536 x h 1670 - 1685mm	
CRMS2/16/4S12D	4 x 100mm, 12 x 150mm Trays	70.2kg
CRMS2/16/6S12D	6 x 100mm, 12 x 150mm Trays	73.2kg
CRMS2/16/8S10D	8 x 100mm, 10 x 150mm Trays	72.3kg
CRMS2/16/10S10D	10 x 100mm, 10 x 150mm Trays	75.3kg
CRMS2/16/12S8D	12 x 100mm, 8 x 150mm Trays	74.6kg
CRMS2/16/16D	16 x 150mm Trays	71.5kg



Options

Work Flap

Features

- Fold down flap (provides work area 300 x 455mm)
- Machined from solid grade laminate allowing for easy cleaning
- Attaches to upper carerail
- Maximum load 2kg

Standard Colours

- · Paint Grey White
- Solid Grade Laminate Grey White

Dimensions (w x d x h)

• 325 x 465 x 70mm



Item	Description	Weight
CT/WF	Work Flap	3.2kg

Infusion Pole & Bracket

Features

- Two hook infusion pole & bracket (pole length 1160mm)
- Attaches to upper carerail
- Maximum load on hooks 2kg

Standard Colours

· Paint - Grey White

Dimensions (w x d x h)

• 238 x 104 x 1157mm



Item	Description	Weight
CT/IV	Infusion Pole & Bracket	1.6kg

Distributør i Norge: Gymo AS - www.gymo.no - post@gymo.no - (+47) 4000 3522

Cylinder Holder - C/CD/D Size

Features

- Compatible with C/CD/D size cylinders (max diameter 100mm)
- Cylinder held in position using a chain, whilst rubber pads protect wear points & reduce noise
- · Attaches to upper carerail

Standard Colours

• Paint - Grey White

Dimensions (w x d x h)

• 145 x 135 x 400mm



Oxygen cylinder not supplied

Item	Description	Weight
CT/CYD	Cylinder Holder - C/CD/D Size	1.4kg

Airway Tube Holder

Features

- Internal dimensions of holder 95 x 50 x 510mm
- · Outer sleeve removable for cleaning
- Fitted with a 15mm square hook
- Attaches to upper carerail

Standard Colours Dimensions (w x d x h)

• Paint - Grey White

• 140 x 70 x 555mm



Item	Description	Weight
CT/AIR	Airway Tube Holder	1.4ka

Alcohol Gel Holder

Features

- Accepts bottles with maximum 90 x 90mm square base
- Fitted with a 15mm hook
- Attaches to upper/lower carerails

Standard Colours

Stainless Steel

Dimensions (w x d x h)

• 95 x 110 x 115mm



Alcohol gel not supplied



Item	Description	Weight
CT/AGH	Alcohol Gel Holder	0.3kg

Glove / Mask Box Holder One Box, 15mm Hook

Features

- Accommodates glove & mask boxes up to 140mm in width
- One piece construction
- 15mm hook attaches to upper/lower carerails

Standard Colours

· Paint - Grey White

Dimensions (w x d x h)

• 143 x 122 x 177mm





Item	Description	Weight
GMBH/15	Glove / Mask Box Holder - One Box, 15mm Hook	0.6kg

Sharps Box Holder One Box, 15mm Hook

Features

- · Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a toggle
- 15mm hook attaches to upper/lower carerails

Standard Colours

• Paint - Grey White

Dimensions (w x d x h)

• 127 x 118 x 230mm



Item	Description	Weight
SBH/15	Sharps Box Holder - One Box, 15mm Hook	0.5kg

Sharps Box Bracket - Universal

Features

- Compatible with the following brackets (customer to supply the plastic brackets & fit fixings supplied).
- · Daniels Healthcare
- POUDS® Wall-Bracket (up to 5L) SHARPSGUARD® mini, 1, 2.5, 3.75, 5 & com plus
- POUDS® Wall-Bracket (7L & above) SHARPSGUARD®
 7. 11.5 & 22
- Frontier (Sharpsafe) 2, 3, 4 & 7
- Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS
- 15mm hook attaches to upper/lower carerails

Standard Colours

• Paint - Grey White

Dimensions (w x d x h)

• 75 x 32 x 152mm









Plastic brackets not supplied

Item	Description	Weight
CT/SBB/15	Sharps Box Bracket - Universal	0.2kg

Debris Bin & Bracket

Features

- Plastic debris bin & bracket (internal dimensions 195 x 130 x 95mm)
- 15mm hook attaches to upper/lower carerails

Standard Colours

· Paint - Grey White

Dimensions (w x d x h)

• 210 x 164 x 109mm





Item	Description	Weight
CT/BIN	Debris Bin & Bracket	0.4kg

Catheter Racks

Features

- Internal dimensions of each section 64 x 75 x 400mm Removable clear perspex panel to aid cleaning
- 15mm hook attaches to upper/lower carerails

Standard Colours

• Paint - Grey White

Dimensions (w x d x h)

- CT/CR2 130 x 108 x 495mm
- CT/CR5 325 x 108 x 405mm



Item	Description	Weight
CT/CR/2	Two Section	1.4kg
CT/CR/5	Five Section	3.0kg

Document Holder

Features

- A4 size suitable for A4 ring binders or foolscap notes (pocket size 274 x 63mm)
- Open front design allows for easy patient identification
- 15mm hook attaches to upper/lower carerails

Standard Colours

• Paint - Grey White

Dimensions (w x d x h)

• 275 x 82 x 402mm



lt	tem	Description	Weight
C	CT/DH/A4	Document Holder	1.7kg

Flat Shelf - Removable

Features

- Supported on the runner guides (non pull out)
- Machined from solid grade laminate allowing for easy cleaning

Standard Colours

• Solid Grade Laminate - Grey White

Dimensions (w x d x h)

• 413 x 457 x 4mm





Item	Description	Weight
CT/FS	Flat Shelf - Removable	1.1kg

Laptop Holder - Removable, Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with models CTM1/10, CTM2/10, CTS1/10, CTS2/10
- Laptop holder fitted to the rear right upright of the trolley
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Laptop supported in a cradle fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)
- · Design allows easy access to laptop ports
- Retrofit model available, compatible with existing CT Caretray trolleys

Standard Colours

· Paint - Grey White

Laptop Holder - Non-Removable, Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with models CTM1/10, CTM2/10, CTS1/10, CTS2/10
- Laptop holder fitted to the rear right upright of the trolley
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Laptop supported in a cradle fitted with a steel mounting bracket that can be fixed using either nylock or wing nuts
- Laptop can be secured using a security clamp (supplied)
- Design allows easy access to laptop ports
- Retrofit model available, compatible with existing CT Caretray trolleys

Standard Colours

· Paint - Grey White



Item	Description	Weight
CT/LHR/S	Removable, Swing Arm	4.4kg
R/CT/LHR/S	Removable, Swing Arm (Retro-Fit)	4.4kg

Item	Description	Weight
CT/LHNR/S	Non-Removable, Swing Arm	5.7kg
R/CT/LHNR/S	Non-Removable, Swing Arm (Retro-Fit)	5.7kg

Motion C5 Holder - Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with models CTM1/10, CTM2/10, CTS1/10, CTS2/10
- Motion C5 holder fitted to the rear right upright of the trolley
- Tablet located in a RAM-In-Motion passive dock allowing for easy insertion & removal
- Retrofit model available, compatible with existing CT Caretray trolleys

Standard Colours

· Paint - Grey White



Item Description Weight CT/LH/C5 Swing Arm 3.4kg R/CT/LH/C5 Swing Arm (Retro-Fit) 3.4kg

Tilt Bins

Features

- Available with 4. 5 or 6 tilting compartments
- Clear tilt bins provide visible storage & assists organisation
- Compartments tilt forward to a 45° angle allowing easy access of contents & can be removed to aid refiling/cleaning
- · Label holder for easy identification of the contents
- · Retaining bar secures compartments in place, preventing bins from tipping forward during transit
- Supplied with a mild/stainless steel mounting plate
- Attaches to an overbridge (order separately)

Dimensions (w x d x h)

Compartment Dimensions

- TBIN4 143 x 127 x 190mm
- TBIN5 108 x 95 x 146mm
- TBIN6 89 x 67 x 95mm

Standard Colours

- · Plastic Cream
- · Paint Grev White
- · Stainless Steel Self Colour



Middle: TBIN5

Bottom: TBIN4

Mild Steel	Stainless Steel	Description	Dimensions (w x d x h)	Weight
TBIN4/MS	TBIN4/SS	Four Compartment	600 x 171 x 207mm	2.3kg
TBIN5/MS	TBIN5/SS	Five Compartment	600 x 135 x 164mm	1.6kg
TBIN6/MS	TBIN6/SS	Six Compartment	600 x 93 x 113mm	0.9kg

Overbridge

Features

- · Fixes onto the rear uprights of the trolley
- · Available as either painted mild steel or self-colour stainless steel
- Provides a fixing location for tilt bins (order separately)
- Factory fitted

Standard Colours

- · Paint Grev White
- · Stainless Steel Self Colour

Item	Description	Weight
CT/OB/MS	Overbridge, Mild Steel	2.0kg
CT/OB/SS	Overbridge, Stainless Steel	2.0kg





CTM1/08/3S1D/CTD1/08CA (Single door with CAM Lock)

Locking Options





(Double door with Security Seal)





(Double door with Electronic Push Button Lock)

Caretray[™] Trolley Doors & Back Panels

Features

- Machined solid grade laminate doors allowing for easy cleaning
- Doors fold back against the sides of the trolley when open
- Double door units feature a return to ensure that only one lock is required to secure both doors
- · Three locking options
 - · CAM Lock with Wing Handle
 - Supplied with two keys
 - Wing handle can operate independently of the key
 - Locks keyed to differ, mastered as standard
 - · Keved alike on request
 - Security Seal
 - · Accepts standard plastic security seals
 - Security seals not included with door (see page 281)
 - Code Lock Electronic Push Button
- Units fitted with a door also supplied with a removable back panel to ensure contents of the trolley are secure
- Factory fitted

Code Lock Electronic Push Button

- Uses a 10 button keypad and a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- Four operating codes with differing levels of authority
 - · Master code
 - Sub-master code
 - User code
 - Technician code
- · User codes are 4 digits
- · Master & sub-master codes are 8 digits
- 15.000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code



Description	CAM Lock with Wing Handle	Security Seal	Code Lock Electronic Push Button	Trolley Compatibility
Single Door & Back Panel	CTD1/08CA	CTD1/08SS	CTD1/08CL	CTM1/08, CTS1/08
Double Door & Back Panel	CTD2/08CA	CTD2/08SS	CTD2/08CL	CTM2/08, CTS2/08
Single Door & Back Panel	CTD1/10CA	CTD1/10SS	CTD1/10CL	CTM1/10, CTS1/10
Double Door & Back Panel	CTD2/10CA	CTD2/10SS	CTD2/10CL	CTM2/10, CTS2/10
Single Door & Back Panel	CTD1/14CA	CTD1/14SS	CTD1/14CL	CTM1/14, CTS1/14
Double Door & Back Panel	CTD2/14CA	CTD2/14SS	CTD2/14CL	CTM2/14, CTS2/14

Additional Trays

Features


- 100/150mm deep translucent high impact plastic trays with integral handle (max load per tray 5kg)
- Each tray is supplied with plastic clip on label cover & four removable dividers (2 x long & 2 x short)
- Trays are supplied without coloured label strips (see CT/LAB)

Internal Dimensions (w x d x h) • CT/TRAY/S - 355 x 405 x 100mm • CT/TRAY/D - 355 x 405 x 150mm

Item	Description	Weight
CT/TRAY/S	100mm Shallow Tray	1.5kg
CT/TRAY/D	150mm Deep Tray	1.8kg

Additional Dividers

Features



Item	Description	Weight
CT/DIV/S	Set of 4 dividers for 100mm Shallow Tray	0.6kg
CT/DIV/D	Set of 4 dividers for 150mm Deep Tray	0.4kg

Additional Label Strips

Features

- · Coloured label strips
 - Red
 - Yellow
 - Green
 - Blue
 - White
- Supplied as a set of 6 sheets

Item	Description
CT/LAB	Coloured Label Strips



Polymer Castors

Features

- Offers reduced rolling resistance making the trolley easier to manoeuvre, lower noise levels & smooth surfaces to assist cleaning
- Supplied as a set of 4 castors
 - 2 x braking / 2 x non-braking
 - 4 x braking
 - 4 x non-braking
- · Factory fitted



Item	Description	Weight
C004	100mm, Grey White Polymer, 2 x Braking / 2 x Non-Braking (set of 4)	2.0kg
C005	100mm, Grey White Polymer, 4 x Braking (set of 4)	2.0kg
C006	100mm, Grey White Polymer, 4 x Non-Braking (set of 4)	2.0kg