

OPTIONAL MODIFICATIONS



URGENT with
FIXED
mattress



URGENT with
REMOVABLE
mattress



Mattress for HPL patient area
200x72 cm - thickness: 5 cm, 7 cm or 10 cm




SIDE RAILS

Tilt side rails
stainless steel



Foldable side rails with
SAFE GAPS SYSTEM




MODIFICATIONS



Foldable transport handle
at head part



Transport tray for monitor
or medical equipment




MODIFICATIONS




Oxygen bottle holder
horizontal



Horizontal IV pole holder




Paper roll holder



Pair of stainless steel
bars for accessories



MODIFICATIONS



IV pole holder
at head part



IV pole
steel-chromed / aluminum



4x pedals of central brake



Safety belts
3 pcs



Colour sampler of leatherette



ver. 2016/01

Note: The colour patterns may slightly vary from the actual colours.

PROMA REHA, s.r.o.
Riegrova 342, 552 03 Česká Skalice, Czech Republic, EU
Tel.: +420 491 112 233, e-mail: info@promareha.cz

EUROPEAN UNION
European Regional Development Fund
Operational Programme Enterprise
and Innovations for Competitiveness

www.promareha.com

TRANSPORT STRETCHER URGENT



Transport stretcher URGENT is designed for quick and safe transport of patients in healthcare facilities. Ergonomic and intuitive controls allow instant use of functions. All functions are designed for maximum safety of patient and medical staff. Side rails of the stretcher are equipped with SAFE GAPS SYSTEM which eliminates the risk of pinching or trapping during folding the side rails down. Large wheels and lightweight construction ensure great mobility and handling. Comfortable width allows transport even for larger patients. Thanks to the high variability you can customized the URGENT to affordable version for transfer of patients at regular departments or to the maximally equipped stretcher for use in emergency. URGENT can be fitted with a wide range of accessories.

- ➡ SIDE RAILS WITH SAFE GAPS SYSTEM
- ➡ SUITABLE FOR PATIENTS UP TO 230KG
- ➡ INTUITIVE CONTROLS
- ➡ HIGH VARIABILITY OF MODIFICATIONS
- ➡ UNIVERSAL USE

TRANSPORT STRETCHER URGENT



COMFORTABLE ADJUSTMENT

Adjustment of the back part is assisted by a gas spring, which allows instant effortless adjustment to the desired position. The nursing staff can easily set individual positions by using handles located on both sides of the trolley.

SAFETY PROTECTORS

Bumper wheels at the corners and stainless steel bumpers on the sides protect stretcher and surrounding objects from unwanted damage.

SAFE WORKING LOAD 230 kg

Solid and sturdy construction with safe working load up to 230 kg create a secure support for patient care. Extremely high durability is ensured thanks to used materials and technologies.

SAFE GAPS SYSTEM

Keep your hands
safe with NEW type
of side rails



FOLDABLE TRANSPORT HANDLE

Transport handle at patient head allows fast folding down in case of emergency.

SIDE RAILS WITH SAFE GAPS SYSTEM

Fast folding of the side rail ensures immediate access to the patient. Thanks to the SAFE GAPS SYSTEM there is not risk of pinching or catching a hand during folding.

PAPER ROLL HOLDER

Holder under patient area for attachment of hygienic paper roll.

SPACE FOR OXYGEN BOTTLE

Oxygen bottle can be safely transported together with the stretcher thanks to the holder underneath patient area.

EASY HEIGHT SETTING

Lift of the patient area is controlled by a foot pedal accessible from both sides. Wide range of vertical positioning ensures compatibility with most hospital beds.

PLASTIC UNDERCARRIAGE COVER

It facilitates and speeds up daily cleaning of the stretcher. The cover protects the mechanisms of the undercarriage thereby extends the life of and reliability of the stretcher.

UNIVERSAL WHEELS

High-quality medical wheels make transport and handling of the bed very easily. The wheels are made of permanent non-marking material. Plastic wheel covers make cleaning easier, protect the bearings from dirt and significantly extend their lifetime.

OPTIMAL WIDTH

Width of the stretcher ensures comfortable transport even for larger patients and also provides perfect access of medical staff to the patient.

FIXED OR REMOVABLE MATTRESS

The stretcher can be equipped with fixed upholstery covered with sanitary leatherette or removable mattress with handles.

CENTRAL BRAKE

Pedals accessible from both sides of the stretcher allow safe braking. Locking one of the wheels in the straight position ensures easy transport by one person.

SPACE FOR BELONGINGS

Patient belongings or necessary medical equipment can be stored in space underneath the patient area.

ACCESSORIES


Compatibility with accessories is ensured by universal holders located at the head part. They allow the use of the most common accessories such as infusion pole. If the IV pole will not be in use, it can be placed in a holder under the patient area.

TECHNICAL SPECIFICATIONS


Type of mattress	Fixed / Removable
Overall dimensions	87 x 208 cm
Dimensions of patient area	72 x 200 cm
Height of patient area	56 – 94 cm / 51 - 89 cm (without mattress)
Height adjustment	hydraulic
Angle of back part	0 - 85°
Back part adjustment	gas spring
Wheels	200 mm
Safe working load	230 kg
Weight	82 kg / 78 kg
Central brake pedals	2x
Bumper wheels	4x
Plastic undercarriage cover	✓
Fixed transport handle	1x

OPTIONAL MODIFICATIONS


Removable mattress (5 cm, 7 cm, 10 cm)
Foldable side rails with SAFE GAPS SYSTEM
Tilt side rails
Foldable transport handle at head part
Paper roll holder
IV pole holder at head part
IV pole holder under patient area
IV pole - aluminum
IV pole - steel-chromed
Oxygen bottle holder
Transport tray for monitor or medical equipment
Pair of stainless steel bars for accessories
Safety belts




Back part adjustment (0° - 85°)



Height adjustment (56 - 94 cm)



Folding side rails with SAFE GAP SYSTEM



Foldable handle